

VIDAS® B•R•A•H•M•S PCT™ for Procalcitonin (PCT) Guided Antibiotic Therapy

November 10, 2016

bioMérieux

Meeting of the Microbiology Devices Panel

Urgent Need for Biomarkers for Appropriate Antibiotic Use and Reducing Antibiotic Resistance

Mark Miller, MD

Chief Medical Officer

bioMérieux

Serious and Growing Threat to U.S. and Global Public Health

- Antibiotic overuse, inappropriate initiation and prolonged use
- Safety risk to patients and contributed to rise of antibiotic resistance
- 2 million illnesses and ~23,000 deaths per year in U.S.¹

Diagnosing Bacterial Infection that Requires Antibiotics is Difficult



- Bacterial cultures can take 2-3 days to perform
- Often provide false positive or false negative results
- Faster, more accurate indicators of infection needed to make critical antibiotic decisions
- Diagnostic tests can support antimicrobial stewardship

50% of the 69 Million Antibiotic Prescriptions for Acute Respiratory Conditions are Unnecessary

34.3 Million

Antibiotic
prescriptions
unnecessary

34.6 Million

Antibiotic
prescriptions
appropriate



When Used Inappropriately, Antibiotics Carry Risks without Benefit

- Side effects and toxicity¹
 - Allergic reactions, diarrhea, other GI issues
- Collateral damage²
 - Such as *Clostridium difficile* infections
- Antibiotic resistance³
- Spread of drug-resistant pathogens⁴

1. NHS. *Antibiotics – side effects*

2. Lessa et al. *NEJM* 2015;372:825-34

3. Laxminarayan et al. *Lancet Infect Dis* 2013;13:1057-98

4. WHO Fact Sheet: www.who.int/mediacentre/factsheets/fs194/en/

Increasing Number of Drug-Resistant Infections

- More serious illness and disability
- Higher death rate
- Need for more complex therapeutic options
- Extended hospitalizations

PCT Assay to Help Guide Appropriate Antibiotic Use in LRTI and Sepsis

- Lower respiratory tract infection
 - Whether to **start** antibiotics
 - When to **stop** antibiotics
- Sepsis
 - When to **stop** antibiotics

LRTI: Whether to Start Antibiotics Difficult to Determine

- 78-year-old male in ER
 - Fever, cough, chest pain
- Examination: Fever, abnormal chest sounds, elevated white blood cell count, abnormal chest X-ray
- Diagnosis: Community acquired pneumonia?
 - Not definitively bacterial or viral

***Physicians will often prescribe antibiotics
“just in case”***

Sepsis: When to Discontinue Antibiotics Difficult to Determine

- 50-year-old female in ICU
 - Shock of undetermined origin
- Examination: Fever, multiple lung abnormalities on X-ray and CT scan of chest
- Diagnosis: Sepsis from pneumonia?
 - Started on antibiotics, ventilation

Unclear when/if she stops benefitting from continued antibiotic treatment

Current Use: VIDAS B•R•A•H•M•S PCT Assay

Two current FDA-cleared uses

1. To aid in risk assessment for sepsis and septic shock on first day of ICU admission
2. To aid in assessing the cumulative 28-day risk of all-cause mortality for patients with severe sepsis or septic shock in the ICU using change of PCT levels over time

Two New Proposed Indications: VIDAS B•R•A•H•M•S PCT Assay

Used in conjunction with other laboratory findings and clinical assessments, VIDAS B•R•A•H•M•S PCT is intended for use as follows:

- *To aid in decision making on antibiotic therapy for inpatients or outpatients, with suspected or confirmed lower respiratory tract infections (LRTI) defined as community-acquired pneumonia (CAP), acute bronchitis and acute exacerbation of chronic obstructive pulmonary disease (AECOPD);*
- *To aid in decision making on antibiotic discontinuation for patients with suspected or confirmed sepsis.*

Existing Evidence from Randomized Controlled Trials Supports Indications for Use

- Nearly two dozen RCTs have evaluated safety and effectiveness of PCT guidance for antibiotic therapy for LRTI and sepsis
- Plan for systematic reviews and meta-analyses discussed through interactive discussions with FDA
- Goal: Determine whether antibiotic use could be safely reduced

Presentation Overview

- Biology of PCT and role as a host-response biomarker stimulated by bacterial infection
- Results of meta-analyses: PCT-guided treatment decisions can reduce antibiotic use without compromising safety
- VIDAS B•R•A•H•M•S PCT can improve the appropriate use of antibiotics with no added risk

Agenda

Diagnostic Utility of PCT

Sam Bozzette, MD, PhD

Vice President, Medical Affairs – Americas
bioMérieux

Meta-Analysis Methodology

Noam Kirson, PhD

Vice President
Analysis Group, Inc.

Meta-Analysis Results

Philipp Schuetz, MD, MPH

Chief Physician of Endocrinology and
Internal Medicine
University of Basel

Benefit-Risk Profile

Mark Miller, MD

Chief Medical Officer
bioMérieux

Diagnostic Utility of Procalcitonin

Sam Bozzette, MD, PhD

Vice President, Medical Affairs – Americas
bioMérieux

PCT is a Host-Response Biomarker Stimulated by Bacterial Infection


- PCT is a precursor of calcitonin
- C cells of thyroid produce PCT, convert to calcitonin
- PCT also produced by neuroendocrine cells
 - Normal serum concentrations <0.05 ng/mL
- Bacterial infection stimulates PCT production¹⁻³
 - Concentrations typically >0.25 ng/mL
 - Can rise 1000-fold in septic shock

1. Linscheid P et al. *Crit Care Med* 2004;32:1715-21


2. Linscheid P et al. *Endocrinology* 2003;144:5578-84

3. Linscheid P et al. *Endocrinology* 2005;146:2699-708


PCT Levels Significantly Higher in Patients with Bacterial Co-Infection




PCT Properties Favorable for Antibiotic Decision Making




PCT Levels Correlate with Disease Severity



Changes in PCT Levels Have Prognostic Significance in Sepsis Patients



PCT with Clinical Judgment Superior to Clinical Judgment Alone in Sepsis



PCT Levels Have a High Negative Predictive Value in LRTI


- NPV = probability condition is absent given negative test

| | Endpoint (Prevalence) | Sensitivity | Specificity | PPV | NPV |
|------------------------|--|-------------|-------------|-----|-----|
| Rodriguez ¹ | Confirmed bacterial co-infection (20%) | 90% | 31% | 25% | 92% |
| Stolz ² | Need for antibiotics (24%) | 84% | 98% | 93% | 94% |


1. Rodriguez et al. *J Infect* 2016;72:143-51

2. Stolz et al. *Swiss Med Wkly* 2006;136:434-40


bioMérieux Antibiotic Management Algorithms: Two Components



LRTI: bioMérieux Proposed Antibiotic Initiation Algorithm



LRTI and Sepsis: bioMérieux Proposed Antibiotic Discontinuation Algorithms



PCT Guidance is Widely Used

- PCT in use clinically for 10 years
- ~50% of U.S. hospitals use PCT¹
- PCT guidance incorporated in several treatment guidelines, for example:
 - *International Guidelines for the Management of Severe Sepsis and Septic Shock*
 - *European Respiratory Society Guidelines*
 - *IDSA Guidelines for Antibiotic Stewardship Programs*
- >36 million PCT determinations worldwide in 2015²

1. Kadri et al. Epidemiology of Procalcitonin Use in United States Hospitals. *Submitted*.

2. Thermo-Fisher

VIDAS B•R•A•H•M•S PCT

- Automated, individual test
- Assay takes 20 minutes
- Two components:
 - VIDAS B•R•A•H•M•S PCT kit
 - VIDAS family of instruments

VIDAS B•R•A•H•M•S PCT Produces Accurate Results

- Performed on serum or lithium heparinate plasma
- Detection limit: 0.03 ng/mL
- Quantitative linear range: 0.05-200 ng/mL

PCT is a Useful, Reliable Biomarker for Detecting Bacterial Infection

- PCT levels
 - Rise rapidly after bacterial insult (4-6 hours)
 - Decline rapidly with control of infection (half-life of ~24 hours)
 - Correlate with severity of illness
- Aids in assessing risk and prognosis

PCT is Useful Tool for Antibacterial Decision Making

- High NPV for bacterial infection
- Adds to clinical judgment
- Widely accepted and cited in international treatment guidelines
- Measured in the VIDAS system with accuracy and precision

Meta-Analysis Methodology

Noam Kirson, PhD

Vice President

Analysis Group, Inc.

Meta-Analysis Evaluated Safety and Effectiveness of PCT-Guided Approach

- Worked with FDA to develop meta-analysis approach
- Compared PCT guided therapy to standard of care for LRTI and sepsis
- Two primary goals:
 - Evaluate reduction in antibiotic use
 - Evaluate impact of PCT guidance on safety

Meta-Analyses Summarize the State of Published Literature

- Provide greater precision in estimating effects than individual trials
- Different RCTs reflect various clinical factors
- Leverages variation to best capture effect
- Averages effect across multiple RCTs
 - Give different weights to different studies

Random Effects Model Accounts for Variation in Effect Size

- Different study characteristics can result in differences in treatment effects
 - Researchers, settings, patients, etc.
- Confidence intervals appropriately account for variation

Design of RCTs Included in Meta-Analysis

Patients with LRTI or sepsis randomized to:

- Control group
 - Treatment based on clinical judgment under standard of care
- PCT group
 - Treatment based on clinical judgment plus results of PCT assay

Two Types of Meta-Analyses

Study-Level

- Combines aggregate study-level information (e.g., mean differences, odds ratios)
- Limited ability to account for patient-level characteristics (e.g., age, gender)

Patient-Level

- Combines individual-level information
- Access to raw datasets
- Greater flexibility to address heterogeneity in patient characteristics

Established Predefined Methodology: Systematic Literature Review



**Study-
Level**

**LRTI
Sepsis**

**Patient-
Level**

**LRTI
Sepsis**

Study-Level Search Procedures

- Databases searched
 - PubMed
 - Cochrane Database of Systematic Reviews
- Keywords prospectively identified
- Publications selected on pre-defined criteria
- 2 reviewers independently screened all articles

Patient-Level Search Procedures

- Literature identified in prior systematic review conducted in 2011 on PCT and acute respiratory infection¹
- Databases searched
 - Cochrane Controlled Trials Registry
 - MEDLINE
 - Embase
- Publications selected on pre-defined criteria
- 2 reviewers independently screened all articles

Study-Level Literature Review Results

| | LRTI | Sepsis |
|---|---------------|---------------|
| Articles screened based on abstracts | 263 | 340 |
| No original data (e.g., review, editorial) | (136) | (148) |
| PCT not measured for intended purpose | (55) | (100) |
| PCT not measured in target population | (20) | (57) |
| Other | (41) | (25) |
| Retained for meta-analysis | 11 | 10 |
| Total patients in analyses | N=4090 | N=3489 |

Patient-Level Literature Review Results

| | | Articles | |
|----------------------------|--|----------|--------|
| Search of database | | 327 | |
| Not RCTs | | (283) | |
| Ongoing trials | | (14) | |
| Duplicates | | (4) | |
| Other | | (12) | |
| Retained for meta-analysis | | 14 | |
| | | LRTI | Sepsis |
| RCTs in final analyses | | 13 | 5 |
| Total patients in analyses | | N=3142 | N=598 |

Established Predefined Methodology: Data Extraction




Study-Level Analysis: Data Extraction

- Data extraction completed independently by 2 reviewers
- Study characteristics and outcomes were extracted into standardized data form
- Discrepancies between reviewers resolved by consensus

Patient-Level Analysis: Data Extraction

- All data from original acute respiratory infection dataset
- LRTI dataset included patients with:
 - Community-acquired pneumonia (CAP)
 - Acute bronchitis
 - Acute exacerbation of chronic obstructive pulmonary disease (AECOPD)
- Sepsis dataset included patients in ICU with pulmonary infection

Established Predefined Methodology: Meta-Analysis



Effectiveness Endpoints Evaluated in Meta-Analyses

| | LRTI | | Sepsis | |
|-------------------------|-------------|---------------|-------------|---------------|
| | Study-Level | Patient-Level | Study-Level | Patient-Level |
| Antibiotic Initiation | ✓ | ✓ | | |
| Duration of Therapy | ✓ | ✓ | ✓ | ✓ |
| Exposure to Antibiotics | ✓ | ✓ | ✓ | ✓ |

- **Duration**: Time on antibiotics among patients initiating treatment
- **Exposure**: Time on antibiotics among all patients

Example of Duration vs. Exposure

| Patient | Time on Antibiotics (days) | Duration (days) | Exposure (days) |
|---------|----------------------------|-----------------|-----------------|
| #1 | 0 | n/a | 0 |
| #2 | 0 | n/a | 0 |
| #3 | 4 | 4 | 4 |
| #4 | 5 | 5 | 5 |
| #5 | 6 | 6 | 6 |
| | <i>Average</i> | 5 days | 3 days |

- Duration more relevant to patient benefit
- Exposure more relevant to public health benefit


Safety Endpoints Evaluated in Meta-Analyses

| | LRTI | | Sepsis | |
|-------------------------|-------------|---------------|-------------|---------------|
| | Study-Level | Patient-Level | Study-Level | Patient-Level |
| Mortality | ✓ | ✓ | ✓ | ✓ |
| Complications | | ✓ | | |
| Hospital Length of Stay | ✓ | ✓ | | ✓ |
| ICU Length of Stay | | | ✓ | ✓ |

Statistical Methods and Bias Assessment


- Random-effects models
 - Study-level: No adjustment for covariates
 - Patient-level: Adjustment for age
- Subgroup analyses and stratifications performed
 - Study-level: LRTI type, risk of bias, level of adherence
 - Patient-level: LRTI type, PCT level, age, gender, inpatient vs. outpatient setting
- Risk of bias assessed using Cochrane Handbook

Definition of Adherence in Study-Level Meta-Analyses



- Example of non-adherence:
 - LRTI patient with PCT=0.12 ng/mL
 - AB treatment initiated

Meta-Analyses Used Best Practices



Results of Meta-Analyses

Philipp Schuetz, MD, MPH

Chief Physician of Endocrinology and
Internal Medicine

University of Basel

Outline for Meta-Analysis Results

LRTI

Baseline Characteristics

Effectiveness: Initiation of Antibiotics

Effectiveness: Antibiotic Duration/Exposure

Safety

Sepsis

Baseline Characteristics

Effectiveness: Antibiotic Duration/Exposure

Safety

Baseline Characteristics Similar Between Groups in LRTI Patient-Level Meta-Analysis

| Characteristic | PCT Group (N=1536) | Control Group (N=1606) |
|---|-----------------------|---------------------------|
| Age (years), median (IQR) | 66 (50, 79) | 66 (49, 78) |
| Female, n (%) | 671 (44%) | 744 (46%) |
| Diagnosis, n (%) | | |
| CAP | 999 (65%) | 1028 (64%) |
| Acute Bronchitis | 249 (16%) | 282 (18%) |
| AECOPD | 288 (19%) | 296 (18%) |
| PCT value at initiation (ng/mL), median (IQR) | 0.23 (0.10, 0.96) | 0.21 (0.09, 1.04) |

Outline for Meta-Analysis Results

LRTI

Baseline Characteristics

Effectiveness: Initiation of Antibiotics

Effectiveness: Antibiotic Duration/Exposure

Safety


Sepsis

Baseline Characteristics

Effectiveness: Antibiotic Duration/Exposure

Safety

Study-Level Findings for Antibiotic Initiation in LRTI



Outline for Meta-Analysis Results

LRTI

Baseline Characteristics

Effectiveness: Initiation of Antibiotics

Effectiveness: Antibiotic Duration/Exposure

Safety


Sepsis

Baseline Characteristics

Effectiveness: Antibiotic Duration/Exposure

Safety


PCT Guidance Significantly Reduced Duration and Exposure to Antibiotics



Large Reductions in Antibiotic Exposure in Inpatient and Outpatient Settings


Patients on Antibiotics (%)

Inpatient



Patients on Antibiotics (%)

Outpatient




Days from Enrollment

■ PCT
■ Control

Fewer Patients with CAP in PCT Group on Antibiotics Throughout Follow-up

Community-Acquired Pneumonia


■ Control ■ PCT



Fewer Patients with Acute Bronchitis in PCT Group on Antibiotics Throughout Follow-up

Acute Bronchitis


Control PCT



Fewer Patients with AECOPD in PCT Group on Antibiotics Throughout Follow-Up

Acute Exacerbation of COPD

Control PCT



Outline for Meta-Analysis Results

LRTI

Baseline Characteristics

Effectiveness: Initiation of Antibiotics

Effectiveness: Antibiotic Duration/Exposure

Safety


Sepsis

Baseline Characteristics


Effectiveness: Antibiotic Duration/Exposure

Safety

No Adverse Mortality Signal Associated with PCT Guidance for LRTI




Similar Rates of Mortality in LRTI Patient-Level Meta-Analysis



Number at Risk

| | | | | | | | |
|----------------|------|------|------|------|------|------|------|
| PCT | 1536 | 1477 | 1404 | 1330 | 1281 | 1248 | 1020 |
| Control | 1606 | 1553 | 1477 | 1395 | 1338 | 1297 | 1058 |

PCT-Guided Therapy Associated with Lower Incidence of Complications



Complications:

- Death
- ICU admission
- Hospitalization or re-hospitalization
- Acute respiratory infection-specific complications
- Recurrent or worsening infection


30-Day Mortality and Complication Profile of PCT Guidance Consistent Across Subgroups

| Endpoint | PCT (N=1536) | Control (N=1606) | Odds Ratio (95% CI) | Odds Ratio (95% CI) p-value |
|---|-----------------|---------------------|------------------------|-----------------------------------|
| | n (%) | n (%) | | |
| 30-day mortality (Patient-Level) | | | | |
| CAP | 92 (9.2%) | 111 (10.8%) | | 0.92 (0.74, 1.15) p=0.47 |
| Acute Bronchitis | 2 (0.8%) | 0 (0.0%) | | Could not be estimated |
| AECOPD | 9 (3.1%) | 8 (2.7%) | | 1.15 (0.46, 2.89) p=0.76 |
| Complications (Patient-Level) | | | | |
| CAP | 190 (19.0%) | 240 (23.3%) | | 0.77 (0.62, 0.96) p=0.02 |
| Acute Bronchitis | 51 (20.5%) | 54 (19.1%) | | 1.09 (0.70, 1.70) p=0.71 |
| AECOPD | 35 (12.2%) | 45 (15.2%) | | 0.75 (0.46, 1.22) p=0.25 |

0.1 1 10

← Favours PCT Favours Control →

Similar Lengths of Hospital Stay in Both Groups Overall and by LRTI Type



Summary of Findings from LRTI Meta-Analyses

- PCT-guided antibiotic treatment for LRTI associated, on average, with:
 - ~19% relative reduction in antibiotic initiation
 - 1 to 3-day reduction in duration of therapy
 - 3 to 4-day reduction in total exposure
- Across LRTI types, 37-65% reduction in exposure
- No signal for increase in short-term mortality, complications, or length of hospital stay
- Safety and effectiveness findings consistent across subgroups

Outline for Meta-Analysis Results

LRTI

Baseline Characteristics

Effectiveness: Initiation of Antibiotics

Effectiveness: Antibiotic Duration/Exposure

Safety

Sepsis

Baseline Characteristics

Effectiveness: Antibiotic Duration/Exposure

Safety

Baseline Characteristics Similar Between Groups in Sepsis Patient-Level Meta-Analysis

| Characteristic | PCT Group (N=287) | Control Group (N=311) |
|---|----------------------|--------------------------|
| Age (years), median (IQR) | 62 (50, 74) | 65 (53, 75) |
| Female, n (%) | 79 (28%) | 95 (31%) |
| PCT value at initiation (ng/mL), median (IQR) | 1.42 (0.39, 5.78) | 1.20 (0.34, 4.74) |

Outline for Meta-Analysis Results

LRTI

Baseline Characteristics

Effectiveness: Initiation of Antibiotics

Effectiveness: Antibiotic Duration/Exposure

Safety


Sepsis

Baseline Characteristics


Effectiveness: Antibiotic Duration/Exposure

Safety

PCT Guidance Significantly Reduced Duration and Exposure to Antibiotics for Patients with Sepsis



Patients with Sepsis in PCT Group Discontinued Antibiotics Earlier than Control Group



Outline for Meta-Analysis Results

LRTI

Baseline Characteristics

Effectiveness: Initiation of Antibiotics

Effectiveness: Antibiotic Duration/Exposure

Safety


Sepsis

Baseline Characteristics


Effectiveness: Antibiotic Duration/Exposure

Safety

No Adverse Mortality Signal Associated with PCT Guidance for Sepsis




Similar Rates of Mortality in Sepsis Patient-Level Meta-Analysis



Number at Risk

| | | | | | | | |
|----------------|-----|-----|-----|-----|-----|-----|-----|
| PCT | 287 | 278 | 257 | 244 | 229 | 221 | 164 |
| Control | 311 | 305 | 288 | 274 | 253 | 239 | 183 |

Length of ICU Stay and Hospital Length of Stay Similar Between Groups



Summary of Findings from Sepsis Meta-Analyses

- PCT-guided antibiotic treatment for sepsis associated, on average, with:
 - ~1.5 to 3-day reduction in antibiotic use
 - 24% reduction in antibiotic exposure
- Similar mortality rate compared to control
- Similar length of ICU stay and hospital stay compared to control

Summary and Clinical Benefit-Risk Assessment

Mark Miller, MD

Chief Medical Officer

bioMérieux

A Global Public Health Emergency



HIGH-LEVEL MEETING ON ANTIMICROBIAL RESISTANCE



21 SEPTEMBER 2016, UN HEADQUARTERS, NEW YORK



Review on Antimicrobial Resistance

Tackling drug-resistant infections globally

**TACKLING DRUG-RESISTANT
INFECTIONS GLOBALLY:
FINAL REPORT AND
RECOMMENDATIONS**

NATIONAL ACTION PLAN FOR COMBATING ANTIBIOTIC-RESISTANT BACTERIA

MARCH 2015



PCT Levels Better Inform Decisions on Antibiotic Use for LRTI and Sepsis

- 78-year-old male: Possible CAP
 - Initial PCT low, 0.11 ng/mL
 - CT scan found pulmonary emboli
- 50-year-old female: Suspected sepsis
 - Initial PCT low, 0.15 ng/mL
 - Consider diagnosis of heart failure
 - Day 4: PCT 0.05 ng/mL
 - Discontinued antibiotic, good outcome

Rapid PCT Detection Can Help Guide Clinicians on Antibiotic Use

- VIDAS B•R•A•H•M•S PCT guides clinicians to significantly, safely reduce antibiotic use
 - Used with other clinical assessments
 - Can inform decision on initiation or continuation of therapy

Leveraged Knowledge Base of RCTs Evaluating PCT-Guided Approach

- 23 RCTs and more than 7,800 patients
- Evaluated safety and effectiveness of PCT-guided antibiotic decision-making
- Meta-analysis study design is robust
 - Provides greater external validity and generalizability than single RCT
 - Responsive to urgent need for improved antimicrobial stewardship

Summary of Results for LRTI: PCT Guidance Effective and Safe

| LRTI | Study-Level | Patient-Level |
|-------------------------|---------------------------|------------------------|
| Effectiveness | | |
| Initiation | OR = 0.26, $p < 0.001$ | 19% relative reduction |
| Duration | 1.3-day mean reduction | 2.9-day mean reduction |
| Exposure | 2.8-day mean reduction | 3.6-day mean reduction |
| Safety | | |
| Mortality | No safety signal detected | |
| Complications | <i>not evaluated</i> | 15% relative reduction |
| Length of hospital stay | No safety signal detected | |

- Consistent safety profile by age, gender, initial PCT level, inpatient vs. outpatient setting, and type of LRTI

Summary of Results for Sepsis: PCT Guidance Effective and Safe

| Sepsis | Study-Level | Patient-Level |
|-------------------------|---------------------------|---------------------------|
| Effectiveness | | |
| Duration/Exposure | 1.5-day mean reduction | 3.2-day mean reduction |
| Safety | | |
| Mortality | No safety signal detected | |
| Length of ICU stay | No safety signal detected | |
| Length of hospital stay | <i>not evaluated</i> | No safety signal detected |

- Consistent safety profile by age, gender, and initial PCT level


Antimicrobial Stewardship Goal: Right Patient, Right Time, Right Duration

- VIDAS B•R•A•H•M•S PCT can help improve antimicrobial stewardship
- PCT-guided therapy
 - Reduce antibiotic use in LRTI patients with no increased risk
 - Decrease antibiotic exposure in sepsis patients with no impact on safety

Reductions in Antibiotic Exposure with PCT Guidance Exceeds National Goals


Inpatient

National Goal: 20% Reduction



Outpatient

National Goal: 50% Reduction



VIDAS B•R•A•H•M•S PCT: Better, Evidence-Based Decisions

- PCT assay and other clinical information helps
 - Prevent and slow emergence of resistant bacteria
 - Avoid side effects of unnecessary antimicrobials
- Benefits patients and healthcare system

VIDAS® B·R·A·H·M·S PCT™ for Procalcitonin (PCT) Guided Antibiotic Therapy


November 10, 2016

bioMérieux

Meeting of the Microbiology Devices Panel


Backup Slides

Overlap between Study-Level and Patient-Level MA for LRTI



Total: 16 unique RCTs - 4494 unique patients

Overlap Between Study-Level and Patient-Level MA for Sepsis




Total: 11 unique RCTs - 3590 unique patients


LRTI Antibiotic Initiation: Cut-Offs (ng/mL) Used in Retained RCTs

| Study | Antibiotics Strongly Discouraged | Antibiotics Discouraged | Antibiotics Encouraged | Antibiotics Strongly Encouraged |
|----------------------|----------------------------------|-------------------------|------------------------|---------------------------------|
| Branche (2015) | ≤ 0.1 | < 0.25 | ≥ 0.25 | ≥ 0.5 |
| Briel (2008) | < 0.1 | ≤ 0.25 | > 0.25 | - |
| Burkhardt (2010) | - | < 0.25 | ≥ 0.25 | - |
| Christ-Crain (2004) | ≤ 0.1 | ≤ 0.25 | > 0.25 | ≥ 0.5 |
| Christ-Crain (2006) | < 0.1 | ≤ 0.25 | > 0.25 | > 0.5 |
| Corti (2016) | ≤ 0.15 | ≤ 0.25 | > 0.25 | - |
| Kristoffersen (2009) | - | < 0.25 | ≥ 0.25 | > 0.5 |
| Long (2009) | - | < 0.25 | ≥ 0.25 | - |
| Long (2011) | < 0.1 | ≤ 0.25 | > 0.25 | - |
| Schuetz (2009) | < 0.1 | ≤ 0.25 | > 0.25 | > 0.5 |
| Stolz (2007) | < 0.1 | ≤ 0.25 | > 0.25 | - |
| BMX Proposal | < 0.1 | ≤ 0.25 | > 0.25 | > 0.5 |

Funnel Plot – LRTI Mortality



Overlap Between Study-Level and Patient-Level MA for LRTI and Sepsis



Total: 23 unique RCTs - 7807 unique patients

VIDAS B.R.A.H.M.S PCT is a Precise Assay at the Cut-Offs

| Sample | Mean (ng/mL) | Repeatability CV* | Total Precision CV* |
|--------|--------------|-------------------|---------------------|
| 1 | 0.10 | 9.7% | 15.9% |
| 2 | 0.22 | 4.3% | 6.8% |
| 3 | 0.52 | 3.8% | 6.1% |
| 4 | 2.06 | 2.0% | 4.9% |

*Coefficient of Variation

Performed on VIDAS 3
According to CLSI Guidelines

Number of Studies Using VIDAS

- Study-level for LRTI (N=2 studies)
 - Corti, Branche: VIDAS exclusively

- Study-level for Sepsis: (N=3 studies)
 - Deliberato: VIDAS exclusively
 - De Jong, Shehabi : Multiple assays (VIDAS, Kryptor, Roche)

Concordance Study: VIDAS vs KRYPTOR

Overall Agreement Ranging 87-99%

| | Overall Agreement (95% CI) | Kappa |
|-------------------|---------------------------------------|---------------|
| 0.10 ng/mL | 86.7% (81.2, 91.0) | 0.7309 |
| 0.25 ng/mL | 97.5% (94.3, 99.2) | 0.9380 |
| 0.50 ng/mL | 99.0% (96.5, 99.9) | 0.9710 |
| 2.00 ng/mL | 97.5% (94.3, 99.2) | 0.8702 |

PCT and Renal Dysfunction

- PCT half-life-time is not significantly altered during renal dysfunction (26.1-33.1 h, vs. 22.3-28.9 h in normal population)¹
- Renal elimination is not a major mechanism for PCT removal from blood²
- Renal failure occurs often in sepsis but PCT algorithm still effective


1. Meisner et al., Intens Care Med 2000

2. Meisner et al., Eur J Anaesthesiol 2001


Reported Adherence in Study-Level Meta-Analyses

| Author | Adherence | |
|---------------------|-----------|----------------|
| LRTI | | |
| Kristoffersen | 59% | } Median = 82% |
| Corti | 61% | |
| Branche | 64% | |
| Christ-Crain (2004) | 78% | |
| Briel | 85% | |
| Burkhardt | 87% | |
| Christ-Crain (2006) | 87% | |
| Schuetz | 91% | |
| Sepsis | | |
| Bouadma | 47% | } Median = 72% |
| Annane | 63% | |
| Nobre | 81% | |
| de Jong | 93% | |

LRTI: Antibiotic Duration by PCT Adherence




LRTI: Mortality by PCT Adherence




Extensive Training and Materials will Guide Appropriate Use of PCT Assay

- Laboratory training
- Physician and other healthcare provider in-service education (ISE)
- Interpretation of results materials
- Physician to Physician education
- Webinars and other ISE as requested by each healthcare system
- ISE by Nurses trained in Critical Care medicine and PCT


LRTI: Total Antibiotic Exposure in Patient-Level Analyses by Type of LRTI Stratified by Inpatient vs Outpatient



LRTI: Mortality in Patient-Level Analyses by Type of LRTI Stratified by Inpatient vs Outpatient




Possible Harm of Over-Treatment




PCT Stability

- The sera or plasma separated from the clot can be stored at 2-8°C in stoppered tubes for up to 48 hours
- >48 hours longer storage is required, freeze at $-25 \pm 6^{\circ}\text{C}$
- Six-month storage of frozen samples does not affect the quality of results
- Three freeze/thaw cycles were validated.


LRTI - Antibiotic Exposure by Initial PCT Level: Patient-Level Meta-Analysis



Sepsis - Antibiotic Exposure by Initial PCT Level: Patient-Level Meta-Analysis




LRTI: Antibiotic Exposure/Duration (Including Pediatrics)



**PEDIATRIC DATA
NOT PREVIOUSLY REVIEWED BY FDA**

LRTI: Hospital Length of Stay (Including Pediatrics)



**PEDIATRIC DATA
NOT PREVIOUSLY REVIEWED BY FDA**