prismaflex

HF 20 SET


eng	Instructions for use	rus	Инструкция по применению
fre	Notice d'utilisation	hrv	Upute za uporabu
ger	Gebrauchsanweisung	lav	Lietošanas instrukcijas
spa	Folleto de utilización	cze	Návod k použití
ita	Istruzioni per l'uso	slv	Navodila za uporabo
swe	Handhavandeninstruktioner	srp	Uputstvo za upotrebu
dut	Gebruiksaanwijzing	slk	Návod na použitie
por	Instruções de use	hun	Használati utasítás
nor	Bruksanvisning	est	Kasutusjuhend
dan	Bruksanvisning	lit	Naudojimo instrukcijos
fin	Käyttöä opas	bul	Инструкции за употреба
pol	Instrukcja stosowania	rum	Instrucțiuni de utilizare
tur	Kullanma talimatları	vie	Hướng dẫn sử dụng
gre	Οδηγίες χρήσης	ukr	Інструкція з використання

REF 109841


Made in France

CE 2797

Date of Revision: 2019-12-01



Baxter

Blood warmer connection (blue)
Connexion réchauffeur sang (bleue)
Blutwärmeranschluss (blau)
Conexión del calentador de sangre (azul)
Collegamento del riscaldatore ematico (blu)
Anslutning för blodvärmare (blå)
Aansluiting bloedverwarmer (blauw)
Conexão para aquecedor de sangue (azul)
Blodvarmerkobling (blå)
Blodvarmertilslutning (blå)
Verenlämmittimen liitäntä (sininen)
Złącze podgrzewacza krwi (niebieskie)
Kan ısıtıcısı bağlantısı (mavi)
Σύνδεση θερμαντήρα αίματος (μπλε)
Соединение для нагревателя крови (синее)
Priključak grijača za krv (plavo)
Asins silditāja savienojums (zils)
Spojka ohřívače krve (modrá)
Priključek za grelnik krvi (moder)
Konekcija grejača krvi (plava)
Pripojenie ohrievača krvi (modré)
Vérmelegítő csatlakozása (kék)
Veresoojendaja ühendus (sinine)
Kraujo šildytuvo jungtis (mėlyna)
Връзка за нагревател на кръв (синя)
Conexiune încălzitor sânge (albastru)
Kết nối máu ấm hơn (màu xanh)
З'єднання нагрівача крові (синій колір)

eng	English	14
fre	Français	16
ger	Deutsch	18
spa	Español.....	20
ita	Italiano	22
swe	Svenska.....	24
dut	Nederlands	26
por	Português	28
nor	Norsk	30
dan	Dansk	32
fin	Suomi	34
pol	Polski.....	36
tur	Türkçe.....	38
gre	Ελληνικά	40
rus	Русский	42
hrv	Hrvatski	44
lav	Latviešu	46
cze	Český.....	48
slv	Slovenščina	50
srp	Srpski	52
slk	Slovenčina.....	54
hun	Hrvatski	56
est	Eesti.....	58
lit	Lietuvių k.	60
bul	български	62
rum	Română	64
vie	Tiếng Việt.....	66
ukr	Українська.....	68
	General characteristics.....	70
	Installation guide	79

DEFINITION OF SYMBOLS

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The Prismaflex HF20 set is manufactured by GAMBRO Industries,
7 Avenue Lionel Terray, BP 126, 69883 MEYZIEU CEDEX, FRANCE.

⚠ **Caution:** Federal law (USA) restricts this device to sale by or on the order of a physician.

DEFINITION OF EXPRESSIONS USED IN THIS MANUAL

In this document :

⚠ **"Warning"** indicates a hazardous situation which, if not avoided, could result in death or serious injury.

⚠ **"Caution"** indicates a hazardous situation which, if not avoided, could result in minor or moderate injury.

"Note" to give additional information.

SCUF: Slow Continuous UltraFiltration.

CVVH: Continuous Venous-Hemofiltration.

CVVHD: Continuous Venous-Hemodialysis.

CVVHDF: Continuous Venous-Hemodiafiltration.

Predilution: addition of replacement fluid to the blood stream upstream to the filter.

Postdilution: addition of replacement fluid to the blood stream downstream to the filter.

"Control unit" refers to the **PrismaFlex** control unit, or to the **PrisMax** control unit (in countries where **PrisMax** is cleared or registered).

PRODUCT DESCRIPTION

- The **Prismaflex** HF20 set is a disposable, extracorporeal circuit for use with the **PrismaFlex** system or with the **PrisMax** system (in countries where **PrisMax** is cleared or registered).
- The **Prismaflex** HF20 set consists of a PAES hollow fiber hemofilter/dialyzer* and tubing lines; refer to the control unit operator manual drawings for details.
- The **Prismaflex** HF20 set is compatible with the sleeve blood warmers, which should be installed on the return line. Refer to the specific instructions for use.
- All line connectors are compatible with the ISO594-1 and ISO594-2 international standards concerning conical fittings.
- The fluid pathways of the **Prismaflex** set are guaranteed sterile and non pyrogenic.
- The **Prismaflex** HF20 set is sterilized by ethylene oxide (EtO). Deaeration is such that EtO residuals comply with the standards in ISO 10993.
- Expiration date: please refer to product label.

* In this document the hemofilter/dialyzer will be referred to as "filter".

INTENDED USE / INDICATIONS

The **Prismaflex** set is indicated for use only with the **PrismaFlex** control unit or with the **PrisMax** control unit (in countries where **PrisMax** is cleared or registered) in providing continuous fluid management and renal replacement therapies. The system is intended for patients who have acute renal failure, fluid overload, or both.

This set is intended for use in the following veno-venous therapies: SCUF; CVVH; CVVHD; CVVHDF.

All treatments administered with the **Prismaflex** set must be prescribed by a physician. The size, weight, metabolic and fluid balance, cardiac status, and general clinical condition of the patient must be carefully evaluated by the prescribing physician before each treatment.

CONTRAINDICATIONS

There are no known absolute contraindications to continuous renal replacement therapies.

For the following conditions a careful assessment of the individual risk/benefit ratio has to be made by the treating physician (relative contraindications):

- inability to establish vascular access,
- severe hemodynamic instability,
- known hypersensitivity to any component of the **Prismaflex** set.

CAUTIONS AND WARNINGS

Note: refer to the control unit user interface and operator's manual for additional cautions and warnings.

⚠ Cautions

1. Particular attention must be paid to extra corporeal blood volume with respect to patient size. Consider the sum of the **Prismaflex** set blood volume (refer to "Specifications") plus the blood volume of any accessory or device if used. The **Prismaflex** HF20 set should be restricted to patients with a body weight greater than 8kg (18lb).
2. When not using the pre blood pump infusion line, it is recommended to clamp this line close to its connection to the access blood tubing line; this will prevent the sedimentation of blood into the pre-blood pump infusion line.
3. If citrated blood is used for the priming of the extracorporeal circuit before patient connection, it is recommended that the pH of this blood be verified and buffered at a value between 7.3 and 7.5. The ionized calcium should be adjusted at a value > or = 1.0.
4. If the patient is not immediately connected to the **Prismaflex** set after priming is complete, flush the set with at least 500 mL priming solution [saline or alkaline solution (pH ≥ 7.3), with or without heparin added according to usual institutional practice] prior to connecting the patient. This requires use of a new bag of priming solution.
5. Using the **Prismaflex** set with blood flow rates lower than the recommended minimum values (see "Operating Parameters" section) may impair filter performance due to hemoconcentration, or to increased risk of coagulation.
6. Since drugs can be removed by the membrane of the filter, the dosage of associated drug treatments may need to be adjusted for patients on continuous renal replacement therapy. Monitoring of blood drug levels of relevant compounds should be performed. The removal of other water-soluble compounds (e.g. vitamins, trace elements) during therapy also requires clinical consideration.

⚠ Warnings

1. Carefully read these instructions for use and the control unit operator's manual before using this product.
2. The use of operating procedures other than those published by the manufacturer or the use of accessory devices not recommended by the manufacturer can result in patient injury or death.
3. Store the **Prismaflex** set in a dry place, between 0° C (32° F) and 30° C (86° F).
4. Do not use this set if the packaging is damaged, if the sterilization caps are missing or loose, or if any of the lines in the set are kinked.
5. To prevent contamination, this **Prismaflex** set must be used as soon as its packaging and sterilization caps are removed.
6. Do not try to remove the filter from the cartridge plate.
7. Use aseptic techniques when handling all blood and fluid lines in the set.
8. **Prismaflex** sets are compatible with the usual disinfection agents used for aseptic setup; however solvents and other chemicals, if used in contact with the product, could damage the set.
9. During priming and operation, observe closely for leakage at joints within the set, and connections to other approved accessories and bags. Leakage can cause blood loss, fluid imbalance or air embolism. If a leakage is detected at a Luer connection and cannot be stopped by tightening the connections, or if leakage occurs at any other location, replace the set.
10. Tightening Luer connections with an excessive force can damage the connectors.
11. In the case of patients with a high risk of hemorrhage it is recommended not to add heparin to the priming solution.
12. Do not allow air to enter the blood compartment of the filter after priming is started. If a large amount of air enters, the set must be replaced.
13. Should acute allergic reactions (first-use syndrome) occur in patients receiving treatment, **immediately stop the treatment** and administer appropriate intervention.
14. Use a 21-gauge or smaller needle to obtain blood/fluid samples or remove trapped air from the **Prismaflex** set. Use of larger needles can cause holes in the sample sites, resulting in external leak or air intake.
15. External blood leakage may not be immediately identified by monitoring equipment and could result in significant blood loss. Check the filter and all connections of the disposable tubings during treatment to minimize the risk of leakage.
16. To assure adequate filter performance, it is recommended that the set be changed every 24 hours of use. However, the set must be changed after 3 days (72 hours). Continued use beyond this limit could result in rupture of the pump segments, with risk of patient injury or death.

17. Destroy this set after single use, using aseptic technique for potentially contaminated equipment and following local regulation for disposal. Do not re-sterilize. The **Prismaflex** set is intended for single use only. Re-using the **Prismaflex** set may cause serious damage to the product resulting in patient injury or death.
18. Use only drugs compatible with plastics listed in the specifications section. Some plastics can be incompatible with drugs when in contact with solutions with pH > 10.

SPECIFICATIONS

See Tables at the end of the document.

SET MATERIALS

PAES hollow fiber	: PolyarylEtherSulfone
Housing and headers	: Polycarbonate
Potting compound	: Polyurethane
Tubing material	: Plasticized polyvinyl chloride (PVC)
Cartridge	: Polyethyleneterephthalate Glycol

Note: the following information is available from the manufacturer upon request:

- information about test methods used to obtain performance characteristics,
- the number and range of particles in the effluent from the dialyzer prepared as recommended for clinical use,
- the types and amounts of residue from the sterilization process.

Note: the **Prismaflex** set is not made with natural rubber latex.

Note: all fluid pathways in direct or indirect blood contact are DEHP-free.

INSTRUCTIONS FOR USE

Note: use the set by following the detailed on-line instructions provided by the control unit. Additional information is available in the control unit operator's manual.

Load Set

Install the set onto the control unit using the photographs on the inside cover as a guide - the same procedure applies for both **PrismaFlex** and **PrisMax** control units (in countries where **PrisMax** is cleared or registered).

Prepare and Connect Solutions

Hang bag of priming solution saline or alkaline solution (pH ≥ 7.3) with or without heparin added (5000IU heparin/liter) according to usual institutional practice, on priming hook. Connect access (red)/effluent (yellow) Y-line to priming solution bag.

Anticoagulation considerations

According to the literature* continuous heparinization at a rate of 10 to 20.IU/kg/h ensures proper operation of the extracorporeal circuit when performing treatment with patients having normal coagulation status.

Depending on the patient's condition, however, heparinization can be lowered to less than 5 IU/kg/h.

Heparinization can be controlled, for instance, by partial thromboplastin time (PTT) measurements: in this case, PTT could be maintained at 20 to 30 seconds over baseline.

* "Continuous renal replacement therapy in critically ill neonates"; G.ZOBEL & Al., *Kidney International, Vol.53, Suppl. 66 (1998)*.

SPECIAL PROCEDURES IN CASE OF COMPLICATION

External Blood Leaks

Note: see Warning no. 15.

If an external blood leakage is observed, immediately stop the blood pump. Initiate corrective action by securing connections or replacing the **Prismaflex** set.

If necessary, administer adequate replacement solution to the patient to compensate for blood loss.

Hypersensitivity Reactions

Note: see Warning no. 13.

Should acute allergic reactions (first use syndrome) occur within the first few minutes of the treatment, it is important to react immediately by discontinuing the session and administering appropriate treatment.

Adverse reactions may occur due to the complex interaction between blood and the artificial surfaces of the entire extracorporeal circuit. These reactions may also be precipitated and/or exacerbated by other external factors involved with the individual patient's specific disease process and the treatment of renal insufficiency. Certain types of adverse reaction may occur due to operational factors associated with the treatment. Therefore, proper management of the fluid removal, electrolyte balance, anticoagulation and blood flow rate as well as monitoring of the overall treatment parameters are essential to avoid side-effects which may be associated with hemodialysis/hemofiltration therapies.

Hypersensitivity reactions have been observed during dialysis. Symptoms of a hypersensitivity reaction may be gastrointestinal, mucocutaneous, respiratory, cardiovascular or systemic in nature and range from very mild to severe. Such symptoms have been described as anaphylactic-like reactions within the first few minutes. Manifestations include nausea, malaise, weakness, a sensation of burning or heat throughout the body, profuse perspiration, respiratory distress and in some instances hypotension and cardiopulmonary arrest. Should a combination of such symptoms appear, particularly at the start of the treatment session, it is important to react immediately by discontinuing the session and administering appropriate treatment. Blood in the extracorporeal circuit must not be returned to the patient.

Extra care must be taken when treating patients who have exhibited possible hypersensitivity symptoms during previous treatments, or patients who have a history of being highly sensitive and allergic to a variety of substances. A physician must be consulted to evaluate the risk and prescribe the appropriate precautions if a possible sensitivity is suspected.

The following factors are considered essential to minimize the risk of hypersensitivity reaction and other side effects:

- Strict adherence to the set-up, priming and rinsing procedures detailed in the manufacturer's instructions for use.
- Setting up and monitoring the treatment operating parameters according to the manufacturer's recommendations specified for each type of **Prismaflex** set and to the patient's needs and tolerance.
- Strict adherence to all WARNINGS and CAUTIONS given by the manufacturer in the instructions for use.

WARRANTY AND LIMITATION OF LIABILITY

a) The manufacturer warrants that the **Prismaflex** set has been manufactured in accordance with its specifications and in compliance with good manufacturing practices, other applicable industry standards and regulatory requirements.

If provided with the lot/serial number of the defective product, the manufacturer will, by replacement or credit, remedy manufacturing defects in the **Prismaflex** set becoming apparent before the expiration date.

b) The warranty under paragraph a) above is in lieu of, and to the exclusion of, any other warranty, whether written or oral, express or implied, statutory or otherwise, and there are no warranties of merchantability or other warranties, which extend beyond those described in paragraph a) above. The remedy set out above for manufacturing defects is the sole remedy available to any person due to defects in the **Prismaflex** set and the manufacturer shall not be liable for any consequential or incidental loss, damage, injury or expense arising directly or indirectly from the use of the **Prismaflex** set, whether as a result of any defect therein or otherwise.

c) The manufacturer shall not be liable for any misuse, improper handling, non-compliance with warnings and instructions, damage arising from events after the manufacturer's release of the **Prismaflex** set, failure or omission to inspect the **Prismaflex** set before use in order to ensure that the **Prismaflex** set is in proper condition, or any warranty given by independent distributors or dealers.

d) The manufacturer is GAMBRO Industries, 7 avenue Lionel Terray, BP 126, 69883 MEYZIEU CEDEX, FRANCE.

MALAYSIA:

Authorized Representative: Baxter Healthcare (Malaysia) Sdn. Bhd.
B-21-3A, Level 21, The Ascent,
Paradigm, No. 1, Jalan SS7/26A,
Kelana Jaya, 47301 Petaling Jaya,
Selangor Darul Ehsan,
Malaysia Registration number: GC25696881418

INDONESIA:

Distributor : PT. Tawada Healthcare
Rukan Permata Senayan Blok A No. 18-19, Jl. Tentara
Pelajar No. 5 Grogol Utara, Jakarta Selatan - Indonesia
KEMENKES RI AKL 20805714921

PHILIPPINES:

Imported and distributed by Baxter Healthcare Philippines, Inc.
Silangan Industrial Estate, Brgy. Canlubang, Calamba City, Laguna, Philippines
MDR-01243



Prismaflex HF20 SET	
PHYSICAL CHARACTERISTICS ⁽¹⁾	
Membrane effective surface area	0.2 m ²
Fiber internal diameter (wet)	215 µm
Fiber wall thickness	50 µm
Blood volume in set	58 mL
Overall dimensions	
• Length	27 cm
• Width	22 cm
• Height	9 cm
Weight	677 g
OPERATING PARAMETERS	
Maximum TMP	500 mm Hg 66.6 kPa
Maximum blood pressure	500 mm Hg 66.6 kPa
Minimum blood flow rate	20 mL/min
Maximum blood flow rate	100 mL/min
PERFORMANCE SPECIFICATIONS ⁽²⁾	
Maximum ultrafiltration rate (mL/min) ⁽³⁾ (bovine blood, Hct 32%, Cp 60 g/L, 37°C)	
QB (mL/min)	20 50 100
Max.QUF (± 20%)	9 17 24
Sieving coefficient (bovine plasma, Cp 60 g/L, 37°C) QB = 50 mL/min, QUF = 10 mL/min	
• Urea	1.00
• Vitamin B12	1.00
• Inulin	0.92
• Albumin	<0.01
Clearance (mL/min) (saline solution; 37°C) Parameters: QB/QS QUF QD (L/h) QD (mL/min)	
	20 mL/min 0 mL/min
	50 mL/min 0 mL/min
	100 mL/min 0 mL/min
	0.5 1 1.5 2 2.5 0.5 1 1.5 2 2.5 0.5 1 1.5 2 2.5 8 17 25 33 42 8 17 25 33 42 8 17 25 33 42
Urea (± 10%)	8 12.5 14.5 15.7 16.5 8.3 15.5 20.9 25.1 28.7 8.3 16.4 23.6 30.1 36.3
Creatinin (± 10%)	7.8 11.7 13.6 14.7 15.6 8.3 15.0 19.7 23.1 26.0 8.3 16.2 22.9 28.5 33.5
Vitamin B12 (± 20%)	6.5 8.6 9.5 10.0 10.3 7.6 11.7 13.9 15.4 16.5 8.1 13.9 17.6 20.3 22.5
Inulin (± 20%)	5.6 7.0 7.5 7.8 7.9 6.9 9.8 11.2 12.1 12.7 7.9 12.6 15.2 16.9 18.1




⁽¹⁾ Nominal values - given for indication

⁽²⁾ Typical mean values obtained from laboratory testing of post-sterilization sample lots.
Results may vary depending on patient and clinical conditions.

⁽³⁾ Ultrafiltration is controlled by the control unit and is independent of the ultrafiltration coefficient (KUF).

INSTALL SET DIRECTLY ONTO THE CONTROL UNIT

Patented ready to use PrismaFlex set packaging.

BAXTER CONFIDENTIAL - INTERNAL USE ONLY

Part Number: 1000014767	Date: 08-JAN-2020	Proofread No.: 1
Designer: CDS	Page: 1 of 1	
Colour Reference: BLACK		

prismaflex

◆ **EXTRA CORPOREAL CIRCUIT WITH FILTER** / Blood access line (red-striped) / Blood return line (blue-striped) / Dialysate line (green-striped) / Effluent line (yellow-striped) / Replacement solution line (purple-striped) / Pre-blood pump line (white-striped) / 5-liter effluent collection bag. Sterile and non-pyrogenic fluid pathways. **Caution:** Do not use if protective caps are not in place or if packaging is damaged. **Caution:** Federal law (USA) restricts this device to sale by or on the order of a physician.

Indonesia - Distributor: PT. Tawada Healthcare, Rukan Permata Senayan Blok A No. 18-19, Jl. Tentara Pelajar No. 5 Grogol Utara, Jakarta Selatan - Indonesia. KE-MENKES RI AKL 20805711114 (Prismaflex M sets) / 20805714919 (Prismaflex HF1000/ HF1400 sets) / 20805714921 (Prismaflex HF20 set) / 20805713569 (Prismaflex ST sets).
Philippines - Imported and distributed by Baxter Healthcare Philippines, Inc. Silangan Industrial Estate, Brgy. Canlubang, Calamba City, Laguna, Philippines. MDR-01243 (Prismaflex HF20 set) / MDR-07338 (Prismaflex M sets).
Malaysia - Authorised Representative Baxter Healthcare (Malaysia) Sdn. Bhd. B-21-3A, Level 21, The Ascent, Paradigm, No. 1, Jalan SS7/26A, Kelana Jaya, 47301 Petaling Jaya, Selangor Darul Ehsan, Malaysia. Registration number: GC99624912318 (Prismaflex M sets) / GC2566881418 (Prismaflex HF20 set).

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Baxter and Prismaflex are trademarks of Baxter International Inc. or its subsidiaries.