

UNITED STATES OF AMERICA  
DEPARTMENT OF HEALTH AND HUMAN SERVICES  
FOOD AND DRUG ADMINISTRATION

+ + +

CENTER FOR DEVICES AND RADIOLOGICAL HEALTH  
MEDICAL DEVICES ADVISORY COMMITTEE

+ + +

CIRCULATORY DEVICES PANEL

+ + +

October 27, 2020  
9:00 a.m.

Via Zoom

PANEL MEMBERS:

RICHARD A. LANGE, M.D., M.B.A.	Chair
GEORGE W. VETROVEC, M.D., MACC, MSCAI	Voting Member
RANDALL C. STARLING, M.D., M.P.H.	Voting Member
JASON T. CONNOR, Ph.D.	Voting Member
RALPH G. BRINDIS, M.D., M.P.H., MACC, FSCAI	Voting Member
JANET WITTES, Ph.D.	Temporary Voting Member
PRAMOD BONDE, M.D.	Temporary Voting Member
RICHARD PAGE, M.D.	Temporary Voting Member
DAVID YUH, M.D.	Temporary Voting Member
ROBERT W. YEH, M.D., M.Sc, M.B.A.	Temporary Voting Member
KEITH B. ALLEN, M.D.	Temporary Voting Member
WAYNE BATCHELOR, M.D.	Temporary Voting Member
ERIK MAGNUS OHMAN, M.D., FACC	Temporary Voting Member
JOAQUIN E. CIGARROA, M.D.	Temporary Voting Member
JOHN HIRSHFELD, M.D.	Temporary Voting Member
VERGHESE MATHEW, M.D., FACC	Temporary Voting Member
BERNARD GERSH, M.D.	Temporary Voting Member
JOHN SOMBERG, M.D.	Temporary Voting Member
JEFFREY BORER, M.D.	Temporary Voting Member

Free State Reporting, Inc.  
1378 Cape Saint Claire Road  
Annapolis, MD 21409  
(410) 974-0947

JACQUELINE ALIKHAANI  
DEBRA DUNN

Consumer Representative  
Patient Representative

ADEN ASEFA, M.P.H.

Designated Federal Officer

Free State Reporting, Inc.  
1378 Cape Saint Claire Road  
Annapolis, MD 21409  
(410) 974-0947

## FDA REPRESENTATIVES:

BRAM ZUCKERMAN, M.D.  
Director, Office of Health Technology 2 (OHT 2: Cardiovascular Devices)  
Office of Product Evaluation and Quality

LYDIA GLAW, Ph.D.  
Assistant Director, Coronary Interventional Devices  
Division of Health Technology 2 C (Coronary and Peripheral Intervention Devices)  
Office of Health Technology 2 (OHT 2: Cardiovascular Devices)  
Office of Product Evaluation and Quality

JAMES SWINK  
Medical Devices Advisory Committee Panel Coordinator

LINDSEY O' KEEFE  
Press Contact

## FDA PRESENTERS:

SAMUEL RABEN, Ph.D.  
Mechanical Engineer/Team Leader, Coronary Interventional Devices Team  
Office of Health Technology 2 (OHT 2: Cardiovascular Devices)  
Division of Health Technology 2 C (Coronary and Peripheral Intervention Devices)  
Office of Product Evaluation and Quality

ANNABELLE CRUSAN, D.V.M., M.S.  
Animal Study Reviewer  
Office of Health Technology 2 (OHT 2: Cardiovascular Devices)  
Office of Product Evaluation and Quality

RONG (RONA) TANG, Ph.D.  
Statistical Reviewer  
Division of Clinical Evidence and Analysis 2 (Biostatistics)  
Office of Clinical Evidence and Analysis

TARA A. RYAN, M.D.  
Medical Officer  
Office of Health Technology 2 (OHT 2: Cardiovascular Devices)  
Office of Product Evaluation and Quality

## SPONSOR PRESENTERS:

VICKI BEBEAU  
Vice President, Clinical and Regulatory Affairs  
Neovasc, Inc.

TIMOTHY D. HENRY, M.D.  
Medical Director  
Carl and Edyth Lindner Center for Research  
The Christ Hospital, Cincinnati, Ohio

SHMUEL BANAI, M.D.  
Director, Division of Cardiology  
Tel Aviv Medical Center  
Professor of Cardiology/Head, Center for Cardiac Research  
Tel Aviv University  
Medical Director  
Neovasc, Inc.

VICTORIA HAMPSHIRE, V.M.D.  
Capital Preclinical Scientific Research Consultants, LLC

SERGE ROUSSELLE, D.V.M., DACVP  
Director, Imaging and Digital Pathology  
StageBio

GREGG W. STONE, M.D., FACC, MSCAI  
Director, Academic Affairs  
Mount Sinai Heart Health System  
Professor of Medicine (Cardiology) and Population Health Sciences Policy  
The Zena and Michael A. Wiener Cardiovascular Institute  
Icahn School of Medicine at Mount Sinai

## SPONSOR ADVISORS:

CHRIS MULLIN, M.S.  
Biostatistician  
NAMSA

## OPEN PUBLIC HEARING SPEAKERS:

MEG SEYMOUR, Ph.D.  
Senior Fellow  
National Center for Health Research

AMIR LERMAN, M.D.  
Department of Cardiovascular Disease  
Mayo Clinic

GERALD KOENING, M.D., Ph.D.  
Henry Ford Health System

RYAN GINDI, M.D.  
Henry Ford Health System

FREDERICK CASCIANO  
Patient

ANNETTE CASCIANO  
Wife of Patient

MARK SOBERANO  
Patient

STEVEN SUMMERS  
Patient

DONNA SUMMERS  
Wife of Patient

LAURIE VANDENBOSSCHE  
Patient

DONALD SCOTT  
Patient

TAMMY HOPKINS  
Patient

CLYDE HART  
Patient

## INDEX

	PAGE
CALL TO ORDER - Richard A. Lange, M.D., M.B.A.	9
PANEL INTRODUCTIONS	10
CONFLICT OF INTEREST STATEMENT - Aden Asefa, M.P.H.	15
TEMPORARY VOTING MEMBER STATUS STATEMENT - Aden Asefa, M.P.H.	17
GENERAL ANNOUNCEMENTS - Aden Asefa, M.P.H.	17
SPONSOR PRESENTATION	
Introduction - Vicki Bebeau	18
Unmet Need - Timothy D. Henry, M.D.	24
Reducer Mechanism of Action - Shmuel Banai, M.D.	31
Animal Studies - Victoria Hampshire, V.M.D.	33
Midsection Reducer Endothelialization - Serge Rousselle, D.V.M., DACVP	37
Study Design and Results - Gregg W. Stone, M.D., FACC, FSCAI	38
Clinical Perspective - Timothy D. Henry, M.D.	55
PANEL QUESTIONS FOR SPONSOR	58
FDA PRESENTATION	
Regulatory History and Considerations, Device Description, Principle Of Operation, and Nonclinical Evaluations - Samuel Raben, Ph.D.	74
Animal Studies - Annabelle Crusan, D.V.M., M.S.	79
Study Design and Statistical Analysis Plan - Rong (Rona) Tang, Ph.D.	83
COSIRA Clinical Study Results - Tara A. Ryan, M.D.	85
Post-Approval Study and Summary - Samuel Raben, Ph.D.	97
PANEL QUESTIONS FOR FDA	99

## INDEX

	PAGE
OPEN PUBLIC HEARING	
Meg Seymour, Ph.D.	115
Amir Lerman, M.D.	117
Gerald Koenig, M.D., Ph.D.	118
Ryan Gindi, M.D.	121
Frederick Casciano	122
Annette Casciano	123
Mark Soberano	124
Steven Summers	125
Donna Summers	126
Laurie Vandebossche	127
Donald Scott	129
Tammy Hopkins	129
Clyde Hart	130
SPONSOR RESPONSE TO PANEL QUESTIONS	132
FDA RESPONSE TO PANEL QUESTIONS	179
PANEL DELIBERATIONS/FDA QUESTIONS	
Question 1	184
Question 2	189
Question 2a	189
Question 2b	194
Question 3	198
Question 3a	198
Question 3b	204

## INDEX

	PAGE
PANEL DELIBERATIONS/FDA QUESTIONS (cont.)	
Question 4	208
Question 4a	208
Question 4b	212
Question 5	214
Question 5a	215
Question 5b	229
Question 5c	233
Question 6	236
Question 7	241
SUMMATIONS	
FDA - Samuel Raben, Ph.D.	243
Sponsor - Vicki Bebeau, Gregg W. Stone, M.D., FACC, FSCAI and Timothy D. Henry, M.D.	244
FINAL COMMENTS	
Jacqueline Alikhaani, Consumer Representative	248
PANEL VOTE	248
ADJOURNMENT	269



MEETING

(9:00 a.m.)

1  
2  
3 DR. LANGE: Great, it's 9:00 a.m., October 27th. I'd like to call this meeting of the  
4 Circulatory Devices Panel to order. I'm Richard Lange, the chairperson of this Panel. I'm  
5 president of Texas Tech University Health Sciences Center in El Paso, where I'm also dean of  
6 the Paul L. Foster School of Medicine, and previously was an interventional cardiologist.

7 I want to thank everybody for participating on this Panel. These are obviously  
8 challenging times and extenuating circumstances. I want to express my appreciation to the  
9 FDA and our audiovisual personnel who have worked beforehand to try to make this go as  
10 seamlessly as possible.

11 I do want to comment, though. If you have any difficulties, we shared three e-mails  
12 with you, you can e-mail Aden, you can e-mail Jim, or you can e-mail me. Now, if you e-mail  
13 me for audiovisual help, I just want to remind you that, between the third and fifth grade, I  
14 did replace a bulb in the carousel and ran the 8 mm film. So I'll be your last resort.

15 Again, I want to thank everybody for participating and especially those on the West  
16 Coast, Pacific time, where it's 6:00 a.m. So thank you, everybody.

17 I note for the record that the voting members present constitute a quorum as  
18 required by 21 C.F.R. Part 14. I would also like to add that the Panel members participating  
19 in today's meeting have received training in FDA device law and regulations.

20 For today's agenda, the Panel will discuss, make recommendations, and vote on  
21 information regarding the premarket approval application, that is the PMA, for the Neovasc  
22 Reducer System sponsored by Neovasc, Inc.

23 Before we begin, I would like to ask our distinguished Committee members and FDA  
24 attending virtually to introduce themselves. Committee members, to do so, I'll remind you  
25 to please turn on your video monitors if you have not already done so, and you need to

Free State Reporting, Inc.  
1378 Cape Saint Claire Road  
Annapolis, MD 21409  
(410) 974-0947

1 unmute your phone before you speak. I will call your name. At that time if you'll please  
2 state your area of expertise, your position, and your affiliation.

3 And I'll start with Dr. George Vetrovec. Please unmute yourself.

4 DR. VETROVEC: -- after a long interventional career, beginning with the first stent in  
5 1970 -- first PCI in 1977 -- committees at the University and do some teaching as well as  
6 writing and information for the ACC.org. I'm excited about today's meeting.

7 DR. LANGE: Thank you for joining us, George.

8 Dr. Randall Starling.

9 DR. STARLING: Good morning, my name is Randy Starling. I have, of course, an M.D.  
10 and master's in public health and epidemiology. My specialty in cardiology is heart failure  
11 and transplant. I've been at the Cleveland Clinic for 25 years. I'm Professor of Medicine at  
12 the Cleveland Clinic Lerner College of Medicine of Case Western Reserve University. I have  
13 served in an ad hoc basis on previous FDA panels and have recently become a member of  
14 this particular circulatory assist devices panel. So I'm looking forward to working with this  
15 esteemed group of colleagues.

16 DR. LANGE: Thank you, Dr. Starling.

17 Dr. Jason Connor.

18 DR. CONNOR: Hi, I'm Jason Connor, founder and biostatistician at ConfluenceStat  
19 and Assistant Professor of Medical Education at the University of Central Florida College of  
20 Medicine.

21 DR. LANGE: Thank you, Jason.

22 Dr. Ralph Brindis.

23 DR. BRINDIS: Good morning. Ralph Brindis. I'm a Clinical Professor of Medicine at  
24 UCSF at the Philip R. Lee Institute of Health Policy Studies. I'm a cardiologist by training,  
25 experience in cardiovascular outcomes research, and also served as the chief medical

1 officer/external affairs at the National Cardiovascular Data Registry.

2 DR. LANGE: Thank you, Ralph.

3 DR. BRINDIS: And I have no new conflicts to disclose.

4 DR. LANGE: Great, thank you.

5 Dr. Janet Wittes.

6 DR. WITTES: Hi, I'm Janet Wittes. I'm a statistician and president of Statistics  
7 Collaborative, and I was a member of this Panel about 15 years ago, so it's good to be back  
8 and I'm looking forward to today.

9 DR. LANGE: Thank you, Janet, for joining us.

10 Dr. Pramod Bonde.

11 DR. BONDE: Hi, I'm Pramod Bonde. I'm Associate Professor of Surgery at Yale  
12 University. I specialize in adult cardiac surgery with a subspecialty in mechanical circulatory  
13 support.

14 DR. LANGE: Thank you, Pramod.

15 Dr. Richard Page.

16 DR. PAGE: Good morning, my name is Rick Page. I'm Professor of Medicine and  
17 dean of the Larner College of Medicine at the University of Vermont. I am a practicing  
18 cardiac electrophysiologist, although I'm noninvasive these days, and I'm the past chair of  
19 this Panel.

20 DR. LANGE: Thank you, Rick, good to have you on board again.

21 Dr. David Yuh.

22 DR. YUH: Good morning. David Yuh. I'm the Chairman of Surgery at Stamford  
23 Hospital in Connecticut, with an academic appointment at Columbia. I'm a clinically active  
24 adult cardiac surgeon primarily in coronary bypass surgery and minimally invasive valve  
25 repair.

1 DR. LANGE: Thank you, David, for joining us.

2 DR. YUH: You're welcome.

3 DR. LANGE: Dr. Robert Yeh.

4 DR. YEH: Hi, everyone. Robert Yeh. I'm an Associate Professor of Medicine at  
5 Harvard Medical School. I direct our complex coronary intervention program here at Beth  
6 Israel Deaconess Medical Center, and I direct the Smith Center for Outcomes Research in  
7 Cardiology also at the Beth Israel.

8 DR. LANGE: Thank you, Robert.

9 Dr. Keith Allen.

10 DR. ALLEN: Hi there, good morning. Keith Allen. I am at the Mid America Heart  
11 Institute in Kansas City, Missouri. I am a cardiac as well as vascular surgeon, director of  
12 surgical research, and the surgical director of our structural heart program.

13 DR. LANGE: Great. Thank you, Keith, for joining us.

14 Dr. Wayne Batchelor.

15 (No response.)

16 DR. LANGE: All right.

17 Dr. Magnus Ohman.

18 DR. OHMAN: Good morning, everybody. My name is Magnus Ohman, I'm a  
19 Professor of Medicine at Duke University in North Carolina, and I am an interventional  
20 cardiologist. I'm also the director of the Duke program for advanced coronary disease, so  
21 I'm very familiar with these types of patients. I'm looking forward to our discussion today.

22 DR. LANGE: Magnus, thanks for joining us.

23 Dr. Joaquin Cigarroa.

24 DR. CIGARROA: Good morning, I'm Joaquin Cigarroa and I'm Professor of Medicine  
25 at OHSU, clinical chief of the Knight Cardiovascular Institute, and division head of

1 cardiology. I'm a general cardiologist with added training in capabilities in interventional  
2 cardiology.

3 DR. LANGE: Great. Thank you, Joaquin.

4 Dr. John Hirshfeld. John, I'll ask you to unmute yourself.

5 DR. HIRSHFELD: Good morning, I'm John Hirshfeld. I'm a Professor of Medicine at  
6 the University of Pennsylvania School of Medicine. I've been a member of this Panel for a  
7 number of years and chaired it. I'm an interventional cardiologist who practiced  
8 interventional cardiology for 45 years and stopped actually doing procedures 2 years ago.

9 DR. LANGE: Great. John, thanks for joining us.

10 Dr. Verghese Mathew.

11 DR. MATHEW: Good morning, my name is Verghese Mathew and I'm the director of  
12 the Worldwide Network for Innovation in Clinical Education Research, or WNICER, which is  
13 a cardiovascular think tank and research collaborative; also continuing to be a practicing  
14 interventional cardiologist in a regional healthcare system in northeast Wisconsin now, and  
15 spent the last 22 years in academic cardiology including most recently as chair of cardiology  
16 at Loyola University in Illinois, as well as the Mayo Clinic in Minnesota. Thank you.

17 DR. LANGE: Thanks, Verghese.

18 Dr. Bernie Gersh. Bernie, you have to unmute. There you go.

19 DR. GERSH: Good morning, I'm Bernard Gersh. I'm a cardiologist with a wide range  
20 of interests at the Mayo Clinic. I'm Professor of Medicine at the Mayo Clinic College of  
21 Medicine, and this is my first meeting of this group, I'm looking forward to it very much.

22 DR. LANGE: Great. Thank you, Bernie, for joining us.

23 Dr. John Somberg.

24 DR. SOMBERG: Good morning, I'm John Somberg. I'm a Professor Emeritus of  
25 Cardiology, Pharmacology, and Medicine at Rush University. I'm also the editor of

1 *Cardiology Research* and practice cardiology in Lake Bluff, Illinois at the moment, and I have  
2 been on this Panel for over 35 years. So that's a long time and I'm looking forward to our  
3 discussions today.

4 DR. LANGE: Thanks for serving again.

5 Jeff Borer. Jeff, I don't see your picture up here, are you on?

6 (No response.)

7 DR. LANGE: All right, we'll have to see if we can get you back on.

8 Our Consumer Representative, Jacqueline Alikhaani. And, Jacqueline, you'll need to  
9 unmute.

10 MS. ALIKHAANI: Good morning, I'm Jacqueline Alikhaani and I am a heart survivor  
11 and volunteer patient advocate for the American Heart Association here in Los Angeles and  
12 for across the country, and I'm also an ambassador for the Patient-Centered Outcomes  
13 Research Institute, and I also volunteer with a number of other community healthcare  
14 stakeholder organizations like AARP and others.

15 DR. LANGE: Jacqueline, thank you for serving today.

16 Patient Representative Debra Dunn. Debra, are you on?

17 MS. ASEFA: Debra just e-mailed me and she's actually on her way to her office, so  
18 she should pop in soon, hopefully.

19 DR. LANGE: Thank you very much, Aden.

20 Is Gary Jarvis on?

21 MS. ASEFA: Gary Jarvis is also experiencing some IT issues, so he will not be able to  
22 be on until maybe 11:00 a.m. this morning.

23 DR. LANGE: Okay. Thank you, Aden.

24 From the FDA, Dr. Bram Zuckerman.

25 DR. ZUCKERMAN: Good morning, my name is Bram Zuckerman. I'm the director,

Free State Reporting, Inc.  
1378 Cape Saint Claire Road  
Annapolis, MD 21409  
(410) 974-0947

1 FDA Office of Cardiovascular Devices and a cardiologist by background. Thank you.

2 DR. LANGE: Thank you, Bram.

3 And then, Aden, if you'll introduce yourself.

4 MS. ASEFA: Hi, I'm Aden. I'm the DFO for this Committee, so if you have any  
5 questions, please let me know. Thank you.

6 DR. LANGE: All right. Probably the last remark I need to make is -- and I'm sorry I'm  
7 not wearing my coat. Besides being the area most affected by COVID in the United States  
8 right now, we had an ice storm last night, so my suit coat is sitting in the kitchen and my  
9 ex-arcuric (ph.) coat is right here. I put on a hoodie, but I've just decided to go coatless, so  
10 my apologies. It's no disrespect to the FDA or to the Sponsor, either one.

11 So with that, Aden Asefa, the Designated Federal Officer for today's Circulatory  
12 Devices Panel, will make some introductory remarks. Take it away, Aden.

13 MS. ASEFA: Good Morning. I will now read the Conflict of Interest Statement.

14 The Food and Drug Administration is convening today's meeting of the Circulatory  
15 Devices Panel of the Medical Devices Advisory Committee under the authority of the  
16 Federal Advisory Committee Act (FACA) of 1972. With the exception of the industry  
17 representative, all members and consultants of the Panel are special Government  
18 employees or regular Federal employees from other agencies and are subject to Federal  
19 conflict of interest laws and regulations.

20 The following information on the status of this Panel's compliance with Federal  
21 ethics and conflict of interest laws covered by, but not limited to, those found at 18 U.S.C.  
22 Section 208 are being provided to participants in today's meeting and to the public.

23 FDA has determined that members and consultants of this Panel are in compliance  
24 with Federal ethics and conflict of interest laws. Under 18 U.S.C. Section 208, Congress has  
25 authorized FDA to grant waivers to special Government employees and regular Federal

Free State Reporting, Inc.  
1378 Cape Saint Claire Road  
Annapolis, MD 21409  
(410) 974-0947

1 employees who have financial conflicts when it is determined that the Agency's need for a  
2 particular individual's services outweighs his or her potential financial conflict of interest.

3 Related to the discussions of today's meeting, members and consultants of this Panel  
4 who are special Government employees or regular Federal employees have been screened  
5 for potential financial conflicts of interest of their own as well as those imputed to them,  
6 including those of their spouses or minor children and, for purposes of 18 U.S.C. Section  
7 208, their employers. These interests may include investments; consulting; expert witness  
8 testimony; contracts/grants/CRADAs; teaching/speaking/writing; patents and royalties; and  
9 primary employment.

10 For today's agenda, the Panel will discuss, make recommendations, and vote on  
11 information regarding the premarket application (PMA) for Neovasc Reducer System, by  
12 Neovasc, Inc., indicated for patients suffering from refractory angina pectoris despite  
13 guideline-directed medical therapy, who are unsuitable for revascularization by coronary  
14 artery bypass grafting or by percutaneous coronary intervention.

15 Based on the agenda for today's meeting and all financial interests reported by the  
16 Panel members and consultants, no conflict of interest waivers have been issued in  
17 accordance with 18 U.S.C. Section 208.

18 Mr. Gary Jarvis is serving as the Industry Representative, acting on behalf of all  
19 related industry. He is employed by Alpha Medical.

20 We would like to remind members and consultants that if the discussions involve any  
21 other products or firms not already on the agenda for which an FDA participant has a  
22 personal or imputed financial interest, the participants need to exclude themselves from  
23 such involvement and their exclusions will be noted for the record.

24 FDA encourages all other participants to advise the Panel of any financial  
25 relationships that they may have with any firms at issue.

Free State Reporting, Inc.  
1378 Cape Saint Claire Road  
Annapolis, MD 21409  
(410) 974-0947



1 A copy of this statement will be available for review and included as part of the  
2 official transcript. Thank you.

3 And now I will read the Appointment to Temporary Voting Status Statement.

4 Pursuant to the authority granted under the Medical Devices Advisory Committee  
5 Charter of the Center for Devices and Radiological Health, dated October 27th, 1990, and as  
6 amended August 18th, 2006, I appoint the following individuals as voting members of the  
7 Circulatory Devices Panel for the duration of this meeting on October 27th, 2020:

8 Wayne Batchelor, Pramod Bonde, Jeffrey Borer, M.D., Joaquin Cigarroa,  
9 Bernard Gersh, John Hirshfeld, Verghese Mathew, Erik Magnus Ohman, Richard Page,  
10 John Somberg, George Vetovec, Janet Wittes, David Yuh, M.D.

11 For the record, these individuals are special Government employees or regular  
12 Government employees who have undergone the customary conflict of interest review and  
13 have reviewed the material to be considered at this meeting.

14 This was signed by Jeffrey Shuren, Director of Center for Devices and Radiological  
15 Health, on October 15th, 2020.

16 A copy of this statement will be available for review and will be included as part of  
17 the official transcript. Thank you.

18 FDA encourages all other participants to advise the Panel of any financial  
19 relationships they may have with any firms at issue.

20 Before I turn the meeting back over to Dr. Lange, I would like to make a few general  
21 announcements.

22 In order to help the transcriber identify who is speaking, please be sure to identify  
23 yourself each and every time that you speak.

24 Transcripts of today's meeting will be available from Free State Court Reporting, Inc.

25 The press contact for today's meeting is Lindsey O'Keefe.

Free State Reporting, Inc.  
1378 Cape Saint Claire Road  
Annapolis, MD 21409  
(410) 974-0947

1 And for the record, FDA has received one written comment. Thank you very much.

2 And Dr. Lange, please continue.

3 DR. LANGE: Thank you very much, Aden.

4 We will now proceed to the Sponsor's presentation. I would like to invite the  
5 Sponsor begin.

6 I will remind the public observers at this meeting that while this meeting is open for  
7 public observation, public attendees may not participate except at the specific request of  
8 the Panel Chair, and that would be me.

9 The Sponsor will have 90 minutes to present. And please begin your presentation.

10 MS. BEBEAU: Good morning, Mr. Chairman, members of the Committee, and  
11 members of the Food and Drug Administration. My name is Vicki Bebeau, I'm the Vice  
12 President of Clinical and Regulatory Affairs at Neovasc. We're pleased to be here today to  
13 present data on our Reducer device for patients who, despite receiving optimal medical  
14 treatment, continue to suffer from refractory angina pectoris.

15 Results from clinical trials and real-world experience provide strong evidence that  
16 the Reducer is safe and effective in providing symptom relief and improving quality of life in  
17 these patients who have failed other therapies and continue to suffer from this disabling  
18 condition.

19 Refractory angina pectoris severely impacts patients' quality of life, leaving them  
20 with pain and disability. These are common symptoms experienced by patients with  
21 coronary artery disease. More than 1.8 million people in the United States suffer from  
22 refractory angina pectoris. About 26,000 to 52,000 have Class III or IV angina, resulting in  
23 severe limitation in their ability to perform activities of daily living.

24 Despite advances in new drug and device therapies for treating coronary artery  
25 disease, many patients with refractory angina remain severely disabled by this condition.

1 The population of patients suffering with angina is rapidly growing, as available  
2 improvements in cardiovascular care have continued to extend life expectancy without the  
3 ability to treat refractory angina symptoms.

4 The Reducer is a therapy aimed at attenuating the disabling symptoms of angina and  
5 subsequently improving quality of life. It was designed to replicate key aspects of a surgical  
6 procedure developed by Dr. Claude Beck in 1950 to narrow the coronary sinus. This  
7 procedure was done before coronary angioplasty and bypass surgery were even available to  
8 treat patients. Dr. Beck's work demonstrated that the success of this procedure was likely  
9 driven by elevated coronary sinus pressure, which triggered protective mechanisms that  
10 improved perfusion of ischemic territories.

11 Beck's studies have been duplicated by other surgeons with equally positive results.  
12 For example, in a comparative study of 45 patients at Uppsala, 75% of patients treated with  
13 Beck's procedure reported complete or significant relief of their symptoms. Additionally,  
14 the portion of patients fit for work almost tripled in the treatment arm from 20 to 58%.  
15 However, this open chest procedure is too high risk for patients who can't undergo for the  
16 revascularization.

17 We designed the Reducer to replicate Beck's narrowing of the coronary sinus,  
18 creating a slight increase in back pressure in the coronary sinus and the small blood vessels  
19 of the heart using a minimally invasive catheterization procedure that can be performed on  
20 this patient population. This back pressure improves oxygenated blood flow into ischemic  
21 areas of the heart.

22 The Reducer is a balloon-expandable stainless steel device. Inflation of the balloon  
23 expands the Reducer to a final profile determined by the inflation pressure of the balloon.  
24 The Reducer is an over-the-wire catheter with a unique hourglass shaped balloon. The  
25 proximal and distal portions of the balloon have differing diameters to conform to the taper

1 typically encountered in the coronary sinus. It is available in one single model designed to  
2 fit the range of anatomies encountered in most patients. It is compatible with coronary  
3 sinus dimensions from 9.5 to 13 mm in diameter at the proximal implant site.

4 The Reducer is implanted in the coronary sinus, the large vein that drains blood from  
5 the heart muscle. The Reducer is implanted 2 to 4 cm distal to the ostium.

6 The procedure is done under local anesthesia so the patient can be discharged the  
7 same day. This procedure is done using a catheter inserted through a small needle  
8 puncture in the internal jugular vein, then advanced into the right atrium, then into the  
9 coronary sinus.

10 The Reducer is intended to create a narrowing in the coronary sinus to increase back  
11 pressure and redistribute blood flow from non-ischemic myocardium to more ischemic  
12 territories of the myocardium 4 to 6 weeks after implant when all but the center portion of  
13 the metal mesh is covered with tissue in-growth. This has been shown to relieve symptoms  
14 in refractory angina patients and significantly improve their quality of life.

15 Based on the data we will present today, we are proposing the following indication  
16 for use. The Reducer system is intended for patients suffering from refractory angina  
17 pectoris despite guideline-directed medical therapy, who are unsuitable for  
18 revascularization by coronary artery bypass grafting or by percutaneous coronary  
19 intervention. This proposed indication met the criteria for Breakthrough Device  
20 Designation. As noted by the Agency at that time, the final indication for our device will  
21 depend on the review of the marketing application, and we look forward to receiving any  
22 feedback from the Panel today regarding any changes to this proposed indication.

23 We started approval for the COSIRA study in both the United States and Europe in  
24 2010, which, considering the time, many companies were seeking study approvals outside  
25 the United States. We were able to gain approval very quickly in Europe, where we chose

Free State Reporting, Inc.  
1378 Cape Saint Claire Road  
Annapolis, MD 21409  
(410) 974-0947

1 to initiate the study. The Reducer achieved CE mark in 2011 based on first-in-man and  
2 initial COSIRA data.

3 Although approved, we chose not to launch the Reducer until the COSIRA was  
4 completed and we had further safety and effectiveness data. The COSIRA final report was  
5 completed in November 2014 and we began a limited launch of the Reducer in 2015 in eight  
6 countries. Today, the Reducer is distributed in 18 countries. We continue to control  
7 distribution due to resources and the challenges of obtaining reimbursement coverage for a  
8 new device in Europe.

9 As part of our postmarket surveillance, we initiated the REDUCER-I study, which is  
10 monitored and has an independent CEC adjudication committee and planned 5-year follow-  
11 up. In 2017 we received IDE approval for the COSIRA-II study, which is now the basis of our  
12 proposed postmarket study.

13 In October 2018, the Reducer also received right to designation by FDA. This  
14 prompted us to consider submitting the device for PMA approval last December based on  
15 the volume of the data we had, especially given the focus of the program to provide more  
16 timely access of innovative medical devices for patients with unmet needs.

17 The goal of the Breakthrough Devices Program is to provide patients and healthcare  
18 providers timely access to medical devices by speeding up their development, assessment,  
19 and review. The program provides sponsors with a high level of interaction with FDA,  
20 prioritized review of submissions, and opportunities for flexible clinical study design.  
21 Another important feature is the balance of pre- and postmarket data collection allowing  
22 for postmarket trials to reduce any uncertainty regarding the device's clinical benefit.

23 Neovasc has a considerable premarket dataset that consists of a sham-controlled,  
24 double-blind, randomized clinical trial that reached statistical significance for its pre-  
25 specified primary endpoint. We will also present data from our postmarket observational

1 study that confirms the effects seen in COSIRA with current follow-up up to 5 years. This is  
2 supported with follow-up of first-in-human patients out to 12 years.

3 Additionally, there is considerable published real-world experience from  
4 investigator-initiated studies that further confirm the results seen in COSIRA. We will be  
5 presenting data selected from these 50 publications based on reported sample sizes of 40  
6 or more patients.

7 FDA designated Reducer as a breakthrough device because it has the potential to  
8 meet an unmet need in this "no option" patient population.

9 Consistent with this regulatory strategy to balance data collection for novel  
10 breakthrough devices, Neovasc has proposed a robust post-approval study to further  
11 confirm the safety and effectiveness of the Reducer. REDUCER-II is a multicenter,  
12 randomized, 1:1 ratio, double-blind, sham-controlled clinical trial with up to 500 patients at  
13 up to 25 investigational centers in North America; will be randomized and followed at  
14 baseline procedure, discharge, 30 days, 90 days, 6 months, 12 months, and annually  
15 through 5 years. This will include sites in jurisdictions where the product is not yet  
16 approved, as well as select U.S. sites. If requested, patients may be unblinded to their  
17 treatment assignment after completion of their 6-month follow-up visit. Patients  
18 randomized to the sham-control arm will be allowed, but not required, to cross over to the  
19 treatment arm at the 6-month follow-up time point after completion of the study visit,  
20 provided they continue to meet all of the inclusion/exclusion criteria and are reevaluated  
21 and approved by the central screening eligibility committee.

22 FDA has published multiple guidance documents addressing innovative breakthrough  
23 medical devices and the factors that should be considered to ensure timely access to safe  
24 and effective medical devices. This includes the factors to consider when assessing the  
25 benefit-risk profile for PMA approval for novel innovative devices.

Free State Reporting, Inc.  
1378 Cape Saint Claire Road  
Annapolis, MD 21409  
(410) 974-0947

1 Neovasc has a considerable premarket dataset that consists of a sham-controlled,  
2 double-blind, randomized clinical trial that reached statistical significance. We also have an  
3 observational study that confirms the effects seen in COSIRA with current follow-up to 5  
4 years. That is further supported with follow-up of first-in-human patients out to 12 years.  
5 Additionally, there is considerable published real-world experience from investigator-  
6 initiated studies that further confirm the results seen in COSIRA.

7 In summary, we believe the Reducer has a favorable benefit-risk profile for the  
8 following reasons: We met our primary endpoint for effectiveness, established a safe use  
9 profile, and are intending to treat a "no option" patient population with a very poor quality  
10 of life with this designated breakthrough device.

11 In FDA's published guidance on the considerations of uncertainty in determining  
12 benefit-risk, they provide examples that are intended to assist in the interpretation in the  
13 context of the statutory standard. Example 1(a) discusses a breakthrough device intended  
14 to treat a treatment resistant condition. Like that example, Reducer is a breakthrough  
15 device that is also intended to treat a treatment resistant condition. In that example, they  
16 discuss three scenarios of uncertainty. The greater the level of uncertainty, the greater the  
17 reliance on postmarket data collection.

18 As mentioned previously, Neovasc has proposed a robust postmarket study, the  
19 REDUCER-II study. This is also consistent with the regulatory strategy published in 2015 by  
20 FDA in a separate guidance entitled "Balancing premarket and postmarket data collection  
21 for devices subjected to premarket approval."

22 We will spend the next 90 minutes presenting the clinical data to support our  
23 position that the Reducer has a reasonable assurance of safety and effectiveness, and the  
24 benefits clearly outweigh the risks in this "no option" patient population.

25 Here now is the agenda for the remainder of our presentation. We are very

Free State Reporting, Inc.  
1378 Cape Saint Claire Road  
Annapolis, MD 21409  
(410) 974-0947

1 fortunate today to have this esteemed group of physicians presenting. All external  
2 presenters have been compensated for their time. I will now ask Dr. Henry to present on  
3 the unmet need in this refractory angina population.

4 DR. HENRY: Good morning, everyone. I am delighted to be here today. My name is  
5 Tim Henry and I'm an interventional cardiologist and currently the director of the Lindner  
6 Research Center at The Christ Hospital in Cincinnati. Prior to this, I was director of research  
7 at the Minneapolis Heart Institute in Minneapolis, and then chief of cardiology at Cedars-  
8 Sinai Medical Center in Los Angeles for the 5 years before coming here.

9 My presentation today will focus on refractory angina, specifically to describe key  
10 aspects of the patient population and importantly, the unmet clinical need. I have  
11 personally spent more than 25 years taking care of refractory angina patients and, in fact,  
12 I've had a dedicated refractory angina clinic for a long time, almost 2 decades.

13 As we all know, angina is common in the United States. The American Heart  
14 Association, in fact, estimates that there are 8 to 12 million patients with angina. But I want  
15 to specifically emphasize that our focus today is on refractory angina, a very specific subset  
16 of angina that has a major effect on patients' quality of life and leads to increased  
17 healthcare utilization. In fact, frequently, these patients have been referred to as "no  
18 option" patients. We talked about patients who have persistent angina and ischemia  
19 despite optimal medical management or medical intolerance in patients who are unsuitable  
20 for further percutaneous or surgical revascularization.

21 The most commonly used definition in the literature comes from the European  
22 Society of Cardiology, which is quite dated, actually, but it defines refractory angina as at  
23 least 3 months of angina due to demonstrated coronary insufficiency that persists despite  
24 optimal medical therapy in patients who are no longer amenable to percutaneous or  
25 surgical revascularization.



1           We recently published a manuscript describing four specific phenotypes in refractory  
2 angina and I think it's worth reviewing these to just get an idea of the patients we're talking  
3 about.

4           The first patient population in the upper left-hand corner is not the focus of today's  
5 presentation. These patients have microvascular angina, so they have minimal coronary  
6 disease but still severe angina due to microvascular dysfunction. And certainly, those  
7 patients can have refractory angina, but that's not who we're talking about and the ESC  
8 definition defined it that you have coronary insufficiency, that these are patients with  
9 coronary disease.

10           So the three phenotypes that are the focus of today are, first, in the left lower  
11 quadrant, these are patients who have one or two specific territories at risk that are  
12 collateral dependent. So this is the classic CTO patient and many, but not all, may be  
13 amenable to revascularization with advanced CTO technology, but these procedures have  
14 some risk and 15 to 20% are unsuccessful.

15           The third patient population in the upper right-panel describes patients who have  
16 diffuse thread-like coronary atherosclerosis. This is typical Type 1 or Type 2 diabetics that  
17 have diffuse distal and side branch disease that frequently is suboptimal for  
18 revascularization or, even if part of the coronary anatomy is amenable, there are territories  
19 that because of the diffuse disease are not amenable to revascularization.

20           And then finally, a fourth phenotype describes patients who have end-stage or  
21 advanced coronary disease. These are patients that we all know that have had one or more  
22 bypass surgeries, multiple PCIs, often 10 or 12 PCIs, and these are truly end-stage coronary  
23 disease patients that no longer have revascularization options and typically, they have areas  
24 of ischemia and multiple distributions.

25           I think reviewing these phenotypes is helpful because it illustrates the heterogeneity

1 of the refractory angina population.

2           There are many ways to grade angina, but the most common is the Canadian  
3 Cardiovascular Society scale. Most of you are familiar with this, but it's worth reviewing  
4 just to focus our attention today.

5           Class I is angina only with strenuous or prolonged exertion at work or recreation.

6           Class II is a slight limitation of ordinary activity. These are patients who can walk  
7 more than two blocks or more than a flight of stairs and are asymptomatic, but with  
8 extensive exertion or maybe with an additional -- after a meal or in cold weather might  
9 have angina. But those are not who we're talking about.

10           Our focus today really starts with Class III patients and it's a very important  
11 distinction. These are patients that, despite medical therapy and despite revascularization,  
12 have marked limitation of ordinary activity, so walking one or two blocks on the level or  
13 walking one flight of stairs precipitates angina. This is a severe limitation in activity.

14           And finally, Class IV patients are unable to do almost any activity without discomfort  
15 and frequently have angina at rest.

16           I want to emphasize, our focus is on Class III and Class IV angina for the discussion  
17 today. And these are patients who frequently have 20, 30 or even more episodes of angina  
18 per week.

19           The prevalence of refractory angina is another important issue and it's very difficult  
20 to determine and that's because there are no specific codes and there's limited registry  
21 data. So this is actually a figure from a manuscript that we have in press in *Circulation:*  
22 *Cardiovascular Intervention* that should be out, I think, in December, using data from both  
23 cardiac catheterization as well as angina populations.

24           If you start with the overall population of angina in the United States, it's 8 to 10  
25 million. If you start with the cardiac cath lab, roughly 6 to 14% of patients undergoing

1 cardiac catheterization have coronary disease that's not amenable to revascularization.

2 Now, if you only include Class III and Class IV angina, it's estimated that this patient  
3 population is about 600,000 to 1.3 million. But I will point out that of the Class III/Class IV,  
4 most of those patients have an acute coronary syndrome or are not amenable for  
5 revascularization.

6 So finally, when you get to this point where you have stable coronary disease, not  
7 candidates for revascularization, Class III and Class IV angina and then with the additional  
8 greater than seven episodes a week, which is frequently used in clinical trials for enrollment  
9 purposes, this is a very specific population that we estimate to be between 26,000 and  
10 52,000 patients a year, which is basically an orphan population.

11 So it's important to reiterate, there's no specific codes for refractory angina and only  
12 a limited number of registries. Therefore, it is challenging to know the total number. But  
13 clearly this is one of the major challenges in caring for these patients, that it's really difficult  
14 to know the total number.

15 In regards to the natural history of the disease, I'll use data from our dedicated  
16 refractory angina clinic that we published in the *European Heart Journal*. It's now a  
17 prospective registry that includes more than 2,000 patients.

18 In contrast to the previous historical data, we published this data in the *European*  
19 *Heart Journal* that indicates the actual mortality for the population is about 3 to 4% per  
20 year. It's a little higher when you include only Class III and Class IV, but even for Class III  
21 and Class IV patients, it's only about 5 to 6% per year. So this contemporary natural history  
22 data is not much higher in patients with stable coronary disease or age-matched controls.

23 So the issue is not mortality so much as it is these patients are left with severe  
24 limitations in their quality of life and have a significant increase in healthcare utilization. So  
25 from my perspective, it's very important to practice a patient-centered health assessment

1 approach. Measuring success of the therapy needs to focus on relieving symptoms,  
2 improving quality of life and functional capacity.

3 This is more data from our dedicated refractory angina clinic to provide some  
4 insights into the baseline characteristics of refractory angina patients. Whereas many of  
5 these patients have been smokers in the past, the majority have stopped smoking and  
6 current smoking is only roughly 10 to 15%. These are typically patients who are extremely  
7 well treated with both secondary prevention medications and antianginal medications. And  
8 in fact, the group we're talking about, by definition, has to have failed standard medical  
9 therapy.

10 About 30 to 40% of the population is diabetic. This has been consistent in clinical  
11 trials and registries. About 30% will have a history of heart failure. About 75% have a  
12 history of a previous myocardial infarction.

13 And in terms of revascularization, 70 to 80% of patients will have previous CABG or  
14 previous PCI. And if you consider any revascularization, it's over 90%. Of course, this  
15 depends a little bit on whether you include microvascular angina phenotype in your  
16 registry, but even if you only include patients with coronary disease, there are a few  
17 patients, for example, severe diabetics that might have never had revascularizations.  
18 Clearly, for sure, more than 90% will have had previous revascularization and the majority  
19 will have multiple procedures.

20 If you look at mortality, and I think I alluded to it, but to really emphasize again, the  
21 mortality at 5 years was only 17% and two-thirds of those deaths were cardiovascular in  
22 nature, so that means one-third were non-cardiovascular. And again, as we noted before,  
23 this is only slightly higher than stable coronary disease. So clearly, while we've made  
24 progress in reducing mortality, the challenge is to manage the debilitating symptoms.

25 Standard therapy for angina in the United States includes antianginals, beta blockers

1 and calcium channel blockers, short-acting nitrates, and ranolazine. These are frequently  
2 used in a stepwise manner, but often they're not well tolerated either due to side effects or  
3 relative hypotension, and it's rare that someone actually tolerates all four. And the average  
4 in clinical trials is about two and a half or so.

5 Another key factor in the medical treatment is secondary prevention with  
6 antiplatelet therapy, lipid-lowering agents, and blood pressure control, and it's a very  
7 important aspect of treatment of these patients.

8 And then finally, we have revascularization, which includes both percutaneous and  
9 surgical revascularization, both of which have improved over time, including new aggressive  
10 approaches to chronic total occlusion and certainly that's helped our therapy.

11 So these are the standard therapies for angina, which actually work quite well, and  
12 we've clearly made progress overall in the treatment of angina.

13 So it's important to emphasize again, the patients we're talking about today have  
14 Class III and Class IV angina despite secondary prevention, despite optimal antianginal  
15 therapy, and despite revascularization. They are suboptimal or unsuitable for  
16 revascularization for a variety of reasons. These are truly "no option" patients.

17 So what is available for these patients in the United States? Really, only two  
18 treatments. EECF, most of you are familiar with EECF, these are leg cuffs that we put on  
19 that sequentially inflate during diastole and then deflate during systole. And this is  
20 physiology similar to an intra-aortic balloon pump, it was approved more than 20 years ago,  
21 and the only randomized trial was the MUST-EECF trial, which was -- did include a sham  
22 procedure, it enrolled 137 patients and, in fact, only 116 were in the per-protocol analysis.

23 The problem, of course, with EECF is availability. It requires thirty-five 1-hour  
24 sessions 5 days a week for 7 weeks, so it's a significant time commitment and frequently, in  
25 my experience, it's logistically difficult because people have difficulty driving every day to

1 get to where an EECp therapy is. So there's limited availability for many patients in the  
2 United States. In addition, most of the patients that we're talking about have already had  
3 one or more EECp courses of therapy.

4 The second approved therapy for refractory angina in the United States is  
5 transmyocardial laser revascularization, or TMR, so this is laser treatment to create small  
6 channels in the myocardium to promote blood supply from the LV cavity. So TMR is  
7 infrequently used currently because of concerns over morbidity and mortality in the patient  
8 population. And again, there's limited availability and very little TMR is actually being done  
9 in the United States, at least standalone TMR is being done at this time.

10 So both of these therapies have limitations and again, I'll emphasize it's been more  
11 than 20 years since there's been a therapy approved for the treatment of refractory angina.

12 So today Dr. Stone will present data on a one and two CCS class improvement. So I'd  
13 like to at least spend a moment on the clinical relevance of this level of improvement.

14 Certainly, a two-class improvement in angina, going from Class IV to Class II or from  
15 Class III to Class I, is a significant improvement in quality of life. Really, an improvement of  
16 two classes enables the patient to be able to go back to work, to care for themselves, to do  
17 their activities of daily living, spend time with their children, it's really a major improvement  
18 in the quality of life. And I think our goal, when you take care of refractory angina patients,  
19 is to actually ideally get back to Class I or Class II where they can have a near normal or  
20 active life and it's unlikely that we'll ever really make patients completely angina free, but  
21 our goal is to actually improve by at least one class and ideally, two classes. So I think this is  
22 a clinically meaningful definition in improvement.

23 In summary, patients who have refractory angina or "no option" patients are a very  
24 specific subset of angina with severe limitations in quality of life and high healthcare  
25 utilization. These patients, unfortunately, have very limited treatment options in the

1 United States at this time; they are severely disabled and are a very challenging patient  
2 population to care for.

3 Thank you very much for your time and next, Dr. Shmuel Banai will discuss the  
4 mechanism of action for the Reducer.

5 DR. BANAI: Good morning, my name is Shmuel Banai and I am the medical director  
6 for Neovasc. I'm also the director of the Division of Cardiology at the Tel Aviv Medical  
7 Center, Professor of Cardiology at the Tel Aviv University School of Medicine. Additionally, I  
8 am the head of the Center for Cardiac Research and Biomedical Engineering at the Tel Aviv  
9 University. Today I will share with you the Reducer mechanism of action.

10 This is a schematic illustration of the full thickness of the left ventricular wall. At the  
11 top is the epicardial coronary artery and the epicardium, and at the bottom, the endocardial  
12 layers of the myocardium and the adjacent left ventricular cavity. The arrows on the left  
13 indicate blood flow into the myocardium.

14 In the normal healthy heart, blood flow in the subendocardial layers of the  
15 myocardium is higher than in the subepicardial layers of the myocardium. During exercise  
16 and increased demand, the physiologic compensatory mechanism causes selective  
17 sympathetically mediated vessel constriction and increased resistance to flow in the  
18 subepicardial blood vessels, favoring subendocardial perfusion and allowing an appropriate  
19 augmented contractility. This increased demand and augmented blood flow is represented  
20 by the larger arrow on the left of the illustration. As a result, in the healthy heart the ratio  
21 of endocardial to epicardial blood flow is approximately 1.2:1 due to preferential dilation of  
22 the subendocardial arterioles, causing a larger diastolic flow through the subendocardium.  
23 In other words, in the normal heart, the subendocardium receives approximately 20% of  
24 blood flow.

25 In contrast, in the ischemic heart, as represented here with yellow obstructive

1 plaque in the epicardial coronary artery, the situation changes. The subendocardium is  
2 more vulnerable to ischemia than the subepicardium. The compensatory mechanism I told  
3 you about, which normally preserves preferential blood flow through the subendocardial  
4 myocardium, becomes dysfunctional and the transmural perfusion is redistributed  
5 towards the subepicardial layers of the ventricle, represented here by the larger arrow  
6 directing into the subepicardium. Resistance to flow in the subendocardial layers of the  
7 myocardium is increased and the endocardial to epicardial blood flow ratio becomes  
8 abnormal and falls from 1.2 at baseline to 0.5.

9 As a result, the perfusion of the subendocardium during exercise becomes  
10 compromised, causing ischemia, impaired contractility, and elevated left ventricular end-  
11 diastolic pressure, which compresses further on the subendocardial blood vessels, thereby  
12 exacerbating ischemia. This creates a vicious cycle of ischemia. Exercise-induced ischemia  
13 leads to anginal symptoms. Reduced left ventricular contractility with increased left  
14 ventricular end-diastolic pressure results in shortness of breath.

15 The mechanism of action presented here is based on animal research using ischemic  
16 heart model with and without elevating coronary sinus pressure. When the Reducer is  
17 implanted in the coronary sinus of an ischemic heart and after the device is fully covered  
18 with tissue in-growth, it elevates backward venous pressure which results in redistribution  
19 of blood into the ischemic subendocardium, breaking this vicious cycle of ischemia.  
20 Elevating backward pressure results in a slight dilation of the diameter of arterioles and  
21 leads to a reduction of vascular resistance in the subendocardium with a subsequent  
22 normalization of the endocardial to epicardial blood flow ratio back to 1.2 that drives blood  
23 flow back into the ischemic subendocardium.

24 Consequently, blood flow in the ischemic subendocardial layers of the myocardium is  
25 enhanced, contractility improves, and left ventricular end-diastolic pressure decreases. This



1 breaks the vicious cycle of ischemia, improves left ventricular contractility with reduction of  
2 left ventricular end-diastolic pressure, and improves exercise-induced shortness of breath.  
3 Improvement of subendocardial blood flow leads to improvement of angina symptoms.

4 In summary, by increasing coronary sinus pressure, the Reducer restores the  
5 endocardial to epicardial blood flow ratio to normal, improves blood flow to the ischemic  
6 subendocardium, increases left ventricular contractility, and reduces left ventricular end-  
7 diastolic pressure.

8 Thank you. I will now turn the presentation over to Dr. Victoria Hampshire.

9 DR. HAMPSHIRE: Thank you, Dr. Banai. My name is Victoria Hampshire, I'm a  
10 veterinarian by training and I've spent the last 30 years in animal research support in the  
11 evaluation of medical devices. For the first 10 years of my career I was on the front lines of  
12 research support at the NIH for cardiovascular and other animal models coming through the  
13 research pipeline. I then transferred to the FDA and worked there for 15 years. Five years  
14 were spent at the Center for Veterinary Medicine as the chief safety officer overseeing  
15 post-marketed veterinary products, and the last 10 at the Center for Devices and  
16 Radiological Health in the Division of Cardiovascular Devices. In 2015 I left the federal  
17 government to start a preclinical consulting company, where I've been busy helping dozens  
18 of companies at all phases of animal study design, conduct, and reporting. I've been  
19 working with Neovasc for the past 5 years to help communicate the totality of their animal  
20 data to the FDA. I'm very happy to be here today to speak to the preclinical animal data for  
21 the Reducer.

22 From 2002 to 2009, Neovasc conducted three animal studies using pig models. Two  
23 were pilot animal studies that evaluated the Reducer in 39 pigs. The pivotal study  
24 evaluated the Reducer in 13 pigs. Four pigs were studied acutely and two pigs were studied  
25 at 24 hours. The remaining seven pigs were studied chronically between 57 days and as

Free State Reporting, Inc.  
1378 Cape Saint Claire Road  
Annapolis, MD 21409  
(410) 974-0947

1 late as 140 days. This subset of pigs had survived the index procedure of the Reducer  
2 implant and went on to their respective endpoint where terminal end-life assessments and  
3 pathology from their tissues was obtained for evaluation.

4 The evaluators of these studies included an internationally recognized cardiologist;  
5 two highly experienced interventionalists; and three board certified pathologists, one  
6 human and two veterinary. Each evaluator independently reported similar findings.

7 Although not strictly GLP, these studies were undertaken using protocols and  
8 standard operating procedures that included modern methods of humane animal  
9 experimentation. They were conducted at experienced and respected international animal  
10 facilities by trained support personnel, and they were consistent with the level of preclinical  
11 research typically conducted for first-in-man or breakthrough technology.

12 Importantly, formalin-fixed heart and histologic assessments critical to the questions  
13 raised by FDA were provided by independent examiners at GLP pathology laboratories in  
14 the United States.

15 The in-life and afterlife methods include sophisticated preoperative preparation,  
16 interventional procedures and monitoring, imaging, echocardiographic and angiographic  
17 assessment, performance and handling assessment, and gross and histological processing.

18 To begin summarizing, critical to some of the questions today is an understanding of  
19 some key anatomic differences of quadrupeds that are available for an implant the size of  
20 the Reducer, including pigs.

21 Here on the left is a human image showing the human right coronary artery  
22 circulation as compared to the latex-perfused circulation in the pig heart on the right.  
23 Viewers can see that the pig and also, by the way, other quadrupeds have impressive azygos  
24 vein contribution to the right heart circulation which is not present in humans. This vein  
25 drains the four legs.

1 Please remember today that any pressure reports for swine outcomes are useful but  
2 not necessarily translational to humans because of these anatomic differences. The  
3 6-month outcome of the pressure gradient across the coronary sinus in the pig is not  
4 concerning and in point of fact, was a bit of a favorable finding considering this anatomic  
5 difference.

6 Now, turning to the data that FDA did receive, the FDA was provided with final study  
7 reports and a final pathology report. For the pivotal animal study, the FDA was also  
8 provided source records containing data and information from implants, Reducer delivery  
9 system performance, and a final pathology report. This report contained macroscopic and  
10 microscopic data from the chronic animal subset of seven pigs. The Reducers were  
11 sectioned at three levels, proximal, mid, and distal, and were evaluated by light microscopy  
12 and histomorphometry.

13 This slide shows a variety of imaging data provided in the animal study report to the  
14 FDA. Before I describe this slide, I would like to clarify that the animal pathologist calls the  
15 inflow proximal because it is deeper inside the heart. The clinician, however, calls the  
16 outflow of the heart that is towards the right atrium proximal because it is closer to the  
17 ostium of the coronary sinus. The wider end of the Reducer, which is closer to the ostium,  
18 is proximal.

19 Neovasc did submit a variety of other imaging information to the FDA from the final  
20 pivotal animal study and in the final study report, beginning with angiographic  
21 documentation in-life and implant showing the Reducer deployed within the midsection,  
22 prominently demonstrative of the luminal restriction, and then again by macroscopic means  
23 at the terminal endpoint shown here in the middle, also taken from the final study report.  
24 In turn, there were also faxitron images to show preservation of stent integrity, and below  
25 the explanted preserved imprint of the Reducer from the pathology report.

1           The summary of the findings concluded that there was acceptable regional tissue-  
2     contacting gross pathology and histology with absence of migration and perforation. There  
3     was a hundred percent coronary sinus lumen endothelialization and tissue proliferation at  
4     proximal and distal ends of the Reducer. There was incomplete coverage of the midsection  
5     struts and preservation of lumen from proximal to distal. This includes luminal preservation  
6     through the midsection with very low levels of inflammation at all time points and in all  
7     sections, and evidence of 10% stenosis at the mid-plane, consistent with the intended effect  
8     of the Reducer geometry. Notably, these findings do support the mechanism of action of  
9     the Reducer.

10           Implantation of the Reducer is performed with an intentionally approximate 10%  
11     over-sizing of both wide ends of the device. Over-sizing is important and helps to achieve  
12     two goals. First, to anchor the elastic vessel wall to help prevent migration and trigger a  
13     process of tissue proliferation which is observable at 57 to 140 days after implantation, and  
14     covers the gaps between the middle struts. This aligns the data from preclinical feasibility  
15     studies.

16           Second, since the narrow central part of the device is not in direct contact with the  
17     vessel wall, there is no trigger for significant tissue overgrowth at this point and therefore  
18     the vessel lumen at the center of the device remains patent. In fact, the proliferative  
19     process is less intense and sometimes absent in the midsection of the Reducer. Any  
20     neointimal coverage of the midsection is the result of encroachment from either end,  
21     proximal or distal. Thus, embedding both wide ends of the Reducer in tissue is sufficient to  
22     create a tube-shaped narrowing even when the center struts are not fully endothelialized.  
23     The tube-shaped narrowing also establishes a mechanism for a pressure gradient across the  
24     narrow center of the device, which is consistent with the mechanism of action reported in  
25     the feasibility studies and in the literature.

1           And the totality of the histologic data is consistent with the postulated mechanism  
2 of action. Overall, the preclinical studies demonstrate acceptable in vivo performance for  
3 introduction, location, target anatomy, deployment, angiographic visualization, and delivery  
4 system of the device. Additionally, the histology of the device is acceptable for all  
5 inflammation and injury criteria.

6           Thank you for your attention to this important preclinical data. I am now going to  
7 turn over a few minutes of time to Dr. Rousselle, the study pathologist, to go over with the  
8 Panel the details about endothelialization so that everyone who is not a pathologist will  
9 understand the information before them on that particular topic.

10           Dr. Rousselle.

11           DR. ROUSSELLE: Thank you, Dr. Hampshire. My name is Serge Rousselle and I'm a  
12 board certified veterinary pathologist and have done medical device safety assessment for  
13 the last 20 years. I was involved in evaluations in some of the studies that are discussed  
14 here today.

15           On this slide we're showing the Reducer midsection and the final pathologist's report  
16 makes a distinction between coronary sinus luminal coverage, which was 100%  
17 endothelialized, and Reducer midsection struts, which were no more than 60%  
18 endothelialized. In the image on the slide there are two profiles visible, one delineating the  
19 sinus profile in pale blue and the other delineating the Reducer profile. Both outlines have  
20 importance in this model and the scoring of these two features for endothelial coverage  
21 was intended to provide a higher level of data granularity. The data demonstrated that the  
22 sinus lumen tissue was fully endothelialized and the device was also endothelialized  
23 wherever it is in contact with the host tissue, as expected.

24           This slide represents two cases of the most advanced tissue integration of the device  
25 for Animal 1111 and the least advanced integration at the mid-level in Animal 1035. The

1 top example clearly demonstrates that as the device becomes fully integrated, there is  
2 growing neointimal formation, visually estimated at about 30% in this case, which reduces  
3 the lumen commensurately.

4 In the lower case, the sinus lumen remains open. Nonetheless, reduction of the  
5 lumen is evident for the combination of slight neointimal formation and frame of structure,  
6 visually estimated here at 5 to 10%.

7 The negative restenosis values that were presented in the pathology report resulted  
8 from an artifact of calculation, i.e., a reduction of device lumen relative to sinus lumen. The  
9 calculation of sinus lumen restenosis, so sinus lumen compared to original sinus profile,  
10 would deliver positive values at all mid-levels regardless of -- studies.

11 Thank you for your time. I'll now turn the presentation over to Dr. Stone.

12 DR. STONE: Well, hello, everybody. I'm Gregg Stone from the Mount Sinai Heart  
13 Health System, and it's a pleasure to be able to speak to you today. I'm going to go over the  
14 detailed study design for the Reducer studies that we completed to date, including the  
15 sham-controlled, randomized COSIRA data, and go over the principal safety and  
16 effectiveness outcomes. Specially, I'm going to focus primarily on three different datasets.  
17 First, a small first-in-human study, which was a prospective, open-label, multicenter  
18 feasibility study of 15 patients at three sites, and I'll use that mostly to look at the very long-  
19 term data. I'm going to spend most of the time emphasizing the results from the COSIRA  
20 trial. This was a prospective, randomized, double-blind, sham-controlled, multicenter trial  
21 in 104 patients at 11 sites. And as, of course, this Committee well knows, it's relatively  
22 unusual to have sham-controlled trials in interventional cardiology, so I'll be emphasizing all  
23 aspects of that important nature of the study. I'll also talk to you about the REDUCER-I  
24 study, this was a single-arm postmarket observational study, to date, in 241 patients at 20  
25 sites with enrollment ongoing.

1           So let's really start with the COSIRA sham-controlled randomized trial. This was a  
2 study that began with 104 patients with refractory angina, and I'll describe the inclusion and  
3 exclusion criteria. At baseline, once patients passed screening, they had detailed  
4 assessments, of course, of angina status, their demographics, EKG, exercise tolerance was  
5 assessed, they had a Seattle Angina Questionnaire performed in addition to their  
6 assessment of the severity of angina by the Canadian Cardiovascular Society classification.  
7 They had a SPECT thallium performed, they had a dobutamine stress echo, and there were  
8 other imaging subset studies, as well.

9           Patients then underwent their final qualification during cardiac catheterization, at  
10 which time coronary sinus angiography was performed, and if the coronary sinus met the  
11 criteria for inclusion in the trial, which I will also describe, then they were randomized in a  
12 blinded manner 1:1 to either Reducer implantation or control. And I'll describe, again, the  
13 procedures that were put in place for the sham control.

14           Patients, after the procedure, were then followed up at 30 days, 3 months, and at 6  
15 months. The study concluded at 6 months, at which time there was a comprehensive  
16 reevaluation done in addition to clinical outcomes, as well as repeat exercise testing. The  
17 study was assessed over time by independent blinded physicians who performed the pre-  
18 and post-procedure angina assessments both for the CCS scale and Seattle Angina  
19 Questionnaire, and I'll describe some of the blinding processes.

20           So patients were enrolled in COSIRA who were 18 years or greater age, who had  
21 symptomatic coronary disease with chronic refractory angina pectoris defined as Canadian  
22 Cardiovascular Class III or IV angina, despite attempted optimal medical therapy for at least  
23 30 days prior to screening. These patients also had to have had limited treatment options  
24 for revascularization by either bypass surgery or PCI, as determined by the investigator,  
25 with evidence of reversible ischemia attributed to the left coronary artery system by

1 dobutamine echo.

2 The key exclusion criteria were unstable angina during the 30 days prior to the  
3 screening period, acute coronary syndrome within 3 months, a successful PCI or CABG  
4 within 6 months, decompensated heart failure or hospitalization for heart failure within 3  
5 months, and arrhythmia that would require a defibrillator or a pacemaker or a defibrillator  
6 or a pacemaker electrode in their right atrium, right ventricle, or coronary sinus.

7 On the table, as I mentioned, coronary sinus venography was performed and  
8 patients were excluded if they had anomalous or an abnormal coronary sinus; for example,  
9 if the coronary sinus was excessively tortuous, had aberrant branches, if it was a left-sided  
10 SVC, etc., or more importantly, if the coronary sinus diameter at the site of the planned  
11 Reducer implantation was either less than 9.5 mm or greater than 13 mm in diameter.

12 Now, the primary effectiveness endpoint of COSIRA was the improvement in two or  
13 more CCS angina classes from baseline to the 6-month post-procedure evaluation. So let  
14 me go over the CCS angina class and explain a little bit why this was chosen as the primary  
15 endpoint.

16 So CCS angina Class IV is severe limitation, either angina at rest or basically angina  
17 with any activities at all. CCS III is marked limitation with symptoms of angina with  
18 everyday living activities, so most activities would cause angina, but not at rest. Class II  
19 would be slight limitation with angina during ordinary activities, so that might be, for  
20 example, during washing dishes or walking briskly. And Class I would be angina only during  
21 strenuous or prolonged physical activities, so almost no angina, but you could induce  
22 angina.

23 So there obviously can be some variability in the assessment of the determination of  
24 what is Class III or Class IV. So for the primary endpoint, this is a very clinically relevant  
25 classification system, the most widely used classification system still today for angina



1 assessment, and a two-class or greater reduction in angina was chosen because this is really  
2 a marked reduction in angina for which there really can't be that much question.

3 So for example, if the patient -- as you'll see, only patients with Class III and IV  
4 angina were allowed to be enrolled in the trial. So if the patient had Class IV angina where  
5 they were basically having rest angina or with any exertion at all, they had to reduce their  
6 angina to this fact where they are having almost no angina or only to the extent where they  
7 could still perform their daily activities. If the patients had Class III angina at baseline, they  
8 had to be reduced to almost no angina. So it's really a marked reduction in angina for the  
9 primary endpoint of this study.

10 There were multiple secondary effectiveness endpoints and these were not powered  
11 endpoints. The primary endpoint was powered, the secondary effectiveness endpoints  
12 were not. This was to look at improvement of one or more CCS angina class, exercise  
13 tolerance testing parameters changes, dobutamine echo wall motion score index, and  
14 Seattle Angina Questionnaire five-domain improvements. And again, the trial importantly  
15 was not powered to demonstrate improvement in any of these objective evidence of  
16 myocardial ischemia, it would've required greater than 200 patients for such a reduction.  
17 So really, you're just looking to see if there were trends present and these were really  
18 hypothesis-generating endpoints.

19 There were two major safety endpoints that were pre-specified, these were not  
20 powered endpoints, either, but periprocedural serious adverse events defined as the  
21 composite of death, myocardial infarction, cardiac tamponade, life-threatening ventricular  
22 arrhythmias, and respiratory failure through 30 days. This would also include a clinically  
23 driven re-dilation of a failed Reducer in the Reducer arm.

24 And then there were major adverse events, which someone called MACE, and that's  
25 the composite of cardiac death, major stroke or myocardial infarction through hospital

1 discharge and then at 30 days, 3 months, and at 6 months.

2           So in COSIRA, all patients were -- in addition to their maximally tolerated antianginal  
3 agents, all patients were pretreated with dual antiplatelet therapy, which was continued for  
4 6 months. The Reducer patients, in addition, received either heparin or bivalirudin in the  
5 cath lab, where that was not administered to the placebo-treated -- or the sham-control  
6 treated patients.

7           So these were the measures that were put in place to ensure blinding. The  
8 physicians that were performing the implant or the sham procedure, of course, they were  
9 unblinded, but they were instructed before the procedure to behave in the exact same  
10 manner in both groups with a script that was discussed ahead of time.

11           The Sponsor was present for all these cases and ensured that the patients were  
12 either wearing headphones playing music or that the patients were heavily sedated, that is  
13 essentially unconscious, throughout the procedure.

14           The randomizations were opened -- the randomization envelopes were opened in  
15 the control room, not the catheterization laboratory, so if there was any discussion it was  
16 only in the control room.

17           The patients were aseptically draped and the way that the Reducer is put in required  
18 draping the patient's face and so the patients actually couldn't see the monitor throughout  
19 the entire procedure, again, even if the patient was at all awake.

20           The study personnel performing the follow-up Canadian Cardiovascular Society class  
21 and Seattle Angina Questionnaire assessments were blinded to randomization. The  
22 dobutamine echo and stress test core laboratories were blinded to randomization, and the  
23 CEC and the DSMB were blinded to randomization.

24           So let's talk about now the effectiveness results, and I'm going to talk about  
25 effectiveness of COSIRA and then other studies and then we'll transition over to safety.

1           So overall, 166 patients were consented. Of that, 107 met the clinical inclusion and  
2 exclusion criteria. Three patients did not meet the angiographic criteria and so 104 patients  
3 were ultimately randomized, 52 to the Reducer and 52 to no Reducer. Six-month follow-up  
4 was complete in essentially all patients, 104 patients. One patient in the control arm died,  
5 which is why it says 103 here, but it was essentially completed on all patients.

6           So this is what these patients with refractory angina looked like. Their mean age was  
7 about 68 years. Approximately 80% of the patients were male, about 85% of the patients  
8 were Caucasian, and they were mildly obese. Approximately 80% of the patients had Class  
9 III angina, where approximately 20% had Class IV angina. As would be expected, these  
10 patients had complex coronary artery disease, more than half of the patients had a prior  
11 myocardial infarction, and almost all the patients has prior PCIs, often multiple prior PCIs,  
12 and prior bypass surgery.

13           Diabetes was present in about 45% of the patients. About half the patients were  
14 smokers. Almost all had hypercholesteremia. Most had hypertension and a very high  
15 incidence of family history of cardiovascular disease.

16           Looking at the antianginal medications, you can see that there were high usage rates  
17 of beta blockers, calcium blockers, nitrates. Again, the goal was the patients were to be  
18 treated with as many antianginal medications that they could tolerate at maximal doses.  
19 There was also scattered use of nicroandil, ivabradine, and ranolazine. You can see the  
20 number of antianginal medications that patients were taking and approximately 70% of the  
21 patients were taking two or more antianginal medications.

22           So there were high rates of technical and procedural success achieved with the  
23 COSIRA implant, with the Reducer implant in COSIRA. Fifty of the 52 were successfully  
24 implanted, that's a 96.2% rate. In the two that could not be implanted, there was  
25 inadequate guide catheter support, so the patients were treated medically and no

1 complications occurred in those patients.

2 So the procedural successes were among those that had implants, that had no acute  
3 need for clinically driven intervention to address an adverse event or a serious adverse  
4 device-related event prior to discharge, and that was 100% or 50 of 50 implanted patients.

5 So this is the primary endpoint in COSIRA, which is an improvement in two or more  
6 CCS angina classes from baseline to 6 months, and we show you here on the left the  
7 intention-to-treat population, 52 versus 52 patients, and on the right the as-treated  
8 population, which is 50 Reducer implants versus 54 no Reducer implants.

9 You can see by intention to treat that 35% of the Reducer patients had a two or  
10 more CCS angina class improvement over 6 months compared to 15% of the sham-control  
11 group. That's a difference of approximately 19%. You can see that was a statistically  
12 significant difference. The number needed to treat to improve one patient by two or more  
13 CCS angina classes was about five. You can see the as-treated results were similar, if  
14 anything, slightly stronger.

15 Now, this is an important slide which shows the time course of the improvement in  
16 two or more CCS angina classes, and this not only speaks to the mechanism of the Reducer,  
17 but it also speaks to blinding and placebo effects. As you heard from Dr. Banai, we expect  
18 the Reducer to become active, if you will, in about 6 to 8 weeks as it endothelializes. At 30  
19 days, there was no significant difference in the percentage of patients that had a two or  
20 greater CCS angina class improvement, 13.5% versus 11.8%. So I would most likely attribute  
21 this to a placebo effect that was present in both arms, with no significant difference. And  
22 again, we wouldn't expect there to be complete endothelialization at 30 days yet.

23 But at 3 months and at 6 months you can see that there are now differences that  
24 favor the Reducer arm, which becomes statistically significant by 6 months. You can also  
25 see that the control arm basically preserves the rate of improvement of two or more CCS

1 angina classes, whereas the Reducer arm continues to increase.

2           So to me, again, I think this is consistent with both the timing of the mechanism of  
3 action of the Reducer, that is when it endothelializes between 6 and 8 weeks. It's also  
4 consistent with the early modest placebo effect in both groups. It also shows that that  
5 placebo effect in the control group was likely continued over 6 months, suggesting that  
6 those patients were not unblinded, otherwise I would've expected that they would lose that  
7 placebo effect. So I think this important slide again points to both the mechanism and the  
8 timing of the effect, as well as the success in the sham-controlled randomized trial.

9           Now, the study was only 104 patients, so we don't have a lot of power to look at  
10 subgroups, but these were five pre-specified subgroups based on LVEF, CABG, diabetes, sex,  
11 and age, and if you just look basically over to the right you can see that there were no  
12 significant interactions between any of the subgroups and the primary endpoint  
13 effectiveness outcome.

14           Looking at now some of the hypothesis-generating secondary endpoints,  
15 improvement in one or more CCS angina class was also substantially better in the Reducer  
16 arm than in control. And here you can see approximately 71 to 72% of the Reducer patients  
17 improved by one or more CCS angina class compared to 42% of control patients. The  
18 number needed to treat therefore to incrementally benefit one treated patient with a  
19 reduction in angina was only about three and a half patients.

20           Here are all the different angina classes. You can see red is IV, yellow is III, blue is II,  
21 green is I, gray is no angina. All the patients had Class III or IV angina at baseline, and if you  
22 just look at kind of the blue and the green, you can see the improvement at 6-month follow-  
23 up.

24           So this slide shows, again, overall changes in angina class, and if you look to the right  
25 you can see we've got the Reducer treatment in the dark bars and the sham-control group

1 in the light bars. You can see the greater than or equal to two-class improvement, those  
2 are the data that I already showed you. The third set of bars are exactly one-class  
3 improvement and that favored the Reducer group. Looking at no change, that favored the  
4 Reducer group in a positive way. Fewer patients didn't change. And if you looked at one-  
5 class worsening, that's the set of bars over to the left, you could see that fewer patients had  
6 worsening angina that were treated with the Reducer compared to the sham control.

7 This slide shows the Seattle Angina Questionnaire and, as you know, there's five  
8 domains that were assessed here, quality of life, angina stability, angina frequency, physical  
9 limitations, and treatment satisfaction. And this slide is showing the change in each of  
10 those domains from baseline to 6 months.

11 Again, we were not powered to show differences in these domains, so these are  
12 hypothesized -- hypothesis-generating data. But you can see there actually was a  
13 statistically significant improvement in quality of life in the Reducer-treated patients  
14 compared to the sham-control treated patients, and there were directional improvements  
15 in angina stability and angina frequency, consistent with the primary endpoint results.  
16 Physical limitation, again, in the right direction, but the improvements were smaller. And  
17 treatment satisfactions were actually relatively high in both groups beforehand and there is  
18 no difference between the two groups at 6 months.

19 Now we go to some of the, again, hypothesis-generating underpowered exercise test  
20 data. You can see here -- and for time, I'm just going to show you some selected data, but  
21 here if you can look at the exercise testing, which was bicycle ergometry done by the ACIP  
22 protocol. You can see the Reducer from baseline to 6 months, those patients had on  
23 average a 65-second improvement in exercise time, whereas the control patients basically  
24 had absolutely no difference in exercise time from baseline to follow-up, and this trend did  
25 not quite reach statistical significance.

1           Here we look at the dobutamine stress wall motion score index, a negative number  
2 here is a better number. This is looking at the change from baseline to 6 months and you  
3 can see the Reducer treatment, in the dark bars, tended to have a greater change in wall  
4 motion score compared to the sham control. If you look at the right set of bars, that's  
5 probably the most relevant set because that is just in the left coronary artery distribution,  
6 which is really the major blood flow egress the Reducer would impact. Here you can see  
7 there tended to be a doubling of the improvement in the wall motion score, but again,  
8 these numbers aren't close to what was needed for statistical significance and so these  
9 values are NS but it is, again, suggestive.

10           So now let me go and look at the post-approval REDUCER-I study and look at the  
11 effectiveness outcomes here. There were actually three different cohorts of patients that  
12 were included in REDUCER-I. This whole experience would be up to 400 subjects, it's  
13 currently enrolling. Two hundred forty-one patients were enrolled as of March 12th for  
14 which we have the data on. The prospective arm, if you will, of REDUCER-I is called Arm 1  
15 and these are new Reducer patients that were included in this formal study. These patients  
16 had similar inclusion and exclusion criteria as the COSIRA randomized trial, but it did also  
17 allow CCS Class II patients to be enrolled.

18           The Arm 2 patients, there were 11 patients where the patients in COSIRA were  
19 actually consented and were followed long term in this registry because this is a long-term  
20 registry.

21           And then there was Arm 3, which were patients implanted under CE mark before this  
22 formal REDUCER-I registry started. Those patients were also consented and added to the  
23 REDUCER-I observational study.

24           So overall, 241 patients; and these patients underwent assessment at baseline for  
25 the procedural outcomes, 30 days, 6 months, and then annual follow-up through 5 years.

1 And all their adverse events were adjudicated by a clinical events committee.

2 So if we compare the REDUCER-I patients in the middle set of data to the COSIRA  
3 Reducer arm patients, these 241 patients, you can actually see their demographics were  
4 very similar in terms of their age, their sex, about 20% women. Again, half prior MI. Almost  
5 all had prior revascularization procedures. You can see the biggest difference here is that  
6 we now have about 29% of the patients in REDUCER-I had Class II angina, whereas in  
7 COSIRA they all had Class III or IV angina. But again, diabetes, smoking,  
8 hypercholesteremia, hypertension, very similar.

9 Now, these are the outcomes data, this is just now the REDUCER-I study. These are  
10 patients that had at least one CCS class angina improvement, and these are paired data.  
11 Given the data that's available from baseline to 6 months, baseline to 12 months, and  
12 baseline to 2 years, you can see the number of patients. And if we just focus here on -- let's  
13 start with a one-class or more improvement in angina, you can see that it starts at 70% at 6  
14 months and it slightly increases to 82% at 2 years. If you look at greater than or equal to  
15 two-class improvement, it's 24% at 6 months and it slightly increases to 31% at 2 years.

16 Now, if you look overall here at comparing these data to COSIRA, you can see we've  
17 got baseline and 6 months and the upper set is COSIRA and the lower set is the REDUCER-I.  
18 But for here we've just done the REDUCER-I and Class III or IV patients. And so if you focus  
19 on the dark set of bars, that might be the best way to go, now you're looking at Reducer to  
20 Reducer, you can see here that looking at the Class III/IV patients only for REDUCER-I at  
21 baseline, 93% of those were Class III, 7% were Class IV, and you can see that the Class I, II,  
22 III, IV at 6-month follow-up, 30% and 44% respectively were Class I and Class II, and these  
23 are really very similar to what we saw in the Reducer arm in the sham-controlled COSIRA  
24 trial. Slightly less Class IVs in Reducer compared to COSIRA. That might be because there  
25 were less Class IVs in REDUCER-I compared to COSIRA.



1           And again, these are all now the angina data on all angina classes in just REDUCER-I,  
2 all the way out through 5 years now. You can see, once we get to 3, 4, and 5 years, the  
3 number of patients start to get small, but these are all the data. And if you focus again on  
4 the -- I think the blue and the green segments of the bars, that's the Class I and Class II  
5 angina, you can see that the vast majority of the patients have Class I and Class II angina  
6 and you can get a sense that there seems to be a persistent effect all the way through  
7 5-year follow-up.

8           This slide looks at, again, REDUCER-I but only the Class III or IV angina patients at  
9 baseline. This is the group, of course, that we're looking for labeling for, the ones who can  
10 really benefit the most from the Reducer. So if you look at these patients at baseline, of  
11 course, they all had Class III or IV angina and now you can see in Class I all the way out  
12 through 5 years approximately 10 to 20% of the patients had Class III or IV angina. So  
13 approximately 80% of these patients no longer have Class III or IV angina after Reducer  
14 implant. And of course, this is open label compared to COSIRA, that is its limitation, but  
15 nonetheless it looks quite effective.

16           Looking at the Seattle Angina Questionnaire assessments that were made in  
17 Reducer, this is the prospective arm only, they're the ones who got the Seattle Angina  
18 Questionnaire and you can see the improvements in -- now using the patient as their own  
19 control from baseline to 6 months, you can see substantial improvements in all five  
20 domains, quality of life, angina stability and angina frequency, with somewhat less  
21 improvements, but still highly statistically significant, in physical limitations and treatment  
22 satisfaction.

23           Going now just to the Class III or IV angina patients and thus comparing the COSIRA  
24 SAQ scores to the REDUCER-I SAQ scores, this is all Class III or IV in both groups, so it is  
25 apples to apples, you can see the data look actually very similar. You can see, remember in

1 COSIRA, there was a statistically significant improvement in quality of life. You can see a  
2 substantial improvement in quality of life in REDUCER-I from baseline to 6 months with the  
3 Reducer, as well, and you can see in all of the domains the outcomes tend to be greater in  
4 both studies with the Reducer. Perhaps the one exception is treatment satisfaction, which  
5 again tends to be high. So the patients were basically reasonably treated, that their doctors  
6 were doing a good job, as good as could be done at baseline; hard to improve that, but  
7 there was a little bit of improvement with Reducer treatment in both groups.

8 An interesting metric that was collected in the prospective arm of REDUCER-I is to  
9 look to see if there was a reduction in emergency department visits and you can see here,  
10 before the Reducer was implanted in these 113 patients, that they've now reached the  
11 1-year post-Reducer implant, they were having an average of 0.7 emergency department  
12 visits per year and that was reduced to 0.2 within 1 year after the Reducer, a highly  
13 statistically significant number. So basically, it was 42% of the patients had visited the ER  
14 before the Reducer and now it's reduced to 13% after the Reducer, a total of 78 ED visits  
15 compared to 22 after the Reducer.

16 In addition to this data that I've shown you, there have been other investigator  
17 sponsored separate studies that have been reported and I show you now unselected data,  
18 all the studies that have been published with 39 or more patients in which they reported  
19 CCS angina class. If you look at the dark bars, the first set is COSIRA treatment with the  
20 Reducer, looking at a two or more greater CCS angina class at 6 months. Remember, that  
21 was 35%. In REDUCER-I, it was 24%. But again, probably a little less because there were  
22 Class II patients enrolled. If you look, though, at these other studies at 6 months, 14  
23 months, 18 months, and 2 years respectively, you can see 36 to 50% improvements in two  
24 CCS angina classes.

25 This is the same data, same studies, now looking at one or more CCS angina class

1 improvement and you can see that's about 70% in both the COSIRA treatment group and  
2 the REDUCER-I implant group studies, and that Class I is about 76 to 85% in the four  
3 investigator separately reported studies. So again, consistent effectiveness.

4 And you see the same thing if you look at the Seattle Angina Questionnaire  
5 improvements. You can see the improvements in both the treatment arm of the COSIRA  
6 trial and the REDUCER-I study and similar or slightly greater improvements in the  
7 investigator reported studies.

8 So if I was going to summarize the results of effectiveness about the Reducer  
9 implant, most importantly, the sham-controlled COSIRA trial met its primary pre-specified  
10 powered endpoint. A two-class or greater improvement in CCS angina class from baseline  
11 to 6 months was present in 35% of Reducer patients versus 15% of sham-control patients,  
12 statistically significant with a number needed to treat of about five patients, which is quite  
13 a low number, as you know.

14 The secondary non-powered endpoints were directionally supportive of the primary  
15 effectiveness results. The open-label REDUCER-I study demonstrated effectiveness with  
16 sustained response during long-term follow-up and the results were consistent in the  
17 COSIRA, REDUCER-I, and independent investigator-led studies.

18 So now let's look at the safety results. I remind you about the two pre-specified  
19 COSIRA safety endpoints, I won't read them again. These were not powered endpoints but  
20 were pre-specified.

21 So if we look to start at all serious adverse events within 6 months that were  
22 reported from the COSIRA study, in the Reducer group there were six patients or 12% that  
23 experienced one or more SAEs, a total of 10 SAEs. In the no Reducer group there were 11  
24 patients or 20.4% of patients that experienced a total of 24 SAEs. So this trend towards  
25 fewer SAEs in the Reducer group was not statistically significant, but you could see it was

1 more than a 50% reduction in the number of SAEs with the Reducer compared to no  
2 Reducer.

3 If we look at the timing of the SAE reports, the first pre-specified safety endpoint  
4 was these SAEs within the first 30 days. You could see, prior to discharge, there were no  
5 SAEs in either group. In the discharge, the 30-day period, each group had two SAEs, so  
6 obviously very similar. Looking at 30-day to 3-month and 3-month to 6-month, this is when  
7 you see fewer SAEs after the Reducer compared to the sham-control group.

8 This lists all the different SAEs. You can see the 10 total SAEs in six patients in the  
9 Reducer group, the 24 total SAEs in 11 patients in the not Reducer group, and you can see  
10 the biggest differences here were less angina, chest pain, and unstable angina that were  
11 considered serious adverse events in the no Reducer group compared to the Reducer group  
12 -- maybe two versus zero acute coronary syndrome.

13 So this slide shows the second pre-specified safety endpoint from the COSIRA sham-  
14 controlled randomized trial, and that's the rate of major adverse events or also called major  
15 adverse cardiac events after 30 days and up to 6 months. And here you can see in the 50  
16 Reducer-implanted patients there was only one MACE event or a 2% rate in the no Reducer  
17 patients, of which there were 54. There were four MACE events or 7.4%. And as I  
18 mentioned before, of the two Reducer patients that were initially assigned to the Reducer  
19 that did not have an implant, none of those patients had a MACE event. So the four MACE  
20 events were all in the 52 original patients assigned to the sham-control group. The one  
21 event in the Reducer group was a myocardial infarction between 30 days and 6 months.  
22 The four events in the patients who did not get a Reducer were cardiac death in one patient  
23 and three myocardial infarctions.

24 Now going to the REDUCER-I study where we can look at that same endpoint, MACE  
25 within 5 years, you can see now that there was a 15.7% rate of MACE within 5 years.

1 Looking at cardiac deaths, the cardiac death rate within 5 years was 5.2%, so about a 1%  
2 death rate per year in these 228 patients by Kaplan-Meier estimates and that's pretty  
3 remarkable given this extensive coronary artery disease and "no option" patients, 1%  
4 cardiac death rate per year. There were four total strokes and 16 patients had a myocardial  
5 infarction, approximately again a 2% or 2.4% per patient myocardial infarction rate per  
6 year. Only one of these events were considered as possibly being procedure related and it  
7 was a myocardial infarction that it was really adjudicated to be unknown whether or not it  
8 was related to the device or not.

9 Thus, looking at the timing of the MACE over time in the REDUCER-I, if you look, I  
10 think, most importantly in the periprocedural period, the first 30 days, there were no  
11 cardiac deaths, no strokes, and there were four myocardial infarctions out of the 241  
12 patients. And then looking all the way out to 3 years you can see that a few patients are  
13 having events basically each year of each type, with myocardial infarction being the most  
14 common, but even these were very infrequent.

15 Now looking at the 6-month rates of MACE across all the Reducer studies, so this is 6  
16 months looking at, to put it in perspective, COSIRA, the sham-control arm in the first data  
17 column; then after the vertical line, the COSIRA Reducer arm 6-month outcome; the first in  
18 human, 15 patients, 6-month outcome, of which there were no events; the REDUCER-I  
19 6-month outcome excluding the 11 patients that were in the COSIRA group. So you have a  
20 total of 295 patients in yellow and at 6 months you can see there was no cardiac death,  
21 there was one stroke, there were eight MIs, and two procedure-related adverse events.

22 The independent studies also support Reducer safety and again, here I show you on  
23 the five lines here, all the safety data that have been reported with all the individual  
24 investigator sponsored and reported studies. You could see here, from Konigstein, Giannini,  
25 Gallone, D'Amico, and Ponticelli, that the rates of cardiac death, stroke, and myocardial

1 infarction have been quite low in these patients; in general, were consistent with what we  
2 saw in COSIRA and REDUCER-I.

3         As we mentioned, these truly are "no option" patients; they really have, despite  
4 having severe angina, no revascularization option, which of course would be the first way  
5 that we would want to treat them. You could see in COSIRA, looking at the Reducer arm,  
6 the no Reducer arm, not a single Reducer patient had a revascularization within 6 months.  
7 Two no Reducer patients did have an attempted revascularization procedure, a PCI, because  
8 of just severe refractory angina. In the REDUCER-I study, out to 5 years you could see that  
9 there's been about a 12% rate of attempted PCI or about 2% per year and no patient had a  
10 bypass surgery.

11         And then finally, looking at the coronary sinus related events, there is a learning  
12 curve, of course, to using the device. So we show you here the premarket outcomes in  
13 COSIRA in 52 attempted patients with the Reducer, the postmarket REDUCER-I study in 204  
14 patients, and the postmarket commercial experience in 1840 patients. I will mention, the  
15 1840 patients, this is a voluntary reporting system of these 1840 distributed devices, so it's  
16 up to the sites to report whether complications occurred. Although most of the  
17 complications on this slide are quite serious and related directly to the Reducer, so I would  
18 anticipate that most of them would've been reported, but we don't know for sure, there  
19 could be a component of underreporting here.

20         You could see here that the incidence of Reducer malposition, migration, or  
21 embolization with two patients in COSIRA, two patients in REDUCER-I, and there have been  
22 10 patients or 0.5% incidence reported from the commercial experience.

23         Looking at coronary sinus dissections and perforations and fractures, dissections and  
24 perforations have been very rare. There have been only two cases of tamponade that have  
25 been reported. There have been no fractures reported, no coronary sinus thrombosis or

1 occlusion reported, and major bleeding events have been uncommon.

2           And then finally, we have the first-in-human follow-up safety data that goes all the  
3 way out to 12 years, that's the longest experience with the Reducer, and specifically repeat  
4 CTAs were performed in seven of these patients at 12 years in whom follow-up was  
5 completed, and all of these devices were patent, no thromboses were detected, no device  
6 fractures were observed, and no migrations were evident.

7           So in summary, while I could go down the litany of everything I just presented, I  
8 think I can basically conclude the totality of the evidence in more than 2,000 distributed  
9 Reducer devices has demonstrated that the Reducer implant is safe when used as indicated.

10           So if I were to try to summarize all of the data that I've shown you, the effectiveness  
11 and the safety, there's no doubt that it would be wonderful if we had a 500-patient sham-  
12 controlled trial with definitive evidence of reduced myocardial ischemia with more of a  
13 contemporary kind of clinical trial metrics that we're using in the last few years. But again, I  
14 do want to emphasize that what's important is to try to place these data in perspective in  
15 terms of, you know, have we demonstrated reasonable evidence of effectiveness and  
16 safety, in particular in the perspective of the unmet clinical need. And these really are  
17 patients who are desperate and would have been through multiple revascularization  
18 procedures, who are taking, in many cases, multiple medications and have just a miserable  
19 life with no therapeutic alternatives. And I believe that the totality of the data, in sum, has  
20 shown that we have a reasonable confidence that this is an effective device and it's quite  
21 safe for these patients with such an unmet clinical need.

22           So I'm going to turn it back over to Tim Henry now, who is going to give you his  
23 clinical perspectives.

24           DR. HENRY: Thanks, Gregg. And for the clinical perspective, I will keep my  
25 comments brief.

1           As noted before, I've spent more than 20 years taking care of patients with  
2 refractory angina, including the development of a specific refractory angina clinic, not only  
3 in Minneapolis and Los Angeles, but now at The Christ Hospital in Cincinnati. This includes a  
4 prospective database with more than 2,000 people. I will tell you that this is a particularly  
5 challenging patient population and they really have extremely limited options. I think the  
6 comment of "no option" patients is really appropriate for them.

7           I've also had extensive clinical trial experience with refractory angina, including  
8 novel pharmaceutical agents, EECPP, spinal cord stimulation, therapeutic angiogenesis with  
9 protein, gene and cell therapy, and then novel interventional techniques including CTO,  
10 shockwave therapy, as well as the Reducer. I've been a national PI or steering committee  
11 for nearly 20 different trials and I will tell you that doing these clinical trials is very  
12 challenging. It's challenging because when you carefully define refractory angina, as I tried  
13 to do earlier, it becomes a relatively small patient population and it's very difficult to find  
14 sites that are good at recruiting these patients.

15           For example, if you take Class III and Class IV patients who are not candidates for  
16 revascularization and who have more than seven episodes of angina per week, you're left  
17 with a very, very specific patient population that's very difficult to enroll. And this is  
18 reflected in the fact that it's been more than 20 years since a new therapy was approved for  
19 this challenging patient population.

20           So let me review the key points. Although angina is common, refractory angina  
21 defined as Class III or Class IV is not and potentially should be considered an orphan  
22 population.

23           By definition, "no option" patients have failed medical therapy and they have failed  
24 current revascularization techniques, both percutaneous and surgical.

25           We clearly have improved mortality. Realistically, it's 3 to 5% per year with



1 contemporary therapy. Slightly higher, maybe 5 to 6% percent per year if you have Class III  
2 and Class IV patients.

3 Importantly, these patients then are left with severe debilitating symptoms which  
4 interfere with their job, interfere with their family life, interfere with their ability to do  
5 activities of daily living, and as you will hear later today during the open mike, these  
6 patients are miserable and they're desperate.

7 Finally, for "no option" patients in the United States, the only therapeutic options  
8 currently available are EECF and most patients have either done that already or it's not  
9 logistically possible for them to do it, and TMR, which is seldom used in the United States  
10 for the reason we talked about before. So this really clearly leaves a majority of patients  
11 with absolutely no options at all. Again, truly "no option" patients.

12 So refractory angina patients in the United States are desperately in need of new  
13 options and I can tell you, I deal with this every week in clinic. It's important that these  
14 options are safe because we know the mortality rate is low and so you want to avoid risky  
15 therapies. Therapy also has to be effective in improving symptoms of refractory angina and  
16 it has to work. So from my perspective, the coronary sinus reducer meets those criteria. In  
17 particular, amazing, outstanding safety now in over 2,000 people that have been treated  
18 around the world and with selected but significant efficacy as you've seen from COSIRA-I.  
19 These results have been supported by ongoing registry data and in COSIRA-I, if you  
20 remember, there is clear significant improvement in angina and quality of life, and the  
21 subsequent extensive registry data and publications are very consistent in both safety and  
22 efficacy.

23 So I personally believe that the coronary sinus reducer is an option that should be  
24 available for "no option" patients in the United States who are truly, truly desperate. So I  
25 appreciate and I thank you for your time and attention.

Free State Reporting, Inc.  
1378 Cape Saint Claire Road  
Annapolis, MD 21409  
(410) 974-0947

1 DR. LANGE: Great. I believe that concludes the Sponsor's presentation and I would  
2 like to thank the Sponsor's representatives for their presentation.

3 And this is an opportunity for the Panel to ask any brief clarifying question. We have  
4 about 25 minutes dedicated to this before our break. I will look for hands to go up and I'll  
5 be multitasking. I see several. I see Bram, Dr. Somberg, Keith Allen, Rick Page, David Yuh.  
6 I'll try to do those in order and I'll just remind you to unmute.

7 So first, Bram, to you.

8 DR. ZUCKERMAN: Yes, thank you, Dr. Lange.

9 First of all, I want to thank the Sponsor for an excellent presentation, but I do want  
10 to note for the record that if this had been done in person instead of virtually, I would've  
11 stopped the presentation at several critical points. Specifically, our rules for sponsor  
12 presentations are such that it's standard practice that when data are being presented in a  
13 sponsor presentation that has not previously been reviewed by the FDA, the presenter  
14 notes that just to help the Panel properly interpret the strength of the data.

15 From my notes, I believe that, for the record, Slide 64, 76, 78, 79, 92, 94, and 96  
16 contain data not previously reviewed by FDA. That's fine, but I would ask the Panel  
17 members, in their independent review, to properly consider this fact. Also, I would ask the  
18 Sponsor, in responding to subsequent Panel questions, to clearly demarcate whether or not  
19 the FDA has had a chance to independently review data being referred to. Finally, another  
20 standard practice that I would ask the Sponsor to adhere to is to indicate more clearly  
21 whether or not certain hypotheses were pre-specified and/or whether multiplicity  
22 adjustments were made. While this information is to be found on the slides, it's in very  
23 small type and many of us, including myself, have poor vision at our advanced age.  
24 Anything that the Sponsor can do to more clearly present the science would be greatly  
25 appreciated. Thank you.

1 DR. LANGE: Thank you very much, Bram.

2 And again, the Sponsor has an opportunity -- many of these questions, clarifying  
3 questions, we'll ask the Sponsor if there's additional data to acquire that over the lunch  
4 period or during the public portion, public presentation. But apropos to Dr. Zuckerman's  
5 comments, in the future, any slides referred to or any data, please specify whether the FDA  
6 has approved it and whether it was -- hypothesis was pre-specified and multiplicity  
7 examined.

8 Thank you, Bram.

9 Dr. Somberg, I think I saw your hand first. Go ahead, John, you're unmuted.

10 DR. SOMBERG: Okay, yes, John Somberg here.

11 I failed to mention that my area of expertise is cardiovascular pharmacology and I'd  
12 like to drill down upon that because it's repeatedly stated by the Sponsor that these are "no  
13 option" patients. So can they define how they went through over that 30-day period, a  
14 sequential test if pharmacologic therapy for angina was maximized? Or made maximum, I  
15 should say. Could they also address why a good third of the patients are on blood pressure-  
16 lowering therapies that are not antianginals such as ARBs, ACEs, and diuretics? Could they  
17 also mention why or if the medications were remaining constant throughout the study in  
18 COSIRA?

19 And I'm specifically talking about the randomized sham-controlled study. So can we  
20 drill down on why these are considered "no option" patients, pharmacologically? And also  
21 was there any review of the cardiac catheterization material by a central committee to say  
22 they had no interventional option? I'm trying to define this population, are they really drug  
23 resistant, are they really procedurally excluded? Thank you.

24 DR. LANGE: Great. Dr. Somberg, again, just for the Sponsor, when we come back  
25 from lunch, the things you mentioned, again, to discuss the 30-day sequential

1 pharmacologic maximization. Explain why a third of the patients weren't on antianginals or  
2 were on blood pressure medicine that were not antianginals, whether the medicines remain  
3 constant and whether the dose was maximized and finally, a review of cardiac  
4 catheterization data.

5 Thank you, Dr. Somberg.

6 Dr. Page, I saw your hand, as well.

7 DR. PAGE: Yes, thank you, Dr. Lange. Rick Page, for the record.

8 Two quick questions in terms of the -- I'm trying to get a handle on the mechanism of  
9 action here and the physiology that underpins a procedure that these patients with no hope  
10 are undergoing. As I read the information, I believe that that procedure, in addition to  
11 constricting the coronary sinus, involved pericardial irritation and wonder if there's any  
12 response to the fact that the response in patients was, for sure, related to the coronary  
13 sinus stenosis as opposed to the pericardial procedure.

14 The other thing is looking at Slide 45, and I believe Dr. Hampshire spoke about the  
15 intentional nature of making this hourglass stent, if you will, and I'm really struck by Animal  
16 1111, which has a reverse hourglass in terms of the angio that I would imagine one would  
17 see in that it's more narrow at the area of the wide part of the hourglass stent and it's more  
18 open at the area where it would be narrower. So it's almost a reverse hourglass or a double  
19 hourglass effect, and wonder whether this stent is -- how they reconciled that with the  
20 shape of the stent, was this intentional or is this just observational and does this procedure  
21 require two stenoses? Are we interested in the pressure gradient at one -- at what we'll call  
22 the distal end and then at the drainage to the more proximal coronary sinus? I'm just  
23 interested if they can comment on the physiology behind what the stent is doing to  
24 pressure in the coronary sinus.

25 DR. LANGE: Great. So again, to summarize Dr. Page's two questions, the presenters

1 specified that the Beck surgery was due to coronary sinus occlusion or stenosis, but  
2 obviously there's pericardial irritation, how do they know that it is due -- the benefits are  
3 due to coronary sinus stenosis, and asking about Slide 45 and the reverse hourglass where  
4 there's some -- thank you very much.

5 Dr. Yuh, did you have your hand up, as well?

6 DR. YUH: Yes, thank you very much. David Yuh.

7 Two quick questions. I was struck in Dr. Henry's presentation that one specific group  
8 of anginal patients excluded intentionally were those suffering from microvascular disease  
9 and I was curious, because it seems like the putative mechanism behind this therapy might,  
10 in theory, benefit those patients the most and I wonder if that exclusion was at least  
11 partially responsible for the fact that so few women were enrolled in the study, you know,  
12 based on the preponderance of microvascular disease in that population.

13 And then secondly, in terms of the sham design, if I could get an approximate  
14 average time for deployment of the actual device, was it substantial enough that there was  
15 a significant time discrepancy, a total procedural time discrepancy between the sham and  
16 the study -- the treatment groups, whereby it might be discerned one way or the other  
17 what therapy or if the patient got or did not get the stent? Those are the two questions  
18 that I had.

19 DR. LANGE: Okay. So Dr. Yuh, again just to summarize, the discussion of exclusion of  
20 microvascular disease and whether that led to the low number of women in the study. And  
21 secondly, procedure time, sham versus treatment to see if there's a difference.

22 Okay, thank you very much.

23 Dr. Allen. And then I've got Dr. Cigarroa and then Dr. Wittes. Dr. Verghese, do you  
24 have your hand up? I'll get you next after that.

25 So, Dr. Allen, you're on mute right now. And then Magnus, okay. And then Ralph.

1 DR. ALLEN: Sorry about that. This is Keith Allen. Thank you all for letting me  
2 participate.

3 I've got several questions, they're relatively short. The first one is, were surgeons,  
4 coronary surgeons, involved in deeming these patients no option? I see an investigator who  
5 I assume is a cardiologist who can deem a patient not a candidate for PCI, but I don't ever  
6 see surgeons involved in deeming these patients no option.

7 Second question is to the effectiveness of blinding. While I applaud the team on  
8 doing a sham-controlled trial, I don't ever see any simple logistics confirming that patients  
9 actually stayed blinded. All you would need to do would be to look at a chest X-ray and you  
10 could see the device in the coronary sinus, so how is that done?

11 The third question relates to this and as to who was managing these patients' post-  
12 op medications. And it's important that medications didn't change in one group more  
13 significantly than the other, and I don't see that presented anywhere in the randomized  
14 trial. And it's also important that the person managing the medications was blinded. We  
15 only have to look back at, for example, CardioMEMS and the issues that had with being kind  
16 of overzealous and biased management of the treatment group that led to more favorable  
17 results.

18 And then the final question goes to pressures. I'm astounded by the fact that not  
19 once in the presentation do I actually ever see any validation that the hypothesis of how the  
20 Reducer works was documented, i.e., where are pressure measurements done in any animal  
21 that documents that? And I understand that they mention that quadrupeds have a problem  
22 doing that, but I also don't see any human data. And I would argue, as a vascular surgeon,  
23 that constriction of a vein, unlike an artery, leads to dilation of that vessel and the  
24 constriction of the coronary sinus will behave very differently and may not yield the results  
25 the Sponsor is actually looking at to validate their hypothesis of how this works.

1 DR. LANGE: Great. Okay, to summarize for the Sponsor, were surgeons involved in  
2 evaluation of revascularization? Talk about the confirmation of the effectiveness of  
3 blinding, please. Please discuss who managed the patient post-op and the number of  
4 medication changes made in the sham and control group or sham and treatment group.  
5 And then finally, provide pressure measurements in humans, please.

6 I've got right now Dr. Cigarroa, Dr. Wittes, Dr. Verghese, Dr. Ohman, and Dr. Brindis.  
7 So Dr. Cigarroa. And then Dr. Borer.

8 DR. CIGARROA: Good morning from the West Coast, and thank you to the Sponsors  
9 for the presentations, and to our Chairman for providing me the opportunity for clarifying  
10 questions.

11 I have three. The first question is, were individuals who assessed the functional  
12 capacity, including angina, during the follow-up period blinded to medications, including  
13 whether or not an individual was on dual antiplatelet therapy or not, is my first question.

14 The second question goes to the point that, by my calculations, there were 34 of 52  
15 in the active arm that did not have any change in angina versus 44 over 52. Were there any  
16 differences in diastolic filling period between those who responded and those who did not  
17 respond during this period?

18 And my final question is a follow-up to Dr. Allen's question with regards to pressure  
19 measurements and, as we know, change in pressure across a reduced cross-sectional area is  
20 influenced by flow and therefore understanding whether or not there is a pressure change  
21 across the Reducer at rest versus during increased heart rate would also be important to  
22 understand relative to this biologic plausibility of increasing CS pressure.

23 DR. LANGE: Okay. So to summarize, again, the three questions Dr. Cigarroa asked,  
24 and correct me if I misstate these. Are the people that assessed angina, were they blinded  
25 to medications and particularly to dual antiplatelet therapy? Was there any difference in

1 diastolic filling period between those treated and those who weren't, that had no change in  
2 their angina, no improvement? And then any measurements across the stent either at rest  
3 or with increased heart rate.

4 Did I summarize those, Dr. Cigarroa, adequately? Okay, great.

5 DR. CIGARROA: Yes, sir.

6 DR. LANGE: Great. Dr. Wittes.

7 DR. WITTES: Hi, this is Janet Wittes.

8 I have four questions, they're all pretty short, I think. The first is what countries  
9 were these, the 104 patients from? That's question one.

10 The second question is what percentage were current smokers? You showed 50%,  
11 more than 50% current or former smokers, but I would like to see what -- because this, the  
12 50%, seems really high for current smokers in the U.S.

13 The other question is can you show Class IIIs had to move to Class I, the Class IVs had  
14 to go to II or I? Can you show the movement in both treatment and control, where the IIIs  
15 and IVs went?

16 And then finally, the more complicated question comes from Slide 75, 77, 78, and 81  
17 and there may have been a few others, as well. You show a number of people at 1 year,  
18 some at 2 years, some at 5 years, and then you calculate the percentage of various things as  
19 a function of those denominators, but you didn't show how many people had reached 1  
20 year, 2 years, 5 years. And so my question is did the number and the percentage that  
21 you're showing reflect the entire group of people that reached those landmark number of  
22 years or is it a selected group in some way and there's missing data?

23 DR. LANGE: Okay, thank you.

24 To summarize Dr. Wittes' comments, what countries enrolled a hundred and four  
25 patients in COSIRA, what percent are current smokers, please show the movement of Class



1 III and IV patients in the treatment and control groups. And then in the last Slides, 75, 77,  
2 78, 81, do those patients represent the entire group at 6 months, 1 year, and 2 years?

3 Thank you, Dr. Wittes.

4 Okay, I've got Dr. Verghese, Ohman, Dr. Brindis, Dr. Borer, and then  
5 Dr. Connor and then Dr. Somberg, okay?

6 I think, Dr. Verghese, I'm going to turn it over to you.

7 DR. MATHEW: Dr. Lange, it's actually Dr. Mathew. Verghese is my first name.

8 DR. LANGE: I'm sorry, I'm sorry.

9 DR. MATHEW: That's all right, just for the record.

10 A couple of quick points, one to extend on Dr. Somberg's comment. Is there any  
11 information on -- there was a mention, not necessarily in this morning's presentation, but  
12 the advanced materials, that about 25% of patients were on zero or one antianginal  
13 therapy. That would be somewhat inconsistent with the refractory group, so if you can  
14 expand on that.

15 And then is there information on changes to or adherence to medication therapy  
16 over the course of follow-up in COSIRA at 6 months or the REDUCER trial out to 5 years,  
17 either increases or decreases in medical therapy in the respective groups of patients.

18 Second, I would echo sort of what three panelists have made comments on already  
19 in terms of mechanism. While mechanism may not be a requirement for approval in this  
20 process, I understand, I'm having difficulty drawing the lines between mechanism and the  
21 clinical benefit that's been attributed to that mechanism.

22 The third is the patient population with microvascular angina, which doesn't seem to  
23 have been included in this population. I would just make the statement that the IFU, while  
24 it includes patients who are not candidates for CABG or PCI, one might infer that it doesn't  
25 specifically exclude patients with microvascular angina. So the proposed IFU, I would just

1 say myocardial modification.

2 Then the last point, also in the advanced materials but not mentioned this morning,  
3 there was a categorization of technical mishandling of the device in eight out of the 50  
4 COSIRA patients. If the Sponsor could expand upon on that and whether that should be of  
5 concern in extrapolating this procedure to a broader number of interventionalists and  
6 patients. Thank you.

7 DR. LANGE: Great. So again, 25% on zero or one medication don't appear to be  
8 refractory, please address that. Again, discuss the change in medications, adherence to  
9 medications in both groups in both COSIRA and REDUCER. Please address whether anybody  
10 with microvascular disease was enrolled, i.e., what did the angiography show? And then  
11 finally, please address the 8 over 50 individuals in which there were technical difficulties  
12 handling the device and whether it ought to be a concern.

13 Did that address -- summarize, Dr. Mathew?

14 DR. MATHEW: Indeed. Thank you, Dr. Lange.

15 DR. LANGE: Thank you very much, sir.

16 Dr. Ohman.

17 DR. OHMAN: Good morning, and this is Magnus Ohman for the record. I want to  
18 thank the Sponsor for an excellent presentation.

19 I have three fairly brief questions. In the original protocol of this device there were  
20 124 patients, I believe, to have been randomized but only 104 actually did so. So the  
21 question I have, is this trial subject to early termination bias, i.e., that we're seeing more  
22 positive results because they were terminated early? It might be useful for the Sponsor to  
23 calculate a worst-case scenario, i.e., for the full sample size by applying the placebo rate or  
24 the control group rate to the remaining COSIRA patients, not the full sample size, just to get  
25 a feel for how sensitive this result is due to the smaller sample size than planned.

1           This is an extension, my second question is an extension of Dr. Somberg's question. I  
2 would like to see the COSIRA results with the 25% who were on minimal therapy excluded  
3 to see if there's any remaining treatment effect in this study.

4           And finally, my third question is in going to the Sponsor's Slide 90 and 91, I counted  
5 up two myocardial infarctions versus one. So here on 90 you see acute myocardial  
6 infarction on top and a little bit further down you see myocardial infarction again, and then  
7 on the next slide there's only actually one given, I believe, if I saw the slide correctly. So the  
8 question here, is that accurate? And of course, if we're talking about two out of 50 versus  
9 one out of 50, there's a fairly big difference in this small sample size for the confidence  
10 interval around such observations. Thank you.

11           DR. LANGE: Thank you, Dr. Ohman.

12           Again, to ask the Sponsor for the sensitivity analysis, considering the possibility of  
13 early termination on account of some of the effects. Exclude the 25% of COSIRA patients  
14 that were on zero or no medications to see if there was still an effect. And then to  
15 reconcile the number of MIs in Slides 90 and 91.

16           Okay, thank you, Dr. Ohman.

17           Dr. Brindis.

18           DR. BRINDIS: Thanks. And again, I want to thank the Sponsor for a terrific  
19 presentation.

20           Most of my questions have already been posed by my esteemed colleagues. My last  
21 one is focusing on safety. So the Sponsor's presentation stated that out of approximately  
22 2,000 follow-ups there's been no coronary sinus thrombosis noted to date. In some of the  
23 materials sent to us there was a CT scan or MRI scan evidence of thrombus, not insignificant  
24 within the occluder device. So I was interested in more comments related to that and also  
25 the utilization of antiplatelet therapy long term, dual or -- appreciating, of course, that a lot

1 of these patients with underlying coronary disease would be on antiplatelet therapy.

2 DR. LANGE: Okay. So the two issues, Dr. Brindis, one is the presenter made a  
3 comment of 2,000 follow-ups, there was no CS thromboses, please explain how that was  
4 documented. And then address the issue of dual antiplatelet therapy and how long that  
5 was continued. Great.

6 DR. BRINDIS: Yes, and the issue of presence of thrombus within the device --

7 DR. LANGE: Great.

8 DR. BRINDIS: -- on scanning.

9 DR. LANGE: Okay. Now I'm going to take the last two or three questions and then  
10 we'll have more opportunity to do this. And I'll get to you, Dr. Yeh, before we finish out and  
11 we'll have other opportunities afterwards, but I'll try to keep on time, as well.

12 So I've got Dr. Borer, Dr. Connor, and then Dr. Yeh.

13 Jeff, you're -- that's right, you're muted. Now we've got you. Go ahead, Jeff.

14 DR. BORER: Thank you. The one question I have left that hasn't already been dealt  
15 with relates to Slide 49 and it has to do with what I think is our fundamental issue here.  
16 We're dealing with the capacity of a device to minimize the magnitude of a subjective  
17 symptom in a population that we have to define. In Slide 49, as I understand it, the  
18 judgment as to whether the patients were of limited options was made by the investigator,  
19 who, if I'm understanding correctly, was unblinded and was determining -- outcomes that  
20 were being determined where angina classes improved, which is pretty subjective.

21 So I'd like to have some comment about that. Were the investigators actually  
22 unblinded and would that have any impact on the assessment of the outcomes here? Is it  
23 likely that they would be, because ultimately what we're talking about is assessing the  
24 magnitude of effect of a device on subjective symptoms, so I'd like a comment on that.

25 DR. LANGE: Okay, Jeff. So the question again, were the investigators unblinded with

1 the assessment both of the inclusion of the patient and also the assessment of the  
2 outcome, as well? Great.

3 DR. BORER: Right.

4 DR. LANGE: Okay, Dr. Connor and then Dr. Yeh.

5 DR. CONNOR: Yeah, I'm Jason Connor.

6 Two quick easy questions. So the first, I apologize if I missed, but is there an  
7 explanation as to why the trial stopped short of planned enrollment?

8 And then second question is I would love to see the primary outcomes for the main  
9 trial by site. I know a bunch of sites are small, but you have at least multiple sites with a  
10 decent number of patients to show the sham group, the treatment group by site, just so we  
11 can see the consistency of effect with different providers.

12 DR. LANGE: Great. Okay, Dr. Connor, again to summarize: Why did the trial stop  
13 prematurely, an explanation for that, and then outcomes by site. Thank you.

14 And, Dr. Yeh, you have the final question, at least during this part. And then we'll  
15 take a break.

16 DR. YEH: Thank you. Robert Yeh.

17 Two quick comments/questions. One is I know the trial was conducted, COSIRA was  
18 conducted in a time where maybe modern complex coronary revascularization techniques,  
19 depending on where the trial was conducted and who were the operators that are involved,  
20 may not have necessarily been there. So have the angiograms been reviewed, I'm curious,  
21 by a complex coronary interventionalist or a cardiac surgeon, as Dr. Allen stated earlier, by  
22 2020 standards or by today's standards as opposed to at the time the trial was conducted,  
23 which may be different.

24 The second is a question related to mechanism, which is I know that the trial,  
25 COSIRA, enrolled patients who had ischemia attributable only to the left coronary arterial

1 system, yet there's not a comment about that in the proposed indications for use about  
2 limiting it to only those patients with ischemia attributable to the left coronary system and  
3 how does that -- if the Sponsor can speak a little bit to the mechanism of why that was  
4 limited to that, what should be the expectation based on location of the occluder, of the  
5 placement of the Reducer, etc., related to coronary venous outflow. Just more on  
6 mechanism as it's been previously stated already by many.

7 DR. LANGE: Okay, so that would be the last one. So again -- okay, John, I'll get to  
8 you in just one moment. You'll have the last word, John.

9 Two questions was, were the angiographies reviewed by interventionalists and  
10 surgeons apropos to what's going on? And then finally, mechanism of occluder of the left  
11 versus non-left, to have them address it.

12 So Bram, you'll get the final word, if that's okay. The final final word. So John, you  
13 get the final word. Bram, you get the final final word.

14 So Dr. Somberg.

15 DR. SOMBERG: Well, actually, a question to add is could the Sponsor also discuss  
16 why they decided to use aspirin, a P2Y12 inhibitor. Since it's on the venous side, the  
17 occluder, why wasn't anticoagulation considered for this type of patient? How was it  
18 distinctly different, for instance, someone who has embolizations from the legs and they  
19 put in one of those filters, when you talk about Coumadin, heparin and Coumadin, not  
20 aspirin or P2Y12. Thank you.

21 DR. LANGE: Thank you, Dr. Somberg.

22 And then Dr. Zuckerman.

23 DR. ZUCKERMAN: Yes, I think Dr. Batchelor may have a final question for the  
24 Sponsor. Dr. Batchelor.

25 DR. LANGE: Are you -- there you are, I'm sorry, in the bottom left. I'm sorry I missed

1 you, Wayne. If you'll unmute, you get the last question, Dr. Batchelor. My apologies.  
2 You're muted right now. So, Dr. Batchelor, can you hear me?

3 DR. BATCHELOR: Yeah. Sorry, can you hear me now? Sorry, I was muted.

4 DR. LANGE: We can hear you now, sir. Yes, sir.

5 DR. BATCHELOR: Thanks. And I apologize this morning, with the log-in issues, but  
6 this is Wayne Batchelor from Inova.

7 A couple of quick questions. One that deals with logistical robustness of the results.  
8 You know, a trial like this is very prone to regression to the mean and there were a few  
9 more patients numerically in the Reducer arm that were Class IV angina versus the control  
10 arm. A lot of speakers have discussed the issues around blinding, but also regression to the  
11 mean could play a role here. So I just looked at this quickly and when you look at the  
12 chi-squared tables, in the control arm there were eight who got better out of 52 and that  
13 number is less than 10. So when we can get to these cells that have smaller numbers  
14 you've got to question whether or not to do a Yates correction or other statistical  
15 correction.

16 So my question is can we maybe see if one patient, for example, was misclassified? I  
17 did the calculation and if one was misclassified in each arm, the statistical results would've  
18 been completely removed. And, in fact, if one was misclassified in the Reducer arm, with  
19 the Yates correction also statistical significance would be lost. I think this is important  
20 when you're dealing with a trial that's naturally prone to regression to the mean, placebo  
21 effects, and also blinding concerns. So I'd like them to address that.

22 And then the other question relates to downstream, once this device is in I presume  
23 you can't take it out, so some of these patients are going to need coronary sinus  
24 intervention. So if a patient needs a biventricular AICD in the future, that completely  
25 precludes, I presume, this?

1 DR. LANGE: Okay.

2 DR. BATCHELOR: Thank you.

3 DR. LANGE: Great. So with that, we'll close with those two questions. One is to  
4 have the Sponsor address regression to the mean and just changing the one patient. And  
5 then finally, does this exclude use of CS for future procedures?

6 We're just about 5 minutes behind, so let's take -- instead of a 15-minute break, if  
7 we could take a 10-minute break. We'll call this our bio break and I'll ask everybody to join  
8 -- it's now 9:21, let's join back at 9:31. Thank you.

9 (Off the record at 11:22 a.m.)

10 (On the record at 11:32 a.m.)

11 DR. LANGE: Before we proceed with the FDA presentation, I failed to see two hands  
12 that were up, so I'd like to ask Dr. Starling, to give him the opportunity to pose his questions  
13 to the Sponsor and then Dr. Gersh.

14 So Dr. Starling.

15 DR. STARLING: Yeah, thank you very much, Dr. Lange, and thank you to the Sponsor  
16 for the presentation. This is Randall Starling.

17 And just for the record, I'd like to acknowledge that "no option" patients do have an  
18 option that hasn't been discussed, which is cardiac transplantation and of course, we know  
19 that's a highly restricted option but it is an indicated option for refractory angina. So my  
20 questions fall into three different categories.

21 The first is physiologic mechanisms, which have been addressed already, but my  
22 specific question relates to any pressure measurements as well as coronary sinus oxygen  
23 content, and we know that, and the third physiologic would be coronary vascular  
24 resistance. It seems the mechanism of this is related to coronary flow.

25 Second comment has to do with patient populations for which this is indicated, and



1 in the 2015 *New England Journal* article that was referenced in the appendix, it indicates  
2 that an exclusion is a right atrial pressure over 15 as well as a right atrial or a right  
3 ventricular or a coronary sinus lead. So I would like to know the specific patient  
4 populations with refractory angina that this is going to be applicable to and which are  
5 excluded.

6 And thirdly would be with durability. So we're told how many improve, but any  
7 information related to the durability and movement of the patients, which has already been  
8 alluded to by other questions, would be appreciated. Thank you.

9 DR. LANGE: Great. Thank you for those questions, Dr. Starling.

10 And then Dr. Gersh.

11 DR. GERSH: Thank you, Dr. Lange, I'll keep my comments very short.

12 The first one is to reiterate what others have stated and that is it does bother me  
13 that zero -- I think 25% of patients were on zero or one antianginal agent and how does that  
14 fit with the definition of refractory angina?

15 Number two. Again, it's been pointed out, I'd like to know more about the early  
16 termination. It's hard for me to understand, having served on many DSMBs, why you would  
17 terminate the trial prematurely when 65% of patients in one group didn't respond versus  
18 80% in the other and the endpoint is one symptom, not mortality or morbidity.

19 Thirdly, I think it may have been Dr. Yeh, but I would like to see more -- I'd like to see  
20 the data that show clinical correlations with mechanistic findings. We've seen the animal  
21 studies, and I tried to look them up, and one or two clinical studies have shown  
22 redistribution of coronary blood flow, but it would have been nice to hear a summary of  
23 clinical data in regard to the mechanisms.

24 And then my last question is really just a comment, it's not a criticism, but it is  
25 interesting to me to look at the treatment satisfaction data. Now, I realize that it's a

1 secondary endpoint but, as Dr. Stone pointed out, both groups were pretty satisfied with  
2 their treatment. It is actually a very small difference between the two groups and that just  
3 seems discordant and what it may suggest -- and I'd be interested in their interpretation  
4 and it may be that, you know, with any sham-control trial, 50% of patients will receive the  
5 device and 50% won't and the patients going into this trial know that. So why were they so  
6 happy with their treatment if, in fact, there was quite a marked discrepancy in terms of  
7 angina relief? And I think this just points to some of the difficulties even with a sham-  
8 control trial of ensuring blinding. Thank you.

9 DR. LANGE: Okay. And I think the Sponsor has those.

10 Again, in addition to Dr. Starling's questions, the issue of whether this is a refractory  
11 angina group when a fourth of patients had zero or one antianginal. Why it was terminated  
12 early. What is the correlation between the clinical results and any mechanism. And why  
13 the discordance between angina relief and satisfaction with treatment. So thank you.

14 All right. And again, there will be additional time if we have additional issues that  
15 need to be brought up.

16 At this time I'd like to give the FDA an opportunity to provide their presentation.

17 Again, I would like to remind the public observers at this meeting that while the  
18 meeting is open for public observation, public attendees may not participate except at the  
19 specific request of the Panel Chair.

20 At this time the FDA will also have 90 minutes for their presentation and please  
21 proceed.

22 DR. RABEN: Good morning, my name is Sam Raben and I'm a mechanical  
23 engineer/team leader in the Coronary Interventional Devices Team in the Office of  
24 Cardiovascular Devices in the Center for Devices and Radiological Health at FDA. We  
25 welcome you to the cardiovascular panel of the Medical Devices Advisory Committee

1 meeting to discuss the premarket approval application for the Neovasc Reducer System.  
2 Today's FDA presentation will be presented by myself along with Drs. Anabelle Crusan, Rona  
3 Tang, and Tara Ryan.

4 Briefly, I would like to acknowledge the complete review team for this submission.  
5 These team members helped review all of the information provided by the Sponsor,  
6 including the preclinical and clinical data and analyses.

7 The Sponsor has requested an indication for patients suffering from refractory  
8 angina pectoris despite guideline-directed medical therapy, who are unsuitable for  
9 revascularization by coronary artery bypass grafting or by percutaneous coronary  
10 intervention.

11 Based on this proposed indication for use, the Panel will be asked to comment on  
12 whether a reasonable assurance of safety and effectiveness has been established for the  
13 proposed indication for use based on the totality of the evidence provided herein.  
14 Additionally, the Panel will be asked to comment and make recommendations on whether  
15 the evidence provided adequately defines the patient population proposed to support the  
16 IFU.

17 Here is an outline of FDA's presentation. I will give a short introduction on the  
18 regulatory history and regulatory considerations, a brief description of the device, followed  
19 by an overview of the device principle of operation, as well as the nonclinical evaluations  
20 conducted in support of the submission. Dr. Annabelle Crusan will discuss the animal  
21 studies performed by the Sponsor, while Dr. Rona Tang will discuss the statistical  
22 information regarding the study design and statistical analysis plan. Then Dr. Tara Ryan will  
23 provide a presentation of the clinical study results from the primary COSIRA clinical trial  
24 performed by the Sponsor. Finally, I will provide a brief discussion of the proposed post-  
25 approval study along with a summary of today's presentation.

Free State Reporting, Inc.  
1378 Cape Saint Claire Road  
Annapolis, MD 21409  
(410) 974-0947

1 I would now like to briefly talk about FDA's breakthrough device designation  
2 program. This program is intended to provide the American public with timely access to  
3 new devices with the potential for significant impact by speeding up the development,  
4 assessment, and review while preserving the statutory standards. One of the key benefits  
5 of this program is that it offers increased interaction with FDA through several different  
6 mechanisms to help provide timely feedback in an effort to speed up the developmental  
7 process.

8 The Reducer device received a breakthrough designation from FDA, as it has the  
9 potential to provide benefit to a patient population that currently has limited treatment  
10 options.

11 It is important for the Panel to note that while the Breakthrough Devices Program  
12 offers increased communication and collaboration with FDA, it does not modify or reduce  
13 the statutory requirements for device approval. FDA still requires that the Reducer device  
14 demonstrate a reasonable assurance of safety and effectiveness.

15 In August 2019, FDA published a guidance document that discussed considerations  
16 of uncertainty in making benefit-risk determinations and premarket requests for medical  
17 devices. In this guidance document FDA indicates that it is willing to accept greater  
18 uncertainty of the benefit-risk profile if this uncertainty can be balanced by other factors.  
19 These other factors include the probable benefit to patients from earlier access to the  
20 device, timely postmarket data collection, or transparency and accountability in the device  
21 labeling.

22 Again, it should be noted that while FDA is willing to accept greater uncertainty in  
23 the benefit-risk profile, its statutory requirements of reasonable assurance of safety and  
24 effectiveness still hold. Additionally, FDA considers uncertainty of the benefit-risk profile  
25 based on the novelty and unique characteristics of each device. While breakthrough

1 devices have been acknowledged to have the potential to offer increased benefit, the  
2 breakthrough designation is not required to be considered under this framework.

3           During today's presentation we will highlight some of the key considerations and  
4 limitations FDA believed are present in the information provided by the Sponsor. These  
5 considerations and limitations include the preclinical animal study data's ability to  
6 demonstrate a trans-reducer pressure gradient or coronary sinus narrowing, clinical  
7 uncertainty of the primary study endpoint results for effectiveness due to study design,  
8 execution, and analysis issues. Interpretation of the secondary endpoints present  
9 challenges due to the limited sample size, high levels of missing data, and lack of pre-  
10 specified hypotheses for these endpoints. FDA is concerned about the uncertainty in the  
11 benefit-risk profile for a chronic cardiovascular implant in a population with a low mortality  
12 rate. Breakthrough devices, while offering potential benefit to patients, are still required to  
13 meet the statutory standards of a reasonable assurance of safety and effectiveness. And  
14 finally, the proposed post-approval study attempts to clarify major questions about the  
15 current clinical data; however, execution may be difficult.

16           Now, I would like to briefly discuss the device description, proposed principle of  
17 operation, and regulatory history.

18           The Neovasc Reducer is a metallic stent that is provided pre-mounted on the  
19 delivery balloon. The Reducer is designed to establish a narrowing of the coronary sinus  
20 with the goal of creating an upstream pressure gradient capable of improving perfusion to  
21 the ischemic myocardium.

22           While the COSIRA data was collected with the current device generation, the earlier  
23 first-in-human trial was performed using the prior device model. The difference between  
24 these generations were changes to the delivery system catheter, changes to the implant  
25 design, and changes to supply the device pre-mounted on the delivery balloon and no

1 longer requiring hand-crimping of the device.

2 As discussed, the Reducer device is intended to treat refractory angina. It is  
3 proposed that after implantation, local fluid disruption and vascular reaction will induce a  
4 neointimal proliferative response that will occlude the fenestrations in the metal mesh.  
5 However, while demonstration is not required as part of the PMA approval process, FDA  
6 believes that this proposed mechanism of action remains unclear and contributes to the  
7 overall uncertainty of the data provided. Additional information regarding this uncertainty  
8 will be provided later in the presentation by Dr. Crusan and Dr. Ryan.

9 Regarding the regulatory history of the Reducer device, an original IDE was  
10 submitted in 2010. FDA and the Sponsor could not agree on certain aspects of the study  
11 design and the Sponsor chose to conduct the COSIRA trial outside the U.S. In 2016, the  
12 Sponsor returned to FDA and sought approval for a new IDE seeking to build on the  
13 knowledge obtained during the COSIRA trial. This new trial, referred to here as COSIRA-II,  
14 was again proposed to be a double-blind, sham-controlled trial seeking to enroll refractory  
15 angina patients. COSIRA-II included a primary endpoint with an objective assessment of  
16 ischemic change using a Bruce treadmill exercise test, a larger sample size, powered  
17 secondary endpoints, and included a longer assessment duration with primary endpoints  
18 evaluated at 1 year with consent for follow-up going out to a total of 5 years.

19 However, during the review of the information provided in the IDE application,  
20 which included the COSIRA OUS data, FDA noted a few concerns that resulted in the FDA  
21 conveying to the Sponsor that additional data beyond what was proposed in COSIRA-II may  
22 be required to support a future PMA. FDA continues to believe that due to the novelty of  
23 this device and the current data limitations, a more direct observation of the device's  
24 performance is necessary. Ultimately, the COSIRA-II IDE was approved by FDA and remains  
25 open; however, the Sponsor has chosen to forego initiation of the trial, believing that the

1 current COSIRA data is sufficient to support their PMA.

2           Regarding the nonclinical testing provided as part of the submission, the Sponsor has  
3 provided several nonclinical bench tests in support of their device. The nonclinical  
4 engineering testing, biocompatibility, sterilization, shelf life, and packaging tests were all  
5 found to be acceptable.

6           I would now like to turn the presentation over to Dr. Annabelle Crusan to discuss the  
7 animal study information provided in this submission.

8           DR. CRUSAN: I'm Annabelle Crusan, a veterinary medical officer in the Office of  
9 Cardiovascular Devices. I am a laboratory animal veterinarian with training and experience  
10 in laboratory animal medicine, including cardiovascular device implantation and in  
11 laboratory animal pathology. I will be presenting FDA's review of the Neovasc Reducer non-  
12 GLP porcine studies.

13           Typically, a fundamental objective of an in vivo animal study in device evaluation are  
14 to provide insights into device proof of principle in support of a scientific rationale for  
15 treatment and to provide information on the basic safety of devices and the procedures.  
16 Large animal models are often utilized in cardiovascular studies because their anatomy and  
17 physiology can provide pathophysiologic insights into human device use. The animal  
18 studies that will be discussed today utilized a pig as the large animal model.

19           Good laboratory practices for animal care and study conduct, as specified in 21  
20 C.F.R. Part 68, assures the quality and integrity of animal study data. GLP animal studies are  
21 typically submitted to support premarket approval applications to FDA. The Sponsor  
22 submitted three non-GLP studies to FDA.

23           In Study CA04001 and R01, their first generation Reducer device was implanted in  
24 the coronary sinus of miniature pigs. Study REDACL003 used farm pigs and mini-pigs and  
25 the finalized model of the Reducer device that was implanted in the human COSIRA trial

1 was used. These animal studies were intended to evaluate the safety and efficacy of  
2 coronary sinus narrowing with the Reducer when implanted in the model that simulated  
3 human use.

4 FDA has concerns with the overall quality and integrity of the Reducer animal study  
5 data. Final reports were incomplete and there was a 6-year delay in the final pathology  
6 report submission. There were missing data on findings in some of the study animals. We  
7 were unable to verify that the protocols were followed due to the missing protocol  
8 information and the lack of amendment and deviation records. Minimum quality assurance  
9 requirements to ensure data validity and completeness were not verified due to the lack of  
10 independent audits for all three animal studies.

11 For a postulated principle of operation, the Reducer is intended to create a  
12 functional coronary sinus stenosis by a process that occurs when metallic stents are  
13 implanted in blood vessels. As depicted in the middle figure, intimal smooth muscle cells  
14 migrate and proliferate in areas in which the device is opposed to the vessel wall. As shown  
15 in the bottom figure, the Sponsor postulates that near-intimal growth covers enough of the  
16 Reducer to restrict coronary sinus blood flow to the stenotic central orifice of the device.  
17 Consequently, if the Reducer performs as intended, a pressure gradient is produced across  
18 the device resulting in elevated coronary sinus pressure.

19 It's important to note that the Sponsor performed a limited histological evaluation of  
20 coronary sinus stenosis and Reducer neointimal coverage. Specifically, for each Reducer  
21 device, only three histological transverse sections from the proximal, middle, and distal  
22 portions of the device were evaluated. Coronary sinus stenosis was assessed by light  
23 microscopy, qualitatively guided visually or by a quantitative morphometry in selected  
24 animals. Quantitative morphometry raw data were not entirely provided for FDA review. In  
25 addition, scanning electron microscopy is required to fully assess tissue coverage of



1 intravascular devices. It is notable that scanning electron microscopy was not performed in  
2 the animal studies.

3 An additional important consideration of coronary sinus stenosis and Reducer  
4 neointimal coverage is the limited histology reports provided from pigs with a device  
5 implanted for 2 to 6 months, when healing would be expected to be complete. In the  
6 finalized pathology reports there were a total of 12 animals implanted with a first  
7 generation device and seven animals implanted with the finalized device. Not all animals  
8 were accounted for in the pathology reports and photos were not provided for all animals,  
9 so FDA could not independently review all reported gross and histologic findings. Of note,  
10 photomicrographs were provided only for the chronic animals implanted with the final  
11 device.

12 In the evaluation of coronary sinus stenosis and neointimal growth, limited  
13 assessments showed that the Reducer did not create a significant vessel narrowing; 17 of 19  
14 or 89% of specimens demonstrated a less than or equal to 50% coronary sinus stenosis and  
15 no vessel had a more than 75% stenosis. Qualitatively, the study pathology noted there is  
16 an incomplete endothelial cell coverage of the device, variable to low levels of neointimal  
17 proliferation, and no appreciable coronary sinus luminal stenosis.

18 It is well established that there are limits to smooth muscle cell migration and  
19 proliferation in response to vascular injuries. Here are representative photomicrographs  
20 submitted to FDA from three pigs representing 57 days, 104 days, and 140-day time points  
21 from animal studies that used the finalized version of the Reducer. Note that there is very  
22 limited or very mild to at most moderate in-device neointimal thickening in the proximal  
23 and distal sections, and there is minimal eccentric neointimal or no neointimal in the  
24 midsection in which the Reducer is not opposed to the coronary sinus wall. Most  
25 importantly, there is no significant coronary sinus stenosis along the entire length of the

1 device.

2 In the study CA04001, the Sponsor reported a coronary sinus trans-reducer pressure  
3 gradient measured 15 minutes post-implantation of  $3.71 \pm 1.75$  mm/Hg. The coronary sinus  
4 pressure gradient across the Reducer at 2 to 6 months remained low and similar to  
5 baseline; that is, there was no evident coronary sinus stenosis. Of note, information on  
6 number of animals studied with evaluable data and information on pressure wire  
7 measurement methods were not provided to FDA.

8 In the CA04001 study using the first generation Reducer, there were eight pigs that  
9 developed myocardial ischemia after placement of an ameroid constrictor on the left  
10 complex artery. Four pigs underwent Reducer implantation in the coronary sinus and four  
11 ischemic pigs did not receive Reducer and served as controls. The Sponsor reported  
12 improved left ventricular contractility by dobutamine stress echo and myocardial perfusion  
13 by myocardial contrast echo in pigs treated with the Reducer at 6 weeks and at 6 months  
14 assessment. However, FDA could not confirm improvements in the left ventricular  
15 contractility and myocardial perfusion in Reducer animals due to the lack of data for review  
16 and there is only one control animal that was available for assessment at 6 months.

17 So in summary, in contrast to the postulated device mechanism of action shown in  
18 the upper figure, neointimal proliferation, as depicted in the bottom figure, was not  
19 sufficient to cover the Reducer and restrict coronary sinus blood flow to the central orifice  
20 of the device. Specifically, there was no confirmation of sufficient device coverage by  
21 neointima to restrict coronary sinus blood flow to the narrowed device central orifice, no  
22 confirmation of coronary sinus stenosis or elevated coronary sinus pressure, and limited  
23 evidence of improved myocardial contractility and blood flow.

24 I will now turn FDA's presentation over to Dr. Rona Tang, who will discuss the  
25 statistical aspects of the COSIRA trial.

1 DR. TANG: Good morning, my name is Dr. Rona Tang, FDA statistician. My  
2 presentation focuses on the study design, statistical method, and analysis results of the  
3 COSIRA study.

4 The COSIRA study is a prospective, multicenter, randomized, double-blinded, and  
5 sham-controlled study. It was conducted at 11 OUS sites. Subjects were randomized to two  
6 study arms. The treatment arm received the Reducer device while the control had a sham  
7 procedure. The originally planned sample size was 124 patients with 62 patients per arm.

8 The primary effectiveness endpoint of the COSIRA study is change of two or more  
9 CCS angina grades from the baseline to 6 months. There are several safety endpoints, but  
10 there's no primary safety endpoints nor any statistical hypotheses.

11 The primary effectiveness endpoint was evaluated in a hypothesis where the  
12 proportion of patients with improvement of two or more CCS grades was compared  
13 between the two arms.

14 The pre-specified test was a Pearson chi-squared test with continuity correction at a  
15 two-side alpha level of 5%.

16 One interim analysis was planned after 50% of the cohort completed their 6-month  
17 follow-up. The Lan-DeMets method using O'Brien-Fleming function was used to account for  
18 this interim look. The alpha level allocated to the interim analysis was 0.0031 and the alpha  
19 allocated to the final analysis was 0.0469.

20 Secondary effectiveness endpoints for the COSIRA study include change of one or  
21 more CCS grades, dobutamine echo wall motion score index, Seattle Angina Questionnaire  
22 score, and exercise tolerance testing. There's no pre-specified statistical hypotheses.

23 The study populations of the COSIRA study include the following: the intent-to-treat  
24 population has all randomized patients. The per-protocol population contains all  
25 randomized patients except patients who did not complete the study or treatment group

1 patients who did not receive the Reducer device. The as-treated population is all  
2 randomized patients; however, patients were analyzed based on the actual treatment  
3 received, regardless of the treatment assignment.

4 The enrollment period for the COSIRA study was 2010 to 2013. The originally  
5 planned sample size was 124 and the patient enrollment was stopped after 104 patients.  
6 The ITT population has exactly 52 patients per arm. Two patients randomized to the  
7 treatment arm did not have a Reducer implanted and they were removed from the  
8 treatment group in the per-protocol analysis and these two patients were counted as  
9 control in the as-treated analysis.

10 An interim analysis was performed for the primary outcome after 50% of the  
11 originally planned cohort completed 6-month follow-up. The Sponsor decided to continue  
12 the study after the interim analysis.

13 For the final analysis in the ITT population, the Reducer group has 18 out of 52  
14 subjects showing two or more CCS improvement, while the control group has eight  
15 successes. A p-value of 0.024 was obtained using Pearson chi-squared test without  
16 continuity correction. The per-protocol and as-treated analyses had similar results.

17 For the secondary endpoints, because there is no pre-specified statistical hypothesis,  
18 the descriptive results will be presented and discussed by FDA clinician Dr. Ryan in her  
19 presentation. I will discuss the missing data observed in the ITT cohort for the secondary  
20 endpoints.

21 For the DSE and SAQ endpoints, the control group has higher missing rates  
22 comparing to the treatment group as shown in the graph.

23 For the ETT endpoints, the missing rates in the treated and the control groups are  
24 both high. There are more than 71% missing data for endpoints timed to 1 mm and the  
25 maximum ST depression.

1 For the secondary data analysis, three missing data imputation methods were used,  
2 including the last observation carried forward, multiple imputation, and a tipping point  
3 analysis. Both LOCF and the multiple imputation methods require specific assumptions  
4 about missing pattern. In particular, the LOCF assumes patients are stable with their angina  
5 status during the follow-up period, which may not be true in this situation. As for the  
6 tipping point method, this is usually used when there are pre-specified hypotheses. Overall,  
7 due to high percentage and imbalance of missing data, the robustness of all secondary  
8 endpoint analyses are open to question.

9 In summary, the primary effectiveness endpoint was met for the COSIRA study.  
10 However, this study does not have a primary safety endpoint and there is no statistical  
11 hypothesis for the secondary endpoints. With a large number of secondary endpoints and  
12 no multiplicity adjustment, there is a real potential for false positive findings. Therefore,  
13 statistical inference such as p-value should be interpreted with caution. Additionally, there  
14 is a high percentage of missing data and the potential missing pattern issue in the  
15 secondary endpoints. The interpretation of the secondary endpoint analysis is challenging.

16 This concludes my presentation and I will now hand over to Dr. Ryan.

17 DR. RYAN: Good morning, my name is Dr. Tara Ryan. I am a medical officer in the  
18 Office of Cardiovascular Devices at FDA, and I will be presenting the clinical data on the  
19 Neovasc Reducer device.

20 The proposed indications for use for the Neovasc Reducer device is, "The Reducer  
21 System is intended for patients suffering from refractory angina despite guideline-directed  
22 medical therapy, who are unsuitable for revascularization by coronary artery bypass  
23 grafting or by percutaneous coronary intervention."

24 Refractory angina is a debilitating chronic condition in patients with coronary artery  
25 disease. It's characterized by severe unremitting cardiac pain, significantly limiting the

1 patient's daily activities. The pain in these patients cannot be controlled by optimal drug  
2 therapy, and these patients are not candidates for percutaneous coronary intervention or  
3 bypass surgery.

4 The number of patients with true refractory angina is difficult to quantify, but it's  
5 thought that the incidence and prevalence of refractory angina will continue to rise as  
6 coronary artery disease survival rates improve and the populations age.

7 So this morning I will discuss the clinical data provided to support the PMA  
8 application. The primary dataset in support of this PMA is from the COSIRA study. This is  
9 the study that the Sponsor believes confirms the safety and effectiveness of the Reducer  
10 device. Also provided in the FDA Executive Summary is supportive clinical data from the  
11 first-in-man pilot study and the REDUCER-I observational study. My presentation will focus  
12 on the results from the COSIRA trial.

13 The COSIRA study was a multicenter, randomized, double-blind, sham-controlled  
14 trial conducted at 11 sites outside of the U.S. Of note, this study was not an IDE study  
15 approved by FDA. The results from the COSIRA trial were published in 2015 in the *New*  
16 *England Journal of Medicine*, and a copy of that article has been provided in the panel pack.

17 The objective of the COSIRA trial was to evaluate whether the implantation of the  
18 coronary sinus-reducing device could effectively improve angina symptoms in patients with  
19 obstructive coronary artery disease who had concomitant evidence of reversible myocardial  
20 ischemia and who were not considered to be candidates for revascularization.

21 The data collected in the study was managed and analyzed by a contract research  
22 organization. The trial was overseen by an independent steering committee, a clinical  
23 events committee, and a data and safety monitoring board. Also incorporated in the study  
24 were core labs for the testing outlined in this slide.

25 Key inclusion criteria included adult patients with Canadian Cardiovascular Society

Free State Reporting, Inc.  
1378 Cape Saint Claire Road  
Annapolis, MD 21409  
(410) 974-0947

1 Class III or Class IV angina, despite efforts to control symptoms with medical therapy for a  
2 minimum of 30 days before screening. Medical therapy included beta blockers, calcium  
3 channel blockers, nicroandil, ivabradine, and short-acting and long-acting nitrates used at  
4 maximum tolerated doses.

5 With reference made to the Canadian Cardiovascular Society/Canadian Pain Society  
6 Joint Guidelines, all patients were required to have evidence of reversible myocardial  
7 ischemia and a left ventricular ejection fraction of more than 25%. Only patients who were  
8 not considered to be candidates for coronary revascularization were eligible to participate  
9 in the study. This eligibility was decided by the heart team at each institution.

10 Only patients with coronary sinus anatomy that was suitable for implantation of the  
11 device were eligible to undergo randomization and the relevant anatomical features are  
12 summarized here.

13 Patients were not eligible for participation if they had had a successful PCI or CABG  
14 procedure in the previous 6 months, patients with unstable angina, decompensated  
15 congestive heart failure, severe COPD, or arrhythmias requiring placement of a pacemaker  
16 or internal defibrillator.

17 Other exclusion patients were those who had a pacemaker or defibrillator electrode  
18 in the right atrium, right ventricle, or coronary sinus, patients with chronic renal failure, or  
19 patients who had undergone tricuspid valve surgery.

20 This figure outlines the flow of patients. A hundred and sixty-six patients were  
21 consented and ultimately a hundred and four patients were enrolled in the study. Of note,  
22 the original COSIRA protocol planned to randomize a hundred and twenty-four patients in  
23 order to generate an evaluable cohort of a hundred and twelve patients at 6 months post-  
24 procedure. In 2013, enrollment was stopped. This decision was based on difficulty  
25 enrolling patients as well as a less-than-expected loss to follow-up. It was not safety

1 related.

2 This figure outlines the flow of patients after randomization took place. A total of 52  
3 patients were in the treatment group and 52 patients were in the control group. One  
4 patient died before reaching the 6-month follow-up period leaving a hundred and three  
5 patients available for follow-up.

6 Enrollment in the study began in September of 2010. A total of a hundred and four  
7 patients were enrolled in the study. Enrollment stopped in May of 2013 and 6-month  
8 follow-up was completed on all patients by November of 2013.

9 So very briefly, patient demographics show that 85% of the patients were male and  
10 87% were Caucasian. The average age of patients enrolled was 68. The COSIRA study was  
11 not representative for the minority or female population as evidenced by the gender and  
12 ethnic information summarized on this slide.

13 The Reducer group had numerically fewer diabetics, 40.4% versus 48.1%, but the  
14 two groups were comparable with regard to smoking history.

15 The two groups were also comparable with respect to hypertension, family history of  
16 cardiovascular disease, hypercholesteremia, hypertension, and other relevant medical  
17 history. Of the total 104 subjects, 73.1% had a history of a previous PCI, 76.9% had a  
18 history of a prior CABG, and 54.8% had a history of a previous myocardial infarction.

19 Dobutamine echo testing was done in 50 patients in each group. Eighty-eight  
20 percent of Reducer patients were positive for reversible ischemia and 58% had wall motion  
21 abnormalities at rest. The ejection fraction was 53.5% in the Reducer group and 54.8% in  
22 the control group.

23 SPECT testing was done at baseline in 48 Reducer patients and 45 control patients.  
24 This testing was positive for reversible ischemia in 86% of Reducer patients and 77% of  
25 control patients. Irreversible defects were also noted in a percentage of each group.



1 As noted previously, a CCS score of III or greater was required for patients to be  
2 included in the COSIRA study. Eighty-one percent of Reducer patients were CCS Class III and  
3 19% were Class IV. With respect to control patients, 87% of patients were Class III and 13%  
4 were Class IV.

5 The table presented in this slide shows the number of antianginal medications that  
6 patients in each group were on at the time of enrollment. Thirty-eight percent of patients  
7 in the Reducer group were on three or more antianginal medications and 31% of control  
8 patients were also on three or more antianginal medications. Twenty-seven percent of  
9 Reducer patients were on less than or equal to one antianginal medication and 25% of  
10 control patients were on less than or equal to one antianginal medication.

11 The table presented in this slide shows the baseline cardiovascular medications that  
12 the COSIRA patients were on at the time of enrollment. Of note, 77% of patients in each  
13 group were on beta blockers, 55% were on calcium channel blockers. For a pivotal IDE trial  
14 of a device intended for the refractory angina population, FDA would typically request data  
15 to assure that patients were on maximally tolerated, guideline-directed antianginal  
16 medications, including drug dosages used and the rationale if patients were not on maximal  
17 therapy.

18 The primary endpoint for the COSIRA study was an improvement of two or more CCS  
19 grades from baseline to 6 months post-procedural evaluation in the Reducer group. Of  
20 note, this was an efficacy endpoint. There was not a primary safety endpoint specified, but  
21 periprocedural serious adverse events and major adverse events were recorded.

22 The study was powered to detect a difference in the proportion of patients  
23 improving greater than or equal to two CCS angina classes at 6 months after the  
24 implantation. The hypothesis corresponding to the primary objective as well as  
25 assumptions are summarized on this slide. The statistical analysis of the primary endpoint

1 used the intention-to-treat population. My FDA colleague, statistician Dr. Rona Tang, has  
2 provided additional details and the FDA perspective regarding the statistical methodology  
3 and analyses that were completed, including an interim analysis that was performed on the  
4 primary effectiveness endpoint.

5 So, for the primary endpoint in the treatment Reducer group, 18 out of 52 patients,  
6 that is 34.6%, experienced a decrease in two or more CCS grades from baseline to 6 months  
7 versus 8 out of 52 patients or 15.4% in the control group. This corresponds to a p-value of  
8 0.024. The primary endpoint was therefore met. Of note, while 34.6% of the Reducer  
9 group demonstrated success using the protocol definition, 65.4% did not. Also of note, the  
10 response/success rate in the control group was almost half that noted in the treatment  
11 group, highlighting the potential for a placebo effect, which I will discuss later in my  
12 presentation.

13 Secondary endpoint success in the COSIRA trial included a CCS classification decrease  
14 of one or more CCS grades, dobutamine echo wall motion score index, Seattle Angina  
15 Questionnaire scores, and exercise tolerance testing which included total exercise duration,  
16 time to ST-segment depression, and maximal ST-segment depression. It is important to  
17 note that there were no pre-specified hypotheses for any of these secondary endpoints.

18 So, for the secondary endpoint of improvement of one or more CCS grades from  
19 baseline to 6 months, the proportion of subjects experiencing an improvement of one or  
20 more CCS classes from baseline to 6 months in the intention-to-treat population was 71.2%  
21 in the Reducer group compared to 42.3% in the control group.

22 Based on the clinical data in the submission, FDA performed an analysis of the  
23 individual changes in CCS grades and this table provides an overview of the analysis. This  
24 table illustrates the number of patients who did not show a response at all or got worse.  
25 More specifically 13 out of 52, or 25%, of Reducer patients showed no change from baseline

1 and 24 out of 51, or 47.1%, of control patients showed no change. With respect to an  
2 increase in CCS grade at 6 months, 2 out of 52, or 3.8%, of Reducer patients got worse,  
3 while 5 out of 51, or 9.8%, of control patients got worse.

4 Dobutamine echo wall motion index was calculated using a total of 16 segments.  
5 The scores for each segment were summed and the total was divided by the number of  
6 segments analyzed. Additionally, since the Reducer is placed in the coronary sinus distal to  
7 the right coronary artery, venous drainage, and modified left coronary artery, WMSI was  
8 calculated using only the 11 segments attributed to the LCA system. The WMSI and the  
9 modified LCA WMSI were calculated both at rest and during stress.

10 Testing was conducted using an acquisition protocol established by the core lab.  
11 Changes from baseline to 6 months in the wall motion indices are outlined in this slide.  
12 Given that there were no pre-specified hypotheses, any inferences must be drawn with  
13 caution. Also of note is a significant amount of missing data, which I will touch upon again  
14 later.

15 The Seattle Angina Questionnaire quantifies the physical and emotional effects of  
16 coronary artery disease. The questionnaire is a 19-item self-administered tool and provides  
17 results in five scales that measure clinically important dimensions of coronary artery  
18 disease, physical limitation, anginal stability, anginal frequency, treatment satisfaction, and  
19 disease perception/quality of life.

20 The current slide shows the data for the intention-to-treat analysis using paired data  
21 for physical limitations, anginal stability, and anginal frequency. Results for both anginal  
22 stability and anginal frequency favor the Reducer device. As mentioned previously, there  
23 was no hypothesis specified and any inferences should be drawn with caution.

24 And this slide shows summary data for the treatment satisfaction and quality of life  
25 component of the Seattle Angina Questionnaire. Once again, this is ITT paired data. A

1 change in the score for quality of life from baseline to 6 months was 18.6 in the Reducer  
2 group and 8.7 in the control group, favoring the Reducer group in improving quality of life  
3 scores.

4 A bicycle ergometry stress test was performed by the ETT core lab using standard  
5 operating procedures. Multiple parameters were intended to be recorded at baseline and  
6 6-month follow-up for a comparative analysis. This included total exercise duration, time to  
7 1 mm ST-segment depression, and maximal ST-segment depression.

8 The table on the current slide shows that the paired data analysis trended in favor of  
9 the Reducer group with the total exercise duration mean increasing by 65 seconds in the  
10 Reducer group versus a mean increase in the control group of 4.31 seconds.

11 However, it is important that we look at the issue of missing data. While all subjects  
12 were expected to have dobutamine stress echo testing, Seattle Angina Questionnaire data,  
13 and exercise tolerance testing as part of their baseline and 6-month follow-up, some of  
14 these analyses had a substantial amount of missing data, particularly related to the ETT  
15 analysis. This table shows an overview of the missing information. It's noted that there's  
16 more missing data in the control group as compared to the treatment group in the  
17 dobutamine stress echocardiography data and the Seattle Angina Questionnaire data. The  
18 missing data is problematic with respect to being able to draw conclusions regarding  
19 efficacy of the Reducer device in the improvement of functional ischemia.

20 There were also measures that were considered to be observational in the COSIRA  
21 study, including SPECT segmental analyses and CT angiographic analysis. The current slide  
22 shows differences in SPECT data. There were no significant trends noted.

23 With respect to CT angiographic follow-up, 27 patients in the Reducer group  
24 underwent angiographic follow-up. There were no signs of migration of the device.  
25 Contrast flow was noted in 94% of patients who underwent imaging. Thrombus was

1 reported in 43.2% of patients.

2 A total of 62 devices were prepared for use in the COSIRA study. Fifty-two were  
3 used in patients randomized to the treatment group. Ten were discarded due to what the  
4 Sponsor terms a device "deficiency" or a device malfunction.

5 For the sake of consistency with respect to characterization of device performance,  
6 this slide summarizes the events as device malfunctions. Eight of these were listed as  
7 operator mishandling events. In one case, the operator had to use a second device due to  
8 clots that were noted in the guide catheter and on the wire. In one patient, the Reducer  
9 device slipped on the balloon while advancing. The physician did not notice this until the  
10 device was inflated and deployed without essential narrowing. A peripheral balloon was  
11 inserted and expanded within the mal-deployed Reducer.

12 Technical success was assessed in the Reducer group only and defined as successful  
13 delivery and deployment of the Reducer device to the intended site as assessed by the  
14 investigator. In the 52 patients who were randomized to the Reducer group, 50 out of 52,  
15 or 96.2%, had a Reducer successfully implanted.

16 Periprocedural serious adverse events were assessed in both the Reducer and  
17 control groups. In the control group, periprocedural serious adverse events were defined as  
18 a composite of death; MI; life-threatening arrhythmias, that is ventricular tachycardia or  
19 ventricular fibrillation; and respiratory failure through 30 days post-procedure as  
20 adjudicated by the CEC.

21 In the Reducer group, periprocedural serious adverse events were defined as the  
22 same, with the addition of cardiac tamponade and clinically driven re-dilation of a failed  
23 Reducer as qualifying events. In the current slide, we see that there was one documented  
24 case of a non-STEMI. Of note, cardiac enzyme data was not captured in all patients post-  
25 procedure.

1 Major adverse events were defined as a composite of cardiac death, major stroke,  
2 and MI in the Reducer and control groups through hospital discharge and at 30 days, 3  
3 months, and 6 months post-procedure evaluation. There were a total of five major adverse  
4 events as adjudicated by the CEC. There was one myocardial infarction in the Reducer  
5 group and three myocardial infarctions and a cardiac death in the control group. None of  
6 the five events occurring after 30 days post-procedure were attributed to the procedure or  
7 the investigational device. Only one MI was considered by the CEC to be related to a study  
8 specific assessment as it occurred during the study-required dobutamine stress echo at the  
9 6-month follow-up.

10 As noted several times in my discussion, the primary dataset that has been  
11 presented as supportive of safety and effectiveness is the COSIRA study. Additional clinical  
12 data has been provided in the PMA application, including a first-in-man study using a first  
13 generation Reducer device at three sites in India and Germany. A total of 15 patients were  
14 treated.

15 Also ongoing is an observational study called the REDUCER-I and this slide presents  
16 summary information on this study. First, I would like to emphasize that the REDUCER-I is  
17 an observational study. There is no control group and therefore FDA is unable to interpret  
18 the data and it is difficult to draw conclusions with respect to effectiveness.

19 The study is currently continuing to enroll patients. The goal is to enroll 400  
20 patients. There are three patient groups in the REDUCER-I study. FDA has received interim  
21 reports, the most recent having been provided in May of 2020. As of March 12th, 2020,  
22 there had been 241 subjects enrolled in the study.

23 With respect to adverse events in REDUCER-I, not all events had been adjudicated by  
24 the CEC as of the last data lock in March of 2020. There had been 13 deaths of study  
25 subjects reported. Ten events have been adjudicated by the CEC as unrelated to the device

1 and/or the procedure. Three events are pending adjudication. There's been one reported  
2 case of cardiac tamponade. An MI rate of 6.6% has been reported. There have been four  
3 stroke events, corresponding to a rate of 1.7% reported in the REDUCER-I study.

4 So to conclude my presentation, I will discuss limitations of the COSIRA study with  
5 respect to consideration as a pivotal study to demonstrate safety and effectiveness, and  
6 also comment on the adjunctive clinical data presented. First, we would like to relate to  
7 Panel members that study limitations have been discussed with the Sponsor, both in writing  
8 and during multiple in-person and teleconference interactions. FDA worked with the  
9 Sponsor to develop a protocol called COSIRA-II, a study protocol that was approved under  
10 an IDE application in November of 2017.

11 That being said, I would like to begin by commenting on the placebo effect. There's  
12 substantial clinical trial evidence for the placebo effect, especially for the relief of pain.  
13 Angina is known to be a placebo-responsive condition as evidenced in previous blinded  
14 studies in the field of refractory angina, including myocardial revascularization studies.  
15 We've also seen the placebo effects highlighted in other recent cardiac trials such as  
16 ORBITA.

17 The importance of exercise tolerance testing as a measure of functional capacity and  
18 predictable ischemic threshold has resulted in its use as a primary efficacy endpoint in  
19 clinical trials that have evaluated anti-ischemic treatments. And in appropriate control  
20 studies, it remains a valid endpoint with the ability to discriminate among effective and  
21 likely ineffective therapies. The COSIRA trial utilized a subjective endpoint of angina.

22 The relatively small sample size for the COSIRA study is problematic. Small sample  
23 size limits our ability to draw conclusions regarding the effectiveness of the coronary sinus  
24 Reducer device. A small positive trial does not negate of the possibility of negative results  
25 in a larger adequately powered study.

1 With respect to blinding in the sham-controlled trial, methods were incorporated  
2 during the procedure to help ensure that the patient and family remained blinded as to  
3 what group the patient was randomized to. However, there was no formal assessment or  
4 analysis regarding the success of blinding throughout the entire study.

5 As noted previously, the COSIRA study was not statistically powered to detect an  
6 improvement in ischemia by objective measures such as stress testing or wall motion index.  
7 Objective criteria of exercise time and stress wall motion could not be shown to improve  
8 significantly given the underpowered nature of the study.

9 As illustrated in my next slide, there was a significant amount of missing data for the  
10 secondary endpoints. Outlined previously in 32, we see there was a significant amount of  
11 missing data for the secondary endpoints. For the DSE and SAQ data, missing data were  
12 approximately twice as common in the control group as compared to the Reducer group.  
13 For the exercise tolerance testing, the data for time to 1 mm depression and maximal ST  
14 depression were mostly missing.

15 With respect to the statistical methodologies employed for data analysis, you heard  
16 from FDA statistician Dr. Tang regarding some of the limitations involving the analyses.  
17 Because optimal medical therapy was not clearly defined, it's difficult to know whether all  
18 patients were comparable with respect to being truly "no option," that is medication  
19 compliance, maximally tolerated doses of medications, and patient-specific conditions  
20 precluding revascularization.

21 Long-term safety data in the controlled COSIRA trial is limited. Seven patients from  
22 the first-in-man study had been reported at a site in India and 20 patients had follow-up  
23 data at 5 years in the REDUCER-I observational study. Given that the Reducer is a  
24 permanent implant, the long-term risks versus benefits cannot be determined. The long-  
25 term angiographic data collected does not assess the device stenosis and/or the resulting



1 increase in coronary sinus pressure.

2 Another limitation is the fact that the demographics of the study, who have been  
3 implanted with the Reducer device, is not representative of the U.S. coronary artery disease  
4 population. More specifically, no African American or Hispanic patients were enrolled in  
5 the COSIRA trial and the majority of patients enrolled were male.

6 Finally, the high non-responder rate noted in the COSIRA trial highlights the need to  
7 be able to identify patients who would likely benefit from receiving this permanent implant.  
8 Baldetti and others in *EuroIntervention* discuss that the presence of alternative venous  
9 drainage systems to the coronary sinus may be crucial in determining whether a patient is  
10 likely to benefit from the Reducer. Other tools may prove to be important. This is an item  
11 that we would like the Panel to consider.

12 So this concludes my presentation regarding the clinical data for the Neovasc  
13 Reducer device. Dr. Raben will continue now with comments regarding the Sponsor's  
14 proposed post-approval study and then concluding FDA remarks. Thank you for your  
15 attention.

16 DR. RABEN: Now I would like to briefly discuss the Sponsor's proposed post-  
17 approval study along with providing some concluding remarks for the presentation.

18 It is common for new and novel PMA devices to have a post-approval study  
19 identified as a condition of approval. These post-approval studies are conducted to provide  
20 additional information required to further support labeling for these new devices and can  
21 be used as a way to address remaining uncertainties identified during the review. Please  
22 note that the discussion of a post-approval study should not be interpreted to mean that  
23 FDA has made a decision or is making a recommendation on the approvability of this PMA.

24 For this device, the Sponsor is proposing a randomized, double-blind, sham-  
25 controlled trial to be performed in a country where the Reducer device has not been

1 approved. Given the concerns identified regarding the statistical and clinical data as  
2 discussed earlier in the presentation, we intend to ask the Panel if a post-approval study  
3 that is randomized and sham controlled is appropriate and feasible.

4 In conclusion, today's presentation has discussed the data and information provided  
5 by the Sponsor to support the Reducer system. We have discussed the regulatory history,  
6 the COSIRA trial along with COSIRA-II. We have discussed some of the regulatory  
7 considerations regarding the Breakthrough Devices Program and FDA's uncertainty  
8 guidance. FDA provided a brief overview of the nonclinical bench testing provided by the  
9 Sponsor.

10 Dr. Crusan provided the discussion of the animal study data collected in support of this  
11 device. Drs. Tang and Ryan provided a detailed discussion of the statistical and clinical  
12 results regarding the COSIRA trial that the Sponsor believes demonstrates sufficient  
13 performance of their device. Dr. Ryan also discussed the currently ongoing REDUCER-I  
14 observational trial and some of the information that has been provided from the study.  
15 And lastly, we discussed the Sponsor's proposed post-approval study.

16 Additionally, as part of FDA's presentation, we have also highlighted some of the key  
17 considerations and limitations FDA believes are present in this information. These  
18 considerations and limitations include significant limitations regarding the preclinical animal  
19 studies' data's ability to demonstrate a trans-reducer pressure gradient in coronary sinus  
20 narrowing; substantial clinical uncertainty of the primary study endpoint results or  
21 effectiveness due to study design execution and analysis issues.

22 Interpretation of the secondary endpoints present major challenges due to limited  
23 sample size, high levels of missing data, and a lack of the pre-specified hypotheses for these  
24 endpoints; uncertainty in benefit-risk profile for a chronic cardiovascular implant in a  
25 population with a low mortality rate. And while this device is a breakthrough device, it is

1 still required to meet the statutory requirements for a reasonable assurance of safety and  
2 effectiveness.

3 And finally, the post-approval study provided by the Sponsor attempts to clarify  
4 major questions about the current clinical data; however, execution of this trial may be  
5 difficult. We believe these issues should be considered as you begin your deliberations this  
6 afternoon.

7 This concludes FDA's presentation and we thank you for your attention and are  
8 happy to answer any questions you may have.

9 DR. LANGE: Great. I would like to thank the FDA speakers for their presentation and  
10 this portion of our Panel deliberation is for us to ask any clarifying questions. And again, I'll  
11 ask you to raise your hand and I'll try to notify. I think Dr. Cigarroa I saw first and if I'm not  
12 mistaken, then Dr. Somberg and Dr. Allen.

13 DR. CIGARROA: Thank you, Dr. Lange.

14 Two questions for FDA, one on Slide 22, point 3, with regards to the statement that  
15 FDA could not confirm improvements in LV contractility nor perfusion and I just want a  
16 clarification from FDA if the data was reviewed and they could not confirm it, if primary  
17 data was reviewed, I'm just a little confused on that point.

18 The second question is a statistical question and would like the perspective from  
19 FDA on the probability of Type 1 error given the number of patients enrolled at a hundred  
20 and four versus the planned statistical analysis at the outset of a hundred and twenty-four.  
21 Thank you.

22 DR. LANGE: Okay. And as we pose questions, there will be some questions we may  
23 ask the Sponsor to address because they've got the primary data and this may be one of  
24 those and so Sponsor, if you'll -- I do think the FDA needs to explain why data could not be  
25 confirmed, so I'll ask them to address that and then I'll ask the Sponsor to address the

1 probability of a Type 1 error in a hundred and four versus a hundred and twenty-four  
2 patients.

3 Dr. Somberg.

4 DR. SOMBERG: A couple of questions to the FDA. One is they talked about missing  
5 data, but the primary endpoint of the COSIRA study is the CCS data. Was there missing data  
6 in that? I didn't see that on the table they listed.

7 Two is, there's been a lot of talk of were these patients refractory angina. The FDA  
8 pointed out a table of -- with about 25% not receiving anginal therapy. Is it the FDA's  
9 conclusion that this was not a refractory angina patient population?

10 And my third question is about the COSIRA-II trial, which my question is, is there a  
11 procedure in that trial that would ensure that the patients are a "no option" resistant  
12 population for both drugs and intervention?

13 DR. LANGE: Okay, I'm going to -- so I'll ask Dr. Zuckerman to address any one of  
14 those three, whether missing CCS data -- I'm not sure the FDA is going to render an opinion  
15 about whether they believe there was refractory angina or not, other than just presenting  
16 the data, but -- and then COSIRA-II.

17 So Dr. Zuckerman.

18 I see your hand, Dr. Vetrovec, as well.

19 So Dr. Zuckerman.

20 DR. ZUCKERMAN: Yes, Dr. Lange, thank you for your points.

21 With regard to Dr. Somberg's question number two, I would ask that you direct the  
22 Sponsor to present their interpretation regarding that question after lunch. The data that  
23 FDA has received on medications is extremely difficult for us to decipher and that was the  
24 best that we could do with that slide, Dr. Somberg.

25 DR. LANGE: Great.

1 DR. SOMBERG: What about the missing data? Was there any missing data on the  
2 primary endpoint for the COSIRA --

3 DR. ZUCKERMAN: We'll be able to answer that, but the Sponsor also can quickly  
4 answer that question for sake of efficiency.

5 DR. LANGE: Great. I've got Dr. Allen, Dr. Vetovec, and Dr. Yeh and then  
6 Dr. Mathew. Verghese, pardon me.

7 So Keith, go.

8 DR. ALLEN: Yeah, hi. Keith Allen.

9 Can the FDA maybe comment on -- normally they would make suggestions on how to  
10 manage missing data and it may be since the secondary endpoints, which is where most of  
11 the missing data was, weren't pre-specified, the FDA just didn't go there. But can the FDA  
12 make any recommendations on how they could analyze the missing data in those key  
13 endpoints and shed light on different imputations, worst-case scenario, for example?

14 DR. LANGE: I see Dr. Zuckerman's hand. That may be a question we'll ask  
15 Dr. Connor to address. And Dr. Wittes.

16 DR. ZUCKERMAN: Yes. In general, missing data are extremely bad. I think that  
17 Drs. Connor and Wittes, in the afternoon, can talk about the missing data problem,  
18 especially when there's a sense that the missing data may not be missing at random as  
19 reflected in several of the datasets, such as the DSE and SAQ datasets.

20 DR. LANGE: Great. I've got Dr. Vetovec, Dr. Yuh, Dr. Mathew and then Dr. Starling.  
21 Thank you, Randy.

22 So Dr. Vetovec.

23 DR. VETROVEC: Yes, I have a quick question and Bram may be able to clear this up.  
24 It's pretty clear that the mechanism of action is not clear, but is that a requirement for  
25 approval? What is the definition of that?

1 DR. ZUCKERMAN: Yes. Underlined in the FDA presentation, a clear understanding of  
2 the mechanism of action is not required to conclude that a device is safe and effective.  
3 However, it may be helpful sometimes to know something about the mechanism of action  
4 when there is uncertainty regarding some of the other variables where there is uncertainty.  
5 But you are absolutely correct, Dr. Vetrovec.

6 DR. LANGE: Great. I've got Dr. Yuh. Go ahead.

7 DR. YUH: Yes, one question. FDA made a point of emphasizing the potential placebo  
8 effect in this patient population based -- was there any review of the literature to indicate  
9 the magnitude of the placebo effect we would expect from this type of a study? I mean,  
10 15% positive response in the control group, at least, seems to me qualitative to be high, but  
11 I just want to scale it or calibrate it to what FDA may have gleaned in terms of what would  
12 be the expected placebo rate.

13 DR. LANGE: And I'll hold that for after lunch if the FDA would like to address that  
14 and if not, either myself or other panelists may want to address that, as well.

15 I see Dr. Mathew, Dr. Starling. And then I think, Dr. Connor, I think you had your  
16 hand up, as well.

17 Okay, so Dr. Mathew.

18 DR. MATHEW: Thank you, Dr. Lange.

19 A couple questions. One is the post-approval study that's being discussed is a  
20 randomized trial, and my question is are we married to that particular version if this goes to  
21 that degree or are there alternate models, such as a well-defined pre-specified registry  
22 after approval, that might be considered in our discussion?

23 The second point had to do with missingness of data, but based on Dr. Zuckerman's  
24 comment, I presume we're going to have to further discussion on that this afternoon, so  
25 that's fine.

1 Thirdly is just really just a statement, that part of the inclusion criteria here is  
2 demonstration of ischemia by dobutamine stress echocardiogram, which may exclude some  
3 patients in real life from a pragmatic perspective for this sort of therapy, and then  
4 demonstrating improved ischemia similarly over the course of the trial may also exclude  
5 some patients that may symptomatically benefit. Again, that might be heretical when we're  
6 talking about approval here, but I think in the real world, since almost exclusively this -- the  
7 goal here is improvement of symptoms, not really ischemia, I think, in the patient  
8 population that we see, it's just something to consider as we talk about this later this  
9 afternoon. Thank you.

10 DR. LANGE: Thank you, Dr. Mathew. And with regard to a post-approval study, we'll  
11 have the opportunity to provide our opinions and recommendations to the FDA during our  
12 discussion.

13 DR. MATHEW: Great, thank you.

14 DR. LANGE: Yeah, it's not a fait accompli. Actually, we'll be interested in this.

15 Dr. Zuckerman from the FDA.

16 DR. ZUCKERMAN: Yeah, I'd really like to underline the point that Dr. Lange just  
17 made, we're married to no specific post-approval study design. The point is that the  
18 Sponsor believes that the uncertainty present in the current data can be adequately  
19 addressed with a post-approval randomized controlled trial called COSIRA-II. The Panel, in  
20 the afternoon, will be needed to address that particular point under Dr. Lange's direction.

21 DR. LANGE: Great. I've got Dr. Starling, Dr. Connor, Dr. Cigarroa, and Dr. Page.

22 So Dr. Starling.

23 DR. STARLING: Thank you, Dr. Lange. Randall Starling speaking.

24 During the FDA presentation, when Slide 45 was up, a comment was made regarding  
25 the heart team. My question is do we have any specific information on what the heart

1 team was composed of?

2 My second question for the FDA is, does the FDA have information that they have  
3 received from the Sponsor regarding anginal relief beyond the 6-month primary endpoint?

4 DR. LANGE: Okay. So I'll ask the Sponsor to address the heart team issue and that is  
5 who decided -- and we asked this before, who adjudicated that the individual was not  
6 eligible to be revascularized? And then I don't believe there's any other data, but Bram, if  
7 you want to address that.

8 DR. ZUCKERMAN: Well, I think the Sponsor did give us this morning some longer-  
9 term data, Dr. Starling, but the subscript on the slide indicated -- on most of those slides,  
10 FDA hadn't reviewed the data. So perhaps the Sponsor can again review those data with  
11 you, but begin with the proviso of what has and has not been reviewed by FDA for your  
12 determinations.

13 DR. LANGE: Apropos to the question, I think, is there longer-term data on COSIRA  
14 and other than those 11 patients that may be in the Reducer study, there are no additional  
15 data available. Okay.

16 Dr. Connor.

17 DR. CONNOR: Yeah, Jason Connor here.

18 My question or comment is kind of a follow-up to Dr. Yuh's when he asked for maybe  
19 a lit review or comment on placebo effect. So I would basically ask the same question but  
20 ask it about natural variability. I'm hesitant to call many things placebo effects. So really, I  
21 assume patients just live with the disease. When they're worse is when they go to the  
22 doctor and when they go to the doctor is when they get invited to a clinical trial. So just in  
23 the variability of their disease they only enter here, which means they can only get better,  
24 or natural variability means they're likely to get better than worse.

25 So what we're seeing in the sham responders is really just regression to their own



1 mean. So that's what I would like to understand, is there any history on these patients or  
2 other patients in the literature just so we can understand, on the scale for the primary  
3 outcome, how likely patients are to change by two or more steps over a 6-month period,  
4 keeping in mind that, by design, they enrolled patients when they were at a random high.

5 DR. LANGE: Okay. So again, we'll ask the Sponsor to address that, that is the natural  
6 variability of the regression to mean either among these patients or history of life (ph.)  
7 patients. Is that correct, Dr. Connor?

8 DR. CONNOR: That's right, thanks.

9 DR. LANGE: All right, very good.

10 Dr. Cigarroa, I think you had your hand up next.

11 DR. CIGARROA: Thank you, Dr. Lange. For the record, this is Joaquin Cigarroa.

12 A question with regards to potential safety concerns for FDA. Did FDA consider  
13 and/or were they concerned about the potential for an elevation in coronary sinus pressure  
14 potentially having a deleterious effect on diastolic function? Early data from the 1970s and  
15 early '80s suggested that as a possibility.

16 DR. LANGE: All right. And I'll throw that back to the Sponsor and that is, were they  
17 concerned about it, did they assess it in any way, can they provide any information about  
18 diastolic function with CS? Is that fair, Dr. Cigarroa?

19 DR. CIGARROA: Yes, sir.

20 DR. LANGE: Okay, great.

21 Dr. Page.

22 DR. PAGE: Thank you, Dr. Lange. For the record, Richard Page.

23 My question is actually to Dr. Zuckerman with regard to his comment about  
24 understanding the mechanism of action of devices. A number of us have been through a  
25 number of these panels where we see statistics reaching barely significance or just missing

1 significance and while we don't necessarily need to understand fully the mechanism of a  
2 device's action, aren't we obligated to consider scientific plausibility as we assess the data  
3 overall and come to a conclusion as to a reasonable assurance of effectiveness?

4 DR. ZUCKERMAN: I think you've summarized some key points well, Dr. Page, that  
5 would be worthwhile to further explain this afternoon, especially since there's been  
6 mention of the recent FDA uncertainty guidance and how we're directed to look at and  
7 make a final determination of safety and effectiveness. But the basic scientific and logical  
8 points you're pointing to are correct. We generally want to look at the totality of data such  
9 that we can, at the end of the day, conclude that there's a reasonable benefit-risk profile.

10 Now, while it is not absolutely required that a mechanism of action be shown  
11 definitively, it can be helpful in that determination. On the other hand, though, I think you  
12 can fully appreciate other cases whereby you might be looking at a thousand-patient  
13 cardiovascular trial where a substantial decrement in mortality is shown at a p-value less  
14 than 0.001 and the mechanism of action is unknown and still, because you believe the  
15 clinical trial was well executed, that you would conclude safety and effectiveness. So the  
16 good common sense in scientific principles that you just enumerated are the ones that the  
17 Panel will be utilizing this afternoon in a thorough discussion.

18 DR. LANGE: Great. I've got Dr. Vetovec, Dr. Ohman, and Dr. Gersh and then  
19 Dr. Allen again.

20 So Dr. Vetovec.

21 DR. VETROVEC: Mine was an accidental hand raise.

22 DR. LANGE: I'll wave back at you, George.

23 All right, Magnus.

24 DR. OHMAN: Magnus Ohman here.

25 I have a question, if the FDA has looked at the discrepancy between the change in

1 the Canadian classification and the SAQ frequency of angina questions. In other words, it  
2 looks to me that these two are not quite aligned and I would be interested in the FDA's  
3 position of why they might not be aligned.

4 DR. LANGE: I'm happy to allow them to talk about that over lunch and come back  
5 and so this is asking their opinion. Bram, are you comfortable with that?

6 DR. ZUCKERMAN: That's an excellent question the Sponsor may want to opine also.

7 DR. LANGE: We'll ask both parties to do so. Thank you. Thank you, Bram.

8 Dr. Gersh.

9 DR. GERSH: For the record, Bernard Gersh.

10 A quick comment and a question. The comment relates to what Dr. Connor said  
11 about regression to the mean. You know, one of the pivotal, I think, sham-controlled trials  
12 is SIMPLICITY -- and after the trial was published, we did an analysis with Dr. Stewart  
13 Pocock, the statistician, and showed unequivocally and published in *JACC* that there clearly  
14 was regression to the mean and it occurred in both groups. So I think it's -- I wouldn't be at  
15 all surprised if there's regression to the mean. At least in SIMPLICITY it was quite easy to  
16 show it statistically.

17 The other comment is I know the measurement -- this comment about contractility  
18 increased with dobutamine stress echo, I think we've been looking for measurements of  
19 contractility for about the last 30 to 40 years and no one has really come up with a clear  
20 measure, but one thing is clear to me, in someone with an ejection fraction that is normal,  
21 which it is in the vast majority of patients here, I don't see how you can assess contractility  
22 with a dobutamine stress echo in someone beginning with an EF of 55 to 60%. It doesn't  
23 make sense to me.

24 DR. LANGE: Thank you, Dr. Gersh.

25 I've got Dr. Allen and then Dr. Wittes and then Dr. Connor.

1 So Dr. Allen.

2 DR. ALLEN: Thank you. Keith Allen.

3 Dr. Zuckerman can probably answer this very easily, but I understand that the  
4 COSIRA trial is not an IDE trial, the company chose to conduct this trial outside the United  
5 States and the FDA either wasn't privy or certainly didn't agree to how the COSIRA trial was  
6 designed. But educate me a little bit, Dr. Zuckerman. I'm struck by -- for an intravascular  
7 permanent device, safety at 6 months means quite short. What would be the normal  
8 standard that the FDA would hold a company to for a permanent intravascular device  
9 looking at safety?

10 DR. ZUCKERMAN: I think that's an excellent point that you're making, Dr. Allen. It  
11 would certainly be longer than 6 months and, in fact, we'll be asking the Panel that question  
12 this afternoon and what should be done to ensure that safety is verified.

13 DR. LANGE: Great. I've got Dr. Wittes and then Dr. Connor.

14 I'm sorry, you're on mute.

15 DR. WITTES: Yeah, Janet Wittes.

16 So my question is actually related a little bit to Dr. Allen's and then I have another  
17 one. The question related to Dr. Allen's is what aspects of this study did you not -- did the  
18 FDA not approve going forward? So that's a question for the FDA.

19 The other question I have, and I'm struggling a lot with the kind of -- the numerically  
20 marginal results. If I do a p-value in a different way, if I do a Fisher's exact test, that 0.024  
21 becomes a 0.04. If I look at the secondaries, they're not really very strong.

22 And then there's the issue about the 104 and the 124. And this question goes to the  
23 Sponsor. What did the DMC tell you? They had an interim analysis, they said to go on,  
24 which meant that you had some sense that there was probable benefit. What did you know  
25 when you decided to stop at 104?

1 DR. LANGE: Great. So we had asked that before and that is, why stop prematurely?  
2 I'll ask the Sponsor again to address that.

3 Dr. Zuckerman.

4 DR. ZUCKERMAN: Yes. So Dr. Wittes has asked a number of good questions and it's  
5 important to remember that initial negotiations for this trial began around 2010. So for the  
6 sake of being truly transparent and factual, in addition to the FDA explaining some of the  
7 negotiations around 2010, I think Dr. Wittes would benefit also from the Sponsor  
8 recounting what their interpretation of events was.

9 DR. LANGE: Okay, we'll do that.

10 Dr. Wittes, just to save the Sponsor some time, you asked some questions before  
11 that the FDA had in their presentation and if you've got a pen there, you asked about how  
12 many countries were involved, Slide 25 has that. Number of current smokers, Slide 53. And  
13 the movement is Slide 64 from the FDA. I'm sorry, you're on mute, Janet.

14 DR. WITTES: The movement was a little hard to -- I can figure it out. It would be so  
15 much nicer if it was just a table that showed the cross-classification, but I can figure this  
16 out.

17 DR. LANGE: Okay. Well, the Sponsor may still want to do that, but you may want to  
18 look at Slide 64 from the FDA, they had some of that information. But we'll still ask the  
19 Sponsor to address that to your satisfaction.

20 Dr. Connor.

21 DR. CONNOR: Yeah. And I had asked why they stopped. Janet, Dr. Wittes, had asked  
22 what they knew at the interim, which I think are two subtly but very importantly different  
23 questions. So I just wanted to point that out.

24 But then my question was actually it seems like the missing data issue is largely an  
25 issue with secondary analyses, but I assume that the primary outcome and the secondaries

1 were all collected at the same visit, you know, at the end of the trial, like the 6-month visit.  
2 So my question is one to the Sponsor, which is operationally, how do you have all this  
3 missing data in secondaries but not in the primary?

4 DR. LANGE: Okay, excellent point, so let me -- again, for the Sponsor, what did they  
5 know at the interim that may or may not have influenced their decision to stop  
6 prematurely? I appreciate that nuance. And then how do we have missing data when it  
7 was an in-person interview at 6 months and we have CCS but we don't have SAQ, which is a  
8 question.

9 DR. CONNOR: Yeah.

10 DR. LANGE: Okay, we'll ask the Sponsor to address that.

11 All right, Dr. Batchelor, I see your hand. And then I see Dr. Hirshfeld.

12 DR. BATCHELOR: Thank you, Dr. Lange.

13 So, you know, we've discussed the fact that we may not have to absolutely have a  
14 clear mechanistic explanation for this device and for any device for that matter, as long as  
15 safety and efficacy reaches certain thresholds. The animal data that was presented by the  
16 FDA was obviously inadequate. I think that was nicely presented, there was missing data  
17 and not all the data was presented in a fashion that the FDA is satisfied with. Furthermore,  
18 the gradients across this device were relatively small, 2-3 mm/Hg, which I suppose we  
19 would expect with this device. We don't have any sense for what would be a normal  
20 response.

21 Two quick questions. Number one, is myocardial information -- you know, they just  
22 presented data on histopathologic assessments related to the coronary sinus. If one of  
23 these things occluded and you're going to have healing in a neointimal fashion, it's going to  
24 be sort of a bell-shaped curve and some are going to heal, I'm sure, worse than others and  
25 some will even develop thrombus, and this study is too small to detect this stuff. I'm just

1 thinking, worst-case scenario, is the FDA satisfied with having very little animal data and  
2 just putting all of their eggs in the basket of clinical trial information and further studies,  
3 whether they're -- it's an observational registry or randomized trial? Or is it important to  
4 have sound animal data to come up with some sense for how this device actually, really,  
5 truly heals and whether or not it affects the myocardium in any way, shape or form? If you  
6 look at case reports of coronary sinus occlusion, it can lead to myocardial edema,  
7 pericardial effusions, and death in the extreme cases. So I'm trying to understand and  
8 understand how this device really heals in the absence -- to Keith Allen's point, in the  
9 absence of safety data that goes long enough in enough patients enrolled.

10 DR. LANGE: So, Dr. Batchelor, I'm going to turn that back around for a second. I  
11 think to ask the question is FDA satisfied with the animal data, I think the answer to that is  
12 no. The next question is do they need additional animal data, and I think they're going to  
13 direct that question to us and say are we satisfied and do we want additional data, and I  
14 know we're going to be talking about that in our discussion. The FDA is specifically not  
15 offering an opinion but actually soliciting ours for this particular thing.

16 Am I speaking -- I'm speaking on behalf of the FDA, Dr. Zuckerman. Dr. Zuckerman,  
17 did I speak correctly?

18 DR. ZUCKERMAN: I just want to underline the points that Dr. Lange has made.  
19 Today is not a time for FDA discussion. These are the questions on the table, Dr. Batchelor,  
20 among others, and we and the Sponsor look forward to a robust discussion by this very  
21 august, independent advisory panel.

22 DR. LANGE: Speaking of august, that goes to Dr. Hirshfeld. You're on mute, John.

23 DR. HIRSHFELD: So I guess that's my signal to leave.

24 DR. LANGE: No.

25 (Laughter.)

1 DR. HIRSHFELD: So I've been thinking a lot about the efficacy issue of this and the  
2 mechanism and I'd like clarification. It would seem logical that as part of this study one  
3 would get systematic imaging of the device at follow-up to determine to what degree the  
4 device has actually done what it's intended to do. A CT angio at 6 months would seem like a  
5 very logical way to figure this out and I don't recall seeing that in any of the panel material  
6 and I just want to be certain that that data is not available, because it would seem like such  
7 logical data to include.

8 DR. LANGE: Okay. So again, thank you, Dr. Hirshfeld, for that insight. We'll ask the  
9 Sponsor if they collected CT angio data on the patients in the COSIRA study and if so, to  
10 present that. Great.

11 Any other hands up? Any other -- all right. Let me propose, we're 10 minutes early.  
12 I think the robustness -- the discussion is going to be very robust and I'd like to preserve as  
13 much time for that. So with your permission, instead of a 45-minute lunch, could we --  
14 would you all -- I'm going to ask you to vote on a 45-minute lunch or a 30-minute lunch and  
15 I'll ask for a raise of hands and of course, I'm only going to see the ones that I'm voting for,  
16 but if you would prefer a 45-minute lunch at this point, can I see your hands?

17 (Show of hands.)

18 DR. LANGE: One, two, three. If you would prefer a 30-minute lunch.

19 (Show of hands.)

20 DR. LANGE: Okay, there are more, so it means that three of you will come back in 45  
21 minutes, the rest of you will come back in 30. All right.

22 George, do you want to say something?

23 DR. VETROVEC: Yes. The question is, does that allow enough time for the Sponsor  
24 to answer the questions and for the FDA? That seems to me to be an important question.  
25 There were a lot of questions raised. We'd all like a quicker lunch, but is that going to



1 interfere with the information we have asked for?

2 DR. LANGE: And I think I appreciate you bringing that up, Dr. Vetrovec.

3 With the FDA, will a shorter lunch be okay for the FDA? There were very few  
4 questions posed to the FDA.

5 DR. ZUCKERMAN: The FDA is okay with the shorter lunch, but we need to hear from  
6 the Sponsor.

7 DR. LANGE: And I'm going to ask them next.

8 Vicki --

9 DR. BEBEAU: Yes.

10 DR. LANGE: -- what would be your preference?

11 DR. BEBEAU: In lieu of the number of questions that we have, I'd really appreciate  
12 the 45 minutes.

13 DR. LANGE: We'd be happy to do that, okay.

14 DR. BEBEAU: Thank you so much.

15 DR. LANGE: I'm going to say it's 11:10, so if we could convene at 11:55, would that  
16 be okay with everybody? I'm sorry, it's 11:55 my time, but those of you on the East Coast,  
17 you're probably getting hunger pains right now because it's 1:05, so let's reconvene in 45  
18 minutes at -- yeah. Fifty-five, 1:55, 12:55, 11:55, and 10:55. All right. See you guys in 45  
19 minutes. And please turn off -- and I'll remind you not to discuss this among yourselves or  
20 with anybody else, and I'll also ask you to turn off your video and mute yourself and we'll  
21 see you in about 45 minutes. Thank you.

22 (Whereupon, at 1:09 p.m. a lunch recess was taken.)

23

24

25

1  
2  
3  
4  
5  
6  
7  
8  
9  
10  
11  
12  
13  
14  
15  
16  
17  
18  
19  
20  
21  
22  
23  
24  
25

AFTERNOON SESSION

(1:58 p.m.)

DR. LANGE: Great, it's just now a hair past 11:55 Mountain Time or 1:15 Eastern Time, 12:15 Central Time, and 10:15 -- 10:55 Pacific Time. My apologies.

We'd like to proceed with the Open Public Hearing portion of the meeting. Public attendees are given an opportunity to address the Panel to present data, information, or views relevant to the meeting agenda. These are all prerecorded, but before then, Aden will read the Open Public Hearing Disclosure Process Statement.

MS. ASEFA: Both the Food and Drug Administration and the public believe in a transparent process for information gathering and decision making. To ensure such transparency at the Open Public Hearing session of the Advisory Committee meeting, FDA believes that it is important to understand the context of an individual's presentation.

For this reason, FDA encourages you, the Open Public Hearing speaker, at the beginning of your written or oral statement, to advise the Committee of any financial relationship that you may have with any company or group that may be affected by the topic of this meeting. For example, this financial information may include a company's or a group's payment of your travel, lodging, or other expenses in connection with your attendance at the meeting. Likewise, FDA encourages you, at the beginning of your statement, to advise the Committee if you do not have any financial relationships. If you choose not to address this issue of financial relationships at the beginning of your statement, it will not preclude you from speaking.

DR. LANGE: Great. Thank you very much, Ms. Asefa.

FDA has received 13 requests to speak prior to the final date published in the *Federal Register*. Each speaker will be given 3 to 5 minutes to speak.

Free State Reporting, Inc.  
1378 Cape Saint Claire Road  
Annapolis, MD 21409  
(410) 974-0947

1 The first speaker is Meg Seymour from the National Center for Health Research  
2 Cancer Prevention and Treatment Fund. Erroneously in the agenda it lists her as being --  
3 this is a nonprofit organization, I just want to state that for the record. And with that, we'll  
4 ask 3D Communications to roll with the public hearing portion of these deliberations.

5 DR. SEYMOUR: Thank you for the opportunity to speak today on behalf of the  
6 National Center for Health Research. I am Dr. Meg Seymour, a senior fellow at the center.  
7 Our center analyzes scientific and medical data to provide objective health information to  
8 patients, health professionals, and policymakers. We do not accept funding from drug or  
9 medical device companies, so I have no conflicts of interest.

10 This Committee is being asked numerous questions regarding whether the data  
11 presented is sufficient to support the safety and effectiveness of the Reducer device and if  
12 so, for whom. The device in question is permanent and therefore should be held to a high  
13 standard when considering approval for a permanent device.

14 The FDA Executive Summary points out that there's a large placebo effect that has  
15 been shown in previous blinded studies in the field of refractory angina, including 42% of  
16 the sham group in the COSIRA trial. Given this strong placebo effect, a permanent device  
17 should not be implanted in patients without demonstrating meaningful improvements for  
18 patients compared to an effective sham-controlled group. And as FDA pointed out, it is  
19 important to know if the sham control is effective in blinding most of the patients where  
20 they did not know if they had the Reducer device. Moreover, studies must be conducted on  
21 a representative sample of patients.

22 Let's start with the shortcomings of the patients enrolled in the study. The patients  
23 studied are not a representative sample of refractory angina patients. Number one, more  
24 than 25% of patients in the study were on either one or no antianginal medications. In  
25 other words, they did not have refractory angina.

1           Number two, there was no information about medication compliance.

2           Number three, zero black or Hispanic patients in the COSIRA trial and more than 80%  
3 of the patients were white, male, or both. White males do not comprise most patients, so it  
4 would not be ethical to claim safety and effectiveness for a device that was not tested on  
5 most of the types of patients who will use it. And yet I don't think the Sponsor would be  
6 satisfied if the Reducer device was approved only for white men.

7           In addition to enrollment, the data provided by the Sponsor is severely lacking.

8           Number one, although there's a statistically significant improvement in the primary efficacy  
9 endpoint, we agree with the FDA that the sample is too small to provide reasonable  
10 assurance that the results reflect the outcomes for most patients. This raises questions  
11 about whether the improvement occurred by chance.

12           Number two, only one-third of the patients achieved the level of effectiveness  
13 required by the primary endpoint, and almost 30% of the patients experienced no change at  
14 all from the implant.

15           Number three, the 6-month follow-up period was too short to evaluate safety or  
16 effectiveness for a permanent implant.

17           Number four, as I previously noted, there was a large placebo effect which,  
18 according to the FDA, presents challenges for interpreting the data given the limited sample  
19 size.

20           The bottom line is that the true effectiveness of the device is unknown. We agree  
21 with the FDA review that additional data are needed to identify the patient population most  
22 likely to have a clinically meaningful benefit with the Reducer device, considering risk for  
23 the implant procedure and the device.

24           Now let's talk a bit more about the safety of the Reducer device. I'm sure you share  
25 my concerns and that of the FDA that the Sponsor did not define the pre-specified primary

1 endpoint for safety and instead relied on secondary safety endpoints, which is the  
2 successful implant of the device without a need for intervention.

3 To address some of the shortcomings in the data, the FDA has asked you to discuss  
4 and make recommendations regarding the Sponsor's proposal to perform a post-approval  
5 randomized sham-controlled trial in a country where Reducer is not approved. That would  
6 basically mean that U.S. patients would be the experimental test subjects for a product  
7 whose risks and benefits are not quite proven for them, especially for women and non-  
8 whites. Instead, the FDA should require the Sponsor applying for approval to do a better  
9 job of recruiting patients with refractory angina, as well as patients who are female and/or  
10 people of color. If they didn't achieve that when they were trying to get approval, we know  
11 the chances are even lower in a postmarket study.

12 In summary, we urge the Committee to require additional premarket data from a  
13 randomized, properly blinded, sham-controlled study with a representative patient  
14 population. Such data are necessary to determine the effectiveness of the device and  
15 whether it is worthy of approval and if so, for whom. The data it provided thus far simply  
16 do not provide enough data evidence that there are benefits outweighing the risks for this  
17 permanent device. Additional data are necessary before such consequential conclusions  
18 can be drawn. Thank you.

19 DR. LERMAN: My name is Dr. Amir Lerman, I'm Professor of Medicine in the  
20 Department of Cardiovascular Disease at the Mayo Clinic in Rochester, Minnesota. I serve  
21 as the vice chair for the department, but also the director of the Chest Pain and Coronary  
22 Physiology Clinic at the Mayo Clinic. I would like to talk to you today about the urgent need  
23 for a new therapy for patients with intractable angina.

24 In our chest pain clinic, in the last more than a decade we focused on patients that  
25 have intractable angina that has several phenotypes. One of them is the group of patients

1 that had that severe underlying coronary disease and already underwent several  
2 procedures and are not technically amenable to any new procedure. The other group of  
3 patients are patients that have intractable angina not responsive to treatment and these  
4 are the patients that have non-obstructive coronary disease and fall under the category of  
5 microvascular angina. Both of these patients have increased evidence of cardiovascular  
6 events, decreased quality of life, and increased cost of care. We are coming to your  
7 meeting in the need of therapy for this group of patients.

8         Current medical therapy is not designed and currently there are no guidelines for  
9 treating of patients with intractable angina. The current device, the coronary sinus reducer,  
10 in my opinion, what I've seen from the data and what I've seen of the device, it's a safe,  
11 easy-to-use and insert device that can bring us with solution to treat these patients. We  
12 currently have requested an FDA approval and have an IDE approval to use this device with  
13 patients with microvascular angina and we are planning to start this study in the next few  
14 weeks.

15         Thus, I think this is the right time for the approval of this relatively noninvasive  
16 device to treat this highly prevalent number of patients.

17         DR. KOENING: Good afternoon, my name is Gerald Koenig. I'm a general and  
18 interventional cardiologist at Henry Ford Hospital, which is a large tertiary medical center in  
19 Detroit, Michigan with a wide referral base in the state. I've been in practice for nearly 15  
20 years and a good portion of my patient population has either coronary or peripheral arterial  
21 disease, which I medically treat and manage and, when clinically appropriate, also treat  
22 with invasive interventions.

23         Since medical school at Harvard and my residency and fellowship at the University of  
24 Michigan, I've had the opportunity of encountering a countless number of patients with a  
25 diagnosis of refractory stable angina. With the prevalence in the United States of nearly 1.8

1 million patients with such a diagnosis, and an annual incidence of nearly 75,000 patients,  
2 that would be not surprising.

3           Currently, my arsenal of treatments for those no longer having a percutaneous or  
4 even a surgical option for revascularization, it's typically medications. Only a general  
5 category of four medication groups have shown clinical efficacy. But this option is not  
6 uncommonly limited; a patient's ability to tolerate the medications, all the pharmacy issues  
7 with interactions with other medications, and the overall effectiveness of the treatments.

8           Additional treatment options exist, such as the external enhanced counterpulsation,  
9 stem cell therapy, extracorporeal shockwave treatment, or previously excimer laser  
10 therapy. But they've been fraught with limited clinical success or cost effectiveness over  
11 prior years, as well as the limitations of the ability for patients to tolerate these regimens.  
12 At that point I find myself and my patients to be significantly debilitated by their limitations.  
13 This frequently leads them to become physically de-conditioned and in a worsened clinical  
14 state for their overall cardiovascular health. In addition, they frequently have a significant  
15 decline in their mental well-being, being more depressed, anxious, and oftentimes  
16 withdrawn.

17           Fortunately, more recently I've had the opportunity to be one of the only U.S.  
18 cardiologists to actually use the coronary sinus Reducer device by Neovasc. This was done  
19 through the FDA and our institutional IRB approval on a compassionate use basis for two  
20 patients that had severe limitations and markedly reduced quality of life with refractory and  
21 persistent stable angina. Both had residual ischemia demonstrated on stress testing and  
22 had been exhausted of current therapeutic options. They were quite miserable and  
23 depressed as a result.

24           This therapeutic option was initially brought to my attention already back in 2015  
25 when the COSIRA clinical trial results were published in the *New England Journal of*

1 *Medicine* and carried out through about 11 centers in the world. However, we've yet to  
2 access these therapeutic technologies in the United States despite having CE mark approval  
3 in Europe and regular success in treating patients based upon registry data.

4         These two patients, which I believe will also be providing their own testaments, have  
5 had quite a significant response in reducing their angina classification category and have  
6 much improved quality of life. Their success has brought a new ray of hope, really, to  
7 potentially treating a number of my other patients, and this also opened up referrals for my  
8 institutional colleagues as well as others at institutions across the street.

9         Again, having had the experience of implanting this device, I can attest to the  
10 relatively low-risk nature of the procedure and really any significant complications, and the  
11 relatively limited training and proctoring that's needed for a high rate of successful  
12 implantation.

13         Patient tolerance of the procedure is exceptional and much like that of a standard  
14 right heart catheterization procedure, and certainly much less than that of an arterial  
15 coronary or peripheral stent implantation. Recovery is rapid and results begin to become  
16 evident after about 6 weeks when the device has had a chance to endothelialize and  
17 actually form the narrowing or stenosis within the coronary sinus.

18         From my experience as a general cardiologist in both medical management to these  
19 challenging medically refractory stable angina patients, to an interventional cardiologist  
20 perspective in performing such a low-risk, high-yield procedure, I definitely feel this  
21 technology would make a significant impact in the quality of life in a countless number of  
22 patients in the U.S., and also be relatively cost effective, especially in the backdrop of the  
23 high healthcare cost associated with treating these patients with relatively few effective  
24 options. Thanks for the opportunity to speak with you today, and I'm always available to  
25 answer any further questions you may have.

Free State Reporting, Inc.  
1378 Cape Saint Claire Road  
Annapolis, MD 21409  
(410) 974-0947



1 DR. GINDI: My name is Dr. Ryan Gindi and I'm a clinical and noninvasive cardiologist  
2 at Henry Ford Hospital. I completed my cardiology fellowship a little over 2 years ago and  
3 during my training, I frequently encountered patients with angina refractory to medical  
4 therapy. The most common and best evidence-based treatment we had at the time was  
5 ranolazine, which in my experience has been marginally efficacious. There were multiple  
6 mechanical attempts to help these patients, such as mechanical counterpulsation and  
7 transmyocardial laser revascularization. This became a subject of interest to me and I had  
8 read about the coronary sinus Reducer, but it wasn't available in the U.S.

9 In my last year of fellowship I cathered a patient who had severe coronary disease and  
10 multiple prior interventions and he was suffering from severe anginal symptoms refractory  
11 to medical therapy including ranolazine. He had no targets for revascularization. It  
12 occurred to me that this would be an ideal patient for a coronary sinus Reducer and I had  
13 mentioned this to my attending during the cath. I was then challenged with obtaining FDA  
14 approval for compassionate use.

15 To make a long story short, I succeeded, and this patient improved from having chest  
16 pain while walking to his mailbox, to 3 months later walking around Washington, D.C. for  
17 several hours without any symptoms. He recently rode a bicycle 35 miles and was symptom  
18 free.

19 I have since obtained FDA approval on a total of six patients. However, we only  
20 implanted one additional patient who had similar positive results. The other five patients  
21 were postponed due to COVID-19. These patients are all suffering and all would have been  
22 implanted over 1 year ago if it wasn't for delays related to obtaining FDA and IRB approval.

23 Following my experience, there are cardiologists throughout our organization and  
24 even the entire country who have reached out for help with their patients suffering from  
25 refractory angina. Angina refractory to medical therapy is a common dilemma we face in

1 cardiology. There have been multiple mechanical interventions proposed; however, none  
2 have panned out over the past several decades and some have only posed higher risks to  
3 our patients. I am baffled that a high-risk CTO intervention that carries significant mortality  
4 risk is acceptable for treating refractory angina and a low-risk coronary sinus Reducer  
5 implantation that has quality data supporting its use and is in the European Society of  
6 Cardiology guidelines, remains unavailable to the majority of our patients. I strongly urge  
7 the FDA to consider approval of this device, which has the potential to help many patients  
8 across the country with a very low risk of even minor adverse outcomes.

9 MR. CASCIANO: My name is Fred Casciano. I'm 58 years old, I live in Traverse City,  
10 Michigan, and was diagnosed with heart disease at the age of 38, with five-way bypass at  
11 41. I had angioplasty with three stents placed within 6 weeks of my bypass surgery. I was  
12 then free of angina until 51.

13 Prior to angina, I was very active, I worked 60 to 70 hours a week, started my own  
14 painting business with a crew of four to six employees. I built five homes, I maintained a  
15 home and a cabin, mowed an acre and a half of lawn with a walk-behind mower and snow-  
16 blown 300 feet of driveway with a walk-behind blower. I cut and stacked firewood,  
17 enjoyed snowmobiling, bowling, and waterskiing.

18 But my life disappeared after angina diagnosis. I suffered severe depression not  
19 knowing from day to day if I was going to live or die. I constantly worried about dying. I  
20 worried about my wife if I was to die. I felt guilty for stress and responsibility put on my  
21 wife. I worried about the cost of medications and medical bills. I also had anxiety for each  
22 time I went to the hospital for a procedure or surgeries, that I would not wake up. Unable  
23 to sleep more than an hour at a time, I was afraid I would not wake up. I was unable to  
24 walk 300 feet to the mailbox, unable to snowmobile, snow-blow, or mow the lawn. I  
25 consistently was reminded of my heart disease with every episode of angina.

1 Under the care of Dr. William O'Neill, Dr. Ryan Gindi, Dr. Gerald Koenig at Henry  
2 Ford Hospital in Detroit, I was informed there was a new procedure that might help people  
3 with my condition. I had the procedure January 2018 and since then my quality of life has  
4 returned. I know that I still have heart disease, but since I received the Reducer I am able to  
5 get back to walking. Six months after the surgery, I went for a 35-mile bike ride. I am able  
6 to get back to doing things around the house, I have been able to paint rooms at my home  
7 and for my mother. I feel my depression is decreasing and I am enjoying life again. I feel  
8 very lucky I received this device, it has given me hope again. Please consider my story as  
9 you make your decision today. There are many other people out there like me suffering  
10 from angina who need this kind of hope, too.

11 MS. CASCIANO: My name is Annette Casciano. My husband is Fred Casciano. We  
12 have been together for 33 years, married for 18 years. Fred was diagnosed with heart  
13 disease at the age of 38. By 41, he had survived several heart attacks and a five-way  
14 bypass. About 6 weeks after his bypass surgery, he had to undergo angioplasty with three  
15 additional stents. Fred was always very active, extremely driven to succeed in his goals.  
16 We always enjoyed boating, snowmobiling, spending time with family and friends. Fred  
17 started his own painting business with several employees and we built five houses.

18 When Fred began to suffer with chronic angina at age 51, you could see him decline  
19 physically and mentally. He became angry that he was unable to do things around the  
20 house. The simple task of walking 300 feet to the mailbox was a struggle. He felt like he  
21 was a burden to me. He worried about the cost of medications, medical insurance, and  
22 procedures. While he was worrying about these things, I worried about him constantly. I  
23 watched the toughest guy I know giving up on life at 57. I was always afraid I would get a  
24 phone call or I would come home to find him suffering a massive heart attack or find that he  
25 had passed away.

1           When we were told about the opportunity for Fred to receive the coronary sinus  
2 Reducer, we were given something we haven't had in a long time, which was hope. We  
3 understand that Fred's cardiac disease would still be there, but having the hope that Fred  
4 could possibly be able to do things he used to do and we could get back to a life that we  
5 had before, one without so much fear.

6           Since the coronary sinus Reducer, Fred's energy has improved, he's become more  
7 confident, he is able to do things without the constant reminder of his heart disease. His  
8 depression has decreased, he's laughing again. Family and friends have noticed a big  
9 improvement. I feel that I have my husband back and the stress that I had been living with  
10 all this time is also becoming less. I don't worry about losing him on a daily basis. I am so  
11 thankful Fred was given this opportunity for this device. It has definitely made an  
12 improvement in our lives and brought us renewed hope for the future. I hope that you can  
13 find it in your wisdom today to bring this kind of improvement and hope to patients and  
14 their families that struggle with angina. Thank you for your time.

15           MR. SOBERANO: My name is Mark Soberano, I worked as a nurse anesthetist for  
16 special needs families at Cincinnati Children's Hospital for 20 years. Now my life is far  
17 different. On January 18th, 2016, after ignoring 18 hours of chest pain, I had a cardiac  
18 arrest. I was revived by my wife and children and taken to a regional hospital where I  
19 received one stent.

20           Since then, I have continued to have 8 out of 10 chest pain with any kind of activity  
21 and continually need nitroglycerin. I am on disability. I cannot walk for more than 50 feet  
22 without having 6 out of 10 chest pain. I take nitroglycerin at least five times a day.  
23 Although I'd like to take more, the side effects preclude that. I am a husband, father of  
24 seven children with ages 28 to 11. But as time goes by, my condition gets worse. I feel like  
25 my family's life is moving forward while mine is standing still. My five sons are all in

Free State Reporting, Inc.  
1378 Cape Saint Claire Road  
Annapolis, MD 21409  
(410) 974-0947

1 scouting and we used to go camping with the scout troop at least once a month. Now I am  
2 unable to join them for even one night. As a family, we would take hikes in the parks near  
3 our home. Now I get chest pain when I walk anywhere outdoors. I have had trouble  
4 walking my service dog. My daughter's a gymnast and I cannot attend a full meet without  
5 being exhausted.

6 This treatment will not change my diagnosis. I will still have complex microvascular  
7 cardiac disease for the rest of my life, but it will give me back some quality in life, the kind  
8 of life that I can only dream of where I actually used to be part of my family, like seeing my  
9 kids participate in a favorite sport and enjoying the outdoors. I beg you to remember that  
10 this treatment has a high potential to change my life so that I can be an active participant in  
11 my world and the world of my family, something that you take for granted, but I cannot.  
12 Thank you.

13 MR. SUMMERS: My name is Steven Summers, I'm 58 years old and I live in Macomb,  
14 Michigan. My first heart attack was in 2012 at the age of 50. I experienced chest pain and  
15 angina as well as chest pressure. I believe at that time I had two stents placed during the  
16 catheterization procedure. The year prior to the first heart attack, I had started seminary  
17 with the goal of becoming a full-time minister. As a result of the heart attack and  
18 subsequent periods of angina and continued procedures, I was unable to continue the  
19 pursuit of a master's degree.

20 My cardiologist has, I believe, run out of options, that is, room for any additional  
21 stents for controlling the angina. I am also not a candidate for bypass surgery. I'm  
22 currently on disability as a result of a mild stroke in 2018, combined with the continuing  
23 heart issues. Prior to the cardiac sinus Reducer procedure, I relied on the nitroglycerin  
24 tablets to ease any chest pain. Any exertion, including taking the garbage to the curb,  
25 mowing the lawn, going for a walk around the block, going up stairs, would have caused

1 chest pain. My cardiologist mentioned the sinus Reducer as a possibility in 2017 and got  
2 approval in 2018. The first 6 months following the procedure, I was still experiencing  
3 angina pretty regularly. The past year and a half have been nearly angina free. I still use  
4 the nitroglycerin occasionally. Rarely, actually, two to three times a month, but the pain is  
5 not nearly as intense and is rare.

6 The angina I experienced prior to the sinus Reducer procedure was tough on my  
7 wife, especially. But now the angina has subsided considerably as a result of the sinus  
8 Reducer procedure. My wife and I were recently on a weekend and we're able to bike ride,  
9 horseback ride, and walk with only one instance of needing a nitro pill. It's actually pretty  
10 rare now that I even think about angina, but when I do, I'm reminded of the life-changing  
11 sinus Reducer procedure and the remarkable improvement of my life. Please remember my  
12 story as you make a decision. Thank you.

13 MS. SUMMERS: Hello, my name is Donna Summers and I have been married to  
14 Steven Summers for 33 years. Steven had his first heart attack at 50 years of age, 8 years  
15 ago, and was diagnosed with diffuse narrowing of his cardiac artery. Before his diagnosis,  
16 we lived active and a busy life with work, church, family, activities, and vacations. And all of  
17 that changed 8 years ago.

18 Living with a husband with angina is worrisome, constantly aware of the need for  
19 potential medical care. I've lost count of the number of stents Steven has had placed and  
20 the number of hospitalizations. Our life changed from being able to just pick up and go to  
21 needing to remind Steven as we leave the house on meds, insulin, nitro. The anxiety and  
22 fear of the nitro not working, the fear of the limiting impact on our family was our constant  
23 companion.

24 Steven took a lot of medications and any activity would trigger the chest pain.  
25 Simple things lots of people would take for granted, like enjoying a walk with the grandkids,

1 that was gone. And other things could trigger it, as well, that were out of our control.  
2 Basically, our family and our church activities revolved silently around Steven's limitation.  
3 We'd always try, but we could tell he was struggling and because he could not work any  
4 longer, he needed to apply for disability, something we never thought would happen. The  
5 sadness that our life was becoming was so limited by this disability and it was profound for  
6 Steven. He struggled with depression and feeling useless and feeling old at a young age.

7         Thankfully, we were told about this procedure, a device that could help Steven's  
8 quality of life and so we talked it over at length with the team and decided it would help.  
9 For the most part now, Steven can do daily work without chest pain and care, it's very rare  
10 for him to have to and he can take walks with the dog and play with the grandkids and  
11 we've gone on two vacations recently and gone on long walks without any chest pain. He is  
12 more hopeful and so am I, as we've just had a 1-year anniversary without any  
13 hospitalization.

14         We know that his cardiac condition is serious, but are appreciative of every moment  
15 we get together. But with this procedure, our quality of life has improved exponentially  
16 since it has allowed us some semblance of normality and enjoyment of basic abilities.  
17 Without this procedure, we were becoming very limited to our home. This procedure gave  
18 us hope for a life maximized instead of a life limited by his heart issues. This would not  
19 have been possible without the procedure and I hope other families are given this  
20 opportunity, as well. This illness impacts the health and well-being of the spouse, too.

21         MS. VANDENBOSSCHE: Good afternoon, my name is Laurie Vandebossche. I live in  
22 Port Huron, Michigan, and I'm 59 years old. In the past two and a half years I've had two  
23 heart attacks and was diagnosed with angina. In those two and a half years I've been  
24 hospitalized several times and had six heart catheterizations, angioplasty, and stents. In  
25 spite of these health issues, I've had an absolutely wonderfully blessed active life with my

1 husband and children. Unfortunately, since suffering from angina, my life has changed  
2 dramatically. I truly hope by sharing some details you will gain a better understanding of  
3 what it's like to live with angina.

4 I had planned on working for another 5 years, but had to take an early retirement  
5 from a very rewarding job that I loved, teaching struggling first graders how to read. In my  
6 retirement, I was hoping to be involved in volunteer work, but I'm afraid to commit for fear  
7 I won't be physically able to fulfill the duties expected of me. Although I've always been a  
8 positive person, I now struggle with anxiety and depression. I would love to wake up in the  
9 morning and not have the feeling of being stressed about whether I'll have angina today on  
10 my morning walk or throughout the day or after I eat or doing yard work or household  
11 chores or even running errands like grocery shopping.

12 I was very social, but now I'm often withdrawn and fearful. I've always enjoyed  
13 walking with friends, but now I make excuses because one of my angina triggers is the  
14 feeling of indigestion and belching, which is not only uncomfortable but embarrassing. I'm  
15 fearful to make any travel plans because of the frequency of my hospitalizations.

16 Last year while on vacation in Chicago, I wound up spending 4 days at Rush Medical  
17 University Hospital. It breaks my heart that our family vacations have had to drastically  
18 change because of my angina.

19 Although I realize this procedure will not cure my cardiovascular disease, I'm very  
20 hopeful it will provide some relief from the deep anxiety I experience every day and help  
21 restore some normality to my life and the lives of my family. After consulting with my  
22 cardiologist, I know this may be my only option to enhance and improve my quality of life.  
23 It has the potential of being a life-changing experience for all those suffering from this  
24 debilitating disease. Thank you for allowing me to speak to you today. I hope and pray you  
25 seriously consider approving this procedure.

Free State Reporting, Inc.  
1378 Cape Saint Claire Road  
Annapolis, MD 21409  
(410) 974-0947



1 MR. SCOTT: Hello, my name is Don Scott and I am here to talk to you today about  
2 my experience living with angina. I am 76 years old and was working full time until I was 71.  
3 I started cutting grass as a teenager, worked as a union teamster on a dock at a trucking  
4 company loading and unloading trucks, and then I spent the last 30 years of my career  
5 working as maintenance and grounds manager at a retirement apartment complex. In all  
6 the time I worked, I carried things, fixed things, and was able to help people.

7 Since retiring for medical reasons, I am no longer able to do even some of the most  
8 basic things that I have done my entire life. My wife cuts the grass at our home and I can  
9 barely walk up or down a flight of steps without being winded or having to take a break.  
10 Even walking from my car to the baseball field where my grandson plays is a challenge. I  
11 really miss the quality of life that I used to have and I'd love to be able to help take care of  
12 my own home and my own family and pick up my grandchildren without fear of getting light  
13 headed or dizzy.

14 I believe I am a good candidate for this procedure and hope you approve it so I can  
15 bring back some of those things that I used to be able to do. Please remember my story as  
16 you make your decision today.

17 MS. HOPKINS: Hi, my name is Tammy Hopkins. I am a 42-year-old female that lives  
18 in Cincinnati, Ohio with my husband. I was diagnosed with angina at the age of 29 and had  
19 a massive heart attack 6 months later that required me to undergo a CABG, which led to  
20 having a stent placed several months later and every year after that, for 4 years, with  
21 continuous hospitalization every 3 to 4 months in between having my stent reopened and  
22 countless surgeries, or just having pain from angina with excruciating pain and being very  
23 short of breath.

24 Before having angina, I used to swim all the time and enjoy life with my dogs and my  
25 family, but now my quality of life is limited because I'm short of breath all the time, never

1 know when my angina is going to come on or how long it's going to last. It limits me to do a  
2 lot of activities like going to Six Flags, riding rides, playing and running with my grand babies  
3 or just normal day walks with my husband and my dogs. I am a nurse and love what I do,  
4 but I have a lot -- I have to give up my passion and life because my everyday normal life of  
5 living with angina is too much for me. Just regular everyday breathing and watching TV or  
6 walking devastates me because of my excruciating pain.

7 I realize that a treatment like this won't change the fact that I have a cardiac  
8 condition, but it can give me a better quality of life that's vibrant, joyful, and worth living.  
9 My mother and sister died at a young age of a massive heart attack, as with myself with no  
10 other options available to me, I live with the same fear. So when I heard of this procedure, I  
11 thought it could be a helpful option of what it could bring to me and my family, as well as  
12 others like me in the same space. I always say if there is a will, there's a way. So please  
13 take into consideration of what this would mean to a lot of people like me and our loved  
14 ones because all we ask for is a chance for a better life. Thank you.

15 MR. HART: Good afternoon and thank you for allowing me to be one of the speakers  
16 today. My name is Clyde Hart and I'm a 74-year-old heart patient with stable angina and  
17 severe breathing issues, and I know I have the best cardiology team trying to guide me  
18 through these issues and help to improve my quality of life. I have had one open heart  
19 surgery in 2014, numerous tests, 14 stents, a triple A surgery, high blood pressure, and  
20 diabetes. My doctors tell me I am not really a viable candidate for any re-vascular  
21 procedures.

22 My cardiologist heard of this procedure and began to research the possibility of  
23 getting this done. He knew it had not been approved in the United States but was being  
24 used overseas and had some promising results. He asked me if he could apply to our FDA  
25 for compassion use approval. He sent in a request and after a short period of time, we

1 were approved. I am not sure how many others across the U.S. were approved, but I do  
2 know there were several. In this process the COVID pandemic came up and a doctor from  
3 Israel has not been able to come, he is going to proctor the procedure with Dr. Balarn (ph.)  
4 and therefore, even though approved, we will have to wait for the infection rate to diminish  
5 prior to getting the device implanted.

6 The quality of life for me has been lessened in the past 3 months, mostly in the  
7 ability to breathe well, and also the angina issue. I think I'm a pretty tough old guy who can  
8 get along with the angina, but the breathing issue really hinders my ability to get around. I  
9 have always been a very active person and honestly, at this point, I would rate my quality of  
10 life about a two. I realize I am not supposed to feel like I did when I was 25, but not being  
11 able to do things with my grandson, not being able to walk into my church and having to  
12 take three times as long for my wife if I have to go anywhere, I would love to be able to go  
13 my grandson's ball games, but would not be able to walk from the parking lot to the field  
14 and it would be too much of a burden on the rest of the family.

15 I would urge you to help with getting this device approved and our doctors trained  
16 to help patients like myself to hopefully have a little better quality of life. Again, thanks for  
17 trying to give me a reason to go on with my life.

18 DR. LANGE: I want to thank each of the 13 speakers. Your perspective and your  
19 experience is valuable, and so thank you for taking the time to share that with this Panel. I  
20 now pronounce the Open Public Hearing to be officially closed and we will proceed with  
21 today's agenda.

22 We'll now begin the Panel deliberations. Although this portion is also open to public  
23 observers, public attendees may not participate except at the specific request of the Panel  
24 Chair. Additionally, we request that all persons who are asked to speak identify themselves  
25 each time. This helps the transcriptionist identify the speakers.

1           During the next approximately hour and a half or so, we will open the floor to  
2 questions for both the Sponsor and the FDA. We have a number of questions that have  
3 already been posed and what I'll do is take this opportunity to repeat the question from  
4 each of the panelists and allow the Sponsor an opportunity to respond to those questions.

5           So is the Sponsor prepared, are they ready? Great.

6           DR. BEBEAU: This is Vicki Bebeau.

7           Yes, we are. I'm actually going to turn it over to Dr. Stone, who will moderate for us.

8           DR. LANGE: Okay. All right, Dr. Stone.

9           DR. STONE: Well, thank you very much. We really appreciate the effort and the  
10 time you're taking today to hear all of the details about the Reducer and I do think that the  
11 patient summaries were actually very powerful.

12           But we have spent the entire lunch period and we will try to address every single one  
13 of your questions. What we're going to do is kind of -- some of them overlap each other, so  
14 we're kind of bunching them in groups. So we may not specifically say who addressed the  
15 question, but we'll try to address every single one and if there are any we have not  
16 addressed at the end, please ask us.

17           DR. LANGE: All right. So, Gregg, let me say, I'm going to do this the opposite. I'm  
18 going to go by the speakers and pose the question and allow you to answer and that way --

19           DR. STONE: I see.

20           DR. LANGE: Yeah, but I agree, that's great to have additional clarifications. But I  
21 thank you very much.

22           DR. STONE: Okay, so that's fine. So yes, we are ready, and myself and Ms. Bebeau  
23 and Dr. Shmuel Banai will be answering the questions, as well as the statistician.

24           DR. LANGE: Terrific. Let me start with Dr. Somberg's question, he asked about what  
25 you can tell us about the pharmacologic maximization for the 30 days prior to the

1 procedure.

2 DR. STONE: Yes. So most of these patients had had chronic refractory angina. The  
3 inclusion criteria were, yes, just greater than 30 days and I will say it was not collected in  
4 the case report form as to the number of months or years that they had refractory angina,  
5 but they had to be on stable medical regimen for at least 30 days, and that was also  
6 reflected in the Seattle Angina Questionnaire forms. So we do have a lot of stability data on  
7 medications post-randomization which we can show you, if you'd like to see that now or  
8 wait until later.

9 DR. LANGE: No, please, now would be a good time.

10 DR. STONE: Okay, so let us -- let's pull up those. Let's see, let me look here. So first  
11 let me show you this slide. We've got, I think, three separate slides I can show you. So this  
12 slide shows you the antianginal medications at baseline and at 6 months, and you can see  
13 the Reducer arm compared to control at baseline, Reducer arm compared to control. This is  
14 just a number of antianginal medications. And basically, what you can see that is in both  
15 arms, just being on a class of medication was fairly stable. Perhaps one exception was  
16 nitrates, there tended to be perhaps slightly greater control group patients taking nitrates  
17 at follow-up compared to the Reducer patients, but pretty much it shows stability in the  
18 number of angina medications and the classes of angina medications. So now I think that  
19 we can --

20 DR. LANGE: Can you go back to that slide? Dr. Somberg, you had asked for this  
21 information. Can you go back to that slide, please? Dr. Somberg.

22 DR. STONE: Yeah. Can you pull up that slide again, please?

23 DR. SOMBERG: Can you hear me?

24 DR. STONE: Yes, John.

25 DR. SOMBERG: Okay. Well, two points. One is, so you're saying that there was

1 really no attempt at maximization of therapy in the pre-randomization and --

2 DR. STONE: Well, no, no. No, no, I'll get to that. That actually was not the case.

3 DR. SOMBERG: Okay, so then -- okay, you'll get to it. And the second thing is do we  
4 have any idea -- I mean, we know the class of drugs, but we don't know the dosages and the  
5 -- you know, was there any up-titration over the course of the study?

6 DR. STONE: Both of those two issues will be in the next two slides. So let me show  
7 you these next two slides and then we can address whether there are any other questions.  
8 So first, here I'll show you this slide, which is the mean antianginal medication dosages at  
9 baseline and 6 months.

10 And before even talking to this slide, let me just state that all the sites in the  
11 investigational protocol and the discussions with the sites, before the patients were  
12 enrolled, was that they and the general cardiologist, the interventional cardiologist and the  
13 cardiac surgeons, had to optimize the medications that these patients were on and they had  
14 to optimize all the revascularization options.

15 So as you can see, most of the patients had had -- almost all the patients had prior  
16 revascularization, most of them had had multiple revascularization procedures. All of these  
17 patients were being evaluated constantly for years in their lives by interventional  
18 cardiologists and cardiac surgeons and they still had Class III or IV angina. All of them were  
19 also constantly being evaluated for maximal medical therapy.

20 Now, when we get to it, and I'll show this in the next slide after this one, the reason  
21 that many of these patients were not on what you might consider four or five drugs is  
22 multiple reasons. As you know, the more drugs you're on and these patients are on not  
23 only antianginal drugs but drugs for diabetes, drugs for hypertension, drugs for  
24 hyperlipidemia and multiple other conditions, it's very difficult to take 8 or 9 or 10 drugs a  
25 day. Number one. Number two, these drugs, and in particular, beta blockers and nitrates,

1 have a lot of side effects and the patients were almost uniformly tried on these drugs, but  
2 many of them could not tolerate the drugs. And as I'll show you, the usage of the drugs  
3 were very similar to what we've seen in every other refractory angina study. So let me  
4 show you these two slides and then I can address any other questions.

5 First of all, this slide that you're seeing right now speaks to the constancy of the drug  
6 dose over time and you can see, basically, this is just a very simple counting up milligrams of  
7 beta blockers at baseline and in the Reducer baseline and then at 6 months.

8 DR. SOMBERG: What does this mean?

9 DR. LANGE: Yeah, I mean, the people are on 10 mg of a beta blocker? What does  
10 that mean?

11 DR. STONE: So as you know, there are different beta blockers and -- you know, from  
12 bisoprolol to metoprolol, etc., these are the means and there -- and I'll be the first one to  
13 acknowledge that there are tables where you avert the average daily doses from one to the  
14 other, that was not --

15 DR. SOMBERG: That's well established.

16 DR. STONE: And that's well established, but that was not done here. This was to try  
17 to get you whatever information was available at short notice and here you can see what  
18 we're able to show you.

19 DR. LANGE: Okay. Gregg, that's fine. Let's go to the next slide. In other words --

20 DR. STONE: Okay.

21 DR. LANGE: -- here's what I'm saying, we have a lot of questions and I want to get  
22 through them.

23 DR. STONE: Yes.

24 DR. LANGE: And so I want to keep our answers short and succinct.

25 DR. STONE: I appreciate that. At least what this shows you is similar doses in both

1 groups and constancy over 6 months, at least with the manner that these data were  
2 calculated.

3 Now let me show you this next slide, and this is a slide that I don't believe FDA has  
4 seen or been able to adequately evaluate, but what this shows you is a literature review of  
5 the principal or the largest refractory angina trials that have been done, RENEW, ACT34,  
6 AGENT-3 and 4, the OPTIMIST registry, and COSIRA, and these are all refractory angina  
7 patients. And here you can see just the -- at baseline, the medications that were used,  
8 nitrates, calcium channel blockers, and beta blockers.

9 Now, in those other four studies, those studies were actually performed before  
10 there were other antianginal drugs available with the exception in RENEW of ranolazine, so  
11 those asterisks basically are drugs not available. You can also see in COSIRA that there are  
12 patients on other antianginal drugs, as well.

13 So I think when you add all of this up, it's approximately similar to all refractory  
14 angina trials, the drugs that patients are taking and able to tolerate, given side effects,  
15 hypotension, etc.

16 DR. LANGE: So Dr. Somberg, do you have enough -- I think you --

17 DR. SOMBERG: Well, I appreciate the Sponsor going to this effort, but I do not think  
18 either of these slides answered the question.

19 DR. LANGE: Okay, all right. Sufficient.

20 Dr. Page had asked about the mechanism of action with regard to the Beck, knowing  
21 if there was both CS narrowing at the procedure and also asbestos applied to the  
22 pericardium and occluder rods, but the comments from the Sponsor that it was -- relief of  
23 angina was due to CS, is there any data that the Sponsor can offer to document that?

24 DR. STONE: Dr. Banai will address that.

25 DR. BANAI: Hi, Shmuel Banai here. Thank you for the opportunity to answer a few



1 of your questions.

2 The Beck operation included several different operations. The Beck 1 operation in  
3 human included both narrowing of the coronary sinus, as you see here, with a metal stylet  
4 in the diameter of 3 mm, within the ligature, around the coronary sinus and the metal stylet  
5 and then taking out the metal stylet and leaving out with the ligature that narrowing the  
6 coronary sinus for 3 mm. This was done in human with -- in some of the patients scrubbing  
7 with asbestos the pericardium, and the thought was to promote angiogenesis inflammatory  
8 response and more blood vessels coming from the chest wall into the myocardium.

9 But Claude Beck also performed experiments in animal, in dog model, of severed  
10 circumflex artery and he showed that only with increasing the pressure in the coronary  
11 sinus by his operation he can double the amount of backflow coming from the distal cut of  
12 the severed circumflex artery by increasing the pressure or relieving the pressure in the  
13 coronary sinus. So even without asbestos, scrubbing of the pericardium, just the narrowing  
14 of the coronary sinus, he's shown in animal models that it preserves myocardium and  
15 improves ischemia.

16 DR. LANGE: Dr. Page also had asked about the narrowing in the coronary sinus, it  
17 looked like a reverse hourglass and was that -- that's on Slide 45. Was that the intended  
18 result of the procedure or something unexpected?

19 DR. BANAI: Yeah, let me take you through the mechanism of action and what really  
20 happens when you implant the Reducer, and I'll also address the question that was about  
21 the pressure gradient across the Reducer.

22 So when you implant the Reducer in the coronary sinus, you intentionally oversize it  
23 by 10%. The intention of over-sizing has two reasons. One reason is to uncork (ph.) the  
24 device in the coronary sinus and prevent migration and the other reason is that both ends,  
25 both wide ends of the Reducer are engaged into the vessel wall and the interaction

1 between that embedment of the struts into the vessel wall start a process of tissue  
2 proliferation, what we call injury induced smooth muscle cell proliferation, tissue  
3 proliferation. And this process takes a few weeks, sometimes 8 or 6 to 8 weeks and this  
4 tissue proliferation grows over the struts to cover the struts and the gaps between the  
5 struts, and in the histology cuts that you see below from an animal model, you see the  
6 proximal end, the distal end, and the midsection.

7 Now, sometimes the midsection is not covered with tissue, which is fine. The  
8 narrowing is caused by tissue proliferation mainly at both ends and what you end up in this  
9 specific example is a tube-like narrowing. But really it doesn't matter, it can be both sides  
10 throughout the Reducer, only one spot. The idea is to get the pressure gradient across the  
11 Reducer, which does the work.

12 And this is an example of a pressure gradient measured in animals 6 months after  
13 implantation of the Reducer, you see almost 3 mm/Hg pressure gradient. Remember that  
14 in the pig there is azygos, big azygos vein, so it might not be representative of what we have  
15 in human. In human, we haven't measured the coronary sinus pressure because we felt it is  
16 unethical to expose this patient for another right heart catheterization 2 months after the  
17 procedure, but we do have CT angiography. May I have Slide AE-10, please? We have CT  
18 angiography done in more than 44 patients, both from the COSIRA trial and from the first-  
19 in-man trial. All of the CT angiography were analyzed by an independent core lab.

20 What you can see here is a CT image showing a contrast flow through the coronary  
21 sinus, coming from 12 o'clock and draining into the right atrium, which is at 6 o'clock. The  
22 central narrowing of the Reducer is clearly seen and contrast is flowing only through the  
23 center. You don't see contrast around the struts or outside the center and you clearly can  
24 see that the Reducer is not occluded.

25 Now, out of the 44 CT angiography analyzed, there was no device fracture, all

1 devices were patent, there was no thrombosis at all in none of them. Also, there was no  
2 distinguishable differences in the CT angio findings between the COSIRA Reducer and the 12  
3 years follow-up of the first-in-man Reducer. Additionally, there were 11 CT angio from the  
4 first-in-human that was performed at 3 years follow-up with similar results.

5 One more thing I wanted to show you, this is another example of 3-D reconstruction  
6 of the heart, the back of the heart, showing the Reducer and a cross-sectional area. And  
7 maybe I'll show another example talking about the tissue growth and the neointimal growth  
8 in the Reducer.

9 If you will pay attention to the bottom picture on the left, if we look at the most left  
10 image, you will be able to see -- I don't know if you see it large enough -- you'll be able to  
11 see a transverse cut of the proximal segment of the Reducer and you will see the white  
12 struts, which is the outer circle which is whitish. At the center you will see a contrast of the  
13 true or residual lumen and in between the struts, the white strut and the gray lumen where  
14 we see the contrast, you see a black circle aligning in between and separating the true  
15 residual lumen and the struts. This is actually the tissue growth that grow in response to  
16 the interaction between the stretching of the vessel wall by the Reducer and the vessel  
17 wall. So this tissue growth caused the narrowing along the Reducer.

18 DR. LANGE: Okay.

19 DR. BANAI: May I have also --

20 DR. LANGE: Dr. Banai. Dr. Banai, I'll get to those things. I just want to make sure  
21 that for each of the questioners, their questions were answered and so we'll talk a little bit  
22 more about this, the CT angio, in just a moment. But Dr. Page had asked previous questions  
23 about the Beck and the CS narrowing.

24 Dr. Page, have they adequately --

25 DR. PAGE: Well, yeah, the Beck wasn't as much a question as a comment, but I

1 appreciate the response and the confirmation was that the human procedure included the  
2 pericardial asbestos and whatever happened in pigs or dogs or whatever it is, I haven't  
3 looked at those data, but the humans all received both procedures is what my  
4 understanding was and I just heard confirmed.

5 The issue of the blood flow I thought was really interesting from the pig, but as was  
6 pointed out, that has significant drainage to the coronary sinus that humans don't enjoy.

7 And then finally, the CTA, I think, is interesting. I'm not an angiographer. I'm an  
8 electrician, not a plumber, so I'll leave it to my plumbing colleagues to comment on what  
9 they're seeing there. I do wonder whether that perfect hourglass is typical of the 44 angios  
10 that were obtained or whether that is, as most of us would do, a really good example of one  
11 case where the perfect hourglass is provided, and I'll leave my comments for the later  
12 discussion.

13 DR. LANGE: Okay, so keep those comments because we will come back, there's  
14 more issues about the CT and the thrombosis that was reported in the Reducer study, so  
15 we'll get around to that. So thank you, Dr. Page.

16 But Dr. Banai, thank you for answering his questions regarding --

17 DR. BANAI: Just a moment on the comment.

18 DR. LANGE: No, no. I mean, there will be additional time, in other words, we'll make  
19 sure that you have -- I just want to address it by each speaker so we don't leave anything  
20 out and so you'll have an opportunity to speak, I promise.

21 Dr. Allen asked, in terms of the evaluation of whether people were revascularization  
22 candidates or not, did a surgeon -- was involved with that decision, that the patient had no  
23 other options? And/or an interventional cardiologist?

24 DR. STONE: This is Gregg.

25 So in general, the answer is yes. Again, the instructions that were given to the

1 investigation team is that these patients had to truly be refractory to all revascularization  
2 options, both PCI and surgery. These patients had been evaluated for surgery multiple  
3 times in their lives and as you saw, a majority of patients had actually had prior surgery,  
4 some had had two surgeries and many of the patients had multiple PCIs. So while this study  
5 was done before the concept of a heart team was really formalized in the world, in general,  
6 yes, these patients were seen by a general cardiologist and a comprehensive care team also  
7 consisting of an interventionalist and a surgeon before they were felt to be eligible for  
8 enrollment in the trial.

9 DR. LANGE: Gregg, let me go -- let me ask a question. That's documented or it's just  
10 presumed that because they got in the trial, that's what happened?

11 DR. STONE: Yes, just presumed. There is not detailed documentation about it.

12 DR. LANGE: Thank you, sir.

13 There was a question about how the blinding was confirmed. And both the FDA and  
14 the Sponsor talked about what measures were taken to try to do it. Were there any  
15 measures taken to confirm the blinding was successful?

16 DR. STONE: Yeah, that's a great question. You know, this is one of the first sham-  
17 controlled randomized trials that have been done in interventional cardiology and in fact, I  
18 think, as you know, there's actually even to date been very few of them that have been  
19 done and the procedures that were put in place when this trial was performed are not as  
20 robust today, and they did not include a blinding or what I would even call a perception  
21 questionnaire.

22 So we are left to look back at the robustness of the study procedures, but I think that  
23 even more important than that, we do have evidence that there was effective blinding, and  
24 I'll just bring this up again because it is one thing to have somebody say yes, I know or I  
25 don't know, but I do think we actually have pathophysiologic evidence that there indeed

1 was blinding. I think you got these 30-day results that do show a combination of placebo  
2 effect and yes, I would certainly say there may be some regression to the mean  
3 phenomenon here. And importantly, I think you see that that was balanced between the  
4 Reducer group and the sham-control group. Then I think you see, over 3 months and 6  
5 months, you've got a separation of the groups of patients here and I think very important is  
6 actually the control group. Because they had the control group become unblinded, I don't  
7 think you would see the same maintenance of their improvement. And the same thing  
8 happened for just one-class improvement in angina, as well. This goes along with the  
9 mechanism of the Reducer, that that neointima develops over about a 2-month period,  
10 which is when you would expect the pressure gradient to be created.

11 DR. LANGE: And I appreciate that.

12 DR. STONE: Thank you.

13 DR. LANGE: But the answer is -- and I agree, the study was done differently, it's not  
14 as robust, but the answer to the question is we have no -- I mean, we can surmise based  
15 upon this data but the answer to, I guess, Dr. Allen's question is that we have no  
16 confirmation.

17 DR. STONE: There was no formal blinding questionnaire.

18 DR. LANGE: All right, thank you.

19 Dr. Zuckerman.

20 DR. ZUCKERMAN: Yes. Important data are being discussed and as I mentioned after  
21 the Panel's initial Sponsor presentation, the Sponsor does have the ability to show data that  
22 has not been independently reviewed by the FDA. But I would again ask the Sponsor, when  
23 they are showing data that has not been independently reviewed by the FDA, such as the  
24 last important slide discussion, that this be mentioned such that panelists are made aware  
25 and also that these comments are appropriately put into the official government record.

1 Thank you.

2 DR. LANGE: Thank you. And if we fail to do so, please remind us, Dr. Zuckerman.

3 Thank you, I appreciate that.

4 Dr. Allen also asked who was managing the post-op medications and any change in  
5 medications and was that person blinded to the treatment.

6 DR. STONE: Yes, all the post-op evaluations and management were performed by  
7 blinded investigators.

8 DR. LANGE: Terrific. And then there was a question about pressure measurements  
9 in animals. Dr. Banai, you showed one graph. How many animals actually had pressure  
10 measurements done?

11 DR. BANAI: In this graph that I showed, this is the mean of six animals to 6 months  
12 after implantation of the Reducer.

13 DR. LANGE: Terrific, terrific.

14 So Dr. Allen, does that address the questions that you had?

15 DR. ALLEN: Yes.

16 DR. LANGE: Great. So thank you for the responsiveness of the Sponsor, I appreciate  
17 it.

18 Dr. Yuh asked about microvascular disease in a broader general -- about the  
19 presence of coronary disease because none of that data was shown for either group, the  
20 extent of coronary disease, or there was microvascular disease and did this exclude women.  
21 So let me turn it over to you, Gregg, to address that.

22 DR. STONE: Yeah, thank you.

23 To get into the trial, the patients had to have severe coronary artery disease and  
24 there was no core lab analysis and there was no formal evaluation of the extent of coronary  
25 disease. The patients with microvascular angina or INOCA, ischemia with normal coronary

1 arteries, were not enrolled in this study and I can tell you, speaking for myself and  
2 Dr. Henry, we certainly believe the trial results speak to the patients that were enrolled,  
3 which were patients with severe coronary disease, again, almost all of whom had multiple  
4 coronary revascularization procedures, there's evidence of that. And while we think this  
5 should be effective for patients with microvascular disease, that would need to be the topic  
6 of a future study and I don't believe the Sponsor is asking for approval in that patient  
7 subgroup.

8 DR. LANGE: Great. Dr. Yuh had also asked about the procedure times for the sham  
9 procedures and the control -- I mean the sham procedures and the device procedures. Do  
10 we have that information?

11 DR. STONE: Yeah. As we sit here, I don't have that information, but these are  
12 relatively short procedures, usually on the order of 30 to 45 minutes. It's basically  
13 implanting a stent in a normal vein, which, for any interventionist with a little bit of  
14 training, is quite easy. And of course these patients, again, were mostly heavily sedated  
15 and I can't imagine that there was a major difference in the procedures, but I don't have  
16 that data.

17 DR. LANGE: All right, great.

18 So Dr. Yuh, to the extent possible, did the Sponsor address both of your questions?

19 DR. YUH: For the most part. The first part was more out of curiosity in terms of the  
20 putative mechanism of this intervention, it would seem to favor those patients that had  
21 microvascular disease with a redistribution of blood flow. So I was just curious as to why  
22 that population was not included at all. There must have been some conscious reason why  
23 you wouldn't include that and I was just curious as to what that was.

24 DR. STONE: No, I mean, I think it's a great question and I think those patients should  
25 respond well. I will show this one additional slide because you also asked about women. It



1 is true that women have a higher incidence of INOCA and the fact that those patients were  
2 excluded is probably one of the reasons why it was a male-predominant population.

3 So initially it was decided to study the Reducer in the population of patients with,  
4 again, obliterative coronary artery disease with no revascularization options, especially left  
5 coronary distribution ischemia, and that led to the population that we studied. You can see  
6 here regarding the number of women, and we've put here the major refractory angina  
7 studies. And again, all of these studies did require coronary artery disease, as well, and you  
8 can see it was 19% women in the Reducer and actually 12 to 16% women in the other  
9 studies. So this seems to be quite typical in terms of both race, sex, and age for refractory  
10 angina populations due to coronary disease.

11 DR. LANGE: Great.

12 DR. YUH: Thank you.

13 DR. LANGE: Dr. Cigarroa had asked about whether the individuals that were  
14 adjusting meds and those that were assessing angina, were they independent, were they  
15 blinded, were they interactive in any way?

16 DR. STONE: No, all of the investigators and personnel who were seeing the patients  
17 during follow-up, who were assessing their angina class, who were adjusting their  
18 medications, who were administering their exercise tests, were blinded.

19 DR. LANGE: Right, but were the people who changed the medications blinded to  
20 how they responded to the SAQ or the CCS?

21 DR. STONE: Oh, you mean were the same people that were taking the CCS, were  
22 they the ones who were adjusting medications?

23 DR. LANGE: I believe that's the question.

24 DR. STONE: So I don't know the answer to that, but I would imagine so, because you  
25 have to ask the patient about their angina before you decide whether to change their

1 medications over time.

2 DR. LANGE: All right, thank you. He had also asked about whether, in the 34 of 52  
3 treated patients that had no change, no improvement in their angina, and the 44 of 52 in  
4 the sham group, was there any difference or change in the diastolic filling period?

5 DR. STONE: Dr. Banai will address this.

6 DR. BANAI: Hi, Shmuel Banai here.

7 May I have Slide MO-11, please? The issue of diastolic function was -- is a question  
8 we ask our self because when you narrow the coronary sinus and you increase coronary  
9 sinus pressure, one might suspect that you can cause interstitial edema and affect diastolic  
10 dysfunction. This was not tested during the COSIRA trial, but it was evaluated later by two  
11 clinical trials.

12 The first one that I'm showing you here is an echo study which tested the diastolic  
13 function in 24 patients with severe grade III and IV angina, who were treated with a Reducer  
14 and diastolic function, was evaluated 6 months later and actually, it showed that the  
15 diastolic function improves at 6 months. And the idea why does it improve is because you  
16 improve perfusion and when you improve perfusion, you improve diastolic function.

17 Lately, it was also tested, and this is a slide that was represented in a manuscript  
18 that was published in May, so the FDA might have not the chance to see it. Nevertheless,  
19 this study was done by perfusion MRI and by perfusion MRI, they took 20 patients with  
20 severe angina, implanted the Reducer and repeated that perfusion MRI at 4 months and  
21 they showed that at 4 months the Reducer implantation improved myocardial contractility,  
22 ejection fraction improved, it improved myocardial longitudinal and circumferential strain  
23 and maybe most importantly with the fact of diastolic function, it is improved ischemic  
24 burden and improved intramural perfusion balance and redistribution of blood into the  
25 ischemic subendocardium, which goes along with animal data that we have previously.

1 Thank you.

2 DR. LANGE: Dr. Banai, thank you for sharing that and thank you for identifying that  
3 the FDA hadn't seen that yet. Thank you.

4 Dr. Cigarroa, does that address your questions or additional clarification to your  
5 questions?

6 DR. CIGARROA: So that addressed two out of the three. The one that I had  
7 specifically asked the Sponsors to clarify were between the responders and non-responders,  
8 was there a difference in the diastolic filling period?

9 DR. STONE: So as measured by heart rate, I presume.

10 DR. CIGARROA: That would be --

11 DR. STONE: Yeah.

12 DR. CIGARROA: -- one surrogate for it. Yes, sir.

13 DR. STONE: As we sit here, I don't have that information. The baseline heart rates  
14 were 60 to 65 beats per minute. At baseline, the beta blockers and ivabradine dosage  
15 didn't change over time, but we'll see if before the next break we can get you any of the  
16 follow-up heart rate information, specifically in the responders and non-responders.

17 DR. CIGARROA: And the reason I'd like clarification of that goes back to the impact  
18 of shorter diastolic filling periods on the ratio of epicardial to endocardial blood flow.

19 DR. STONE: No, I understand, that's a great question.

20 DR. LANGE: Thank you. Thank you for addressing that. Dr. Banai and Dr. Stone,  
21 thank you. Thank you, Joaquin.

22 Dr. Wittes had asked about what countries were involved, 11 countries, and I think  
23 you have that information now, Janet, and the current smoking rate I think we have as well,  
24 the FDA. You did ask about showing movement between treatment and controls and I don't  
25 know if the Sponsor would like to address that.

1 DR. STONE: Yes, we can. We can show you the individual patient data and from  
2 baseline to 6 months. If we can pull up those slides, I can show it to you either in the  
3 standard point by point, patient by patient, or in tabular form. I'm waiting for those slides  
4 to come up.

5 (Pause.)

6 DR. STONE: Okay, here we go. Now we've got it. Okay. So here you see point by  
7 point, and the Reducer patients on the left and the control patients on the right, the Class  
8 IV and Class III patients, at which they all started at and then Class 0, I, II, III, IV at 6 months.  
9 And again, I could -- if you want to see this in tabular form, I can --

10 DR. WITTES: Yeah, tabular form is much easier to see.

11 DR. STONE: Okay, let me bring that up.

12 DR. WITTES: Thanks.

13 DR. STONE: Okay, so here you can see the baseline, the Reducer on the left, control  
14 on the right, baseline Class III/IV are the two columns and then the follow-up 6-month  
15 classes are shown in the rows.

16 DR. WITTES: Okay, that doesn't really show it, but the combination of that and what  
17 the FDA showed, one can reconstruct it.

18 DR. LANGE: Okay, great. Great.

19 DR. CONNOR: Dr. Lange, can I ask a question here, a follow-up to that? This is  
20 Jason.

21 DR. LANGE: Yes, sir.

22 DR. STONE: And it's possible that the FDA did not see those slides before.

23 DR. LANGE: Okay. Thank you, Gregg.

24 Go ahead, Dr. Connor.

25 DR. CONNOR: Yeah, on the last slide I was worried because you kind of showed that

1 quickly -- if you can put the last slide up -- during the core presentation, that there were  
2 actually two control patients who went, I guess, from Class III all the way to no angina at 6  
3 months, which is amazing and great for those two patients. I was wondering if you had any  
4 of those in your observational study, I guess, the REDUCER trial, your second trial that you  
5 had a bigger  $n$  on it. Did you have anyone go from like a III to a zero in that trial, either?

6 DR. STONE: So go from a III to a zero. In the REDUCER-I there were only -- only four  
7 classes were captured, so the zeroes were put into Class I. So we certainly had many  
8 patients that went from Class III to Class 0/I, but I can't break it out from Class 0 and Class I  
9 in REDUCER.

10 DR. CONNOR: Okay. Yeah, that kind of went fast during the core and I was just  
11 wondering about that, given that the only two patients we have who had none whatsoever  
12 were controls. Okay, thank you.

13 DR. WITTES: And let me just say that it really is easier if you see the number of  
14 people who were in III who went into I, II, III, and IV rather than having these two cross-  
15 sectionals, but we can figure it out.

16 DR. LANGE: All right. And they may be able to get that during a break, Dr. Wittes.  
17 And I appreciate trying to read our minds to know exactly what we want, so if they can  
18 provide that information, that will be great.

19 And also, Dr. Wittes had asked, the last slides, 75, 77, 78, 81, that were shown, it  
20 talked about people that were -- results at 6 months, 1 year and 2 years, were those all the  
21 patients? In other words, were there --

22 DR. STONE: Yeah.

23 DR. LANGE: You reported on a 181. Were there only a 181 at 6 months or were  
24 there a larger group and only 181 reported?

25 DR. STONE: Yeah. No, we looked at that and I can show you those data. It's

1 approximately at each time period all the way out to 5 years. It was 90 to 95% of all the  
2 patients represented in the CCS angina class. So I can show you these data here, these are  
3 the -- I'm not sure that this has been seen, too, this has not been seen by FDA before, but  
4 here are -- this is the patient visit compliance, all the patients eligible and all the ones that  
5 actually had their visit, and then if you actually cross-correlate this to that actual slide that  
6 showed the number of data with CCS angina classes, you could see that it's basically 90 to  
7 95% at each time period.

8 DR. LANGE: So go back to that one slide just so that we can hit it for Janet.

9 DR. STONE: So here, I'll show you now the original slide. You know, so for example,  
10 if you look at -- what would you like? If you look at --

11 DR. LANGE: Just the slide you showed previous to this.

12 DR. STONE: Yeah, okay. Yeah, I was going to help you. It's easier if there's two  
13 slides up at once, but here's the slide we showed previously.

14 DR. LANGE: So at 6 months there were -- go back where you were, please.

15 DR. STONE: Right. So go back and forth between the two slides.

16 DR. LANGE: Right there, stay right there. There were 199, Janet, at 6 months. They  
17 have data on 181. There are 166 at 12 months, they're reporting data on 140. And there  
18 are 117 at 2 years and they're reporting data on 98.

19 DR. WITTES: Yeah. So I have to say that bothers me because by not showing what's  
20 missing, we interpret the data as if it's not -- the data are not missing. But thanks for  
21 showing.

22 DR. HIRSHFELD: So actually before we leave this, it looks like these percentages add  
23 up to more than a hundred percent.

24 (Audio feedback.)

25 DR. LANGE: The greater than one class improvement includes those with greater

1 than one and greater than two.

2 DR. HIRSHFELD: Okay. All right, got you.

3 DR. LANGE: All right. So Janet, does that answer the questions that you had? And  
4 we can talk about --

5 DR. WITTES: Yes.

6 DR. LANGE: -- the limitations, okay. Great, thank you.

7 Dr. Mathew had wanted to -- and as did many Panel members, wanted to address  
8 the fact that these were people that were -- had refractory chest pain but 25% were on zero  
9 or one medication.

10 DR. STONE: Right. So as I showed you, again, there's a lot of reasons why these  
11 patients who do truly have refractory angina are on only zero or one medication. It's not by  
12 choice, if you will, because they are really having severe angina and they are quite limited.  
13 But in general, it's because there were side effects from those medications, there was  
14 hypotension, dizziness, renal insufficiency, etc.. They're taking many other medications and  
15 as I showed you, I can show it again if you'd like, the angina, number of angina medications  
16 in the classes are really pretty consistent with what you've seen in every truly refractory  
17 angina trial before. I can put it up quickly if you'd like to see it again because it went fast  
18 the last time.

19 DR. LANGE: So Gregg, again, my question that I'll pose to you all, were these things  
20 that were documented or do we just assume that because a patient was only on one  
21 medication, that's all they could tolerate?

22 DR. STONE: Right. So there was not a detailed questionnaire in the case report form  
23 saying for each possible antianginal class, of which there are about seven of them, why was  
24 the patient on it or not on it, or only on that dose and what other side effects in the last 1,  
25 10, 20 years that they had from those medications, that data was not collected in COSIRA.

1 DR. LANGE: Okay. And the other question that Dr. Mathew had asked was what do  
2 we know about the adherence to medical therapy in both COSIRA and the REDUCER trial?

3 DR. STONE: Right. So as I showed you from the other prior two slides the other  
4 time, the classes of medications were consistent in both groups, at least as reported by the  
5 patients. There certainly were no pill counts or anything that were done in COSIRA, that's  
6 really been done in any study. But from patient-reported compliance of medication, I can  
7 put it up, and with the quick and dirty dose equivalence that was assessed, it also was  
8 consistent and similar in both groups. So this just shows you, again, whether the patients  
9 were taking them at baseline and at 6 months in both groups and you can see, in general, I  
10 would say it's pretty similar between the two groups. And this slide may not have been  
11 seen by FDA before.

12 DR. LANGE: Great. And so I guess what you're saying on that is we don't have any  
13 information regarding adherence or compliance.

14 DR. STONE: Right, I think within the limitation reported, right, we didn't have pill  
15 counts or bottle openings from electronic bottles or anything like that.

16 DR. LANGE: The last question Dr. Mathew asked, and we'll make sure he has -- he  
17 asked about that there were 8 of 50 devices that were technically -- there was some  
18 difficulty handling these with technical difficulties and I'll allow the Sponsor to address that.

19 DR. STONE: Well, I will let Dr. Bebeau address it, but I will say that as with any new  
20 procedure, there's always a learning curve. There's a learning curve for interventionalists  
21 working in the coronary sinus, working with this device, and I think as you saw as we went  
22 from COSIRA to REDUCER-I to the commercial experience, the outcomes in terms of  
23 coronary sinus events, etc., were going down in general. So I think some of this is due to  
24 the early learning curve, some are operator errors, some are training experience. As you  
25 can see, the procedures were quite uncomplicated, though, and the Reducer was implanted



1 in 50 of 52. The only times it wasn't implanted was because of poor guide catheter support  
2 and that can be improved by future guides.

3 Ms. Bebeau, I don't know if you want to further address that. You might be on mute.

4 MS. BEBEAU: So sorry, I'm trying to click on here and I apologize. Vicki Bebeau here.

5 For the technical success rate, as you spoke about, we had the 50 of 52 patients with  
6 successful implant. I would have to go to another slide as far as the mishandling, but  
7 typically in the protector that's on the device, sometimes in the handling, I think there was  
8 a case where it caught the -- the strut caught on some gauze. I mean, there's been just a  
9 learning curve with some of the technical successes. But the procedure successes, as we  
10 talked about, it was just the two patients that they had difficulty with the guide catheter.

11 DR. ALLEN: Can I ask a number-counting question because it doesn't make sense.  
12 I've heard you say 50 of 52 were successful, but you had a 3.9% embolization rate. So  
13 something, either you didn't count -- I'm not sure where your numbers are coming from.  
14 You had a 3.9% embolization rate, but yet you state you had 50 out of 52 were successfully  
15 implanted.

16 DR. STONE: So Dr. Banai was at the procedures and I think can address this.

17 DR. ALLEN: I just want to know the -- no, I don't want it addressed, I just want to  
18 know the numbers, so they don't add up.

19 DR. STONE: Well, they were ultimately -- the procedures were ultimately  
20 successfully implanted. So I mean, it's not uncommon at all that in interventional  
21 cardiology you try one thing, it doesn't work, you take the stent out, you try something else.  
22 Some of these --

23 DR. ALLEN: I guess I have a problem, though, with an embolization of a device being  
24 counted as a successful procedure.

25 DR. STONE: I understand. I think it all depends on the definition. In terms of the

1 definition -- you know, interestingly, for those non-interventional cardiologists, believe it or  
2 not, most embolizations don't cause a problem, either the stent can be retrieved or even  
3 when they're not retrieved, it's almost never caused an issue, which is interesting. It's not  
4 desirable and of course, we would try to avoid it and the rate of that has gone down over  
5 time with greater learning with the Reducer. That said, ultimately, the procedures were  
6 successful in 50 of 52, but I share your comment and your concern.

7 DR. BANAI: Shmuel Banai here.

8 They're called embolization when the Reducer moved either on the balloon  
9 backwards, dislodged or it was on the wire. In every case, the Reducer was taken out as it  
10 was on the wire and a new Reducer was implanted successfully. So there was no  
11 embolization to another organ.

12 DR. LANGE: So I'm sorry, I'm looking at this slide and it says all but one occurred  
13 during a procedure and one embolization was reported 3 hours post-procedure and it was a  
14 pulmonary artery. So did I miss that?

15 DR. ALLEN: You missed the same thing I missed.

16 DR. BANAI: This is in the REDUCER-I or COSIRA? Are we talking about COSIRA or  
17 REDUCER-I?

18 DR. LANGE: Yeah, COSIRA.

19 DR. STONE: I think if we put up -- we can put up this slide. Here, I think you can see  
20 COSIRA is at the bottom and there were these two cases where the device was dislodged  
21 from the balloon catheter but did not embolize. The pulmonary artery cases you see were  
22 in the commercial experience of several thousand cases.

23 DR. LANGE: Okay. Thank you, Gregg. Thanks for clarifying that. Okay.

24 Dr. Mathew, did they address the questions that you had?

25 DR. MATHEW: Yeah, I think within the limits of what's possible and I think that

1 addresses it, yes. Thanks, Gregg.

2 Just one other point that I brought up or asked about, maybe this is an  
3 administrative question, but the IFU, proposed IFU, as written, it said patients not eligible  
4 for CABG or PCI, and while it's clear that microvascular angina is not included in this subset  
5 of patients nor included for the approval, I guess maybe I'm being a little overly cautious,  
6 but could somebody interpret that proposed IFU as including somebody with microvascular  
7 angina who obviously wouldn't be a candidate for CABG or PCI? So it would just be a  
8 matter of clarifying the phraseology if it comes to that.

9 DR. STONE: Yeah, I can speak quickly about indication and then Ms. Bebeau can also  
10 speak about why it turns out that was the IFU. I think there's the full expectation that  
11 based on the advice from this Panel and in working with the FDA, that that proposed IFU  
12 and central indication would be changed. I'm a firm believer that the results of any trial  
13 apply to the patients that were enrolled, which were patients here with severe coronary  
14 artery disease who did not have coronary revascularization options with Class III or IV  
15 angina. That's what I would expect. Ms. Bebeau can speak to as to why that ended up as  
16 the initial submission.

17 MS. BEBEAU: Yeah, Vicki Bebeau here.

18 As I mentioned earlier in my presentation, we submitted that indication for use  
19 based on the fact that that was the indication we were provided for the breakthrough  
20 designation, fully understanding, as it was clear from the Agency and the guidance  
21 document, that it would likely require to have some further discussions with FDA. We are  
22 open to any changes that the Panel or FDA would recommend, if the device would be  
23 approved, to add to the IFU and amend it.

24 DR. MATHEW: Okay, thank you.

25 DR. LANGE: Dr. Ohman had some questions regarding the early termination of the

Free State Reporting, Inc.  
1378 Cape Saint Claire Road  
Annapolis, MD 21409  
(410) 974-0947

1 study and whether that introduces some bias and could you provide a sensitivity analysis of  
2 how that might influence the trial.

3 DR. STONE: Yes. Let me ask Chris to address that.

4 MR. MULLIN: Good afternoon. Chris Mullin, I'm a statistician with NAMSA. So the  
5 trial did have the pre-specified analysis plan for potential early stopping with an alpha  
6 spending function. The trial was not stopped at that point, so the decision was made to  
7 continue. Subsequently, though, there was the decision made to end enrollment, but that  
8 was based on logistical concerns, not safety concerns and not crossing the pre-specified  
9 boundary. So really, the effect was reduced power relative to the originally planned sample  
10 size.

11 DR. LANGE: Chris, regarding premature or early termination and the possibility of  
12 bias towards a positive result.

13 MR. MULLIN: Yeah. So because it's not stopped due to awareness of the efficacy  
14 data, I would argue that there's no issue with that. The pre-specified stopping would've  
15 handled that if that boundary had been crossed, but that wasn't the case.

16 DR. STONE: I can also address, importantly, that the Sponsor -- and Ms. Bebeau can  
17 verify this -- had no knowledge of the outcomes data when the decision was made to stop  
18 the trial. The data remained blinded and they just -- they stopped because the enrollment  
19 had slowed down at that point and they were getting better follow-up than anticipated, so  
20 they felt they would likely have enough data for a robust and hopefully a positive  
21 conclusion to the study, if it was going to be positive, so they elected to stop.

22 DR. LANGE: Dr. Connor.

23 DR. CONNOR: And I would just say, from the stats perspective, if it truly is just  
24 logistically you're not enrolling and you stop short, that you're not increasing Type 1 error  
25 because you're not stopping because you have a success. So in theory, it's not increasing

1 Type 1 error if they stop for that reason.

2 DR. WITTES: Can I ask a question? This is Janet.

3 DR. LANGE: Please, Janet. Go ahead.

4 DR. WITTES: And that's related to the slight difference between the question that's  
5 being addressed and the question that I asked. Or that I tried to ask. What was conveyed  
6 by the DMC to the Sponsor and did the Sponsor, when the Sponsor said we're going to stop  
7 the trial because we can't recruit, did they ask the DMC whether that was okay or did they  
8 just do it?

9 DR. STONE: I'll let Ms. Bebeau and Dr. Banai answer that.

10 MS. BEBEAU: Yes, Vicki Bebeau here.

11 The decision to stop the study was made by the board of directors and the board of  
12 directors had no knowledge of the analysis that was -- any of the data, they had no  
13 exposure to the data. It was a decision based off of the difficulty in recruiting and how long  
14 the recruitment had gone on, as well as the small number of lost to follow-up patients.

15 DR. WITTES: But did they ask the DMC for whether that was okay with the DMC?

16 MS. BEBEAU: I don't know the answer to that, so I honestly can't answer that  
17 statement.

18 DR. WITTES: All right, so let me just tell you why that's important, because if they  
19 didn't, then I'm fine, I agree totally with what everybody says, it has no effect on Type 1  
20 error rate. If they asked and got some kind of a signal back, then you wonder whether  
21 there was some bias.

22 DR. LANGE: Thank you, Doctor.

23 Dr. Ohman had also asked about a sensitivity analysis and I don't know, Dr. Mullin, if  
24 you had an opportunity to do such or thought it was worthwhile.

25 MR. MULLIN: Sure, thank you. Chris Mullin again.

1           So I think, if I recall the question, that was the one if the enrollment had continued  
2 to the full originally planned sample size, how sensitive the results would be and of course,  
3 since we did these calculations over lunch, FDA hasn't had a chance to verify them. But by  
4 our calculations it would be still significant, the p-value is about 0.02 if you just take the  
5 observed results and apply that proportionately to the full sample size. And then the  
6 sensitivity, if you switch two patients, two treatment patients, from successes to failures,  
7 that would flip the results so the p-value would be over 0.05. Likewise, two control patients  
8 going the other way or one in each group going the other way.

9           DR. LANGE: And then Dr. Ohman had also asked, if we took the COSIRA data and  
10 took out the individuals that were on minimal medications, did the results change, were  
11 they still affected?

12           DR. STONE: Yeah, and I can show you data in that regard. We can pull up that  
13 subgroup slide. While you're waiting, the answer is there was no interaction between the  
14 number of baseline medications and the primary endpoint results. Okay, yeah. Here we go,  
15 terrific. And I don't believe FDA has seen this slide.

16           So here you can see the primary endpoint on top of two or more CCS angina class  
17 improvement and the bottom half of the slide is one or more CCS angina class improvement  
18 and you see all patients, which is the first and fourth treatment differences in the 95%  
19 confidence interval, and then you see, for those patients who were on zero to one versus  
20 two or more antianginal medications, you can see the response rate and the 95%  
21 confidence intervals. The interaction p-values are 0.27 and 0.35, respectively.

22           I would point out, of course, this is a modest sized study and of course, all subgroups  
23 are at best hypothesis generating and there's certainly no reason to try to pull any -- too  
24 much out of these small subgroups. It is interesting, though, I will say, that the zero to one  
25 antianginal medication group, if anything, you could see the difference was a little bit more

1 positive, which I actually take as a good thing because these are the patients that just could  
2 not tolerate antianginal medications, who had great clinical need and they had the good  
3 response. But again, I would say there's consistency effect across the number of baseline  
4 medications.

5 DR. LANGE: Great.

6 DR. SOMBERG: This is John Somberg.

7 You can just take the other point of view, is that the less medications, the less  
8 treated they are, therefore they have more significant ischemia and therefore the process  
9 works. But if you took a really resistant population that was fully medicated, there would  
10 be little difference.

11 DR. LANGE: Dr. Ohman.

12 DR. STONE: Well, you could take that, but again, these patients really had severe  
13 coronary disease and even those patients on zero and one medications had a high number  
14 of prior PCIs, bypass surgeries, myocardial infarctions. I don't think there's any reason to  
15 believe that they just were fresh patients with Class III or IV angina.

16 DR. LANGE: Gregg, unless you're willing to show that data that shows the zero and  
17 one patient and the three and four, let's not make those comments. Unless we can put  
18 them side by side and show that's the case.

19 DR. STONE: Sure, we can try to get you that after the break.

20 DR. LANGE: Magnus, you had your hand raised.

21 DR. OHMAN: And I think the last slide that Dr. Stone showed is very telling. I mean,  
22 the therapeutic efficacy in the zero to one is twice as much, which makes me wonder how  
23 well the patient population is characterized. And obviously, if you took those out, out of  
24 the equation, this study will not be statistically significant. So basically, it's not reliable or  
25 the findings are not stable enough.

1 DR. LANGE: Okay, all right.

2 DR. STONE: If I can respond.

3 DR. LANGE: Gregg. Gregg, I'm sorry.

4 DR. STONE: Okay.

5 DR. LANGE: Yeah. I mean, with all due respect --

6 DR. STONE: No, no. That's fine, Richard. That's fine.

7 DR. LANGE: I just want to get through the questions and then the Panel will have a  
8 chance after the break --

9 DR. STONE: Absolutely.

10 DR. LANGE: -- and we will call the Sponsor and the FDA back for questions. And  
11 then Magnus also had a question on Slides 90 and 91 about the MI, just number of MIs, and  
12 I know that you all can address that.

13 DR. STONE: Yes, we can pull that up. Magnus is absolutely right, there were two  
14 myocardial infarctions. The other slide that Magnus was looking at was looking between 30  
15 days and 6 months, so there was one MI before 30 days and one MI after 30 days. So I can  
16 show you this slide, Magnus, is all the events. So this is within 6 months and now I'll show  
17 you the other slide that you were looking at before. This slide was after 30 days to 6  
18 months. We had another slide that showed all the outcomes in total and there were two  
19 myocardial infarctions, both non-STEMIs.

20 DR. OHMAN: Thank you.

21 DR. LANGE: Thank you, Gregg.

22 Dr. Brindis had asked about safety. The comment made is that there were over  
23 2,000 individuals that had the Reducer implanted and there is no evidence of coronary sinus  
24 thrombosis and I think he just asked for how that was determined.

25 DR. STONE: Yes. So there's been no reported instances of coronary sinus



1 thrombosis and as we'll show you, the CT data is also consistent with that, and I know there  
2 was a question about that from PROMISE on the CT data, so if you'd like it to be now or  
3 later, but Ms. Bebeau can address that issue.

4 DR. LANGE: That would be great because Dr. Brindis had actually asked about the  
5 presence of thrombosis.

6 DR. STONE: Perfect.

7 DR. LANGE: And I think it had to do with the CT report in the REDUCER trial where it  
8 looked like 16 of 37, or 43%, reported a thrombus was present by the CT angio core.

9 DR. STONE: So Ms. Bebeau.

10 MS. BEBEAU: Yes, thank you.

11 I will pull up a slide here. I apologize. This might be a slide that FDA hasn't seen  
12 before, but when we were working through our deficiencies, when filed for our IDE, we had  
13 met with FDA and talked about concerns that they had potentially with device  
14 patency/occlusion, emboli or thrombus, fracture/durability, and in the COSIRA study there  
15 was data on the 37 CT scans that were done.

16 We'd been approached by the physician from India to look at CT data with his  
17 patients, so we were speaking with FDA, we talked about having the same core lab that  
18 read the CTAs in COSIRA to read the seven patients' 12-year follow-up. In talking with FDA,  
19 they did agree with us that we could reanalyze. We had Dr. Wiseman from MedStar  
20 reanalyze the CT data from the COSIRA study, which he actually then did. Let me just share  
21 another slide here. And again, FDA may not have seen this slide, also.

22 So when Dr. Wiseman re-reviewed the data, again, it was consistent that there were  
23 no fractures but that previously, the thrombus that had been identified was really areas of  
24 low CT density and the device coded as thrombus and the CT core lab indicated that low  
25 density area seen could represent other diagnoses and not possible to differentiate

1 between thrombosis and fibrosis. But in the reanalysis, the thrombosis was defined as  
2 hypoattenuation.

3 Working with FDA, we also had put together -- and I'll show this slide. And the FDA  
4 has seen this during our discussions with them. This was done every single CT scan, both  
5 for the eight patients from India and the 12-year follow-up in addition to the 37 patients  
6 that were analyzed by the same CT core lab during the COSIRA study. And this is an  
7 example of a patient that had been previously reported with thrombus, and I'd like to turn  
8 this over to Dr. Banai to explain further.

9 DR. BANAI: Thank you. Shmuel Banai here.

10 This is a very nice example. If you again look at the bottom pictures, and the one  
11 most to the left and the second one from the left, this is a proximal aspect of the Reducer  
12 and at 4 o'clock you can see hypoattenuation, a black spot, which was previously suspected  
13 to be a thrombus, but that is typical for neointimal formation for tissue growth occupying  
14 the area and encroaching on the lumen. There was no thrombosis seen in none of the  
15 animal studies and none of the CT angios and none reported as adverse events in all the  
16 clinical cases we had in Europe. Maybe I'll elaborate on that a little bit.

17 The question was asked why don't we give anticoagulation like we do to patients  
18 with venous thrombosis, and the reason is that here, although it is a low-pressure system, a  
19 venous system with a venous pressure, because of the narrowing we have a high-flow  
20 system. The funnel shape of the coronary sinus after Reducer is implanted and the  
21 narrowing is established, the funnel shape causes acceleration of blood flow through the  
22 narrowing and because of the high flow, we don't see thrombosis and dual antiplatelet  
23 therapy is probably sufficient. Only in patients with chronic arterial fibrillation we give  
24 anticoagulation and single antiplatelet therapy.

25 DR. SOMBERG: Can I ask why you give --

Free State Reporting, Inc.  
1378 Cape Saint Claire Road  
Annapolis, MD 21409  
(410) 974-0947

1 DR. LANGE: Go ahead, John. Go ahead, John.

2 DR. SOMBERG: No, I just asked why they give dual antiplatelet therapy for the  
3 venous side.

4 DR. BANAI: We decided to give that because we -- exactly like we do in bare-metal  
5 stents in the coronaries, we did the same protocol and then only single antiplatelet therapy  
6 after 6 months.

7 DR. ALLEN: Can I -- yeah, it's totally different on the venous side, but you call this  
8 hypoattenuation, and I guess I'm very confused because we use the same term in structural  
9 heart and this, on a transcatheter valve, would be clot not fibrosis. And it looks like clot on  
10 a CT. And I'm confused because -- and I guess I would want to have clarity from the FDA  
11 and have them weigh in at this point in this process because it's really important we get --  
12 initially it's called thrombus and then you get somebody else to review it who then says no,  
13 it's not thrombus, and a hundred percent of them it goes away. I guess it's pretty confusing  
14 to me.

15 DR. LANGE: And I guess if you go back -- I think probably the best they can do is call  
16 it hypodense and they really can't tell what it is.

17 DR. BANAI: I can remark on this. When you do CT angiography to patients with  
18 stent in the coronary arteries, this is exactly the picture you get for diagnosis of restenosis.  
19 This is the way you diagnose restenosis in a stent, it's exactly the same picture, and that's  
20 what you see here. It's tissue proliferating into the lumen of the artery.

21 DR. LANGE: So Dr. Brindis, you had had some questions. And I'll get to you,  
22 Dr. Batchelor.

23 Dr. Brindis, does this address your questions?

24 DR. BRINDIS: It does.

25 DR. LANGE: And we'll have conversations about it. The difference, I just want the

1 Sponsor to answer and then we'll get involved in conversations about the value, the  
2 limitations, what we believe and don't believe.

3 Dr. Batchelor, you had your hand up, sir.

4 DR. BATCHELOR: Just a quick question. On paper the investigators adjudicating the  
5 CCS class and the outcome were supposed to be blinded, but if only the patients who  
6 received the device are on aspirin and Plavix, how were they blinded? How did you account  
7 for that?

8 DR. LANGE: They actually all got aspirin at once.

9 DR. STONE: All patients received dual antiplatelet therapy.

10 DR. BATCHELOR: Okay, thank you.

11 DR. LANGE: Actually, Dr. Cigarroa.

12 DR. CIGARROA: This is Joaquin Cigarroa. Thank you, Dr. Lange.

13 I just wanted a point of clarification on a statement that Dr. Banai made when  
14 talking about the antiplatelet versus antithrombotic. The comment regarding a high-flow  
15 state induced by the hourglass and the endothelialization is theoretical. There are no  
16 measurements where you actually measure velocity with any Doppler techniques or actual  
17 flow, correct?

18 DR. BANAI: Correct.

19 DR. CIGARROA: And at the narrowest point, there is -- that's the area where there is  
20 lack of complete endothelialization, correct?

21 DR. BANAI: Correct.

22 DR. CIGARROA: And so we actually don't know, in the narrowest point where you  
23 have endothelialization on the inlet and the outlet, whether or not there are eddies there  
24 and errors of spaces in the context of that incomplete endothelialization and what the  
25 velocity may actually be in that segment, correct?

1 DR. BANAI: Correct. Nevertheless we know that in all CT angiography, we see blood  
2 flow through the Reducer.

3 DR. CIGARROA: Understood. Thank you very much.

4 DR. LANGE: Thank you.

5 Dr. Connor, again, you had asked why they stopped the trial prematurely and they  
6 had described logistical issues and if you have a follow-up question to that and also --

7 DR. CONNOR: No.

8 DR. LANGE: Okay. And then he had asked for the outcomes by site and I don't know  
9 if the Sponsor had an opportunity to --

10 DR. STONE: I don't know. Chris Mullin, I don't know if you have that. We may not  
11 have that right now. We might be able to get back to you after the next break on that one.

12 DR. LANGE: Thank you, Dr. Stone, appreciate that.

13 Dr. Somberg asked about antiplatelet agent and anticoagulation and that's been  
14 addressed.

15 Dr. Batchelor had asked about regression to the mean and again, if a patient was  
16 assigned to one group or another. I think, Dr. Batchelor, they had addressed that. Any  
17 other clarifying questions on that? Are you comfortable with the response?

18 DR. BATCHELOR: No, Dr. Mullin nicely addressed to the effects of having one or two  
19 patients misclassified on the statistical significance, so thank you.

20 DR. LANGE: Great, thank you.

21 And then Dr. Batchelor also asked about the use of coronary sinus for other  
22 procedures in the future, like a bi-v lead and the recommendations from the Sponsor  
23 regarding that.

24 DR. STONE: Yes, I can address that. One of the nice things about the Reducer is that  
25 it's actually very easy to re-expand that central hourglass. All one has to do at any time is

1 just put a regular size balloon back into the Reducer and you can re-expand. Here, I'll just  
2 show you this picture. But that can easily be done. So the coronary sinus, you can have a  
3 future lead placed in the coronary sinus through the Reducer quite easily in the future.

4 DR. BANAI: Shmuel Banai here.

5 I want to comment on this. The reason we made the Reducer from stainless steel  
6 was exactly that, because you can always get rid of the narrowing. But since we have some  
7 experience now, I can tell you that it is very easy to implant a CRT-D lead into the coronary  
8 sinus through the Reducer without disrupting the Reducer.

9 DR. LANGE: Great.

10 DR. BANAI: But you can do it if needed.

11 DR. LANGE: Thank you, Dr. Banai.

12 Dr. Starling had asked, in humans, has there been any pressure measurements, any  
13 measurement of CSO<sub>2</sub> content or resistance, coronary vascular resistance.

14 DR. STONE: There has not been -- I mean, I'll just clarify that to do that, we expect  
15 the Reducer to again have its maximum amount of neointimal hyperplasia at around 2  
16 months. So it would require bringing patients back for an additional cardiac catheterization  
17 procedure and that has not been done, it's not been felt to be ethical.

18 DR. LANGE: Okay, great. And he had also asked about one of the exclusion criteria  
19 was an RA pressure greater than 15 mm and the explanation or rationale.

20 DR. STONE: I'll let Dr. Banai address that.

21 DR. BANAI: Yeah, this was a decision made when we were not sure that the device  
22 would work or not and we thought that the patient with very high right atrial pressure,  
23 mean pressure of 15 in the right atrium in the heart, they might also probably have high  
24 pressure in the coronary sinus that might -- you know, it might not be effective. That's why  
25 we decided to not include this patient, but we don't have any data to support that.

1 DR. LANGE: Just a question, Dr. Banai, in terms of physiology. Obviously, the  
2 Reducer is meant to increase coronary sinus pressure to improve the endocardial blood  
3 flow and obviously, an RA pressure of 15 increases CS pressure. Why that doesn't change  
4 blood flow and why it doesn't -- why that's not anti-ischemic.

5 DR. BANAI: It does, but patient with very high right arterial or right-sided pressure,  
6 they have probably other problems. They have right-sided failure, so the problem of angina  
7 is one of their problems. We, in the COSIRA, tried to take pure patients with obstructive  
8 coronary artery disease, preserved left ventricular function with no severe valvular disease  
9 or other cardiac disease that can mask the effect, the pure effect of the Reducer on angina  
10 and ischemia.

11 DR. LANGE: Okay, thank you.

12 And then lastly, Dr. Starling asked about the durability and any information  
13 regarding movement of patients over time in the COSIRA or other trials.

14 DR. STONE: So we don't have any data from COSIRA the past 6 months. We do have  
15 the data that we've shown you from REDUCER-I up to 5 years and from individual  
16 investigator sponsored studies, about four or five of them, in which the overall angina class  
17 appears to be approximately stable. But I agree, there are some limitations with  
18 assessment in a certain proportion of patients in all of those studies. But I don't think that  
19 we have any data that suggests that the effect is not long lasting.

20 DR. LANGE: Dr. Starling, does that address your questions sufficiently?

21 DR. STARLING: Just one last clarification. So a patient with an RA pacemaker lead or  
22 a right ventricular lead, are they eligible for the device? And currently, if the patient has a  
23 CRT device, is it not feasible?

24 DR. STONE: Go ahead.

25 DR. BANAI: Yeah. Pacemaker lead in the right atrium, in the right ventricle, are fine.

1 We allow the device as long as the pacemaker was implanted more than 3 months before,  
2 so there would not be a risk of moving mainly the atrial lead.

3           Regarding CRT-D leads in the coronary sinus, we don't allow it for two reasons. One,  
4 these patients usually have low ejection fraction, they have heart failure, and we don't  
5 know that this device is effective in heart failure patients with very low ejection fraction.  
6 The other thing, we're not sure that the Reducer implanted on a lead will not accidentally  
7 break the lead. That's why until now we did not allow implantation of the Reducer in  
8 patients with a coronary sinus pacemaker lead.

9           DR. STARLING: Thank you.

10           Thank you, Dr. Lange, that answers my questions.

11           DR. LANGE: And then Dr. Gersh had asked about is there any correlation of the  
12 clinical results with any of the -- any mechanisms, that is, associating coronary sinus  
13 reduction with improvement in anginal class.

14           DR. STONE: Well, I guess the answer would be no, we don't have any such evidence  
15 because as you saw from the CTs, the hourglass shape is formed in every single patient,  
16 we've seen that in every single CT. So certainly you saw about 70% of the patients who get  
17 the Reducer do have some reduction in angina. I don't think there's been a detailed  
18 analysis done to try to understand who the non-responders are, but we don't have detailed  
19 follow-up coronary sinus data in those patients because of the need for a repeat cath for  
20 human anatomic measurements, etc. So I don't think we can address that.

21           DR. LANGE: Okay, great.

22           DR. BANAI: Shmuel Banai here.

23           There are data in animals showing relationship between coronary sinus pressure and  
24 regional myocardial blood flow in animals. Not in human.

25           DR. LANGE: Okay. Thank you, Dr. Banai.



1 Dr. Allen.

2 DR. STONE: Thank you. I assume that's --

3 DR. ALLEN: Dr. Banai, can you help me understand, as a vascular surgeon, if I put an  
4 AV fistula in and I narrow a vein, that vein dilates. But in every one of your CTAs, and you  
5 talk about this hourglass, I don't ever see dilated coronary sinus proximal to where the  
6 higher pressure would be. Can you explain why not?

7 DR. BANAI: Shmuel Banai here.

8 I can speculate. First of all, we do not cause edema of the heart, we do not cause  
9 very high pressure gradient going backwards. We are doing a small pressure gradient  
10 elevation just enough to dilate some arterioles in the subendocardium to reduce resistance  
11 to flow. There are reports that even if you completely occlude the coronary sinus, patients  
12 do not die because of possible alternative drainage, venous drainage of the heart to the  
13 right ventricle into the right atrium. But --

14 DR. ALLEN: And I appreciate the explanation, but it's not really answering my  
15 question. Because your CTAs are very impressive, if you believe the CTA, you get a very  
16 impressive narrowing of the coronary sinus and to not show any dilatation of that coronary  
17 sinus before that seems unusual, which makes me suspect of the mechanism.

18 DR. BANAI: Yeah, but that's a fact that in all 44 CT angios, these are the pictures we  
19 see.

20 DR. LANGE: Thank you.

21 And then Dr. Gersh's last question was the discrepancy between what appeared to  
22 be an improvement in CCS class, but people that got sham treatment had the same  
23 satisfaction, angina -- and discordant results.

24 DR. STONE: Yeah, I can address that. The treatment satisfaction question is very  
25 different than all the other questions in the Seattle Angina Questionnaire. It basically is

1 asking are you satisfied with your doctors taking care of you, are you satisfied that they're  
2 trying to do everything that they can, are you satisfied that they are explaining to you what  
3 the situation is. It's not a direct question, are you satisfied that you have no angina.

4 And as a result, when you look at the baseline treatment satisfaction scores, they're  
5 around 80. This a 0 to 100 scale where a hundred is perfect, and the baseline treatment  
6 satisfaction scores are usually about 80 to 85 and that's in every refractory angina trial,  
7 whereas the baseline quality of angina or amount of angina scores are usually around 40 or  
8 45. So there's very little room to improve your treatment satisfaction in either group  
9 because they believe their doctors are doing a good job.

10 DR. GERSH: Thank you, that answers it.

11 DR. LANGE: Dr. Cigarroa had asked, the FDA showed a slide, 2020 -- or Slide 22  
12 looking -- discussing LV contractility, perfusion, and the FDA made a comment that they  
13 were unable to confirm the results and he had inquired as to why that was.

14 DR. STONE: I would leave that to Dr. Banai or Ms. Bebeau.

15 DR. BANAI: Shmuel Banai.

16 Are we talking about animals or humans?

17 DR. LANGE: Animals.

18 DR. BANAI: Okay, so the animal study --

19 (Crosstalk.)

20 DR. LANGE: -- FDA slide, FDA Slide 22.

21 DR. BANAI: Yeah. May I have Slide AD-20 for the animal studies?

22 DR. LANGE: No, let me show you Slide 22 in just a second so you know what he's  
23 talking about, Dr. Banai.

24 (Pause.)

25 DR. LANGE: Jim, can you show FDA Slide 22?

1 MR. SWINK: Yeah. Sam, can you go ahead and share your screen? Thank you.

2 DR. LANGE: Yeah. So Dr. Banai, Dr. Cigarroa asked about the last statement.

3 DR. BANAI: Yeah. Okay, the animal study that we did with the Reducer -- can you  
4 see my slide now? Can you share my screen, please? Okay. The animal study that was  
5 done with the Reducer are ischemic pig model, pigs were implanted with an ameroid  
6 constrictor around the proximal circumflex artery and the ameroid constrictor caused  
7 gradual narrowing until complete occlusion of the circumflex artery.

8 This is a very difficult animal model because the pigs, unlike the dogs, do not grow  
9 collaterals and most of the pigs will either die from this procedure or develop a full-blown  
10 myocardial infarction and a scar. Those pigs who survive the procedure and did not form a  
11 scar had reversible ischemia induced by dobutamine stress and perfusion of contrast, echo  
12 contrast, into the root of the aorta. So there were four pigs in each group and you can see  
13 here -- and then after ischemia was confirmed by dobutamine echo and perfusion, the pigs  
14 were divided into two groups, one group received the treatment with a Reducer and one  
15 did not, and they were reevaluated after 6 weeks and 6 months.

16 What you can see here on the left are examples of two pigs under dobutamine stress  
17 with perfusion and you can see the ischemia territory; actually, black area of the lateral wall  
18 of the pigs and on the right, same pigs under same dobutamine stress, even a slightly higher  
19 heart rate. Six weeks after Reducer implantation there is marked reduction in the ischemic  
20 territory, there's actually perfusion into this area of the circumflex artery. More than that,  
21 we didn't measure contractility here, so the pressure wasn't about contractility, we showed  
22 improved perfusion. Three out of the four control pigs died of pulmonary edema and  
23 ventricular fibrillation and none of the Reducer pigs died at 6 months.

24 Now, we do have data about improved contractility, systolic performance in humans,  
25 both by dobutamine echo wall motion score index and by MRI, the MRI data I showed you,

1 and also data -- if I can show you, if you want, about the improvement in dobutamine wall  
2 motion score index in humans.

3 DR. LANGE: I think we've seen that.

4 So Dr. Cigarroa, does this answer your question?

5 (Off microphone response.)

6 DR. LANGE: It does. Okay, great.

7 Dr. Somberg had asked if there were missing data with regard to the Seattle  
8 questionnaire and exercise stress test and dobutamine stress echo. He asked whether  
9 there were any missing data regarding the CCS data.

10 DR. STONE: So in COSIRA, we have CCS data in every patient.

11 DR. LANGE: Okay, great.

12 Dr. Somberg, does that address your question?

13 DR. SOMBERG: Yes.

14 DR. LANGE: Great. And, Gregg, I think you had addressed this. Dr. Starling had  
15 asked about who the heart team was that evaluated these patients to determine whether  
16 they were refractory or not and I think what you had said is there really wasn't a heart  
17 team, they just -- they were included in the study because of the refractory angina, is that  
18 correct?

19 DR. STONE: Correct. This was before the formal concept of a heart team, but all  
20 these patients, again, were cared for by general cardiologists and had been seeing  
21 interventionalists and surgeons for possible relief, to see if they had a revascularization  
22 option and they didn't, which is how they qualified.

23 DR. LANGE: Okay. Dr. Yuh.

24 DR. YUH: Thank you. David Yuh.

25 Along those same lines, were there any guidelines distributed to these evaluation

1 teams or individual surgeons versus interventionalists, as to what they should be using as  
2 criteria for deeming somebody as a "no option"? And the reason I'm asking is that even  
3 amongst cardiac surgeons looking at an angiogram, there may be differences in opinion  
4 between cardiac surgeons as to what's revascularizable versus not, and I'm just worried that  
5 with such small sample size that any variation in selection mechanism could confound the  
6 data. That's what kind of is nagging me about this -- about the definition you have of "no  
7 option" patients.

8 DR. STONE: Well, these were all very experienced sites that the Sponsor went to. I  
9 mean, that was a criteria for them being in the trial and these are sites that were very used  
10 to taking care of patients with complex coronary artery disease. I think, again, you could  
11 see by the fact that there was bypass surgery in a large percentage of the patients already,  
12 PCIs in a large percentage of the patients already, it's not like these patients were being  
13 ignored in terms of revascularization options.

14 That being said, it is a real-world study and you're absolutely right, there are  
15 differences in terms of threshold for revascularization among interventionalists as well as  
16 surgeons, and it's possible that there may have been more or less aggressive  
17 interventionalists or surgeons at different sites. Nonetheless, this is a multicenter trial and I  
18 think the results do reflect a multicenter experience with revascularization. I don't think  
19 surgical revascularization has changed much over the last decade in terms of what you can  
20 and cannot revascularize. There may be a slight improvement over the last decade for  
21 some interventionalists with advanced CTO techniques but, for the most part, this is  
22 coronary revascularization as reflected in these patients.

23 DR. YUH: David Yuh.

24 Just a follow-up. So as an example, in looking at an angiogram of a potential patient  
25 that had bypass surgery, for example, were grafts evaluated? And different surgeons, I

1 think, would take -- may take a different approach as to whether or not that patient would  
2 likely succeed with a re-operative coronary bypass. Even though the notion of coronary  
3 bypass surgery hasn't changed much there's still a lot of variation, I think, between  
4 surgeons as to what they deem as doable versus not or too high risk.

5 DR. STONE: No. I mean, I can certainly say that in general, yes, all of these patients  
6 had repeat angiograms to see that there were no revascularization options. And I will say,  
7 one thing I've learned from prior studies is that while we do sometimes evaluate these  
8 patients in eligibility committees, we've learned it's never a good idea for a remote  
9 committee to tell a surgeon or an interventionalist you need to operate on that patient or  
10 you need to fix that.

11 So as you know, surgeons like to operate and interventionalists like to stent, and if  
12 there was anything that they could do where they thought they might help the patient, my  
13 presupposition would be that they would have offered that to the patient before putting  
14 them in a trial such as this, especially given the number of procedures they had previously  
15 had.

16 DR. YUH: But don't you think that there would be more consistency if you had, for  
17 example, a core group evaluating, not necessarily dictating --

18 DR. ZUCKERMAN: Dr. Yuh, I don't mean to completely cut you off, but --

19 DR. YUH: Sure.

20 DR. ZUCKERMAN: -- our goal here is just to get some responses from the Sponsor  
21 such that we can then allow you and the other panelists to debate these important topics.  
22 Your point is extremely critical and an additional point is that this study was completely  
23 done outside the United States, so the Sponsor and FDA will be looking for you and other  
24 expert panelists to further develop this topic among yourselves, but if we can continue with  
25 just quickly addressing these questions to the Sponsor.

1 DR. LANGE: Great. So I think I tried to cover everybody's questions. Was there a  
2 question previously asked that hasn't been answered? If not, I want to ask Jacqueline  
3 Alikhaani, who's got a question. So is there any question that the Sponsor did not address?  
4 I see two hands, Dr. Bonde, Dr. Batchelor. And then Jackie, I'll come to you.

5 So Dr. Bonde. And I'll ask you to unmute.

6 DR. BONDE: Hi, Pramod Bonde here.

7 One of the facts that we learned in cardiac surgery is that if you perfuse retrogradely  
8 a coronary sinus, you tend to protect the left side better than the right side, and there have  
9 been instances and studies published of poor protection of the right heart because  
10 perfusion of the right ventricle is not as optimal as that of the left ventricle. In your studies,  
11 were you able to see a difference in the function of the right ventricle versus the left  
12 ventricle? That's one question.

13 And the second question is that how will that impact patients who have a single  
14 territory non-revascularizable disease, for example, a right coronary artery? Those are my  
15 two questions.

16 DR. BANAI: Shmuel Banai here.

17 Thank you, this is a great question and the reason, the thought was that the Reducer  
18 will not affect patient with ischemia or angina due to the right coronary artery was the  
19 anatomy. The vein that drains the internal wall in the right coronary artery territory drains  
20 right at the ostium and when you do retrograde cardioplegia or when we implant the  
21 Reducer, it's further deep in by 1 or 2 cm, it's deeper, so you do not increase pressure in the  
22 vein of draining right coronary artery. That's why in the COSIRA and at the beginning of our  
23 experience we said this is good only for patient with left-sided ischemia. Having said that,  
24 recently there was data published from Europe, and I don't know if FDA ever seen this or  
25 not, but this is a study that was done on patient with CTO and they --

1 DR. LANGE: Dr. Banai, is this a Sponsor study? If not, we're going to move on.

2 DR. BANAI: No.

3 DR. LANGE: Okay.

4 DR. BANAI: No, it's not Sponsor study. It's a study coming from Europe, from several  
5 centers in Italy, the Netherlands, and Belgium, I think, and if you look at the left side, they  
6 tried to evaluate the effect of the Reducer on patients with CTO, that they couldn't open  
7 the CTO and evaluate the effect of the Reducer. And surprisingly enough, patient with  
8 isolated RCA occlusion responded to the Reducer.

9 DR. LANGE: Thank you, Dr. Banai. Thank you. We're here to discuss the Sponsor  
10 and the Sponsor's studies and the results, so we're going to move on rather than talking  
11 about other studies that have been published.

12 So Dr. Batchelor, you had a question. You're on mute, sir.

13 DR. BATCHELOR: My question was related to the RCA, as well, so we're covered.  
14 Thank you.

15 DR. LANGE: Got it. Okay, I think we addressed it.

16 Ms. Alikhaani, I want to give you an opportunity to ask any questions and/or make  
17 comments. And if you'll please unmute.

18 MS. ALIKHAANI: Yes. I did not hear much about the patient voice, particularly  
19 American patients, throughout the study process and I think it is very important, very useful  
20 to have a team of patients and family members and caregivers and nurses to help advise  
21 throughout the conduct of the study. And I just can't -- you know, also, I just can't really  
22 understand why we can't just do the study here in the United States, where we're going to  
23 be using these therapies to help patients who are American patients here in this country.  
24 And as a heart patient, I'm a heart patient, healthcare consumer, volunteer patient  
25 advocate for heart health, and my biggest concern is always safety and effectiveness for



1 heart patients. I know what it's like to live with chest pain. As a woman and as an African  
2 American, I know what it's like to experience disparities in care and I just feel like -- I mean,  
3 there's so many questions and issues that are so legitimate that have been raised here  
4 today during this discussion and I just resonate with all of them. Very, very good discussion.

5 And it seems to me that a lot of these issues should've been taken care of before the  
6 hearing. I think that there's so many questions that bother me, like incomplete data and all  
7 the discussion that we've had trying to flesh all these issues out, incomplete data, missing  
8 data, early termination issues, exclusion of patients, lack of significant diversity in the study.  
9 These are all really critical issues.

10 And it seems like there was, you know, just such a missed opportunity to really do  
11 some good follow-up collaboration with the FDA that wasn't done because I'm seeing that  
12 here in the hearing we're having information presented that the FDA has never even seen  
13 before. And so that's kind of bothersome for me and I just wish that we could've done a  
14 study right here in the United States. I mean, we have such a diverse community of  
15 patients and consumers here in the U.S. and we have women experiencing disparities in  
16 care, we have African Americans, other ethnic communities, traditionally underserved  
17 groups. And so these are some of the concerns I have.

18 DR. LANGE: That's a great, great segue, I think, to go into the Panel deliberations.  
19 Jackie, I appreciate your comments and I think those will be prescient as we head into  
20 these.

21 Before we do that, any last remaining questions, clarifying questions for the  
22 Sponsor? Dr. Ohman. After that we'll go into our deliberations.

23 Dr. Ohman.

24 DR. OHMAN: Yeah, I think I did ask on the FDA side of the equation, but I had a  
25 question about how you reconcile the SAQ questionnaire of angina frequency and the

1 response in Canadian classification improvement that we're seeing. So I'm trying to  
2 reconcile a neutral effect versus a beneficial effect.

3 DR. STONE: So thanks, Magnus.

4 So in general, first, I think that they were quite consistent. You know, the Seattle  
5 Angina Questionnaire, first of all, it's not an optimal questionnaire for refractory angina  
6 patients in particular because of the wide standard deviation and given the standard  
7 deviation, we would've needed about 200 to 230 patients depending on the scale to be  
8 powered and as you know, there was only 104 randomized patients.

9 That being said, in addition to the two CCS angina class reduction, there also was a  
10 statistically significant reduction in angina quality and there were relatively substantial  
11 directional differences in angina frequency and angina stability, suggesting that that was  
12 Type 2 error.

13 So I actually think that it was relatively consistent, and had more patients been  
14 enrolled it would likely have been positive. But I do think that to remove some of that  
15 uncertainty, that is one of the reasons that Dr. Henry and I have been very vocal that we  
16 should perform an additional post-approval randomized trial, allowing this therapy to be  
17 brought to the public because it is such a major unmet clinical need, but at the same time  
18 we will generate more data regarding safety, effectiveness, and mechanism of action.

19 DR. LANGE: So with that, we are about, I'd say, 10 minutes ahead of where we still  
20 got a lot of work to do. I think this is a natural breaking point before we go into our  
21 deliberations. First of all, I want to thank the Sponsor. We asked a number of questions,  
22 some easy, some difficult, but there were a number of them and in a short period of time. I  
23 think you guys have done a terrific job of addressing those, so I want to thank you for that.  
24 And for the FDA, as well.

25 I'm going to recommend we take a break of 10 minutes and then let's get back, let's

1 begin our Panel deliberations. And probably the best way to do that is to begin working  
2 through the questions.

3 So Dr. Zuckerman, your thoughts.

4 DR. ZUCKERMAN: Yes, that's an excellent plan, Dr. Lange, but prior to that you'll  
5 probably want to hear the FDA response to questions which can be posed by the Panel.  
6 That can be done very quickly when we get back from the break. Dr. Lydia Glaw, Assistant  
7 Director of the Coronary Interventional Team, will be the quarterback to help you out on  
8 that.

9 DR. LANGE: Great. And Dr. Zuckerman, I have no preference. Would you prefer to  
10 do those now and then take a break or take a break and then do those?

11 DR. ZUCKERMAN: It depends really on you and the Panel members. It can be done  
12 either way.

13 DR. LANGE: Let's hear the FDA response. Let's do that now and then let's take a  
14 break.

15 So Lydia, let me hand it over to you. And there were some questions specifically  
16 directed toward the FDA. Do you want to address those or do you want me to --

17 (Audio feedback.)

18 DR. GLAW: Yes. Thank you, Dr. Lange. I think I have -- here and I can go through all  
19 the appropriate --

20 DR. LANGE: So Lydia, hold on a second. I think there's some audio -- it sounds like  
21 she's on Halloween voice. Lydia, go ahead and talk again and see if we've got it corrected.

22 (Audio feedback.)

23 DR. LANGE: No.

24 DR. GLAW: Okay.

25 MR. SWINK: Lydia, maybe stop your video, mute your mike and then kind of

1 re-launch it again.

2 DR. GLAW: Hello? Is that any better? Still bad.

3 (Audio feedback.)

4 DR. GLAW: Okay, I'm going to turn it over to someone else.

5 DR. LANGE: I'm sorry, Lydia, I couldn't hear what you said.

6 (Audio feedback.)

7 DR. GLAW: I'm going to see if someone else can take it on --

8 DR. LANGE: Okay.

9 MR. SWINK: We have Sam.

10 DR. RABEN: Yes, I'm happy to do that. Lydia, if you want to try and call back in, we  
11 can see if we can address that but I can handle that in the meantime.

12 Dr. Lange, can you give me the first question and then I can help facilitate the  
13 address?

14 DR. LANGE: Sure. Dr. Allen asked how does FDA handle missing or recommend  
15 handling missing data.

16 DR. RABEN: For that, Dr. Rona Tang will answer that question.

17 DR. TANG: Sure. This is Dr. Rona Tang from the FDA.

18 So for the missing data question, it is challenging because there are so many, there's  
19 a large number of missing data and the missing pattern for the secondary endpoint is not  
20 clear. So actually, we look to the Panel statisticians for recommendations on this issue.

21 DR. LANGE: Okay. And then Dr. Connor had asked whether we have -- because of  
22 the regression to the mean, do we have a history of information regarding these patients or  
23 similar patients?

24 DR. RABEN: Dr. Ryan can answer this question.

25 DR. RYAN: We don't have, with respect to -- are you asking about the placebo effect

1 or a statistical regression to the mean?

2 DR. LANGE: Dr. Connor?

3 DR. CONNOR: Well, I think the question was everyone was throwing around placebo  
4 effect and I'm calling it regression to the mean.

5 DR. RYAN: Okay.

6 DR. CONNOR: So I think the question is any historical evidence that says it's more  
7 one than the other. I would bet regression to the mean, but I always do.

8 DR. RYAN: I got you, okay. So yeah, we did do a literature search about the data  
9 that's available, specifically with respect to the placebo effect and angina. And I mean,  
10 certainly, the placebo effect can vary. We found out that there's a lot of range. Particularly  
11 with patients that are in desperate states and stress levels that are really high, the placebo  
12 effect can definitely be very high. As I mentioned before, there's been various studies that  
13 have shown that the placebo effect can vary in these -- in TMR studies, particularly, was an  
14 example that I brought up. So I don't think we're able to give you a number, we don't have  
15 a threshold that we would expect, for example, this percentage of patients we'd expect the  
16 placebo effect. So we can't really give you an upper or lower limit on that.

17 DR. CONNOR: Okay, thank you.

18 DR. LANGE: Dr. Wittes had asked what did the FDA not approve going forward.  
19 Janet, do you want to expand on that? Yeah, it looks like she's frozen.

20 DR. RYAN: I think I understand -- oh, go ahead. I think Dr. Wittes was asking about  
21 what our concerns were about the COSIRA study and how that was reflected going forward.  
22 So I think, basically, for the COSIRA-II study, and I think the protocol was included or  
23 information was included in the panel pack, one of the pushes we had about COSIRA was  
24 the use of a subjective endpoint rather than objective endpoint, and the COSIRA-II study  
25 protocol would have addressed that issue. So there were other issues, as well, but without

1 getting into a lot of specifics, that was the major concern that we had.

2 DR. LANGE: Okay, great. A couple more chat questions, again, not to the FDA but to  
3 the Sponsor, and why there was so much missing data on the secondaries but not the  
4 primary outcome, again, because they were measured at the same visit. The SAQ and the  
5 CCS were supposed to be done in person at 6 months but were missing SAQ data.

6 DR. STONE: So there's actually relatively little SAQ data that is missing.

7 DR. LANGE: About 15 to 20%, Gregg?

8 DR. STONE: No, not that much. Let me put up the slide.

9 DR. LANGE: I guess, Gregg, why is any missing?

10 DR. STONE: Yeah. Yeah, that's a good question. So here is the data. You can see it's  
11 probably about 5 or 7%. So the SAQ is achieved differently than a CCS angina class. CCS  
12 angina class, you're talking to the patient, you're listening to nuances, etc. The Seattle  
13 Angina Questionnaire is a set series of questions that are on a written form and the patient  
14 is given the form and asked to go into a separate room and then to fill out the form. We  
15 don't have the reasons why in about 5% of those patients. You can see it's approximately 5,  
16 maybe 7% overall. Why the patients did not fill out the form, they may not have  
17 understood the questions. I don't want to speculate, but it's a relatively small percentage.

18 DR. LANGE: Okay, thanks. So we all know.

19 And then finally, again, there was information regarding CT angio. During the break,  
20 if the Sponsor could provide the summary statistics of the CTA findings rather than showing  
21 a picture. There were 44 patients that had CT angiograms done. If we could see the  
22 minimal luminal diameters, persistent flow of those signs, I think that would be lovely.

23 DR. RYAN: We have a slide on that, too, if you'd like to see that. We prepared a  
24 backup slide.

25 DR. LANGE: Go ahead and let's see it and see if --

1 DR. RYAN: Sam, can you share that? I think it's 154 or 155.

2 (Pause.)

3 DR. RYAN: Okay, there it is. There we go. So you see there's a total of 44 CTAs, 37  
4 from COSIRA and then seven from the first-in-man study in India, and the minimal luminal  
5 diameters, both proximally and distally, are on this slide for your review.

6 DR. LANGE: Great.

7 Dr. Hirshfeld, does this answer the question that you have?

8 DR. HIRSHFELD: Well, if I interpret this correctly, I think we'd be interested in two  
9 things, what's the minimal lumen diameter in the neck of the device as opposed to the  
10 proximal and distal ends? And secondly, we saw evidence of, in the animal studies, where  
11 the neck of the device was not endothelialized and there was evidence that there would be  
12 flow around the neck of the device. And so I think the prevalence of that phenomenon in  
13 people would be answered by the CTA data and I haven't seen that, you know, specifically  
14 tabulated.

15 DR. LANGE: Could the Sponsor provide that during the break?

16 DR. STONE: We can certainly look at the break and see what we have. We showed  
17 you in a previous slide it's a hundred percent patency. I don't know what we have, as I sit  
18 here, in terms of the NLD at the neck versus elsewhere. I've not seen a case where there's  
19 flow outside the neck, but I understand you want quantitative tabulation which, of course,  
20 you're right in asking for. I'll see what we have at the break.

21 DR. LANGE: That sounds great.

22 And with that, I think all the questions that -- and a lot of them were addressed or  
23 we've offered and have been addressed. So with that, I'm going to ask that we take a  
24 15-minute break. Let's convene at 4:45 Eastern Standard Time and we'll be ready to delve  
25 into the questions and use these answers as a platform to do that. All right, thank you

1 everybody. See you in 15 minutes.

2 (Off the record at 4:29 p.m.)

3 (On the record at 4:46 p.m.)

4 DR. LANGE: Great. I'd like to reconvene the Panel.

5 At this time let's focus our discussion on the FDA questions. Panel members, copies  
6 of the questions are in your folders. I would ask that each Panel member identify him or  
7 herself each time he or she speaks to facilitate the transcription. In just a second I'm going  
8 to ask the FDA to show the first question. Although there are six or seven questions, and  
9 there's about 11 embedded in there, we have a large panel, I think about 19 of us, so what  
10 we'll try to do is get a consensus about the opinion without getting all 19 to speak on each  
11 question. I'll try to summarize. Obviously, if there are any dissenting views or any  
12 additional views you want me to relay to the FDA, we'll do so. With that, let me turn it over  
13 to the FDA to start with our first question.

14 DR. RABEN: Thank you. Over 90% of the COSIRA subjects were taking at least one  
15 antianginal medication and 36.5 were taking three or more at baseline. However, in a trial  
16 intended for a refractory angina population, greater than 25% of subjects were on only  
17 none or one anginal medication. Additionally, at baseline, approximately 75% and 50% of  
18 subjects were taking beta blockers or calcium channel blockers, respectively. No  
19 justification was provided regarding the proportion of patients prescribed beta blockers,  
20 nitrates, and calcium channel blockers in a refractory angina population. Also, no  
21 information was provided about medication compliance, or whether patients were on  
22 therapeutic or maximally tolerated doses.

23 So FDA's Question Number 1: When determining an acceptable indication for use  
24 statement, FDA must consider if the data provided supports a reasonable assurance of  
25 safety and effectiveness for a defined patient population. Please discuss whether the

Free State Reporting, Inc.  
1378 Cape Saint Claire Road  
Annapolis, MD 21409  
(410) 974-0947



1 COSIRA trial identified and enrolled a defined patient population with refractory angina  
2 (despite optimal medical therapy).

3 DR. LANGE: I'll either look for a hand or I'll call on somebody. Dr. Batchelor and  
4 then Dr. Somberg.

5 DR. BATCHELOR: Yeah, so I think the Sponsor and the previous speakers tried to  
6 address this, but I think it's still an issue. I think the way that this is going to probably be  
7 best handled, in my opinion, is in the way in which the future study is handled, and I know  
8 we're going to get into that, if there's a trial that's going to be done, how is it structured,  
9 designed, etc., but I think it's going to be extremely important that this issue is addressed in  
10 another study because I think it was woefully not handled effectively in this trial. And I  
11 might add that historically, that's been an issue in prior studies, as well.

12 So any study moving forward should really identify what is the definition for  
13 refractory angina, who governs the medications, what's the escalation protocol for  
14 optimizing medical therapy, etc., so we've got a really good even spot to start the trial out,  
15 so I'll just make those comments and leave it at that.

16 DR. LANGE: Great. Dr. Somberg?

17 DR. SOMBERG: What I would add to that is I think it's critical that we do have a  
18 defined patient population and the patient population be resistant to medical therapy.  
19 We're dealing with an interventional procedure and a device which is semi-permanent. As  
20 we put in more, some of them may actually be permanent, so I think it's very critical that  
21 we define this population. A quarter of them were on no meds. Well, maybe they couldn't  
22 tolerate any meds, but that's very, very unusual and it's not just the number of meds, it's  
23 the pushing the meds to their maximum dose and their combinations. And we said well,  
24 some of them are symptomatic, they may be dizzy, they may be unable to tolerate the  
25 medication, but look at the large -- almost half of the people were on ACE inhibitors, ARBs,

1 diuretics, and we know when we treat angina patients we diminish other anti -- drugs that  
2 will reduce blood pressure and give the antianginal drugs a greater chance, and they often  
3 control the blood pressure, as well.

4 So I'm very disappointed in the pharmacologic determination that these are resistant  
5 patients, but I'm also saying that we need a better way to make sure that they are not  
6 interventionally treated, able to be treated by other therapies, as well. So if we're going to  
7 define a last-resort population, we've got to do it medically and we have to do it, also,  
8 interventionally.

9 DR. LANGE: Okay. I'm going to call on Dr. Brindis and then Dr. Mathew.

10 DR. BRINDIS: I agree with both comments, this is Ralph Brindis.

11 The only thing I would add is, to build on John's point, is we need dosages both at  
12 baseline and follow-up because that also tells an additional story.

13 DR. LANGE: Thank you, Ralph.

14 Dr. Mathew and then Dr. Cigarroa, and I think we'll close it at that.

15 DR. MATHEW: I agree with all the comments. I would just sort of reiterate the point  
16 that in real life as well as other clinical trials, it's been very difficult to optimize medical  
17 therapy for reasons that have been mentioned. While this proportion of patients on little  
18 to no meds is probably higher than I might have anticipated in anecdotal clinical experience  
19 in other trials, it may be that this is what it is in a small number of patients and the error of  
20 variation there, as well. So I think in a very few trials such as COURAGE, where there was  
21 dogmatic emphasis on medical therapy, or SYMPLICITY, where there was a roll-in period of  
22 aggressive medical therapy for hypertension that could achieve it, but I think this study is  
23 not necessarily an outlier in terms of aggressiveness of medical therapy.

24 DR. LANGE: Dr. Cigarroa. And I saw Dr. Yeh and Dr. Vetrovec, and then if we have a  
25 consensus, we'll move on, but Dr. Cigarroa.

1 DR. CIGARROA: My answer to Question 1 is no for two reasons. First of all, the  
2 patients studied were primarily Caucasian men, which is not reflective of the demographics  
3 of our country, and secondarily, the percentage of individuals who were intolerant to  
4 medicines without insight as to the reasons left a -- you know, which accounted for about  
5 25%, challenges the ability to interpret these results. So for both those reasons I would say  
6 that we have not yet identified a defined patient population.

7 DR. LANGE: Okay. Dr. Vetovec and Dr. Yeh, anything else to add?

8 DR. VETROVEC: Very quickly. While I think it's -- the number of patients who were  
9 inadequately treated is exceptionally high for this study, I would suggest the idea that some  
10 patients, say, who aren't taking the medicines and won't take the medicines and they are  
11 refractory. And so some percentage of that I don't think is a crime, but I do think it's  
12 important to have a much better balance of good medical therapy.

13 DR. LANGE: Okay. And then Dr. Yeh, you get the last comment here.

14 DR. YEH: I was going to say that I'm sort of playing devil's advocate in my own mind,  
15 but while I agree that it may not be clear that all of these patients have truly refractory  
16 angina despite a very optimal medical regimen, I do think at the same time we don't -- it's  
17 quite clear that all these patients are very symptomatic, they all had Class III or IV CCS  
18 symptoms and very low Seattle Angina Questionnaire scores. And so that being said, I  
19 wonder if, although it didn't necessarily find a refractory angina population, it does find a  
20 population with severe symptomatic coronary disease, mostly on medical therapy. So it  
21 may not be the patient population that they explicitly named here, but it is a clinically  
22 relevant patient population, so I'm just wondering how that plays out and how we think  
23 about that as a group.

24 DR. LANGE: Okay. So I'm going to summarize, Bram.

25 It seems that there is -- the feeling is that the Sponsor has identified people with

1 angina, severe angina, but there's insufficient documentation to define whether they  
2 actually have refractory angina or not. The data that's missing are angiographic data that's  
3 missing, medication, dosing, whether on max dose, whether patients are compliant,  
4 whether there was medication adversity or intolerance, and whether the patients were  
5 truly not amenable to revascularization, and future studies need to address this. The Panel  
6 recognized the difficulty, but feels like the Sponsor needs to make a better effort if they're  
7 going to target refractory angina as opposed to severe angina.

8 Dr. Zuckerman, does this address the question satisfactorily? Additional comments?

9 DR. ZUCKERMAN: Yes. And can you --

10 DR. GERSH: I'd like to make an additional comment.

11 DR. ZUCKERMAN: I guess, Dr. Gersh, did you want to make a comment? And then  
12 I'd like to make one.

13 DR. GERSH: I agree with Dr. Yeh completely. I think these are patients with severe  
14 angina. Are they truly refractory? No.

15 DR. LANGE: Okay.

16 UNIDENTIFIED SPEAKER: I think he muted himself by mistake.

17 DR. LANGE: I think you muted yourself, Bernie.

18 DR. GERSH: But in terms of revascularization, at that time, as Gregg Stone pointed  
19 out, there was no formal heart team and I believe -- and I think, moving ahead, if there's  
20 another study I would design it differently in some ways. But I think I'm comfortable that  
21 these patients -- there was a consensus that they were not revascularized. A large number  
22 had bypass surgery, even repeat bypass surgery, multiple stents, and I'm comfortable with  
23 that. You know, 25% zero to one drugs, that is surprising in terms of antianginal therapy.  
24 But I do think they've answered the question that this is a group with severe symptoms and  
25 probably not revascularizable and we have to remember, the formal heart team was not in

1 existence at that time.

2 DR. LANGE: Great. Dr. Zuckerman.

3 DR. ZUCKERMAN: Okay, but the question is we recognize that this is a patient with  
4 significant angina, but is it well enough defined, Dr. Lange, with respect to all the other key  
5 variables that would allow us to write an indications for use? And am I correct in noting  
6 that the majority of the Panel believes not?

7 DR. LANGE: I think that's correct and that is, it's severe angina. Is it refractory  
8 angina? I don't believe that the Panel has enough information to ascertain that's true, so  
9 they're uncomfortable with that.

10 DR. ZUCKERMAN: Thank you.

11 DR. LANGE: Great. Next question from the FDA.

12 DR. RABEN: So regarding blinding and the role of placebo effect. Although subjects  
13 were blinded to their treatment group, there was no assessment of blinding success, such  
14 as a questionnaire asking subjects to identify the study arm to which they believed they  
15 were assigned. Additionally, the rate of missing data for dobutamine stress  
16 echocardiography at 6-months follow up was notably higher in the control group, which  
17 may indicate problems with the blinding. A notable placebo effect was also observed in the  
18 COSIRA control group, which presents challenges for interpreting the data given the limited  
19 sample size.

20 So Question 2a is please discuss the robustness of the trial results given the lack of a  
21 blinding assessment throughout the course of the study and the limited sample size.

22 DR. LANGE: Magnus, you had asked about blinding and confirmation, so do you  
23 want to address, answer this first? Then we'll let other people chime in.

24 DR. OHMAN: Yeah, so it's never easy to do successful blinding in sham-operated  
25 studies, so I think the study group should be congratulated for attempting, but I think some

1 of the key aspects to this, namely how was it assessed, in other words, was there an  
2 objective way of assessing how well the patient -- given that this is a study of symptomatic  
3 behavior rather than objective ischemia.

4 So I think the study design is very challenged and interpretation, therefore, is very  
5 difficult because it's on a symptomatic scale where sham-operated -- but how well blinded  
6 patients are included, it's very, very hard. Had they been combined with objective evidence  
7 of ischemia, then the story could have been told, but the lack of the latter makes it almost  
8 impossible to understand it.

9 DR. LANGE: Dr. Somberg.

10 (Pause.)

11 DR. LANGE: You're on mute. John, you're on mute.

12 DR. SOMBERG: What I was saying was given all the drawbacks of the trial, I think the  
13 blinding issue is one of the better aspects. They tried hard to blind it, and there's nothing  
14 to suggest that it wasn't blinded. So I'm less concerned of blinding, but I'm most concerned  
15 of the missing data. It wasn't in the Canadian questionnaire, so that's good, but everything  
16 else seems to be Swiss cheese.

17 DR. LANGE: Dr. Batchelor and then Jason, I'm going to come to you because it looks  
18 like there's twice as much missing data in one side than another and if it's blinded, what's  
19 that likelihood?

20 So Dr. Batchelor, let me to come to you.

21 DR. BATCHELOR: Just to stand on the blinding issue, you know, I've been involved in  
22 sham trials and as Magnus mentioned, they're very difficult to do in the cath lab. However,  
23 it's crucial that the manner in which the case is dictated and documented has to be done in  
24 a way you can't figure it out. Remember, all of these folks, people are going to look at the  
25 procedure report and follow-up and if the interventionalist doesn't dictate it in a

1 completely generic fashion, you can quickly figure out who got what. So we didn't have a  
2 chance to ask that question, but it's a very important point moving forward in terms of trial  
3 design for the future. We have to make sure that the interventionalists and all the source  
4 documentation is completely obscure in terms of treatment assignment.

5 DR. LANGE: Great. Jason?

6 DR. CONNOR: Yeah, agree on the blinding and it seems like they did a good job to  
7 attempt it. There are so many one-arm trials we see in devices that I think this is higher  
8 quality evidence than we usually see. At the same time, I think the missing data is very  
9 problematic especially given this idea that the primary is there but the secondaries aren't  
10 and it's imbalanced. Like I guess that if people are sent into a room to fill out an instrument  
11 and they miss a few items in the instrument, you can't get a complete score. It sounds like,  
12 from Dr. Stone's explanation, that's kind of what happened, but I think that's an  
13 unsatisfying answer because in theory, people just are missing questions and you can't get a  
14 whole score for the instrument. There's no reason that should be this different between  
15 groups.

16 So I was really hoping the answer on (a) why is there missing data in the secondary  
17 but not the primary and why is it imbalanced would be, you know, a little more thoughtful  
18 or well-studied and I think without that, that it's really hard to answer how effectively it was  
19 blind or why this imbalance, but I'd love to hear if Dr. Wittes has any insight there, too.

20 DR. LANGE: Dr. Brindis and Dr. Allen, I see both hands up.

21 DR. BRINDIS: So I want to also chime in about the missing data, I think it's a  
22 significant problem, particularly it being unbalanced between the two groups. I wanted to  
23 get back a little bit because I think we glossed over the issue that Joaquin brought up, which  
24 is the patient demographics. Because this study was so imbalanced in terms of ethnicity  
25 and gender, this is a huge step backwards and angina is not the same by gender and

1 ethnicity and that's sorely lacking in the study, it needs to be emphasized.

2 DR. LANGE: Great.

3 DR. ZUCKERMAN: Okay, can I pause the Panel a moment? We are going to get into  
4 the demographics question extensively in Question 5, but if we could try to address the  
5 questions that Dr. Lange is posing to us, that is most helpful to the FDA.

6 DR. LANGE: Great. Dr. Allen and then Dr. Starling.

7 DR. ALLEN: Yeah, Keith Allen.

8 You know, I think the Sponsor is to really to be congratulated on doing a sham  
9 procedure. As somebody who was in Gregg Stone's position in 1998 presenting to the FDA  
10 on TMR data, I understand the constraints and concerns of doing studies like this.  
11 Nonetheless, the investigators and the consultants for the Sponsor have all been involved in  
12 sham trials and they clearly understand what's needed to make that successful.

13 I am really concerned about the imbalance and the missing data for the simple  
14 reason that I'm concerned that the patients did know what they got and usually if the  
15 imbalance was on the treated side, I would be less concerned. But normally the group that  
16 doesn't get the therapy, they lose the flip of the coin or so they perceive they lose the flip  
17 of the coin, are less likely to participate in follow-up studies, particularly ones that require  
18 time and effort, like a treadmill test or a thallium study or extended answering of  
19 questionnaires. So it does really concern me.

20 DR. LANGE: Great. So what I'm hearing, and I'm going to summarize, if there are  
21 additional comments, please add them.

22 Dr. Zuckerman, what I'm hearing is that the blinding is absolutely essential, that the  
23 panelists applaud the Sponsor for attempting to do so. Unfortunately, we don't know how  
24 successful they were at doing that because there was no objective assessment of that and  
25 for some studies it may be easy to figure out, and the missing data imbalance makes



1 individuals concerned that despite the Sponsor's attempt is that the patient and/or the  
2 caring physicians may have figured out which treatment they received and as Dr. Brindis  
3 said, having evidence of ischemia, which is much more objective as opposed to a subjective  
4 symptom, would've been essential to help in this particular trial.

5 But Dr. Brindis, an additional comment to that?

6 DR. BRINDIS: Yeah, the only other point which I thought was excellent, and I can't  
7 remember which panelist made it, the ability of actually picking up the presence of the  
8 device on a chest X-ray, and Keith gets the credit for that --

9 (Crosstalk.)

10 DR. BRINDIS: -- and I know that would be in a X-ray report or anybody who's  
11 following the patient would certainly be looking at that chest X-ray potentially with interest.

12 DR. LANGE: Or an echo that shows a metal device in the coronary sinus because  
13 they all had DSEs, so it would be very difficult. So again --

14 DR. ZUCKERMAN: Thank you, that's an excellent response.

15 DR. STARLING: Dr. Lange, could I make a comment, please? He's muted. You're  
16 muted.

17 DR. PAGE: You're muted, Rick.

18 DR. LANGE: I'm sorry. Someone asked if they can make a comment, I didn't see who  
19 made that -- Dr. Starling, I'm sorry.

20 DR. STARLING: Please. I just wanted to reiterate, having participated in at least a  
21 half dozen sham-controlled trials going back to the ACORN era to the most recent trial with  
22 an interatrial shunt, there are numerous safeguards that can be put in place including  
23 having the investigator sign a document at every visit to verify that they have not been  
24 unblinded to the treatment. So yes, this trial was a bit dated from a few years back, but I  
25 don't think we have any assurances as to appropriate blinding in addition to the other

1 issues related to an outside U.S. study with missing data.

2 DR. LANGE: Very nice. So if a future study is done, having that objective evidence of  
3 blinding documented, that will be quite important. Thank you.

4 DR. GERSH: Richard, I think it's important to bring up the interatrial shunt study  
5 because this is an example of where blinding was verified at the time of the procedure, but  
6 it has been verified throughout the follow-up. I think it's critically important.

7 DR. LANGE: Thank you. And both times, you're right, Dr. Gersh.

8 DR. GERSH: I think we just have to think about it in future trial design, I mean, this  
9 was several years ago.

10 DR. LANGE: Yes, thank you.

11 Next question from the FDA.

12 DR. RABEN: Question 2b is, given that some patients did not appear to receive any  
13 benefit from treatment (only 34.6 achieved a primary endpoint success of a change in CCS  
14 of 2 or more and 28.8 demonstrated no change in CCS from baseline), we would like the  
15 Panel to discuss whether patients who are more likely to receive a significant clinical benefit  
16 can be identified prior to implantation of the Reducer device.

17 DR. LANGE: Dr. Yuh.

18 DR. YUH: Thank you, David Yuh.

19 I think this is one question where we're somewhat handicapped by, at least for me,  
20 not being clear on the mechanism of action. I think that mechanism can speak towards who  
21 would benefit and what the limitations of the procedure are. For example, if the putative  
22 mechanism is really at the micro -- redistribution at the microvascular level, you would  
23 think that patients with primarily microvascular disease would stand to benefit the most  
24 and I think for future studies or future investigation, that that sort of patient population be  
25 included as a subset, at least, in addition to the subsets that were used already.

1 DR. LANGE: Great. Dr. Starling.

2 DR. STARLING: Yeah, I think that we don't know the answer to this question and  
3 ironically, the data that we were shown indicated that those that had the lower anginal  
4 class were most likely to respond. Dr. Stone attempted to provide us with an explanation  
5 for that, but I think the answer to this question is we don't know who will derive the most  
6 benefit.

7 DR. LANGE: Okay. Before I ask for any other hands, does anybody feel like they  
8 know, from the data we have, who would receive a significant clinical benefit?

9 (No response.)

10 DR. LANGE: Does anybody -- okay, I don't think -- I'll go to Dr. Somberg, Dr. Page,  
11 and Dr. Ohman.

12 Dr. Somberg.

13 DR. SOMBERG: Yes, I think there's one anatomic/physiologic characteristic that  
14 we're just not focusing in on and that is that the coronary sinus, for the most part, only  
15 does about half the drainage of the myocardium, the total myocardium, and I think we  
16 should assume that the benefit of this device has something to do with its implantation in  
17 the coronary sinus increasing pressure and therefore more proximal perfusion.

18 So I think it would be useful, and I'm surprised preliminary work wasn't done on this  
19 to try to define where the ischemia is coming from in a stress thallium or wall motion study  
20 and then trying to see if that is a plausible area, that the coronary sinus in that patient will  
21 drain because, for instance, the -- what was it, the Bayesian drainage has anything from 30  
22 to 60% of alternative venous drainage. So I think we're not correlating the plausible  
23 intervention with the physiology of the heart and maybe some studies looking at the  
24 distribution of the drainage compared to where the ischemia's coming from could be used  
25 to try to identify patients who would benefit.

1 DR. LANGE: So at this point, you propose that that's a possible thing to look at but  
2 you don't -- you're not saying at this particular point you know which distribution --

3 DR. SOMBERG: Of course not, nobody knows because it's not been explored, but I  
4 think it's a hypothesis to generate and to try to test instead of putting it to everybody who  
5 has something to do with the circumflex in the LAD and only getting a fraction of those  
6 responding.

7 DR. LANGE: Dr. Page and then Dr. Ohman.

8 DR. PAGE: Yeah, I don't think I have a good enough understanding of the nature of  
9 the patients and the mechanism of action to really say that we can predict who will benefit  
10 from this. I'm always struck when I hear testimonials from patients, and this time their  
11 significant others, and the fact that they're looking for something, but at best, 34% would  
12 have a positive response as defined by two levels of improvement in class and almost half  
13 that many who got a sham procedure had a similar response. So I think the answer to the  
14 question is we can't predict which of those patients are going to respond favorably.

15 DR. LANGE: Okay. Dr. Ohman, I think you have the last word on this.

16 DR. OHMAN: Yeah. I want to echo what others have said, I think it's unclear. But  
17 unfortunate for this study and the very limited sample size, the best therapeutic data seem  
18 to be in those who took the minimal amount of medication, zero or one, which is rather  
19 paradoxical, so there may be some clue there to the future but it may actually not address  
20 the refractory angina.

21 DR. LANGE: Okay. And Dr. Borer, I'll let you have the last word on this, sir. Unmute  
22 yourself, Jeff.

23 DR. BORER: Okay, can you hear me?

24 DR. LANGE: Yes, sir.

25 DR. BORER: Great. Okay, I'm disturbed by what I think is a fundamental problem

1 here. Some patients did not appear to be receiving benefit from the treatment, only a third  
2 achieved a primary endpoint success of a change in CCS of at least two grades. This is an  
3 angina study. For 50 years the standard for angina for success or magnitude of success in  
4 minimizing angina, preventing angina, has been a quantitative exercise test, how many  
5 minutes will a patient be able to go on a quantitative exercise test, and that's not what  
6 we're talking about here, we're talking about a subjective endpoint and I don't -- I think that  
7 that's a real problem. It just magnifies all the other problems that have been discussed  
8 already. I think in developing a new study there has to be a much greater emphasis on  
9 quantitative assessment.

10 DR. LANGE: Okay. So I'm going to summarize.

11 Bram, I think the consensus of the Committee is no, we don't know the answer to  
12 this, we don't know who would benefit. And some information that's lacking is the  
13 mechanism, the correlation of the mechanism with any clinical benefit. There's a lot about  
14 the native patients, we don't know demographics, nature, distribution of coronary anatomy,  
15 so that's ill defined. And then finally, as Dr. Borer mentioned, we don't have any  
16 quantitative data to ascertain which patients received it and so all of those things are  
17 limitations that prevent us from saying who would benefit from this.

18 DR. ZUCKERMAN: Okay, so Dr. Lange, that's an excellent summary after a great  
19 Panel discussion. Just a quick comment to respond to Dr. Borer. In 2010, when this study  
20 was being discussed with FDA, not unexpectedly, Dr. Borer, FDA did recommend the  
21 treadmill endpoint that has been used as a standard. The Sponsor disagreed and was one  
22 of the reasons why the study was done overseas, the Sponsor was confident that this was  
23 an objective enough endpoint. Certainly with COSIRA-II, again, the FDA has really tried to  
24 urge the Sponsor to use the exercise treadmill endpoint as the primary endpoint and when  
25 we get further into the questions, I'd like you to further expand with other Panel members

1 on how objective some of these angina endpoints are. Thank you.

2 DR. LANGE: All right. I think we can move on to the next question. Thanks to the  
3 Panel for opining on that.

4 DR. RABEN: Thank you, Dr. Lange.

5 So the next question is regarding the primary effectiveness endpoint. A two or  
6 greater CCS grade change at 6 months was the primary effectiveness endpoint of the  
7 COSIRA trial. Primary endpoint success was observed in 34.6% of subjects treated with the  
8 Reducer device, while 15.4% of subjects achieved success in the control group. In 28.8% of  
9 the Reducer group and 57.7% of the control group, no change in CCS was observed.  
10 However, angina can be a placebo-responsive condition. Exercise tolerance testing  
11 provides an objective measure of functional capacity in myocardial ischemia. Other clinical  
12 trials evaluating anti-ischemic treatments have used exercise testing results as a primary  
13 effectiveness endpoint.

14 So Question 3a: Please discuss and comment on the subjective assessment of angina  
15 (change in CCS grade) as a clinically meaningful correlate of ischemia to support a  
16 reasonable assurance of Reducer device effectiveness.

17 DR. LANGE: Dr. Borer. And again, unmute yourself, Jeff.

18 DR. BORER: Yeah, a change in CCS grade is certainly not an unreasonable measure of  
19 whether an intervention does something or doesn't, but it's not the best evidence. The  
20 best evidence would be a quantitative assessment that includes an electrocardiographic  
21 measure of ischemia, as well as an antianginal effect. In fact, although I don't know how  
22 the rules are written for devices for this kind of thing, for drugs, for years, an antianginal  
23 drug has to prevent angina by an anti-ischemic mechanism, there has to be a reduction in  
24 ischemia. So I would say that a reduction in CCS grade is nice, but it's not sufficient.

25 DR. LANGE: So I'm going to ask all Panel members to -- a yes or no vote on this, that

1 is going forward, either another study or a postmarket study, how many of you would --  
2 Ralph's shaking his finger at me. I'm not assuming either, okay, I'm not assuming either, but  
3 how many would recommend that the Sponsor use only CCS grade change and would say  
4 that that would be a reasonable assurance of effectiveness? If you say yes, let me see your  
5 hand.

6 (Show of hands.)

7 DR. LANGE: If you say no, let me see your hand.

8 (Show of hands.)

9 DR. LANGE: Then I'll call on Dr. Brindis, okay?

10 DR. BRINDIS: One comment that hasn't been made yet is a publication showing the  
11 lack of correlation between physician reporting of CCS and SAQ scores. I might assume it  
12 might be more aligned in this study because it was blinded, but it is interesting that patient-  
13 reported outcomes is what clinicians like and we haven't mentioned that before.

14 DR. LANGE: Thank you, Ralph.

15 DR. SOMBERG: Dr. Lange, can you tell us the results of the hands? I didn't see  
16 everybody.

17 DR. LANGE: Yes, there were no hands up with reasonable assurance and all hands  
18 went up with they would not recommend that as the sole measure of reasonable assurance.

19 DR. SOMBERG: Thank you.

20 DR. LANGE: It was unanimous.

21 Dr. Yeh.

22 DR. YEH: I mean, I'm thinking about what is the more important endpoint here and  
23 for this category of patients, I just do want to raise the possibility that the goal is the  
24 reduction of angina, not the reduction of ischemia, and one can imagine therapies, not is  
25 the case, per se, but one can imagine future therapies that treat patients with angina not

1 through reduction of ischemia but via some other mechanism. And so this question as  
2 framed, I think, sort of implies that ischemia is the thing that is the gold standard for what  
3 we're trying to reduce here, but I mean, I think, arguably for this patient population, they  
4 don't care about their objective measurement of ischemia, what they care about is what  
5 they're feeling, their angina. And so among these two things about whether or not I sort of  
6 think that the exercise tolerance test is a better study versus angina class, I'm sort of more  
7 concerned about the blinding issue that we talked about previously for this particular issue.

8 If there were perfect blinding, I would be, I think, much more comfortable and I  
9 don't have, I think, the historical viewpoint of this that Dr. Borer does, but I'm more  
10 comfortable. We use Seattle Angina Questionnaire very frequently and I think what it's  
11 subject -- it's a very stable estimate over time and what I think I'm more concerned about  
12 when we use those types of endpoints is the blinding issue, but also exercise tolerance test,  
13 it certainly also could be a placebo-responsive condition. I just want to say that people's  
14 exercise effort may be related to also those things. I don't disagree with anything that's  
15 been said, I just wanted to put that point out there.

16 DR. LANGE: I appreciate it.

17 I've got the following hands up in this order: Dr. Cigarroa, Dr. Mathew, Dr. Borer,  
18 and Dr. Somberg.

19 Dr. Cigarroa.

20 DR. CIGARROA: So I certainly agree that there are approaches that can impact the  
21 perception of ischemia, whether it's dysthymia or classic chest discomfort through means  
22 other than mitigating ischemia. However, in this particular situation where the putative  
23 mechanism is a primary means of impacting ischemia, I think assessing ischemia in addition  
24 to the measurements of angina is critically important and I think especially even more so  
25 given the small sample size that we are looking at in this trial in which flipping two patients



1 in either direction on how one classifies would change the statistical significance. So not  
2 knowing the mechanism of action, having the purported mechanism of action improving  
3 symptoms via reduction in ischemia, I think one has to have a quantitative assessment of  
4 ischemia.

5 DR. LANGE: Okay. Dr. Mathew.

6 DR. MATHEW: Yeah, Verghese Mathew. Thank you, Dr. Lange.

7 So I just want to convey my agreement with what Dr. Yeh brought up. From a  
8 patient perspective, as we see these people in real life, I mean, they're really concerned  
9 about their symptoms and they're not coming in to sort of say we're wanting our ischemia  
10 reduced, they want their symptoms reduced, and I think it goes to the patient populations  
11 that we're testing here and what the goals are, and I think that a larger denominator  
12 would've been helpful in that regard.

13 But I think the limitations of measuring ischemia and the variableness of ischemia, I  
14 think there's a scientific discomfort on my part, and maybe our collective part, that not  
15 having an objective measure like that is, again, discomforting. I don't know that I would  
16 solely rely on Canadian Cardiovascular Society class, but some complement that satisfies  
17 our scientific curiosity, I think, would be important moving forward in another trial.

18 DR. LANGE: Thank you, Dr. Mathew.

19 Dr. Borer.

20 DR. BORER: Yeah, thank you.

21 There are a couple of issues here. First of all, I agree completely with the statement  
22 about blinding that I think Dr. Yeh made a few minutes ago. But the primary reason I think  
23 that the demand is made that reduction in angina as measured on an exercise tolerance  
24 test is accompanied by evidence of a reduction of ischemia is safety. If you were able, with  
25 any kind of intervention, to allow people to walk longer on a treadmill before they develop

1 limiting angina and you weren't limiting their ischemia, they could die while they're walking  
2 on the treadmill. That would be a bad thing and I think we'd all agree with that.

3 So I think that there really must be some assessment of ischemia associated and not  
4 just for purposes of completeness. In Europe the primary criterion for approving an  
5 intervention as an antianginal is evidence of a reduction in ischemia, not angina, so I think  
6 that it's very important to do both, but to quantify both.

7 DR. LANGE: Okay.

8 DR. GERSH: Richard?

9 DR. LANGE: Dr. Somberg, anything?

10 And then, Bernie, I'll come to you.

11 DR. SOMBERG: I would just agree with Jeff's statement, we want to do both, and I  
12 think it's very concerning that we would just put all our emphasis on CCS because we can  
13 give an analgesic and you can take an opioid and you wouldn't feel pain but you would have  
14 significant physiologic consequences. So I think we want to have something to substantiate  
15 that it's just more than a pain reliever, that it's actually an ischemic reliever, as well. So we  
16 might know the absolute mechanism of this device but we want to know that it's doing  
17 something to the ischemic process.

18 DR. LANGE: Okay. And then Dr. Gersh, Allen, and then I'll summarize.

19 DR. GERSH: I totally agree that we need both, I think both measurements of  
20 ischemia or assessments of ischemia and severity of angina have their limitation and one  
21 great example of that is the ISCHEMIA trial that's just being published. I mean, the  
22 rationale for this trial was that people with moderate to severe degrees of ischemia would  
23 benefit from revascularization and although not yet published, it's been presented. In fact,  
24 there was no relationship between the severity of ischemia and the outcome of the  
25 ISCHEMIA trial. So it has its limitations, too, as does the assessment of angina, but -- clinical

1 assessment, but I prefer to see both.

2 DR. LANGE: Great. I'm going to summarize, but -- and then Keith or Randy, if you  
3 have anything to add to it, raise your hand and I'll let you do so.

4 I think the consensus is, Bram, is that the Committee is recommending both. The  
5 CCS alone is not an effective way to reasonably assure effectiveness.

6 So Keith, if you want to go ahead. You're on mute. This is Keith Allen, Dr. Keith  
7 Allen.

8 DR. ALLEN: Yeah, Keith Allen.

9 So I think Joaquin really summarized, and I do think it's very important that  
10 mechanism play a role in when you're choosing your -- how you're going to judge your  
11 primary endpoint, and I just want to give you two examples to maybe make it clear for  
12 people listening.

13 You know, in this case the mechanism is improvement in profusion, it clearly is  
14 related to ischemia. Picking an endpoint like ETT makes a lot of sense. On the other hand,  
15 let's say the mechanism was de-innervation (ph.), that's how you are fixing people's angina,  
16 then improvement in profusion and reduction in ischemia wouldn't be a reasonable  
17 endpoint and perhaps CCS class in a well-blinded study is a much better endpoint. I don't  
18 think you can separate the two, but for this particular example, CCS class is not the right  
19 endpoint.

20 DR. LANGE: All right, great.

21 Dr. Starling, do you want to add anything to that, sir?

22 DR. STARLING: Yeah, I just wanted to say, and this is coming as being a heart failure  
23 specialist, this discussion is somewhat analogous to NYHA class and cardiopulmonary  
24 exercise testing or ejection fraction as far as gaining some objectivity, and this is sort of a  
25 patient-reported outcome, although they did have a blinded assessor, but we were never

1 really told, at least in COSIRA, we were never really told what the rigor was as far as that  
2 blinded assessment other than it was a blinded assessment. But fully agree, we need  
3 something more than CCS.

4 DR. LANGE: Okay. Bram, do you need any other additional comments from the  
5 Panel or do you feel like you have enough information?

6 DR. ZUCKERMAN: This was an excellent discussion, please go on to the next  
7 question.

8 DR. LANGE: Thank you.

9 Go ahead, Sam.

10 DR. RABEN: Thank you, Dr. Lange.

11 Question 3b, we would like the Panel to discuss and comment on the overall primary  
12 effectiveness rate of 34.6%, given the permanent implant nature of this device and  
13 vulnerability of this no-option patient population.

14 DR. LANGE: Dr. Borer.

15 DR. BORER: You know, this is a question that can only be answered in a relative way.  
16 The overall primary effectiveness rate of 34.6% using the endpoint that was used may be  
17 wonderful. Do we know anything that's better? I don't. You know, in people with  
18 refractory angina, we don't have a comparator. So I think that that's reasonable enough.  
19 Yes, the patient population is relatively vulnerable, but the goal of the therapy is to improve  
20 their quality of life, to make them feel better, and if the effectiveness endpoint is achieved,  
21 the measure that was used is achieved in a the third of the patients, that may be pretty  
22 good. So I don't think that's really an important issue. An important question, I think the  
23 important question is do most -- do people generally feel better, because that's what you're  
24 trying to get them to do. We're not talking about survival with this trial.

25 DR. LANGE: I've got Dr. Allen and Dr. Batchelor.

1 DR. ALLEN: So the 34, 35% rate compared to 15% in the control, I think that's an  
2 okay number. If you had a 34% cure rate of a Stage IV cancer drug, you would have a real  
3 winner. I think at least what concerns me with it's not so much the number, but the  
4 veracity of the number and whether the clinical trial was done in such a way with all  
5 considerations such as missing data and things like that, whether I believe that they got to  
6 that number and it was powerful. So the 34% I'm okay with, I just don't know that the trial  
7 makes me comfortable that they actually got that number.

8 DR. LANGE: Okay. Dr. Batchelor, Dr. Cigarroa, Dr. Brindis, and Dr. Somberg.

9 DR. BATCHELOR: Just to springboard off of what Keith just mentioned, I actually  
10 think the response rate is pretty good if the patient population was truly defined as truly  
11 refractory angina patients. I'm sure Dr. Gersh, Dr. Borer, and others who have a lot of  
12 wealth of experience in this thing can chime in, and I think Dr. Gersh mentioned that this is  
13 actually not that bad of a response rate. My issue is it's not -- the absolute risk reduction  
14 and the relative risk reduction is actually pretty good, it's the confidence intervals around  
15 this result and the sample size that we have concerns about where the -- you know, the  
16 confidence intervals are fairly wide. In just a little of the sample size, it only takes one or  
17 two patients to shift the statistical significance. So I think we're talking about two different  
18 things. The delta, the difference, it's a pretty good delta but the confidence around that  
19 reduction is questionable.

20 DR. LANGE: Great. And I've got Dr. Cigarroa, Dr. Brindis, Dr. Ohman, and  
21 Dr. Somberg before that, go ahead.

22 Dr. Cigarroa.

23 DR. CIGARROA: My points are along the same line as Dr. Batchelor and Dr. Allen's.  
24 The number itself is a reasonable effectiveness rate; however, the believability of that  
25 number, given the issues that we discussed about the trial and the sample size and blinding,

1 calls into question whether I believe that number.

2 DR. LANGE: Thank you --

3 (Crosstalk.)

4 DR. BRINDIS: Yeah, so my comment is if we believe the number, NNT is just fine for  
5 me assuming that the device is safe, you know, because that's the tradeoff. So that would  
6 be my comment.

7 DR. LANGE: Thank you.

8 Dr. Somberg and then Dr. Ohman, and then I'll try to summarize.

9 DR. SOMBERG: I guess I'm going to be in a minority position here that I don't think  
10 that number -- I mean, we're talking about 34% response in the intervention versus 15 in  
11 the control, you know, you're going to benefit 15 out of 100 people for a permanent device.  
12 I think that's very marginal.

13 DR. LANGE: Okay. I appreciate even -- especially dissenting opinions.

14 Dr. Ohman.

15 DR. OHMAN: Yes, so I'd like to turn it around and think of it in a different term.  
16 Nearly half the patients did not respond. Now, I can live with that if we showed a reduction  
17 in ischemia because of the variability in angina, but the fact of the matter, that I put in a  
18 device and about half of them are not going to have any benefit at all using this  
19 classification worries me, and that's why I think the importance of angina plus ischemia  
20 becomes critical.

21 DR. LANGE: So I'm going to -- and Dr. Borer.

22 DR. BORER: Yeah, I just want to come back to something John Somberg said a little  
23 bit earlier, he said you could use an opioid and relieve the pain, and you could. Janet Wittes  
24 and Bram will remember, the first meeting at which Janet and I participated for this  
25 Committee was when partial laser revascularization was put forward and there were no

1 exercise tests, there was no measure of ischemia, and I remember making a suggestion that  
2 they could've given everybody morphine. And the trial was done in Detroit at Wayne State  
3 and the progenitors of the trial were up in arms, you can't use opiates for treatment of  
4 coronary disease, that's a terrible idea. Well, no, it's not. If the idea is just to make them  
5 feel better, it may be a good idea. The problem is that the opiate won't relieve ischemia, so  
6 we have no reason to believe it would, we'd have to test it. So I think that's a key issue to  
7 remember. Chest pain isn't enough.

8 DR. LANGE: Keith, the last comment.

9 DR. WITTES: Can I jump in?

10 DR. ALLEN: Yeah. So let me just make -- I just want to --

11 (Crosstalk.)

12 DR. ALLEN: -- because we're not really here to debate TMR and I was involved in a  
13 lot of the TMR studies, but I would beg the issue with what was just stated because Dan  
14 Burkhoff in 1999 presented the ATLANTIC trial which was published in the *Lancet* and their  
15 FDA-mandated primary endpoint was ETT and patients that got TMR had a 111-second  
16 improvement over the control group following surgical sole therapy, TMR. So we do have  
17 objective evidence in other antianginal devices so we're not here to debate TMR but I just  
18 think we need to be careful about statements that maybe aren't correct.

19 DR. LANGE: All right. I'm going to summarize. Some Panel members felt like a 34%  
20 rate was -- the 34, 35% was satisfactory, others did not. Bram, that's obviously balanced  
21 with the safety and the placebo effect --

22 (Audio feedback.)

23 DR. LANGE: -- confidence in the reliability of this information and also of the patient  
24 population. I think there was less concern about the number and more concern about --  
25 especially if there's a marginal benefit, more concern about the reliability of the patient

1 population and the assessment. Does that answer your question?

2 DR. ZUCKERMAN: Yes, it does. It's a very good summary of a very good discussion.

3 DR. LANGE: Thank you.

4 All right, Sam.

5 DR. RABEN: Great, thank you.

6 The next question, I would like to discuss secondary effectiveness analysis. In a  
7 secondary effectiveness analysis in COSIRA, the ETTs, specifically the bicycle ergometry and  
8 dobutamine stress echocardiography, were used to objectively assess ischemia. Subjects in  
9 the Reducer group had numerically longer exercise durations (mean increase of 64.7  
10 seconds versus a mean increase of 4.3 seconds) and time to ST-segment depression versus  
11 control patients, which was 76.3 seconds versus 33.8. However, the study was (1)  
12 underpowered to detect improvements in functional ischemia between treatment groups,  
13 and (2) there was a substantial amount of missing information. For DSE data, missing data  
14 was noted in roughly 15% of Reducer subjects, while 30% was missing for control subjects.  
15 Total exercise duration testing was missing in about 25% of all patients, and ST depression  
16 data was missing from 70 to 88% of patients. These two factors impact the conclusions that  
17 may be drawn from these ischemia data.

18 Question 4a is, please discuss overall Reducer device effectiveness observed in the  
19 COSIRA trial, considering the small sample size (underpowered study for ischemia  
20 endpoints), high control group response rate, significant amounts of missing data for  
21 objective ischemia assessments, and lack of pre-specified hypotheses tests for objective  
22 ischemia assessments.

23 DR. LANGE: Dr. Gersh and Dr. Allen and Dr. Somberg.

24 DR. GERSH: Thank you. I want to go back to the *New England Journal* publication of  
25 this paper and the editorial that was written by Christopher Granger and myself, and

Free State Reporting, Inc.  
1378 Cape Saint Claire Road  
Annapolis, MD 21409  
(410) 974-0947



1 basically, what we see in there, I think -- and that is the bar that was set for the magnitude  
2 of benefit was quite on, I mean, two Canadian CCS classes, as was pointed out by Gregg  
3 Stone. If you go from four to two, that's a huge change, and if you go from three to one,  
4 that's a large change, but what we did say, the study is underpowered. If other studies  
5 show this then we may have a highly effective treatment for a very difficult group of  
6 patients, but it is an underpowered study. I do not understand still why it was stopped.

7 And thirdly, I think that the -- in terms of the secondary endpoints, the discrepancy  
8 and the lack of missing data and the discrepancy between the two groups totally invalidates  
9 that. I don't know what to make of the secondary endpoints at all, so I think if this degree  
10 of benefit can be substantiated by the studies, and I think we have a very interesting and  
11 effective form of treatment, and that's what we said a couple of years back.

12 DR. LANGE: Great. I've got Dr. Allen, Dr. Somberg, and then Dr. Wittes.

13 Keith, you're muted. Sorry about that.

14 DR. ALLEN: My fault, excuse me. Keith Allen.

15 Now, I have a hard time, quite honestly, putting much faith in the secondary  
16 effective analysis and obviously the company and the Sponsor didn't either because they  
17 didn't do a pre-specified hypothesis, they didn't power the study for any of these, so they  
18 clearly weren't confident they might meet these secondary effectiveness endpoints, so I'm  
19 not sure why I should be empowered to think they're positive.

20 DR. LANGE: Okay. Dr. Somberg and then Dr. Wittes.

21 DR. SOMBERG: With that degree of missing data in the control group of 30%, the  
22 data really will not be informative. I think any study that's presented to the FDA that has a  
23 critical component with that much missing data really should not be brought to the Panel.

24 DR. LANGE: Okay. Janet.

25 DR. WITTES: Yeah, so Dr. Zuckerman asked us about how Jason and I would handle

1 missing data, and I agree that when you have 30% in one group and 15% in the other group  
2 missing, there's no --

3 DR. LANGE: Janet, we just lost you.

4 DR. ALLEN: Janet needs to pay her Wi-Fi bills better.

5 UNIDENTIFIED SPEAKER: Yeah, she might want to stop sharing her video, that might  
6 help her.

7 DR. LANGE: Yeah. Janet, if you can hear me, turn your video off, leave your audio  
8 on.

9 DR. WITTES: My video is off and my audio is on, can you hear me?

10 DR. LANGE: We can, so I'm sorry, we lost you just at the -- go ahead.

11 DR. WITTES: Okay. So what I was going to say is --

12 (Audio feedback.)

13 DR. LANGE: All right. Janet, we're having trouble with your audio, as well.

14 DR. WITTES: Okay. Never mind.

15 DR. LANGE: I'm so sorry. I know what you're saying is important. As a non-  
16 statistician, I can't even help you out.

17 (Audio feedback.)

18 UNIDENTIFIED SPEAKER: -- format.

19 DR. LANGE: Yeah. So Bram, I'm going to summarize this.

20 So Dr. Zuckerman, I think everybody recognizes and the Sponsor did, too, it's  
21 underpowered, lacks a hypothesis for any of the secondary endpoints. We also  
22 acknowledge that it's missing data, as well, and I think the Panel is actually putting --

23 (Audio cut out.)

24 DR. ZUCKERMAN: I think the Sponsor was partially aware of some of these problems  
25 and that's why the Sponsor presented to the Panel the REDUCER-I observational data which

1 are nonrandomized data where the Sponsor believes that the same trends are being shown  
2 and that these are important data. Could someone like Dr. Allen just comment on this  
3 strategy? Do the REDUCER-I data help you at all, Dr. Allen?

4 DR. LANGE: And then after that, Jason, I'm going to come to you because you had  
5 your hand raised.

6 So Dr. Allen.

7 DR. ALLEN: Thanks, Dr. Zuckerman.

8 The short answer and very short answer is no, and I've really discounted the  
9 REDUCER-I data altogether. Those were observational studies; a lot of these are single  
10 center, some of the data they presented were single center, obviously none of it is  
11 randomized. Where the placebo effect concerning any data that's not randomized and  
12 appropriately blinded comes into suspect. Commingling randomized data like COSIRA with  
13 REDUCER-I data actually only confuses the issue, in my mind, and it draws into question a  
14 lot of the issues that we've already discussed.

15 DR. LANGE: Dr. Connor, you had your hand up and if you want to make a comment  
16 there, and also whether you think that the REDUCER-I data overcomes the shortcomings  
17 that we've identified.

18 DR. CONNOR: Yeah, my comment was just going to be regarding pre-specification of  
19 secondary endpoints, pre-specification of missing data, things like that. I had requested the  
20 SAP because I always request the SAP when I'm on this, and didn't receive it. I received the  
21 protocol but I've not seen any SAP that may have gotten into, you know, alpha spending,  
22 how secondary endpoints would be in alliance from a prospective standpoint, what order  
23 they focus on, things like that. So we don't know that, either.

24 Regarding the supportive study, I mean, I think it's -- you know, it just doesn't, as far  
25 as I can tell, also have at least pre-specified performance goals and things like that. If it

1 were a one-armed trial it would at least have those and they put a stake in the ground and  
2 see what -- so I mean, we've seen the data, but I think that's it.

3 DR. LANGE: Okay. Dr. Starling.

4 DR. STARLING: Thank you, Dr. Lange. Randall Starling.

5 I'd like to respond to Dr. Zuckerman's question just by saying that I think the concern  
6 raised is we have an inadequate COSIRA trial and then an open-label, non-randomized  
7 REDUCER, and I would raise the word cardiomyoplasty and just say that the placebo effect, I  
8 think, is an overwhelming consideration that we have to reflect on.

9 DR. LANGE: Okay. So I think what you're hearing, Dr. Zuckerman, and I appreciate  
10 you advancing the question, but I'm not -- I don't think the Panel feels like the information  
11 or data derived from REDUCER-I overcomes the shortcomings of the secondary  
12 effectiveness analysis.

13 DR. ZUCKERMAN: Thank you, Dr. Lange, a very helpful discussion.

14 DR. LANGE: Thank you.

15 Sam, next question, sir.

16 DR. RABEN: Thank you.

17 For Question 4b, please also discuss if additional premarket objective ischemia  
18 assessment data are needed to support the Reducer effectiveness (e.g., primary endpoint of  
19 the COSIRA-II trial, which was change in total exercise duration in modified Bruce treadmill  
20 exercise tolerance testing at 6 months).

21 DR. LANGE: And I think addressing the previous question, we've addressed this, as  
22 well. The Panel feels like obviously anginal information and relief of angina, change in angio  
23 class is important, but it has to be combined with some objective evidence.

24 MS. ALIKHAANI: That's not what he asked.

25 DR. LANGE: Please discuss if additional premarket objectives are needed to support

1 Reducer effectiveness. Okay, I'm sorry, I misunderstood. In other words, whether the  
2 current data are sufficient or whether there's additional objective ischemia assessment.

3 Dr. Somberg. Unmute.

4 DR. SOMBERG: Since we haven't established the appropriate patient population, I  
5 should say -- resistant population, COSIRA-II needs to do that, and then you need to show  
6 the ischemia measurement using exercise testing or some other, dobutamine study, echo or  
7 what have you. So I think COSIRA-II needs to define the population and then show that you  
8 have a reduction in angina and that is related to ischemia, they're all tied together.

9 DR. LANGE: Okay. Dr. Page.

10 DR. PAGE: Yeah, I'm not exactly sure what the question is, whether we are asked if  
11 we think premarket objective assessment is needed or if it were done, whether we would  
12 need some objective evaluation. I think we've answered the second part, but it just seems  
13 like we're getting ahead of ourselves right now if we're saying we need further premarket  
14 data because it gets to whether we're including whether it's effective or not and we can tip  
15 our hand now or not.

16 DR. LANGE: Can we get some clarification from the FDA? What would you like the  
17 Panel to opine on?

18 DR. ZUCKERMAN: Sam, do you want to begin on that?

19 DR. RABEN: Although I think the Panel has done a great job of discussing the value  
20 of an ischemia assessment, I think it would also be valuable for us to understand if the  
21 current totality of the data is sufficient or if additional data is needed to further support the  
22 proposed indication.

23 DR. PAGE: Okay, go back to my question, so if you substitute "whether" instead of  
24 "if," we're getting into effectiveness. If you're saying it's a hypothetical if we were doing  
25 more would we need objective evidence in addition to the subjective, I think we've

1 answered that. So if you're changing this middle, this first "if" in the question to "whether,"  
2 then we're getting into effectiveness, we can do that, but it seems like we're getting ahead  
3 of ourselves.

4 DR. LANGE: I agree. So can we hold off on that, Sam, until we talk about  
5 effectiveness?

6 DR. RABEN: I think that's reasonable. Bram, do you concur?

7 DR. ZUCKERMAN: Yes. Thank you, Dr. Lange and Dr. Page.

8 DR. LANGE: Okay, I want to have something entered into the record. Dr. Janet  
9 Wittes typed in a response to the previous question and I just want to read it so it's entered  
10 in the record.

11 She said, "Sorry you couldn't hear me. I find the fact that secondaries are so weak  
12 and have so much missing data really bothers me. The combination of a primary outcome  
13 that is subjective with non-convincing data on the secondaries makes the whole study  
14 uninterpretable coupled with the fact that this population is so unlike the U.S. is another  
15 problematic issue. Dr. Zuckerman asked whether we two statisticians should suggest what  
16 we need to do with the missing data, my response is that too much data are missing to be  
17 able to do any statistical fix." I just read that for entering into the record and thank  
18 Dr. Wittes for her comments.

19 All right, Sam, let's go to the next question.

20 DR. RABEN: Great. So now I'd like to discuss some of the study limitations. As  
21 discussed in FDA's executive summary, there are limitations to the currently provided  
22 dataset. These limitations include, but are not limited to:

- 23 • Lack of a non-exercise primary effectiveness endpoint and no pre-specified  
24 hypothesis tests for objective secondary endpoints;  
25 • Small sample size;

- 1 • Significant missing secondary endpoint information;
- 2 • Lack of a formal assessment for coronary sinus stenosis or severity;
- 3 • Lack of evidence of a coronary sinus pressure gradient across the device;
- 4 • High placebo response rate;
- 5 • Trial cohort demographics are not representative of the U.S. population.

6 In addition, the Reducer device is intended to create a coronary sinus stenosis  
7 resulting in a functionally significant increase in coronary sinus pressure gradient that may  
8 reduce myocardial ischemia by redistributing subepicardial blood flow to the  
9 subendocardium. However, in vivo animal studies were not sufficient to confirm tissue  
10 coverage to restrict coronary sinus blood flow to the Reducer's central orifice. Further,  
11 neither in vivo animal nor clinical data were provided to show that the Reducer device  
12 performed as intended, because there were no adequate studies that assessed:

- 13 • The presence of severity of coronary sinus;
- 14 • A coronary sinus pressure gradient across the device;
- 15 • The association of a coronary sinus stenosis or a coronary sinus pressure gradient  
16 with reduced angina or ischemia.

17 So Question 5a: Please discuss and make recommendations whether additional  
18 premarket data from a randomized sham-controlled clinical study are needed to support  
19 the safety and effectiveness of the Neovasc Reducer System given the concerns and  
20 limitations with the current available data.

21 DR. LANGE: Okay, and the FDA is going to be really interested in this particular  
22 question.

23 Dr. Page.

24 DR. PAGE: Yeah, now we're getting to the "whether" question and I don't -- a  
25 number of us have been doing this for quite a while and I don't know when I've seen a litany

1 of problems this long and this significant, and I think I could add to them, but I won't. I  
2 think it's -- from my standpoint it's absolutely clear that we need another trial that is better  
3 performed, in the U.S., and provides both objective as well as an arguably important  
4 subjective measure of angina and ischemia.

5 DR. LANGE: Great. Dr. Allen. You're muted, Keith. I'm sorry. You're still muted.

6 DR. ALLEN: Sorry, I'm a heart surgeon, I'm a slow learner. So I agree with what  
7 Richard said completely. I find it very odd that the very postmarket study that the company  
8 is proposing actually should be their pivotal IDE trial and I find it unusual that we're asked  
9 to essentially approve a product and then do the study that we would all like to see done,  
10 and so we absolutely need more and a better study.

11 DR. LANGE: Dr. Batchelor and then Dr. Zuckerman.

12 DR. BATCHELOR: So to follow up with that question, perhaps it's the result of this  
13 being submitted as breakthrough technology, and so I would like to hear the FDA's response  
14 on is the reason for the -- was the answer to the question that Keith has proposed really  
15 related to the breakthrough technology. And to just answer the Question 5a a little bit  
16 more, I think most of us, if not all of us, feel that the current data so far is not good enough  
17 to provide any answer in terms of effectiveness that we're confident on.

18 However, I must say I think it's a safe procedure, I think it's relatively safe. I think  
19 we've seen 2,000 procedures; it's not 100% safe but within the realm of interventional  
20 procedures I think it's quite safe. So I think what's outstanding is effectiveness, and I think  
21 a new trial designed properly in the United States with adequate site selection that would  
22 represent the U.S. demographic better would be obviously strongly advised. But maybe we  
23 can ask -- have the FDA answer the question of is this all because of breakthrough  
24 technology in terms of Keith's question.

25 DR. LANGE: Dr. Zuckerman.



1 DR. ZUCKERMAN: Yes. Thank you, Drs. Batchelor and Allen for asking the key  
2 question which is, why are we here today? And to a large extent, it revolves around a  
3 disagreement between the Sponsor and FDA regarding what is implied currently for the  
4 breakthrough device designation for this particular device.

5 At the beginning of this Panel meeting today, in the Sponsor's Slides 12 and 13 and  
6 in subsequent FDA slides, we all recognize that this device has a breakthrough designation  
7 as granted by FDA and therefore has the potential to be an important treatment modality  
8 and it benefits from increased interaction between FDA and the Sponsor.

9 As a result of the breakthrough designation and other recent guidance documents  
10 regarding uncertainty, there are several other factors to consider. Most importantly, the  
11 standard for a reasonable assurance of safety and effectiveness for a PMA device is not  
12 altered. However, in our consideration of whether this particular device is safe and  
13 effective, we do do our usual due diligence and we look at the benefit variables and the risk  
14 variables.

15 Now, by virtue of the breakthrough designation, it allows FDA and the Sponsor to  
16 deal with the uncertainty in these particular variables in a very flexible manner when  
17 appropriate. For example, by virtue of the breakthrough designation and guidance  
18 documents that go with it, there may be certain instances where the uncertainty in the  
19 benefit with the remaining uncertainty in the benefit and risk variables could be  
20 appropriately handled in the postmarket arena, and that's where the Sponsor is right now.

21 On the other hand, the FDA believes that given the uncertainty here for the benefit  
22 and risk variables, and the context with which we're dealing, meaning the patient  
23 population, known features of the device, etc., that this type of uncertainty would be better  
24 cleared up with a premarket randomized trial that has been referred to as COSIRA-II.

25 Regardless of what the Sponsor and FDA believe at this point, it's very important

1 that we hear now, independently from this Advisory Panel, what the Advisory Panel would  
2 recommend given that we have a breakthrough device with the potential for treatment of a  
3 very important population, and we have the COSIRA-I study.

4 So I hope, Dr. Batchelor, that with that preface you can help us out here.

5 DR. LANGE: So I've got Dr. Mathew, Dr. Yuh, and Dr. Vetovec and Dr. Yeh and  
6 Dr. Ohman.

7 DR. MATHEW: Verghese Mathew. Thank you, Dr. Lange.

8 I think we all agree that the pivotal trial hasn't been done, but that's what we'd like  
9 to see. Now, for whatever reason that couldn't be done the first time around in terms of  
10 sample size and the robustness of data collection and the difficulty in enrollment, so as I put  
11 on the hat of bringing therapies to patients that don't have good options, I guess my  
12 question is how likely is it going to succeed the second time around if that's the direction  
13 we're going, is that we really need that pivotal trial, whether it's premarket or postmarket.

14 And I would say that if it was postmarket, it would be virtually impossible, at least in  
15 my mind, to enroll in a trial where you approve a device and then ask people to randomize  
16 yet again. I think the path of least resistance is patients that would fit the bill for inclusion  
17 would probably just get this as a clinical device. So if there are other pathways to explore  
18 this in a registry-type fashion postmarket, I think that would be reasonable to discuss. But I  
19 think randomization after approval will be less successful than the first go-around.

20 Thank you.

21 DR. LANGE: I've got Dr. Yuh, Dr. Vetovec, Dr. Yeh, and Dr. Ohman.

22 DR. YUH: Thanks, Dr. Lange. David Yuh.

23 You know, I understand the good intentions for this breakthrough device rubric, but  
24 I'm concerned in a general sense that this is kind of a slippery slope where you can do a  
25 suboptimal study and then in a way work closely with the FDA to try to smooth out the

1 edges in an effort to get it approved. I just worry that we might be setting a precedent that  
2 could get out of hand in terms of new devices.

3 Many devices can be, you know, depending on how you define them, considered a  
4 breakthrough device, but that's just an overall concern I have at looking at the litany of  
5 flaws in this particular study. Even though that they probably represent some sort of  
6 positive signal for efficacy, it's still the quality of the study that bothers me and setting a  
7 precedent for accepting something like that through this mechanism.

8 DR. ZUCKERMAN: Dr. Yuh, that's why the FDA guidance document on uncertainty is  
9 very carefully written with the breakthrough devices guidance document so that if the  
10 probability of success of accomplishing something in the post-approval setting is not high,  
11 then that is not a recommended strategy. So again, I believe the guidances are written very  
12 consistent with good clinical judgment and as you've noted, we're looking for good clinical  
13 input here to solve a difficult situation.

14 DR. LANGE: I've got Dr. Vetovec, Dr. Yeh, Dr. Ohman, Dr. Gersh, and Dr. Hirshfeld,  
15 and then I'll close it after that.

16 DR. VETROVEC: The concern here, Dr. Zuckerman, I think this was very helpful for  
17 you to outline this, I'm not sure I've quite caught the nuances of this. I think we could say  
18 that this device seems to have a relative benefit, but the better approach, we might not  
19 have to randomize if we could come up with a really definitive mechanism of action so that  
20 we could measure that and show there was a difference with that, and it correlated with  
21 the symptom relief seen in a registry. So I don't know how flexible these decisions can be.  
22 So I think there's a small pathway to work our way through this. But the question then  
23 becomes if you do the postmarketing study and you do all I say and it turns out negative,  
24 are you going to take it off of the market?

25 DR. ZUCKERMAN: You're asking the FDA?

1 DR. VETROVEC: Yes.

2 DR. ZUCKERMAN: Again, I would first of all ask you to think about this from a clinical  
3 perspective. Do you have a high enough uncertainty right now that what would you  
4 recommend if that case actually occurred and how would it affect your decision making  
5 because that should be a factor, that benefit-risk paradigm that you just pointed out and  
6 was included in one of Dr. Raben's introductory slides this morning, so I must put it back to  
7 you as a Panel member.

8 DR. LANGE: In other words, it would come back to Panel at that point, George.

9 DR. VETROVEC: Okay. I guess my sense would be that I think there's a very narrow  
10 path that might --

11 (Audio feedback.)

12 DR. VETROVEC: -- despite all the limits, but I'm just not quite sure how that all plays  
13 out. But I think one could devise it but it would be very important to have some better  
14 understanding of how it affects ischemia or whether there's a gradient relationship or  
15 something else that you could correlate with the results.

16 DR. LANGE: Dr. Yeh.

17 DR. YEH: Thanks, Dr. Lange. Robert Yeh.

18 I was going to -- you know, when I put on my sort of clinical hat and I think about it, I  
19 think this is a really, really important patient population for which I am often sent patients  
20 and we struggle to treat these patients. So I understand the breakthrough designation, this  
21 is a really important patient population that has limited therapies, and when I think about  
22 some of the things that we do offer to these patients currently in clinical practice, I must  
23 say that many of the things that we offer these patients have -- are on equally thin evidence  
24 bases and I'm sure it would be riskier, I think. I direct our CTO and PCI program and I have  
25 to agree with one of the comments which is I think that we're on thinner evidence even

1 there and certainly those procedures are risky. So when I think about it and if this were to  
2 be approved, would I offer it to some of those patients? I think that I would, based on  
3 what's been presented.

4 At the same time, when I put my scientific hat on, I say do I feel confident in the  
5 effect size that's been demonstrated? I don't. Do I feel confident in a mechanism? I don't  
6 at all feel confident in the mechanism based on what's been presented. You know, do I feel  
7 good about safety? I do feel good about the safety and I don't think any additional  
8 premarket information is going to alleviate any further concerns about safety that were  
9 raised about intracoronary sinus dilation, that will all be -- need really long term follow-up  
10 and sort of registry types of data, also.

11 So at the end, I'm sort of mixed here because I think what the science has presented  
12 does not adequately convince me that this is a very efficacious device. At the same time,  
13 the patient population is a high-need population where we are grasping at different  
14 treatments. That's the sort of mixed feelings I have about this.

15 DR. LANGE: I'm going to ask Dr. Ohman, Dr. Gersh, Dr. Hirshfeld, and Dr. Connor.  
16 Dr. Ohman.

17 DR. OHMAN: Yes, thank you.

18 I think this helps a lot, Dr. Zuckerman, that you sort of laid it out in this regard, and  
19 being one who sees these patients on a very regular basis, I can concur that they seek, they  
20 desperately seek, things that work for them. But we also know that treating angina, there's  
21 a fair bit of placebo effect and many of these patients are convinced that certain things  
22 work and other things don't work. I think what we have to safeguard is that we understand  
23 the science and that we can actually back it up with objective evidence that show that it  
24 does indeed reduce ischemia. Otherwise, it will be very hard. If there would be a multitude  
25 of devices available through mechanisms that are not entirely straightforward, we would

1 face much harder clinical decision making than we do right now.

2 Finally, I would say on the breakthrough aspect of this, while these patients are  
3 suffering and they're having terrible quality of life, they're actually, as Dr. Henry pointed out  
4 in one of his earlier slides on mortality, yearly mortality is actually quite low and it's  
5 comparable to other patients with angina. So in a way, I can understand where we'd want  
6 to care and take care of the quality of life of our patients, I think it becomes less urgent as  
7 patients live and while they have poor quality of life, I think it's important for us to establish  
8 what really works.

9 DR. LANGE: Thank you.

10 Dr. Gersh, Dr. Hirshfeld, and then Dr. Connor. And then I'll summarize.

11 DR. GERSH: Thank you, thank you.

12 Basically, I'm in agreement. I agree, that I don't think safety is a real concern. I feel  
13 fairly comfortable with the data that have already been presented, but it is very difficult to  
14 me to make a decision about efficacy on the basis of a hundred and two randomized  
15 patients and many of the other issues we've heard. So I think we need a carefully designed  
16 randomized trial performed in the United States, blinded with a sham, and the endpoints to  
17 be determined as has been discussed earlier today.

18 DR. LANGE: Great. Thank you, sir.

19 Dr. Hirshfeld.

20 DR. HIRSHFELD: Yeah, I just wanted to mention I'm having a déjà vu moment here  
21 because it's just about 2 years ago that this group was here with the Impulse Dynamics  
22 Optimizer PMA and that was the first breakthrough device that this Panel had to deal with.  
23 And I recall very well that we were really concerned about the fact that this device did not  
24 have any objective data in terms of physiologic data that supported it and yet, because of  
25 the breakthrough device categorization, the Panel ultimately voted to approve it and it's on

1 the market. I think we're in a similar situation here that we have, at best, a small effect size  
2 demonstrated clinically and we have inadequate functional data to support that. And so I  
3 think, in terms of the direction in which we go, are we going to relax criteria for  
4 breakthrough devices or are we going to stick with the more rigorous for efficacy?

5 DR. LANGE: Dr. Connor.

6 DR. HIRSHFELD: And I've got Bram's attention.

7 DR. LANGE: I'm sorry. Dr. Connor and then Dr. Zuckerman.

8 DR. CONNOR: Yeah, this may be a question for my clinical colleagues, I almost asked  
9 this earlier and didn't, but given what Dr. Yeh said, I wonder at times like these, for  
10 instance, this is supposed to mimic the Beck procedure noninvasively, you know, we can't --  
11 FDA can't regulate the practice of medicine, it can only regulate what tools are on the  
12 market. So yeah, so I was wondering, given what Dr. Yeh said, kind of about what people  
13 do without this and is it -- you know, they do more risky things, they try anything because  
14 patients are willing to take risks to feel better and live their lives, like we've heard. So I was  
15 curious, like is the Beck procedure actually done these days? I guess that it can't be done  
16 for a lot of these patients, but I assume it can on some. So do people actually do that?

17 DR. LANGE: Not anymore, not anymore.

18 DR. CONNOR: Okay.

19 UNIDENTIFIED SPEAKER: No, that's not done.

20 DR. LANGE: Yeah. Dr. Starling.

21 DR. CONNOR: Good, thank you.

22 DR. LANGE: Dr. Starling and then I'll try to summarize so we can move on.

23 Dr. Starling. And then Dr. Zuckerman.

24 DR. STARLING: Thank you, Dr. Lange. Randall Starling.

25 And, you know, as a non-interventionalist, the point that I have to make is that as far

1 as breakthrough is concerned, I don't think if I can establish efficacy in my mind, I think  
2 safety, yes, but not efficacy. And what I worry about when we -- if we release something  
3 like this is that patients won't get proven therapies and it's very easy to default to  
4 something new versus pushing the envelope with what is already established. So this is  
5 what my concern is with respect to approval as breakthrough and then moving to a pivotal  
6 study. I'm in the camp that we need pivotal data now.

7 DR. LANGE: Okay. Dr. Zuckerman and then our patient rep would like to have a  
8 comment, as well.

9 So Dr. Zuckerman.

10 DR. ZUCKERMAN: Yes. I want to, first of all, thank Dr. Hirshfeld for his comments  
11 because, as he points out, we're trying to balance a multitude of factors. Number one, we  
12 would like, if at all possible, U.S. patients to have timely access to important new therapies.  
13 On the other hand, by virtue of our regulations, we do want to make sure that any  
14 approved devices have appropriate safety and effectiveness. This often can be a difficult  
15 situation and that's why we need an expert advisory panel, as we have assembled here  
16 today, to help us make this decision by giving us their best advice possible.

17 The best advice, though, that I would give to this Panel is found in a paragraph on  
18 page 16 of our uncertainty guidance, which is the paragraph right before the slide example  
19 that the Sponsor gave this morning. And let me read directly from the paragraph because I  
20 think it will help Dr. Hirshfeld and others.

21 "Finally, as noted above, the decision as to whether or not a device provides a  
22 reasonable assurance of safety and effectiveness is based on the totality of the valid  
23 scientific evidence, including clinical studies and nonclinical testing. The appropriate extent  
24 of uncertainty of benefits and risks in a given case will depend on considerations of the  
25 factors set forth in Section IV (e.g., the disease or condition at issue, the availability of



1 alternative products, and risk mitigations) and other relevant information concerning the  
2 device. We anticipate that the greatest extent of uncertainty discussed in the examples  
3 below would only be appropriate under rare circumstances and, in any case, the sponsor  
4 must show, among other things, that the totality of the valid scientific evidence provides a  
5 reasonable assurance of safety and effectiveness of the device."

6 So I would like to ask Dr. Hirshfeld and any other Panel member if this helps them  
7 put into appropriate context how to employ the fact that the device has been designated as  
8 a breakthrough device, but yet the regulatory bar still remains a reasonable assurance of  
9 safety and effectiveness.

10 Dr. Hirshfeld.

11 DR. HIRSHFELD: I agree. I think we should be very careful not to allow that  
12 designation to take our eye off the safety and effectiveness ball.

13 DR. LANGE: Thank you.

14 Jacqueline Alikhaani, you have some thoughts about this, and then I'll go to  
15 Dr. Batchelor. Jackie, you're on mute right now.

16 MS. ALIKHAANI: Right. I think that, as a consumer, as a patient, given the great  
17 discussion we've had here today, there's so much missing information and discrepancies  
18 related to that information that it just makes me feel uncomfortable, I mean, as a patient.  
19 Patients, when they get the consent form, they don't know about all this discussion we're  
20 having, you know, relative to breakthrough devices and whatnot. I still think we have to  
21 keep our eye on the ball. I think we need that new trial, we need it in the United States, we  
22 need to get it right so we can have all the evidence and there won't be any need for --

23 (Audio feedback.)

24 MS. ALIKHAANI: -- pressure from the doctors to make these decisions and trying to  
25 make the best safest decision for their patient with a need. Just let's just do the trial, that's

1 what I'd like to see.

2 DR. LANGE: Thank you.

3 MS. ALIKHAANI: And get the trial right this time.

4 DR. LANGE: Thank you, Jackie. Thanks for those comments.

5 Dr. Batchelor.

6 DR. BATCHELOR: I want to just start off by saying those are very sage comments,  
7 Jackie, and I'm actually going to really pound home kind of the same points. The way I look  
8 at this is what is the risk to society and our patients in not doing -- in letting this go through  
9 now and get approved and then having a real challenge doing a randomized trial, because I  
10 think to Dr. Mathew's point, it's going to very hard to do a randomized trial that's going to  
11 be really effective after this, after approval, and versus the risk of delay.

12 If we delay, I just don't think we're putting patients at risk in terms of death or  
13 serious complications other than just a lot of angina, but they've had angina for a long time  
14 prior to this and honestly, I look at this a little bit differently in terms of breakthrough  
15 technology. So if this was a breakthrough technology that was reducing myocardial  
16 infarction or death, it would be looked at, in my opinion, in a different light. But this is a  
17 breakthrough technology that's going to make people feel better but their annual mortality  
18 is still going to be somewhere between 2 and 4 to 5%, which is reasonable.

19 So if we delay, to me, it's not a huge jeopardy that we're exposing patients to and I  
20 think it's worthwhile because the science we'll get out of it will be much better and we'll be  
21 able to look Jacqueline in the eye and our patients in the eye and say we know this works  
22 and we're not going down a slippery slope where we're letting this creep in and then we're  
23 going to have a really hard time with a randomized trial later. So my advice is get it right  
24 and to hold off. I'm an interventional cardiologist, I'd love to do these devices, but those  
25 are my thoughts.

1 DR. WITTES: May I say something? This is Janet. Can you hear me?

2 DR. LANGE: Janet?

3 DR. WITTES: Yeah. Can I say something?

4 DR. LANGE: Who is Janet?

5 DR. HIRSHFELD: Janet Wittes.

6 DR. WITTES: Janet Wittes.

7 DR. LANGE: I'm sorry?

8 DR. WITTES: Never mind, okay.

9 DR. LANGE: I didn't see your voice.

10 (Crosstalk.)

11 DR. WITTES: That's all right, forget it.

12 DR. LANGE: Go ahead, Dr. Wittes, my apologies.

13 DR. WITTES: Okay. What I was going to say is if I believe if the data were convincing  
14 to me that this works for quality of life and symptoms, I would be very strongly favoring  
15 recommending approval of it. The fact that the mortality rate is low is not, to me,  
16 dispositive. What makes a difference to me is an inability, all the other things we've talked  
17 about that make it not clear that what we think we're seeing in the primary outcome is  
18 actually reproducible.

19 DR. LANGE: Great. All right, now I'm going to bring this to -- I'm going to bring this  
20 to a close, we have other questions.

21 Debra, I will ask for your opinion and Gary's at the very end before we vote, okay?  
22 But just because we've had a lot of discussion and we still have more to talk about, I think I  
23 can give you a show hands right now based upon the totality of information. The question  
24 is do we need additional premarket data from a randomized sham-controlled clinical study?  
25 And I'm just going to ask for a show of hands. Those that say yes, raise your hands.

1 (Show of hands.)

2 DR. LANGE: Those who say no, raise your hands.

3 (Show of hands.)

4 DR. LANGE: And there's one no and the rest are affirmative. Now, we can have  
5 additional discussion, we can talk about safety and efficacy, but Dr. Zuckerman, does this  
6 give you the direction and the opinion of the Panel right now?

7 DR. ZUCKERMAN: Yes, it does. Thank you.

8 DR. LANGE: Okay, all right. I appreciate that really robust discussion. Normally,  
9 right now I'd be looking around the room because you all are trying to get out to catch your  
10 flights, but I'm looking around, there's nobody in an airport, so I got you for a little bit  
11 longer, okay?

12 MS. ASEFA: I think Debra Dunn said she had a work commitment to keep and that's  
13 kind of -- that's why she wanted to speak --

14 DR. LANGE: I'm sorry. Debra, unmute yourself. At the very end of these questions  
15 I'll obviously call upon you and Jackie and Gary, but let's get your opinion right now based  
16 upon what you've heard.

17 MS. DUNN: Hear me?

18 DR. LANGE: Yes, ma'am.

19 MS. DUNN: Can you hear me? My name is Debra Dunn, I'm a patient advocate, and  
20 everything that we've been talking about here today really hits home. I'm a heart failure  
21 patient and during a simple device switch out, device number seven, my biventricular ICD, I  
22 got a staph infection in the pocket. So I had lead extraction and I had a cardiac tamponade.  
23 I was left with a septal opening of one inch and I received the Amplatzer device after the  
24 staph infection was taken care of, and I have been an expert witness for the FDA on the  
25 Amplatzer, which recycled back around, I'm not sure if we talked about that that was a

1 device that obviously came before a panel. I was not on that panel but I was an expert  
2 witness when it came back around because it was fragmenting and so I did have to present  
3 for that.

4 Everything here has been taken very seriously. I do live with angina, I'm a serious  
5 heart failure patient, and I respect the Panel, that they would like to take this back under  
6 advisement, get more clinical data to make a wise decision, but I can tell you living with  
7 angina every day as a patient, it's very unsettling and it's very scary.

8 But I guess when it all comes down to it, as a patient, we want to make sure that  
9 we're safe and whatever's being put in our body. My Amplatzer, who knows what will  
10 happen? But it was something that at the time it saved my life. So I'm 18 years out from  
11 my first heart attack and I'm doing great, I'm getting device number seven switched out and  
12 a new one put in real soon, so thank you for all your consideration.

13 DR. LANGE: Debra, thanks for that heartfelt comment, thanks for the perspective,  
14 we appreciate that. I appreciate you joining the Panel. And I can tell you, everybody on the  
15 Panel wishes you well with your next device and I hope you do, so thanks so much for  
16 sharing with us.

17 MS. DUNN: Thank you very much.

18 DR. LANGE: Sam.

19 DR. RABEN: So going to the next question, regarding patient demographics and  
20 other baseline characteristics, the populations were similar between the treatment groups.  
21 The average age of subjects was 67.8 years and ranged from 35 to 87 years old. The  
22 majority of subjects (80.8%) were male and white (86.5%). The groups had comparable  
23 heart rates and blood pressure. However, the study included a limited number of female  
24 (19.2%) and minority (5.8%) patients.

25 Question 5b is, the demographics of the patients enrolled in the COSIRA trial had

1 differences compared to the U.S. refractory angina population (i.e. no black or Hispanic  
2 patients enrolled and under-representation of females). Please discuss the applicability of  
3 the study results to the U.S. refractory angina population and whether there's a need for  
4 additional clinical data of the safety and effectiveness of the Neovasc Reducer device in a  
5 more demographically representative population.

6 DR. LANGE: I'm going to summarize that because we've talked about that, several  
7 individuals have, and Jackie, I'm going to give you credit, I don't think anybody could've said  
8 it any better than you. I think as difficult as it may be to get a more demographically  
9 diverse representative population, I think the feeling of the counsel of the Panel is that's  
10 absolutely necessary. If there's any disagreement -- John, how could you disagree with me?  
11 We've been working together for 15 years.

12 (Crosstalk.)

13 DR. SOMBERG: -- disagree. No, no, no.

14 DR. LANGE: Go ahead, John. Then I think --

15 DR. SOMBERG: I just want to say that it's advisable to have a heterogeneous  
16 population, but it's not been shown that Hispanics or people of African American descent  
17 have a more physiologic difference here. I mean, I think the important thing would be to  
18 include women and to have more people -- you know, not exclude microvascular disease. I  
19 think that is the clinically relevant discussion.

20 DR. LANGE: Okay. Dr. Borer and then Dr. Cigarroa and then Dr. Batchelor.

21 Jeff, you're on mute.

22 DR. BORER: I agree with John. I think the issue is we have other more important  
23 concerns than whether we have enough black or Hispanic patients or even enough women.  
24 The fact is that we have what we have. If the study were adequate otherwise, which it's  
25 not, then we could write a label that says these data apply to the population that was

1 studied, which was deficient in black or Hispanic patients and females. But if another trial is  
2 to be done, every effort has to be made to enhance the representation of black, Hispanic,  
3 and women patients. You can't be sure you're going to be able to find them. They couldn't  
4 in the first trial. So I think that this is very important, but it's not a drop-dead issue.

5 DR. LANGE: I'll go to Dr. Cigarroa and Dr. Batchelor. But Jeff, I'll remind you that the  
6 trial was done in Europe.

7 DR. BORER: Yeah, that's fine.

8 DR. LANGE: So Dr. Cigarroa, Dr. Batchelor, and Dr. Zuckerman.

9 DR. CIGARROA: I think this is a -- sorry, this is Joaquin Cigarroa for the record.

10 So I think the applicability of the study results to our refractory angina population is  
11 challenging to extrapolate. I do believe that there's need for additional clinical data for  
12 effectiveness in a demographically representative population.

13 If one looks at African Americans as an example and looks at LV mass and capillary  
14 density, which is really important in impacting, again, reversal of flow from the  
15 subendocardial to the epicardial, especially as left ventricular diastolic pressure increases,  
16 that difference in LV mass in a hypertensive population may be substantial and may  
17 adversely impact the efficacy of this, the purported mechanism the Sponsor has stated. So  
18 those differences are important.

19 We can say the same for women. Just because one has epicardial coronary artery  
20 disease does not decrease, per se, the commingling of microvascular disease and in fact,  
21 they often commingle. So I would state that this is an important issue that must be  
22 addressed prospectively.

23 DR. LANGE: Great. Dr. Batchelor.

24 Thank you, Dr. Cigarroa.

25 Dr. Batchelor.

1 DR. BATCHELOR: Thank you. Again, this is off of what was just said by Dr. Cigarroa  
2 very nicely. There are symptomatic differences in the way that patients present with angina  
3 and some African Americans are more likely to have shortness of breath and probably  
4 because of the superimposed hypertension and LV diastolic function issue, women  
5 obviously tend to present older and are more likely to present with atypical symptoms.

6 So I think Dr. Cigarroa also talked about some physiologic differences that might  
7 exist. I think we have to get really careful when you state there's never been a difference in  
8 the efficacy of a medication and/or an intervention across racial or ethnic groups because it  
9 really hasn't been studied. So to say that there's no difference doesn't mean that there's no  
10 difference, it's just that we don't know. We do know that there are -- we do tend to come  
11 to the table and to the clinician with slightly different characteristics depending on many  
12 factors including race, age, and gender so I would just caution us on that assumption,  
13 number one.

14 And then finally, remember, with the Food and Drug Administration Safety and  
15 Innovation Act from 2012 and what's come about from that, I think there's a responsibility  
16 now, I think we've got to start into a new phase of conversation where it's been signed into  
17 law and we have a responsibility to try to at least do our best to represent various ethnic  
18 and racial groups and women in clinical trials, otherwise we're just speaking into the air  
19 without really having anyone listen to us.

20 DR. LANGE: So Dr. Zuckerman, I think what you're hearing is we're not certain  
21 whether there are mechanistically different reasons, we're not sure whether different  
22 gender or different racial or ethnicities will respond differently, but we will never know  
23 unless they're studied and a study outside of the U.S. is a population that is different from  
24 us may not answer that question. So I think that the Panel overall would say -- I don't think  
25 anybody on the Panel would say we need less diversity in this study.



1 DR. ZUCKERMAN: So this has been an excellent discussion and I think it's also  
2 important, too, for the panelists to understand what the FDA regulations are. Number one,  
3 America is a heterogeneous country and our FDA regulations specifically state that if we're  
4 going to use outside U.S. data to make a decision about safety and effectiveness for a U.S.  
5 population, then that OUS population data really needs to be extrapolate-able to the U.S.  
6 situation and I think we've heard in this discussion some questions about that.

7 DR. LANGE: Great, thank you.

8 All right, Sam, may I have you move on?

9 DR. RABEN: Great. Question 5c. Acknowledging that an understanding of the  
10 Reducer's mechanism of action is not required for PMA approval, please discuss the  
11 principal data supporting the intended clinical benefit in your assessment of the strengths  
12 and limitations of the data supporting device effectiveness

13 If you recommend additional premarket data to support a reasonable assurance of  
14 safety and effectiveness of the Reducer device, please describe the types of studies (e.g.,  
15 animal or human) that would be most helpful. Please comment on and make  
16 recommendations regarding whether the recommended data could be obtained using a  
17 protocol similar to the COSIRA-II trial.

18 DR. LANGE: Dr. Allen.

19 DR. ALLEN: So I think, as Dr. Zuckerman and some of his esteemed colleagues at the  
20 FDA pointed out earlier, determining a mechanism doesn't -- it may enhance the benefit of  
21 getting an approved device, but it's not mandatory and I'm not as concerned about not  
22 having great data on mechanism as I am on the clinical trial design and how it was executed,  
23 that bothers me a lot more. If this had been a great trial, for example, this had been  
24 COSIRA-II, I probably wouldn't have minded not knowing what the mechanism is. I don't  
25 really how know aspirin works but I take it every day. So that would be my take on this.

1 DR. LANGE: Dr. Somberg.

2 DR. SOMBERG: If the Sponsor chose a feel-better primary endpoint, it would be  
3 important to have a secondary endpoint of ischemia measured either by treadmill,  
4 dobutamine echo or something else. The mechanism isn't critical. But to know that we're  
5 dealing with ischemia as an underlying improvement of symptoms is, I believe, very  
6 important.

7 DR. LANGE: Dr. Page and then Dr. Gersh.

8 DR. PAGE: Yeah, I agree with what was said about if it clearly works and it's clearly  
9 safe and we don't fully understand why, I can live with that. When we're on the cusp, I  
10 would like to have a better -- a fuller understanding of the scientific plausibility behind this  
11 and would hope that in the upcoming pivotal study, if that were to be conducted, we would  
12 have a better understanding of the anatomy, the shape of the human coronary sinus with  
13 this mature device in place, and a better understanding about truly the physiology  
14 underlying whatever response we're getting.

15 DR. LANGE: Great. Dr. Gersh.

16 Thank you, Dr. Page.

17 DR. GERSH: I think the pivotal study, if it's going to be done, should be in people  
18 with obstructive coronary disease and where I think mechanisms may be important is, I  
19 mean, in many ways the concept if this could work for microvascular angina is quite  
20 attractive. But in order to look at that, long before we do a randomized trial we have to  
21 identify the people with endothelial function studies and physiological hemodynamic  
22 studies at cath and then look at it in that group of patients completely where the benefits  
23 and the mechanisms may be different, but first of all we'll have to identify them and then  
24 think about mechanisms. But I think we mustn't lump together microvascular angina with  
25 microvascular dysfunction that is probably present in most people with obstructive disease.

1 I think of them as two very different entities and it may be that the mechanisms of benefit,  
2 if they're there, may vary between those who have non-obstructive and obstructive disease.

3 DR. LANGE: Dr. Brindis, do you have an opinion?

4 (Off microphone response.)

5 DR. LANGE: No. Dr. Mathew?

6 DR. MATHEW: I'd agree that mechanistic studies may be important if we're going to  
7 do a pivotal trial because I think we had the conversation earlier about subsets of patients  
8 that are going to be more likely to benefit and I think we still may have that issue to  
9 contend with if there's a second pivotal trial. So I think that the trouble of going through  
10 another well-designed and well-executed trial is going to go forward, I think mechanism  
11 should be some subset of that population.

12 DR. LANGE: Okay, thank you.

13 Dr. Starling.

14 DR. STARLING: Yeah, I think a mechanistic is an option that should be a secondary  
15 and it's going to get down to funding and resources. I think the hard endpoint that's  
16 described from the Sponsor for COSIRA-II of an exercise treadmill and an increase in 60  
17 seconds in exercise time is an appropriate primary endpoint.

18 DR. LANGE: Okay. Dr. Yuh. You're muted, David. I'm sorry.

19 DR. YUH: I don't really have anything useful to add to the previous comments.

20 DR. LANGE: Okay. And so I'll ask -- I'll call the last person, Dr. Allen. And then I'll call  
21 on Dr. Batchelor.

22 Keith, you need to unmute yourself. Any comments?

23 DR. ALLEN: I really don't have anything to add.

24 DR. LANGE: Okay. Dr. Batchelor, any opinion on this?

25 DR. BATCHELOR: I think the previous speakers have really laid this out nicely.

1 Noting else.

2 DR. LANGE: So what I'm hearing, Bram, is it certainly needs to be a robust clinical  
3 study and the values of animal studies would be helping to identify subsets or providing  
4 mechanisms when the results of the pivotal clinical study are "on the bubble," as Dr. Page  
5 mentioned, in which case the mechanism study would kind of push them over the finish  
6 line.

7 DR. ZUCKERMAN: Thank you.

8 DR. LANGE: Great. Thanks for everybody's opinion on that. Sam.

9 DR. RABEN: Okay. Question Number 6 is regarding the benefit and risk. Given the  
10 totality of the evidence regarding the effectiveness and safety profile for the device, please  
11 comment on the benefit-risk profile for this device.

12 DR. LANGE: Do we want to do that or do we just want to go straight to safety and  
13 effectiveness at this point? I think we've talked a lot about this already previously or do you  
14 want to -- do you feel comfortable with that, Dr. Zuckerman, or do you want us to address  
15 this?

16 DR. ZUCKERMAN: If there are any brief Panel comments regarding the benefit-risk  
17 ratio that Panel members would like to make for the open public record, this is the time  
18 because subsequently they will need to vote on these questions. So if they'd like to bring  
19 up any points, comments, it's appropriate at this point.

20 DR. LANGE: Okay. And I see Dr. Page would like to, Dr. Mathew would like to.

21 DR. PAGE: Yeah, we talk about this being a relatively low-risk procedure and as  
22 invasive cardiac procedures go, it is. But it's not a zero risk. The embolization in real-world  
23 is reported. And the other thing I think we should keep in mind that, just as the very  
24 compelling statement from that patient who's been through various procedures including  
25 lead extractions and the like, some of these patients are going to need a biventricular

1 pacemaker and we don't have any data as to accessing the coronary sinus with a  
2 biventricular lead. We could do an intervention to open it up that undoes the effect of this  
3 device, perhaps.

4 So I think risk is low, risk would be acceptable for a device that was clearly showing  
5 effectiveness, but it's not zero risk. Any of the invasive people on the line know any time  
6 you get invasive something bad can happen.

7 DR. LANGE: Okay. Dr. Mathew.

8 DR. MATHEW: Yeah, I generally concur. I mean, I think this is relatively safe, it  
9 sounds like there's a learning curve from the data they presented but it's relatively safe, I  
10 presume, if you proctor people through this for a case or two. So I think the denominator's  
11 low, I think the benefit is modest at best and by that function, the ratio probably is  
12 favorable but I think we all agree that the benefit is modest at best based on the data we've  
13 seen so far.

14 DR. LANGE: Okay. Dr. Allen. And then Dr. Cigarroa.

15 DR. ALLEN: So I guess I agree with what mostly has been said, but I really think we  
16 focus more on the acute risk of the procedure and while it's an invasive procedure, I think  
17 the acute risk is low. I really don't know what the long-term risk of this implantable device  
18 is and 6 months seems to me a very short window for an implanted intravascular device to  
19 determine its safety even in the intermediate term.

20 DR. LANGE: Okay. Dr. Cigarroa and then Dr. Yuh and then Dr. Somberg.

21 DR. CIGARROA: This is Joaquin Cigarroa for the record.

22 So as I think about the totality of the evidence, effectiveness and safety profile, the  
23 benefit-risk profile to me is impacted by what is the endpoint, is it death, is it infarct or  
24 stroke? Or is it, in this case, a symptom which clearly is important but is not death or  
25 myocardial infarction?

1           And therefore when I think about safety, I also then say well, is this effective and as  
2 we have all been discussing, the degree of uncertainty regarding effectiveness, at least in  
3 my mind, is on a low, medium or high, is substantially high. And therefore even though the  
4 short-term safety appears reasonable, it's hard for me to believe that the benefit-risk  
5 profile with the data we know and the endpoint that has been chosen is favorable without  
6 additional insights and without additional high-quality data.

7           DR. LANGE: Dr. Yuh.

8           DR. YUH: Thank you, Dr. Lange. David Yuh.

9           Yeah, I think looking at the profile of the adverse events for this particular device,  
10 I'm pretty satisfied given certain factors and the patient population we're looking at in  
11 terms of characterizing the risk of the device, at least in the short term, and I think I feel  
12 that I would feel confident in giving informed consent to the patients with all the caveats,  
13 patients that are truly looking for a last-ditch alternative for a very debilitating condition, so  
14 that's my take on the safety profile. Under other circumstances I might be more -- have  
15 more stringent requirements, but I think I'm satisfied with the overall picture of safety for  
16 this device.

17          DR. LANGE: Dr. Somberg.

18          DR. SOMBERG: I just want to reiterate what I said earlier about the FDA should  
19 question the use of antiplatelet agents with this device, it sits on the venous side, the  
20 question is why not antithrombotic agents and also the possibility of late stent thrombosis.  
21 One of the problems with coronary sinus in this situation is if you had a coronary sinus  
22 thrombosis, it often presents as angina or myocardial infarction and that would be on the  
23 continuum of the natural progression of these patients' disease. So how do you  
24 differentiate that? I think some thought should be given, an appropriate therapy  
25 introduced early, maybe to anti-coagulate patients for a period of time while they

1 endothelialize the stent.

2 DR. LANGE: Our Consumer Rep, Jackie Alikhaani, has asked to speak. Jackie. And  
3 then I'll got to Dr. Yeh. You're muted, Jackie. Thank you.

4 MS. ALIKHAANI: Yes. Yes, like was already said, I think that every procedure has a  
5 risk, we all know that, and I think that patients -- I mean, I don't think it's possible to do a  
6 really good benefit-risk assessment without having all the information we need. We don't  
7 have all the evidence because we don't -- you know, it wasn't there in the trial, that's what  
8 we've been talking about for hours.

9 So patients, even the patients who are having a lot of quality of life issues and a lot  
10 of pain, I think people are depending on us to do really due diligence in terms of getting all  
11 that evidence, getting everything that we can and just lay that out for patients to know  
12 about and they really can make informed decisions based on the evidence that is given to  
13 them so people won't be making decisions in the dark and later on, when something goes  
14 wrong, then who's going to be liable for that? I think we should do everything we can to  
15 get all the information we can that -- it's just worth it.

16 DR. LANGE: Thank you, Ms. Alikhaani.

17 Dr. Yeh.

18 DR. YEH: Thank you, Robert Yeh for the record.

19 I think that, you know, there was a comment about the short-term safety of the  
20 device and 6 months being too short, I do think that the safety evaluation probably is the  
21 one thing that will be based on single-arm type of data and for those data, I do feel  
22 confident that those long-term data that were presented from -- the one thing that maybe  
23 the Reducer single-arm study might be good for is to show that weren't sort of a large  
24 number of late complications that were unanticipated at, say, the multi-year -- you know, 3  
25 to 5-year period.

1 But when I think about, you know, if this were the only data that we were ever going  
2 to get, it feels to me that the device is safe and if I were betting on it from just a pure  
3 probability Bayesian standpoint, I would say that the evidence suggests that the device has  
4 probably a greater chance than not of being effective. And so based on just a point  
5 estimate, I think that effectiveness probably exceeds -- the benefit probably exceeds the  
6 risk, but I don't think we're choosing on the basis, I don't think we make this decision on the  
7 basis of point estimates especially when there are wide confidence intervals around all of  
8 these assessments.

9 And so even though I think that the benefit-risk profile does favor the device, I think  
10 there's enough uncertainty here and that we're probably making decisions not based on  
11 point estimates, but we probably should make the decision based on how confident we are  
12 in those point estimates. And I think there's not enough certainty here to say definitively  
13 that the benefit-risk profile is greater than zero.

14 DR. LANGE: Great. And I'll add one final comment because I'm not -- my comments  
15 have been consistent with what everybody else has said. I don't feel like we have adequate  
16 safety data, there's no secondary endpoints in a small group. We have, again, 2000 patients  
17 in long-term follow-up, we assess they don't have coronary sinus thrombosis because no  
18 patient reported it. We don't have really adequate long, good follow-up and it's not been  
19 very rigorous, so I can't really can't assess -- I can't say it's safe. I'm not saying it's harmful, I  
20 just can't -- I don't think there's enough data to suggest it's safe. That's just my opinion.

21 So we've addressed Question 6, I think. Bram, you've heard from the Committee, in  
22 toto they don't feel like the benefit-risk profile of the device is favorable primarily because  
23 of lack of confidence in the benefit. Most Committee members feel like the risk is low, I'm  
24 probably on the outside of that bubble.

25 DR. ZUCKERMAN: Thank you.



1 DR. LANGE: Great.

2 DR. RABEN: So next, to discuss the proposed post-approval study. In response to  
3 the concerns identified by FDA during the initial rounds of review, the Sponsor has  
4 proposed the following potential post-approval study: Neovasc has committed to do a  
5 post-approval randomized, double-blind, sham-controlled study in a country where the  
6 Reducer is not approved to allow the collection of data to reduce the amount of remaining  
7 uncertainty FDA may have.

8 Note: This requested discussion item related to the proposed post-approval study  
9 should not be interpreted to mean that FDA has made a decision or is making a  
10 recommendation on the approvability of this PMA.

11 FDA may require a post-approval study (or studies) at the time of approval of a PMA  
12 to provide information on the continued safety and effectiveness of the approved device.  
13 These studies are not intended to provide a reasonable assurance of safety and  
14 effectiveness, as the determination must be established prior to device approval, and are  
15 typically not randomized. FDA is concerned that after making a determination that there is  
16 a reasonable assurance that the device is safe and effective in a limited option patient  
17 population, it may not be appropriate or feasible, due to the lack of clinical equipoise, to  
18 mandate a new trial as a condition of approval in which patients who lack alternative  
19 treatments would be randomized to a sham control.

20 So for Question 7, please discuss and make recommendations regarding the  
21 Sponsor's proposal to perform a post-approval randomized sham-controlled trial. Please  
22 also discuss what alternative postmarket approval studies could provide the data needed to  
23 support this device.

24 DR. LANGE: So Bram, I think that the Panel addressed this previously and that is all  
25 except for one member thought that a premarket approval study was necessary to move

1 forward with approval. And so based upon that, I'll speak for the Panel, anybody's willing to  
2 disagree, there's no postmarket study approval, again, except with one member voting  
3 otherwise that would give the assurance of effectiveness and safety that would warrant us  
4 to recommend FDA approval.

5 DR. ZUCKERMAN: Thank you, Dr. Lange.

6 DR. LANGE: Dr. Brindis.

7 DR. BRINDIS: Yeah, so you know, I agree with everything that you just said, Richard,  
8 but if -- I'm still puzzled by the issue of a breakthrough technology and despite reading the  
9 definition just put up and also Bram's sage comments on how we should look at it. So I'm  
10 wondering if you could have a setup with the company, a preapproval study which allows  
11 this so-called breakthrough technology that allows at some point the patient being given  
12 the option to roll over to the active group, sort of I think how it was implied in the trial set  
13 up -- proposed by the company. And I would not -- in terms of this post-approval study, for  
14 all the reasons, particularly well articulated by Wayne, I would like to have that here, not in  
15 another country.

16 DR. ZUCKERMAN: Okay. Dr. Brindis, the question is really centered on there is  
17 uncertainty right now and the need for additional data and is it realistic to expect that these  
18 additional data to decrease uncertainty or, as Dr. Yeh pointed out, to decrease the width of  
19 the confidence interval, should it be -- could it be realistically done in the post-approval  
20 setting or in the premarket setting. And I think, as Dr. Lange has summarized, unless other  
21 panelists disagree, that realistically it needs to be done in the premarket setting.

22 DR. BRINDIS: Right.

23 DR. ZUCKERMAN: What you brought up is a different question that the Agency can  
24 discuss with the Sponsor at a different time, which is after a certain -- if a randomized,  
25 controlled trial is being conducted, are there certain criteria where patients who might be

1 randomized to the control group and are failing standard treatment be crossed over, and  
2 that's always an important consideration that can be discussed with the Sponsor, you  
3 pointed out some reasons to consider it, but it's really independent from this question.

4 DR. BRINDIS: I agree.

5 DR. ZUCKERMAN: Thank you.

6 DR. LANGE: I think that may be the last question prior to moving to safety and  
7 effectiveness, is that correct?

8 DR. RABEN: That is the last question, but then we will move to summations.

9 DR. LANGE: Absolutely. So at this time, the Panel will hear summations, comments  
10 or clarifications from the FDA. The FDA will have 10 minutes. Following that, the Panel will  
11 hear summaries, comments, and clarifications from the Sponsor and the Sponsor will also  
12 have 10 minutes. So with that, let me turn over to the FDA first.

13 DR. RABEN: Great, thank you.

14 Today's Panel meeting has focused on the clinical data and supplementary  
15 information provided by the Sponsor to support the Neovasc Reducer System. We  
16 discussed the regulatory history of the Reducer device and the regulatory considerations  
17 regarding FDA's Breakthrough Devices Program and uncertainty guidance. FDA explained  
18 that the breakthrough device program is intended to provide the American public with  
19 timely access to new beneficial devices with the potential for significant impact but does  
20 not modify or reduce the statutory requirements for device approval. That is, the Reducer  
21 device must demonstrate a reasonable assurance of safety and effectiveness.

22 We particularly appreciate the Panel's view on the strengths and limitations of the  
23 data including whether the COSIRA patients had refractory angina despite an optimal  
24 antianginal therapy, study blinding, the amount and imbalance of missing data,  
25 identification of patients most likely to benefit, the device's mechanism of action, the

1 importance of anginal symptoms and objective ischemia endpoints, treatment effect size,  
2 and lastly, enrolled study subjects' demographics distributions, compared to the U.S. We  
3 thank the Panel for their attention and discussion and for your careful consideration of this  
4 PMA.

5 DR. LANGE: Great. Thank you very much from the FDA both for your summary and  
6 for an excellent presentation at the beginning of the meeting.

7 And now for the Sponsor, 10 minutes as well for comments.

8 MS. BEBEAU: Thank you, Dr. Lange. I would like Dr. Stone to provide some  
9 commentary, as well as Dr. Henry.

10 DR. STONE: Well, thank you very much, this has been a long day and an excellent  
11 discussion. I think most of the points that have been made are actually points that we  
12 realize. You know, I've led over 120, maybe 130 randomized trials in my career and I was  
13 not involved in COSIRA, so I guess you might ask so why I am sitting here today helping the  
14 Sponsor with this data with my interpretation and analysis? And in that regard, I want to  
15 bring it back to the patient and this is really a desperate patient cohort. There's plenty of  
16 patients that, as you've heard from testimonials, including the patient advocate, are really  
17 leading unhappy lives at best and sometimes desperate lives at worst because of severe  
18 ongoing angina with minimal activities or no activities.

19 And so to have a therapy that may benefit these patients, I think myself as a  
20 cardiologist who's taken care of these patients for 30 years, is very important. When you  
21 look at, again, the options for patients who have exhausted coronary revascularization and  
22 antianginal medications, there are few. In fact, if you look at the randomized trial at EECF,  
23 it wasn't studied in Class IV patients and two-thirds of the patients had Class I or II. So then  
24 we're really left with surgical TMR or heart transplantation, which are really both therapies  
25 much more beyond that do have substantial safety concerns.

1           So for me, when I've got a patient in front of me that really needs something, the  
2 first thing I'm going to ask is okay, is it a safe therapy and I really do think that there's a lot  
3 of data to suggest that there's safety with this device. The procedure itself for an  
4 experienced interventional cardiologist is not that much more than putting in a stent in a  
5 normal vein. It's through a central line approach which residents do every single day, let  
6 alone cardiologists, and it really can be trained quite easily. And when you look at the bulk  
7 of the safety data, I mean, the SAEs in the COSIRA trial tended to favor the device arm, the  
8 MACE events and may (ph.) events tended to favor the device arm and you saw a training  
9 effect over time from experience to experience. So this really does, to me, seem like a safe  
10 device.

11           So it really does come down to effectiveness, and I understand the limitations of the  
12 dataset that you've seen, but I would say that angina is what matters most to these  
13 patients, I think it was Dr. Yeh that said it. It's not ischemia, it's not exercise duration of a  
14 few extra seconds, it's angina. And a two class or more reduction in CCS angina class is a  
15 robust reduction. I mean, there's really no misinterpreting that. And the Seattle Angina  
16 Questionnaire tracked along perfectly with that. In fact, there was a statistically significant  
17 improvement in Seattle angina quality.

18           The p-value, if you will, the alpha for the primary pre-specified endpoint was 0.2, so  
19 there's only a 2% chance that that's a false positive finding. It's not from 0.01, but it's  
20 pretty good for the pre-specified primary endpoint.

21           The number needed to treat was only about three to five which, for a relatively safe  
22 device, especially to get a two-class reduction in CCS angina is quite low and much better  
23 than antianginal drugs.

24           So then the last point I'll make is yes, the REDUCER-I data, of course, is open label  
25 and the five or six investigator-sponsor studies are open label and without concurrent

1 control arm. So yes, there's much more limited implications you can draw from that data.  
2 But the fact that they are consistent with the COSIRA study in terms of angina improvement  
3 and Seattle Angina Questionnaire improvement, I do think is some supportive data in terms  
4 of effectiveness.

5 So the question I guess I would ask is, if you have a patient in front of you, for the  
6 cardiologists here, or if you have a family member that really had disabling angina and you  
7 were to have this potentially relatively safe therapy available to them, would you want to  
8 offer it to them?

9 I do believe that we can do a post-approval trial. Not all patients in the United  
10 States will get reimbursed for this device, we will allow crossovers at 6 months, we will  
11 enroll patients outside the United States, we will ensure that there's diversity, which I do  
12 think is a very important issue and I take all those comments to heart. We will ensure that  
13 the blind is maintained. Although in that regard the one thing I'll say is you cannot see this  
14 device on a chest X-ray or an echo, you can't see it, okay, it's too small. You have to have a  
15 CT scan to be able to see it. So I do believe that the blinding, although we don't have a  
16 questionnaire to prove it, I do believe that the blinding in COSIRA was reasonable.

17 So bottom line, we don't have a perfect dataset by any means and as I said in my  
18 early comments, do I wish we had had a 500-patient sham-controlled trial, do I wish we  
19 were a hundred percent confident in mechanism and so in ischemia reduction? Yes, I do.  
20 But I just keep coming back to the unmet clinical need, the fact that this was recognized by  
21 the FDA by giving it a breakthrough designation, everyone on the Panel knows about these  
22 patients and while they don't deny if we were to start tomorrow a randomized trial, it  
23 would probably take 4 to 5 years before we would have that answer, especially in the  
24 COVID era. Again, not an excuse, but just the reality. Thank you very much.

25 MS. BEBEAU: Dr. Henry.

Free State Reporting, Inc.  
1378 Cape Saint Claire Road  
Annapolis, MD 21409  
(410) 974-0947

1 DR. HENRY: Thank you. I appreciate the opportunity to take a few moments and  
2 really, I'll tell you, I agree a hundred percent with what Gregg just summarized, including  
3 that we wish we had more data and I do believe that doing a postmarket trial will be  
4 doable. But I want to add a couple key points.

5 First of all, this is an incredibly challenging patient population to do clinical trials and  
6 in particular, they are miserable. And currently in the United States, they have no options,  
7 and currently what we're providing to those patients are riskier than what we saw today  
8 with the safety. So I do believe that the safety profiles are excellent based on the data that  
9 you saw with more than 2,000 patients and I do believe, taking care of these patients every  
10 day, that the risk-benefit profile is good. And I do believe that we need to do more data  
11 with -- and like COSIRA-II in terms of the postmarketing trial.

12 The last point I will make is that I would love to have us say that it's different in the  
13 United States, but the reality of this trial is the baseline characteristics and the medications  
14 are basically identical to the clinical trials that have been done in the United States  
15 including MUST-EECP, AGENT 3 and 4, RENEW --

16 (Audio cut out.)

17 DR. HENRY: -- I will end and I -- did I lose me?

18 DR. LANGE: We lost you for just in your last sentence, Tim, you're back on.

19 DR. HENRY: Yeah, so I think I'll stop with that but -- because I think Gregg's summary  
20 was outstanding and I want to thank the Panel because it really is an outstanding Panel with  
21 extensive experience taking care of these patients and so this discussion has been  
22 outstanding today.

23 MS. BEBEAU: Thank you, Dr. Henry.

24 From Neovasc, I want to thank the Panel for their time and for the careful  
25 considerations. And I would add again that the company is definitely committed to doing a

1 postmarket REDUCER-II study. Thank you for your time.

2 DR. LANGE: Great. I want to express my appreciation to the Sponsor, Dr. Stone,  
3 Dr. Henry, Chris Mullin, as well. Your presentations were excellent, you all were very  
4 responsive to the questions, very informative, and very respectful, as well. So to our  
5 colleagues, Dr. Stone, Dr. Henry, to Neovasc, I want to express my appreciation. Thank you.

6 We're now ready to move on to vote on the Panel's recommendations to the FDA.  
7 Prior to doing that, I'd like to ask Jacqueline Alikhaani, as our Consumer Rep, if she has any  
8 last comments. She'll have the last comment before we vote. Unfortunately, she's not a  
9 voting member. So Jackie.

10 MS. ALIKHAANI: Yes. Great discussion today, it was a really very informative and  
11 educational process and I appreciate the opportunity to be a part of this Panel. I think that  
12 we need to include the patient voice in more and more of these kind of discussions so  
13 patients can be able to make -- patients want to be empowered to make informed decisions  
14 and they have to do that in collaboration with the doctor and that includes the doctor  
15 having the best evidence from research.

16 And you can't get the best evidence without the best designed trials that show the  
17 best outcomes to inform everybody about what's the best decision. And I think when we do  
18 that, I think everybody feels good when they know you've done your best. And I know the  
19 best is not always easy to get to, I know the trials are complicated. But thank you so much  
20 for allowing me to be a part of this and speak out.

21 DR. LANGE: Ms. Alikhaani, thank you for participating and staying with us to the  
22 end. I really appreciate your perspectives and your opinions, thank you for sharing.

23 MS. ALIKHAANI: You're welcome.

24 DR. LANGE: We're now ready to vote on the Panel's recommendations to the FDA  
25 for Neovasc Reducer System. I'm happy to report I got no pre-election votes and we're not



1 mailing it into a ballot. So we will be conducting via e-mail which is a little unusual, we've  
2 not done that before. So Aden will give you some instructions for that. The Panel is  
3 expected to respond to three questions relating to safety, effectiveness, and benefit versus  
4 risk. And Aden Asefa will now read two definitions to assist in the voting process.

5 MS. ASEFA: Thank you, Dr. Lange.

6 Before we go into the vote, I just want to read a COI that I did not mention into the  
7 record. For the duration of the Circulatory System Devices Panel meeting on October 27th,  
8 2020, Ms. Jacqueline Alikhaani has been appointed to serve as a temporary nonvoting  
9 Consumer Representative. For the record, Ms. Alikhaani serves as a consumer rep for the  
10 Cardiovascular and Renal Drugs Advisory Committee at the Center for Drug Evaluation and  
11 Research. Ms. Alikhaani is a special Government employee who has undergone the  
12 customary conflict of interest review and has reviewed the materials to be considered at  
13 the meeting. The appointment was authorized by Russell Fortney, Director of the Advisory  
14 Committee Oversight and Management Staff, on September 22nd, 2020.

15 DR. LANGE: All right, Jackie, you're official now.

16 MS. ASEFA: Okay, Jim.

17 The Medical Device Amendments to the Federal Food, Drug and Cosmetic Act, as  
18 amended by the Safe Medical Devices Act of 1990, allow the Food and Drug Administration  
19 to obtain a recommendation from an expert advisory panel on designated medical device  
20 premarket applications that are filed with the Agency. The PMA must stand on its own  
21 merits and your recommendation must be supported by safety and effectiveness data in the  
22 application or by applicable publicly available information. The definitions of safety and  
23 effectiveness are as follows:

24 Safety as defined in 21 C.F.R. Section 860.7(d)(1) - There is reasonable assurance that  
25 a device is safe when it can be determined, based upon valid scientific evidence, that the

1 probable benefits to health from use of the device for its intended uses and conditions of  
2 use, when accompanied by adequate directions and warnings against unsafe use, outweigh  
3 any probable risks.

4 Effectiveness as defined in 21 C.F.R. Section 860.7(e)(1) - There is a reasonable  
5 assurance that a device is effective when it can be determined, based upon valid scientific  
6 evidence, that in a significant portion of the target population, the use of the device for its  
7 intended uses and conditions of use, when accompanied by adequate directions for use and  
8 warnings against unsafe use, will provide clinically significant results.

9 Panel members, we will now begin the voting process. I will read the proposed  
10 indication for use and the three voting questions and send each of the voting members an  
11 e-mail to respond to. Once I read all the three questions, we will tally the votes and read  
12 them into the record.

13 The indications for use is as follows: The Reducer System is intended for patients  
14 suffering from refractory angina pectoris despite guideline directed medical therapy, who  
15 are unsuitable for revascularization by coronary artery bypass grafting (CABG) or by  
16 percutaneous coronary intervention (PCI).

17 Voting Question Number 1: Is there a reasonable assurance that the Neovasc  
18 Reducer System is safe for patients who meet the criteria specified in the proposed  
19 indication? Please vote now yes, no or abstain.

20 (Panel vote.)

21 MS. ASEFA: Voting Question Number 2 reads as follows: Is there a reasonable  
22 assurance that the Neovasc Reducer System is effective for use in patients who meet the  
23 criteria specified in the proposed indication? Please vote now yes, no or abstain.

24 (Panel vote.)

25 MS. ASEFA: The third and final voting question reads as follows: Do the benefits of

1 the Neovasc Reducer System outweigh the risks for use in patients who meet the criteria  
2 specified in the proposed indication? Please vote now yes, no or abstain.

3 (Panel vote.)

4 DR. ALLEN: We have gotten an e-mail by now from you?

5 (Crosstalk.)

6 MS. ASEFA: -- tally and verify the official votes.

7 DR. LANGE: I'm sorry, if you've --

8 DR. ALLEN: I've not gotten any e-mails.

9 DR. LANGE: Yeah, I haven't gotten an e-mail, either.

10 MS. ASEFA: Oh, Dr. Lange, you wouldn't have gotten an e-mail. So I sent the e-mail  
11 at 6:20 p.m.

12 DR. PAGE: I received an e-mail and voted.

13 DR. SOMBERG: Is this an interactive e-mail or you just write (1) yes, no; (2) yes, no?

14 MS. ASEFA: Yes.

15 DR. ALLEN: If I just check it and go to the next thing, it doesn't hold.

16 MS. ASEFA: Oh, it's just an e-mail, it's just a reply e-mail.

17 DR. SOMBERG: So you just reply --

18 MS. ASEFA: Yeah.

19 DR. SOMBERG: -- and give your answers?

20 MS. ASEFA: Yeah.

21 DR. SOMBERG: Okay, it's not like a Google --

22 MS. ASEFA: Yeah, I know.

23 DR. SOMBERG: -- form or something? No, okay.

24 DR. ALLEN: I haven't gotten -- I haven't gotten an e-mail.

25 (Crosstalk.)

1 DR. LANGE: Okay.

2 UNIDENTIFIED SPEAKER: I voted.

3 DR. LANGE: If you've not gotten an email, would you raise your hand?

4 (Show of hands.)

5 DR. LANGE: Dr. Allen did not. Dr. Brindis did not. Dr. Vetovec did not. Dr. Connor,  
6 Dr. Wittes, Dr. Bonde. I guess we knew how you guys were going to vote and we decided  
7 not to extend one.

8 DR. BRINDIS: I have your e-mail --

9 (Crosstalk.)

10 DR. BRINDIS: I can send it on to you, but I've not received --

11 DR. WITTES: Can we just tell you how we're going to vote?

12 DR. LANGE: No, we got to -- no, we have to do it anonymously first, Janet, and then  
13 we tally them. But thank you --

14 DR. WITTES: All right, well.

15 DR. LANGE: Thank you for asking.

16 DR. SOMBERG: Well, it's another example of social media blocking the election.

17 (Laughter.)

18 DR. LANGE: That's fake news, John.

19 DR. SOMBERG: Richard, that's my line usually.

20 DR. LANGE: Okay, again let me see -- if you've not gotten it, let me see hands again.

21 (Show of hands.)

22 DR. LANGE: All right. And so Jim Swink, let me -- keep them up for a second until I  
23 call and then put your hand out. So Keith Allen has not, Ralph Brindis has not, Randall  
24 Starling has not, George Vetovec has not, Jason Connor has not, Janet Wittes has not, and  
25 Pramod Bonde has not.

1 DR. BRINDIS: You know, Richard, one option is we could vote privately, chat to her if  
2 that may not be working. I just got it, it looks like.

3 DR. LANGE: Yeah. We thought about that but sometimes it becomes less private  
4 than you'd like.

5 DR. STARLING: I have the e-mail, Rich.

6 DR. LANGE: Okay, great. I'm going to wait for 30 --

7 DR. CONNOR: No, I don't. But Aden, I have your e-mail. Aden, can I just send you an  
8 e-mail that's not in response to yours?

9 MS. ASEFA: I mean, that's fine. It's just basically saying yes, no or abstain to the  
10 three voting questions.

11 DR. LANGE: One, two, three, okay. Has anybody else not gotten the e-mail at this  
12 point?

13 (Audio feedback.)

14 DR. LANGE: I still got Dr. Vetovec, Dr. Bonde, and Dr. Allen still.

15 DR. ALLEN: Yeah, I just got -- I just got the e-mail. I'm just going to write one, two,  
16 and three and put my vote.

17 DR. LANGE: That would be great.

18 DR. ALLEN: It's just a blank e-mail to respond to.

19 DR. LANGE: That sounds great.

20 Dr. Vetovec.

21 DR. VETROVEC: I got mine.

22 DR. LANGE: Super, super.

23 Dr. Bonde, do you have yours yet? You're muted. So Dr. Bonde, do you have Aden's  
24 e-mail address?

25 DR. BONDE: Yeah.

1 DR. LANGE: If you'll just put one, two, and three and vote for each of those, safety,  
2 effectiveness, and benefit-risk. And then we'll make sure that Aden has all the votes.

3 I see, Bernie, you're trying to vote three times. Don't do that, Bernie. Jason asked  
4 what e-mail are you using, to Jim Swink.

5 DR. CONNOR: No, you can ignore that, I thought that was just to James.

6 DR. LANGE: Okay.

7 DR. CONNOR: I sent Aden my e-mail.

8 DR. LANGE: Got it, super.

9 (Tally of votes.)

10 DR. LANGE: All right.

11 DR. VETROVEC: How many times can we vote?

12 (Laughter.)

13 DR. BATCHELOR: As many times as you need to get it right, George.

14 (Laughter.)

15 DR. VETROVEC: I lived in Chicago as a boy, so I'm familiar with this.

16 DR. SOMBERG: Then there's no limit if you're from Chicago.

17 (Crosstalk.)

18 DR. ZUCKERMAN: Could we pause a moment?

19 Dr. Lange, do you know if Dr. Wittes is having a technical problem? I thought I heard  
20 her trying to speak again.

21 DR. LANGE: Janet, can you hear me?

22 (Audio feedback.)

23 DR. LANGE: So Aden, does Janet -- has Janet voted? Has Dr. Wittes voted?

24 MS. ASEFA: Yes, I can see her vote.

25 DR. LANGE: Okay. And Dr. Wittes, your audio -- your audio sounds like it's

1 somewhere between Mars and Jupiter. If it gets to that point where we can't hear you, I  
2 may have you type your response in terms of any no responses so we can enter it for the  
3 record.

4 But Aden, have you got everybody's vote?

5 MS. ASEFA: I think so. I'm still tallying them up.

6 DR. LANGE: While she's doing that, let me thank the Panel and Jim and Dino sitting  
7 in audiovisual at the end who has spent hours, countless hours, before the meeting making  
8 sure it works and with few exceptions, Dr. Wittes being one of them, it's worked pretty  
9 well. And I want to thank you, all the Panel members, this kind of a meeting goes a little bit  
10 longer because we're not sitting in the same room. But on the down -- or the upside is you  
11 don't have to fly back home or drive back home.

12 MS. ASEFA: Janet, did you vote? I did get an e-mail, but you didn't vote.

13 DR. WITTES: I'm trying to vote, I'm trying to.

14 MS. ASEFA: Okay.

15 DR. SOMBERG: Do we get paid by the hour for this or --

16 (Laughter.)

17 DR. LANGE: You do and you're on Mountain Time right now.

18 DR. SOMBERG: Illinois we have very strict rules about payment. Fifteen dollars an  
19 hour, even.

20 (Pause.)

21 DR. SOMBERG: Bram, you look so worried. The polls are in already.

22 (Pause.)

23 MS. ASEFA: Okay. Dr. Page, did I get your vote? Oh, yeah, I got it.

24 DR. PAGE: I sent one.

25 MS. ASEFA: Yes, I got it. Yeah. And then Dr. Wittes, Janet Wittes. Just waiting on

1 her --

2 DR. WITTES: No, but have you gotten -- I sent several -- several times I sent. You  
3 didn't get any of them?

4 MS. ASEFA: I didn't get the results, but I did get an e-mail from you. Hold on one  
5 second.

6 DR. LANGE: Okay, let me ask a question. Has everybody else voted?

7 MS. ASEFA: Yeah, everyone has voted.

8 DR. LANGE: Okay. So Janet, if you'd like to give your vote verbally right now, we'll  
9 take it. We'd ultimately reveal that anyway. Do you want --

10 DR. WITTES: No, no, no.

11 DR. LANGE: No? Okay.

12 DR. WITTES: Yeah, I voted no, no, no.

13 DR. LANGE: No, no, no. Okay. So let's enter those.

14 DR. WITTES: Sadly.

15 DR. LANGE: I thought that was your response to me when I said do you want to give  
16 a verbal thing and you said no, no, no, I thought --

17 DR. SOMBERG: And she said yes, yes, yes to that.

18 DR. WITTES: No, but with sadness I voted no, no, no.

19 DR. LANGE: Okay, we'll get the final tally here in just a moment.

20 (Pause.)

21 DR. LANGE: He's working on the hanging chads right now, we should have those  
22 resolved pretty quickly.

23 MS. ASEFA: Sorry about this.

24 DR. SOMBERG: Just imagine if the U.S. post office was involved. Three months.

25 (Pause.)



1 MS. ASEFA: Sorry, I have people double-checking my work.

2 DR. LANGE: Thank you, Aden. It's important to get it right, so thank you.

3 DR. VETROVEC: Does this go to the Supreme Court if there are any questions?

4 (Laughter.)

5 MS. ASEFA: Sure.

6 DR. ZUCKERMAN: Yes, it can, George. That's why it's really important just to give  
7 Aden a few minutes. I know that everyone wants to eat dinner and so forth, but it's really  
8 important.

9 DR. VETROVEC: No, I agree. I just thought it was funny.

10 DR. ZUCKERMAN: And as a preview, after Aden summates the votes, Dr. Lange will  
11 be asking Panel members to specify why they voted in a certain way. So right now, during  
12 this downtime, if you could, think about your responses so that that section can move along  
13 efficiently.

14 (Pause.)

15 MS. ASEFA: Okay. So the votes have been captured and I will now read the votes  
16 into the record.

17 On Question 1, the Panel voted 14 yes, 4 no, and 0 abstain the data shows  
18 reasonable assurance that the Neovasc Reducer System is safe for use in patients who meet  
19 the criteria specified in the proposed indication.

20 On Question 2, the Panel voted 1 yes, 17 no, and 0 abstain that there is a reasonable  
21 assurance that the Neovasc Reducer System is effective for use in patients who meet the  
22 criteria specified in the proposed indication.

23 And on Question 3, the Panel voted 3 yes, 13 no, and 2 abstain that the benefits of  
24 the Neovasc Reducer System outweigh the risks for use in patients who meet the criteria  
25 specified in the proposed indication.

1 The three voting questions are now complete.

2 DR. LANGE: Aden, thank you for taking the time to get it right. I appreciate it,  
3 thanks.

4 MS. ASEFA: Thank you.

5 DR. LANGE: I will now ask the Panel members to discuss their votes and please feel  
6 free to opine about why you voted yes or no on Questions 1, 2, and 3. And I'll run down the  
7 list.

8 Dr. Vetovec, you first.

9 DR. VETROVEC: I think the whole question really turns on how refractory patients or  
10 a patient is. I voted that there was safety, yes. I voted no to proving that it is effective.  
11 And I voted yes for risk-benefit, but I qualified that, which I can't do, I realize, but my verbal  
12 qualification is it depends on how symptomatic the patient is and where they are because  
13 it's not a terribly risky device. And so if the patient is symptomatic enough, I think I'd be  
14 willing to take it or anybody would be willing to take it, so I think that's the dilemma.

15 DR. LANGE: Thank you, Dr. Vetovec.

16 Dr. Starling.

17 DR. STARLING: Yes, thank you. Randall Starling.

18 So I'm very sympathetic to this patient group and frequently see these patients at  
19 the request of my interventional colleagues. As far as safety, I voted yes, based on the  
20 totality of data I was comfortable. The second question as far as effectiveness, I voted no,  
21 being uncomfortable with the data presented from the randomized trial and the sample of  
22 just over a hundred patients and the lack of the Reducer information, pushing it over for the  
23 reasons that we discussed during the meeting. As far as risk-benefit, I again voted no and  
24 my no vote was based -- primarily based on lack of any convincing efficacy.

25 DR. LANGE: Dr. Connor.

1 DR. CONNOR: Yes, Jason Connor here.

2 So I voted yes to safety, no to efficacy, and abstain for the risk-benefit. I think  
3 safety, given what we've seen in this trial, what we've seen in the open label, and then just  
4 the experience in Europe has all been positive with no big safety concerns.

5 The efficacy, for all the reasons discussed, I won't reiterate, is weak. You know, the  
6 company got FDA feedback, chose not to accept it and ran the trial elsewhere, they have an  
7 open IDE with FDA feedback for a trial they chose not to initiate and I think, you know, they  
8 had plenty of chances to run and obtain the efficacy that we all would like to have seen.

9 When it comes down to risk-benefit, I really appreciate hearing from the patients we  
10 heard from after lunch, I understand how debilitating this can be. And so my abstain vote  
11 there was because not as a clinician seeing these patients, I don't exactly feel qualified to  
12 weigh the benefit and the risk given it seems the risk is very low.

13 DR. LANGE: Thank you, Dr. Connor.

14 Dr. Brindis.

15 DR. BRINDIS: I voted yes for safety, no for efficacy, and no for risk-benefit. Again, I  
16 want to thank you, Richard, for running this historic meeting and but most of all, I do feel a  
17 lot of concern for the -- empathy for the patients involved. I'm really sorry that the  
18 company opted not to take on this trial when the FDA offered it a couple years ago, we'd be  
19 further along in terms of potentially helping our patients with a device that potentially  
20 could be of value.

21 DR. LANGE: Thank you, Dr. Brindis.

22 Dr. Allen.

23 DR. ALLEN: Yes, I voted yes for safety. Although I do question the long-term safety, I  
24 think overall the totality pushed me to vote yes. I voted no on efficacy. There are just too  
25 many issues with the trial. And I also have great conflicts. The FDA is, I think, doing a great

1 job at working with companies proactively to get trials designed correctly and the company  
2 had a chance to do this the right way and chose not to and the data just isn't good enough  
3 for me to vote yes. And risk over benefit, we have a zero in the denominator and you can't  
4 have a doable number, so I voted no on risk-benefit.

5 DR. LANGE: Dr. Wittes, don't surprise me you changed your vote now.

6 (Audio feedback.)

7 DR. LANGE: And unfortunately we can't get the sound issues worked out. Janet, you  
8 might want to turn your video off and see if that helps.

9 DR. WITTES: No, I'll chat.

10 DR. LANGE: That's good.

11 DR. WITTES: Good that I'm chatting or good that you can hear me?

12 DR. LANGE: We can hear you well now, go ahead.

13 DR. WITTES: Okay.

14 (Audio feedback.)

15 DR. LANGE: And now we can't, Janet, I'm so sorry. My apologies. If you chat, I'll  
16 make sure I read it to get it into the record.

17 I'll move on to Dr. Page at this point.

18 DR. PAGE: Thank you. I voted yes in terms of safety from the perspective that if it  
19 were an effective device, the safety, while it's not perfect, it would be acceptable. In terms  
20 of reasonable assurance of effectiveness for the reasons that I and others have  
21 enumerated, it just does not clear the bar. And then the issue of overall, I had to vote no,  
22 as well. I must say, as others have said, I was moved by hearing the patients and their  
23 families. The concept of no option, the concept of feeling desperate, the people talked  
24 about the urgent need and the patients need hope, but we shouldn't give them false hope.  
25 We need to be providing -- even with this vulnerable population, we need to provide them

1 something that we truly believe is going to be effective, and that was not proven today.

2 DR. LANGE: Well said. Thank you, Dr. Page.

3 Dr. Bonde.

4 DR. BONDE: I feel very sympathetic and I truly feel empathy for this patient  
5 population. I voted yes for the safety, yes, it can be safely put into these patients.  
6 However, I voted no for the other two questions because I don't think there's enough data  
7 to support either effectiveness or to prove a risk versus benefit ratio for these patients. So  
8 though we would like to help these patients, we would not like to have a device that is not  
9 having a proven efficacy and benefit.

10 DR. LANGE: Thank you, Dr. Bonde.

11 For the record, I'll read Dr. Wittes' comments. She voted no for safety because there  
12 is not enough long-term safety and I think 52 patients is too small a sample size to assess  
13 safety even in the short term. No for efficacy for all the reasons we all discussed. It made  
14 me very sad to vote no because I wanted to be convinced. I have no for risk-benefit  
15 because of the no for efficacy.

16 Thank you, Dr. Wittes.

17 Dr. Yuh.

18 DR. YUH: Thank you, Dr. Lange.

19 I voted yes for safety. I think, as I mentioned before, the totality of the safety profile  
20 data gave me confidence in the safety of the device. I voted no on the efficacy. You know, I  
21 don't want to enumerate all the flaws, but none of them was a deal killer to me or a fatal  
22 flaw but in totality, I think -- I just didn't have the confidence to meet the standard. The  
23 standard is the standard for efficacy and I just wasn't confident, and therefore I voted  
24 abstain on the ratio only because I think there may very well be a positive signal there but  
25 just couldn't be convinced with this particular trial and that's why I had to abstain.

1 DR. LANGE: Thank you, Dr. Yuh.

2 Dr. Yeh.

3 DR. YEH: I voted a yes for safety for all the reasons that have been enumerated. I  
4 voted no for efficacy because I didn't feel confident enough in that assessment and maybe I  
5 should -- there was an unknown for risk-benefit or uncertain I would've voted for that one,  
6 but I just -- but I voted no for risk-benefit just on the basis of although I think there may be  
7 a benefit, very similar to what Dr. Yuh just said, I'm just not confident enough that that's  
8 truly there.

9 DR. LANGE: Thank you, Dr. Yeh.

10 Dr. Batchelor.

11 DR. BATCHELOR: Yeah, I was really, really torn and I thought the statements from  
12 the patients were quite compelling and some comments at the end I thought were actually  
13 quite helpful. I was looking at this through the lens of it being a breakthrough designation,  
14 as well, and as an interventionalist who treats a lot of these patients who are really  
15 desperate for therapy and who right now are getting things like EECF and chronic total  
16 occlusion interventions which are probably less efficacious than for EECF and more  
17 dangerous. So I think that this is a safe procedure, I was convinced by the over 2,000  
18 patients treated, but I think we would see something come out if it was an unsafe  
19 procedure.

20 In terms of the effectiveness, I switched this last minute and I'm going to be an  
21 outlier because I thought about this. I initially voted no and I switched because I thought all  
22 the things that we discussed regarding the lack of -- the concerns about the placebo effect.  
23 If anything, that would bias toward the null and create less ability for us to actually  
24 distinguish the treatment effect between the two arms and that would, in fact, slightly  
25 disadvantage the device. I think the device is probably -- probably does work for angina, I

1 don't know that it works for ischemia, but as a patient-centered outcome for simply angina,  
2 I think that there's a good chance that this device does work, therefore I actually switched  
3 my vote last minute and said yes, and then because the risk-benefit profile is based on  
4 those two answers, I had to actually vote yes for the risk-benefit relationship. And this has  
5 all been in the context of absolutely requiring more data and that's where I think that's --  
6 that's where I feel that -- you know, this would be completely switched around if this was  
7 an end decision, but linked to the requirement for a randomized trial with the structure that  
8 we described, I'm okay making these -- I'm standing firmly by these answers.

9 DR. LANGE: Thank you, Dr. Batchelor.

10 Dr. Ohman.

11 DR. OHMAN: Yes, thank you, Dr. Lange.

12 I'm going to thank the Sponsor, Neovasc, for actually bringing forward a great new  
13 technology, and I also want to thank the FDA for doing very careful analysis of the very  
14 limited database that was presented to them.

15 I voted no on safety mainly because I really don't know the safety long term beyond  
16 the 6 months. It was fairly good data up to about a year, but after that it drops off quite  
17 dramatically and a lot of the long-term issues such as the impact that this might have on  
18 diastolic heart failure or other things were not available in the dataset, which I think is very  
19 important because many of these patients do get that.

20 I voted no for efficacy for the reasons that I couldn't link ischemia, which is  
21 objective, with a device that I would put in permanently when over about nearly half the  
22 patients had very little treatment benefit, so that makes it very hard. And then the final  
23 aspect is, of course, if I voted no on both I couldn't actually say that there was any risk-  
24 benefit equation that would be favorable. Thank you.

25 DR. LANGE: Thank you, Dr. Ohman.

1 Dr. Cigarroa.

2 DR. CIGARROA: So I, too, would like to thank the FDA and the Sponsor for their  
3 thoughtful approach to presenting the data and analyzing it and framing critical discussions  
4 for us today. These patients are very challenging, having cared for them on a frequent basis  
5 for over 25 years now, and because of this unmet need that adversely impacts quality of  
6 life, I wanted to vote yes on each of them. However, the conversations from today and the  
7 data presented, to me, provided clarity even though I didn't like the answers I was provided  
8 as it relates to the vote because of the unmet need for patients.

9 As it relates to safety, although I have concerns about the long-term components, I  
10 think the overall data presented made me vote yes despite the limitations of that data. As  
11 it relates to efficacy, I would say that the primarily male Caucasian European and non-U.S.  
12 approach made it very difficult for me to extrapolate our heterogeneous population, and  
13 then the challenges of a small sample size with the flaws that we discussed further  
14 cemented my response for voting no, which led to a vote of no on the third.

15 DR. LANGE: Thank you, Dr. Cigarroa.

16 Dr. Hirshfeld.

17 DR. HIRSHFELD: I'd like to also really compliment FDA on an outstanding and  
18 rigorous presentation and keeping us at a very high level through this whole meeting.

19 I voted no, no, and no, and for many reasons that were already outlined. I thought  
20 that the 6-month safety data was too short for a permanently implanted device. For as far  
21 as efficacy was concerned, I thought that the subjective signal of efficacy was there, but  
22 weak, and the objective signal of anti-ischemia is lacking and I thought we needed some  
23 more strength in each of those in order to vote for efficacy and that led me to, obviously,  
24 no for risk-benefit.

25 I had two other thoughts that I'd like to share. The first is I believe we heard data



1 today that suggested that in the countries in which this device is approved there's been  
2 somewhere around 2,000 devices implanted and that's over 9 years since the CE mark was  
3 granted. So I wonder if the clinicians in the countries where it's approved for the CE mark  
4 are telling us something if this patient population is as large as it is and there's been only  
5 that many implanted.

6 The second thing is that I think that if FDA were to approve this device now on the  
7 premise that a randomized pivotal trial is going to be done subsequently, there's a certain  
8 hypocrisy in that course of action, and so I think that it would be inconsistent to say we're  
9 going to approve it but you haven't -- to go on and demonstrate efficacy, so those are all  
10 reasons why I was a no, no, no person.

11 DR. LANGE: Thank you, Dr. Hirshfeld.

12 Dr. Mathew.

13 DR. MATHEW: Yeah, Verg Mathew here, Verghese Mathew. Thank you, Dr. Lange.  
14 Thanks for the opportunity to be a part of this very comprehensive discussion, I respect the  
15 opinion and contribution of every single person on the call here, so thank you for that.

16 I voted yes for safety for the reasons that have been enumerated. I voted no for  
17 efficacy. I was willing to buy the primary endpoint was met and I think that's important to  
18 say. I would've been more compelled if a secondary endpoint that was objective would've  
19 been appropriately powered and met, and I was concerned about the missingness of data in  
20 that regard, as well.

21 In terms of benefit-risk, I voted yes for similar reasons that Dr. Batchelor outlined, I  
22 think there's probably a signal for effect here. I think the relative risk is low. I think more  
23 data is important and I think how one formulates that, whether it's a randomized trial  
24 preapproval, I think for reasons I mentioned before, I think a randomized trial post-approval  
25 is challenging but whether there would be a pathway forward for a meticulous agreement

1 on a mandatory registry-type thing post-approval would be another consideration that I  
2 mentioned earlier with David. Those are my thoughts and thanks again for the opportunity  
3 to participate.

4 DR. LANGE: Thank you for participating, Dr. Mathew, appreciate it.

5 Dr. Gersh.

6 DR. GERSH: Thank you, Dr. Lange.

7 Well, this is my first meeting, it's really been a great experience for me. I'm very  
8 pleased to be part of the Panel and thank the FDA for inviting me. I really, really enjoyed  
9 the day and I have been very impressed with what I've heard and the level of discussion.

10 I voted yes for safety, but I accept that longer-term data would be beneficial. In  
11 terms of efficacy, I voted no for all the reasons you've heard and I think it's just too few  
12 patients and the signal, there's a signal there but given the other limitations to the data, I  
13 could not -- I didn't feel I could vote for efficacy.

14 My third was no really based upon Number 2. Until we have some proof of efficacy,  
15 I'm not sure I feel comfortable assessing risk-benefit. If there's a randomized trial that  
16 shows even a mild to moderate but significant signal, then I think I would support the risk-  
17 benefit ratio. Thank you.

18 DR. LANGE: Thank you, Dr. Gersh. And welcome to the Committee, we're glad to  
19 have you for the first time.

20 Dr. Somberg, this is not your first ride.

21 DR. SOMBERG: You should've said rodeo.

22 DR. LANGE: Go, John.

23 DR. SOMBERG: Okay. Well, I voted no for safety because I felt the long-term data  
24 was really inadequate for a permanently implanted device. I do not think the risk of  
25 anticoagulation has been evaluated. Because it's on the venous side I think there's going to

1 be a need for an antithrombotic, warfarin, for instance, and that might increase the bleeds.  
2 And I also think there's -- I was disturbed by the discussion of stent thrombosis and whether  
3 there's a thrombosis or not, and I think that needs to be worked out pre-clinically so we  
4 know if there is an incidence of thrombosis because if you do have a coronary sinus  
5 thrombosis it may present with angina or MI and it's going to be a question if it's disease  
6 progression or not.

7         Okay, I also voted no for efficacy because I think the population was not defined and  
8 the data was very interesting in that those people with the mildest or the least drugs, which  
9 we will say the least resistant patients, benefitted the most. So I'm not sure this drug will  
10 be effective in a drug-resistant end-stage population. I think it may be an antianginal but  
11 then again it's going to be a risk-benefit ratio whether you want to take a medicine that can  
12 be easily stopped versus you want to have a permanent implant with the risk that that  
13 ensues. So not knowing the safety, not knowing the efficacy, obviously I voted no for the  
14 third question.

15         DR. LANGE: Thank you, Dr. Somberg.

16         DR. SOMBERG: Thank you.

17         DR. LANGE: And Dr. Borer. And again, Jeff, if you'll unmute. Thank you.

18         (No response.)

19         DR. LANGE: Can't hear you yet. Find that unmute button one more time. Should be  
20 in the top right, I think.

21         DR. CONNOR: And mine's lower left if you're on a Mac.

22         DR. SOMBERG: Yeah, lower left.

23         DR. LANGE: I'll just ask you to unmute. Yeah, so Jeff, are you able to see that mute  
24 button down on the lower left?

25         DR. CONNOR: Is your phone on mute?

1 (No response.)

2 DR. LANGE: I'm going to ask you to unmute. And I'm sorry, neither I nor the  
3 audiovisual host is able to unmute anybody. Wait for just a minute.

4 DR. SOMBERG: Maybe he can chat in.

5 DR. LANGE: Yeah. So Jeff, would you like to just type in on the chat and I can enter  
6 it into the record?

7 (No response.)

8 DR. LANGE: You don't see the chat, either. Okay.

9 MS. ALIKHAANI: Maybe he can just send an e-mail.

10 DR. LANGE: Yeah. Jeff, you're there. Yeah, you're there. Don't touch anything.

11 DR. BORER: Okay.

12 DR. LANGE: Just talk.

13 DR. BORER: With regard to safety, I voted yes.

14 DR. LANGE: Okay.

15 DR. BORER: With trepidation because I agree totally with John and others that the  
16 long-term data aren't there, that the issues about the need for antithrombotic therapy  
17 hasn't been worked out. I mean, there are a lot of problems with it, but from the data that  
18 exist, it seems to me unlikely that this device is going to cause an important safety problem,  
19 so I voted yes.

20 With regard to efficacy, I had to vote no, first because the data aren't there and  
21 second and perhaps more importantly, because the population was poorly defined. I have  
22 no idea who we would want to give this device to. It's people with refractory angina, but  
23 that was never defined, so I couldn't possibly vote yes for that. And if I couldn't vote yes for  
24 effectiveness, then how could I vote yes for the relationship between effectiveness and  
25 safety, so I voted no for that, also. So that's why I voted the way I did.

1 DR. LANGE: Great. And I don't vote unless there's a tie, but I always appreciated the  
2 fact and respected Dr. Page who always told us how he voted regardless, and I would've  
3 voted no on all three for the reasons that you all discussed already.

4 I would like to thank the Panel, it really is an august panel, very accomplished  
5 individuals on here. Thank you for those who I have served with before and those of you  
6 who are new serving on the Panel. I want to thank the FDA, excellent presentations,  
7 excellent analysis, excellent discussion, and excellent preparation. The same to the  
8 Sponsor, I think the presentations were outstanding and again, I appreciate you addressing  
9 the issues. I look forward to you working with the FDA, and I hope that instructions and  
10 observations that we provided will help you to work with the FDA to get this to a point  
11 where we can reevaluate it and assess its safety and efficacy and hopefully, approve it for  
12 this very needy population. And I want to thank our public hearing speakers for their  
13 contributions and Jackie Alikhaani, our Consumer Rep, as well.

14 Dr. Zuckerman, do you have any final comments?

15 DR. ZUCKERMAN: Yes. I especially want to thank Dr. Lange for running an extremely  
16 successful virtual advisory panel meeting. This is the first devices advisory panel meeting  
17 that has been run in the COVID era virtually and Dr. Lange has certainly set the standard. So  
18 thank you. Also, I'd like to thank all the Panel members for the extremely excellent work  
19 that was done this afternoon. The FDA and more importantly, the American people owe  
20 you a debt of gratitude. So thank you and Godspeed.

21 DR. LANGE: With those comments, this meeting of the Circulatory Devices Panel is  
22 now adjourned. Again, thank you to everybody. Have a good night.

23 (Whereupon, at 8:01 p.m., the meeting was adjourned.)

24

25

C E R T I F I C A T E


This is to certify that the attached proceedings in the matter of:

CIRCULATORY DEVICES PANEL

October 27, 2020

Via Zoom

were held as herein appears, and that this is the original transcription thereof for the files of the Food and Drug Administration, Center for Devices and Radiological Health, Medical Devices Advisory Committee.

A handwritten signature in black ink that reads "Tom Bowman". The signature is written in a cursive style with a horizontal line underneath it.

TOM BOWMAN

Official Reporter