



Premarket Approval Application (PMA) for the Lutonix 014 Drug Coated Balloon (Lutonix 014 DCB)

Cardiovascular Devices Advisory Panel Meeting
February 17, 2021

*Presentations By: Eleni Whatley, Ph.D.,
Rong (Rona) Tang, Ph.D., and Donna Buckley, M.D., M.S.
Office of Cardiovascular Devices (OCVD)/Office of Clinical Evidence and Analysis (OCEA)
Office of Product Evaluation and Quality (OPEQ)
Center for Devices and Radiological Health (CDRH)*



Review Team

Eleni Whatley, PhD	- Lead and Engineering Reviewer
Donna Buckley, MD, MS	- Clinical Reviewer
Rong (Rona) Tang, PhD	- Statistical Reviewer
Karen Manhart, VMD	- Animal Studies Reviewer
Hajira Ahmad, PhD	- Biocompatibility Reviewer
Michelle Johnson, PhD	- Chemistry Reviewer
Ron Brown, PhD	- Toxicology Reviewer
Parnali Chatterjee, PhD	- CDER PK reviewer
Monica Cooper, PhD	- CDER CMC reviewer
Steve Hertz, PE	- CDER cGMP reviewer
Kimberly Kopecki, MS, RAC	- Lead and GMP reviewer
Megha Reddy	- Sterility Reviewer

Presentation Outline



- Introduction
- Clinical Background
- Device Description and Indications for Use
- Regulatory History
- IDE Clinical Trial Design and Statistical Analysis Plan
- Clinical Data Review
- Additional Data Evaluations
- Post-Approval Study Proposal
- Benefit/Risk Summary and Conclusions



Introduction, Clinical Background, Device Description, and Regulatory History

Eleni Whatley, Ph.D.
Biomedical Engineer
CDRH/OPEQ/OCVD/PMDT




Meeting Purpose

- To request Panel feedback on whether the Lutonix 014 DCB provides a reasonable assurance of safety and effectiveness and a favorable benefit/risk profile for the treatment of below-the-knee (BTK) critical limb ischemia (CLI).
- The panel will be asked if:
 - For *safety*, whether the Lutonix 014 DCB is associated with “**the absence of unreasonable risk** of illness or injury associated with the use of the device for its intended uses and conditions of use.” (21 CFR 860.7(d)(1))
 - For *effectiveness*, if it has been determined that “based upon valid scientific evidence, that in a significant portion of the target population, the use of the device for its intended uses... will provide **clinically significant results**.” (21 CFR 860.7(e)(1))

Clinical Background



- Peripheral arterial disease (PAD):
 - Narrowing of the blood vessels in the lower extremities
 - Results in inadequate blood flow to downstream tissues
 - Clinical signs and symptoms: Leg pain when walking (claudication)
 - Categorized: Rutherford Classification (RC) 1 to 3
- Critical Limb ischemia (CLI):
 - Clinical signs and symptoms: Limb pain at rest, tissue loss, non-healing wounds, gangrene, and potential amputation
 - Categorized: RC 4 to 6



Clinical Background






- Available treatment strategies for BTK CLI include:

- Surgical

- Bypass

- Intravascular

- PTA
- Atherectomy
- Off-label device use




- If approved, the Lutonix 014 DCB would be the first-of-a-kind device for the treatment of BTK CLI in the US beyond vessel preparation strategies (PTA and atherectomy).

Device Description



- The Lutonix 014 DCB is an over-the-wire PTA catheter with a paclitaxel drug coating on the surface of the balloon.
- The primary mode of action for the Lutonix 014 DCB is achieved through the mechanical dilatation of the vessel during the balloon inflation. Drug delivery during the dilatation is designed to provide the ancillary benefit of preventing restenosis.





Device Description

- The Lutonix drug coating contains the anti-proliferative drug Paclitaxel as well as various excipients and solvents.
- The balloon is coated with a constant $2 \mu\text{g}/\text{mm}^2$ of paclitaxel; thus the total drug dosage is correlated to the balloon surface area.
- Available balloon diameters 2-4 mm, and available balloon lengths 40-150 mm.

Balloon Diameter	Total Dosage (mg) per Respective Balloon Length					
	40 mm	60 mm	80 mm	100 mm	120 mm	150 mm
2.0 mm	0.5	0.8	1.0	1.3	1.5	1.9
2.5 mm	0.6	0.9	1.3	1.6	1.9	2.4
3.0 mm	0.8	1.1	1.5	1.9	2.3	2.8
3.5 mm	0.9	1.3	1.8	2.2	2.6	3.3
4.0 mm	1.0	1.5	2.0	2.5	3.0	3.8

Proposed Indications for Use



The Lutonix 014 Drug Coated Balloon PTA catheter is indicated for patients with critical limb ischemia who have obstructive de novo or non-stented restenotic lesions in native popliteal, tibial, and peroneal arteries up to 320 mm in length and 2.0 to 4.0 mm in diameter.

Regulatory History – IDE Approval



Original Investigational Device Exemption (IDE) approved May 30, 2013 for:

- 480 subjects
- Randomized 2:1 to the Lutonix 014 DCB or PTA
- Primary safety endpoint: Freedom from below-the-knee (BTK) major adverse limb event (MALE) + perioperative death (POD) at 30 days
- Primary effectiveness endpoint: A composite of limb salvage and primary patency at 12 months

Regulatory History



- After IDE approval in 2013, the study enrolled for 5 years
- Multiple protocol changes were enacted, some to assist in enrollment
 - 34 IDE supplements for various changes, including to the design and manufacturing of the device, clinical protocol modifications, and other administrative changes
 - The two most significant changes were to the primary endpoint analysis and primary endpoint assessment timepoint (to be discussed)
- Trial terminated early and data were submitted to FDA
 - Predefined study success criteria not met
 - FDA concluded that a robust primary endpoint treatment effect was not present and we found no greater benefit of the Lutonix DCB vs. PTA for secondary effectiveness endpoints

Revision to FD&C Act, July 2012



- FDA will only disapprove an IDE due to concerns related to subject safety and protections
- FDA will not disapprove an IDE because we do not believe that the investigation will support approval of a device
- Feedback or concerns related to the study design can be provided as study design considerations

Regulatory History – IDE Changes



- Major protocol changes made during course of investigation include:
 - Removal of hemodynamic inclusion criteria
 - Rutherford Classification 3 patient enrollment permitted
 - Increased sample size and addition of interim analyses
 - Shortened primary effectiveness endpoint from 12 months to 6 months
 - Added co-primary effectiveness endpoint for “proximal segment flow pathway”
 - Added secondary effectiveness endpoint excluding early mechanical recoil
 - Early study termination

Primary Effectiveness Endpoint Assessment Change From 12 months to 6 Months




- Change occurred:
 - Approximately 3 years after original IDE approval
 - After enrollment of 325 patients
- Sponsor rationale: Clinical meaningfulness of a 6-month assessment due to aggressive nature of the disease.
- FDA feedback:
 - Acknowledged that an improvement at 6 months may be clinically meaningful
 - **Durability of the treatment effect was also valued by patients and physicians and would be important in demonstrating a reasonable assurance of effectiveness and a favorable benefit-risk profile for the Lutonix 014 DCB.**

Added Proximal Segment Co-Primary Effectiveness Endpoint



- Change occurred:
 - Approximately 4.5 years after the original IDE approval
 - After enrollment of 440 patients
- Background:
 - “Flow pathway”: analysis unit for the primary effectiveness endpoint
 - Includes the following arteries: popliteal, tibioperoneal, anterior tibial, posterior tibial, and peroneal
 - Multiple lesions may be present in a single flow pathway
 - “Proximal segment flow pathway”: Vessels that are entirely within the proximal 2/3 segment of the target flow pathway boundary or are split across the 2/3 segment cut-off



Proximal Segment Analysis - Background



- The co-primary effectiveness endpoint included two assessments:
 - Assessment 1: For “full flow pathway”
 - Assessment 2: If the “full flow pathway” assessment did not show superiority, the analysis was repeated for the “proximal segment flow pathway”
- Success = either of the analyses reach statistical significance
 - Due to this change, and addition of interim analyses (to be discussed), the alpha level to reach statistical significance was reduced to 0.0085 for both co-primary endpoints in order to control overall study type 1 error rate at 0.025
- Only the full flow pathways analysis is presented

Regulatory History – PMA Timeline



- PMA submission after early IDE termination
- Major Deficiency Letter – questions regarding effectiveness
- Agency Directed Assignment (ADA) - requested input on effectiveness
- Two Not Approvable (NOAP) letters issued
- Sponsor request for an Advisory Committee Meeting





IDE Study Design and Statistical Analysis Plan (SAP)

Rong (Rona) Tang, Ph.D.
Statistical Reviewer
CDRH/OCEA/DCEAII

Pivotal Clinical Trial Design



- Prospective, multicenter, 2:1 randomized, controlled trial
- Compared the Lutonix 014 DCB (test group) vs. standard balloon angioplasty (control group) for treatment of BTK arteries
- Sample Size: Approved for up to 1000 subjects to obtain 840 flow pathways evaluable for the primary effectiveness endpoint
- Sites: 75 Global Sites

Primary Safety Endpoint



Definition: Freedom from both BTK MALE (major adverse limb event) + POD (peri-operative death) at 30 days

- MALE:
 - Above ankle amputation; or
 - Major reintervention (new bypass graft, jump/interposition graft revision, or thrombectomy/thrombolysis) of the index limb involving a BTK artery
- POD: All-cause mortality

Primary Safety Endpoint Statistical Hypothesis and Analysis



- Hypothesis:

$$H_0: p_{\text{Control}} - p_{\text{DCB}} \geq \delta \text{ and } H_1: p_{\text{Control}} - p_{\text{DCB}} < \delta$$

where p is the success rate in each treatment group arm through 30 days and δ is the non-inferiority margin of 12%

- Pre-Specified Statistical Analysis:

- Primary analysis population: Intent-to-Treat
- Analysis unit: Per-Patient
- Farrington and Manning test with one-sided alpha = 0.025

Primary Effectiveness Endpoint



Definition: Limb Salvage + Primary Patency at 6 months

- Limb Salvage: Freedom from above ankle amputation in index limb
- Primary Patency: Freedom from target lesion occlusion and freedom from clinically driven target lesion revascularization (CD-TLR)*

*CD-TLR defined as CEC-adjudicated reintervention due to worsening Rutherford class, stagnant or worsening wound healing, or a new or recurrent wound

Co-Primary Effectiveness Endpoint Statistical Hypothesis



- Hypothesis :

$$H_0: \beta \leq 0 \text{ vs } H_1: \beta > 0$$

Where β is the increment in log odds of success in the treatment group

- Up to two hypotheses may be tested
 - Full flow pathway analysis
 - If not successful, proximal segment flow pathway population

Co-Primary Effectiveness Endpoint Statistical Analysis



Pre-Specified Statistical Analysis


- Primary analysis population: intent-to-treat
- Analysis unit: per-flow pathway
- Two hypotheses each tested at an alpha = 0.0085
 - First test the full flow pathway
 - If the full flow pathway test is not successful, retest the hypothesis for the proximal segment flow pathway
- Logistic regression with random effects

Interim Analyses




- Bayesian decision process was used to adjust the sample size to a maximum of 840 flow pathways.
- May occur when 400, 500, 600, and 700 randomized flow pathways are treated.

Adaptive by Design Process for Primary Effectiveness Endpoint




Adaptive by Design Process for Primary Effectiveness Endpoint



Interim Analysis and False Positive Rate Inflation




- Without interim analysis:
 - Comparing p-value to 0.025
- With interim analyses to allow for early stopping for effectiveness:
 - Increased chance of type I error (false positive rate)
 - Reduce significance level required to maintain an overall 0.025 false positive rate



Alpha Level Determination

Pre-specified alpha allocation for the primary effectiveness endpoint

Two final hypothesis tests with
significance level = 0.0085 each



Interim analysis with early stopping
for potential efficacy (up to 4 analyses)



Overall one-sided type I error rate
controlled at 0.025 (via simulation)

Pre-Specified Hypothesis Testing for Secondary Endpoints



- Superiority may be tested at 6 months for the following:
 - Primary patency with exclusion of early mechanical recoil
 - Primary patency
 - Freedom from clinically-driven TLR
 - Composite of freedom from above-ankle amputation, unhealed wound, ischemic rest pain, target vessel occlusion, and clinically-driven target vessel revascularization (TVR)
- Sequential testing performed only if the study showed success for the primary endpoints

Since the primary effectiveness endpoint was not met, secondary endpoint hypotheses were not tested



IDE Study Design and Clinical Data Review

Donna Buckley, M.D., M.S.
Interventional Radiologist/Medical Officer
CDRH/OPEQ/OCVD/PIDT



Key Selection Criteria

- Rutherford Classification
 - Rutherford 5 (minor tissue loss)
 - Rutherford 4 (ischemic rest pain)
 - Rutherford 3 (severe claudication)
 - Added to increase enrollment; accounted for ~10% of enrolled patients
- Stenosis >70% in one or two BTK arteries
- Lesion Dimensions
 - Total length of all target lesions <320 mm
 - 2-4 mm diameter
- Reconstitution of flow at or above the ankle
- Hemodynamic criteria were removed to increase enrollment



Design Elements

- Third Party Review
 - Data Monitoring Committee (DMC)
 - Clinical Events Committee (CEC)
 - Core Laboratories: Duplex Ultrasound (DUS) and angiography
- Blinding
 - Patient was blinded to 6 months
 - CEC, Core Labs, and Sponsor were blinded
 - Physicians performing the index procedure and follow-up clinical evaluations were *not blinded*



Select Secondary Endpoints

- Primary Safety Endpoint Evaluated Longer Term
- Primary Effectiveness Endpoint Evaluated Longer Term
- Wound Healing
- Amputation: Major (above ankle) and minor (below ankle)
- Freedom from clinically-driven target lesion revascularization (CD-TLR)
- Primary Patency (both including and excluding early mechanical recoil)
- Secondary Patency
- Hemodynamic outcomes (ABI, TBI)
- Change in Rutherford Classification
- Quality of Life (EQ-5D; Change in Walking Impairment Questionnaire (WIQ))
- All-Cause Mortality

Enrollment and Accountability



SAFETY

442 patients randomized (2:1)

Per Patient
(287 DCB, 155 PTA)

Evaluable for 30-Day Safety
(Per-Patient)

286 DCB Patients
(99.7%)

155 PTA Patients
(100%)

6-Month Follow-Up
(Per-Patient)

272 DCB Patients
(94.8%)

146 PTA Patients
(94.2%)

EFFECTIVENESS

507 flow pathways enrolled

Per Flow Pathway
(323 DCB, 184 PTA)

Evaluable for 30-Day Safety
(Per-Flow Pathway)

294 DCB Flow paths
(91.0%)

156 PTA Flow paths
(84.8%)

Evaluable for 6-Month Effectiveness
(Per-Flow Pathway)

269 DCB Flow paths
(83.3%)

137 PTA Flow paths
(74.5%)

Baseline Demographics and Lesion Characteristics



Baseline Demographics

	DCB Subjects (N=287)	PTA Subjects (N=155)
Mean Age, Years	72.9	72.9
Gender, n/N (%)		
Male	202/287 (70.4%)	104/155 (67.1%)
Female	85/287 (29.6%)	51/155 (32.9%)
Race, n/N (%)		
American Indian/Alaska Native	1/287 (0.3%)	0/155 (0.0%)
Asian	25/287 (8.7%)	15/155 (9.7%)
Black or African American	33/287 (11.5%)	12/155 (7.7%)
White	226/287 (78.7%)	127/155 (81.9%)
Other	2/287 (0.7%)	1/155 (0.6%)
Rutherford Category, n (%)		
3	26 (9.1%)	16 (10.3%)
4	100 (34.8%)	52 (33.5%)
5	161 (56.1%)	87 (56.1%)
History of Risk Factors, n/N (%)	285 / 287 (99.3%)	155 / 155 (100.0%)
Diabetes	204 / 287 (71.1%)	106 / 155 (68.4%)
Dyslipidemia	225 / 287 (78.4%)	116 / 155 (74.8%)
Hypertension	264 / 287 (92.0%)	148 / 155 (95.5%)
Cigarette Smoking	170 / 287 (59.2%)	89 / 155 (57.4%)
Current	43 / 287 (15.0%)	19 / 155 (12.3%)
Former	127 / 287 (44.3%)	70 / 155 (45.2%)

Lesion Characteristics

	DCB Lesions (N=352)	PTA Lesions (N=212)
Lesion Type, n (%)		
Distal 1/3	17 (4.8%)	14 (6.6%)
Proximal 2/3	194 (55.1%)	121 (57.1%)
Split across 2/3 reference line	126 (35.8%)	70 (33.0%)
Unknown/NA	15 (4.3%)	7 (3.3%)
Target Lesion Length (mm)		
Mean (SD)	111.8 (92.64)	94.7 (85.36)
Any Calcification, n/N (%)	211/352 (59.9%)	115/212 (54.2%)
Severe Calcification, n/N (%)	53/352 (15.1%)	28/212 (13.2%)
TASC Lesion Type, n (%)		
A	182 (51.9%)	131 (62.7%)
B	61 (17.4%)	32 (15.3%)
C	62 (17.7%)	28 (13.4%)
D	46 (13.1%)	18 (8.6%)

Similar

- Demographics & Medical Conditions
- Lesion Characteristics
- Medication Usage



Primary and Long-Term Safety Outcomes



Primary Safety Endpoint Results


Freedom from both BTK MALE (major adverse limb event) + POD (peri-operative death) at 30 days (ITT Per-Patient)

	DCB Subjects (N=287)	PTA Subjects (N=155)	Difference	P value
Free from Primary Safety Events through 30 Days	284 / 286 (99.3%) (97.5%, 99.9%)	154 / 155 (99.4%) (96.5%, 100.0%)	-0.1% (-3.9%, 3.8%)	< 0.0001
Primary Safety Events				Non-inferiority margin = 12%
Above Ankle Amputation ≤ Day 30	0	0		SAFETY ENDPOINT MET
Major Re-intervention ≤ Day 30	1	0		
Death ≤ Day 30	1	1		



Long-Term (36-Months) Safety Results

**Freedom from BTK MALE (Major Adverse Limb Event)
+ POD (Peri-Operative Death) through 36 Months**
(ITT Per-Patient; Binary Outcome)




Visit	DCB Subjects (N=287) Response Rate	PTA Subjects (N=155) Response Rate	Difference (95% CI)
30 Days	284 / 286 (99.3%)	154 / 155 (99.4%)	-0.1% (-1.6%, 1.5%)
6 Months	265 / 272 (97.4%)	139 / 146 (95.2%)	2.2% (-1.7%, 6.2%)
12 Months	242 / 253 (95.7%)	123 / 131 (93.9%)	1.8% (-3.1%, 6.6%)
24 Months	202 / 218 (92.7%)	100 / 110 (90.9%)	1.8% (-4.6%, 8.1%)
36 Months	146 / 162 (90.1%)	66 / 77 (85.7%)	4.4% (-4.7%, 13.5%)

Long-Term (36-Months) Safety Results



**Freedom from BTK MALE (Major Adverse Limb Event)
+ POD (Peri-Operative Death) through 36 Months
(ITT Per-Patient; KM Estimates)**





Primary and Long-Term Effectiveness Outcomes



Primary Effectiveness Endpoint Results

Limb Salvage + Primary Patency at 6 months (ITT Per-Flow Pathway)


		DCB n/N (%) (95% CI)	PTA n/N (%) (95% CI)	Difference in Response (95% CI)	P- value
Free from Primary Effectiveness Failure at 6 Months		201/269 (74.7%) (69.1%, 79.8%)	88/137 (64.2%) (55.6%, 72.2%)	10.5% (0.3%, 18.8%)	0.0222 P-value <u>not</u> <0.0085 for superiority
Composite Endpoint Failure Events ≤ Day 210, n					EFFECTIVENESS ENDPOINT <u>NOT</u> MET
Patients with major amputation		4 (1.4%)	3 (2.0%)		
Pathways with clinically driven TLR		28 (10.4%)	30 (21.9%)		
Pathways with total occlusion on imaging		42 (15.6%)	21 (15.3%)		

	DCB n/N (%)	PTA n/N (%)
All evaluable pathways with clinically driven TLR	28/303 (9.2%)	30/172 (17.1%)

Long-Term (36-Months) Effectiveness Results



Limb Salvage + Primary Patency through 36 Months (ITT Per-Flow Pathway; Binary Outcome)




Visit	DCB Pathways (N=323) Response Rate	PTA Pathways (N=184) Response Rate	Difference (95% CI)
30 Days	283 / 294 (96.3%)	144 / 156 (92.3%)	4.0% (-1.0%, 7.9%)
6 Months	201 / 269 (74.7%)	88 / 137 (64.2%)	10.5% (0.3%, 18.7%)
12 Months	128 / 251 (51.0%)	75 / 132 (56.8%)	-5.8% (-17.0%, 5.2%)
24 Months	84 / 228 (36.8%)	54 / 123 (43.9%)	-7.1% (-17.5%, 4.5%)
36 Months	58 / 210 (27.6%)	29 / 100 (29.0%)	-1.4% (-11.6%, 11.3%)

Long-Term (36-Months) Effectiveness Results




Limb Salvage + Primary Patency through 36 Months (ITT Per-Flow Pathway; KM Estimates)





Missing Effectiveness Outcome Data

Missing Data (ITT Effectiveness)



- Missing data higher for the PTA group at all time points
- Even at the 6 month time point, there is substantial amount of missing data, notably LTFU in the PTA group
- High rates of missing data present across all follow-up time points

High rates of missing data add uncertainty regarding study conclusions




Secondary Endpoints Wounds

Wound Healing



Percentage of Healed Wounds (ITT Per Patient)



Limitations


- No uniform wound assessment or healing scale used
- Wounds assessed by the unblinded physicians performing treatment
- Wound photographs did not undergo third-party independent review
- No photo required if a wound was deemed healed

Wound Healing




Percentage of Healed Wounds (ITT Per Patient)

Infection Present



Gangrene Present






Secondary Endpoints

Freedom From Amputation


Amputation



**Freedom From Major Amputation
(ITT Per-Patient KM Estimates)**



**Freedom From Minor Amputation
(ITT Per-Patient KM Estimates)**



Freedom from Major Amputation (ITT Per-Patient Binary Outcomes)

Visit	DCB Response Rate	PTA Response Rate	Difference 95% CI
30 Days	286/286 (100.0%)	154/155 (99.4%)	0.6% (-0.6%, 1.9%)
6 Months	267/271 (98.5%)	142/145 (97.9%)	0.6% (-2.1%, 3.3%)
12 Months	244/251 (97.2%)	127/130 (97.7%)	-0.5% (-3.8%, 2.8%)
24 Months	204/215 (94.9%)	103/109 (94.5%)	0.4% (-4.8%, 5.6%)
36 Months	148/159 (93.1%)	67/ 74 (90.5%)	2.5% (-5.2%, 10.3%)

Freedom from Minor Amputation (ITT Per-Patient Binary Outcomes)

Visit	DCB Response Rate	PTA Response Rate	Difference (95% CI)
30 Days	272 / 286 (95.1%)	143 / 155 (92.3%)	-2.8% (-7.7%, 2.0%)
6 Months	241 / 274 (88.0%)	128 / 147 (87.1%)	-0.9% (-7.5%, 5.8%)
12 Months	219 / 257 (85.2%)	110 / 135 (81.5%)	-3.7% (-11.6%, 4.1%)
24 Months	183 / 226 (81.0%)	88 / 117 (75.2%)	-5.8% (-15.1%, 3.6%)
36 Months	135 / 178 (75.8%)	53 / 83 (63.9%)	-11.9% (-24.1%, 0.1%)




Secondary Endpoints
Freedom From Clinically-Driven
Target Lesion Revascularization (CD-TLR)

CD-TLR



Freedom From CD-TLR

(ITT Per-Flow Pathway; KM Estimates)



Freedom From CD-TLR

(ITT Per Plow-Pathway; Binary Outcomes)

Visit	DCB (N=323) Response Rate	PTA (N=184) Response Rate	Difference (95% CI)
30 Days	317/321 (98.8%)	179/184 (97.3%)	1.5% (-2.0%, 4.1%)
6 Months	275/303 (90.8%)	142/172 (82.6%)	8.2% (1.5%, 13.3%)
12 Months	216/281 (76.9%)	116/152 (76.3%)	0.6% (-9.8%, 8.5%)
24 Months	169/249 (67.9%)	85/130 (65.4%)	2.5% (-9.2%, 12.4%)
36 Months	115/203 (56.7%)	52/99 (52.5%)	4.1% (-9.4%, 16.8%)




Secondary Endpoints

Secondary Patency

Secondary Patency



Definition: Freedom from total occlusion independent of whether or not patency was re-established via an endovascular procedure







Secondary Endpoints

Ankle Brachial Index (ABI) and Toe Brachial Index (TBI)

ABI and TBI



Change in Hemodynamic Outcomes (ITT Per-Patient)






Secondary Endpoints Rutherford Classification



Rutherford Classification

Rutherford Classification Through 36 Months (ITT Per-Patient)






Secondary Endpoints Quality of Life (QoL) Assessment




EQ-5D and WIQ

EQ-5D Change from Baseline



Walking Impairment Questionnaire (WIQ)






Secondary Endpoints

All-Cause Mortality



Mortality

All-Cause Mortality Through 36 Months (ITT Per-Patient)



Visit	DCB Pathways (N=287) Response Rate	PTA Pathways (N=155) Response Rate	Difference (95% CI)
30 Days	3 / 286 (1.0%)	1 / 155 (0.6%)	0.4% (-1.3%, 2.1%)
6 Months	14 / 280 (5.0%)	6 / 150 (4.0%)	1.0% (-3.0%, 5.0%)
12 Months	23 / 270 (8.5%)	11 / 139 (7.9%)	0.6% (-5.0%, 6.2%)
24 Months	38 / 247 (15.4%)	16 / 124 (12.9%)	2.5% (-4.9%, 9.9%)
36 Months	47 / 200 (23.5%)	23 / 94 (24.5%)	-1.0% (-11.5%, 9.5%)



Subgroup Analyses

Subgroup Analyses



- For most subgroups, there was no statistically significant interaction effect (p-value threshold of 0.15) for key variables
- No interaction for:
 - Gender
 - Geographic & Site Location (US vs. non-US; EU vs. Japan vs. OUS)
 - Rutherford Category (3 vs. 4 vs. 5)
- FDA did not identify baseline patient or lesion characteristic subgroups in which the Lutonix DCB demonstrated a differential benefit vs PTA



Pivotal Study Summary

Pivotal Trial Summary (1)



- Primary Safety Endpoint - **Met**
 - Freedom from BTK MALE (major adverse limb event) + POD (peri-operative death) at 30 days (ITT Per-Patient) was 99.3% for DCB vs. 99.4% for PTA meeting the 12% non-inferiority margin ($p < 0.0001$).
- Primary Effectiveness Endpoint - **Not Met**
 - Limb salvage + primary patency at 6 months (ITT per-flow pathway) was 74.7% for DCB vs. 64.2% for PTA ($p = 0.0222$; $p \leq 0.0085$ required for superiority).
 - 10.5% effect size in favor of the DCB

Pivotal Trial Summary (2)



- Longer-term primary effectiveness endpoint rates favored the PTA group at 12 months through 36 months
- FDA did not observe trends in favor of the Lutonix DCB vs. PTA for:
 - Freedom from CD-TLR at 12 months and beyond
 - Mortality
 - Major amputations
 - Wound healing and infections
 - Hemodynamic parameters: ABI & TBI
 - Rutherford Classification
 - Quality of life and walking impairment

Pivotal Trial Summary (3)



Study execution and design limitations

- Trial terminated early prior to enrolling the required sample size
- Multiple protocol modifications during the ongoing study
- High rate of missing primary effectiveness endpoint data
- Wound healing assessments flawed by important methodologic limitations

Pivotal Study Conclusions



The Lutonix 014 DCB was associated with:

- A short-term lower rate of target lesion revascularization vs. PTA, which is of uncertain clinical significance
- No clear evidence to FDA for additional clinical benefits vs. PTA for a durable reduced repeat revascularization rate, target lesion patency, arterial hemodynamics, major amputations, wound healing, Rutherford classification, or quality of life



Summary of Adjunctive Data Sources

**Eleni Whatley, Ph.D.
Biomedical Engineer
CDRH/OPEQ/OCVD/PMDT**

Adjunctive Data Sources



- The Sponsor provided additional datasets intended to support the safety and effectiveness of the Lutonix 014 DCB:
 - Global BTK Real-World Registry
 - Pooled analysis of the BTK IDE Pivotal Trial and the Global BTK Real-World Registry
 - Real-world data from the Society of Vascular Surgery (SVS) Vascular Quality Initiative (VQI) Database
 - Japan HD RCT
- The sponsor provided literature sources from small, single center studies.



Global BTK Real-World Registry

- Primary safety endpoint: Composite all-cause death, above-ankle amputation, or major re-intervention at 30 days
- Primary effectiveness endpoint: Freedom from CD-TLR at 6 months
- Results:

Measure	BTK Registry (N=371)	
	Freedom from primary safety events % (n/N)	95% CI ¹
Primary Safety Endpoint	98.3% (354/360)	96.4%, 99.4%

¹ Exact binomial confidence interval

Measure	BTK Registry (N=371)	
	Success % (n/N)	95% CI ¹
TLR-Free at 12 Months	79.9% (239/299)	74.9%, 84.3%
TLR-Free at 24 Months	74.2% (187/252)	68.3%, 79.5%

¹ Exact binomial confidence interval. All confidence intervals are based on nominal levels and not adjusted for multiple comparisons.

Measure	BTK Registry (N=371)	
	TLR-Free % (n/N)	95% CI ¹
6-Month Primary Effectiveness Endpoint	90.0% (289/321)	86.2%, 93.1%

¹ Exact binomial confidence interval

Global BTK Real-World Registry



- Strengths: Large sample size, real world subjects, prospectively designed, CEC Adjudicated
- Limitations: No active control, limited objective evaluations (e.g., no imaging with patency being determined by investigator), and high rates of missing data

Pooled analysis - IDE Trial and Global BTK Registry



Pooled DCB data (BTK IDE
Trial + BTK Global Registry)

vs.

PTA data (BTK IDE Trial)


- The pooled dataset included 658 DCB subjects and 155 PTA subjects
- Primary endpoints:
 - Safety: Composite all-cause death, above-ankle amputation, or major re-intervention at 30 days
 - Effectiveness: Composite of Limb Salvage and Primary Patency at 6 months

Pooled analysis - IDE Trial and Global BTK Registry



	Freedom from Primary Effectiveness Event	
Time Point	DCB (95% CI)	PTA (95% CI)
30 Days	95.7% (93.7%, 97.0%)	86.3% (79.8%, 91.0%)
6 Months	76.6% (72.9%, 79.9%)	58.9% (50.5%, 66.9%)
12 Months	58.1% (53.8%, 62.2%)	52.3% (43.8%, 60.7%)
24 Months	36.8% (32.1%, 41.7%)	39.8% (31.7%, 48.6%)

Limb Salvage and Primary Patency



Pooled analysis - IDE Trial and Global BTK Registry



- Strengths:
 - Larger sample size that includes real world subjects
 - PTA control group from the IDE trial
- Limitations:
 - Retrospective analysis of known data from IDE trial and global registry
 - Important differences in study oversight, eligibility criteria, and primary endpoint definitions
 - Concerns with propensity score methodology
 - Final SAP did not contain sufficient details for the propensity score method
 - Insufficient overlap present in the propensity score distribution between treatment groups
 - Subjects with missing outcome data excluded




Real-world data from the SVS VQI Database

- Real-world and off-label use of the 4 mm diameter SFA product that is approved in the US utilizing the VQI Peripheral Vascular Reintervention database
- Included 167 consecutive Lutonix DCB cases and 397 consecutive PTA cases who were propensity matched
- The primary endpoint assessment was freedom from TLR at 6 months.

SVS VQI Data



Day	DCB TLR Free Estimate (95% CI)	PTA TLR Free Estimate (95% CI)
0	100.0% (100.0%, 100.0%)	100.0% (100.0%, 100.0%)
30	98.9% (94.7%, 99.8%)	97.2% (94.9%, 98.5%)
180	96.1% (91.0%, 98.4%)	95.2% (92.5%, 97.0%)
365	91.8% (85.1%, 95.6%)	88.6% (84.3%, 91.7%)



- Strengths: Real world subjects, non-randomized comparator
- Limitations: Off-label use of one device size, protocol concerns, and statistical issues (e.g., ambiguity regarding patient selection, protocol changes to primary analysis method, and issues with how missing data were handled)

Adjunctive Data Sources



- The robustness of these studies were hindered by one or more of the following limitations:
 - Small sample sizes
 - No active controls
 - Retrospective analyses of known data
 - Unresolved statistical concerns
 - Short term follow up

FDA considered these additional studies in the totality of the data, but they did not result in any different findings or new conclusions, and they did not provide additional support beyond the pivotal IDE trial



Post Approval Study Proposal

Regulatory Purpose



- The requirement for a post-approval (PAS) study does not alter the requirements for premarket approval
- Pre-market data must reach the threshold for providing a reasonable assurance of safety and effectiveness
- Discussion regarding the proposed PAS does not indicate that FDA has made a decision or is making a recommendation on the approvability of this PMA



Post Approval Study (PAS) Proposal

- Continued follow-up PAS for their pivotal and global cohorts
- New Enrollment:
 - Vascular Quality Initiative (VQI) prospective data collection
 - Approximately 200 subjects
 - 1 year primary endpoint of TLR and vital status through 5 years
 - Success criteria and questions to be addressed are unclear



Benefit/Risk Considerations and Conclusions

Benefit/Risk Considerations



FDA determines whether the totality of the data provide a “reasonable assurance of safety and effectiveness” by weighing any probable benefit to health from the use of the device against any probable risk of injury or illness from such use.

- A reasonable assurance of safety occurs when “it can be determined, based upon valid scientific evidence, that the probable benefits... outweigh any probable risks,” and can be demonstrated by establishing “**the absence of unreasonable risk of illness or injury** associated with the use of the device for its intended uses and conditions of use.” (21 CFR 860.7(d)(1))
- A reasonable assurance of effectiveness occurs when “it can be determined, based upon valid scientific evidence, that in a significant portion of the target population, the use of the device for its intended uses... will provide **clinically significant results.**” (21 CFR 860.7(e)(1))

Benefit/Risk Considerations



Some factors that should be considered when evaluating benefits and risk of a medical device include:

- Benefits
 - Type of benefit
 - Magnitude of benefit
 - Likelihood of patients experiencing one or more benefit
 - Duration of effects
 - Patient perspective on benefit
 - Medical Necessity
- Risks
 - Severity of harm
 - Likelihood of risk
 - Duration of exposure to population
 - Patient tolerance of risk

Uncertainty in Benefit-Risk Determinations



- The appropriate extent of uncertainty regarding the benefits and risks in a given case will depend on consideration of multiple factors and other relevant information concerning the device, including:
 - The extent of the probable benefits and risks of the device
 - The extent of the public health need
 - The probable benefits of earlier patient access to the device
 - The ability to reduce uncertainty post-market
 - The likely effectiveness of mitigations, such as labeling, and other tools

Discussion of Benefits



- A 10.5% numerical improvement for the Lutonix 014 DCB was observed for the primary effectiveness endpoint at 6 months
 - This improvement was no longer evident by 12 months and the event rate numerically favored PTA at 12 months and beyond
 - FDA did not observe any evidence of added benefit of the Lutonix 014 DCB vs. PTA from pre-specified and post hoc secondary effectiveness endpoint evaluations including for:
 - Freedom from CD-TLR at 12 months and beyond
 - Major amputations
 - Wound healing and infections
 - ABI & TBI
 - Rutherford Classification
 - Quality of life and walking impairment

Discussion of Risks



- No additional safety risks associated with the Lutonix DCB vs. PTA identified in the current study data
- A safety signal for increased late mortality has been observed in patients treated with paclitaxel-coated devices in the superficial femoral artery
 - This trend was not seen in the current BTK study
- Potential concerns associated with paclitaxel-coated devices, including mortality, should be considered if there is uncertainty regarding clinical benefits

Benefit Risk Considerations for the Panel



FDA is requesting Panel input on whether:

- The data demonstrate a reasonable assurance of effectiveness for the Lutonix 014 DCB
 - A reasonable assurance of effectiveness occurs when “it can be determined, based upon valid scientific evidence, that in a significant portion of the target population, the use of the device for its intended uses... will provide **clinically significant results.**”
- The type, magnitude, and duration of the probable benefits of the Lutonix 014 DCB outweigh the known and probable risks



Thank you!

Questions?