

OCS™ Heart System for the Resuscitation, Preservation, and Assessment of Donor Hearts

April 6, 2021

Circulatory System Devices Panel



Introduction

Waleed Hassanein, MD

President and CEO

TransMedics

TransMedics Introduction



Founded in 1998 to develop Organ Care System (OCS™) technology to increase donor organ utilization for transplantation and improve post-transplant clinical outcomes

Clinically driven organization that pioneered concept of extracorporeal perfusion of donor hearts, lungs, and livers for transplantation


- Sponsored 8 US FDA pivotal trials

The OCS is developed and manufactured in US

- OCS Lung FDA approved
- OCS Liver under review by FDA
- OCS Heart approved internationally and > 1,000 cases transplanted to date worldwide



Only ~30% of Donor Hearts Are Used for Heart Transplants




Significant underutilization of deceased donor hearts

Limiting access to patients in need for heart transplantation

Patients waiting are not guaranteed a heart




Longer Wait Times Associated with Higher Mortality



Cold Storage Limits Utilization of Donor Organs and Shown to Negatively Impact Post-Transplant Outcomes



Cold Storage

-  Severe time-dependent injury (ischemia)
-  No organ optimization capabilities
-  No assessment of organ viability



Only 3 out of 10
DBD hearts used¹



~15-31%²
PGD
Primary Graft
Dysfunction (PGD)



OCS Heart System: Integrated Portable Platform Designed to Address Limitations of Cold Storage



OCS™ Heart Console



OCS™ Heart Perfusion Set




OCS™ Heart Solutions

VIDEO




Evolution of OCS Heart – 3 FDA Pivotal Trials for 3 Different Clinical Indications

OCS PROCEED II Trial –
Standard-Criteria Heart Donors




Primary Dataset for this PMA

OCS Heart EXPAND + CAP Trials –
Extended-Criteria Heart Donors



OCS Heart DCD + CAP Trials –
DCD Heart Donors





Proposed Indication for Use Consistent with Study Criteria

The TransMedics® Organ Care System (OCS™) Heart System is a portable extracorporeal heart perfusion and monitoring system indicated for the resuscitation, preservation, and assessment of donor hearts intended for a potential transplant recipient in a near-physiologic, normothermic and beating state. OCS Heart is indicated for donor hearts with one or more of the following characteristics:

- Expected cross-clamp or ischemic time ≥ 4 hours due to donor or recipient characteristics (e.g., donor-recipient geographical distance, expected recipient surgical time); or
- Expected cross-clamp or ischemic time ≥ 2 hours AND one or more of the following:
 - Donor Age ≥ 55 years; or
 - Donors with history of cardiac arrest and downtime ≥ 20 minutes; or
 - Donor history of alcoholism; or
 - Donor history of diabetes; or
 - Donor Left Ventricle Ejection Fraction (LVEF) $\leq 50\%$ but $\geq 40\%$; or
 - Donor history of Left Ventricular Hypertrophy (LVH) (septal or posterior wall thickness of $> 12 \leq 16$ mm); or
 - Donor angiogram with luminal irregularities but no significant coronary artery disease (CAD)



Overview of Heart EXPAND + CAP Results

Primary effectiveness endpoint met in EXPAND ($p < 0.0001$)

84% of extended-criteria donor hearts (**refused an average of 60 times**) were successfully transplanted using the OCS Heart System

8% ISHLT severe PGD well below rates reported in literature

97% all-cause patient survival at 30-days post-transplant is comparable to routine heart transplant outcomes (96%; Colvin et al, 2020)

92% & 87% all-cause patient survival at 6 & 12 months, respectively

96% cardiac-related patient survival at both 6 & 12 months



Agenda

Clinical Need to Expand Donor Heart Utilization

Maryjane Farr, MD

Associate Professor of Medicine, Division of Cardiology
Medical Director, Adult Heart Transplant Program
Columbia University Medical Center / New York-Presbyterian Hospital

Heart EXPAND & EXPAND CAP Trials

Jacob Schroder, MD

Assistant Professor of Surgery
Surgical Director, Heart Transplantation Program
Duke University School of Medicine

PROCEED II Trial Summary

Waleed Hassanein, MD

Statistical Considerations for Long-Term Survival Modeling

Chris Mullin, MS

Director, Product Development Strategy, NAMSA

TransMedics Position on FDA Questions Training & Post-Approval Programs

Waleed Hassanein, MD

Clinical Perspective & Benefit-Risk Assessment

Ashish Shah, MD

Professor of Cardiac Surgery
Alfred Blalock Endowed Director and Chairman, Department of Cardiac Surgery
Vanderbilt University Medical Center



Additional Experts

Pathology

Anthony J Demetris, MD

Starzl Professor of Liver and Transplant Pathology
University of Pittsburgh

Regulatory

Miriam Provost, PhD

Vice President, Global Regulatory Affairs
TransMedics, Inc.

Medical Monitor

John Wallwork, CBE, FRCS, FMedSci

Emeritus Professor, Cardiothoracic Surgery
Royal Papworth Hospital, Cambridge University, UK



Clinical Need to Expand Donor Heart Utilization

Maryjane Farr, MD

Associate Professor of Medicine, Division of Cardiology
Medical Director, Adult Heart Transplant Program

Columbia University Irving Medical Center /
New York-Presbyterian Hospital



End-Stage Heart Failure – Major Public Health Issue

6.5 million adults in U.S. have heart failure¹

46% increase in heart failure prevalence estimated by 2030¹


5-10% of patients with heart failure are “end-stage” or “advanced”¹

40-60% 1-year mortality rate for patients with end-stage heart failure²

Durable Left Ventricular Assist Devices Used as Bridge to Heart Transplant

Ventricular Assist Devices (VADs)


- Beneficial for selected patients
 - Bridge to transplant
 - Destination therapy if LV failure only
- 2-year survival HM3 79%
- Complications
 - Stroke 10%
 - Bleeding 24.5%
 - Infection 58.3%
 - Right ventricular failure 34.2%



Cardiac Transplant is Gold Standard

Cardiac Transplant

- Definitive replacement therapy
 - Requires intensive long-term care
- 88-92% survival at 1 yr, 72-80% at 5 yrs
- Improves
 - Functional status
 - Health-related quality of life





Heart Transplant Challenges Due to Supply

12,588 deceased organ donors in 2020 in U.S.¹


3,658 heart transplants in 2020 in U.S.¹



Only ~3 hearts out of every 10 donated hearts are used for transplant



Oversold Transplant Waiting List



Still Dependent on Cooler Despite Medical Advances



COLD STORAGE

Limitations of Cold Storage

- Donor hearts can only safely be preserved for ~4 hours
- Time-dependent ischemia
- No ability to provide therapeutic intervention to organ
- No ability to assess heart function
- Limited to standard-criteria hearts

Recent Federal Changes Mandate that Organs Be Allocated to Sickest Patients First

- Jan 2020 – UNOS with HHS determined local allocation not consistent with fair access
- Now able to transplant to sickest patients first
 - Longer travel times
 - Longer ischemic times
- Need new technologies to mitigate adverse effects of long ischemic times



Growing Pool of Lost Hearts to DCD Donation

- 2019 – 1,543 DCD donors between ages 18-49 years old
 - ≥ 1 organ (liver, kidney, lung) transplanted
 - Only 7 hearts able to be utilized
 - All successfully transplanted using OCS Heart System in clinical trial (IDE G180272; pivotal trial ongoing)



Unmet Need to Address Limitations of Cold Storage

Heart transplant is the gold standard therapy for end-stage heart failure

Cold storage is the only available option for preservation despite severe limitations that restrict utilization to ~3 out 10 donor hearts

Significant unmet need for new heart preservation technologies to address limitations of cold storage and increase number of life-saving heart transplants



Heart EXPAND & EXPAND CAP Trials

Jacob Schroder, MD

Assistant Professor of Surgery

Surgical Director, Heart Transplantation Program

Duke University School of Medicine



Heart EXPAND Trial Overview

Objective

Evaluate ability of OCS Heart System to significantly increase utilization of donor hearts that are rarely transplanted today due to limitations of cold storage

Target Enrollment

Extended-criteria donor hearts in typical heart recipients

Trial Design

Single-arm, multicenter U.S. clinical trial

Why a Single-Arm Trial

Unethical to randomize extended-criteria donor hearts to cold storage

Concurrent controls would be of limited clinical value as this is a different donor heart population

Donor Heart Eligibility Criteria

- Expected total cross-clamp time \geq 4 hours; OR
- Expected total cross-clamp time \geq 2 hours plus \geq 1 additional risk factor
 - Donor age 45-55 years with no coronary catheterization data
 - Donor age \geq 55 years
 - Left ventricular septal or posterior wall thickness of > 12 and ≤ 16 mm
 - Reported down time of ≥ 20 min with stable hemodynamics at final assessment
 - Left heart ejection fraction 40-50%
 - Donor angiogram with luminal irregularities with no significant CAD
 - History of carbon monoxide poisoning with good cardiac function at time of donor assessment
 - Social history of alcoholism with good cardiac function at time of donor assessment
 - History of diabetes with negative coronary angiogram for CAD



Donor Heart on OCS Transplant Acceptance Criteria

- OCS arterial lactate levels < 5 mmol/L at end of OCS perfusion period with stable lactate trend
- Recommended ranges for OCS heart perfusion parameters after stabilization
 - Coronary flow (CF) 400-900 mL/min
 - Aortic pressure (AOP) 40-100 mmHg
- Clinically satisfied with donor heart evaluation on OCS

**Recipient's Surgical Procedure Should Not Be Initiated
Until Donor Heart Accepted on OCS**

Role of OCS Parameters and Clinical Judgment in Donor Heart Acceptance

- Arterial lactate is an important biomarker for myocardial ischemia
- Prospective analysis of OCS perfusion parameters to predict graft failure¹
 - First 49 patients transplanted on OCS with standard-criteria hearts
 - Lactate trend, rate of change, and ending lactate had high sensitivity and specificity
 - 5 mmol/L limit for standard-criteria hearts (baseline lactate \leq 1 mmol/L)
- In EXPAND, a stable lactate trend added as qualifier because extended-criteria hearts may have higher starting lactate values (eg, 3-4 mmol/L)
 - Some hearts turned down in EXPAND because lactate increased despite attempts to maximize perfusion and AOP (certain cases where absolute lactate $<$ 5 mmol/L)
- Other OCS parameters and clinical judgement also used to determine transplantability
- Providing additional data to inform clinical judgment is benefit of OCS not a risk



Recipient Inclusion and Exclusion Criteria

Recipient Inclusion Criteria

- Registered male or female primary heart transplant candidate
- Age \geq 18 years
- Provided informed consent

Recipient Exclusion Criteria

- Prior solid organ or bone marrow transplant
- Chronic use of hemodialysis or diagnosis of chronic renal insufficiency
- Multi-organ transplant required

Primary Effectiveness and Safety Endpoints

- Primary effectiveness composite endpoint
 - Survival at Day 30
 - Absence of ISHLT severe PGD in first 24 hours
 - Performance goal of 65% (assumed rate of 80%)
- Primary safety endpoint
 - Incidence of heart graft-related SAEs (HGRSAEs) in first 30 days



Rationale for Performance Goals in EXPAND

No published literature on rate of PGD in extended-criteria donor heart transplants

Higher range of published PGD of ~30% added to the 5% rate of 30-day mortality for standard-criteria hearts to derive the 65% PG

With sample size of 75 patients, success rate assumed to be 80% to meet 65% PG

No prior data on moderate and severe ISHLT PGD through 30 days post transplant
Primary effectiveness endpoint includes safety components (survival and PGD)

Secondary Endpoints

- Patient survival at Day 30
- Incidence of severe PGD in first 24 hours
- Rate of donor heart utilization for transplantation


PGD Assessment and Adjudication

- Moderate or severe PGD defined in protocol according to ISHLT consensus recommendations
- Clinical events independently adjudicated by medical monitor
 - Dr. John Wallwork
 - Founding Member and Past President, ISHLT*
 - Emeritus Professor of Cardiothoracic Surgery, Papworth Hospital and Cambridge University*



Donor Heart Risk Factors for EXPAND Trial Inclusion – 52% of Donors Had 2 or More Risk Factors

93 Donor Hearts Included






UNOS Heart Transplant Database Validated EXPAND Criteria Are Seldomly Utilized in the U.S.


| Donor Risk Factors / Inclusion Criteria | OCS EXPAND (N=93) | UNOS SRTR Data* (N=10,426) | p-value |
|---|----------------------|-------------------------------|----------|
| Expected cross-clamp time \geq 4 hours | 37 (39.8%) | 1,607 (15.4%) | < 0.0001 |
| Downtime \geq 20 minutes | 33 (35.5%) | 240 (2.3%) | < 0.0001 |
| LVEF \geq 40% and \leq 50% | 24 (25.8%) | 481 (4.6%) | < 0.0001 |
| Donor age \geq 55 years | 11 (11.8%) | 295 (2.8%) | < 0.0001 |
| LVH $>12 \leq 16$ mm | 18 (19.4%) | Not collected | |
| 2 or More Risk Factors / Inclusion Criteria | | | |
| Cross-clamp \geq 4 hours plus other factors | 13 (14.0%) | 464 (4.5%) | 0.0003 |
| Downtime \geq 20 minutes plus other factors | 9 (9.7%) | 58 (0.6%) | < 0.0001 |
| Donor age \geq 55 years plus other factors | 7 (7.5%) | 104 (1.0%) | < 0.0001 |



EXPAND Enrolled Donor Hearts with Significantly More Match Run Refusals




81% of Donor Hearts Utilized for Transplantation



Reasons for Turn Down on OCS

- Continuous rising lactate and final lactate ≥ 5 mmol/L (n=8)
- Continuous rising lactate (n=7)
- Continuous rising lactate and RV dysfunction (n=2)
- Continuous rising lactate and inability to wean off pacing (n=1)

Arterial Lactate Trend Served as a Key Indicator for Assessment of Donor Hearts on OCS



Demographics and Baseline Characteristics of Transplanted Recipient Population

| Recipient Characteristics | (N=75) |
|---|-----------------|
| Baseline Characteristics | |
| Age (years), mean \pm SD | 55.5 \pm 12.6 |
| Male, n (%) | 61 (81%) |
| BMI (kg/m ²), mean \pm SD | 27.7 \pm 4.7 |
| Risk Factors, n (%) | |
| Age > 65 | 18 (24%) |
| History of mechanical circulatory support | 48 (64%) |
| Female donor to male recipient mismatch | 12 (16%) |
| Renal dysfunction | 11 (15%) |
| Status, n (%) | |
| Status 1A | 52 (69%) |
| Status 1B | 22 (29%) |
| Status 2 | 1 (1%) |


Donor Heart Preservation Characteristics

| Parameter | Transplanted Donor Hearts (N=75) |
|-----------------------------------|-------------------------------------|
| Cross-clamp time (hours) | |
| Mean \pm SD | 6.3 \pm 1.6 |
| Min – Max | 2.9 – 11.4 |
| Cold ischemic time (hours) | |
| Mean \pm SD | 1.7 \pm 0.4 |
| Min – Max | 1.1 – 2.8 |

85% of hearts transplanted in US < 4 hour maximum with cold storage*



Primary Effectiveness Endpoint Met






Favorable Results on Secondary Endpoints

| Secondary Endpoints (Components of Primary Composite Endpoint) | (N=75) |
|--|-----------------------|
| Patient survival at day 30 post transplantation | |
| Proportion (n/N) | 94.6% (70/74*) |
| Severe PGD (left or right ventricle) in first 24 hours post transplantation | |
| Proportion (n/N) | 10.7% (8/75) |

*One patient had graft failure requiring re-transplant



PGD with OCS Heart Similar or Lower than Other Studies




Important Considerations for Assessing Mortality in Heart Transplantation

- Most patients undergoing heart transplant are not otherwise healthy
 - Majority previously on VADs with associated complications
 - On long-term immunosuppressives
- Mortality in initial post-transplant period likely related to transplant procedure or cardiac graft
- After initial post-transplant period, recipients subject to competing risks for non-cardiac causes of death
- Cardiac-related survival is a clinically appropriate endpoint to assess preservation technology




All-Cause and Cardiac-Related Long-Term Survival



| Month | 0 | 1 | 6 | 12 | 18 | 24 |
|-----------------|----|----|----|----|----|----|
| Overall | 75 | 70 | 65 | 59 | 50 | 30 |
| Cardiac-related | 75 | 70 | 65 | 59 | 50 | 30 |

Causes of Death in EXPAND Trial



* This death was outside of the adjudication window; however, the original death note was obtained from the site



Low Incidence of Primary Safety Endpoint Events


| | (N=75) |
|--|----------------|
| Primary safety endpoint | |
| Mean \pm SD | 0.2 \pm 0.37 |
| HGRSAEs by type | |
| Moderate or severe PGD, n/N (%) | 11/75 (14.7%) |
| Primary graft failure requiring re-transplantation | 1/75 (1.3%) |



EXPAND CAP & Pooled EXPAND + CAP Analyses



Findings from EXPAND Continued Access Protocol Further Support Effectiveness and Safety of OCS Heart System



41 Patients with \geq 30-day follow-up


91% Utilization for transplantation

100% 30-day survival

2.4% Incidence of severe PGD



Pooled Results from EXPAND + CAP Support Effectiveness and Safety of OCS Heart System




116 Patients

84% Utilization for transplantation

97% 30-day survival

8% Incidence of severe PGD

Pooled EXPAND + CAP Long-Term Survival



| Month | 0 | 1 | 2 | 3 | 4 | 5 | 6 | 12 |
|------------------------|-----|-----|-----|----|----|----|----|----|
| Overall | 116 | 111 | 100 | 98 | 97 | 96 | 85 | 65 |
| Cardiac-related | 116 | 111 | 100 | 98 | 97 | 96 | 85 | 65 |



Data from EXPAND + CAP Trials Provide Substantial Evidence of Safety and Effectiveness of OCS Heart System

Primary effectiveness endpoint met in EXPAND ($p < 0.0001$)

84% of extended-criteria donor hearts (**refused an average of 60 times**) were successfully transplanted using the OCS Heart System

8% ISHLT severe PGD well below rates lower reported in literature

97% all-cause patient survival at 30-days post-transplant is comparable to routine heart transplant outcomes (96%; Colvin et al, 2020)

92% & 87% all-cause patient survival at 6 & 12 months, respectively

96% cardiac-related patient survival at both 6 & 12 months



PROCEED II Trial Summary

Waleed Hassanein, MD

President and CEO

TransMedics



PROCEED II – First Trial Designed for Any Ex-Vivo Organ Perfusion Technology and for OCS Heart System


Designed based on 2001 Celsior cold preservation solution trial design: non-inferiority RCT with 30-day follow-up for standard-criteria hearts

Met primary effectiveness and safety endpoints

Unplanned, post-hoc UNOS follow-up revealed increased overall mortality; however comparable cardiac-related mortality in OCS arm vs control

Substantial differences in device design, use model, and donor and recipient populations limit the utility of PROCEED II data for this PMA

PROCEED II Enrolled Standard-Criteria Donors That Are Substantially Different from EXPAND & CAP Population




EXPAND and CAP Directly Support the Proposed Indication for Use, Not PROCEED II




Post-Hoc Long-Term Survival Trend from UNOS Registry

Overall Survival



Cardiac-related Survival




| Month | 0 | 6 | 12 | 24 | 36 | 48 | 60 |
|----------------|----|----|----|----|----|----|----|
| OCS | 56 | 49 | 45 | 41 | 36 | 33 | 32 |
| Control | 62 | 59 | 58 | 54 | 51 | 48 | 48 |


| Month | 0 | 6 | 12 | 24 | 36 | 48 | 60 |
|----------------|----|----|----|----|----|----|----|
| OCS | 56 | 49 | 45 | 41 | 36 | 33 | 32 |
| Control | 62 | 59 | 58 | 54 | 51 | 48 | 48 |



Causes of Mortality in PROCEED II Trial ≤ 60 Days Post Transplant





Causes of Mortality in PROCEED II Trial > 60 Days Post Transplant: UNOS Registry





Published International Studies Show Favorable OCS Long-Term Survival (N=165 OCS Cases)





Statistical Considerations for Long-Term Survival

Chris Mullin, MS

Director, Product Development Strategy
NAMSA

Background on FDA Models to Extrapolate Long-Term Survival

- FDA developed model based on preliminary PROCEED II data (2015)
- FDA used this model to extrapolate long-term survival results in:
 - PROCEED II
 - EXPAND

Concerns About Underlying Assumptions of FDA Model

- Parametric models use strong assumptions for underlying hazard rates
- Specific concerns with piecewise-exponential model
 - Model choice and cut points all post hoc
 - Data-driven cut points can lead to bias / error inflation
- No clinical justification provided for cut points selected by FDA
 - Calls into question validity of model



FDA Model Underestimates Long-Term Survival

Difference Between Kaplan-Meier and Piecewise Exponential Survival Estimates (Percentage Points)

| Time Post Transplantation | PROCEED II OCS Heart | PROCEED II SOC | EXPAND OCS Heart |
|---------------------------|----------------------|----------------|------------------|
| 1 Year | 1.0 | 0.1 | 0.0 |
| 2 Years | 0.3 | 1.0 | -0.1 |
| 3 Years | 1.1 | -1.4 | 0.3 |
| 4 Years | -1.4 | -3.2 | ? |
| 5 Years | -3.5 | -8.3 | ? |

Summary of Concerns About FDA's Model to Extrapolate Long-Term Survival

- Choice of model for extrapolation is questionable
- FDA's model is inaccurate at predicting long-term survival
- Results show FDA's model appears neither valid nor reliable for extrapolation of long-term data in heart transplantation



TransMedics Positions On FDA Questions

Waleed Hassanein, MD

President and CEO

TransMedics



Question 1: EXPAND Design & Conduct

EXPAND is Well Designed, Clinically and Statistically Robust. UNOS Registry Data Analysis Validated the Inclusion Criteria of Extended-Criteria Donor Hearts in the U.S.

Study Design

- Randomization of extended-criteria hearts to cold storage **is not ethical**
- Concurrent controls **could be obtained from UNOS database for U.S. transplants**
- PG predicated on observed success rate of **80%**
- Donor heart inclusion criteria **independently defined by experts and proven to be significantly different than standard U.S. heart transplants based on UNOS database analysis**



Question 1: EXPAND Design & Conduct

The EXPAND & CAP Trial Were Conducted According to Highest Clinical Standards and the Adjudication Process Strictly Adhered to Trial Protocol

Trial Conduct

- **There were no revisions to donor heart inclusion criteria** – additional criteria recorded in original source documents were retabulated to provide the complete picture of risk factors in EXPAND donor population
- **There were no changes to severe PGD definition** – Medical Monitor consistently and strictly followed ISHLT criteria throughout adjudication
- **Sensitivity analysis demonstrated robustness of the endpoint regardless of adjudication**





Question 2: EXPAND Inclusion Criteria

EXPAND & CAP Enrolled Donor Hearts That Are Seldomly Used for Transplants in the U.S. Based on UNOS Data Analysis

| Donor Characteristics | EXPAND + CAP (N=138) | UNOS/SRTR* (N=10,873) | p-value |
|--|-------------------------|--------------------------|----------|
| Age (years) – mean \pm SD | 36.4 \pm 12.1 | 32.1 \pm 11.0 | < 0.0001 |
| Age \geq 55 years | 13 (9%) | 309 (3%) | 0.0002 |
| LVH $>12 \leq 16$ mm | 18 (19%) | Not collected | |
| Cross-clamp time \geq 4 hours (Expected) | 66 (48%) | 1730 (16%) | < 0.0001 |
| Cross-clamp time \geq 4 hours (Actual) | 113 (97%) | 1730 (16%) | < 0.0001 |
| LVEF between 40% - 50% | 30 (22%) | 500 (5%) | < 0.0001 |
| Downtime \geq 20 minutes | 43 (31%) | 255 (2%) | < 0.0001 |
| Cross-clamp \geq 4 hours plus \geq 1 risk factor | 23 (17%) | 500 (5%) | < 0.0001 |
| Downtime \geq 20 minutes plus \geq 1 risk factor | 10 (7%) | 61 (1%) | < 0.0001 |

EXPAND Evaluated Donor Hearts That Are Not Routinely Transplanted in the US Today





Question 3: Transplantability

Lactate Level and Trend Are Useful Guides to Managing Perfusion of Donor Hearts on OCS in Conjunction with AOP, CF and Overall Clinical Judgment

Prior clinical data have demonstrated that **lactate is a sensitive (63%) and highly specific (98%) biomarker for graft dysfunction¹**



Lactate level/trend is **NOT** the only parameter or consideration for transplantability
Used in conjunction with AOP, CF, and clinical judgment of the overall clinical test condition of the donor heart on OCS Heart System

This OCS Heart use model has been successfully used internationally with excellent clinical outcomes in DBD and DCD donor heart transplants



Question 4a & 4b: PROCEED II and EXPAND Study Analysis

Favorable Long-Term Survival for EXPAND + CAP Support the Approval of OCS for the Proposed Indication for Extended-Criteria Donor Hearts






Question 4c: Donor Hearts with ≥ 4 Hours Cross-Clamp Time

OCS Heart is Safe and Effective for Donor Hearts with ≥ 4 Hours of Expected Cross-Clamp Time

≥ 4 -Hour Cross-clamp Criteria in EXPAND & CAP Population (N=33)



All 4 Deaths Were Not Cardiac-related

- Vasoplegia leading multiorgan failure
- Complications of pre-existing liver cirrhosis
- Non-recoverable CVA
- Motor vehicle accident (14 months)

| Month | 0 | 1 | 2 | 3 | 4 | 5 | 6 | 12 |
|-----------------|----|----|----|----|----|----|----|----|
| Overall | 33 | 31 | 26 | 26 | 26 | 26 | 21 | 15 |
| Cardiac-related | 33 | 31 | 26 | 26 | 26 | 26 | 21 | 15 |



Question 5: Pathophysiology and Pathology

No Definitive Clinical Evidence of OCS-Related Injury of Donor Hearts. OCS Ability to Turn Down Potentially Bad Extended-Criteria Donor Hearts is a Clinical Benefit Not a Risk

- Donor hearts studied in EXPAND & CAP had significant risk factors making them highly unlikely to be used for transplantation – **many of these factors could contribute to pathological findings**
- **Brain death associated with significant physiologic changes** could show as pathological findings of a donor heart on histological examination of myocardium
- To our knowledge **there has never been any published or presented literature linking OCS Heart System to myocardial injury during perfusion**
- **84%** successful utilization of extended-criteria hearts is a significant clinical benefit



Question 6: Indications for Use


EXPAND & CAP Enrolled Donor Hearts That Are Seldomly Used for Transplants in the U.S. Based on UNOS Data Analysis and Match Run Refusals

The TransMedics® Organ Care System (OCS™) Heart System is a portable extracorporeal heart perfusion and monitoring system indicated for the resuscitation, preservation, and assessment of donor hearts intended for a potential transplant recipient in a near-physiologic, normothermic and beating state. OCS Heart is indicated for donor hearts with one or more of the following characteristics:

- Expected cross-clamp or ischemic time ≥ 4 hours due to donor or recipient characteristics (e.g., donor-recipient geographical distance, expected recipient surgical time); or
- Expected cross-clamp or ischemic time ≥ 2 hours AND one or more of the following:
 - Donor Age ≥ 55 years; or
 - Donors with history of cardiac arrest and downtime ≥ 20 minutes; or
 - Donor history of alcoholism; or
 - Donor history of diabetes; or
 - Donor Left Ventricle Ejection Fraction (LVEF) $\leq 50\%$ but $\geq 40\%$; or
 - Donor history of Left Ventricular Hypertrophy (LVH) (septal or posterior wall thickness of $> 12 \leq 16$ mm); or
 - Donor angiogram with luminal irregularities but no significant coronary artery disease (CAD)

Question 7: Benefit / Risk

84% Utilization of Extended-Criteria Donor Hearts Seldomly Used in the U.S. to Increase Heart Transplants Is a Significant Clinical and Public Health Benefit



Risk Profile of Turned-down Extended-criteria Hearts

- **72** mean match run refusals
- **55%** (12/22) had multiple donor heart risk factors



Training & Post-Approval Programs

Waleed Hassanein, MD

President and CEO

TransMedics

Dedicated Clinical Training Infrastructure




15,000 square-foot facility equipped with latest surgical and diagnostics equipment to replicate a retrieval environment



> 400 HCPs trained from 90 global academic and clinical institutions

Multi-Faceted Training Course



Classroom training for didactics and lessons learned



Demonstration and hands-on training for instrumentation and clinical management of the donor heart on OCS



Summary of OCS Heart Training Program



TRAINING

Initial hands-on clinical training and certification of every new clinical center starting OCS Heart program



SUPPORT

TransMedics provides 24/7 retrieval support via phone, messaging and email



TECHNOLOGY

Dedicated OCS Heart iPad® training and support application

Post Approval Program – FDA Question 8

New Enrollment Post-Approval Registry

- **175-patient**, single-arm, prospective post-approval registry
- Primary effectiveness endpoint: 12-month cardiac-related survival
 - PG of **86%**
 - Assumed success rate of **93%** required to meet the PG with 175 patients
- **12 months** of primary follow-up, then UNOS registry follow-up through 5 years

Continued Follow-up of EXPAND Patients

- UNOS registry follow-up of all EXPAND patients through 5 years
- Assess cardiac-related and all-cause survival



Clinical Perspective & Benefit-Risk Assessment

Ashish Shah, MD

Professor of Cardiac Surgery
Alfred Blalock Endowed Director and Chairman
Department of Cardiac Surgery
Vanderbilt University Medical Center



OCS Heart System Proposed Indication

The TransMedics® Organ Care System (OCS™) Heart System is a portable extracorporeal heart perfusion and monitoring system indicated for the resuscitation, preservation, and assessment of donor hearts intended for a potential transplant recipient in a near-physiologic, normothermic and beating state.


OCS Heart is indicated for donor hearts with one or more of the following characteristics:

- Expected cross-clamp or ischemic time ≥ 4 hours due to donor or recipient characteristics (e.g., donor-recipient geographical distance, expected recipient surgical time); or
- Expected cross-clamp or ischemic time ≥ 2 hours AND one or more of the following:
 - Donor Age ≥ 55 years; or
 - Donors with history of cardiac arrest and downtime ≥ 20 minutes; or
 - Donor history of alcoholism; or
 - Donor history of diabetes; or
 - Donor Left Ventricle Ejection Fraction (LVEF) $\leq 50\%$ but $\geq 40\%$; or
 - Donor history of Left Ventricular Hypertrophy (LVH) (septal or posterior wall thickness of $> 12 \leq 16$ mm); or
 - Donor angiogram with luminal irregularities but no significant coronary artery disease (CAD)




EXPAND Targeted Donor Hearts Not Typically Used for Transplantation in the U.S.

| Donor Characteristics | EXPAND + CAP (N=138) | UNOS/SRTR* (N=10,873) | p-value |
|--|-------------------------|--------------------------|----------|
| Age (years) – mean \pm SD | 36.4 \pm 12.1 | 32.1 \pm 11.0 | < 0.0001 |
| Age \geq 55 years | 13 (9%) | 309 (3%) | 0.0002 |
| LVH >12 \leq 16 mm | 18 (19%) | Not collected | |
| Cross-clamp time \geq 4 hours (Expected) | 66 (48%) | 1730 (16%) | < 0.0001 |
| Cross-clamp time \geq 4 hours (Actual) | 113 (97%) | 1730 (16%) | < 0.0001 |
| LVEF between 40% - 50% | 30 (22%) | 500 (5%) | < 0.0001 |
| Downtime \geq 20 minutes | 43 (31%) | 255 (2%) | < 0.0001 |
| Cross-clamp \geq 4 hours plus \geq 1 risk factor | 23 (17%) | 500 (5%) | < 0.0001 |
| Downtime \geq 20 minutes plus \geq 1 risk factor | 10 (7%) | 61 (1%) | < 0.0001 |



Increased Donor Heart Utilization for Transplantation – Huge Clinical Benefit for Heart Failure Patients in U.S.




- This high rate of donor heart utilization could double the number of transplants
- Significant clinical and public health benefit to increase heart transplant procedures in the U.S.



OCS Provides Additional Data to Inform Clinical Decision-Making on Transplantability vs Flying Blind with Cold Storage



Favorable Long-Term Clinical Outcomes in EXPAND + CAP





| Month | 0 | 1 | 2 | 3 | 4 | 5 | 6 | 12 |
|------------------------|-----|-----|-----|----|----|----|----|----|
| Overall | 116 | 111 | 100 | 98 | 97 | 96 | 85 | 65 |
| Cardiac-related | 116 | 111 | 100 | 98 | 97 | 96 | 85 | 65 |

*Stehlik et al, Circ Heart Fail 2017



OCS Heart Resulted in Distant Procurement of Donor Hearts That Could Not Be Achieved by Cold Storage





OCS May Enable National Sharing of Donor Hearts – Pushing Historical Boundaries of Cold Storage




*Recent case from EXPAND CAP. Data not yet submitted to the FDA.

OCS Heart Enables Use of Donor Hearts That Cannot Be Safely Preserved with Cold Storage

- 38-year-old donor who died from cerebrovascular hemorrhage
 - Distance from recipient hospital > **1,000 miles**
 - Heart turned down **327 times** by other transplant centers before acceptance in EXPAND study
- 60+-year-old recipient with cardiomyopathy
 - Blood type O
 - Status 1A
 - On LVAD for ~1 year prior to transplant
- Patient transplanted, discharged within 2 weeks, doing well 4 years post transplant



What Approval of OCS Heart Could Mean for Heart Transplant in the United States



OCS™ Heart System for the Resuscitation, Preservation, and Assessment of Donor Hearts

April 6, 2021

Circulatory System Devices Panel




Back up Slides Shown




Profile of Turned Down Hearts in EXPAND + CAP

Match Run Refusals



Inclusion Criteria





Summary of OCS Heart Training Program



TRAINING

Initial hands-on clinical training and certification of every new clinical center starting OCS Heart program



SUPPORT

TransMedics provides 24/7 retrieval support via phone, messaging and email



TECHNOLOGY

Dedicated OCS Heart iPad® training and support application



OCS May Enable National Sharing of Donor Hearts – Pushing Historical Boundaries of Cold Storage





1. Baran et al, Circ Heart Fail 2019

*Recent case from EXPAND CAP. Data not yet submitted to the FDA.




OCS Heart Resulted in Distant Procurement of Donor Hearts That Could Not Be Achieved by Cold Storage






Post-Hoc Analysis of Overall Survival in EXPAND + CAP vs SRTR Standard-Criteria Heart Transplant Recipients at Same Sites During Same Time Frame



| Month | 0 | 1 | 2 | 3 | 4 | 5 | 6 | 12 |
|-------------------|-------|-------|-------|-------|-------|-------|-------|-------|
| EXPAND + CAP OCS | 116 | 111 | 100 | 98 | 97 | 96 | 85 | 65 |
| SRTR Cold Storage | 1,813 | 1,752 | 1,724 | 1,707 | 1,695 | 1,687 | 1,626 | 1,284 |



Post-Hoc Analysis of Overall Survival in EXPAND + CAP vs SRTR Standard-Criteria Heart Transplant Recipients During Same Time Frame




| Month | 0 | 1 | 2 | 3 | 4 | 5 | 6 | 12 |
|-------------------|--------|--------|--------|--------|--------|--------|-------|-------|
| EXPAND + CAP OCS | 116 | 111 | 100 | 98 | 97 | 96 | 85 | 65 |
| SRTR Cold Storage | 10,873 | 10,484 | 10,330 | 10,240 | 10,172 | 10,116 | 9,689 | 7,855 |



Question 4c: Donor Hearts with ≥ 4 Hours Cross-Clamp Time

OCS Heart is Safe and Effective for Donor Hearts with ≥ 4 Hours of Expected Cross-Clamp Time

≥ 4 -Hour Cross-clamp Criteria in EXPAND & CAP Population (N=33)



All 4 Deaths Were Not Cardiac-related

- Vasoplegia leading multiorgan failure
- Complications of pre-existing liver cirrhosis
- Non-recoverable CVA
- Motor vehicle accident (14 months)


| Month | 0 | 1 | 2 | 3 | 4 | 5 | 6 | 12 |
|-----------------|----|----|----|----|----|----|----|----|
| Overall | 33 | 31 | 26 | 26 | 26 | 26 | 21 | 15 |
| Cardiac-related | 33 | 31 | 26 | 26 | 26 | 26 | 21 | 15 |

Data are Poolable: Similar Primary Effectiveness Endpoint Results Across Sites

| Site | Primary Effectiveness Endpoint | |
|---------|--------------------------------|-----|
| | n/N | % |
| Overall | 106/116 | 91 |
| C-12 | 2/2 | 100 |
| 06/C-01 | 48/53 | 91 |
| 08 | 1/1 | 100 |
| 09 | 7/7 | 100 |
| C-05 | 4/4 | 100 |
| 10 | 10/12 | 83 |
| 02 | 7/7 | 100 |
| C-03 | 5/5 | 100 |
| 04 | 1/1 | 100 |
| 05 | 1/2 | 50 |
| 03/C-11 | 12/14 | 86 |
| C-06 | 8/8 | 100 |

**P-value for heterogeneity
across sites = 0.91**
(sites with < 5 patients pooled)

EXPAND Evaluated Donor Hearts That Are Not Routinely Transplanted in the US Today





Proposed Indication for Use Consistent with Study Criteria


The TransMedics® Organ Care System (OCS™) Heart System is a portable extracorporeal heart perfusion and monitoring system indicated for the resuscitation, preservation, and assessment of donor hearts intended for a potential transplant recipient in a near-physiologic, normothermic and beating state.

OCS Heart is indicated for donor hearts with one or more of the following characteristics:

- Expected cross-clamp or ischemic time ≥ 4 hours due to donor or recipient characteristics (e.g., donor-recipient geographical distance, expected recipient surgical time); or
- Expected cross-clamp or ischemic time ≥ 2 hours AND one or more of the following:
 - Donor Age ≥ 55 years; or
 - Donors with history of cardiac arrest and downtime ≥ 20 minutes; or
 - Donor history of alcoholism; or
 - Donor history of diabetes; or
 - Donor Left Ventricle Ejection Fraction (LVEF) $\leq 50\%$ but $\geq 40\%$; or
 - Donor history of Left Ventricular Hypertrophy (LVH) (septal or posterior wall thickness of $> 12 \leq 16$ mm); or
 - Donor angiogram with luminal irregularities but no significant coronary artery disease (CAD)




Kaplan-Meier Analysis of Overall Survival Expand + CAP OCS Patients with ≥ 4 hours of Ischemic Time and Any Additional Inclusion Criteria



| Month | 0 | 1 | 6 | 12 |
|------------------|----|----|----|----|
| Expand + CAP OCS | 20 | 20 | 16 | 12 |



Geographic Distance Between Donor and Recipient Hospital for EXPAND and Proceed II





Conclusions from Histopathological Assessment of Turned-Down Donor Hearts in EXPAND

- Ischemia-reperfusion injury (IRI) is unavoidable ***regardless of preservation method***
- Most donor hearts already harbor some ischemic-type damage from peri-operative events^{1,2}
 - Histopathological findings do not strongly correlate with function unless damage is extensive and heart is reperfused¹⁻²
- Issue with cold storage: IRI is ***undetectable*** and is not observed ***until the donor heart is transplanted into the recipient***
- OCS offers benefit of allowing for proactive identification, monitoring, and responding to IRI *ex vivo* rather than reacting *in vivo* after transplant
 - Consistent with low rate of PGD in EXPAND and CAP vs literature rates with cold storage

No convincing evidence that OCS is “damaging” hearts – all observations consistent with IRI that would be present under any circumstances

Observed Myocyte Dissolutions Take More Time to Develop than OCS Perfusion Times



Animal Study Pathological Evaluation of OCS Preserved Hearts

- (n=40) were maintained on the OCS in a beating state for approximately 6 hours while being perfused with warm oxygenated blood supplemented with the TransMedics Maintenance Solution.
- Microscopic sections were taken from the ventricles (8), atria (2), pulmonary artery (1), aorta (1) and coronary arteries (2).
- Pathologists blinded to the experimental conditions, each slide was scored for ischemia, hemorrhage and edema on a scale from 0-4, and the scores for the ventricular slides.
- Pathologic evaluation validated the association of OCS enabled assessment parameters (CF, AOP, ending lactate levels) with ischemic damage.

352

Pathology Influences Device Development: The TransMedics Organ Care System for Heart

R.F. Padera¹, P. Lezberg², E. Hansen², D. Sousa², G.L. Winters¹


¹Brigham and Women's Hospital, Boston, MA; ²TransMedics, Inc., Andover, MA

Purpose: Pathologic evaluation was used to optimize the device development of the Organ Care System (OCS) device, its operating parameters and functional assessment parameters. The objective of this analysis was to determine the correlation between OCS perfusion parameters, metabolic measurements and histological assessment.

Methods and Materials: Porcine hearts (n=40) were maintained on the OCS in a beating state for approximately 6 hours while being



Kaplan-Meier Analysis of Overall Survival EXPAND + CAP OCS Patients vs. Proceed II Patients Who Met at Least One Expand Criteria



| Month | 0 | 1 | 6 | 12 | 24 |
|-------------------------|-----|-----|----|----|----|
| Expand + CAP OCS | 116 | 111 | 85 | 65 | 30 |
| Proceed II OCS | 17 | 16 | 15 | 15 | 13 |



Mechanical Support Use in PROCEED II

| | OCS (N=62) | Control (N=66) |
|---|------------------|-------------------|
| Post-Transplant Mechanical Circulatory Support (MCS) – Any type (IABP, ECMO, VAD etc.) | 9 (14.5%) | 7 (10.6%) |
| IABP Only | 3 | 5 |
| VAD Only | 0 | 1 |
| ECMO Only | 2 | 1 |
| ECMO and IABP | 1 | 0 |
| VAD and IABP | 2 | 0 |
| VAD and ECMO | 1 | 0 |
| Patients with MCS + death within 30 days | 3* | 1# |
| Patients who were discharged alive post MCS | 6 (9.7%) | 6 (9.1%) |

* OCS Deaths: 1 for PGD, 1 for Hyperacute Rejection and 1 for Acute Protamine Reaction

Control Deaths: 1 for PGD