

# Ischemic Stroke System for Treatment of Patients with Acute Ischemic Strokes

**December 10, 2021**

Neurologic Devices Panel

BrainsGate

1



## Introduction

**Jeffrey Saver, MD, FAAN, FAHA**

Director, UCLA Comprehensive Stroke and Vascular Neurology Program

SA Vice-Chair and Professor of Neurology, DGSOM

2

## Acute Ischemic Strokes are Devastating Events, Result in Lifetime of Disability and Reduced QoL

- Guidelines recommend timely reperfusion with pharmacologic IV-tPA or Endovascular Thrombectomy
  - Improve neurological outcomes
  - Reduces long-term disability
- Use is time dependent, requires administration within 6 hours
  - Treatment in 6-to-24-hour window limited to select patients

3

## ISS500 is First-of-a-Kind Treatment for Patients with Acute Ischemic Stroke

- Clinical program started in 2006
- 4 clinical trials
  - 2 randomized, sham-controlled
- > 1,400 patients enrolled
- 100 sites globally (11 US)

4

## Totality of Evidence Supports Positive Benefit-Risk for Ischemic Stroke System

- Patients with confirmed cortical involvement, treated 8 to 24 hours after stroke onset, achieved consistent improvements
  - Favorable disability outcomes
  - Improved long-term QoL
- Favorable safety profile
  - Reduced risk of symptomatic ICH
- Final system has efficient and reliable usability

**SPG Stimulation fills treatment gap for many patients who are ineligible or have no access to current reperfusion therapies**

5

## Proposed Indication for Ischemic Stroke System

*“The ISS500 is indicated to increase cerebral blood flow and reduce disability in adult patients with acute ischemic stroke with confirmed cortical involvement in the anterior circulation who are ineligible or have no access to IV-tPA and endovascular thrombectomy. Treatment is to be initiated between 8-24 hours from stroke onset (last known well).”*

6

# FDA Discussion Questions / Concerns Addressed in Presentation

CO-7

Topic	Keywords
Timing of major changes, including CCI subgroup	CCI addition
	Device change
Generalizability to US population	Generalizability US
Use of sliding dichotomous scale	Sliding dichotomy
mITT analysis	mITT analysis
Overall safety and potential AEs	Safety
Implantation skills	Implantation
Real-world selection of CCI patients	Patient selection



FDA Discussion “Keywords”

7

## Ischemic Stroke System

### Tom Devlin, MD, PhD, FSVIN

Director, CHI Memorial Stroke & Neuroscience Center  
 Director, Chattanooga Center for Neurologic Research LLC  
 Professor of Neurology, University of Tennessee Health Science Center

## Unmet Need & Pathophysiology

### Michael Hill, MD, MSc, FRCP

President, Canadian Neurological Sciences Federation  
 Professor, Dept Clinical Neuroscience & Hotchkiss Brain Institute  
 Cumming School of Medicine, University of Calgary & Foothills Medical Centre

## MoA of SPG Stimulation & Effectiveness & Safety Results

### Jeffrey Saver, MD, FAAN, FAHA

Director, UCLA Comprehensive Stroke and Vascular Neurology Program  
 SA Vice-Chair and Professor of Neurology, DGSOM

## Training & Post-Approval Plan

### Eyal Shai, EMBA

Chief Technology Officer  
 BrainsGate

## Clinical Perspective

### Michael Hill, MD, MSc, FRCP

CO-8

8

## Additional Experts

### Scott E. Kasner, MD MSCE

Ruth M. and Tristram C. Colket, Jr. President's Distinguished Professor  
 Vice Chair for Clinical Affairs  
 Department of Neurology, Perelman School of Medicine, University of Pennsylvania  
 Director, Comprehensive Stroke Center, University of Pennsylvania Health System

### Chris Mullin

Statistician  
 Director, Product Development Strategy  
 NAMSA

9



## Ischemic Stroke System

### Tom Devlin, MD, PhD, FSVIN

Director, CHI Memorial Stroke & Neuroscience Center  
 Director, Chattanooga Center for Neurologic Research LLC  
 Professor of Neurology, University of Tennessee Health  
 Science Center

10

# ISS500 System Components



## Neurostimulator

Positioned next to sphenopalatine ganglion (SPG) via natural canal in hard palate



## Navigation System and Injector

Facilitates safe and accurate positioning while minimizing risks or complications



## External Treatment System

Activates neurostimulator to deliver electrical pulses within predefined range

11

# Preparation for Implant Procedure

## 1 Patient Preparation

- Reference marker
- Nose marker


## 2 Implantation CT

## 3 Local Anesthesia

## 4 Planning




Reference Markers



Implanter Marks  
Canal Path on CT

12



13


**Simple, Technically Straightforward Procedure  
Poses Minimal Risk to Patient** CO-14

**Day 1  
Immediately after  
Implant Procedure**

**Day 5  
After Implant Removed**

14

# Current Implant Reduces Procedure Times and Implant Complications




	First Implant	Current Implant
Mechanical structure	Fragile	Flexible, strong
Requires trocars	Yes	No
Stimulation delivery	Identical	


All changes intended to optimize implant placement and improve overall safety; stimulation delivery - identical

FDA Concern: Device Change

# Optimal Dosing Range Identified During Clinical Development



- Final Treatment System delivers stimulation within predefined optimal range



- Stimulation level set based on physiologic signs of SPG activation





## Unmet Need and Pathophysiology


**Michael Hill, MD, MSc, FRCPC**

President, Canadian Neurological Sciences Federation  
Professor, Department Clinical Neuroscience &  
Hotchkiss Brain Institute  
Cumming School of Medicine, University of Calgary &  
Foothills Medical Centre

## Acute Ischemic Stroke is Major Cause of Death and Disability in United States


- > 800,000 patients experience stroke each year in US
  - ~85% of all strokes are ischemic
- 1<sup>st</sup> leading cause of acquired neurological disability
- 2<sup>nd</sup> leading cause of dementia
- 5<sup>th</sup> leading cause of death in US
  - 100,000 acute ischemic stroke deaths per year in US

# Strokes Cause Cascade of Cellular and Molecular Events, Leading to Extensive Damage



19

# ~50% of Patients Have Salvageable Tissue 24 Hours After Stroke




High variability of infarct progression by patient

- **Poor collaterals** “fast progressors”
- **Adequate collaterals** slower progression, more salvageable tissue

Wheeler, 2015

20

## Blood-Brain-Barrier Compromised Due to Stroke



Bernado-Castro, 2020


21

## Two Primary Treatment Options for Patients with Ischemic Stroke



### Pharmacologic IV-tPA

- First-line therapy
- Recommended initiation within 4.5 hours from stroke onset
  - Use restricted to within 3 hours in US
- Select patients may benefit up to 9 hours post-stroke



### Endovascular Thrombectomy (EVT)

- Endovascular approach to remove offending thrombus and restore anterograde perfusion
- Typically < 6 hours from onset
- Use limited to select patients between 6 and 24 hours

IV-tPA = intravenous tissue plasminogen activator

22


# Unmet Need for Patients with Acute Ischemic Stroke < 24 Hours

### IV-tPA Ineligibility

- Outside time window
- On anticoagulant

### EVT Ineligibility

- Clots too distal
- Tortuous anatomy
- Carotid artery chronically occluded
- Damage too extensive



Bahr Hosseini, 2018; Desai 2020



## SPG Stimulation and Data Supporting MoA


**Jeffrey Saver, MD, FAAN, FAHA**

Director, UCLA Comprehensive Stroke and Vascular Neurology Program

SA Vice-Chair and Professor of Neurology, DGSOM

# The Sphenopalatine Ganglion (SPG)


- Anatomy
  - Pterygopalatine fossa
  - Posterior to maxillary sinus
  - 3 mm, 60,000 nerve cells
- Components
  - Parasympathetic cell bodies
  - Traversing sympathetic and sensory
- Functions
  - Dilation of anterior cerebral circulation
  - Dilation of meningeal and dural vessels
  - Secretomotor function to nasopharyngeal mucosa and lacrimal gland




25

# Direct Reperfusion vs Collateral Flow Augmentation

EVT




SPG Stimulation




26

# SPG Stimulation Increases Blood Flow in Ischemic Penumbra




27

# SPG Stimulation Reduces Infarct Volume and Stabilizes Blood Brain Barrier




28

# Cerebral Blood Flow Increases in Inverse U-Shaped Dose Response



Levi H, 2012

# Case Study Shows Greatest SPG Effect on Cortical Region After Stimulation



## Clinical Trial Effectiveness & Safety Findings

31

## Clinical Studies – Over 1,400 Patients in 4 Trials

	ImpACT-1	ImpACT-24A	ImpACT-24B	ImpACT-24M
<b>Aim</b>	Feasibility Tolerability	Supportive	Pivotal	Usability
<b>Design</b>	Single Arm	Phase 2B RCT	Pivotal RCT	Single Arm
<b>Primary Endpoint(s)</b>	Rx Completion SAEs	mRS at 3 Months	mRS at 3 Months	Optimized Delivery CCA Flow/Hand Strength
<b>Dates</b>	2006 – 2008	2009 – 2011	2011 – 2018	2017 – 2018
<b>Patients (N)</b>	98	253	1000	50
<b>Centers</b>	12 OUS	6 US / 35 OUS	7 US / 66 OUS	4 OUS

32



## ImpACT-24B Pivotal Trial

33

## ImpACT-24B: Prospective Randomized, Double Blind, Sham-Controlled, Parallel Arm, Multi-Center Trial

	Study Design
<b>Objective</b>	Efficacy & safety in anterior circulation stroke in 8 – 24h window
<b>Design</b>	Randomized, double-blind, sham-controlled
<b>Primary endpoint</b>	mRS improvement beyond expectations at 3 months (sliding dichotomy)
<b>Enrollment</b>	18 countries, 73 sites, 1,000 mITT patients, June 2011 – March 2018

34

## Key Inclusion / Exclusion Criteria

Inclusion	Criteria	Exclusion	Criteria
Age	Male 40–80 Female 40–85	Imaging	<ul style="list-style-type: none"> <li>▪ ICH</li> <li>▪ Massive (&gt;2/3)</li> <li>▪ Lacunar</li> <li>▪ Posterior circulation</li> </ul>
NIHSS	7–18		
Time from stroke onset	8–24h		
Clinical & radiological	Anterior circulation	Reperfusion therapy	<ul style="list-style-type: none"> <li>▪ IV thrombolysis</li> <li>▪ EVT</li> </ul>

35

## ImpACT-24B Study Flow

Time Period	Activity
Day 1	<ul style="list-style-type: none"> <li>▪ Baseline imaging</li> <li>▪ 1:1 dynamic randomization</li> <li>▪ Neurostimulator / Sham implantation</li> </ul>
Days 1-5	<ul style="list-style-type: none"> <li>▪ Daily SPG / Sham stimulation</li> </ul>
Day 5	<ul style="list-style-type: none"> <li>▪ Follow-up imaging</li> <li>▪ Implant / Sham removal</li> <li>▪ Day 5 mRS, NIHSS</li> </ul>
Follow-up	<ul style="list-style-type: none"> <li>▪ Day 30, 60 mRS, NIHSS</li> <li>▪ Day 90 mRS (blinded observer), NIHSS, SIS-16</li> </ul>

36


## Blinding in ImpACT-24B Was Extensive and Effective

Study Procedure	SPG	Sham Control
Baseline CT	Brain + PGP canal and fossa	Brain
Patient reference / navigation marker	Y	Y
Local anesthesia	Y	Y
Implantation procedure	Mucosa puncture + implant placement	Mucosa puncture
5 days treatment	Stimulation + vibration	Vibration
Transmitter sticker/positioning	Y	Y
Stimulation adaptation	Comfortable tolerance level (stim)	Comfortable tolerance level (vib)
Day 5 follow-up CT	Brain + electrode position	Brain
Implant removal	Mucosa touching + implant removal	Mucosa touching
Assessments by blinded observer	Y	Y

*James Blinding Index: blinding success in patients and assessors*

37


## Individualized Dose Selection Method Had to Be Modified to Maintain Blinding in ImpACT-24B



1. Start low, escalate in steps
2. Reach minimal sensation / lacrimation level, set as comfortable tolerance level (CTL)  
Treat at this intensity


38

# Individualized Dose Selection Method Had to Be Modified to Maintain Blinding in ImpACT-24B



39

# Primary Endpoint – Disability (mRS) at 90d



Chaisinanunkul, 2015

40

## Primary Endpoint: Success Threshold Based on Patient Age, NIHSS, Brain Side (VISTA Model)

### Success Threshold

	Favorable			Unfavorable		
1	0	1	2	3	4	5/6
2	0	1	2	3	4	5/6
3	0	1	2	3	4	5/6
4	0	1	2	3	4	5/6
5	0	1	2	3	4	5/6

Hosmer Lemeshow goodness of fit: Good model performance

FDA Question: Sliding Dichotomy

AbESTT II, 2008; PAIS, 2009; European Stroke Organization Outcomes Working Group, 2012; Stroke Therapy Academic Industry Roundtable, 2013

41

## Two Pre-Specified Primary Analysis Populations

### 1. Modified intention to treat (mITT)

- All patients in SPG and sham-control groups who received at least 1 full stimulation session

### 2. Confirmed cortical involvement (CCI)

- All mITT patients with
  - NIHSS  $\geq 10$
  - At least one cortical ASPECTS region

FDA Concern: Patient Selection

42

## Multiplicity Controlled Using Hochberg Procedure

- Study succeeds if:
  - $p < 0.05$  in both populations
  - or**
  - $p < 0.025$  in one population

Hochberg, 1988; Multiple Endpoints in Clinical Trials: FDA Draft Guidance, 2017; Lees, 2018

43

## CCI Addition Approved By FDA Prior To Unblinding

<b>2014 &amp; 2016</b>	Interim analyses (sample size / futility) <ul style="list-style-type: none"> <li>▪ No subgroup analyses conducted by DSMB</li> <li>▪ <b>Sponsor and steering committee blinded</b></li> </ul>
<b>Jun 2017</b>	External development– DAWN result <ul style="list-style-type: none"> <li>▪ 1<sup>st</sup> study to show benefit in 24-hour window</li> <li>▪ Salvageable tissue identifiable by NIHSS &amp; imaging</li> <li>▪ Post-hoc evaluation on ImpACT-24A</li> </ul>
<b>Jan 2018</b>	CCI added as co-primary analysis population <ul style="list-style-type: none"> <li>▪ FDA approved – Jun 2018</li> </ul>
<b>Jul 2018</b>	<b>ImpACT-24B data unblinded</b>

FDA Question: CCI Addition

44

## Secondary Efficacy Endpoints

- At Day 90
  - Functional independence (mRS 0 – 2 vs 3 – 6)
  - Able to walk + body self-care (mRS 0 – 3 vs 4 – 6)
  - Stroke-related QoL (Stroke Impact Scale [SIS]-16)
  - Disability-related QoL (utility-weighted mRS, post-hoc)
  
- Long-term (at 180 and 360 days)
  - Patient-reported stroke-impact assessment (RIKS-Stroke)

45

## Demographics

	mITT Population		CCI Population	
	SPG N = 481	Control N = 519	SPG N = 244	Control N = 276
Age; median (years)	70	71	70	72
Female	50%	52%	48%	49%
Pre-stroke mRS > 0	8.5%	5.6%	8.6%	6.2%
Hypertension	87%	84%	87%	85%
Diabetes	24%	27%	22%	24%
Atrial Fibrillation	25%	26%	34%	31%

46

## Baseline Characteristics

	mITT Population		CCI Population	
	SPG N = 481	Control N = 519	SPG N = 244	Control N = 276
<b>NIHSS; median (IQR)</b>	<b>12 (9, 14)</b>	<b>12 (9, 14)</b>	<b>13 (12, 15)</b>	<b>13 (11, 15)</b>
<b>Left Brain Stroke</b>	<b>57%</b>	<b>50%</b>	<b>57%</b>	<b>52%</b>
<b>ASPECTS; median (IQR)</b>	<b>7 (6, 9)</b>	<b>7 (6, 9)</b>	<b>7 (5, 8)</b>	<b>7 (5, 8)</b>
<b>Time from last-known-well to first stimulation; median (hours)</b>	<b>19.9</b>	<b>18.7</b>	<b>19.7</b>	<b>18.5</b>

47


## Primary Efficacy Results

	SPG	Control	Absolute Difference	Odds Ratio (95% CI)	p-value	Threshold
<b>mITT (N = 1000)</b>	<b>48.6%</b>	<b>45.5%</b>	<b>3.2%</b>	<b>1.14 (0.89, 1.46)</b>	<b>0.31</b>	<b>0.05</b>
<b>CCI (N = 520)</b>	<b>49.6%</b>	<b>39.9%</b>	<b>9.7%</b>	<b>1.48 (1.05, 2.10)</b>	<b>0.0258</b>	<b>0.025</b>

48




# Consistent Effect Indicated by All Primary and Secondary Efficacy Endpoints




49

# Benefits Persisted Through 1-Year Follow-up Further Confirming the Positive Effects of SPG Stimulation



50

# SPG Effect Magnitude Consistent with Standard of Care for Ischemic Stroke




FDA Question: Sliding Dichotomy

1. Cranston, 2017; 2. IST1+CAST, 2000, POINT, 2018; 3. ECASS, 2008  
 4. ImpACT-24B, 2019; 5. NINDS-tPA Study, Yafeh, 2007

51


# Strong Dose Response Relationship (CCI)



1. Guyat, J Clin Epidemiol 2011  
 2. FDA Guidance for Industry

52

# Greatest Benefit Seen When Stimulation Delivered Within Optimal Final Dose Range



53

# US vs OUS: No Statistical Differences

Sliding Dichotomy	US	OUS	Interaction p-value
<b>Odds ratio</b> (95% CI)	<b>1.11</b> (0.26, 4.72)	<b>1.50</b> (1.05, 2.15)	<b>0.69</b>

- No evidence of treatment effect difference
  - Non-significant interaction for treatment by geography
- Despite randomization, comparisons of geographic subgroups are sensitive to small imbalances in baseline covariates

FDA Question: Generalizability US

54

## Adjusted Analyses Show Large, Consistent Benefit in US and OUS Patients

Sliding Dichotomy	US	OUS	Interaction p-value
<b>Odds ratio</b> (95% CI)	<b>1.62</b> (0.30, 8.63)	<b>1.46</b> (0.98, 2.18)	<b>0.91</b>

- Differences in baseline characteristics of US subgroup between SPG stimulation and sham-control groups
- Patients in all geographic regions received similar standard of care for ischemic stroke per guidelines, including medications

FDA Question: Generalizability US

55

## No Bias Introduced by mITT Analysis

	CCI Patients Allocated SPG included in Primary Analysis N = 520	Patients Allocated SPG not included in Primary Analysis N = 34
<b>Mean Age, years (SD)</b>	<b>70 (10)</b>	<b>71 (10)</b>
<b>Sex (female)</b>	<b>49%</b>	<b>50%</b>
<b>NIHSS (mean)</b>	<b>13.5 (2.5)</b>	<b>13.8 (2.1)</b>
<b>Stroke side (left brain)</b>	<b>55%</b>	<b>47%</b>
<b>Median ASPECTS (IQR)</b>	<b>7 (5, 8)</b>	<b>7 (5, 7)</b>
<b>Time from stroke onset (hours); median (IQR)</b>	<b>16.3 (13.5, 19.4)</b>	<b>15.4 (13.2, 18.3)</b>

FDA Concern: mITT Analysis

56

# ImpACT-24A – “Phase 2B” RCT

## Consistent Findings Supportive of Pivotal Trial

57

# ImpACT-24A Greatest Benefit in CCI Population

Primary mRS Sliding Dichotomy

mITT	SPG	Control	Odds Ratio (95% CI)	p-value
ImpACT-24A (N = 253)	49.7%	40.0%	1.48 (0.89, 2.47)	0.13

CCI	SPG	Control	Odds Ratio (95% CI)	p-value
ImpACT-24A (N = 87)	50.0%	27.0%	2.70 (1.08, 6.73)	0.03

FDA Question: CCI Addition

58

## ImpACT-24A Consistent with ImpACT-24B

### Primary mRS Sliding Dichotomy

mITT	SPG	Control	Odds Ratio (95% CI)	p-value
ImpACT-24A (N = 253)	49.7%	40.0%	1.48 (0.89, 2.47)	0.13
ImpACT-24B (N = 1000)	48.6%	45.5%	1.14 (0.89, 1.46)	0.31

CCI	SPG	Control	Odds Ratio (95% CI)	p-value
ImpACT-24A (N = 87)	50.0%	27.0%	2.70 (1.08, 6.73)	0.03
ImpACT-24B (N = 520)	49.6%	39.9%	1.48 (1.05, 2.10)	0.0258

59

## Efficacy Conclusions Supported by Pooled Analysis


### Primary mRS Sliding Dichotomy

mITT	SPG	Control	Odds Ratio (95% CI)	p-value
ImpACT-24A (N = 253)	49.7%	40.0%	1.48 (0.89, 2.47)	0.13
ImpACT-24B (N = 1000)	48.6%	45.5%	1.14 (0.89, 1.46)	0.31
<b>Pooled mITT (N = 1253)</b>	<b>48.9%</b>	<b>44.6%</b>	<b>1.20 (0.96, 1.49)</b>	<b>0.12</b>

CCI	SPG	Control	Odds Ratio (95% CI)	p-value
ImpACT-24A (N = 87)	50.0%	27.0%	2.70 (1.08, 6.73)	0.03
ImpACT-24B (N = 520)	49.6%	39.9%	1.48 (1.05, 2.10)	0.0258
<b>Pooled CCI (N = 607)</b>	<b>49.7%</b>	<b>38.3%</b>	<b>1.61 (1.16, 2.23)</b>	<b>0.004</b>


60

## Pooled CCI: Consistent Effect in All Endpoint








61

## Pooled CCI: Homogenous Treatment Effect in All Subgroups



62

# Benefits Achieved Regardless of Core Size and Time from Stroke Onset

Primary Endpoint Sliding Dichotomy mRS	SPG	Control		Odds Ratio (95% CI)
<b>Pooled CCI Population</b>	<b>49.7%</b>	<b>38.3%</b>		<b>1.61 (1.16, 2.23)</b>
<b>ASPECTS ≥ 7</b>	<b>54.3%</b>	<b>44.3%</b>		<b>1.49 (0.94, 2.36)</b>
<b>ASPECTS &lt; 7</b>	<b>45.1%</b>	<b>32.9%</b>		<b>1.68 (1.05, 2.67)</b>
<b>TFSO ≤ 18</b>	<b>49.6%</b>	<b>38.6%</b>		<b>1.56 (0.96, 2.55)</b>
<b>TFSO &gt; 18</b>	<b>49.7%</b>	<b>38.1%</b>		<b>1.61 (1.04, 2.47)</b>

0.2

1

5

Favors SPG

TFSO = Time from Stroke Onset

## Safety Results



## Safety Analysis Populations

### 1. All patients

- SPG stimulation and sham-control patients who had mucosal puncture performed

### 2. CCI

- All safety population patients with
  - NIHSS  $\geq 10$
  - At least one cortical ASPECTS region

65

## Safety Endpoints

### Pre-Specified


- SAEs
- Neurological deterioration ( $\geq 4$  NIHSS within 1st 10 days)
- Mortality by day 90
- Stimulation-related SAEs & AEs
- Implantation-related SAEs & AEs
- Failed implantations

### Additional

- Pneumonia SAEs
- Symptomatic Intracranial Hemorrhage (sICH)

66


# All Patients: No Increase in Risk of Mortality, SAEs or Other Common Stroke Complications



**FDA Question: Safety**

67

# CCI: Consistent Favorable Safety Profile with SPG Stimulation



68

## Stimulation Related SAEs in Both Groups

Possibly Related Events*	All – ImpACT-24B	
	SPG N = 536	Control N = 519
<b>Total stimulation SAE</b>	<b>3 (0.6%)</b>	<b>2 (0.4%)</b>
Stroke in evolution	1 (0.2%)	1 (0.2%)
Hemorrhagic transformation stroke	1 (0.2%)	1 (0.2%)
Epileptic seizure	1 (0.2%)	-

\* No events were classified as “definitely” or “probably” related

69

## Non-Serious Stimulation Related AEs (> 1%)

Related Events	All – ImpACT-24B	
	SPG N = 536	Control N = 519
<b>Lacrimation increased</b>	<b>71 (13.2%)</b>	<b>3 (0.6%)</b>
<b>Headache</b>	<b>19 (3.5%)</b>	<b>4 (0.8%)</b>
<b>Pain</b>	<b>118 (22.0%)</b>	<b>4 (0.8%)</b>
<b>Medical device discomfort</b>	<b>5 (0.9%)</b>	<b>6 (1.2%)</b>

- All transient
- Lacrimation: known sign of SPG activation, resolves at end of treatment session
- Headache may be side effect of SPG activation
- Facial pain & discomfort: avoidable by not exceeding CTL

70

## Implantation Related SAEs with ISS500; All Resolved Without Sequela

Preferred Term	ImpACT-24B	
	Current Implant N = 197	First Implant N = 339
<b>Total implant related SAE</b>	<b>1 (0.5%)</b>	<b>2 (0.6%)</b>
Complication of device removal	1 (0.5%)	1 (0.3%)
Device breakage	-	1 (0.3%)

71

## Implantation is Safe

	ImpACT-24B	
	Current Implant N = 197	First Implant N = 339
<b>Skin-to-skin (minutes); median (IQR)</b>	<b>17 (12, 23)</b>	<b>35 (25, 52)</b>
<b>SAE</b>	<b>0.5%</b>	<b>0.6%</b>
<b>AE</b>	<b>7.6%</b>	<b>36.9%</b>
<b>Misplacements</b>	<b>2.0%</b>	<b>8.3%</b>
<b>Incomplete procedures</b>	<b>2.0%</b>	<b>5.0%</b>

72

# Modified Implant Mitigated Implantations Risks


Implantation-Related Non-Serious AEs – FDA’s List of Concerns

	ImpACT-24B	
	Current Implant N = 197	First Implant N = 339
Acute pain	1% (2)	11.5% (39)
Bleeding (implant site hemorrhage)	0	3.8% (13)
Swelling (including Infection, Erythema)	0	1.5% (5)
Chronic neuropathic pain / nerve injury	0	1.5% (5)
<b>Micro-aspiration</b>		
Pneumonia Aspiration	0	0.6% (2)
Bronchopneumonia	0	0.3% (1)
Apnea	0	0.3% (1)
Airway endangerment / Laryngospasm	0	0
Palate laceration	0	0

FDA Question: Safety


73

# All Patients Pooled: Safety Results Further Support Favorable Safety Profile



74

# CCI Pooled: Consistent Safety Results








FDA Question: Safety

75

## Overall Benefit-Risk Assessment

76


# Efficacy Summary



# Safety Conclusion

**Implantation is Safe** CO-72

	ImpACT-24B	
	Current Implant (N=197)	First Implant (N=339)
Skin-to-skin (minutes), median (IQR)	17 (12, 23)	35 (25, 52)
SAE	0.5%	0.6%
AE	7.6%	36.9%
Misplacements	2.0%	8.3%
Incomplete procedures	2.0%	5.0%



## ImpACT-24M Trial

### Usability

79


## ImpACT-24M Objectives

1. Speed and accuracy of implantation with final device
2. Confirm dose of SPG stimulation can be efficiently set by non-noxious physiological effects

80



# Implantation in ImpACT-24M Accurate and Efficient



FDA Question: Implantation

81


# Implantation is Safe and Simple (< 5 minutes)

	ImpACT-24M	ImpACT-24B	
	Final Device N = 50	Current Implant Old Navigation N = 197	First Implant Old Navigation N = 339
Skin-to-skin (minutes); median (IQR)	4 (3, 7)	17 (12, 23)	35 (25, 52)
SAE	0% (0)	0.5% (1)	0.6% (2)
AE	2% (1)	7.6% (15)	36.9% (125)
Misplacements	0% (0)	2.0% (4)	8.3% (28)
Incomplete procedures	0% (0)	2.0% (4)	5.0% (17)

FDA Question: Safety


82

# ImpACT-24M Stimulation Level Set Based on Non-Noxious Physiologic Signs of SPG Activation




83

# 92% of Patients Treated Within Final Dose Range Using CTL Identification



84

# Correct Stimulation Validated by Increased Common Carotid Artery Flow Volume





Ultrasound doppler image of common carotid artery

	Baseline	During Stimulation	Increase	p-value
Diameter (mm)	8.0	8.9	11%	<0.0001
Peak systolic velocity (cm/s)	65.6	76.8	17%	0.0001
Peak systolic flow (cc/s)	32.5	46.9	44%	<0.0001
End diastolic velocity (cm/s)	14.0	17.1	22%	0.0004
End diastolic flow (cc/s)	7.1	10.8	52%	<0.0001

Parameters assessed between 45-60 minutes after initiation of stimulation

85

# Improved Motor Deficit Confirms Correct Stimulation Setting



86

## ImpACT-24M Supports Applicability to Real-World Clinical Practice

- Current device facilitates accurate and simple implantation
  - Shortened procedure time
  - 100% correct placements
  - No SAEs
- Setting stimulation intensity using non-noxious physiologic markers yields intensity levels in optimal range
- Immediate effects on blood flow and motor function

87



### Training and Post-Approval Plan

**Eyal Shai**

Chief Technology Officer  
BrainsGate

88

## Implanter Training Program

- Implantation performed by medical doctor and an assistant
- One day in-person training
  - Didactic session (2 hours)
  - Practice on head model (5 hours)
- 5 Implantations, remote guidance & support
- Qualification - 3 procedures with remote supervision



FDA Question: Implantation

89

## Treatment Training Program

- Online training for healthcare professionals
  - Physiological markers identification
  - System operation
  - Online test

FDA Question: Implantation

90

## Post Approval Study Plans

1. System's performance (Automatic data collection)
  - Guidance system accuracy
  - Implantation procedure time
  - Stimulation level
2. Registry Data collection
  - Clinical outcome
  - Failed implantations
  - Complications (device-related), safety incidents
  - FDA identified theoretical risks

FDA Question: Safety

91




### Clinical Perspective

**Michael Hill, MD, MSc, FRCPC**

President, Canadian Neurological Sciences Federation  
Professor, Department Clinical Neuroscience &  
Hotchkiss Brain Institute  
Cumming School of Medicine, University of Calgary &  
Foothills Medical Centre

92



# Unmet Need for Patients Who are Ineligible or Unable to Receive Current SoC



No effective alternatives to reduce disabling effects of acute ischemic stroke

Bahr Hosseini, 2018; Desai 2020







# > 75% of Patients in Treatment Gap have Worse Outcomes with mRS > 3 at Discharge



Bahr Hosseini, 2017; Desai 2020

LVO = Large Vessel Occlusion; MVO = Medium Vessel Occlusion; SVO = Small Vessel Occlusion

# Clinical Situations Where ISS500 Could Provide Benefit for Patients Presenting within 24 Hours but Ineligible for SOC

<p><b>Distal MVOs</b> (e.g. M2/M3/M4/A1/A2/A3)</p> 	<p><b>Large Core But Salvageable Penumbra Present</b></p> 	<p><b>Thrombectomy-Capable Centers Too Far Away</b></p> 
<p><b>Chronic Cervical ICA Occlusion+ Acute MCA/ACA Thrombus</b></p> 	<p><b>Chronic Cervical or Intracranial Occlusion + Hemodynamic Failure</b></p> 	<p><b>Tortuous Aortic / ICA Anatomy</b></p> 

Menon, 2018; Grossberg, 2018; Chen 2018; Momjian-Mayor, 2005; Saba, 2015; Bahr Hosseini, 2018; Adeoye, 2014






# Identification of CCI Patients Based on Practices Routinely Used Today

- CCI determined by neurological examination using assessments routinely done in clinical practice
  - Combination of NIHSS and ASPECT score
  - Non-contrast CT and other standard imaging modalities


FDA Concern: Patient Selection



## Data Demonstrate that SPG Stimulation is Effective Treatment Option: Addresses Treatment Gap

Primary Endpoint Sliding Dichotomy mRS	SPG	Control		Odds Ratio (95% CI)
<b>Pooled CCI Population</b>	<b>49.7%</b>	<b>38.3%</b>		<b>1.61 (1.16, 2.23)</b>
<b>ASPECTS ≥ 7</b>	<b>54.3%</b>	<b>44.3%</b>		<b>1.49 (0.94, 2.36)</b>
<b>ASPECTS &lt; 7</b>	<b>45.1%</b>	<b>32.9%</b>		<b>1.68 (1.05, 2.67)</b>
<b>TFSO ≤ 18</b>	<b>49.6%</b>	<b>38.6%</b>		<b>1.56 (0.96, 2.55)</b>
<b>TFSO &gt; 18</b>	<b>49.7%</b>	<b>38.1%</b>		<b>1.61 (1.04, 2.47)</b>

0.2                      1                      5

 **Favors SPG**

TFSO = Time from Stroke Onset

97

## Safety Profile Well-Characterized, Demonstrating Minimal Risks Associated with SPG Stimulation

- SPG stimulation associated with mild, transient stimulation AEs
  - Pain, lacrimation and headache
  - Resolved during stimulation, or upon completion of therapy
- No increased risk of mortality
- SAEs, neurological deterioration and pneumonia less common with SPG compared with Control
- SPG significantly reduces rate of symptomatic ICH


98

# Simplifications to Implant Procedure Designed to Ensure Safe Transition to Real-World Use

- Simple bedside procedure
- Navigation system facilitates safe and accurate positioning of the neurostimulator
  - Minimizes risk of implant complications
- Final device ensures stimulation within optimal range

FDA Question: Implantation

# Effectiveness and Safety Results Represent Clinical Meaningful Outcomes for Patients



# Totality of Evidence Supports Positive Benefit-Risk for Ischemic Stroke System

## Unmet Need

- Guidelines recommend reperfusion therapies
- Use is time dependent
- Many patient's ineligible or do not have access to care

## Effectiveness

- Target CCI population achieved consistent improvements
- Benefits regardless of stroke severity and time from stroke onset
- Final device ensures optimal stimulation

## Safety

- Favorable safety profile
- Significantly reduced risk of sICH
- Final device reduces procedure time and implant complications

# Ischemic Stroke System for Treatment of Patients with Acute Ischemic Strokes

**December 10, 2021**

Neurologic Devices Panel

BrainsGate