


September 7th, 2022
Rahul Thakar
Head of Technical & Scientific Management

Characteristics of Pharmaceutical Elastomers in Container Closure Systems

Agenda

- **Introduction to Container Closure Systems**
 - Vials, PFS, Cartridges with a focus on elastomeric components
- **Introduction to Pharmaceutical Elastomers**
 - Selection criteria and key considerations
 - Physical and chemical properties
 - Applications, variations, and functionalities
 - E/L profiles
 - Common Concerns & Industry Trends
- **Pharmaceutical Elastomer Manufacturing Process**
- **Processing of Elastomeric Components**
 - Fundamentals of RFS and RTU components
 - **Importance of siliconization** and selection criteria
 - Basics of camera inspection
 - Sterilization choices and elastomer packaging selection
- **Importance of Manufacturing environment**

Pharmaceutical Elastomers & Seals



Pharmaceutical Elastomers

Customers dealing with different type of primary containers:

- Stoppers and aluminum seals **for Vial applications**
- Plungers and Combiseals **for Cartridge applications**
- Plungers, Tip Caps and Needle Shields **for Prefilled Syringe applications**
- Custom designs for **auto-injectors, drug delivery devices**



Sealing solutions for Vials



Small volume parenteral (SVP)

Large volume parenteral (LVP)

VIAL STOPPERS

	SERUM		IGLOO		LYO	
		OmniFlex™		OmniFlex™	2-LEG	OmniFlex™
13 mm						
Coated	-	Yes	-	Yes	-	Yes
Design	V9024	V9401	V9250	V9402	HPP002	HPP079
Compound	FM140 FM257 FM457	FM457	FM140 FM257 FM457 FM460	FM457	FM140 FM257 FM457 FM460	FM457

	SERUM		IGLOO		LYO	
		OmniFlex™		OmniFlex™	2-LEG	OmniFlex™
20 mm						
Coated	-	Yes	-	Yes	-	Yes
Design	V9025	V9407	V9172	V9397	V9154	V9396
Compound	FM140 FM257 FM457	FM457	FM140 FM257 FM457 FM460	FM457	FM140 FM257 FM457 FM460	FM457

	INFUSION	INFUSION		LYO	
		DIN	ISO	OmniFlex™	4-LEG
29 mm					
32 mm					
Coated	-	-	-	Yes	-
Design	5294	V9003	V9240	V9240	V9262
Compound	FM140 FM257 FM457	FM140 FM257 FM457	FM140 FM257 FM457	FM259	FM140 FM257 FM457 FM460

Sealing solutions for Vials


Key take-away messages:

Stopper types – serum, lyophilization, infusion


Typical sizes – 13, 20, 28, 29, 32 mm or custom designs

Elastomer formulations – typically halobutyls


Blowback – European or American or no blowback



No Blowback



European Blowback



American Blowback

Image courtesy: Schott

Sealing solutions for Vials



Aluminum seals

Flip Caps
(plastic caps and aluminum seals)

STANDARD CAPS	13 mm	20 mm	32 mm
ALUMINIUM HOLE CAP	 AHC13.0005	 AHC20.0004	 AHC32.0009
CENTER TEAR OFF	 CTC13.0006	 CTC20.0010	 CTC32.0009
SCORELINE TEAR OFF		 STC20.0004	 STC32.0006
DOUBLE TEAR OFF		 DTC20.0001	 DTC32.0002
FLIP CAP <small>Also available in flush design</small>	 FC13.0009	 FC20.0016	 FC32.0011
SCORELINE FLIP CAP <small>Also available in flush design</small>	 SFC13.0008	 SFC20.0020	 SFC32.0001
UNIVERSAL TEAR OFF	Flip-Tear UTC20.0007 Tear-Off UTC20.0002	 UTC32.0016	
PULL OFF CAP			 POC32.0001

PrimeCap™


13 mm



20 mm



Sealing solutions for Cartridges



From 1–3 ml

DuraCoat™



Type	BI-LAYER	MONO-LAYER	BI-LAYER	MONO-LAYER	BI-LAYER
Design	ACS0001	ACS0002	ACS0003	ACS0004	ACS0006
Compound	FM257 + H1-7-207 FM457 + H1-7-207	FM257 FM457	FM257 + H1-7-207 FM457 + H1-7-207	FM257 FM457	H5-23-5 + H1-7-207 FM257 + H1-7-207 FM457 + H1-7-207

ø 6.85 mm


ø 8.65 mm

ø 9.65–9.7 mm



Coated	–	–	–	–	Yes
Design	V4528	V9321	V9336	V9478	V9520
Compound	FM257 FM457	FM257 FM457	FM257 FM457	FM257 FM457	FM457

Sealing solutions for Prefilled Syringes



SOFT NEEDLE SHIELD

½"



Coated

-

Design

V6544

Compound

FM27
FM30

RIGID NEEDLE SHIELD

½"



-

V9812

FM27
FM30

TIP CAPS

RIBBED

MUSHROOM



-

V9406

FM27
FM30
FM257

-

V9257

FM27
FM30
FM257

PLUNGERS

0.5 ml

1 ml long

1-3 ml

5 ml

10 ml



Coated

-

Yes

-

Yes

-

Yes

-

-

Design

V9315

V9503

V9283

V9519

V9258

V9517

V9319

V9344

Compound

FM257
FM457

FM457

FM257
FM457

FM257
FM457

FM257
FM457

FM457


FM257
FM457

FM257
FM457

Introduction to Pharmaceutical Elastomers


Elastomers in Pharmaceutical Applications

Vials



- Primary & secondary packaging components for small and large volume parenterals
- Fluoropolymer coated option

Prefilled syringes & cartridges



- High-precision components for injection systems and lined seals
- Fluoropolymer coated option

Elastomers in Pharmaceutical Applications

Blood Collection



- Sealing components and needle sleeves for blood collection applications
- Tailor-made designs

Single-use Syringe



- Sealing components, not made with natural rubber, for single-use syringes

IV Systems



- Various components, not made with natural rubber, for IV systems

Key Considerations for Pharmaceutical Elastomers

- Is the API or formulation absorbed by the elastomer?
- API and formulation stability
- Does the rubber react with the API and leach out impurities? (E/L profile)
- Temperature range at which closure and product are stable?
- Long term effects of storage
- Effects of sterilization on closures
- Seal integrity over the course of shelf life; CCI..?



Common Concerns

Drug Sensitivity

Drug products (especially biotech) are more sensitive to yesterday's packaging solutions:

- Extractables and leachables
- Silicone fluid
- Particles



Common Concerns

Component Functionality


- Stoppers
 - Coring/ fragmentation
 - Needle penetration forces
 - Needle self-sealing capacity
 - Machinability challenges
 - CCI
- Alu Seals
 - Pull-off forces
- Plungers in PFS systems
 - Break loose/ gliding behavior
 - CCI

USP <381> guidance is available for component functionality




Revision of USP <381>

USP <381> Elastomeric Closures for Injections



- USP <381> Elastomeric Components in Injectable Pharmaceutical Product Packaging/Delivery Systems
- USP <1381> Assessment of Elastomeric Components Used in Injectable Pharmaceutical Product Packaging/Delivery Systems
- USP <382> Elastomeric Component Functional Suitability in Parenteral Product Packaging/Delivery Systems
- USP <1382> Assessment of Elastomeric Component Functional Suitability in Parenteral Product Packaging/Delivery Systems

Evolution of Pharmaceutical Elastomers



Key Considerations for Pharmaceutical Elastomers

Relevant Physical Properties:

- Elastomer hardness (ISO 48-4)
- Density (ISO 2781)
- Ash Content (Ph.Eur. 2.4.16)
- Compression Set (ISO 815-1)
- Tensile Strength, Elongation, Modulus (ISO 37)
- Water vapor transmission rate (ASTM F-1249)
- Oxygen transmission rate (ASTM D-3985)



Relevant Chemical Properties and Regulatory Compliance

- Ph. Eur. 3.2.9, ISO 8871-1, USP <381> <382>, JP 7.03
- E/L profile – tested in WFI (basic, neutral, acidic), hexane, isopropyl alcohol

Important to have data packages with physical, chemical and regulatory information for pharmaceutical elastomers

**Aforementioned testing is not a comprehensive list of tests that are needed or performed.*

Pharmaceutical Elastomers

Pharmaceutical elastomers are thermoset rubbers as opposed to thermoplastic materials. Examples of thermoset rubbers:

- Natural rubber (NR)
- Polyisoprene (IR)
- Polychloroprene (CR)
- Styrene butadiene (SBR)
- Nitrile butadiene rubber (NBR)
- Ethylene propylene diene monomer (EPDM)
- Butyl rubber (IIR)
- Silicone rubber (Q)

Pharmaceutical Elastomers

Natural rubber / Polyisoprene

- Natural rubber: latex allergy discussions
- Historically the oldest elastomer type
- Need complex curing systems
- Good elastic properties
- Polyisoprene (synthetic) replaces Natural rubber

SBR (styrene-butadiene rubber)

- Intermediate permeability
- Typically used for pre-assembled EtO sterilized components (e.g. Needle Shields)

Pharmaceutical Elastomers, cont.


Halobutyl (BromoButyl, ChloroButyl)

- Cleanest curing system
- Lowest permeability
- High resistance to ageing
- Regular **butyl** still on the market, and also newer types like **BIMS** (Brominated isobutylene para-methylstyrene)

Pharmaceutical Elastomers, cont.

Nitrile rubber


Typically used for mineral oil based drugs



Silicone rubber

High permeability


Typically not used for parenteral applications



EPDM rubber


For niche applications

Oxygen Transmission Rate




- Data shown for typical SBR (styrene-butadiene) and isoprene compounds
- FM 259 is a fluoropolymer coated compound

Water Vapor Transmission Rate



- Data shown for typical SBR (styrene-butadiene) and isoprene compounds
- FM 259 is a fluoropolymer coated compound

Elastomer Formulations



Elastomer Formulations

Typical ingredients in elastomers (a summary):

- **Elastomer base** – Halobutyl, SBR, Isoprene
- **Fillers** – Clays, mineral powders; provide structure
- **Plasticizers** – processing aids, provide flexibility
- **Curing System**
 - Crosslinking agent – introduce elasticity; co-determine strength and hardness
 - Accelerator – no longer used in modern formulations due to increased E/L
 - Activators – activate crosslinking sites, increase vulcanization efficiency
- **Pigment** – Impart color
- **Antioxidant** – prevent rubber degradation
- **Wax** – processing aids



Elastomer Formulations

Fillers give mechanical strength (stiffness) to a rubber

Attributes physical properties to a rubber compound

- More filler = Harder compound
 - Better for gliding profile plungers
 - Better against stickiness in bulk
 - Worse for stopper piercing (coring!)


Inorganic fillers ('white compounds')

- Aluminum silicate (clay)
- Magnesium silicate (talc)
- Silicate
- [Calcium carbonate]

Carbon black ('black compounds')

- Undesired for cleanliness reasons
- May be associated with polynuclear aromatic hydrocarbons (PAHs)

Elastomer Formulations



- **Cure systems:**
 - Crosslinking agents
 - Activator: gives the onset of vulcanization
 - Accelerator: speeds up the vulcanization, easily extractable organic molecules such as thiurams, sulfenamides, thiazoles, ...
- **Modern cure systems**
 - Aim at giving less extractables
- **Historic cure systems**
 - Used easily extractable organic accelerators

Elastomer Formulations

Inorganic pigments

- Titanium dioxide
- Traces of carbon black
- Oxides of iron

Organic pigments

- Avoided in modern compounds

Halobutyl polymer stabilizers

- Calcium stearate
- Epoxydized soybean oil

Anti-oxidants

- Already present in halobutyl elastomer
- Hindered phenol type anti-oxidants
- Additionally added to improve environmental stability (ageing)

Plasticizer, Waxes, Oils


- High polymeric weight plasticizers, Paraffinic oil
- To tune a formulation (e.g. reduce coring)

Extractables & Leachables

To understand the difference between leachables and extractables

- **Extractables** = The entire set of chemical species which “could” leach from the packaging over a wide pH, temperature, and solvent range.
 - **Leachables** = Real time/game ready. Tested at the end of shelf life under normal conditions.
-
- Request available data packages from vendors which detail broad based analytical testing of packaging materials
 - Understand toxicity and immunogenicity issues with identified leachables
 - There are so many pathways to destabilize your drug product; if it’s not in there it can’t cause problems.


Extractables & Leachables



Extractables & Leachables

Difference in Extractable Results for an **OLD** vs **NEW** rubber

(IPA Extract; GC/MS analysis)



“OLD” Pharmaceutical RUBBER

“NEW” Pharmaceutical RUBBER


Common Concerns

Industry Trends

- Efficient filling/stoppering operations
- Less human intervention in the fill/finish process
- Need for zero defects
- Counterfeit products
- Smaller batch sizes/individualized medicine
- Lyophilization
- Cold chain storage issues
- Novel device requirements





Fluoropolymer Spray Coating Technology



Fluoropolymer Spray Coating Properties


Ultra low extractable: spray coating acts as a barrier to a wide range of chemical species



Sensitive formulations, large molecules, biologics and (even with small molecules)

Elastomer Manufacturing and Component Processing


Elastomer Manufacturing Process



Component Processing

- Fundamentals of RfS and RTU components
- Importance of siliconization and selection criteria
- Basics of camera inspection
- Sterilization choices and elastomer packaging selection

Product Options



Ready for Sterilization (RfS)

Why are rubber components washed?

- Washing removes elastomeric debris, dirt and particles generated during molding and die-trimming
- Bring products in a controlled state of microbiological cleanliness (bioburden and endotoxins)
- Meet regulatory requirements:
 - FDA's guidance for industry – *Sterile Drug Products Produced by Aseptic Processing — Current Good Manufacturing Practice cGMP of 2004*
 - EMA guidance: *Note on quality of water for pharmaceutical use*

Why are components siliconized?

- Elastomers are inherently sticky
- To allow for functional properties, machinability and transportation

Ready for Sterilization (RfS)


What makes components RfS?

- Washed with a validated washing program
 - Purified water, and with Water-for-Injection (WFI) in last rinse
 - Defined specifications on low level of microbiological contamination
- Microbiological specifications apply:
 - for FirstLine™ → Bioburden: ≤ 0.05 CFU/cm²; Endotoxins: ≤ 0.02 EU/cm²
- Packaged in steam sterilizable (RfS) bags
 - Flexible packaging configurations including multiple Tyvek bags, RTP bags
- Visual inspection results with particulate analysis and silicone monitoring per batch



Ready for Sterilization (RfS)

'Easy-to-wash'



'Difficult-to-wash'




Product Design: Impact of Silicone

Lubricious barrier coatings:

Siliconization of elastomeric components


#particles per 10 cm² (ISO 8871-3)




Datwyler provides both silicone options and in varied quantities to match the drug packaging needs.

Product Design: Impact of Silicone

Lubricious barrier coatings:
Ultra-low visible and subvisible particle levels



Product Design: Impact of Silicone



- High viscosity (30,000 cSt) and low viscosity (350 cSt)
- High viscosity silicone produces ~7 times less silicone particulate than low viscosity silicone
- Linear relationship between the quantity of silicone



Introduction to Camera Inspection

Chance of Detection

Perspective



Contrast of the defect



Camera Inspection


REJECTED DEFECTS:

Capabilities: What can we see and remove?


- Particulates (fiber, hair, dust, etc.)
- Discoloration -- Embedded metal
- Dimensional errors
- Incomplete features
- Damages or defects

What are the limits?


- System will capture and remove particulates of area size $<0.03\text{mm}^2$ or $\sim 175 \times 175$ microns in area (threshold)
- Minimum 10% contrast ratio difference in gray values (product color is important)




Deformed product



Foreign Matter



Damaged Coating



Introduction to Sterilization Methods

Gamma Sterilization

- An irradiation method where electromagnetic gamma rays are used to sterilize material
- Due to higher penetrability of gamma rays (relative to beta), frequently used to sterilize dense materials in the bulk
- Typical dosage is from 15-40 KGy, and is material dependent

Steam Sterilization

- A physical sterilization method where material to be sterilized must be in contact with steam at 121°C, 30 min.
- Typically, material is packaged in Tyvek bags to allow steam permeation

Ethylene Oxide (ETO) Sterilization

- A chemical sterilization method that is widely used for disposable medical devices and for sterilization of nested and tubed empty prefilled syringes with assembled tip cap or needle shield

An Assessment – Gamma-RTU and Steam-RTU

Gamma-RTU

- **Flexible packaging conditions** – due to high penetrability of gamma rays, dense materials packaged in bulk and stacked can be sterilized
- **Low cost solution that is easy to implement** – No isolators or clean rooms required for sterilization

Advantages

- **No measurable effect on chemical/physical properties**
- **High compatibility** with elastomeric formulations

Steam-RTU

- **Measurable effect on chemical/physical properties that is dependent on the elastomer formulation and irradiation dosage.** Elastomeric formulations need to be formulated that are compatible with gamma. – For eg. hardness, stickiness, fragmentation may be altered

Disadvantages

- **Residual moisture content monitoring** in elastomer may not be desirable for certain applications – For eg. Lyophilization
- **Exposure to high temperature dry heat** – Subsequent drying step (from dry heat) may adversely affect elastomer or the outside packaging material
- **Tyvek bags** are a source of particulate contamination

Importance of Manufacturing Environment

Pharmaceutical Trends



Molecular Mass

**Small molecule
Low sensitivity**

**Large molecule
High sensitivity**

Should the component manufacturing environment be an extension of the API manufacturing environment?

Thank you!