

OPTIMIZED COLLOIDAL GOLD ASSAY



C-SYNC™

COVID-19 ANTIGEN TEST

FAST | EASY | RELIABLE | RESULTS IN 10 MINUTES

For use under the Emergency Use Authorization (EUA) only.

For *in vitro* diagnostic use.

For prescription use only.

Instructions for Use

Read the complete IFU prior to performing a test



COV9712A

INTENDED USE

The C-Sync™ COVID-19 Antigen Test is a lateral flow immunoassay intended for the qualitative detection of nucleocapsid protein and spike protein antigens from SARS-CoV-2 in direct anterior nasal swab specimens collected from individuals who are suspected of COVID-19 by their healthcare provider within the first seven (7) days of symptoms onset, when tested at least twice over three days with at least 48 hours between tests, and from individuals without symptoms or other epidemiological reasons to suspect COVID-19, when tested at least three times over five days with at least 48 hours between tests.

Testing is limited to laboratories certified under the Clinical Laboratory Improvement Amendments of 1988 (CLIA), 42 U.S.C. §263a, that meet the requirements to perform moderate, high or waived complexity tests. This product is authorized for use at the Point of Care (POC), i.e., in patient care settings operating under a CLIA Certificate of Waiver, Certificate of Compliance, or Certificate of Accreditation.

The C-Sync™ COVID-19 Antigen Test does not differentiate between SARS-CoV or SARS-CoV-2 viruses.

Results are for the identification of SARS-CoV-2 nucleocapsid and spike protein antigens, which are generally detectable in anterior nasal swab specimens during the acute phase of infection. Positive results indicate the presence of viral antigens, but clinical correlation with patient history and other diagnostic information is necessary to fully determine infection status. Positive results do not rule out bacterial infection or coinfection with other viruses. The agent detected may not be the definitive cause of disease.

Laboratories within the United States and its territories are required to report all results to the appropriate public health authorities.

All negative results are presumptive and confirmation with a molecular assay if necessary for patient management, may be performed. Negative results do not rule out SARS-CoV-2 infection and should not be used as the sole basis for treatment or patient management decisions, including infection control measures such as isolating from others and wearing masks. Negative results should be considered in the context of an individual's recent exposures, history, and the presence of clinical signs and symptoms consistent with COVID-19.

The C-Sync™ COVID-19 Antigen Test is intended for use by healthcare professionals or operators who are proficient in performing tests in point of care settings. The C-Sync™ COVID-19 Antigen Test is only for use under the Food and Drug Administration's Emergency Use Authorization. This product has not been FDA cleared or approved.

SUMMARY AND EXPLANATION OF THE TEST

SARS-CoV-2 is a highly contagious virus which emerged at the end of 2019 and is spread primarily through respiratory particles.

The C-Sync™ COVID-19 Antigen Test is a visually read rapid lateral flow immunoassay for the qualitative detection of SARS-CoV-2 directly from nasal swabs, without Viral Transport Media (VTM). Results are available in less than 10 minutes, making it a valuable evaluation tool for use in Point of Care.

PRINCIPLES OF THE TEST

The C-Sync™ COVID-19 Antigen Test is a lateral flow immunochromatographic assay for the detection of nucleocapsid and/or spike protein antigens specific to SARS-CoV-2 in nasal swab specimens, directly collected from individuals with or without symptoms, as described in the intended use.

SARS-CoV-2 specific antibodies and a control protein are immobilized onto a membrane support as two distinct lines and combined with other reagents/pads to construct a test strip encased in a plastic cassette with a sample well.

An anterior nasal swab sample collected from the patient is eluted into a vial with an extraction buffer and four drops are added to the sample well. The sample flows along the membrane. Test results are interpreted at 10 minutes. A red/purple Test (T) line and a Control (C) line appearing on the test strip, indicates that SARS-CoV-2 antigen was detected. The presence of one colored line at the C-line indicates that SARS-CoV-2 antigen was not detected. No appearance of a colored line at the C-line indicates an invalid test.

MATERIALS AND SUPPLIES PROVIDED

The C-Sync™ COVID-19 Antigen test kits are available with 20 tests per kit including the following components:

- 20 Test cassettes individually packaged in a foil pouch with a desiccant
- 20 Extraction buffer tubes with dropper cap
- 20 Sterile Nasal Swabs for sample collection
- Instructions for Use (IFU)
- Quick Reference Instructions

MATERIALS REQUIRED BUT NOT PROVIDED

- External Positive and Negative Quality Controls sold separately - BioSynchronicity Product Catalog Number: COV9712ACS available at www.biosynchronicity.com
- Personal Protective Equipment per local recommendations
- Timer
- Tube rack for specimens

WARNINGS AND PRECAUTIONS

- Read all instructions carefully before performing the test. Failure to follow the instructions may result in inaccurate test results.
- For prescription and *in vitro* diagnostic use only.
- In the USA, the emergency use of this product is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of *in vitro* diagnostic tests for detection and/or diagnosis of COVID-19 under Section 564(b)(1) of the Act, 21 U.S.C. § 360bbb-3(b)(1), unless the declaration is terminated or authorization is revoked sooner.


TEST PROCEDURE

A. TEST PREPARATION

1. Allow the kit components to reach a temperature between 15°C - 30°C (59°F - 86°F) prior to testing.
2. Remove a C-Sync™ test from a pouch. Place the C-Sync™ Test Cassette on a flat and clean surface, for the entire duration of test operation.

B. SPECIMEN COLLECTION PROCEDURE

Remove a nasal swab from the package.



Insert the entire collection tip of the swab provided (usually ½ to ¾ of an inch, or 1 to 1.5 cm) inside the nostril.



Gently but firmly, collect a sample from the nasal wall by slowly rotating the swab into a circular path against the nasal wall at least 4 times for a total of 15 seconds. Be sure to collect any drainage that may be present on the swab. Slowly remove the swab.



Using the same swab, repeat the same sample collection process in the other nostril.



If immediate testing is not possible, and to maintain the best performance and avoid possible contamination, it is highly recommended the nasal swab specimen be placed in a clean, unused plastic tube and capped tightly at room temperature (15°C-30°C) for up to (1) hour prior to testing. Ensure the swab fits securely within the tube and the cap is tightly closed. **If greater than 1 hour delay occurs, dispose of the sample. A new sample must be collected for testing.**

C. SPECIMEN PREPARATION PROCEDURE

Remove the foil seal from the extraction buffer tube.



Carefully insert the swab inside the tube and plunge it to the bottom. Rotate the swab against the wall of the tube for 10 seconds.



Tightly squeeze and release the bottom of the tube against the absorbent part of the swab. This will make the liquid buffer inside the tube rise and fall. Make sure the buffer rises and falls to fully cover and wash the absorbent part of the swab. Twirl the swab while squeezing and washing the swab with the buffer. Repeat this squeezing and twirling technique 5 times.



Stop squeezing the bottom of the tube and rotate the swab against the wall of the tube for 10 seconds.



Remove the swab slowly while squeezing tightly the top part of the tube against the absorbent part of the swab to extract as much liquid as possible from the swab with the goal of drying the swab. Firmly place the dropper cap attached to the vial onto the top of the tube. Make sure the dropper cap is securely in place.



D. SPECIMEN TESTING

Set the timer at 10 minutes, but do not start yet. Invert the extraction tube and squeeze the ridge area of the tube to dispense 4 drops, drop by drop, into the sample well of the C-Sync™ Test Cassette. Tap the device if the specimen fluid appears viscous when dropping the buffer into the sample well to aid in the wicking process.



Apply only 4 drops to the sample well. Erroneous results may occur if less/more than 4 drops are applied.




Start the timer for 10 minutes. A visual observation of the red/purple dye entering the evaluation window of the test within 10 to 15 seconds after the test timer has been started shows the test is progressing. The test must be read within 10 minutes, as indicated below in Reading and Interpretation of Results.

Do not read test results before 10 minutes or after 15 minutes. Erroneous results may occur when reading before 10 minutes or after 15 minutes.


E. READING AND INTERPRETATION OF RESULTS

After dispensing the drops in the sample well, the results will be available as follows:

Positive: If the Control line (C-line) is red/purple within 10 minutes and the Test line (T-line) is visible within 10 minutes, the test is positive. The color intensity of the Test line (T-line) will depend on the amount of SARS-CoV-2 nucleocapsid and/or spike protein antigens in the sample. Any very faint colored Test line (T-line) should be considered positive. **Tests should not be interpreted past 15 minutes**



Negative: If the Control line (C-line) is red/purple at 10 minutes, and the Test line (T-line) is not visible at 10 minutes, the test is negative.



Invalid: If the Control line (C-line) is not visible at 10 minutes, and regardless if the Test line (T-line) is visible or not visible, the test is invalid. Discard the test. A new test should be repeated.



Do not interpret the results before 10 minutes. Results read after 15 minutes may be invalid.

Repeat testing is needed to improve test accuracy. Please follow the table below when interpreting test results for COVID-19.

Status on First Day of Testing	First Result Day 1	Second Result Day 3	Third Result Day 5	Interpretation
With Symptoms	Positive	N/A	N/A	Positive for COVID-19
	Negative	Positive	N/A	Positive for COVID-19
	Negative	Negative	N/A	Negative for COVID-19
Without Symptoms	Positive	N/A	N/A	Positive for COVID-19
	Negative	Positive	N/A	Positive for COVID-19
	Negative	Negative	Positive	Positive for COVID-19
	Negative	Negative	Negative	Negative for COVID-19

Results should be considered in the context of an individual's recent exposures, history, and the presence of clinical signs and symptoms consistent with COVID-19.

POSITIVE RESULT

The presence of both the control line (C) and test line (T) indicates the presence of SARS-CoV-2 antigen. The result suggests current SARS-CoV-2 infection. Samples with low levels of antigen may produce a faint test line. A visible test line is considered positive.

Repeat testing does not need to be performed if patients have a positive result at any time.

A positive test result means that the virus that causes COVID-19 was detected in the sample, and it is very likely the individual has COVID-19 and is contagious. Please contact the patient's doctor/primary care physician (if applicable) and the local health authority immediately and instruct your patient to adhere to the local guidelines regarding self-isolation. There is a very small chance that this test can give a negative result that is incorrect (a false negative).

- This product has not been FDA cleared or approved, but has been authorized by FDA under an EUA for use by authorized laboratories; use by laboratories certified under the Clinical Laboratory Improvement Amendments of 1988 (CLIA), 42 U.S.C §263a, that meet the requirements to perform moderate, high, or waived complexity tests. This product is authorized for use at the Point of Care (POC), i.e., in patient care settings operating under a CLIA Certificate of Waiver, Certificate of Compliance, or Certificate of Accreditation.
- **Serial testing should be performed in individuals with negative results at least twice over three days (with 48 hours between tests) for symptomatic individuals and three times over five days (with at least 48 hours between tests) for asymptomatic individuals. You may need to purchase additional tests to perform this serial (repeat) testing.**
- If the individual has symptoms longer than 7 days, you should consider testing at least three times over five days with at least 48 hours between tests.
- Federal Law restricts this test to sale by or on the order of a licensed practitioner (U.S. only).
- This product has been authorized only for the detection of proteins from SARS-CoV-2, not for any other viruses or pathogens.
- Laboratories within the United States and its territories are required to report all results to the appropriate public health laboratories.
- Do not use the test kit if the cassette pouch is damaged or improperly sealed.
- Do not use the test kit beyond expiration date printed on the outer package.
- Do not open the test cassette pouch until just before use. Once opened, the test card should be used immediately.
- **Do not read test results before 10 minutes or after 15 minutes. Results read before 10 minutes or after 15 minutes may lead to a false positive, false negative, or invalid result.**
- Use only the components supplied with the kit.
- Swabs, vials, and tests are single use only. Do not re-use test components.
- Freshly collected specimens must be tested within an hour.
- Do not interchange the kit components from different lots.
- Inadequate specimen collection can adversely affect results.
- Do not store or test specimens in Viral Transport Media (VTM).
- Temperature extremes and high humidity can adversely affect results.
- Wear appropriate personal protection equipment and gloves when running each test and handling patient specimens. Change gloves between handling of specimens suspected of COVID-19.
- Dispose of test cassette and materials as biohazardous waste in accordance with federal, state, and local requirements.
- Keep testing kit and kit components away from children and pets before and after use. Avoid contact with your skin, eyes, nose, or mouth. Do not ingest any kit components. The reagent solution contains harmful chemicals (see table below). If the solution contacts your skin, eyes, nose, or mouth, flush with large amounts of water. **If irritation persists, seek medical advice:** <https://poisonhelp.hrsa.gov/> or **1-800-222-1222.**

Hazardous Ingredient for the Reagent Solution		
Chemical Name	Harms (GHS) code for each ingredient	Concentration
Triton X-100	H302 Acute oral toxicity H315 Skin irritation H318 Serious eye damage H400 Short-term (acute) aquatic hazard H410 Long-term (chronic) aquatic hazard	1%

For more information on EUAs please visit: <https://www.fda.gov/emergency-preparedness-and-response/mcm-legal-regulatory-and-policy-framework/emergency-use-authorization> For the most up to date information on COVID-19, please visit: <https://www.cdc.gov/coronavirus/2019-ncov/index.html>

STORAGE AND STABILITY

- Store the kit between 2°C to 30°C (36°F to 86°F).
- Avoid direct exposure to sunlight.
- Do not freeze the content of the kit.
- The shelf life of the kit is indicated on the outer package.
- Do not use the test kit beyond its expiration date.
- The test device must remain in the sealed pouch until use.

QUALITY CONTROL

A built-in control in each test cassette acts as an internal procedural control for the test and is needed to assess if the test is working properly. The appearance of a red/purple Control line (C-line) is used to confirm sufficient flow of the sample along the membrane has occurred and the functional integrity of the test has been maintained. If the Control line does not develop within 10 minutes, the test is considered invalid and should be discarded. Retesting with a new cassette is recommended. If the expected results are not obtained, do not perform a test and contact technical support.

EXTERNAL QUALITY CONTROLS

External Quality controls, required to demonstrate the integrity of the cassette and to ensure the reagents and the assay procedure are working properly, are not included in the C-Sync™ COVID-19 Antigen Test kit but are available from BioSynchronicity Corporation. The External Quality Control kit is composed of one Negative Control swab, one SARS-CoV-2 Nucleocapsid Protein Positive Control swab, and one SARS-CoV-2 Spike Protein Positive Control swab. It is recommended that positive and negative external quality controls be run once with every new lot, every shipment, every new user, or in accordance with Local, State, and Federal regulations, or accreditation requirements.

External Quality Controls

for C-Sync™ COVID-19 Antigen Test

For use under the Emergency Use Authorization (EUA) only
 For use with the C-Sync™ COVID-19 Antigen Test only
 For *in vitro* diagnostic use
 For prescription use only

Instructions for Use (IFU)



COV9712ACS

INTENDED USE

The C-Sync™ COVID-19 Antigen Test External Quality Controls are to be used exclusively with the C-Sync™ COVID-19 Antigen Test to monitor the entire assay and provide assurance the test is performing within specifications

SUMMARY AND EXPLANATION OF THE TEST

The C-Sync™ COVID-19 Antigen Test Quality Controls are external quality controls formulated specifically to demonstrate the integrity of the cassette and to ensure the reagents are working properly. The Quality Controls consist of 1 SARS-CoV-2 nucleocapsid protein positive control swab, 1 SARS-CoV-2 spike protein positive control swab, and 1 negative control swab. It is the responsibility of each laboratory or healthcare setting using the C-Sync™ COVID-19 Antigen Test to establish an adequate quality assurance program to ensure the performance of the test kit under its specific locations and conditions of use. The Positive and Negative Controls should be run once with every new lot, every new shipment, and every new user, and quality control requirements should be followed in conformance with local, state, and federal regulations or accreditation requirements and the user laboratory's standard quality control procedures.

MATERIALS PROVIDED

- 1 Set:
- 1 C-Sync™ COVID-19 Antigen Test Non-Infectious Recombinant SARS-CoV-2 Nucleocapsid Protein Positive Control Swab.
 - 1 C-Sync™ COVID-19 Antigen Test Non-Infectious Recombinant SARS-CoV-2 Spike Protein Positive Control Swab.
 - 1 C-Sync™ COVID-19 Antigen Test Negative Control Swab.

TEST PROCEDURE

A. TEST PREPARATION

Wear appropriate personal protective equipment and gloves when running the test.

B. NEGATIVE CONTROL SWAB PROCEDURE


Remove a negative Quality Control swab from the package.



1


Do not open the External Quality Control swab until just before use. Once opened, the External Quality Control swab should be used immediately.

Remove the foil seal from the extraction buffer tube.



2

Carefully insert the swab inside the vial and plunge it to the bottom. Rotate the swab against the wall of the vial for 10 seconds



3

Tightly squeeze and release the bottom of the vial against the absorbent part of the swab. This will make the liquid buffer inside the vial rise and fall. Make sure to cover and wash the absorbent part of the swab. Twirl the swab while squeezing and washing the swab with the buffer. Repeat this squeezing and twirling technique 5 times.



4

Stop squeezing the bottom of the vial and rotate the swab against the wall of the vial for 10 seconds.



5

Remove the swab slowly while squeezing tightly the top part of the vial against the absorbent part of the swab to extract as much liquid as possible from the swab with the goal of drying the swab. Firmly place the dropper cap attached to the vial onto the top of the vial. Make sure the dropper cap is securely in place.



6



Place a C-Sync™ cassette on a flat surface
 Add 4 drops of the extraction buffer vial to the sample well of the cassette

7



Read the results at 10 minutes. The tests should not be interpreted past 15 minutes

The Control line (C-line) must be visible and the Test line (T-line) must not be visible. If the C-line is not visible or the T-line is visible, the results are invalid. Do not conduct tests with patients or report results. **Contact technical support.**

8

C. POSITIVE CONTROL SWAB PROCEDURE for each positive swab


Remove a positive Quality Control swab from the package.



1


Do not open the External Quality Control swab until just before use. Once opened, the External Quality Control swab should be used immediately.

Remove the foil seal from the extraction buffer tube.




2

Carefully insert the swab inside the vial and plunge it to the bottom. Rotate the swab against the wall of the vial for 10 seconds



3

Tightly squeeze and release the bottom of the vial against the absorbent part of the swab. This will make the liquid buffer inside the vial rise and fall. Make sure to cover and wash the absorbent part of the swab. Twirl the swab while squeezing and washing the swab with the buffer. Repeat this squeezing and twirling technique 5 times.



4

Stop squeezing the bottom of the vial and rotate the swab against the wall of the vial for 10 seconds.



5

Remove the swab slowly while squeezing tightly the top part of the vial against the absorbent part of the swab to extract as much liquid as possible from the swab with the goal of drying the swab. Firmly place the dropper cap attached to the vial onto the top of the vial. Make sure the dropper cap is securely in place.



6



Place a C-Sync™ cassette on a flat surface
 Add 4 drops of the extraction buffer vial to the sample well of the cassette

7



Read the results at 10 minutes. The tests should not be interpreted past 15 minutes

The Control line (C-line) must be visible and the Test line (T-line) must not be visible. If the C-line is not visible or the T-line is visible, the results are invalid. Do not conduct tests with patients or report results. **Contact technical support.**

8

EMERGENCY USE AUTHORIZATION


In the USA, this product has not been FDA cleared or approved; but has been authorized by FDA under an EUA for use by authorized laboratories; use by laboratories certified under CLIA, 42 U.S.C. §263a, that meet requirements to perform moderate, high or waived complexity tests. This product is authorized for use at the Point of Care (POC), i.e., in patient care settings operating under a CLIA certificate of waiver, certificate of compliance, or certificate of accreditation. This product has been authorized only (to be used with the C-Sync™ COVID-19 Antigen Test) for the detection of proteins from SARS-CoV-2, not for any other viruses or pathogens. In the USA, the emergency use of this product is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of *in vitro* diagnostics for detection and/or diagnosis of the virus that causes COVID-19 under Section 564(b)(1) of the Federal Food, Drug and Cosmetic Act, 21 U.S.C. § 360bbb-3(b)(1), unless the declaration is terminated or authorization is revoked sooner.

TECHNICAL SUPPORT

Phone: +1 (425) 898-3431
 E-mail: support@biosynchronicity.com
 www.biosynchronicity.com

 Manufacturer

 CE mark


 Contains sufficient for <n> test

 Store between 2°C and 30°C (36°F and 86°F)


 In Vitro Diagnostic device


 Keep dry


 Keep away from sunlight

 Consult Instructions for Use

 For prescription use only

 Do not reuse

 Do not use if package is damaged

 Use by date

 Catalog Number

External Quality Controls for C-Sync™ COVID-19 Antigen Test