IMPORTANT! How To Use This Test

Nano-Check™ COVID-19 Antigen **At-Home Test Quick Reference Instruction**

Catalog No. MD-8150

Serial testing should be performed in all individuals with negative results; individuals with symptoms of COVID-19 and initial negative results should be tested again after 48 hours. Individuals without symptoms of COVID-19, and with initial negative results, should be tested again after 48 hours and, if the 2nd test is also negative, a 3rd time after an additional 48 hours. You may need to purchase additional tests to perform this serial (repeat) testing

- If you test negative but continue to have symptoms of COVID-19, and both your first and second tests are negative, you may not have COVID-19, however you should follow-up with your healthcare provider.
- If your test is positive, then proteins from the virus that causes COVID-19 have been found in your sample and you likely have COVID-19.



Below are photos of actual positive tests Please note that the Ag line may be faint.



based tests such as PCR tests. If the test is negative but COVID-19-like symptoms, e.g., fever, cough, and/or shortness of breath continue, follow up testing for SARS-CoV-2 with a molecular test or testing for other respiratory disease should be considered. If applicable, seek follow up care with the primary health care provider.

- the country and informs public health decisions
- Manufactured for Nano-Ditech Corp. 259 Prospect Plains Road, Building K. Cranbury, NJ 08512, USA 1-855-297-7877 https://www.nanoditech.com/ P/N EP-3438-QR0 (May 2023) Cat. No. ND-MD8150

Nano-Check™ COVID-19 Antigen At-Home Test Quick Reference Instruction

Catalog No. MD-8150

The Nano-Check™ COVID-19 Antigen At-Home Test is only for in vitro diagnostic use under the Food and Drug Administration's Emergency Use Authorization. This product has not been FDA cleared or approved. Please read this instruction for use before using the test. For use with direct anterior nasal swab specimens.

INTENDED USE

The Nano-Check™ COVID-19 Antigen At-Home Test is a lateral flow immunoassay intended for the qualitative detection of nucleocapsid protein antigen from SARS-CoV-2.

This test is authorized for non-prescription home use with self-collected anterior nasal (nares) swab samples from individuals aged 14 years or older or adult collected anterior nasal (nares) swab samples from individuals aged two years or older. This test is authorized for individuals with symptoms of COVID-19 within the first 6 days of symptom onset when tested at least twice over three days with at least 48 hours between tests, and for individuals without symptoms or other epidemiological reasons to suspect COVID-19, when tested at least three times over five days with at least 48 hours between tests.

The Nano-Check™ COVID-19 Antigen At-Home Test does not differentiate between SARS-CoV or SARS-CoV-2.

Results are for the identification of SARS-CoV-2 nucleocapsid protein antigen, which is generally detectable in anterior nasal (nares) samples during the acute phase of infection. Positive results indicate the presence of viral antigens, but clinical correlation with past medical history and other diagnostic information is necessary to determine infection status. Positive results do not rule out bacterial infection or co-infection with other viruses and the agent detected may not be the definitive cause of disease. Individuals who test positive with Nano-Check™ COVID-19 Antigen At-Home Test should self-isolate and seek follow up care with their physician or healthcare provider as additional testing may be necessary.

All negative results are presumptive and confirmation with a molecular assay, if necessary for patient management, may be performed. Negative results do not rule out SARS-CoV-2 infection and should not be used as the sole basis for treatment or patient management decisions, including infection control measures such as isolating from others and wearing masks. Negative results should be considered in the context of an individual's recent exposures, history, and the presence of clinical signs and symptoms consistent with COVID-19.

Individuals who test negative and continue to experience COVID-19 like symptoms of fever, cough and/or shortness of breath may still have SARS-CoV-2 infection and should seek follow up care with their physician or healthcare provider.

Individuals should provide all results obtained with this product to their healthcare provider for public health reporting and to receive appropriate medical care. All healthcare providers will report all test results they receive from individuals who use the authorized product to relevant public health authorities in accordance with local, state, and federal requirements, using appropriate LOINC and SNOMED codes, as defined by the Laboratory In Vitro Diagnostics (LIVD) Test Code Mapping for SARS-CoV-2 Tests provided by CDC.

The Nano-Check™ COVID-19 Antigen At-Home Test is intended for non-prescription self-use and/or as applicable an adult lay user testing another person 2 years of age or older in a non-laboratory setting. The Nano-Check™ COVID-19 Antigen At-Home Test is only for in vitro diagnostic use under the Food and Drug Administration's Emergency Use Authorization. This product has not been FDA cleared or approved.

FREQUENTLY ASKED QUESTIONS

WHAT ARE THE KNOWN AND POTENTIAL RISKS AND BENEFITS OF THE TEST?

Potential risks include:

- Possible disconfort during sample collection.
 Possible incorrect test result (see Warnings and Result Interpretation sections for more information). Potential benefits include:
- The results, along with other information, can help you and your healthcare provider make informed recommendations about your care.
- The results of this test may help limit the potential spread of COVID-19 to your family and others in your community. For more information on EUAs go here: https://www.fda.gov/emergency-preparedness-and-response/mcm-
- legal-regulatory-and-policy-framework/emergency-use-authorization

WHAT IS THE DIFFERENCE BETWEEN AN ANTIGEN AND MOLECULAR TEST?

There are different kinds of tests for the SARS-CoV-2 virus that causes COVID-19. Molecular tests detect genetic material from the virus. Antigen tests, such as the Nano-Check™ COVID-19 Antigen At-Home Test, detect proteins from the virus. Due to the lower sensitivity of antigen tests, there is a higher chance this test will give you a false negative result when you have COVID-19 than a molecular test would.

HOW ACCURATE IS THIS TEST?

Clinical studies have shown that antigen tests more accurately determine whether you are infected with the virus that causes COVID-19 when taken multiple times across several days. Repeat testing improves test accuracy. This serial testing approach is recommended to minimize the risk of incorrect results. For more information on the performance of the test and how the performance may apply to you, please refer to the performance data in the Healthcare Provider Instructions for Use (IFU), available at https://www.nanoditech.

WHAT IF I HAVE A POSITIVE TEST RESULT?

A positive result means that it is very likely you have COVID-19 because proteins from the virus that causes COVID-19 were found in your sample. You should self-isolate from others and contact a healthcare provider for medical advice about your positive result..

WHAT IF I HAVE A NEGATIVE TEST RESULT?

A negative test result indicates that antigens from the virus that causes COVID-19 were not detected in your sample. However, if you have symptoms of COVID-19, and your first test is negative, you should test again in 48 hours since antigen tests are not as sensitive as molecular tests. If you do not have symptoms and received a negative result, you should test at least two more times with 48 hours in between tests for a total of three tests. If you have a negative result, it does not rule out SARS-CoV-2 infection; you may still be infected and you may still infect others. It is important that you work with your healthcare provider to help you understand the next steps you should take.

WHAT DOES AN INVALID TEST RESULT MEAN?

An invalid result means the test was not able to tell if you have COVID-19 or not. If the test is invalid, a new swab should be used to collect a new nasal specimen and you should test again with a new test.

IMPORTANT

Do not use this test as the only guide to manage your illness. Consult your healthcare provider if your symptoms persist or become more severe. Individuals should provide all results obtained with this product to their healthcare provider

WARNING, PRECAUTIONS, AND SAFETY INFORMATION

• Read all instructions carefully before performing the test. Failure to follow the instructions may result in inaccurate test results.

• In the USA, this product has not been FDA cleared or approved, but has been authorized by FDA under an Emergency Use Authorization. This product has been authorized only for the detection of proteins from

SARS-CoV-2, not for any other viruses or pathogens. The emergency use of this product is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of in vitro diagnostics for detection and/or diagnosis of COVID-19 under Section 564(b)(1) of the Federal Food, Drug, and Cosmetic Act, 21 U.S.C. § 360bbb-3(b)(1), unless the declaration is terminated, or authorization is revoked sooner

Serial testing should be performed in individuals with negative results at least twice over three days (with 48 hours between tests) for symptomatic individuals and three times over five days (with at least 48 hours between tests) for asymptomatic individuals. You may need to purchase additional tests to perform this serial (repeat) testing. An anterior nasal (nares) swab sample can be self-collected by an individual age 14 years and older.

Children aged 2 to 13 years should be tested by an adult. • Do not use on anyone under 2 years of age.

- Wear a safety mask or other face-covering when collecting a specimen from a child or another individual.
 Do not use if any of the test kit contents or packaging is damaged.
- Test components are single-use. Do not re-use.
- Do not touch the swab tip.
- Once opened, the test card should be used within 90 minutes. Do not read test results before 15 minutes or after 30 minutes. Results read before 15 minutes or after
- 30 minutes may lead to a false negative, false positive, or invalid result.
 Ensure that there is sufficient lighting for testing and interpretation.
 When collecting an anterior nasal (nares) swab sample, only use the swab provided in the kit.
 Do not use kit past its expiration date. For information about current expiration dates for at-home OTC
- COVID-19 diagnostic tests, visit <u>At-Home OTC COVID-19 Diagnostic Tests</u>. If you have had symptoms longer than 6 days you should consider testing at least three times over five days with at least 48 hours between tests.
- Freshly collected specimens should be processed as soon as possible.
- This product is only for the detection of proteins from SARS-CoV-2, not for any other viruses or pathogens.
 Testing should be performed in an area with good lighting and room temperature condition.
- Dispose of all materials in household waste.
 Wash hands thoroughly or use hand sanitizer before and after the test.
- Failure to follow the instructions for use may adversely affect test performance and/or invalidate the test result. Make sure to swirl and plunge the swab up and down in extraction buffer while squeezing the sides of the tube for 15 times; squeezing the swab head at least once or more in the reagent tube during the swab removal procedure. Insufficient swirling or squeezing of the swab head may produce false negative results.

Keep testing kit and kit components away from children and pets before and after use. Avoid contact with your skin, eyes, nose, or mouth. Do not ingest any kit components. The reagent solution contains harmful chemicals (see table below). If the solution contacts your skin, eyes, nose, or mouth, flush with large amounts of water. If irritation persists, seek medical advice: <u>https://www.poisonhelp.org</u> or 1-800-222-1222

Chemical Name	Harms (GHS Code) for each ingredient	Concentration
Gentamicin	Skin sensitization (H317)	4%
Sodium Azide	Acute Tox. 2 (Oral), H300, Acute Tox. 1 (Dermal), H310	0.09%

For more information on EUAs please visit: <u>https://www.fda.gov/emergency-preparedness-and-</u>

response/mcm-legal-regulatory-and-policy-framework/emergency-use-authorization • For the most up to date information on COVID-19, please visit: www.cdc.gov/COVID19

LIMITATIONS

- The performance of this test was established based on the evaluation of a limited number of clinical specimens collected between June 2022 and July 2022. The clinical performance has not been established for all circulating variants but is anticipated to be reflective of the prevalent variants in circulation at the time and location of the clinical evaluation. Performance at the time of testing may vary depending on the variants circulating, including newly emerging strains of SARS-CoV-2 and their prevalence, which change over time.
- There is a higher chance of false negative results with antigen tests than with laboratory-based molecular tests due to the sensitivity of the test technology. This means that there is a higher chance this test will give a false negative result in an individual with COVID-19 as compared to a molecular test, especially in samples with low viral load.
- All COVID-19 antigen test negative results are presumptive and confirmation with a molecular assay may be necessary.
- If you continue to have symptoms of COVID-19, and both your first and second tests are negative, you may not have COVID-19, however you should follow-up with a healthcare provider.
- If the test is positive, then proteins from the virus that causes COVID-19 have been found in the sample and the individual likely have COVID-19.
- These test results are shown as lines of color. Because these lines can be very faint, users with conditions affecting their vision such as far-sightedness, glaucoma, or color blindness are encouraged to seek assistance to interpret results accurately (e.g., reading glasses, additional light source, or another person). Incorrect test results may occur if a specimen is incorrectly collected or handled.
- This test detects both viable (live) and nonviable SARS-CoV-2. Test performance depends on the amount of virus (antigens) in the sample and may or may not correlate with viral culture results performed on the same sample.
- Test results must be evaluated in conjunction with other clinical data available to the physician.
- Positive test results do not exclude co-infection with other pathogens.
 Negative test results are not indicative of the presence/absence of other viral or bacterial pathogens.
- If the differentiation of specific coronaviruses and strains is needed, additional testing, in consultation with state or local public health departments, is required.
- False negative results may occur in individuals who have indicated or whose clinical status or history would indicate they are currently taking high doses of biotin. There was no interference up to 3,500 ng/mL of biotin in the samples.

STORAGE CONDITIONS

Store the Nano-Check™ COVID-19 Antigen At-Home Test between 36-86°F (2-30°C) in a place out of direct sunlight and out of reach of children. Reagents and devices must be used at room temperature 65-86°F (18-30°C) before use. The unsealed cassette may be used for 1.5 hours. It is recommended to use the test kit immediately after opening. The expiration date assigned at manufacturing is on the outer package box. For information about current expiration dates for Nano-Check™ COVID-19 Antigen At-Home Test, <u>visit At-Home</u> OTC COVID-19 Diagnostic Tests.

Index of Symbols





