

## Appendix 4. Supplemental Clinical Data

### *Primary Effectiveness Endpoint for 3 studies for multiple cohorts*

Table 4.1 shows the change in daytime ASBP for the Intent to treat (ITT), per protocol (PP), and complete ABP (CA) patient populations in each of the three studies.

**Table 4.1. Primary Effectiveness Endpoint for RADIANCE-HTN Studies, Daytime Ambulatory SBP (mmHg) - ITT, PP, CA**

	uRDN			Sham			Baseline Adjusted	
	Baseline	2 months	Change	Baseline	2 months	Change	Mean/Median Difference (95% CI) (uRDN - Sham) <sup>1</sup>	p-value <sup>1</sup>
<b>SOLO</b>								
ITT uRDN = 74 Sham = 72	150.3 ± 7.8 149.1 [134.8, 165.8]	141.9 ± 11.9 140.6 [113.3, 166.7]	-8.5 ± 9.3 -8.5 [-28.9, 13.6] (-10.6, -6.3)	150.0 ± 9.8 150.1 [134.5, 176.7]	147.9 ± 13.3 147.7 [118.1, 177.4]	-2.2 ± 10.0 -1.1 [-34.2, 25.0] (-4.5, 0.2)	-6.3 [-9.4, -3.1]	0.0001
PP uRDN = 64 Sham = 58	150.5 ± 7.3 149.6 [134.8, 165.7]	142.0 ± 11.7 141.0 [113.3, 166.7]	-8.5 ± 9.6 -8.6 [-28.9, 13.6] (-10.9, -6.1)	149.5 ± 9.9 148.8 [134.5, 176.7]	149.4 ± 12.0 148.7 [124.0, 177.4]	-0.1 ± 8.5 -0.8 [-22.2, 25.0] (-2.3, 2.1)	-8.2 [-11.5, -5.0]	<.0001
CA uRDN = 73 Sham = 71	73 150.5 ± 7.7 149.4 [134.8, 165.8]	73 141.9 ± 12.0 140.6 [113.3, 166.7]	73 -8.6 ± 9.3 -8.7 [-28.9, 13.6] (-10.7, -6.4)	71 149.9 ± 9.8 150.1 [134.5, 176.7]	71 147.7 ± 13.3 147.6 [118.1, 177.4]	71 -2.2 ± 10.1 -1.2 [-34.2, 25.0] (-4.6, 0.2)	-6.3 [-9.5, -3.1]	0.0001
<b>TRIO</b>								
ITT uRDN = 69 Sham = 67	150.0 ± 11.9 147.1 [134.5, 179.8]	141.0 ± 16.1 138.6 [105.8, 193.2]	-9.0 ± 14.5 -8.0 [-61.3, 26.0] (-12.5, -5.5)	151.1 ± 12.6 148.8 [133.8, 202.0]	146.3 ± 18.8 143.0 [107.6, 200.7]	-4.8 ± 15.9 -3.0 [-67.9, 47.9] (-8.6, -0.9)	-4.5 [-9.6, 0.6] -4.5 [-8.5, -0.3] <sup>2</sup>	0.0809 (0.0223*)

	uRDN			Sham			Baseline Adjusted	
	Baseline	2 months	Change	Baseline	2 months	Change	Mean/Median Difference (95% CI) (uRDN - Sham) <sup>1</sup>	p-value <sup>1</sup>
PP uRDN = 55 Sham = 57	149.9 ± 12.4 145.7 [134.5, 179.8]	140.8 ± 17.6 138.6 [105.8, 193.2]	-9.1 ± 14.8 -8.7 [-61.3, 26.0] (-13.2, -5.1)	150.2 ± 12.3 147.5 [134.8, 202.0]	145.1 ± 19.0 142.4 [107.6, 200.7]	-5.1 ± 15.9 -3.3 [-67.9, 47.9] (-9.3, -0.9)	-4.1 [-9.8, 1.7] -5.4 [-9.5, -1.3] <sup>2</sup>	0.1614 (0.0109*)
CA uRDN = 63 Sham = 67	150.3 ± 12.3 146.8 [134.5, 179.8]	140.4 ± 16.7 137.6 [105.8, 193.2]	-9.8 ± 14.9 -9.7 [-61.3, 26.0] (-13.6, -6.1)	151.1 ± 12.6 148.8 [133.8, 202.0]	146.3 ± 18.8 143.0 [107.6, 200.7]	-4.8 ± 15.9 -3.0 [-67.9, 47.9] (-8.6, -0.9)	-5.3 [-10.6, -0.0] -5.8 [-9.7, -1.6] <sup>2</sup>	0.0491 (0.0051*)
<b>RADIANCE-II</b>								
ITT uRDN = 145 Sham = 73	150.2 ± 8.6 149.0 [135.0, 169.0]	142.3 ± 13.4 140.0 [104.0, 183.0]	-7.9 ± 11.6 -8.0 [-41.0, 26.0] (-9.8, -6.0)	151.3 ± 9.0 151.0 [136.0, 170.0]	149.5 ± 11.1 149.0 [127.0, 176.0]	-1.8 ± 9.5 -1.0 [-30.0, 26.0] (-4.1, 0.4)	-6.3 [-9.3, -3.2] -7.0 [-9.0, -4.0] <sup>2</sup>	<.0001 (<.0001*)
PP uRDN = 131 Sham = 63	149.8 ± 8.5 148.0 [135.0, 169.0]	141.9 ± 13.4 139.0 [104.0, 183.0]	-7.9 ± 11.7 -9.0 [-41.0, 26.0] (-9.9, -5.9)	150.0 ± 8.3 148.0 [136.0, 168.0]	149.0 ± 10.8 148.0 [127.0, 176.0]	-1.0 ± 9.0 -1.0 [-22.0, 26.0] (-3.3, 1.3)	-6.9 [-10.2, -3.6] -7.0 [-10.0, -4.0] <sup>2</sup>	<.0001 (<.0001*)
CA uRDN = 145 Sham = 73	150.2 ± 8.6 149.0 [135.0, 169.0]	142.3 ± 13.4 140.0 [104.0, 183.0]	-7.9 ± 11.6 -8.0 [-41.0, 26.0] (-9.8, -6.0)	151.3 ± 9.0 151.0 [136.0, 170.0]	149.5 ± 11.1 149.0 [127.0, 176.0]	-1.8 ± 9.5 -1.0 [-30.0, 26.0] (-4.1, 0.4)	-6.3 [-9.4, -3.2] -7.0 [-9.0, -4.0] <sup>2</sup>	<.0001 (<.0001*)

Data displayed as Mean±SD, Median [Range], and 95% CI for change. Change is calculated as 2 months - Baseline

<sup>1</sup>Mean difference with 95% CI and p-value value from ANCOVA, adjusting for baseline value. In the event that the change from baseline in either cohort is non-normal, the p-value (\*) from a baseline adjusted ANCOVA on the ranks is also provided.

<sup>2</sup>Median Difference - Hodges-Lehmann estimate of location shift and 95% asymptotic CI which are not associated with the p-value (\*) via a baseline adjusted ANCOVA on the ranks

Note that all listed p-values are not adjusted for multiplicity.

Secondary Effectiveness Endpoints for TRIO – Completers population

Table 4.2 shows the secondary endpoints at 2 months of 24-hour SBP/SBP, Nighttime SBP/DBP, and Daytime DBP at 2 months for the CA population. The differences between the uRDN and Sham groups are larger than in the ITT population. The median differences in 24-hour and nighttime SBP/DBP are statistically significant, and only 24-hour ASBP were above the clinically meaningful threshold of 5.0 mmHg.

**Table 4.2: Secondary Endpoint BP Change from Baseline to 2 Months (Complete 2M ABPM Population)**

	Renal Denervation (n=69)			Sham Procedure (n=67)			Unadjusted	Baseline Adjusted	
	Baseline	2 month	Difference	Baseline	2 month	Difference	Median Difference (95% CI) (uRDN – Sham) <sup>1</sup>	Mean Difference (95% CI) (uRDN - Sham) <sup>2</sup>	p-value <sup>2</sup>
24 Hour Ambulatory systolic blood pressure (mmHg)	63 144.0 ± 13.9 138.8 [123.3, 180.1]	63 134.5 ± 16.4 132.5 [102.5, 194.7]	63 -9.5 ± 14.3 -9.4 [-55.2, 30.8] (-13.1, -5.9)	67 145.4 ± 14.0 142.4 [125.0, 201.6]	67 140.5 ± 18.7 138.2 [103.1, 197.2]	67 -4.8 ± 16.5 -2.9 [-62.7, 61.4] (-8.9, -0.8)	-5.6 [-9.5, -1.3]	-5.1 (-10.3, 0.1)	0.0523 (0.0043*)
24 Hour Ambulatory diastolic blood pressure (mmHg)	63 89.1 ± 8.5 87.6 [75.6, 113.1]	63 83.3 ± 11.4 81.9 [63.0, 122.5]	63 -5.7 ± 8.9 -5.6 [-34.5, 20.5] (-8.0, -3.5)	67 89.5 ± 9.5 87.2 [76.8, 131.8]	67 85.8 ± 12.0 84.6 [66.5, 125.3]	67 -3.7 ± 10.8 -2.4 [-44.5, 38.4] (-6.4, -1.1)	-2.7 [-5.3, 0.0]	-2.1 (-5.5, 1.2)	0.2141 (0.0423*)
Nighttime Ambulatory systolic blood pressure (mmHg)	63 134.1 ± 18.4 129.6 [104.7, 181.2]	63 125.2 ± 18.5 123.7 [91.2, 196.9]	63 -8.9 ± 16.3 -9.0 [-49.2, 39.2] (-13.0, -4.8)	67 136.4 ± 18.6 132.4 [100.5, 199.4]	67 131.9 ± 20.9 129.9 [89.1, 193.1]	67 -4.5 ± 19.5 -1.8 [-72.1, 77.9] (-9.3, 0.2)	-5.0 [-10.1, 0.5]	-5.3 (-11.0, 0.4)	0.0692 (0.0146*)
Nighttime Ambulatory diastolic blood pressure (mmHg)	63 81.4 ± 11.0 79.7 [60.0, 114.7]	63 75.8 ± 12.5 73.9 [54.9, 121.0]	63 -5.5 ± 10.4 -5.7 [-28.9, 23.8] (-8.2, -2.9)	67 81.3 ± 12.1 80.3 [57.9, 125.5]	67 78.4 ± 13.2 77.0 [49.4, 119.7]	67 -2.8 ± 12.9 -2.0 [-54.2, 48.3] (-6.0, 0.3)	-3.6 [-7.0, -0.1]	-2.7 (-6.4, 1.1)	0.1639 (0.0249*)
Daytime Ambulatory diastolic blood pressure (mmHg)	63 94.0 ± 7.9 91.6 [82.6, 112.6]	63 88.1 ± 12.1 84.1 [66.7, 123.8]	63 -5.9 ± 9.5 -6.0 [-39.3, 18.7] (-8.2, -3.5)	67 94.6 ± 9.1 91.6 [82.9, 136.1]	67 90.7 ± 12.2 89.5 [68.8, 137.1]	67 -3.9 ± 10.5 -2.0 [-41.7, 30.2] (-6.5, -1.3)	-2.5 [-5.4, 0.3]	-2.1 (-5.5, 1.4)	0.2367 (0.0811*)

Data displayed as Mean ± SD, Median [Range], and 95% CI for change. Change is calculated as 2 months - Baseline

<sup>1</sup>Hodges-Lehmann estimate of location shift and 95% asymptotic CI.

<sup>2</sup>Mean difference with 95% CI and p-value value from ANCOVA, adjusting for baseline value. In the event that the change from baseline in either cohort is non-normal, the p-value (\*) from a baseline adjusted ANCOVA on the ranks is also provided.













Note that all listed p-values are not adjusted for multiplicity.

*Patient Preference Study Supplementary Data*

Table 4.3 shows the patient attributes and levels for ReCor’s patient preference study. Table 4.4 shows an example of a discrete choice experiment (DCE).

**Table 4.3 Attributes and Levels**

Attributes		Possible Levels	
Patients who have CV event in next 10 years		Under 55 years old	55 years old and over
		15%	30%
		30%	45%
		45%	60%
Current treatment	Number of pills per day	0 pills a day	
		1 pill a day	
		2 pills a day	
		3 pills a day	
	Minimally invasive procedure	No procedure	
		Procedure	
Patients who have mild-to-moderate side effects requiring more doctor visits		20%	
		30%	
		40%	
		50%	
		70%	
Patients who have serious side effects requiring hospitalization		0%	
		5%	
		10%	
		20%	
Patients who have serious side effects requiring a procedure		0%	
		5%	
		10%	
		No additional treatment	
Future treatments required to manage your blood pressure		Additional pill in 1 year	
		Additional pill in 3 years	
		Additional procedure in 3 years	
		Additional procedure in 7 years	

	Treatment A	Treatment B
Patients who have CV event in next 10 years	 15 out of 100 (15%)	 15 out of 100 (15%)
Current treatment	 No procedure 2 pills per day	 A minimally invasive procedure No pill
Future treatments required to manage your blood pressure	 1 year Additional pill in 1 year (Total 3 pills per day)	 3 years Additional pill in 3 years (Total 1 pill per day)
Patients who have mild-to-moderate side effects requiring more doctor visits	 50 out of 100 (50%)	 30 out of 100 (30%)
Patients who have serious side effects requiring hospitalization	 0 out of 100 treated (0%)	 0 out of 100 (0%)
Patients who have serious side effects requiring a procedure	 0 out of 100 (0%)	 5 out of 100 (5%)

**Figure 4.4. Example Discrete Choice Experiment**