

# **Premarket Approval Application (PMA) for ReCor Medical's Paradise Ultrasound Renal Denervation System**

**Circulatory System Devices Advisory Committee Meeting  
August 22, 2023**



# Introduction, Background, Clinical Study Design

**Paul Warren, PhD**  
**Biomedical Engineer**  
**Office of Cardiovascular Devices**

# Review Team

Lead Reviewer

Clinician (Vascular Surgery)

Clinician (Cardiology)

Clinician (Nephrology)

Statistics

Statistics

Patient Preference Study

Assistant Director, PIDT

Chief Medical Officer, OHT2

Paul Warren, PhD

Robert Lee, MD

Meir Shinnar, MD, PhD

Douglas Silverstein, MD

Adrijo Chakraborty, PhD

Wei-Chen Chen, PhD

David Gebben, PhD

Misti Malone, PhD

Andrew Farb, MD

# Outline

- Background
- Device Description and Proposed Indications for Use
- Clinical Study Design
- Clinical Study Results
  - Safety
  - Effectiveness
  - Patient preference study
- Post-approval Study
- Conclusions

# Clinical Background

- Hypertension (HTN) is a major public health problem
  - Prevalence ~45% of US adults
  - Higher rate among African Americans (57.1%) vs. Caucasians (43.6%) & Hispanics (43.7%) (NHANES, 2017-2018)<sup>1</sup>
- Associated with increased risk of serious conditions including<sup>2</sup>
  - Stroke
  - Heart disease
  - Heart failure
  - Noncardiac vascular disease
  - Renal Disease
- BP medications are the mainstay of HTN therapy, but:
  - BP medication adherence present in approximately 60% of patients<sup>3</sup>
  - Target BP achieved in approximately 45% of patients<sup>4</sup>

1. National Health and Nutrition Examination Survey Fact Sheet. CDC. July 2020. Available at: [https://www.cdc.gov/nchs/data/factsheets/factsheet\\_nhanes.pdf](https://www.cdc.gov/nchs/data/factsheets/factsheet_nhanes.pdf)
2. Carey RM et al. Prevention and Control of Hypertension: JACC Health Promotion Series. *J Am Coll Cardiol* 72(11). 2018.
3. Choudhry NK et al. Medication Adherence and Blood Pressure Control: A Scientific Statement From the American Heart Association. *Hypertension*. 79(1). 2022.
4. Dorans KS et al. Trends in prevalence and control of hypertension according to the 2017 ACC/AHA Guidelines. *J Am Heart Assoc* 7(11). 2018.

# Defining Hypertension (1)

2017 US Societal Guideline Classification of Blood Pressure in Adults<sup>5</sup>

Category	SBP		DBP
Normal	<120 mmHg	AND	<80 mmHg
Elevated	120-129 mmHg	AND	<80 mmHg
Hypertension			
Stage 1	130-139 mmHg	OR	80-89 mmHg
Stage 2	≥140 mmHg	OR	≥90 mmHg

5. Whelton PK et al. ACC/AHA/AAPA/ABC/ACPM/AGS/APhA/ASH/ASPC/NMA/PCNA Guideline for the prevention, detection, evaluation, and management of high blood pressure in adults: A report of the American College of Cardiology / American Heart Association task force. *Circulation* 138(17). 2018.

# Defining Hypertension (2)


- Uncontrolled HTN: Above BP goal
  - Due to non-adherence to treatment; or
  - Despite adherence to treatment
- Resistant HTN: Above BP goal despite the use of 3 HTN medications (including a diuretic) with complementary mechanisms of action

# Role of Renal Physiology in Hypertension

- Renal vasculature innervated by mainly efferent sympathetic nerves
- Stimulation of efferent nerves leading to:
  - Increased reabsorption of Na and water
  - Reduced renal blood flow and GFR (vasoconstriction)
  - Increased activity of the RAAS



**Increased BP**





# Renal Denervation (RDN)

- An approach to reduce renal sympathetic activity by ablating the surrounding nerves
- Early single-arm clinical studies of percutaneous RDN technologies showed large decreases in BP
- However, initial sham-controlled trials did not see same degree of BP reduction, and no significant difference between RDN treatment and sham
- After denervation, some animal studies show re-innervation<sup>6-8</sup>
  - If re-innervation occurs in humans, sustained BP reduction could be impacted

6. Mulder J et al. Renal sensory and sympathetic nerves reinnervate the kidney in a similar time-dependent fashion after renal denervation in rats. *Am J Physiol* 304(8). 2013.  
 7. Booth LC et al. Reinnervation following catheter-based radio-frequency renal denervation. *Exp Physiol* 100(5). 2015.  
 8. Kiuchi MG et al. Renal denervation update from the International Sympathetic Nervous System Summit: JACC State-of-the-Art Review. *J Am Coll Cardiol* 73(23). 2019.

# 2018 FDA Advisory Committee on HTN Devices

- Discussed clinical trial designs to evaluate safety and effectiveness of devices for HTN
- Key Panel recommendations<sup>9</sup>
  - Sham controlled trials
  - Trial designs
    - Medication withdrawal (off-med) study
    - Standardized BP medication (on-standardized med) study
  - Ambulatory BP measurement (ABPM) used as the primary BP assessment method
  - A  $\geq 5$  mmHg difference in BP reduction between active treatment and sham as clinically significant
  - Value of patient preference information

# Paradise Ultrasound Renal Denervation (uRDN) System



Catheter tip with ultrasound transducer and balloon



Ultrasound generator and display screen

## Proposed indications for use

The Paradise uRDN System is indicated to reduce blood pressure in patients with *uncontrolled hypertension, who may be inadequately responsive to, or who are intolerant to anti-hypertensive medications.*

# Breakthrough Devices Program

- Paradise uRDN System granted breakthrough status in Dec 2020 for patients with resistant or uncontrolled hypertension
- Breakthrough devices may provide for more effective treatment or diagnosis of life-threatening or irreversibly debilitating diseases or conditions
- Intended to provide patients with timely access to certain devices by expediting their development, assessment, and review

*Contains Nonbinding Recommendations*

## **Breakthrough Devices Program**

---

### **Guidance for Industry and Food and Drug Administration Staff**

Document issued on December 18, 2018.

The draft of this document was issued on October 25, 2017.

This document supersedes “Expedited Access for Premarket Approval and De Novo Medical Devices Intended for Unmet Medical Need for Life Threatening or Irreversibly Debilitating Diseases or Conditions,” issued on April 13, 2015.

# Breakthrough Device Program

- **Does allow for:**
  - Interactive and timely communication with FDA
  - Prioritized submission review
  - Efficient and flexible clinical study design
  - Expedited manufacturing and quality systems review
  - Pre/Postmarket balance of data collection
- **Does not** alter/reduce the statutory requirement for premarket approval
  - A reasonable assurance of safety and effectiveness

# Pre/Postmarket Balance of Data Collection



- FDA may accept greater uncertainty for a premarket submission along with timely postmarket data collection if the uncertainty is sufficiently balanced
- Benefit/Risk considerations include:
  - Probable benefits from earlier access
  - vs.
  - Probable risk of harm should postmarket data show that the device is ineffective or unsafe



# Nonclinical and Preclinical Device Evaluation


- Catheter Engineering Testing
  - Bench testing
  - Ultrasound output and delivery
- Generator Engineering Testing
  - Electrical safety
  - Software validation
  - Cybersecurity
- System Compatibility
- Biocompatibility
- Sterilization & Packaging
- Preclinical Animal Studies

**No outstanding non-clinical study issues**

# CLINICAL STUDY DESIGN



# Clinical Study Design Overview



## Medication Escalation

- Between 2-6 months to goal ( $\leq 135$  mmHg)

## Blinding:

- SOLO and TRIO: Thru 6 months
- RADIANCE II: Thru 12 months

## Sham crossover to RDN permitted:


- SOLO and TRIO: At 6 months
- RADIANCE II: At 12 months

# Key Enrollment Criteria

- Between 18 and 75 years old
- History of HTN
- Suitable renal anatomy

	<b>SOLO Off-Med (Mild-to-moderate HTN)</b>	<b>RADIANCE-II Off-Med (Stage 2 HTN)</b>	<b>TRIO On-Standardized Med (Resistant HTN)</b>
<b>OBP</b>	$\geq 140/90$ & $< 180/110$ mmHg on 0, 1 or 2 meds; <i>or</i> $\leq 140/90$ on 1 or 2 meds	$\geq 140/90$ & $< 180/120$ mmHg on 0, 1, or 2 meds; <i>and</i> Previous or currently prescribed BP meds	$\geq 140/90$ on $\geq 3$ meds (including a diuretic)
<b>Daytime ABP</b>	$\geq 135/85$ & $< 170/105$ mmHg after washout	$\geq 135/85$ & $< 170/105$ mmHg after washout	$\geq 135/85$ mmHg after stabilization
<b>BP Meds</b>	0, 1, or 2	0, 1, or 2	At least 3
<b>Sample size</b>	146 (1:1 randomization)	224 (2:1 randomization)	136 (1:1 randomization)

**OBP/ABP = Office/ambulatory blood pressure**



# Statistical Analysis Plan

**Wei-Chen Chen, PhD**

**Statistician**

**Office of Clinical Evidence and Analysis**

# Paradise uRDN Clinical Trials

## Differences in statistical analysis plan

### 1. SOLO & TRIO:

- No prespecified primary safety endpoint
- Primary effectiveness endpoint was prespecified and tested

### 2. RADIANCE-II:

- Primary safety endpoint was prespecified and tested
- Primary effectiveness endpoint was prespecified and tested

# Primary Safety Endpoint

## *Analysis Population: RADIANCE-II*

- The primary safety endpoint was defined as the occurrence of at least one of major adverse events (MAE) at 30 days
  - All-cause mortality
  - New onset acute end-stage renal disease
  - Significant embolic event resulting in end-organ damage
  - Renal artery perforation requiring invasive intervention
  - Renal artery dissection requiring an invasive intervention
  - Major vascular complications
  - Hospitalization for hypertensive or hypotensive crisis
  - Hospitalization for major cardiovascular or hemodynamic-related events
  - New stroke or new MI

or

New renal artery stenosis (RAS), defined as a >70% diameter stenosis confirmed by CTA or MRA at 6 months

**Additional analysis for pooled SOLO, TRIO, & RADIANCE-II trials**

# Primary Safety Evaluation

## *Statistical Hypothesis and Analysis*

- **Safety event rate performance goal (PG) = 9.8%**, derived from literature review
- The primary safety null and alternative hypotheses are:
  - $H_0: \pi \geq 9.8\%$
  - $H_a: \pi < 9.8\%$where  $\pi$  is the proportion of subjects who had experienced at least one safety endpoint event.
- The upper limit of one-sided exact 95% confidence interval was used to determine endpoint success (i.e., less than PG).

# Primary Effectiveness Endpoint

For each of three trials (SOLO, TRIO, RADIANCE-II), the primary effectiveness endpoint was defined as the change in daytime ASBP from baseline to two months

–Daytime: 7:00 AM – 10:00 PM

# Primary Effectiveness Endpoint

## *Statistical Hypothesis*



### Hypothesis

$$H_0: \beta_{\text{txt}} = 0$$

$$H_a: \beta_{\text{txt}} \neq 0$$

where  $\beta_{\text{txt}}$  is the regression coefficient from the **Analysis of Covariance (ANCOVA)** model:

$$Y = \beta_0 + \beta_{\text{txt}} * X_{\text{trt}} + \beta_{\text{bl}} * X_{\text{bl}} + \text{residual}$$

and

- $Y$  : the reduction in daytime SBP from baseline to 2 months post procedure
- $\beta_0$  : the intercept
- $\beta_{\text{txt}}$  : the regression coefficient associated with the treatment effect (i.e., parameter of interest)
- $X_{\text{trt}}$  : an indicator variable with a value of 1 for uRDN and 0 for Sham
- $\beta_{\text{bl}}$  : the regression coefficient associated with the baseline daytime SBP
- $X_{\text{bl}}$  : the baseline daytime SBP (i.e., baseline adjustment)



# Primary Effectiveness Endpoint

## *Statistical Analysis*



- Primary analysis based on Intention-to-treat (ITT) population
- The treatment effect ( $\beta_{\text{txt}}$ ) was tested at the two-sided 0.05 alpha level
  - Mean difference between uRDN and Sham adjusting with baseline BP was reported
- Additional prespecified analysis
  - When the assumption of normality for the ANCOVA model may be in question, an **ANCOVA on the ranks** (Quade (1967)) was performed.
- Exploratory analysis
  - Median (Hodges-Lehmann estimate): Summarizes the distribution of BP reduction without adjusting for baseline BP

# Primary Effectiveness Endpoint

## *Missing Data*



### Imputation

- In all three trials (SOLO, TRIO, and RADIANCE-II)
  - For subjects who met “High BP Action” escape criteria, the last BP measurement prior to medication change was used
- In SOLO and TRIO
  - For subjects missing BP values, a value of zero was used for the BP reduction
- In RADIANCE-II
  - Multiple imputation methods
  - Imputation model included age, sex, and baseline ambulatory SBP
  - Twenty (20) imputed datasets were produced and combined via Rubin’s rules for statistical inference

# Secondary and Additional Effectiveness Endpoints

- Secondary endpoints
  - Change from baseline to two months
    - 24-hour SBP & DBP
    - Nighttime SBP & DBP
    - Daytime DBP
    - Office SBP & DBP
    - Home SBP and DBP
- Additional endpoints
  - Incidence of BP reduction ( $\geq 5, 10, 15, 20$  mmHg)
  - Proportion with BP control ( $< 135/85$  mmHg)
  - Medication burden
    - Number of antihypertensive medication
    - Antihypertensive medication load index

$$\sum_{AH\ Meds} \left( class\ weight \frac{prescribed\ dose}{standard\ dose} \right)$$

# Subgroup Analysis

## *Primary Effectiveness Endpoint*

- 1) Age
- 2) Race
- 3) Sex
- 4) Geography: US vs EU/UK
- 5) Baseline daytime ASBP (<median vs. ≥median)
- 6) Baseline office SBP (<median vs. ≥median)
- 7) Baseline home SBP (<median vs. ≥median)
- 8) Average baseline 24-hour ABP heart rate (above and below median)
- 9) eGFR (<60 vs ≥60)
- 10) Abdominal obesity (split for male >102cm and ≤102cm; and for female >88cm and ≤88cm)
- 11) Number of ablations performed (Treatment Group only)
- 12) Presence of untreated accessory arteries (Treatment Group only)
- 13) Balloon size (Treatment Group only)
- 14) Pre/post COVID-19 enrollment

# Statistical Considerations

- Gatekeeping procedure used to control overall type I error rate for secondary endpoints
  - 24-hour SBP & DBP (SOLO, TRIO, RADIANCE-II)
  - Nighttime SBP & DBP (SOLO and TRIO only)
  - Daytime DBP (RADIANCE-II only)
  - Office SBP & DBP (RADIANCE-II only)
  - Home SBP and DBP (RADIANCE-II only)
- P-values in other secondary, additional, and exploratory analyses were not adjusted for multiplicity and should be interpreted with caution.




## Clinical Study Results

Douglas Silverstein, MD

Nephrologist

Office of Gastrorenal, ObGyn, General Hospital,  
and Urology Devices

# Clinical Study Design Overview



## Medication Escalation

- Between 2-6 months to goal ( $\leq 135$  mmHg)


## Blinding:

- SOLO and TRIO: Thru 6 months
- RADIANCE II: Thru 12 months


## Sham crossover to RDN permitted:

- SOLO and TRIO: At 6 months
- RADIANCE II: At 12 months


## SOLO (Off-Med) Accountability



## RADIANCE-II (Off-Med) Accountability



## TRIO (On-Standardized Med) Accountability





## Key Baseline Patient Features

Measure	SOLO (Off-Med)		RADIANCE-II (Off-Med)		TRIO (On-Standardized Med)	
	uRDN (n=74)	Sham (n=72)	uRDN (n=150)	Sham (n=74)	uRDN (n=69)	Sham (n=67)
<b>Sex</b>						
Male	62.2%	54.2%	68.7%	77.0%	81.2%	79.1%
Female	37.8%	45.8%	31.3%	23.0%	18.8%	20.9%
<b>Age</b>	54.4 ± 10.2	53.8 ± 10.0	55.1 ± 9.9	54.9 ± 7.9	52.3 ± 7.5	52.8 ± 9.1
<b>Geography</b>						
US	47.3%	47.2%	66.7%	62.2%	40.6%	37.3%
OUS	52.7%	52.8%	33.3%	37.8%	59.4%	62.7%
<b>Race</b>						
American Indian or Alaska Native	0.0%	0.0%	0.0%	0.0%	0.0%	1 (1.52%)
Asian	1.3%	0.0%	0.0%	1.3%	1.5%	1.5%
Black	16.2%	18.0%	14.0%	20.2%	20.6%	19.7%
Caucasian	81.0%	72.2%	76.0%	75.6%	66.2%	77.3%
Hispanic or Latino	1.3%	5.5%	10.0%	2.7%	7.4%	0.0%
Native Hawaiian or other Pacific Islander	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%
Other/Mixed Race	0.0%	0.0%	10.0%	2.7%	4.41%	0.0%
<b>BMI</b>	29.9 ± 5.9	29.0 ± 5.0	30.1 ± 5.2	30.6 ± 5.2	32.8 ± 5.7	32.6 ± 5.4
<b>Abdominal circumference (cm)</b>	101.5 ± 14.2	98.5 ± 15.1	102.4 ± 12.3	104.3 ± 13.1	109.4 ± 15.5	109.2 ± 12.9
<b>Office SBP (mmHg)</b>	142.6 ± 14.7	144.6 ± 15.9	155.8 ± 11.1	154.3 ± 10.6	161.9 ± 15.5	163.6 ± 16.8
<b>Office DBP (mmHg)</b>	92.3 ± 10.1	93.6 ± 8.3	101.3 ± 6.7	99.1 ± 5.6	105.1 ± 11.6	103.3 ± 12.7

# BP Medications at Screening



Number of BP Meds	Randomized (%)	uRDN (%)	Sham (%)
<b>SOLO (Off-Med)</b>			
0	19%	16%	22%
1	42%	45%	39%
2	38%	38%	38%
3	1%	1%	1%
<b>RADIANCE II (Off-Med)</b>			
0	34%	36%	31%
1	34%	35%	34%
2	31%	29%	34%
3	0%	0%	0%
4	0%	0%	1%
<b>TRIO (On-Standardized Med)</b>			
3	40%	39%	42%
4	32%	29%	36%
5	19%	23%	15%
6+	8%	9%	8%

# **SAFETY ENDPOINT RESULTS**

# Primary Safety Endpoint

## *Analysis Population: RADIANCE-II*

- The primary safety endpoint was defined as the occurrence of at least one of major adverse events (MAE) at 30 days
  - All-cause mortality
  - New onset acute end-stage renal disease
  - Significant embolic event resulting in end-organ damage
  - Renal artery perforation requiring invasive intervention
  - Renal artery dissection requiring an invasive intervention
  - Major vascular complications
  - Hospitalization for hypertensive or hypotensive crisis
  - Hospitalization for major cardiovascular or hemodynamic-related events
  - New stroke or new MI

or

New renal artery stenosis (RAS), defined as a >70% diameter stenosis confirmed by CTA or MRA at 6 months

**Additional analysis for pooled SOLO, TRIO, & RADIANCE-II trials**

# Primary Safety Results: RADIANCE-II and Pooled uRDN Subjects

	N	Rate	95% CI	Performance Goal (PG)	p-value
<b>RADIANCE-II</b>	150	0.0%	0 - 1.63%	9.8%	<0.0001
<b>Pooled uRDN subjects</b>	367	1.1%	0.3% - 2.77%		

# Individual Safety Events for uRDN Subjects - Initial Treatment & Crossover



	SOLO	RADIANCE II	TRIO	Combined
<b>N (Initial + Crossover)</b>	109 (72 + 37)	168 (149 +19)	90 (69 +21)	367
<b>30-day events</b>				
All-cause mortality	0	0	2	2
Major vascular complications	0	0	2	2
Hospitalization for hypertensive or hypotensive crisis	0	0	1	1
Hospitalization for major CV or hemodynamic events	1	0	0	1
<b>6-month RAS</b>	0	0	0	0
<b>Total</b>	<b>1 (0.9%)</b>	<b>0</b>	<b>5 (3.3%)</b>	<b>6 (1.1%)</b>

*Note that endpoints with 0 events are not shown here  
Number of Events (% of Subjects with Events)*

# RADIANCE-II Renal Artery Diameter Stenosis at 6 & 12-Months by CTA/MRA



Follow-up Timepoint	Total Evaluable	1-30% Stenosis	31-50% Stenosis	51-70% Stenosis	>70% Stenosis
<b>6 Months</b>					
uRDN	137	1.5% (2)	0.7% (1)	0.0% (0)	0.0% (0)
Sham	58	0.0% (0)	3.4% (2)	0.0% (0)	0.0% (0)
<b>12 Months</b>					
uRDN	112	2.7% (3)	2.7% (3)	1.8% (2)	0.0% (0)
Sham	14	0.0% (0)	0.0% (0)	0.0% (0)	0.0% (0)

*Data based on core lab analysis*

# Renal Function (GFR) at 12 Months in SOLO

	uRDN (n = 68)			Sham (n = 67)			Baseline Adjusted	
	Baseline	12 months	Change	Baseline	12 months	Change	Mean Difference (95% CI) (RDN - Sham) <sup>1</sup>	p-value
<b>eGFR</b>	84.46 ± 16.47	83.47 ± 14.61	-0.99 ± 11.18	82.23 ± 15.79	85.63 ± 16.99	3.40 ± 11.53	-3.79 (-7.40, -0.19)	0.0391
<b>Serum Creatinine</b>	0.90 ± 0.18	0.90 ± 0.18	0.01 ± 0.11	0.90 ± 0.18	0.87 ± 0.17	-0.03 ± 0.10	0.04 (0.00, 0.07)	0.0326

*Note not all p-values are adjusted for multiplicity.*

*\*Data shown for SOLO study. TRIO study results show statistically similar GFR results in the RDN and Sham arms.*



# **EFFECTIVENESS ENDPOINT RESULTS**

# Primary Effectiveness Endpoint

Mean difference in change of daytime (7 AM-10 PM) ASBP between treatment and sham groups from baseline to 2 months after RDN (or Sham) procedure

- Primary analysis population: ITT cohort

# Primary Effectiveness Results at 2 Months

ITT Cohort	uRDN change (mmHg)	Sham change (mmHg)	Baseline BP Adjusted Mean Treatment Effect (mmHg) <sup>1</sup>	p-value
<b>Off-Med Trials</b>				
SOLO	-8.5 ± 9.3 (74)	-2.2 ± 10.0 (72)	-6.3 (-9.4, -3.1)	0.0001
RADIANCE-II	-7.9 ± 11.6 (145)	-1.8 ± 9.5 (73)	-6.3 (-9.3, -3.2)	<.0001
<b>On-Standardized-Med Trial</b>				
TRIO	-9.0 ± 14.5 (69)	-4.8 ± 15.9 (67)	-4.5 (-9.6, 0.6)	0.0809
<p>Data presented as mean ± SD (n); difference presented as mean (95% CI); p-value via a baseline adjusted ANCOVA (two-sided).</p> <p><sup>1</sup> Negative value favors uRDN.</p>				

# Additional TRIO Results at 2 Months

ITT Cohort	uRDN change (mmHg)	Sham change (mmHg)	Unadjusted Median Treatment Effect (mmHg) <sup>1</sup>	p-value
<b><i>On-Standardized-Med Trial</i></b>				
TRIO (ANCOVA on the ranks)				0.0223
TRIO (median)	-8.0 (69)	-3.0 (67)		
TRIO (Hodges-Lehmann estimate)			-4.5	

Data presented as median (n); p-value via a baseline adjusted ANCOVA on the ranks (two-sided).

<sup>1</sup> Negative value favors uRDN.

# Primary and Secondary Endpoints:


## Mean 24-hour SBP, Daytime SBP, and Office SBP at 2 Months

Note not all p-values are adjusted for multiplicity.

**SOLO ITT**

**Off-Med Studies**


**RADIANCE-II ITT**



# Primary and Secondary Endpoints: Mean 24-hour SBP, Daytime SBP, and Office SBP at 2 Months

**On Standardized Med Study**


**TRIO ITT**



BL Adj. Mean Diff.	<b>-4.3</b>		<b>-4.5</b>		<b>-5.4</b>	
95% CI	<b>(-9.3, 0.7)</b>		<b>(-9.6, 0.6)</b>		<b>(-11.9, 1.1)</b>	
p-value	<b>0.0895</b>		<b>0.0809</b>		<b>0.1042</b>	
N =	69	67	69	67	64	66
BL SBP	143.9	145.4	150.0	151.1	155.6	154.9

*Note not all p-values are adjusted for multiplicity.*

# Proportion of Subjects with $\geq 5$ , $\geq 10$ , & $\geq 15$ mmHg Reduction in Daytime SBP at 2 Months





Note not all p-values are adjusted for multiplicity.

# Off-Med Studies

ASBP  $\leq 135$  at 2 Months as a Function of Baseline ASBP

ASBP Change at 2 Months as a Function of Baseline ASBP






# On Standardized Medication Study


**ASBP  $\leq 135$  at 2 Months as a Function of Baseline ASBP**

**ASBP Change at 2 Months as a Function of Baseline ASBP**




# Off-Med Studies Subgroup Analysis at 2 Months

**SOLO ITT**



Change from baseline daytime ASBP (mmHg)  
Mean difference (uRDN – Sham)

**RADIANCE-II ITT**



Change from baseline daytime ASBP (mmHg)  
Mean difference (uRDN – Sham)


Note not all p-values are adjusted for multiplicity.

Results generally consistent across subgroups in SOLO and RADIANCE-II

# On Standardized Medication Study Subgroups 2 Months



**TRIO ITT**





*Note not all p-values are adjusted for multiplicity.*

**Treatment difference in favor of uRDN more pronounced in US subjects and subjects with baseline ASBP ≥ the median**

# **DURABILITY OF BP REDUCTION**

# Daytime and Nighttime ASBP at Baseline, 2 Months, and 6 Months

## RADIANCE-II (OFF-Med Study)



**BP Meds:**





- Thru 2 months: Added for emergency only thru 2 months
- 2-6 months: Added as needed to optimize BP

**Patients blinded to treatment group thru 6 months**


RADIANCE II  
Mild-Moderate Hypertension

# Off-Med Studies



BP  $\Delta$  from Baseline to Latest Available Follow-up



# On Standardized Med Study



## BP $\Delta$ from Baseline to Latest Available Follow-up



# Medication Burden

## Assessments


- Number of prescribed antihypertensive medications
- Antihypertensive Medication Load Index: Composite based on the dosages of medications

$$\sum_{AH\ Meds} \left( class\ weight \frac{prescribed\ dose}{standard\ dose} \right)$$



# Medication Burden – SOLO

		Change in Daytime ASBP (mmHg)	Change in office SBP (mmHg)	Average number of BP Meds	Med Load Index
2 months	uRDN (N = 74)	-8.5 ± 9.3	-10.8 ± 13	0.1 ± 0.3	0.0 ± 0.1
	Sham (N = 72)	-2.2 ± 10.0	-3.9 ± 17.4	0.2 ± 0.5	0.1 ± 0.3
6 months	uRDN (N = 69)	-18.1 ± 12.2	-18.2 ± 14.2	1.0 ± 0.9	0.5 ± 0.4
	Sham (N=71)	-15.6 ± 13.2	-15.9 ± 17.2	1.3 ± 0.9	0.7 ± 0.5
12 months	uRDN (N = 65)	-16.5 ± 12.9	-18.1 ± 14.9	1.0 ± 0.9	0.5 ± 0.5
	Sham (N = 67)	-15.8 ± 13.1	-13.6 ± 17.2	1.3 ± 0.9	0.7 ± 0.5




# Medication Burden Off-Med Study



## RADIANCE-II


		Daytime ASBP Δ (mmHg)	Office SBP Δ (mmHg)	Average # of BP Meds	Med Load Index
2 months	uRDN	-7.9 ± 11.6 (145)	-11.0 ± 13.5 (137)	0.1 ± 0.3	0.1 ± 0.2
	Sham	-1.8 ± 9.5 (73)	-5.5 ± 12.9 (71)	0.1 ± 0.4	0.1 ± 0.2
6 months	uRDN	-17.5 ± 11.4 (143)	-20.9 ± 14.8 (143)	1.3 ± 1.0	0.7 ± 0.6
	Sham	-17.4 ± 14.0 (67)	-20.2 ± 16.4 (57)	1.5 ± 1.0	0.8 ± 0.6



Note not all p-values are adjusted for multiplicity.

# Medication Burden – TRIO

		Change in Daytime ASBP (mmHg)	Change in office SBP (mmHg)	Average number of BP Meds	Med Load Index
2 months	uRDN	-9.0 ± 14.5 (69)	-8.5 ± 19.1 (64)	3.1 ± 0.5	2.1 ± 0.4
	Sham	-4.8 ± 15.9 (67)	-2.8 ± 20.7 (66)	3.1 ± 0.8	2.1 ± 0.5
6 months	uRDN	-11.8 ± 14.2 (65)	-10.4 ± 16.8 (63)	3.8 ± 1.0	2.3 ± 0.6
	Sham	-12.3 ± 14.2 (64)	-11.2 ± 22.7 (64)	4.1 ± 1.1	2.4 ± 0.6
12 months	uRDN	-12.1 ± 14.1 (59)	-12.6 ± 19.8 (59)	3.7 ± 1.5	2.3 ± 0.9
	Sham	-10.9 ± 18.3 (59)	-7.8 ± 28.9 (59)	4.0 ± 1.2	2.5 ± 0.7



# Blinding Assessment

- Blinding generally maintained across 3 trials using Bang and James Blinding Indices
- In off-med trials, more sham subjects able to guess correct treatment at 6 months vs baseline.



# Patient Preference Study

**David Gebben, PhD**

**Health Economist**

**Office of Strategic Partnerships and  
Technology Innovation**

# Patient Preference Information (PPI)

**CDRH Guidance Document:** *Patient Preference Information – Voluntary Submission, Review in PMAs, HDE Applications, and De Novo Requests and Inclusion in Decision Summaries and Device Labeling. August 2016*

PPI Definition: Qualitative or quantitative assessments of the relative desirability or acceptability to patients of specified alternatives or choices among outcomes or other attributes that differ among alternative health interventions

- Not a patient-reported outcome (PRO) or other clinical trial endpoint or outcome

# Benefit-Risk Determination

- Before considering benefit-risk (B-R), establish a reasonable assurance of safety and effectiveness
- CDRH recognizes the patient preference information can supplement the assessment of benefits and risks
- Patient preference studies consider how patients weigh the benefits and risks of treatment options\*

*\*FDA Guidance: Factors to Consider When Making Benefit-Risk Determinations in Medical Device Premarket Approval and De Novo Classifications Guidance for Industry and Food and Drug Administration Staff (Issued August 30, 2019)*

# Recommended Qualities of Patient Preference Studies\*



Well-designed and conducted patient preference studies can provide valid scientific evidence regarding patients' risk tolerance and perspective on benefit. This may inform FDA's evaluation of a device's benefit-risk profile during the PMA, HDE application, and De Novo request review processes.

## A. All about Patients


- Patient Centeredness
- Sample Representativeness
- Capturing Heterogeneous Patient Preferences
- Comprehension by Study Participants

## B. Good Study Design

- Established Good Research Practices
- Effective Benefit-Risk Communication
- Minimal Cognitive Bias
- Relevance

## C. Good Study Conduct and Analysis

- Study Conduct
- Logical Soundness
- Robustness of Analysis of Results




\*FDA Guidance: Patient Preference Information - Voluntary Submission, Review in Premarket Approval Applications, Humanitarian Device Exemption Applications, and De Novo Requests, and Inclusion in Decision Summaries and Device Labeling (Aug 2016)



# ReCor PPI Study Qualities


- 258 respondents to PPI survey
- Qualities consistent with CDRH PPI Guidance:
  - Follows guidelines for good research practices established by recognized professional organizations
  - Followed good ethical research practices
  - Survey understandable to respondents
- A few attribute levels do not correspond to levels supported by the evidence - may have tilted the patient preference study results toward the uRDN procedure.



Example Choice Task


# PPI Study Results by Attributes and Levels

- Results generally as expected with levels in order of expected preference
- Decrease in “10 years CV risk” was the main driver of preference choices



# PPI Study Results

- 42%\* of patients would choose uRDN procedure versus taking an additional pill all else being equal
- May be an overestimate due to levels
  - “Risk of mild-to-moderate side effects requiring more doctor visits” - Lowest risk presented as 20% - the highest risk for pills 10% based on the literature
  - “Treatment durability” - durability of the uRDN procedure 7 years - longer than currently available clinical data



\*Recor Executive Summary page 43.



# Proposed Post-Approval Study and FDA Conclusions

**Paul Warren, PhD**  
**Biomedical Engineer**  
**Office of Cardiovascular Devices**

# Proposed Post-Approval Study

- Continued follow-up of SOLO and TRIO subjects through 3 years and RADIANCE-II subjects through 5 years
- Global Paradise System Registry (US-GPS) – new enrollment of 500 US subjects at up to 100 sites per the approved indication for use

# Proposed Post-Approval Study Endpoints

## Safety

### *30 days*

- All-cause mortality
- Major vascular complications
- Hospitalization for hypertensive crisis
- Hospitalization for major cardiovascular- or hemodynamic-related events
- Renal artery injury requiring invasive intervention

### *6-month, 12-month, and annually to 5 years*

- All-cause mortality
- New onset renal artery stenosis >70%
- Significant decline in renal function
- Hospitalization for hypertensive crisis
- Hospitalization for major cardiovascular- or hemodynamic-related events

## Effectiveness

- Primary
  - Change in home SBP & DPB at 3 months
  - Change in office SBP & DBP at 3 months
- Secondary
  - Home & office BP, heart rate, and pulse pressure through 5 years
  - Number of BP meds
  - Quality of Life
  - Percentage of subjects with controlled office and home BP

# Proposed Post-Approval Subgroup Analysis

- *Sex*
- *Race (Black versus non-black)*
- *Age*
- *Baseline Office systolic BP*
- *Baseline Home systolic BP*
- *Heart Rate*
- *Abdominal obesity*

- *Body Mass Index*
- *eGFR*
- *Heart Failure (NYHA I, II, III, IV)*
- *Isolated Systolic Hypertension*
- *Diabetes*
- *Number and class of antihypertension medications*

# Conclusions (1)

- Primary safety endpoint met
  - In RADIANCE-II, no MAEs through 30 days and no cases of > 70% RAS through 6 months → primary safety event rate of 0%
  - Pooled safety event rate of 1.1%
- Primary effectiveness endpoint met for Off-med SOLO and RADIANCE-II trials at 2 months
  - Between-group difference in mean daytime ASBP reduction 6.3 mmHg for SOLO and RADIANCE-II in favor of uRDN
- In the On Standardized Meds TRIO trial, the between-group difference in mean daytime ASBP reduction was 4.5 mmHg in favor of uRDN at 2 months



# Conclusions (2)

- Strengths
  - Three sham-controlled randomized trials
  - Studies independently powered for effectiveness
- Limitations
  - Small long-term RCT data sample size for uRDN patients
    - SOLO: 51 SOLO subjects at 3 years
    - TRIO: 51 TRIO subjects at 2 years
  - Challenges in interpreting the clinical significance of BP reduction durability
    - BP medication changes beyond 2 months
    - Subject blinding
    - Cross-over reduced the sample size of the control group
- Patient preference study
  - Some patients may prefer uRDN to an additional BP pill

# PANEL QUESTIONS

# QUESTION 1

# #1. Safety Profile

## Primary Safety Endpoint Results

- Composite of 30-day MAE & 6-month new onset RAS >70% DS
- RADIANCE-II: **0%** (0/150)
  - Upper 95% CI 1.63%, met the 9.8% pre-specified performance goal
- Pooled: **1.1%** (6/367; Upper 95% CI 2.75%)
  - 2 deaths
  - 2 major vascular complications (1 pseudoaneurysm, 1 DVT)
  - 1 hypotensive crisis
  - 1 hospitalization for presyncope

# #1. Safety Profile

## Renal Artery Stenosis

- Evaluable CTA/MRA imaging at 12 months (pooled between SOLO, TRIO, R-II)
- No hemodynamically significant RAS
  - 0.0% (0/238) of >70% DS
  - 0.8% (2/238) of 51-70% DS
  - 2.1% (5/238) of 31-50% DS
  - 1.3% (3/238) of 1-30% DS
- No clinically significant changes in eGFR or serum creatinine

# #1. Safety Profile

- Safety events were low and generally transient in nature
- No significant RAS was observed. Although mild to moderate renal artery narrowing is not associated with a functional reduction in renal blood flow, long-term follow-up data are limited, and renal arterial lesions may progress over time.

**Please discuss the acute and mid-term procedural and device safety profile of uRDN and the clinical significance of renal arterial responses to uRDN treatment.**

# QUESTION 2

## #2. Effectiveness – BP Measurement Method

- Data have been presented using both ambulatory blood pressure measurement (ABPM) and office blood pressure measurement (OBPM)
- Most prior hypertension trials have used OBPM
- ABPM has been shown to have greater prognostic value and was identified as preferable at the 2018 Panel Meeting. This may be due to the large number of measurements made for ABPM that are free from potential biases (e.g., white coat effect)






# #2. Effectiveness – BP Measurement Method

### SOLO (Off-medication)\*

### RADIANCE-II (Off-medication)\*

### TRIO (On standardized medication)\*



BL Adj. Mean Diff.	<b>-4.1</b>		<b>-6.3</b>		<b>-6.5</b>	
95% CI	<b>(-7.1, -1.2)</b>		<b>(-9.4, -3.1)</b>		<b>(-11.3, -1.8)</b>	
p-value	<b>0.0061</b>		<b>0.0001</b>		<b>0.0073</b>	
N =	74	72	74	72	74	72
BL SBP	142.6	143.8	150.3	150.0	154.5	153.6


BL Adj. Mean Diff.	<b>-6.2</b>		<b>-6.3</b>		<b>-5.4</b>	
95% CI	<b>(-9.1, -3.4)</b>		<b>(-9.3, -3.2)</b>		<b>(-9.0, -1.8)</b>	
p-value	<b>&lt;.0001</b>		<b>&lt;.0001</b>		<b>0.0035</b>	
N =	144	72	145	73	137	71
BL SBP	143.2	144.5	150.2	151.3	156.8	156.7

BL Adj. Mean Diff.	<b>-4.3</b>		<b>-4.5</b>		<b>-5.4</b>	
95% CI	<b>(-9.3, 0.7)</b>		<b>(-9.6, 0.6)</b>		<b>(-11.9, 1.1)</b>	
p-value	<b>0.0895</b>		<b>0.0809</b>		<b>0.1042</b>	
N =	69	67	69	67	64	66
BL SBP	143.9	145.4	150.0	151.1	155.6	154.9

Note not all p-values are adjusted for multiplicity.

\*At 2 months

# #2. Effectiveness – BP Measurement Method



*Note not all p-values are adjusted for multiplicity.*

**Off-medication study; at 2 months**

Please discuss the relative value of ABPM and OBPM in assessing changes in blood pressure for purposes of evaluating effectiveness of uRDN.

# QUESTION 3

### #3. Effectiveness – Magnitude of BP Reduction



- 2018 Panel recommendation for a clinically significant therapeutic effect:
  - Difference in SBP reduction between treatment and sham groups should be 5-7 mmHg (ASBP)
- Primary endpoints for the Paradise uRDN System
  - Adjusted difference in mean reduction of daytime ASBP at 2 months favored uRDN over Sham by 4.5 mmHg (on-standardized med) to 6.3 mmHg (off-med)

# #3. Effectiveness – Magnitude of BP Reduction

## Off-Med Studies

### SOLO ITT

### RADIANCE-II ITT



BL Adj. Mean Diff. **-4.1** **-6.3** **-6.5**  
 95% CI **(-7.1, -1.2)** **(-9.4, -3.1)** **(-11.3, -1.8)**  
 p-value **0.0061** **0.0001** **0.0073**

BL Adj. Mean Diff. **-6.2** **-6.3** **-5.4**  
 95% CI **(-9.1, -3.4)** **(-9.3, -3.2)** **(-9.0, -1.8)**  
 p-value **<.0001** **<.0001** **0.0035**

N =	74	72	74	72	74	72
BL SBP	142.6	143.8	150.3	150.0	154.5	153.6


N =	144	72	145	73	137	71
BL SBP	143.2	144.5	150.2	151.3	156.8	156.7

Note not all p-values are adjusted for multiplicity.

# #3. Effectiveness – Magnitude of BP Reduction

## On Standardized Med Study

TRIO ITT



BL Adj. Mean Diff.	<b>-4.3</b>		<b>-4.5</b>		<b>-5.4</b>	
95% CI	<b>(-9.3, 0.7)</b>		<b>(-9.6, 0.6)</b>		<b>(-11.9, 1.1)</b>	
p-value	<b>0.0895</b>		<b>0.0809</b>		<b>0.1042</b>	
N =	69	67	69	67	64	66
BL SBP	143.9	145.4	150.0	151.1	155.6	154.9

Note not all p-values are adjusted for multiplicity.

# #3. Effectiveness – Magnitude of BP Reduction






Please discuss the relative importance of absolute BP reduction in uRDN subjects vs. the difference in BP reduction between uRDN and Sham groups in evaluating the treatment effect for SOLO, TRIO, and RADIANCE-II.

# QUESTION 4



# #4. Effectiveness – Durability of BP Reduction



- 2 months – differences between groups statistically significant
- After 2 months – absolute reductions persisted in both groups
- 6 months and later - differences in mean daytime ASBP reduction not significant

# #4. Effectiveness – Durability of BP Reduction

Challenges in interpreting BP data beyond 6 months:

- **Medication** – From 2-6 months, medication escalation to target BP for all studies (med index only significantly different for SOLO at 6 and 12 mo)
- **Unblinding** – occurred at 6-mo for SOLO and TRIO; 12-mo for RADIANCE-II
- **Crossover** reduced Sham sample size for later timepoints – allowed at 6-mo for SOLO and TRIO; 12-mo for RADIANCE-II
- **Limited long-term data** – RADIANCE-II has limited data beyond 6-mo, OBP available for most patients at 24-mo for SOLO (68%) and TRIO (63%)

## #4. Effectiveness – Durability of BP Reduction

Please discuss the strengths and limitations of longer-term BP data in patients treated with uRDN and what conclusions can be drawn about whether uRDN provides a durable reduction in BP.

# QUESTION 5

# #5 Patient Preference Study

- ReCor conducted a patient preference study with 258 patients to ascertain preferences for uRDN procedure compared to pills only.
- Generally, the study aligned with the CDRH PPI Guidance document.
- Based on the preference weights 42% of respondents would choose the uRDN procedure over an additional pill.
- Two attribute levels do not correspond to the evidence which may have impacted the respondents' choices.

# #5 Patient Preference Study

Please discuss the degree of importance that the patient preference study results should be given when considering supplemental benefit-risk assessment information.

# QUESTION 6

# #6. Indications for Use

	<b>SOLO (off-med)</b>	<b>TRIO (standardized med)</b>	<b>RADIANCE-II (off-med)</b>
Patient Population	Mild-to-Moderate HTN	Resistant HTN	Stage 2 HTN
Sample size	146 (69 uRDN: 37 Sham)	136 (74 uRDN: 72 Sham)	224 (150 uRDN: 74 Sham)
OBP	<ul style="list-style-type: none"> <li>• <math>\geq 140/90</math> &amp; <math>&lt; 180/110</math> mmHg on 0, 1 or 2 meds;</li> <li><i>or</i></li> <li>• <math>\leq 140/90</math> on 1 or 2 meds</li> </ul>	<ul style="list-style-type: none"> <li>• <math>\geq 140/90</math> on <math>\geq 3</math> meds, including a diuretic</li> </ul>	<ul style="list-style-type: none"> <li>• <math>\geq 140/90</math> &amp; <math>&lt; 180/120</math> mmHg on 0, 1, or 2 meds;</li> <li><i>and</i></li> <li>• Previous or currently prescribed AH therapy</li> </ul>
Daytime ABP	$\geq 135/85$ & $< 170/105$ mmHg after washout	$\geq 135/85$ mmHg after stabilization	$\geq 135/85$ & $< 170/105$ mmHg after washout
Antihypertensive Medication	0, 1, or 2	At least 3	0, 1, or 2



## #6. Indications for Use

The proposed indications for use are:

*The Paradise uRDN System is indicated to reduce blood pressure in patients with uncontrolled hypertension, who may be inadequately responsive to, or who are intolerant to anti-hypertensive medications.*

- a. Please discuss whether the available clinical data support the proposed indications for use.**
- b. Please discuss if “inadequately responsive to or intolerance to anti-hypertensive medications” should be further defined in the labeling, and if so, please discuss definitions.**

# QUESTION 7

# #7. Labeling

Please discuss labeling recommendations for post-uRDN renal artery imaging (e.g., imaging modality, follow-up imaging timing and frequency, and site training and expertise).

# QUESTION 8

## #8. Benefit/Risk

Given the totality of the evidence presented regarding the safety and effectiveness of the device, please comment on the benefit-risk profile of this device.

# QUESTION 9

# #9. Post-market Study

ReCor proposed a postmarket registry study that will incorporate uRDN subjects from RADIANCE-II and the continued access study with enrollment of up to 500 new subjects that meet the indications for use (uncontrolled hypertension, who may be inadequately responsive to, or who are intolerant to anti-hypertensive medications). This proposed study will collect office and home measured BP (and not 24-hour ABP).

- a. Please comment on the sample size, proposed endpoints, and BP measurement methods.
- b. Please discuss whether the PAS enrollment should pre-specify more diverse patient subgroups.
- c. Please discuss the strengths and limitations of a single arm study design for the PAS.
- d. No renal arterial imaging follow-up is planned. Therefore, please discuss the need for a pre-specified imaging follow-up protocol to confirm long-term uRDN safety.







# Supplemental Slides

# Study Design Elements

- Third Party Review
  - DSMB
  - CEC
  - Core labs
- Blinding
  - Patients blinded until 6 months (SOLO and TRIO) or 12 months (R-II)
  - Study personnel measuring BP blinded for study duration
  - Sponsor and DSMB blinded to primary effectiveness data

# Follow-up Schedule

	Screening	Baseline	Procedure	1M	2M	6M	12M	24-60M
<b>OBPM</b>	X	X		X	X	X	X	X
<b>HBPM</b>		X		X	X	X	X	
<b>ABPM</b>		X			X	X	X	
<b>Renal DUS<sup>1</sup></b>		X <sup>2</sup>			X	X		X <sup>5</sup>
<b>CTA/MRA</b>			X		X <sup>3</sup>	X <sup>3,4</sup>	X <sup>5</sup>	
<b>Urine chemistry and drug metabolite</b>		X			X	X		
<b>Blood chemistry</b>		X			X	X	X	
<b>Quality of Life</b>	X	X			X	X	X	
<b>Blinding assessment</b>			X (discharge)		X	X		

OBPM/HBPM/ABPM: Office/home/ambulatory blood pressure measurement; CTA: computed tomography angiography; MRA: magnetic resonance angiography

<sup>1</sup> SOLO/TRIO only

<sup>2</sup> Recommended. A recent (within 6 months of consent) good quality renal duplex ultrasound is acceptable

<sup>3</sup> if required in the event of clinical suspicion of renal artery stenosis (RAS)

<sup>4</sup> Required for all RII subjects (Sham and uRDN)

<sup>5</sup> procedure was conducted on uRDN treated subjects

## Key Baseline Patient Features

Measure	SOLO (Off-Med)		RADIANCE-II (Off-Med)		TRIO (On-Standardized Med)	
	uRDN (n=74)	Sham (n=72)	uRDN (n=150)	Sham (n=74)	uRDN (n=69)	Sham (n=67)
<b>Sex</b>						
<b>Male</b>	46 (62.2%)	39 (54.2%)	103 (68.7%)	57 (77.0%)	56 (81.2%)	53 (79.1%)
<b>Female</b>	28 (37.8%)	33 (45.8%)	47 (31.3%)	17 (23.0%)	13 (18.8%)	14 (20.9%)
<b>Age</b>	54.4 ± 10.2	53.8 ± 10.0	55.1 ± 9.9	54.9 ± 7.9	52.3 ± 7.5	52.8 ± 9.1
<b>Geography</b>						
<b>US</b>	35 (47.3%)	34 (47.2%)	100 (66.7%)	46 (62.2%)	28 (40.6%)	25 (37.3%)
<b>OUS</b>	39 (52.7%)	38 (52.8%)	50 (33.3%)	28 (37.8%)	41 (59.4%)	42 (62.7%)
<b>Race</b>						
<b>American Indian or Alaska Native</b>	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.00%)	1 (1.52%)
<b>Asian</b>	1 (1.3%)	0 (0.0%)	0 (0.0%)	1 (1.3%)	1 (1.5%)	1 (1.5%)
<b>Black</b>	12 (16.2%)	13 (18.0%)	21 (14.0%)	15 (20.2%)	14 (20.6%)	13 (19.7%)
<b>Caucasian</b>	60 (81.0%)	52 (72.2%)	114 (76.0%)	56 (75.6%)	45 (66.2%)	51 (77.3%)
<b>Hispanic or Latino</b>	1 (1.3%)	4 (5.5%)	15 (10.0%)	2 (2.7%)	5 (7.4%)	0 (0.0%)
<b>Native Hawaiian or other Pacific Islander</b>	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)
<b>Other/Mixed Race</b>	0 (0.0%)	0 (0.0%)	15 (10.00%)	2 (2.7%)	3 (4.41%)	0 (0.0%)
<b>BMI</b>	29.9 ± 5.9	29.0 ± 5.0	30.1 ± 5.2	30.6 ± 5.2	32.8 ± 5.7	32.6 ± 5.4
<b>Abdominal circumference (cm)</b>	101.5 ± 14.2	98.5 ± 15.1	102.4 ± 12.3	104.3 ± 13.1	109.4 ± 15.5	109.2 ± 12.9
<b>Office SBP (mmHg)</b>	142.6 ± 14.7	144.6 ± 15.9	155.8 ± 11.1	154.3 ± 10.6	161.9 ± 15.5	163.6 ± 16.8
<b>Office DBP (mmHg)</b>	92.3 ± 10.1	93.6 ± 8.3	101.3 ± 6.7	99.1 ± 5.6	105.1 ± 11.6	103.3 ± 12.7

# Procedural Characteristics

Renal Denervation Group	SOLO Study	TRIO Study	RADIANCE-II
	N=74	N = 69	n=150
Procedure time (sheath removal - sheath insertion) (min) <sup>1</sup>	71.9 ± 23.2	83.0	76.7 ± 25.2
Contrast volume (cc)	138.5 ± 66.6	176.9 ± 77.0	135.7 ± 67.4
Fluoroscopy exposure (minutes)	13.7 ± 6.8	19.0 ± 11.5	15.9 ± 8.6
Total Number of Emissions <sup>2</sup>	5.3 ± 1.1	5.8 ± 1.2	5.6 ± 1.0
Number of Subjects with Accessory and/or Proximal Side Branch Emissions	9 / 74 (12.16%)	17 / 69 (24.64%)	30 / 150 (20.00%)
Treatment successfully delivered (minimum 2 emissions bilateral)	71 / 74 (95.95%)	67 / 69 (97.10%)	148 / 150 (98.67%)
Total Emission Time (seconds)	37.4 ± 8.0	40.7 ± 8.1	38.9 ± 7.3

# Primary Safety Endpoint

## *RADIANCE-II*



- The primary safety endpoint was defined as the occurrence of at least one of following major adverse events (MAE):
  - a. 30-day
    - 1) All-cause mortality
    - 2) New onset (acute) end-stage renal disease (eGFR<15 mL/min/m<sup>2</sup> or need for renal replacement therapy)
    - 3) Significant embolic event resulting in end-organ damage (e.g., kidney or bowel infarct, lower extremity ulceration or gangrene, or doubling of serum creatinine)
    - 4) Renal artery perforation requiring invasive intervention
    - 5) Renal artery dissection requiring an invasive intervention
    - 6) Major vascular complications (e.g., clinically significant groin hematoma, arteriovenous fistula, pseudoaneurysm) requiring surgical repair, interventional procedure, thrombin injection, or blood transfusion (>2 units of packed red blood cells within any 24-hr period during the first 7 days post-randomization)
    - 7) Hospitalization for hypertensive or hypotensive crisis
    - 8) Hospitalization for major cardiovascular- or hemodynamic-related events (e.g., HF; MI; stroke)
    - 9) New stroke
    - 10) New MI
  - b. 6 Month: New onset renal artery stenosis (RAS), defined as a >70% stenosis, confirmed by CTA/MRA

**Additional analysis for pooled SOLO, TRIO, & RADIANCE-II trials**

# Primary Safety Endpoint – R-II Imaging Completed



**Table X: RADIANCE-II Imaging Completed  
(ITT subjects in window)**

	uRDN (n=150)	Sham (n=74)
<b>6 Month follow-up</b>		
CTA/MRA	94.5% (138/146)	85.5% (59/69)
CTA	76.7% (112/146)	65.2% (45/69)
MRA	17.8% (26/146)	20.3% (14/69)
<b>12 Month follow-up</b>		
CTA/MRA	94.4% (117/124)	--
CTA	80.6% (100/124)	--
MRA	13.7% (17/124)	--

*Data displayed as % (n/N)*

*Data based on available imaging forms; 6 and 12 month data are ongoing.*

*NOTE: Sham subjects are not required to have imaging at 12 months under the current protocol.*



# Safety - RAS data through 12 months

Study	Total	No measurable stenosis % (n)	1-30% stenosis % (n)	31-50% stenosis % (n)	51-70% stenosis % (n)	71-99% stenosis % (n)	Renal artery occlusion % (n)
RADAIANCE-II	90	94.4% (85)	3.3% (3)	2.2% (2)	0.0% (0)	0.0% (0)	0.0% (0)
SOLO*	64	96.9% (62)	0.0% (0)	1.6% (1)	1.6% (1)	0.0% (0)	0.0% (0)
TRIO	53	100.0% (53)	0.0% (0)	0.0% (0)	0.0% (0)	0.0% (0)	0.0% (0)
RADIANCE-II CO	6	83.3% (5)	0.0% (0)	16.7% (1)	0.0% (0)	0.0% (0)	0.0% (0)
SOLO CO	25	92.0% (23)	0.0% (0)	4.0% (1)	4.0% (1)	0.0% (0)	0.0% (0)
Combined	238	95.8% (228)	1.3% (3)	2.1% (5)	0.8% (2)	0.0% (0)	0.0% (0)

*Data included in listing is based on data available from the Core Lab*

*Not all Crossover (CO) subjects have reached 12M CO f/u.*

*\* One SOLO subject randomized to treatment was not treated, and subsequently crossed over.*

*Their data is included in the SOLO CO summaries.*


# Pooled Renal Artery Diameter Stenosis (uRDN Subjects at 12 Months by CTA or MRA)



Study	N	Renal Artery Diameter Stenosis					
		0%	1-30%	31-50%	51-70%	71-99%	Occluded
SOLO	64	96.9% (62)	0.0% (0)	1.6% (1)	1.6% (1)	0.0% (0)	0.0% (0)
RADIANCE-II	112	92.9% (104)	2.7% (3)	2.7% (3)	1.8% (2)	0.0% (0)	0.0% (0)
TRIO	53	100.0% (53)	0.0% (0)	0.0% (0)	0.0% (0)	0.0% (0)	0.0% (0)
SOLO CO	25	92.0% (23)	0.0% (0)	4.0% (1)	4.0% (1)	0.0% (0)	0.0% (0)
RADIANCE-II CO	6	83.3% (5)	0.0% (0)	16.7% (1)	0.0% (0)	0.0% (0)	0.0% (0)
<b>Total</b>	<b>238</b>	<b>95.8% (228)</b>	<b>1.3% (3)</b>	<b>2.1% (5)</b>	<b>0.8% (2)</b>	<b>0.0% (0)</b>	<b>0.0% (0)</b>

CO = Crossover Subjects

# Means and Medians for TRIO



Note not all p-values are adjusted for multiplicity; p-value is based on ANCOVA; p-value(Rank) is based on ANCOVA on the ranks.

Figure X: SBP at 2 months for TRIO (means and medians)

# Primary Safety Results: Imaging Performed (SOLO, TRIO)

	SOLO		TRIO	
	uRDN (N=74)	Sham (N=72)	uRDN (N=69)	Sham (N=67)
2 Month Follow-up				
Ultrasounds completed	97.3% (72/74)	90.3% (65/72)	93.9% (62/66)	91.0% (61/67)
Ultrasounds not done	2.7% (2/74)	9.7% (7/72)	6.1% (4/66)	9.0% (6/67)
Ultrasounds trigger MRA/CTA	10.8% (8/74)	8.3% (6/72)	3.0% (2/66)	9.0% (6/67)
Of triggered, MRA/CTA completed	50.0% (4/8)	50.0% (3/6)	50.0% (1/2)	66.7% (4/6)
MRA/CTA completed, but not due to Ultrasound trigger	0.0% (0/74)	0.0% (0/72)	0.0% (0/66)	0.0% (0/67)
6 Month Follow-up				
Ultrasounds completed	98.6% (72/73)	95.8% (68/71)	92.4% (61/66)	95.3% (61/64)
Ultrasounds not done	1.4% (1/73)	4.2% (3/71)	7.6% (5/66)	4.7% (3/64)
Ultrasounds trigger MRA/CTA	11.0% (8/73)	8.5% (6/71)	4.5% (3/66)	9.4% (6/64)
Of triggered, MRA/CTA completed	75.0% (6/8)	16.7% (1/6)	66.7% (2/3)	50.0% (3/6)
MRA/CTA completed, but not due to Ultrasound trigger	0.0% (0/73)	0.0% (0/71)	1.5% (1/66)	0.0% (0/64)
24 Month Follow-up				
Ultrasounds completed	91.7% (55/60)	--	88.9% (48/54)	--
Ultrasounds not done	8.3% (5/60)	--	11.1% (6/54)	--
Ultrasounds trigger MRA/CTA	13.3% (8/60)	--	11.1% (6/54)	--
Of triggered, MRA/CTA completed	50.0% (4/8)	--	50.0% (3/6)	--
MRA/CTA completed, but not due to Ultrasound trigger	1.7% (1/60)	--	1.9% (1/54)	--

# Secondary & Observational BP Changes @ 2 months



	SOLO		TRIO		RADIANCE-II	
	Mean Difference (95% CI) (uRDN - Sham) <sup>1</sup>	p-value	Mean Difference (95% CI) (uRDN - Sham) <sup>1</sup>	p-value	Mean Difference (95% CI) (uRDN - Sham) <sup>1</sup>	p-value
Daytime Ambulatory DBP (mmHg)	-2.6 (-4.6, -0.6)	0.0118 (0.0060*)	-1.6 (-4.9, 1.7)	0.3415 (0.1835*)	-3.9 [-5.6, -2.2]	<.0001
24 Hour Ambulatory SBP (mmHg)	-4.1 (-7.1, -1.2)	0.0061	-4.3 (-9.3, 0.7)	0.0895 (0.0162*)	-6.2 [-9.1, -3.4]	<.0001
24 Hour Ambulatory DBP (mmHg)	-1.8 (-3.7, 0.2)	0.0715	-1.7 (-4.9, 1.5)	0.3054 (0.1228*)	-4.1 [-5.7, -2.4]	<.0001
Nighttime Ambulatory SBP** (mmHg)	-2.5 (-6.0, 0.9)	0.1534	-4.4 (-9.9, 1.2)	0.1213 (0.0441*)	-5.8 [-9.0, -2.6]	0.0004 (<.0001*)
Nighttime Ambulatory DBP** (mmHg)	-1.4 (-3.8, 1.0)	0.2492	-2.2 (-5.8, 1.4)	0.2242 (0.0534*)	-4.2 [-6.3, -2.2]	<.0001 (<.0001*)
Home SBP** (mmHg)	-7.1 (-10.4, -3.8)	<.0001 (<.0001*)	-4.3 (-8.6, 0.0)	0.0524	-7.6 [-10.1, -5.0]	<.0001
Home DBP** (mmHg)	-3.6 (-5.6, -1.5)	0.0009 (<.0001*)	-2.6 (-5.2, 0.0)	0.0527	-4.3 [-5.9, -2.8]	<.0001
Office SBP** (mmHg)	-6.5 (-11.3, -1.8)	0.0073 (0.0007*)	-5.4 (-11.9, 1.1)	0.1042 (0.0374*)	-5.4 [-9.0, -1.8]	0.0035
Office DBP** (mmHg)	-4.1 (-7.0, -1.3)	0.0045	-3.2 (-7.5, 1.1)	0.1375 (0.1598*)	-2.3 [-4.9, 0.2]	0.0755

Note not all p-values are adjusted for multiplicity.

# Subgroup analysis of primary effectiveness endpoint at 2 months – US vs OUS



	US		OUS		Interaction p-value
	Mean Difference (95% CI) (uRDN - Sham)	p-value btw uRDN and Sham	Mean Difference (95% CI) (uRDN - Sham)	p-value btw uRDN and Sham	
<b>SOLO</b>	-8.4 (-13.1, -3.7)	0.0006	-4.2 (-8.5, 0.1)	0.0543	0.1905 (0.1796*)
<b>TRIO</b>	-7.7 (-15.7, 0.3) -9.0 (-16.1, -2.3) <sup>1</sup>	0.0593 (0.0048*)	-2.0 (-8.8, 4.7) -1.5 (-6.1, 3.0) <sup>1</sup>	0.5477 (0.5388*)	0.2901 (0.0846*)
<b>RADIANCE-II</b>	-4.7 (-8.6, -0.9)	0.0172	-9.0 (-13.8, -4.1)	0.0005	0.1477 (0.1500*)

Note not all p-values are adjusted for multiplicity.

# Medication Burden

- Medication Burden also assessed via Defined Daily Dose (DDD)
- DDD = sum of the average maintenance dose per day the subject is taking

# Statistical Backup Slides



# Statistical Analysis Populations

- **Intention-to-Treat (ITT)** cohort: subjects according to their randomization assignment
- **Per-Protocol (PP)** cohort: subjects treated per their assigned treatment group without deviation from major enrollment criteria
- **Complete ABPM (CA)** cohort: subjects treated per their assigned treatment group that have ABP values at both baseline and follow-up
- **Crossover (CO)** cohort: subjects who received uRDN after being randomized to Sham.
  - Crossover allowed:
    - After 6-months follow-up in SOLO and TRIO
    - After 12-months follow-up in RADIANCE-II

# Sample Size and Power

- Primary safety endpoint (at 6 months):
  - RADIANCE-II only:
    - The evaluable sample size of 128 uRDN subjects provides about 95% power for the performance goal of 9.8% if the composite MAE rate is expected at 3.0% and analyzed by the upper one-sided exact 95% confidence bound (i.e., one-sided 0.05 alpha level).
- Primary effectiveness endpoint (at 2 months):
  - SOLO and TRIO:
    - Based on a two-sample t-test, for an assumed mean  $\pm$  standard deviation difference of  $6\pm 12$  mmHg with a two-sided 0.05 alpha level, a planned evaluable sample size of 128 subjects provides about 80% power (i.e., 64 subjects per arm).
  - RADIANCE-II:
    - Based on a 2:1 randomization, two-sample t-test, for an assumed mean  $\pm$  standard deviation difference of  $6\pm 12$  mmHg with a two-sided 0.05 alpha level, a planned evaluable sample size of 192 subjects provides about 90% power (i.e., 128 subjects in uRDN and 64 in Sham).


# Proportion of Subjects with $\geq 5$ , $\geq 10$ , & $\geq 15$ mmHg Reduction in Daytime SBP at 2 Months

Daytime SBP Reduction	uRDN	Sham Procedure	p-value
<i>SOLO (Off Med)</i>	(n=74)	(n=72)	
$\geq 5$ mmHg	49/74 (66.2%)	24/72 (33.3%)	<.0001
$\geq 10$ mmHg	32/74 (43.2%)	13/72 (18.1%)	0.0010
$\geq 15$ mmHg	19/74 (25.7%)	8/72 (11.1%)	0.0234
<i>RADIANCE-II (Off Med)</i>	(n=150)	(n=74)	
$\geq 5$ mmHg	64.1% (93/145)	34.2% (25/73)	<.0001
$\geq 10$ mmHg	47.6% (69/145)	16.4% (12/73)	<.0001
$\geq 15$ mmHg	25.5% (37/145)	9.6% (7/73)	0.0057
$\geq 20$ mmHg	11.7% (17/145)	6.8% (5/73)	0.2594
<i>TRIO (On Standardized Med)</i>	(n=69)	(n=67)	
$\geq 5$ mmHg	42/69 (60.9%)	28/67 (41.8%)	0.0260
$\geq 10$ mmHg	29/69 (42.0%)	17/67 (25.4%)	0.0401
$\geq 15$ mmHg	21/69 (30.4%)	10/67 (14.9%)	0.0311


Note not all p-values are adjusted for multiplicity.

# Durability of Effect – SOLO & RADIANCE-II


SOLO ITT Daytime SBP




R-II ITT Daytime SBP





SOLO ITT Office SBP



R-II ITT Office SBP




# Durability of Effect – TRIO




# TRIO, 2 months

## ANCOVA Results (ITT & PP & CA)


**TRIO (ITT)**  
Daytime SBP (chg ~ bp + arm)



**TRIO (PP)**  
Daytime SBP (chg ~ bp + arm)



**TRIO (CA)**  
Daytime SBP (chg ~ bp + arm)



*Note not all p-values are adjusted for multiplicity.*

# TRIO, 2 months, ITT

- Diagnostics for normality assumption of ANCOVA model
- In the middle row and left column, outliers can be seen at each end of the tails
  - Dots deviate away from the theoretical line derived from the normal distribution










Figure 6: TRIO Goodness of fit daytime ASBP 2-month ITT baseline adjusted

# TRIO, 2 months, ITT

- Top: uRDN
- Bottom: Sham
- Left: Baseline
- Middle: 2 Months
- Right: Difference
- x-axis: Daytime SBP (mmHg)
- y-axis: Frequency





# ANCOVA for Adjusting Baseline BP

- FDA Guidance: “Adjusting for Covariates in Randomized Clinical Trials for Drugs and Biological Products Guidance for Industry”
  - *This guidance provides recommendations for the use of covariates in the analysis of randomized, parallel group clinical trials that are applicable to both superiority trials and noninferiority trials.*
  - *The ICH E9 guidance strongly advises prespecification of “the principal features of the eventual statistical analysis,” including “how to account for [covariates] in the analysis to improve precision and to compensate for any lack of balance between treatment groups.”*
  - *The ICH E9 guidance also cautions against adjusting for “covariates measured after randomization because they could be affected by the treatments.”*

# ANCOVA on the ranks

- Quade (1967) JASA, Vol. 62, No. 320, pp. 1187-1200.
  - *The problem is then to test the hypothesis  $H_0$  that the conditional distribution of  $Y$  given  $X$  is the same for each population, where the alternatives of interest are those which imply that some populations tend to have greater values of  $Y$  than others for all fixed values of  $X$ .*
  - Note:
    - $H_0$ : The conditional distribution of BP change given baseline BP is the same for each group.
      - Rejecting this  $H_0$  may be due to the difference of distribution, either in shape or in location and both combined.
      - This estimand is different from the treatment effect obtained by the ANCOVA model.

# Hodges-Lehmann Estimates

- Hodges and Lehmann (1963)
  - The median of all paired differences (one subject from each group)
- Hollander and Wolfe (1999)
  - The asymptotic lower and upper confidence limits are for the location shift.
- Note:
  - Point estimate and confidence limits are neither adjusted for baseline BP nor associated with ANCOVA or ANCOVA on the ranks.
  - Symmetric distribution is preferred when summarizing the population parameter of distribution location (e.g., median of treatment effect).