

Relative sensitivity stratified by days post symptoms onset.

DPSO	RT-PCR Positives	Rapid Antigen Positives
0-1	3	2
2	17	15
3	10	10
4	12	12
5	1	1
6	1	1
All	44	41

Performance stratified by age groups

Age group	RT-PCR Positives	Rapid Antigen Positives
< 14	6	5
14 - 24	8	8
>24 - 64	29	27
≥ 65	1	1
All	44	41

Another prospective clinical study was conducted between January 2021 and May 2022 as a component of the Rapid Acceleration of Diagnostics (RADx®) initiative from the National Institutes of Health (NIH). A total of 7,361 individuals were enrolled via a decentralized clinical study design, with a broad geographical representation of the United States. Per inclusion criteria, all individuals were asymptomatic upon enrollment in the study and at least 14 days prior to it and did not have a SARS-CoV-2 infection in the three months prior to enrollment. Participants were assigned to one of three EUA authorized SARS-CoV-2 OTC rapid antigen tests to conduct serial testing (every 48 hours) for 15 days. If an antigen test was positive, the serial-antigen testing result is considered positive. At each rapid antigen testing time point, study subjects also collected a nasal swab for comparator testing using a home collection kit (using a 15-minute normalization window between swabs). SARS-CoV-2 infection status was determined by a composite comparator method on the day of the first antigen test, using at least two highly sensitive EUA RT-PCRs. If results of the first two molecular test were discordant a third highly sensitive EUA RT-PCR test was performed, and the final test result was based upon the majority rule. Study participants reported symptom status throughout the study using the MyDataHelps app. Two-day serial antigen testing is defined as performing two antigen tests 36-48 hours apart. Three-day serial antigen testing is defined as performing three antigen tests over five days with at least 48 hours between each test. Out of the 7,361 participants enrolled in the study, 5,609 were eligible for analysis. Among eligible participants, 154 tested positive for SARS-CoV-2 infection based on RT-PCR, of which 97 (62%) were asymptomatic on the first day of their infection, whereas 57 (39%) reported symptoms on the first day of infection. Pre-symptomatic subjects were included in the positive percent agreement (PPA) of asymptomatic individuals, if they were asymptomatic on the first day of antigen testing, regardless of whether they developed symptoms at any time after the first day of testing.

Performance of the antigen test with serial testing in individuals is described in the following table.

Data establishing PPA of COVID-19 antigen serial testing compared to the molecular comparator single day testing throughout the course of infection with serial testing. Data is from all antigen tests in study combined:

DAYS AFTER FIRST PCR POSITIVE TEST RESULT	ASYMPTOMATIC ON FIRST DAY OF TESTING			SYMPTOMATIC ON FIRST DAY OF TESTING		
	Second Result Day 3			Performance % PPA		
	Ag Positive / PCR Positive	Ag Positive / PCR Positive	Ag Positive / PCR Positive	Ag Positive / PCR Positive	Ag Positive / PCR Positive	Ag Positive / PCR Positive
	1 Test	2 Test	3 Test	1 Test	2 Test	3 Test
0	9/97 (9.3%)	35/89 (39.3%)	44/78 (56.4%)	34/57 (59.6%)	47/51 (92.2%)	44/47 (93.6%)
2	17/34 (50.0%)	23/34 (67.6%)	25/32 (78.1%)	58/62 (93.5%)	59/60 (98.3%)	43/43 (100%)
4	16/21 (76.2%)	15/20 (75.0%)	13/15 (86.7%)	55/58 (94.8%)	53/54 (98.1%)	39/40 (97.5%)
6	20/28 (71.4%)	21/27 (77.8%)	16/18 (88.9%)	27/34 (79.4%)	26/33 (78.8%)	22/27 (81.5%)
8	13/23 (56.5%)	13/22 (59.1%)	4/11 (36.4%)	12/17 (70.6%)	12/17 (70.6%)	7/11 (63.6%)
10	5/9 (55.6%)	5/8 (62.5%)		4/9 (44.4%)	3/7 (42.9%)	

1 Test = one (1) test performed on the noted days after first PCR positive test result. Day 0 is the first day of documented infection with SARS-CoV-2.
 2 Tests = two (2) tests performed an average of 48 hours apart. The first test performed on the indicated day and the second test performed 48 hours later.
 3 Tests = three (3) tests performance an average of 48 hours apart. The first test performed on the indicated day, the second test performed 48 hours later, and a final test performed 48 hours after the second test.

Analytical performance

Limit of Detection (LoD)

The SARS-CoV-2 positive specimen was prepared by spiking inactivated SARS-CoV-2 (isolate USA-WA1/2020) into pooled negative nasal matrix (PNM) that was confirmed to be negative with PCR. A serial dilution of specimens was tested by applying a sample volume of 50 µL to each nasal swab before elution and sample application was performed according to the test procedure described in the IFU. The LoD was determined to be 1.4 x 10³ TCID₅₀/mL. Based upon the testing procedure for this study the LoD of 1.4 x 10³ TCID₅₀/mL equates to 7.0 x 10³ TCID₅₀/swab

Omicron Testing

The performance of this test device in the detection of the Omicron variant of SARS-CoV-2 was evaluated in a dilution series of clinical specimens which were positive for the Omicron variant. This testing was conducted by the National Institutes of Health (NIH) as a component of the Rapid Acceleration of Diagnostics (RADx®) initiative. Specimen pools were prepared by the RADx team using pooled clinical samples from currently circulating Omicron strains and tested by RADx® to assess performance with the Omicron variant. Results from this dilution series cannot be compared to any devices tested with a different specimen pool and do not indicate that a test will have different clinical performance compared to other EUA authorized tests. Compared to an EUA authorized RT-PCR method, the STANDARD Q COVID-19 Ag Test 2.0 detected 100% of live virus Omicron samples at a Ct-value of 27.7 (n=5). Testing was also compared to two additional EUA authorized OTC antigen tests (Assay #1 and Assay #2). Omicron dilutions at lower viral concentrations (Ct-values higher than 27.7) were not detected by the STANDARD Q COVID-19 Ag Test 2.0 in this study.

Omicron Live Pool 1		Test Assay #1	Test Assay #2	STANDARD Q COVID-19 Ag Test 2.0
Dilutions	Avg N2 Ct (N=9)	Percent Positive N=5	Percent Positive N=5	Percent Positive N=5
Dilution 1	19.4	100	100	100
Dilution 2	20.6	100	100	100
Dilution 3	21.6	100	100	100
Dilution 4	22.4	100	100	100
Dilution 5	23.3	100	100	100
Dilution 6	24.5	0	100	100
Dilution 7	25.6	0	100	100
Dilution 8	26.5	0	0	100
Dilution 9	27.7	0	0	100
Dilution 10	28.5	0	0	0
Dilution 11	29.4	0	0	0
Dilution 12	30.3	0	0	0

Cross-reactivity

No cross-reactivity was observed for the following organisms at the indicated concentrations, except for SARS-coronavirus, which exhibited cross-reactivity when tested at 1.58 x 10⁴ TCID₅₀/mL. A titration of SARS-CoV was performed to determine the concentration at which cross reactivity was no longer observed. Cross reactivity was no longer observed for SARS-CoV at 1.58 x 10¹ TCID₅₀/mL. These results are not unexpected as the STANDARD Q COVID-19 Ag Test 2.0 targets the nucleocapsid protein, which is present in both the SARS-CoV and SARS-CoV-2 viruses and is highly homologous.

Microorganism / Specimen	Concentration Tested for Cross Reactivity	Result
Human coronavirus 229E	1.43 x 10 ⁵ TCID ₅₀ /mL	NEG
Human coronavirus OC43	8.50 x 10 ⁴ TCID ₅₀ /mL	NEG
Human coronavirus NL63	1.17 x 10 ⁵ TCID ₅₀ /mL	NEG
SARS-coronavirus	1.58 x 10 ⁴ TCID ₅₀ /mL	POS
SARS-coronavirus (1:1000)	1.58 x 10 ¹ TCID ₅₀ /mL	NEG
MERS-coronavirus	1.43 x 10 ⁵ TCID ₅₀ /mL	NEG
Adenovirus	1.43 x 10 ⁵ TCID ₅₀ /mL	NEG
Human metapneumovirus 4 Type B2	1.43 x 10 ⁵ TCID ₅₀ /mL	NEG
Parainfluenza virus 1	5.50 x 10 ⁵ TCID ₅₀ /mL	NEG

Parainfluenza virus 2	1.43 x 10 ⁵ TCID ₅₀ /mL	NEG
Parainfluenza virus 3	1.43 x 10 ⁵ TCID ₅₀ /mL	NEG
Parainfluenza virus 4b	1.43 x 10 ⁵ TCID ₅₀ /mL	NEG
Influenza A	1.43 x 10 ⁵ CEID ₅₀ /mL	NEG
Influenza B	1.43 x 10 ⁵ TCID ₅₀ /mL	NEG
Enterovirus 68	1.43 x 10 ⁵ TCID ₅₀ /mL	NEG
Respiratory syncytial virus	1.43 x 10 ⁵ PFU/mL	NEG
Rhinovirus	1.43 x 10 ⁵ TCID ₅₀ /mL	NEG
Haemophilus influenzae	1.00 x 10 ⁶ CFU/mL	NEG
Streptococcus pneumoniae	1.00 x 10 ⁶ CFU/mL	NEG
Streptococcus pyogenes	2.59 x 10 ⁶ CFU/mL	NEG
Candida albicans	1.00 x 10 ⁶ CFU/mL	NEG
Bordetella pertussis	1.00 x 10 ⁶ CFU/mL	NEG
Mycoplasma pneumoniae	2.57 x 10 ⁶ CFU/mL	NEG
Chlamydia pneumoniae	1.00 x 10 ⁶ IFU/mL	NEG
Legionella pneumophila	1.00 x 10 ⁶ CFU/mL	NEG
Mycobacterium tuberculosis	1.00 x 10 ⁶ CFU/mL	NEG
Pneumocystis carinii	1.00 x 10 ⁶ nuclei/mL	NEG
P. jirovecii-S. cerevisiae	8.10 x 10 ⁶ CFU/mL	NEG
Staphylococcus aureus subsp. aureus	1.00 x 10 ⁶ CFU/mL	NEG
Staphylococcus epidermidis	1.00 x 10 ⁶ CFU/mL	NEG
Pooled Negative Matrix	N/A	NEG

MICROBIAL INTERFERENCE

For the microorganisms which did not demonstrate cross-reactivity, additional microbial interference testing with SARS-CoV-2 positive samples spiked into pooled negative nasal matrix were performed and no microbial interference was observed.

ENDOGENOUS / EXOGENOUS CROSS-REACTIVITY STUDY

No cross-reactivity was observed for the following substances at the indicated concentrations. Each substance was spiked into pooled negative nasal matrix for testing.

Potentially Interfering Substance	Concentration	Result
Human Whole Blood (EDTA tube)	4% v/v	NEG
Mucin (Porcine Stomach, type II)	0.5%	NEG
Chloraseptic (Menthol/Benzocaine)	1.5 mg/mL	NEG
Naso GEL (NeilMed)	5% v/v	NEG
Nasal Drops (Phenylephrine)	15% v/v	NEG
Nasal Spray (Oxymetazoline)	15% v/v	NEG
Nasal Spray (Cromolyn)	15% v/v	NEG
Zicam	5% v/v	NEG
Homeopathic (Alkalol)	10% v/v	NEG
Sore Throat Phenol Spray	15% v/v	NEG
Tobramycin	4 µg/mL	NEG
Mupirocin	10 mg/mL	NEG
Tamiflu (Oseltamivir Phosphate)	5 mg/mL	NEG
Fluticasone Propionate	5% v/v	NEG
Body & Hand Lotion (CeraVe)	0.5% w/v	NEG
Body Lotion with 1.2% Dimethicone	0.5% w/v	NEG
Hand Lotion (Eucerin)	5% w/v	NEG
Hand Sanitizer with Aloe, 62% Ethyl Alcohol	5% w/v	NEG
Hand Sanitizer Cream Lotion (Vaseline)	15% v/v	NEG
Hand Sanitizer, 80% Ethanol, Fast Drying	15% v/v	NEG
Hand Soap Liquid Gel (Soft Soap)	10% w/v	NEG

ENDOGENOUS / EXOGENOUS INTERFERENCE SUBSTANCES STUDY

For the substances listed above, additional endogenous / exogenous interference studies were performed with samples containing SARS-CoV-2 in pooled negative nasal matrix. Each substance was spiked into a positive sample and no endogenous / exogenous interference was found, except for hand soap liquid gel, which caused a false negative result at a concentration of 10% w/v and 5% w/v. A positive result (no interference) was observed at a concentration of 1% w/v.

HIGH-DOSE HOOK EFFECT

SARS-CoV-2 cultured virus was spiked into pooled negative nasal matrix. SARS-CoV-2 cultured virus did not show hook effect at the virus stock concentration of 2.80 x 10⁶ TCID₅₀/mL.

BIBLIOGRAPHY

1. Coronaviruses. European Centre for Disease Prevention and Control. <https://www.ecdc.europa.eu/en/covid-19/latest-evidence/coronaviruses>. Accessed 6 Jan 2021.
2. Wu et al. Nature. 2020. 579:265-9.
3. Coronaviridae Study Group of the International Committee on Taxonomy of Viruses. Nat Microbiol. 2020. 5:536-44.
4. <https://www.who.int/emergencies/diseases/novel-coronavirus-2019/technicalguidance/naming-the-coronavirus-disease-%28covid-2019%29-and-the-virus-that-causes-it>
5. WHO. <https://www.who.int/dg/speeches/detail/who-director-general-s-opening-remarks-atthe-media-briefing-on-covid-19---11-march-2020>. Accessed 6 Jan 2021.
6. WHO. <https://www.who.int/publications-detail-redirect/diagnostic-testing-for-sars-cov-2>. Accessed 6 Jan 2021.
7. Centers for Disease Control and Prevention. <https://www.cdc.gov/coronavirus/2019-ncov/symptoms-testing/symptoms.html>. Accessed 6 Jan 2021.
8. World Health Organization (WHO). Antigen-detection in the diagnosis of SARS-CoV-2 infection using rapid immunoassays. Interim guidance. 11 September 2020.



For In Vitro Diagnostic Use.
 For use under an FDA Emergency Use Authorization (EUA) only.
 For Prescription use only.

Distribution in USA by SD BIOSENSOR USA, INC.

www.sdbiosensor.com
 COVID-19 Support +82-80-970-9700

Manufactured by SD Biosensor, Inc.
Head office : C-4th&5th, 16, Deogyeong-daero 1556beon-gil, Yeongtong-gu, Suwon-si, Gyeonggi-do, 16690, REPUBLIC OF KOREA
Manufacturing site : 74, Osongsangmyeong 4-ro, Osong-eup, Heungdeok-gu, Cheongju-si, Chungcheongbuk-do, 28161, REPUBLIC OF KOREA

Please contact us for any complaints/inquiries/suggestions via email (ts@sdbiosensor.com), phone (+82-80-970-9700) or website (www.sdbiosensor.com).

L23COV174FDAR1
 Issue date: 2023.09


STANDARD Q COVID-19 Ag Test 2.0

EN

REF Q-NCOV-08G
Cat No. 09COV174D

For In Vitro Diagnostic Use.
For use under an FDA Emergency Use Authorization (EUA) only.
For Prescription use only.

KIT CONTENTS



QUICK REFERENCE INSTRUCTIONS

Please study the Standard Q COVID-19 Ag Test 2.0 instructions for use thoroughly before using these Quick Reference Instructions or performing a test. This is not a complete product insert.

Intended Use

The STANDARD Q COVID-19 Ag Test 2.0 is a lateral flow immunoassay intended for the qualitative detection of nucleocapsid protein antigen from SARS-CoV-2 in direct anterior nasal swab specimens collected from individuals who are suspected of COVID-19 by their healthcare provider within the first six (6) days of symptom onset when tested at least twice over three days with at least 48 hours between tests, and individuals without symptoms or other epidemiological reasons to suspect COVID-19 when tested at least three times over five days with at least 48 hours between tests. The STANDARD COVID-19 Ag (Positive/Negative) Control swab is intended for use as an external quality control material to monitor the performance of STANDARD Q COVID-19 Ag Test 2.0.

Warning and Precautions

* Read all instructions carefully before performing the test. Failure to follow the instructions may result in inaccurate test results.
* In the USA, this product has not been FDA cleared or approved, but has been authorized by FDA under an Emergency Use Authorization. This product has been authorized only for the detection of proteins from SARS-CoV-2, not for any other viruses or pathogens. The emergency use of this product is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of in vitro diagnostics for detection and/or diagnosis of COVID-19 under Section 564(b)(1) of the Federal Food, Drug, and Cosmetic Act, 21 U.S.C. 360bbb-3(b)(1), unless the declaration is terminated or authorization is revoked sooner.
* Serial testing should be performed in individuals with negative results at least twice over three days (with 48 hours between tests) for symptomatic individuals and three times over five days (with at least 48 hours between tests) for asymptomatic individuals. You may need to purchase additional tests to perform this serial (repeat) testing.


PREPARATION

1 Bring test kit to room temperature (59-86 °F / 15-30 °C).

2 Wash your hands with soap and water, or use hand sanitizer before performing the test. Make sure you rinse thoroughly and your hands are dry before starting.



3 Check test expiry date on the back of the foil pouches. Do not use if the expiry date has passed.



Testing should commence immediately after opening the test device pouch.


4 Open the foil pouch by tearing along the tear-line. Remove the test device and desiccant package from the foil pouch. Place the test device on a flat surface.

5 Ensure that the test device is intact and that there are no green beads in the desiccant package (white and yellow beads are expected). Do not open the desiccant package.



TEST PROCEDURE

1 Remove the swab from the packaging. Ensure that you only touch the handle of the swab and not the soft pad at the tip.




2 Gently insert the swab into the nostril approximately 3/4 of an inch. Rotate the swab at least 5 times, brushing against the inside walls of the nostril, for about 15 seconds.

Do not just spin the swab.

Using the same swab, repeat in the second nostril.

Swab Both Nostrils



Right nostril


Left nostril

Inaccurate test results may occur if the nasal swab specimen is not properly collected. With children, the maximum depth of insertion into the nostril may be less than 3/4 of an inch, and you may need to have a second person to hold the child's head while swabbing.

3 Carefully open the extraction buffer tube avoiding spillage. If any liquid spills, do not use the tube.

4 Insert the swab into the extraction buffer tube until the soft pad is in the liquid. Squeeze the tube at the bottom and stir the swab more than 10 times while squeezing the tube. Remove the swab while squeezing the sides of the tube to extract the liquid from the swab. Dispose of the swab and seal the tube securely with the nozzle cap.

Ensure that the nozzle cap is securely fitted before proceeding to the next step



Failure to squeeze the tube can lead to incorrect results due to excess buffer in the swab.

5 Hold the tube upright above the sample well. Drop 4 drops onto the sample well.



Do not apply the liquid in the rectangular result window.

6 Set the timer and read the test result at 20 minutes. Do not read the result before 20 minutes or after 30 minutes.




Do not move or lift the test device during this time.

CONTROL PREPARATION AND TEST PROCEDURE

Positive/Negative control


Preparation

- Put the positive or negative control swab into an extraction buffer tube. While squeezing the buffer tube, stir the swab more than 10 times.
- Remove the swab while squeezing the sides of the tube to extract the liquid from the swab.
- Press the nozzle cap tightly onto the tube.



Test Procedure

- Apply 4 drops of extraction buffer to the sample well.



Do not apply the liquid in the rectangular result window.

- Read the test result in 20-30 minutes.



Do not move or lift the test device during this time.

INTERPRETATION OF TEST RESULTS



Negative result

If the Control (C) line is visible, but the Test (T) line is not visible, the test is negative.

To increase the chance that the negative result for COVID-19 is accurate, you should:

- Test again in 48 hours if the individual has symptoms on the first day of testing.
- Test 2 more times at least 48 hours apart if the individual does not have symptoms on the first day of testing.



Positive result

If the Control (C) line and the Test (T) line are visible, the test is positive. Any faint visible, colored test (T) line with the control line (C) should be read as positive.

Repeat testing does not need to be performed if patients have a positive result at any time.



Invalid result

If the control (C) line is not visible, the test is invalid. Re-test with a new swab and new test device. If your test result is still invalid, please contact below.

Distribution in USA by SD BIOSENSOR USA, INC.

www.sdbiosensor.com
COVID-19 Support +82-80-970-9700

Manufactured by SD Biosensor, Inc.

Head office : C-4th&5th, 16, Deogyong-daero 1556beon-gil, Yeongtong-gu, Suwon-si, Gyeonggi-do, 16690, REPUBLIC OF KOREA
Manufacturing site : 74, Osongsaengmyeong 4-ro, Osong-eup, Heungdeok-gu, Cheongju-si, Chungcheongbuk-do, 28161, REPUBLIC OF KOREA

Please contact us for any complaints/inquiries/suggestions via email (ts@sdbiosensor.com), phone (+82-80-970-9700) or website (www.sdbiosensor.com).

STANDARD

COVID-19 Ag Control swab

STANDARD™ COVID-19 Ag Control swab

For *in vitro* diagnostic use only.

For prescription use only.

For use under the Emergency Use Authorization (EUA) only.

PLEASE READ INSTRUCTIONS CAREFULLY
BEFORE YOU PERFORM THE TEST


SD BIOSENSOR

INTENDED USE

STANDARD COVID-19 Ag Control swab is intended for use as an external quality control material to monitor the performance of STANDARD Q COVID-19 Ag Test 2.0. For use only with STANDARD Q COVID-19 Ag Test 2.0.

KIT CONTENTS AND ACTIVE INGREDIENTS

- STANDARD COVID-19 Ag Positive Control swab: 5 swabs packed in individual pouches (red colored label). Positive control swabs contain recombinant SARS-CoV-2 nucleocapsid protein, BSA (Bovine Serum Albumin) and excipients. The recombinant nucleocapsid protein works as an artificial antigen and it is used to verify if the system is working properly.
- STANDARD COVID-19 Ag Negative Control swab: 5 swabs packed in individual pouches (blue colored label). Negative control swabs do not contain active ingredients.
- Instructions for Use

**MATERIALS REQUIRED BUT NOT PROVIDED**

- STANDARD Q COVID-19 Ag Test 2.0 (REF No. Q-NCOV-08G)
- Timer
- Disposable gloves
- Biohazard waste container

STORAGE AND STABILITY

Store the STANDARD COVID-19 Ag Control swab at 2-30°C/36-86°F. Kit materials are stable until expiration date printed on the outer box.

WARNING AND PRECAUTIONS

- Bring the kit contents and specimens to operating temperature (15-30 °C/ 59-86 °F) before testing.
- Do not reuse the control swabs.
- Do not use the control swabs if the pouch is damaged or the seal is broken.
- Do not smoke, drink or eat while testing.
- Do not use the control swabs for sample collection from patients.
- If there is evidence of microbial contamination in the reconstituted control in the extraction buffer, discard the control.
- Wear protective clothing, mask, and gloves when handling specimens and reagents. Wash hands thoroughly after the tests are done.
- Clean any spillage by immediately and thoroughly wiping up with a suitable disinfectant such as 1% sodium hypochlorite solution.
- Handle and dispose of all specimens, controls and materials used in testing as though they contain infectious agents.
- Dispose of all samples and materials used to perform the test as biohazard waste. Laboratory chemical and biohazard wastes must be handled and discarded in accordance with all local, state, and national regulations.
- Do not use kit materials if expiration date has passed.
- In the event of damage of packaging, contact the distributor of this product at SD BIOSENSOR Customer Service Center ts@sdbiosensor.com


CIRCUMSTANCES FOR RUNNING QUALITY CONTROL TESTS

It is important to perform quality control tests with positive and negative control materials to ensure your system is working properly. It is recommended that positive and negative controls be run:


- Each new operator prior to performing testing on patient specimens,
- When opening a new test kit lot,
- Whenever a new shipment of test kits is received,
- If the temperature of the test kit storage area falls outside of 2°-30°C (35°-86°F), and
- At periodic intervals as dictated by the user facility, country, state or local regulations and policies :
 - Control tests may be ran prior to performing each serial testing on patient specimens.
 - Serial testing of STNADAR Q COVID-19 Ag Test 2.0 should be performed in individuals with negative results at least twice over three days (with 48 hours between tests) for symptomatic individuals and three times over five days (with at least 48 hours between tests) for asymptomatic individuals.

PREPARING A QUALITY CONTROL TEST

- Bring the STANDARD™ COVID-19 Ag Test 2.0 and the STANDARD™ COVID-19 Ag Control swab to operating temperature (15-30°C / 59-86°F) at least 30 minutes prior to the test.
- Carefully read the Instructions for Use for the STANDARD™ COVID-19 Ag Test 2.0.
- Check the expiration date on the pouches of the control and of the test device. Do not use expired control or test devices.


**PERFORMING A QUALITY CONTROL TEST**

- Insert the positive or negative control swab into an extraction buffer tube which is in the STANDARD Q COVID-19 Ag Test 2.0. Stir the swab at least ten times while squeezing the sides of the buffer tube.




Warning: Failure to squeeze the tube can lead to incorrect results due to excess buffer in the swab.

- Remove the swab while squeezing the sides of the tube to extract the liquid from the swab. Dispose of the used swab in accordance with your biohazard waste disposal protocol.
- Press the nozzle cap tightly onto the tube.



- Apply 4 drops of the prepared control mixture into the specimen well of the test device.
- Read the results in accordance with the Instructions for Use accompanying the STANDARD Q COVID-19 Ag Test 2.0.



Warning: Read the results at 20 minutes. Do not read before 20 minutes or after 30 minutes. Even faint lines should be considered as a valid result.

INTERPRETATION OF QUALITY CONTROL RESULTS

STANDARD COVID-19 Ag Positive Control			
Result		Interpretation	Follow up
Test Line (T)	Control Line (C)		
Positive	Positive	Pass	-
Negative	Positive	Invalid	Retest*
No control line (C)		Invalid	Retest*

STANDARD COVID-19 Ag Negative Control			
Results		Interpretation	Follow up
Test Line (T)	Control Line (C)		
Negative	Positive	Pass	-
Positive	Positive	Invalid	Retest*
No control line (C)		Invalid	Retest*

* Use new test devices and new control for retesting. If the invalid control test result recurs, contact SD BIOSENSOR Customer Service Center: ts@sdbiosensor.com

LIMITATIONS

- For *in vitro* diagnostic use.
- For prescription use only.
- This product has not been FDA cleared or approved; but has been authorized by FDA under an Emergency Use Authorization (EUA).
- The emergency use of this product is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of *in vitro* diagnostics for detection and/or diagnosis of COVID-19 under Section 564(b)(1) of the Federal Food, Drug and Cosmetic Act, 21 U.S.C. § 360bbb-3(b)(1), unless the declaration is terminated or authorization is revoked sooner.
- Do not use kit contents beyond the expiration date printed on the outside of the box.
- The test procedure, precautions and interpretation of results for this test must be followed strictly when testing.
- Failure to follow the test procedure and interpretation of test results may adversely affect test performance and/or produce invalid results.
- This product is provided for quality assurance purposes and must not be used for calibration or as primary reference preparations in any test procedure.
- Adverse storage conditions or use of outdated reagents may produce erroneous results.
- Alterations in physical appearance may indicate instability or deterioration of this product. If there is evidence of microbial contamination in this product, discard it.

FOR OUR CUSTOMERS ONLY-LIMITED WARRANTY:

SD BIOSENSOR warrants that this product will be free from defects in material and workmanship until the expiration date printed on the label. THIS LIMITED WARRANTY IS IN LIEU OF ANY OTHER WARRANTY, EXPRESS OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. IN NO EVENT SHALL SD BIOSENSOR BE LIABLE FOR INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES.

IVD

Distribution in USA by SD BIOSENSOR USA, INC.

www.sdbiosensor.com
COVID-19 Support +82-80-970-9700

Manufactured by SD Biosensor, Inc.

Head office : C-4th&5th, 16, Deogyong-daero 1556beon-gil, Yeongtong-gu, Suwon-si, Gyeonggi-do, 16690, REPUBLIC OF KOREA

Manufacturing site : 74, Osongsaengmyeong 4-ro, Osong-eup, Heungdeok-gu, Cheongju-si, Chungcheongbuk-do, 28161, REPUBLIC OF KOREA

Please contact us for any complaints/inquiries/suggestions via email (ts@sdbiosensor.com), phone (+82-80-970-9700) or website (www.sdbiosensor.com).

L24COVC14JNCENR1
Issue date: 2023.09



Reference number



In vitro Diagnostics



Instructions for Use



Contains Sufficient for n -Tests



Caution



Note



Do not re-use.



Manufacturer



Date of manufacture



Keep away from moisture



Keep away from sunlight



Do not use if packaging is damaged



Use by



Batch code



Temperature limitation



Prescription use only



For use under an FDA Emergency Use Authorization (EUA) only