

cobas[®] SARS-CoV-2 & Influenza A/B v2

**Qualitative assay for use on the cobas[®] 6800/8800
Systems**

For use under Emergency Use Authorization (EUA) only

For in vitro diagnostic use

cobas[®] SARS-CoV-2 & Influenza A/B v2 P/N: 10033401190

cobas[®] SARS-CoV-2 & Influenza A/B Control Kit P/N: 09446133190

cobas[®] Buffer Negative Control Kit P/N: 09051953190

Table of contents

Summary and explanation of the test	4
Reagents and materials	6
cobas®SARSCoV-2 & Influenza A/B v2 reagents and controls.....	6..
cobas®Somni reagents for sample preparation.....	8..
Reagent storage and handling requirements.....	9..
Additional materials required.....	10
Alternative collection kits for swab specimens for use on cobas®6800/8800 Systems.....	10
Instrumentation and software required.....	11.
Precautions and handling requirements	12
Warnings and precautions.....	12
Reagent handling.....	13
Good laboratory practice.....	13
Sample collection, transport, and storage	14
Sample collection.....	14
Nasal (anterior nares) swab specimen collection on a healthcare worker or self-collected on site.....	15
Transport and storage.....	16
Instructions for use	17
Procedural notes.....	17
Running cobas®SARSCoV-2 & Influenza A/B v2.....	17
Specimens collected in cobas®PCR Media, 0.9% physiological saline, UTM® or UVT....	17
Specimens collected using cobas®PCR Media Uni or Dual Swab Sample Kit.....	18
Results	20
Quality control and validity of results.....	20
Interpretation of results on the cobas®6800/8800 Systems.....	20
Procedural limitations.....	24.
Conditions of Authorizations for Labs.....	25

Non-clinical performance evaluation	26
Key performance characteristics.....	26
Analytical sensitivity (Limit of Detection).....	26
Inclusivity	29
Precision (repeatability).....	30
Analytical specificity (cross-reactivity and microbial interference).....	33
Co-infection (competitive interference).....	35
Collection media equivalence.....	35
Clinical performance evaluation	36
Additional information	38
Key test features.....	38
Symbols.....	39
Technical support.....	40
Manufacturer and distributor.....	40
Trademarks and patents.....	40
Copyright.....	40
References.....	41
Document revision.....	42

Summary and explanation of the test

Intended use

cobas® SARS-CoV-2 & Influenza A/B v2 assay for use on the cobas® 800/8800 System (cobas® SARS-CoV-2 & Influenza A/B v2) is an automated multiplex real-time RT-PCR assay intended for simultaneous qualitative detection and differentiation of SARS-CoV-2, influenza A virus, and influenza B virus RNA in healthcare provider-collected nasal and nasopharyngeal swab specimens and self-collected anterior nasal swab specimens (collected in a healthcare setting with instruction by a healthcare provider) from individuals suspected of respiratory viral infection consistent with COVID-19 by their healthcare provider. cobas® SARS-CoV-2 & Influenza A/B v2 is intended for use as an aid in the differential diagnosis of SARS-CoV-2, influenza A, and influenza B in humans and is not intended to detect influenza C. Testing is limited to laboratories certified under the Clinical Laboratory Improvement Amendments of 1988 (CLIA), 42 U.S.C. §263a, that meet requirements to perform moderate to high complexity tests.

RNA from SARS-CoV-2, influenza A, and influenza B is generally detectable in respiratory specimens during the acute phase of infection. Positive results are indicative of the presence of SARS-CoV-2, influenza A, and/or influenza B RNA; clinical correlation with patient history and other diagnostic information is necessary to determine patient infection status. Positive results do not rule out bacterial infection or infection with other viruses. The agent detected may not be the definite cause of disease. Testing facilities within the United States and its territories are required to report all SARS-CoV-2 results to the appropriate public health authorities.

Negative results do not preclude infection from SARS-CoV-2, influenza A, and/or influenza B and should not be used as the sole basis for treatment or other patient management decisions. Negative results must be combined with clinical observations, patient history, and epidemiological information.

cobas® SARS-CoV-2 & Influenza A/B v2 is intended for use by qualified clinical laboratory personnel specifically instructed and trained in the techniques of real-time RT-PCR and on the use of cobas® 800/8800 Systems. In the United States, cobas® SARS-CoV-2 & Influenza A/B v2 is for use under the Food and Drug Administration's Emergency Use Authorization.

Explanation of the test

cobas® SARS-CoV-2 & Influenza A/B v2 is a qualitative test for the use on the cobas® 800 System and cobas® 8800 System for the detection of the 2019 novel coronavirus (SARS-CoV-2), influenza A, and influenza B RNA in both anterior nasal and nasopharyngeal swab samples collected in Copan Universal Transport Medium System (CUTM®) or BD™ Universal Viral Transport System (UVT) and additionally for nasal swab samples collected in PCR Media or 0.9% physiological saline. The RNA Internal Control, used to monitor the entire sample preparation and PCR amplification process, is introduced into each specimen during sample processing. In addition, the test utilizes external controls (a low titer positive control and a negative control).

Principles of the procedure

cobas® SARS-CoV-2 & Influenza A/B v2 is based on fully automated sample preparation (nucleic acid extraction and purification) followed by PCR amplification and detection. The cobas® 800/8800 Systems consist of the sample supply module, the transfer module, the processing module, and the lytia module. Automated data management is performed by the cobas® 800/8800 software, which assigns results for all tests. Results can be reviewed directly on the system screen, and printed as a report.

Nucleic acid from patient samples and added internal control RNA (RNA IC) molecules are simultaneously extracted. Nucleic acid is released by addition of proteinase and lysis reagent to the sample. The released nucleic acid binds to the silica surface of the added magnetic glass particles. Unbound substances and impurities, such as denatured protein, cellular debris and potential PCR inhibitors are removed with subsequent wash steps and purified nucleic acid is eluted from the magnetic glass particles with elution buffer at elevated temperature. External controls (positive and negative) are processed in the same way with cobas® SARS-CoV-2 & Influenza A/Bv2 run.

Selective amplification of SARS-CoV-2 target nucleic acid from the sample is achieved by the use of target-specific forward and reverse primers for ORF1a/b non-structural region that is unique to SARS-CoV-2. Additionally, a conserved region in the structural protein envelope gene was chosen for pan-Sarbecovirus detection. The pan-Sarbecovirus detection set will also detect SARS-CoV-2 virus. For influenza A, selective amplification of target nucleic acid from the sample is achieved by the use of two target-specific sets of forward and reverse primers one for the genomic region encoding matrix proteins 1 and 2 (M1/M2) and one for the gene encoding polymerase basic protein 2 (PB2). For influenza B, selective amplification of target nucleic acid from the sample is achieved by the use of target-specific forward and reverse primers for the nuclear export protein (NEP) / nonstructural protein 1 (NS1) genomic region. Selective amplification of RNA Internal Control is achieved by the use of non-competitive sequence specific forward and reverse primers which have no homology with the coronavirus/influenza genomes. Amplified target is detected by cleavage of fluorescently labeled oligonucleotide probe. A thermostable DNA polymerase enzyme is used for amplification.

The cobas® SARS-CoV-2 & Influenza A/Bv2 master mix contains detection probes which are specific for the coronavirus type SARS-CoV-2, members of the Sarbecovirus subgenus, influenza A virus, influenza B virus and the RNA Internal Control nucleic acid. The coronavirus/influenza A, influenza B and RNA Internal Control detection probes are each labeled with unique fluorescent dyes that act as a reporter. Each probe also has a second dye which acts as a quencher. When not bound to the target sequence, the fluorescent signals of the intact probes are suppressed by quencher dye. During the PCR amplification step, hybridization of the probes to the specific single-stranded DNA template results in cleavage of the probe by the 5' to 3' exonuclease activity of the DNA polymerase resulting in separation of the reporter and quencher dyes and the generation of a fluorescent signal. With each PCR cycle, increasing amounts of cleaved probes are generated and the cumulative signal of the reporter dye increases concomitantly. Each reporter dye is measured at defined wavelengths, which enables simultaneous detection and discrimination of the amplified coronavirus targets, influenza targets and the RNA Internal Control. The master mix includes deoxyuridine triphosphate (dUTP), instead of deoxythymidine triphosphate (dTTP), which is incorporated into the newly synthesized DNA (amplicon). Any contaminating amplicons from previous PCR runs are destroyed by the AmpErase enzyme [uracil-N-glycosylase], which is included in the PCR mix, when heated in the first thermal cycling step. However, newly formed amplicons are not destroyed since the AmpErase enzyme is inactivated once exposed to temperatures above 55°C.

Reagents and materials

The materials provided for cobas® SARS-CoV-2 & Influenza A/B v2 can be found in Table 1. Materials required, but not provided can be found in Table 2, Table 3, Table 4, Table 7, Table 8 and Table 9.

Refer to the **Reagents and materials** section and **Precautions and handling requirements** section for the hazard information for the product.

cobas® SARS-CoV-2 & Influenza A/B v2 reagents and controls

All unopened reagents and controls shall be stored as recommended in Table 1 to Table 4.

Table 1 cobas® SARS-CoV-2 & Influenza A/B v2

(SCoV2-FluA/B v2)

Store at 2-8°C

192 test cassette (P/N1003340119)

Kit components	Reagent ingredients	Quantity per kit 192 tests
Proteinase Solution (PASE)	Tris buffer, < 0.05% EDTA, calcium chloride, calcium acetate, 8% proteinase, glycerol EUH210: Safety data sheet available on request. EUH208: Contains subtilisin from <i>Bacillus subtilis</i> May produce an allergic reaction.	22.3 mL
RNA Internal Control (RNA IC)	Tris buffer, < 0.05% EDTA, 0.001% non-target related armored RNA construct containing primer aptamer specific sequence regions (non-infectious RNA in MS2 bacteriophage), < 0.1% sodium azide	21.2 mL
Elution Buffer (EB)	Tris buffer, 0.2% methyl-4 hydroxybenzoate	21.2 mL
Master Mix Reagent 1 (MMX-R1)	Manganese acetate, potassium hydroxide, < 0.1% sodium azide	7.5 mL
SCoV2-FluA/B v2 Master Mix Reagent 2 (SCoV2-FluA/B v2 MMX-R2)	Tricine buffer, potassium acetate, < 18% dimethyl sulfoxide, glycerol, < 0.1% Tween 20, EDTA, < 0.15% dATP, dCTP, dGTP, dUTPs, < 0.01% upstream and downstream SARS-CoV-2, Sarbecovirus, influenza A and influenza B primers, < 0.01% Internal Control forward and reverse primers, < 0.01% fluorescent-labeled oligonucleotide probes specific for SARS-CoV-2, Sarbecovirus, influenza A, influenza B and the RNA Internal Control, < 0.01% oligonucleotide aptamer, < 0.1% Z05D DNA polymerase, < 0.10% AmpErase (uracil-N-glycosylase) enzyme (microbial), < 0.1% sodium azide	9.7 mL

Table 2 cobas® SARS-CoV-2 & Influenza A/B Control Kit**(SCoV2-FluA/B CTL)**

Store at 2–8°C

(P/N 09446133190)

Kit components	Reagent ingredients	Quantity per kit
SCoV2-FluA/B Positive Control (SCoV2-FluA/B (+) C)	Tris buffer, < 0.05% Sodium azide, < 0.005% EDTA, < 0.003% Poly rA, < 0.01% Non-infectious plasmid DNA (microbial) containing SARS-CoV-2 sequence, < 0.01% Non-infectious plasmid DNA (microbial) containing pan-Sarbecovirus sequence, < 0.01% Non-infectious plasmid DNA (microbial) containing influenza A sequence, < 0.01% Non-infectious plasmid DNA (microbial) containing influenza B sequence	16 mL (16 x 1 mL)

Table 3 cobas® Buffer Negative Control Kit**(BUF (-) C)**


Store at 2-8°C

(P/N 09051953190)

Kit components	Reagent ingredients	Quantity per kit
cobas® Buffer Negative Control (BUF (-) C)	Tris buffer, < 0.1% sodium azide, EDTA, < 0.002% Poly rA RNA (synthetic)	16 mL (16 x 1 mL)

cobas® omni reagents for sample preparation

Table 4 cobas® omni reagents for sample preparation*

Reagents	Reagent ingredients	Quantity per kit	Safety symbol and warning**
cobas® omni MGP Reagent (MGP) Store at 2–8°C (P/N 06997546190)	Magnetic glass particles, Tris buffer, 0.1% methyl-4 hydroxybenzoate, < 0.1% sodium azide	480 tests	Not applicable
cobas® omni Specimen Diluent (SPEC DIL) Store at 2–8°C (P/N 06997511190)	Tris buffer, 0.1% methyl-4 hydroxybenzoate, < 0.1% sodium azide	4 x 875 mL	Not applicable
cobas® omni Lysis Reagent (LYS) Store at 2–8°C (P/N 06997538190)	43% (w/w) guanidine thiocyanate***, 5% (w/v) polydocanol***, 2% (w/v) dithiothreitol***, dihydro sodium citrate	4 x 875 mL	 <p>DANGER</p> <p>H302: Harmful if swallowed. H314: Causes severe skin burns and eye damage. H412: Harmful to aquatic life with lasting effects. EUH032: Contact with acids liberates very toxic gas. EUH071: Corrosive to the respiratory tract. P273: Avoid release to the environment. P280: Wear protective gloves/ protective clothing/ eye protection/ face protection/ hearing protection. P301 + P330 + P331: IF SWALLOWED: Rinse mouth. Do NOT induce vomiting. P303 + P361 + P353: IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water. P304 + P340 + P310: IF INHALED: Remove person to fresh air and keep comfortable for breathing. Immediately call a POISON CENTER/doctor. P305 + P351 + P338 + P310: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER/ doctor. 593-84-0 Guanidinium thiocyanate 9002-92-0 Polidocanol 3483-12-3 (R*,R*)-1,4-dimercaptobutane-2,3-diol</p>
cobas® omni Wash Reagent (WASH) Store at 15–30°C (P/N 06997503190)	Sodium citrate dihydrate, 0.1% methyl-4 hydroxybenzoate	4.2 L	Not applicable

* These reagents are not included in cobas SARS-CoV-2 & Influenza A/B2 test kits. See listing of additional materials required (Table 7, Table 8 and Table 9).

** Product safety labeling primarily follows EU GHS guidance

***Hazardous substance

Reagent storage and handling requirements

Reagents shall be stored and will be handled as specified in Table 5 and Table 6.

When reagents are not loaded on the cobas® 6800/8800 Systems, store them at the corresponding temperature specified in Table 5.

Table 5 Reagent storage (when reagent is not on the system)

Reagent	Storage temperature
cobas® SARS-CoV-2 & Influenza A/B v2	2–8°C
cobas® SARS-CoV-2 & Influenza A/B Control Kit	2–8°C
cobas® Buffer Negative Control Kit	2–8°C
cobas® omni Lysis Reagent	2–8°C
cobas® omni MGP Reagent	2–8°C
cobas® omni Specimen Diluent	2–8°C
cobas® omni Wash Reagent	15–30°C

Reagents loaded onto the cobas® 6800/8800 Systems are stored at appropriate temperatures and their expiration is monitored by the system. The cobas® 6800/8800 Systems allow reagents to be used only if all of the conditions shown in Table 6 are met. The system automatically prevents use of expired reagents. Table 6 allows the user to understand the reagent handling conditions enforced by the cobas® 6800/8800 Systems.

Table 6 Reagent expiry conditions enforced by the cobas® 6800/8800 Systems

Reagent	Kit expiration date	Open-kit stability	Number of runs for which this kit can be used	On-board stability (cumulative time on board outside refrigerator)
cobas® SARS-CoV-2 & Influenza A/B v2	Date not passed	90 days from first usage ^a	Max 40 runs ^a	Max 40 hours ^a
cobas® SARS-CoV-2 & Influenza A/B Control Kit	Date not passed	Not applicable ^b	Not applicable	Max 8 hours
cobas® Buffer Negative Control Kit	Date not passed	Not applicable ^b	Not applicable	Max 10 hours
cobas® omni Lysis Reagent	Date not passed	30 days from loading ^c	Not applicable	Not applicable
cobas® omni MGP Reagent	Date not passed	30 days from loading ^c	Not applicable	Not applicable
cobas® omni Specimen Diluent	Date not passed	30 days from loading ^c	Not applicable	Not applicable
cobas® omni Wash Reagent	Date not passed	30 days from loading ^c	Not applicable	Not applicable

^aThe performance has not been established for suggested use cycles and time, but is based on similar reagents used on the same system.

^bSingle use reagents

^cTime is measured from the first time that reagent is loaded onto the cobas® 6800/8800 Systems.

Additional materials required

Table 7 Materials and consumables for use on cobas® 6800/8800 Systems

Material	P/N
cobas® omni Processing Plate	05534917001
cobas® omni Amplification Plate	05534941001
cobas® omni Pipette Tips	05534925001
cobas® omni Liquid Waste Container	07094388001
cobas® omni Lysis Reagent	06997538190
cobas® omni MGP Reagent	06997546190
cobas® omni Specimen Diluent	06997511190
cobas® omni Wash Reagent	06997503190
Solid Waste Bag and Solid Waste Container or Solid Waste Bag With Insert and Kit Drawer	07435967001 and 07094361001 or 08030073001 and 08387281001
cobas® omni Secondary Tubes 13x75 (optional)	06438776001
cobas® PCR Media Tube Replacement Cap Kit	07958056190
cobas® PCR Media Disposable Tube Stand (optional)	07958064190
MPA RACK 13 or 16 MM ^a	N/A
RD5 RACK – RD Standard rack ^a	N/A

^aPlease contact your local Roche representative for a detailed order list for sample racks, racks for clotting accepted on the instrument and compatible with the assay.

Alternative collection kits for swab specimens for use on the cobas® 6800/8800 Systems

Table 8 Alternative specimen collection kits used with cobas® SARS-CoV-2 & Influenza A/Bv2

Collection Kit	P/N
cobas® PCR Media Uni Swab Sample Kit	07958030190
cobas® PCR Media Dual Swab Sample Kit	07958021190
cobas® PCR Media 100 tube kit	06466281190
cobas® Uni Swab 100 Kit	09205098190

Instrumentation and software required

The cobas® 6800/8800 software and cobas® SARS-CoV-2 & Influenza A/B/2 analysis package (SW cobas® CoV2FluA/B ASAP) must be installed on the instrument. The Instrument Gateway (IG) server will be provided with the system.

Table 9 Instrumentation

Equipment	P/N
cobas® 6800 System (Moveable Platform)	05524245001 and 06379672001
cobas® 6800 System (Fixed Platform)	05524245001 and 06379664001
cobas® 8800 System	05412722001
Sample Supply Module	06301037001
Instrument Gateway	06349595001

For additional information, please refer to the cobas® 6800/8800 Systems – User Assistance and/or User Guide.

Precautions and handling requirements

Warnings and precautions

As with any test procedure, good laboratory practice is essential to the proper performance of this assay. Due to the high sensitivity of this test, care should be taken to keep reagents and amplification mixtures free of contamination.

- For in vitro diagnostic use under the Food and Drug Administration's Emergency Use Authorization (EUA) only.
- For prescription use only.
- cobas® SARS-CoV-2 & Influenza A/Bv2 has not been FDA cleared or approved but has been authorized for emergency use by FDA under an EUA for use by authorized laboratories certified under the Clinical Laboratory Improvement Amendments of 1988 (CLIA), 42 U.S.C. §263a, that meet requirements to perform moderate or high complexity tests.
- cobas® SARS-CoV-2 & Influenza A/Bv2 has been authorized only for the detection and differentiation of nucleic acid from SARS-CoV-2, influenza A, and influenza B, not for any other viruses or pathogens.
- The emergency use of cobas® SARS-CoV-2 & Influenza A/Bv2 is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of in vitro diagnostics for detection and/or diagnosis of COVID-19 under Section 564(b)(1) of the Federal Food, Drug, and Cosmetic Act, 21 U.S.C. § 360b(b)(1), unless the declaration is terminated or authorization is revoked sooner.
- Results (positive and negative) for influenza should be interpreted with caution. If an influenza result is inconsistent with clinical presentation and/or other clinical and epidemiological information, FDA-authorized Influenza NAATs are available for confirmation if clinically indicated.
- Laboratories within the United States and its territories are required to report SARS-CoV-2 results to the appropriate public health authorities.
- All patient samples should be handled as if infectious, using good laboratory procedures as outlined in Biosafety in Microbiological and Biomedical Laboratories and in the CLSI Document M29. Only personnel proficient in handling infectious materials and the use of cobas® SARS-CoV-2 & Influenza A/Bv2 and cobas® 6800/8800 Systems should perform this procedure.
- All human-sourced materials should be considered potentially infectious and should be handled with universal precautions. If spillage occurs, immediately disinfect with a freshly prepared 0.6% sodium hypochlorite distilled or deionized water (dilute household bleach or 1:10) or appropriate site procedures.
- The use of sterile disposable pipettes and nuclease-free pipette tips is recommended by the supplier or specified required consumables to ensure optimal test performance.
- Safety Data Sheets (SDS) are available on request from your local Roche representative.
- Closely follow procedures and guidelines provided to ensure that the test is performed correctly. Any deviation from the procedures and guidelines may affect optimal test performance.
- False positive results may occur if carryover of samples is not adequately controlled during sample handling and processing.

Reagent handling

- Handle all reagents, controls, and samples according to good laboratory practice in order to prevent carryover of samples or controls.
- Before use, visually inspect each reagent cassette, diluent, lysis reagent and wash reagent to ensure that there are no signs of leakage. If there is any evidence of leakage, do not use that material for testing.
- **cobas®omni** Lysis Reagent contains guanidine thiocyanate, a potentially hazardous chemical. Avoid contact of reagents with the skin, eyes or mucous membranes. If contact does occur, immediately wash with generous amounts of water; otherwise burns can occur.
- **cobas®SARS-CoV-2 & Influenza A/Bv2** test kit, **cobas®SARS-CoV-2 & Influenza A/B** Control kit, **cobas®Buffer Negative Control kit**, **cobas®omni** MGP Reagent, and **cobas®omni** Specimen Diluent contain sodium azide as a preservative. Avoid contact of reagents with the skin, eyes, or mucous membranes. If contact does occur, immediately wash with generous amounts of water; otherwise, burns can occur. If these reagents are spilled, dilute with water before wiping dry.
- Do not allow **cobas®omni** Lysis Reagent, which contains guanidine thiocyanate, to contact sodium hypochlorite (bleach) solution. This mixture can produce a highly toxic gas.
- Dispose of all materials that have come in contact with samples and reagents in accordance with country, state, and local regulations.

Good laboratory practice

- Do not pipette by mouth.
- Do not eat, drink, or smoke in designated work areas.
- Wear laboratory gloves, laboratory coats, and eye protection when handling samples and reagents. Gloves must be changed between handling samples and **cobas®SARS-CoV-2 & Influenza A/Bv2** kits, **cobas®SARS-CoV-2 & Influenza A/B** Control kit, **cobas®Buffer Negative Control kit** and **cobas®omni** reagents to prevent contamination. Avoid contaminating gloves when handling samples and controls.
- Wash hands thoroughly after handling samples and kit reagents, and after removing the gloves.
- Thoroughly clean and disinfect all laboratory work surfaces with a freshly prepared solution of 0.6% sodium hypochlorite in distilled or deionized water (dilute household bleach 1:10). Follow by wiping the surface with 70% ethanol.
- If spills occur on the **cobas®6800/8800** instrument, follow the instructions in the **cobas®6800/8800 Systems – User Assistance and/or User Guide** to properly clean and decontaminate the surface of instrument(s).

Sample collection, transport, and storage

Note: Handle all samples and controls as if they are capable of transmitting infectious agents.

Store all samples at specified temperatures.

Sample stability is affected by elevated temperatures.

Always use caution when transferring specimens from a primary collection tube to a secondary tube.

Use pipettes with aerosol-barrier or positive-displacement tips to handle specimens.

Always use a new pipette tip for each specimen.

Ensure samples are equilibrated to room temperature prior to transfer into a secondary tube.

Sample collection

Table 10 summarizes what collection devices can be used with specific sample types.


Table 10 Overview of collection devices and sample types

Collection Device	Nasopharyngeal	Nasal
Copan Universal Transport Media (UTM-RT®)	√	√
BD™ Universal Viral Transport (UVT)	√	√
0.9% Physiological saline		√
cobas® PCR Media Uni Swab Sample Kit		√
cobas® PCR Media Dual Swab Sample Kit		√
cobas® PCR Media Kit (and 100 Tube PCR Media Kit)		√



- Collect nasal and nasopharyngeal specimens according to standard collection technique using flocked or polyester-tipped swabs and immediately place in 3 mL of Copan Universal Transport Medium (UTM-RT®) or BD™ Universal Viral Transport (UVT) or equivalent.
- Collect nasal specimens according to standard collection technique using flocked or polyester-tipped swabs and immediately place into cobas® CR Media tube from cobas® CR Media Kit (P/N 06466281190)
- Collect nasal specimens using cobas® CR Media Uni Swab Sample Kit (P/N 07958030190) or cobas® CR Media Dual Swab Sample Kit (P/N 07958021190) according to instructions below
- Refer to the Instructions for Use of the Collection Devices for hazard information.


Nasal (anterior nares) swab specimen collection - healthcare worker or self-collected on site

WARNING: DO NOT PREWET SWAB IN cobas® PCR MEDIA BEFORE COLLECTION!



OR





- Collect nasal specimens according to standard collection technique using flocked or polyester-tipped swabs and immediately place in 3 mL of 0.9% physiological saline.

Transport and storage

- Transportation of collected specimens must comply with all applicable regulations for the transport of etiologic agents.
- Samples collected in UTM-RT®
 - After collection, specimens can be stored for up to 48 hours at 2-25°C followed by up to 3 days at 2-8°C.
- Samples collected in cobas® PCR Media
 - After collection, specimens can be stored for up to 24 hours at 2-25°C followed by up to 3 days at 2-8°C.
- Samples collected in 0.9% physiological saline
 - After collection, specimens can be stored for up to 48 hours at 2-25°C followed by up to 3 days at 2-8°C.
- If delivery and processing of samples exceeds specified time periods, specimens should be transported in dry ice and once in laboratory frozen at -70°C or colder

Instructions for use

Procedural notes

- Do not use **cobas® SARS-CoV-2 & Influenza A/Bv2** reagents, **cobas® SARS-CoV-2 & Influenza A/B Control Kit**, **cobas® Buffer Negative Control Kit** or **cobas® omni** reagents after their expiry dates.
- Do not reuse consumables. They are for one-time use only.
- Ensure that specimen barcode labels on sample tubes are visible through the openings on the side of the sample racks. Refer to the **cobas® 800/8800 Systems User Guide** for proper barcode specifications and additional information on loading sample tubes.
- Refer to the **cobas® 800/8800 Systems User Assistance and/or User Guide** for proper maintenance of instruments.

Running cobas® SARS-CoV-2 & Influenza A/B v2

cobas® SARS-CoV-2 & Influenza A/Bv2 can be run with a minimum required sample volume of 0.6 mL in the **cobas® omni Secondary Tube** for specimens collected in **Opan Universal Transport Medium (UTM-RT®)**, **BD™ Universal Viral Transport (UVT)**, **cobas® PCR Media** or 0.9% physiological saline. Specimens collected using **cobas® PCR Media Uni Swab Sample Kit**, **cobas® PCR Media Dual Swab Sample Kit** can be run in their primary collection tube with a minimum required sample volume of 1.0 mL.

Specimens collected in cobas® PCR Media, 0.9% physiological saline, UTM-RT® or UVT

Specimens collected in tubes compatible with **cobas® 800/8800 Systems** may be loaded directly onto the **cobas® 6800/8800 Systems**. The swab must be removed from the sample tube prior to direct loading onto the system. Specimens collected in tubes which are not compatible with the **cobas® 800/8800 Systems** must be transferred into a secondary tube prior to processing on the **cobas® 800/8800 Systems**. The **cobas® omni Secondary Tube** is the preferred option. If using frozen samples in secondary tubes, place the samples at room temperature (15-30°C) until completely thawed and then briefly mix (e.g., vortex for 3-5 seconds) and centrifuge to collect all sample volume at the bottom of the tube. Samples should be processed using the sample type selection in the user interface (UI) as described in Table 11. Additional tubes for testing **cobas® SARS-CoV-2 & Influenza A/B v2** are available. Contact your local Roche representative for detailed testing instructions and an order list of primary tubes and secondary tubes compatible with the instruments.

Follow the steps below to transfer patient sample from a primary collection tube into a **cobas® omni Secondary Tube**:

- Unscrew the primary sample tube cap.
- Lift the cap and any attached swab to allow a pipette to be inserted into the sample tube.
- Transfer at least 0.6 mL into the prepared barcoded secondary tube.
- Transfer secondary tube to a rack. Close the primary sample tube cap.

Specimens collected using cobas® PCR Media Uni or Dual Swab Sample Kit

Samples collected using cobas® PCR Media Uni or Dual Swab Sample Kit must be uncapped and can be loaded directly onto racks for processing on cobas® 6800/8800 Systems. Transfer into a secondary tube is not necessary. cobas® PCR Media tubes fit on to the MPA RACK 16 MM and can be processed with the swab remaining in the tube. Samples collected using cobas® PCR Media Uni or Dual Swab Sample Kit should be processed using the cobas® PCR Media swab sample type selection in the user interface (UI) of the cobas® SARS-CoV-2 & Influenza A/Bv2 as described in Table 11.

A properly collected swab specimen should have a single swab with the shaft broken at the scoreline. Swab shafts which are broken above the score line will appear longer than normal and may also be bent over to fit into cobas® PCR Media tube. This may create an obstruction to the pipetting system which might cause the loss of sample, repeat results and/or mechanical damage to the instrument. In the event that a swab specimen has an improperly broken shaft, remove the swab prior to sample processing on cobas® 6800/8800 Systems. Use caution when disposing of specimen swabs; avoid splashing or touching swabs to other surfaces during disposal to prevent contamination.

Incoming cobas® PCR Media primary swab specimen tubes with no swabs or with two swabs have been collected according to the instructions in their respective collection kit IFU and should not be tested. If the sample containing two swabs in the cobas® PCR Media primary tubes must be tested, transfer 0.6 mL into the prepared barcoded secondary tube.


Occasionally, incoming swab specimens contain excessive mucus which may induce a pipetting error (e.g., clot or other obstruction) on the cobas® 6800/8800 Systems. Prior to retesting of specimens that exhibited clots during initial processing, remove and discard the swab, then cap and vortex these specimens for 30 seconds to disperse the excessive mucus. Specimens can be processed twice on cobas® 6800/8800 Systems while the swab is in the collection tube. If additional testing is required, or if the first test fails due to specimen pipetting error (e.g., clot or other obstruction), the swab must be removed and the remaining fluid must have a minimum volume of 1.0 mL. The test procedure is described in detail in the cobas® 6800/8800 Systems- User Assistance and/or User Guide. Figure 1 below summarizes the procedure.

Table 11 Sample type selection in the user interface of the cobas® SARS-CoV-2 & Influenza A/Bv2

Collection kit /Matrix type	Minimum volume (mL) Processing tube	Process as Sample Type
Copan Universal Transport Medium (UTM-RT®) BD™ Universal Viral Transport 0.9% Physiological saline cobas® PCR Media Kit	0.6 mL cobas® omni Secondary tube	VTM
Copan Universal Transport Medium (UTM-RT®) BD™ Universal Viral Transport 0.9% Physiological saline cobas® PCR Media Kit	Compatible tubes without swab inside the tube; for dead volume contact your local Roche representative	VTM
cobas® PCR Media Uni or Dual Swab Sample Kit	1.0 mL Primary tube	cobas® PCR Media swab

The test procedure is described in detail in [cobas® 800/8800 Systems User Assistance and/or User Guide](#). Figure 1 below summarizes the procedure.

Figure 1 cobas® SARSCoV-2 & Influenza A/Bv2 procedure



Results

The cobas® 6800/8800 System automatically detects the SARS-CoV-2, influenza A and influenza B, for each individually processed sample and control, displaying individual target results for samples as well as test validity and overall results for controls.

Quality control and validity of results

- One cobas® Buffer Negative Control [BUF (-) C] and one [SCoV2-FluA/B CTL] are processed with each batch.
- In the cobas® 6800/8800 software and/or report, check for flags and their associated results to ensure the batch validity.
- All flags are described in the cobas® 6800/8800 Systems User Guide.
- The batch is valid if no flags appear for any controls. If the batch is invalid, repeat testing of the entire batch.

Validation of results is performed automatically by the cobas® 6800/8800 software based on negative and positive control performance.

Interpretation of results on the cobas® 6800/8800 Systems

Display examples for cobas® SARS-CoV-2 & Influenza A/B v2 are shown in Figure 2.

Figure 2 Example of cobas® SARS-CoV-2 & Influenza A/B v2 results display on the cobas® 6800/8800 System

Test	Sample ID	Valid*	Flags	Sample type	Overall result*	Target 1	Target 2	Target 3	Target 4
SCoV2-FluA/B 400 µL	Sample_01	NA		VTM	NA	FluA Negative	SCoV2 Negative	PanSarB Negative	FluB Negative
SCoV2-FluA/B 400 µL	Sample_02	NA	Y40T	VTM	NA	Invalid	Invalid	Invalid	Invalid
SCoV2-FluA/B 400 µL	Sample_03	NA		VTM	NA	FluA Positive	SCoV2 Negative	PanSarB Negative	FluB Negative
SCoV2-FluA/B 400 µL	Sample_04	NA		VTM	NA	FluA Negative	SCoV2 Positive	PanSarB Positive	FluB Negative
SCoV2-FluA/B 400 µL	Sample_05	NA		VTM	NA	FluA Negative	SCoV2 Negative	PanSarB Negative	FluB Positive
SCoV2-FluA/B 400 µL	Sample_06	NA		VTM	NA	FluA Negative	SCoV2 Negative	PanSarB Positive	FluB Negative
SCoV2-FluA/B 400 µL	Sample_07	NA	C01H2	VTM	NA	FluA Positive	Invalid	Invalid	Invalid
SCoV2-FluA/B 400 µL	Sample_08	NA	C01H1	VTM	NA	Invalid	SCoV2 Positive	Invalid	FluB Positive
SCoV2-FluA/B	C161420284090428828404	Yes		(-) Ctrl	Valid	Valid	Valid	Valid	Valid
SCoV2-FluA/B	C161420284093009580264	Yes		SCoV2-FluA/B (+) C	Valid	Valid	Valid	Valid	Valid

*The "Valid" and "Overall Result" columns are not applicable to sample **cobas® SARS-CoV-2 & Influenza A/B v2**. Values reported in these columns are not applicable and do not impact the validity of results reported within target result columns. Refer to **Table 12, cobas® SARS-CoV-2 & Influenza A/B v2 results interpretation** for specific instructions on test results interpretation.

Interpretation of results

For a valid batch, check each individual sample for flags in the **cobas® 6800/8800** software and/or report. The result interpretation should be as follows:

- A valid batch may include both valid and invalid sample results.
- Invalid results for one or more target combinations are possible and are reported out specifically for each channel.
- Results of this test should only be interpreted in conjunction with information available from clinical evaluation of the patient and patient history.

Results and their corresponding interpretation for detecting SARS-CoV-2 & Influenza A/B are shown below (Table 12).

Table 12 cobas® SARS-CoV-2 & Influenza A/B v2 results interpretation

Target 1 Influenza A*	Target 2 SARS- CoV-2	Target 3 Pan- Sarbecovirus	Target 4 Influenza B*	Interpretation
Negative	Negative	Negative	Negative	No target RNA Detected
Negative	Negative	Negative	Positive	Influenza B RNA Detected
Positive	Negative	Negative	Negative	Influenza A RNA Detected
Positive	Negative	Negative	Positive	Influenza A and Influenza B RNA Detected

Target 1 Influenza A*	Target 2 SARS- CoV-2	Target 3 Pan- Sarbecovirus	Target 4 Influenza B*	Interpretation
Negative	Negative	Positive	Negative	<p>Presumptive Positive for SARS-CoV-2 RNA. A negative SARS-CoV-2 result and a positive pan-Sarbecovirus results is suggestive of 1) a sample at concentrations near or below the limit of detection of the test, 2) a mutation in the SARS-CoV-2 target region in the oligo binding site, 3) infection with some other Sarbecovirus (e.g., SARS-CoV or some other Sarbecovirus previously unknown to infect humans), or 4) other factors.</p> <p>For samples with a Presumptive Positive result, additional confirmatory testing may be conducted, if it is necessary to differentiate between SARS-CoV-2 and SARS-CoV-1 or other Sarbecovirus currently unknown to infect humans, for epidemiological purposes or clinical management.</p>
Negative	Negative	Positive	Positive	<p>Presumptive Positive for SARS-CoV-2 RNA and influenza B RNA Detected. A negative SARS-CoV-2 result and a positive pan-Sarbecovirus results is suggestive of 1) a sample at concentrations near or below the limit of detection of the test, 2) a mutation in the SARS-CoV-2 target region in the oligo binding site, 3) infection with some other Sarbecovirus (e.g., SARS-CoV or some other Sarbecovirus previously unknown to infect humans), or 4) other factors.</p> <p>For samples with a Presumptive Positive result, additional confirmatory testing may be conducted, if it is necessary to differentiate between SARS-CoV-2 and SARS-CoV-1 or other Sarbecovirus currently unknown to infect humans, for epidemiological purposes or clinical management.</p>
Positive	Negative	Positive	Negative	<p>Influenza A RNA Detected and Presumptive Positive for SARS-CoV-2 RNA. A negative SARS-CoV-2 result and a positive pan-Sarbecovirus results is suggestive of 1) a sample at concentrations near or below the limit of detection of the test, 2) a mutation in the SARS-CoV-2 target region in the oligo binding site, 3) infection with some other Sarbecovirus (e.g., SARS-CoV or some other Sarbecovirus previously unknown to infect humans), or 4) other factors.</p> <p>For samples with a Presumptive Positive result, additional confirmatory testing may be conducted, if it is necessary to differentiate between SARS-CoV-2 and SARS-CoV-1 or other Sarbecovirus currently unknown to infect humans, for epidemiological purposes or clinical management.</p>

Target 1 Influenza A*	Target 2 SARS- CoV-2	Target 3 Pan- Sarbecovirus	Target 4 Influenza B*	Interpretation
Positive	Negative	Positive	Positive	Influenza A RNA Detected, Presumptive Positive for SARS-CoV-2 RNA, and influenza B RNA Detected. A negative SARS-CoV-2 result and a positive pan-Sarbecovirus results is suggestive of 1) a sample at concentrations near or below the limit of detection of the test, 2) a mutation in the SARS-CoV-2 target region in the oligo binding site, 3) infection with some other Sarbecovirus (e.g., SARS-CoV or some other Sarbecovirus previously unknown to infect humans), or 4) other factors. For samples with a Presumptive Positive result, additional confirmatory testing may be conducted, if it is necessary to differentiate between SARS-CoV-2 and SARS-CoV-1 or other Sarbecovirus currently unknown to infect humans, for epidemiological purposes or clinical management.
Negative	Positive	Negative	Negative	SARS-CoV-2 RNA Detected. A positive SARS-CoV-2 result and a negative pan-Sarbecovirus results is suggestive of 1) a sample at concentrations near or below the limit of detection of the test, 2) a mutation in the pan-Sarbecovirus target region, or 3) other factors.
Negative	Positive	Negative	Positive	SARS-CoV-2 RNA and influenza B RNA Detected. A positive SARS-CoV-2 result and a negative pan-Sarbecovirus results is suggestive of 1) a sample at concentrations near or below the limit of detection of the test, 2) a mutation in the pan-Sarbecovirus target region, or 3) other factors.
Positive	Positive	Negative	Negative	Influenza A RNA and SARS-CoV-2 RNA Detected. A positive SARS-CoV-2 result and a negative pan-Sarbecovirus results is suggestive of 1) a sample at concentrations near or below the limit of detection of the test, 2) a mutation in the pan-Sarbecovirus target region, or 3) other factors.
Positive	Positive	Negative	Positive	Influenza A RNA, SARS-CoV-2 RNA, and influenza B RNA Detected. A positive SARS-CoV-2 result and a negative pan-Sarbecovirus results is suggestive of 1) a sample at concentrations near or below the limit of detection of the test, 2) a mutation in the pan-Sarbecovirus target region, or 3) other factors.
Negative	Positive	Positive	Negative	SARS-CoV-2 RNA Detected
Negative	Positive	Positive	Positive	SARS-CoV-2 RNA and influenza B RNA Detected
Positive	Positive	Positive	Negative	Influenza A RNA and SARS-CoV-2 RNA Detected
Positive	Positive	Positive	Positive	Influenza A RNA, SARS-CoV-2 RNA, and influenza B RNA Detected

If any individual target result is invalid, the presence or absence of that individual target cannot be determined. If all initial valid target results can be interpreted as described in the table.

*Results (positive and negative) for influenza should be interpreted with caution. If an influenza result is inconsistent with clinical presentation and/or other clinical and epidemiological information, FDA-cleared influenza NAATs are available for confirmation if clinically indicated.

Procedural limitations

- cobas® SARS-CoV-2 & Influenza A/B v2 has been evaluated only for use in combination with cobas® SARS-CoV-2 & Influenza A/B Control Kit, cobas® Buffer Negative Control, cobas® omni MGP Reagent, cobas® omni Lysis Reagent, cobas® omni Specimen Diluent, and cobas® omni Wash Reagent for use on cobas® 800/8800 Systems.
- Reliable results depend on proper sample collection, storage and handling procedures.
- This test is intended to be used for the detection of SARS-CoV-2, influenza A, and influenza B RNA in nasopharyngeal aspirate swab samples collected in Copan Universal Transport Medium (UTM-RT) or BD™ Universal Viral Transport System (UVT) and nasal swab samples collected in cobas® CR Media and 0.9% physiological saline. Testing of other sample types with cobas® SARS-CoV-2 & Influenza A/B v2 may result in inaccurate results.
- Detection of SARS-CoV-2 and influenza A/B RNA may be affected by sample collection methods, patient factors (e.g., presence of symptoms), and/or stage of infection.
- As with any molecular test, mutations within the target region of cobas® SARS-CoV-2 & Influenza A/B v2 could affect primer and/or probe binding resulting in failure to detect the presence of virus.
- Due to inherent differences between technologies, it is recommended that, prior to switching from one technology to the next, users perform method correlation studies in their laboratory to qualify technology differences. One hundred percent agreement between the results should not be expected due to aforementioned differences between technologies. Users should follow their own specific policies/procedures.
- False negative or invalid results may occur due to interference. The Internal Control is included in cobas® SARS-CoV-2 & Influenza A/B v2 to help identify the specimens containing substances that may interfere with nucleic acid isolation and PCR amplification.
- The addition of AmpErase enzyme into cobas® SARS-CoV-2 & Influenza A/B v2 Master Mix reagent enables selective amplification of target RNA; however, good laboratory practices and careful adherence to the procedures specified in this Instructions For Use document are necessary to avoid contamination of reagents.
- Results (positive and negative) for influenza should be interpreted with caution. If an influenza result is inconsistent with clinical presentation and/or other clinical and epidemiological information, FDA-authorized influenza NAATs are available for confirmation if clinically indicated.
- The performance of this test was established based on the evaluation of a limited number of clinical specimens. The clinical performance has not been established with all circulating variants but is anticipated to be reflective of the prevalent variants in circulation at the time and location of the clinical evaluation. Performance at the time of testing may vary depending on the variants circulating, including newly emerging strains of SARS-CoV-2 and their prevalence, which change over time.

Conditions of Authorizations for Labs

cobas®ARS-CoV-2 & Influenza A/Bv2 Letter of Authorization, along with the authorized Fact Sheet for Healthcare Providers, the authorized Fact Sheet for Patients, and authorized labeling are available on the FDA website: <https://www.fda.gov/medicaldevices/coronavirusdiseases2019/covid-19-emergencyuseauthorizationsmedical-devices/in-vitro-diagnostics-euas>

However, to assist clinical laboratories using cobas®ARS-CoV-2 & Influenza A/B v2 (“this product” in the conditions below), the relevant Conditions of Authorization are listed below:

- A. Authorized laboratories using this product¹ must include with test result reports, all authorized Fact Sheets. Under exigent circumstances, other appropriate methods for disseminating these Fact Sheets may be used, which may include mass media.
- B. Authorized laboratories using this product must use this product as outlined in the authorized labeling. Deviations from the authorized procedures, including the authorized instruments, authorized extraction methods, authorized clinical specimen types, authorized control materials, authorized other ancillary reagents and authorized materials required to use this product are not permitted.
- C. Authorized laboratories that receive this product must notify the relevant public health authorities of their intent to run this product prior to initiating testing.
- D. Authorized laboratories using this product must have a process in place for reporting test results to healthcare providers and relevant public health authorities, as appropriate.
- E. Authorized laboratories must collect information on the performance of this product and report to DMD/OHT7/OPEQ/CDRH (via email: CDRH-EUA-Reporting@fda.hhs.gov) and Roche Diagnostics US Customer Technical Support (via telephone number 1-800-526-1247) any suspected occurrence of false positive or false negative results and significant deviations from the established performance characteristics of this product of which they become aware.
- F. All laboratory personnel using this product must be appropriately trained in real time RT-PCR techniques and instruments used in the cobas®ARS-CoV-2 & Influenza A/B v2, and use appropriate laboratory and personal protective equipment when handling this kit, and use this product in accordance with the authorized labeling.
- G. Roche Molecular Systems, authorized distributor(s) and authorized laboratories using this product must ensure that any records associated with this EUA are maintained until otherwise notified by FDA. Such records must be made available to FDA for inspection upon request.

¹ “This product” refers to cobas®ARS-CoV-2 & Influenza A/B v2. The letter of authorization refers to, “Laboratories certified under the Clinical Laboratory Improvement Amendments of 1988 (CLIA), 42 U.S.C. §263a, that meet requirements to perform moderate or high complexity tests” as “authorized laboratories.”

Non-clinical performance evaluation

cobas® SARS-CoV-2 & Influenza A/B v2 is an updated version of cobas® SARS-CoV-2 & Influenza A/B comprised of an influenza A dual target assay design which improves inclusivity of the test. The influenza B and SARS-CoV-2 assays of cobas® SARS-CoV-2 & Influenza A/B remained unchanged in cobas® SARS-CoV-2 & Influenza A/B v2.

Studies were conducted to demonstrate that the general performance of each target of the assay is unchanged and to demonstrate the effectiveness of the updated design of the influenza A target. The following key performance characteristics data were generated with either cobas® SARS-CoV-2 & Influenza A/B or cobas® SARS-CoV-2 & Influenza A/B v2.

Key performance characteristics

Analytical sensitivity (Limit of Detection)

The Limit of Detection (LoD) study determines the lowest detectable concentration of SARS-CoV-2, influenza A, and influenza B at which greater or equal to 95% of all (true positive) replicates test positive.

To determine the LoD, six cultured viruses—two each of influenza A and influenza B strains as well as the live and the heat-inactivated form of SARS-CoV-2 isolate from a US patient—were serially diluted in simulated clinical matrix to build two co-formulated target panels and three target single-formulated panels with one strain per virus. Seven to eight concentration levels, with two-fold serial dilutions between the levels, were prepared on three days and tested with a total of 63 replicates per concentration across three reagent lots for co-formulated panels and with a total of 21 replicates per concentration using one reagent lot for single-formulated panels. Table 13 to Table 16 summarize the established LoD values.

Table 13 Summary of LoD for influenza A determined with cobas® SARS-CoV-2 & Influenza A/B

Viral Strain	Kit lot	Panel	95% Probit [TCID ₅₀ /mL]	95% CI of Probit [TCID ₅₀ /mL]	Hit rate ≥95% [TCID ₅₀ /mL]	Mean Ct at ≥95% Hit rate
A/Kansas/14/2017 (H3N2)** Cat No 0810586CF Lot 323540	Lot 1	single-formulated	0.050	0.034 – 0.098	0.036	38.2
	Lot 1	co-formulated	0.12	0.073 – 0.28	0.071	36.6
	Lot 2	co-formulated	0.083	0.054 – 0.17	0.14	36.7
	Lot 3	co-formulated	0.062	0.040 – 0.14	0.071	37.0
	Lot 1-3	co-formulated	0.086	0.065 – 0.12	0.071	37.5
A/Brisbane/02/2018 (H1N1)*** Cat No 0810585CF Lot 323771	Lot 1	co-formulated	0.020****	0.013 – 0.048	0.026	37.4
	Lot 2	co-formulated	0.020	0.013 – 0.064	0.026	38.4
	Lot 3	co-formulated	0.025	0.016 – 0.059	0.026	38.1
	Lot 1-3	co-formulated	0.022	0.017 – 0.034	0.026	38.0

* LoD equivalency was demonstrated via performance studies with cobas® SARS-CoV-2 & Influenza A/B v2, refer to Table 17.

** Lot specific factor to convert TCID₅₀ into copy number was determined using NATrol™ Influenza A H3 Stock (Catalog# NATFLUAH3-STQ, Lot: 331079) material. 1 TCID₅₀/mL corresponds to 631 cp/mL.

*** Lot specific factor to convert TCID₅₀ into copy number was determined using NATrol™ Influenza A H1 Stock (Catalog# NATFLUAH1-STQ, Lot: 331080) material. 1 TCID₅₀/mL corresponds to 5811 cp/mL.

**** Claimed LoD was verified testing influenza A H1N1 pdm09 strains containing the C124A (GISAID: EPI_ISL_14387941), and the C124A plus G141A mutations in the M gene (GISAID: EPI_ISL_15803829) with cobas® SARS-CoV-2 & Influenza A/B v2.

Table 14 Summary of LoD for influenza B determined with cobas® SARS-CoV-2 & Influenza A/B*

Viral Strain	Kit lot	Panel	95% Probit [TCID ₅₀ /mL]	95% CI of Probit [TCID ₅₀ /mL]	Hit rate ≥95% [TCID ₅₀ /mL]	Mean Ct at ≥95% Hit rate
B/Phuket/3073/2013 (Yamagata lineage) Cat No 0810515CF Lot 320436	Lot 1	single-formulated	0.011	0.0076 – 0.023	0.017	35.4
	Lot 1	co-formulated	0.019	0.012 – 0.044	0.034	35.1
	Lot 2	co-formulated	0.016	0.0095 – 0.050	0.017	35.4
	Lot 3	co-formulated	0.019	0.010 – 0.084	0.017	35.3
	Lot 1-3	co-formulated	0.017	0.012 – 0.026	0.017	35.3
B/Colorado/06/2017 (Victoria lineage) Cat No 0810573CF Lot 323459	Lot 1	co-formulated	0.027	0.017 – 0.065	0.026	34.9
	Lot 2	co-formulated	0.032	0.019 – 0.084	0.053	34.5
	Lot 3	co-formulated	0.019	0.012 – 0.050	0.026	35.0
	Lot 1-3	co-formulated	0.026	0.019 – 0.040	0.026	34.9

*LoD equivalency was demonstrated via performance studies with cobas® SARS-CoV-2 & Influenza A/B2, refer to table 17.

Table 15 Summary of LoD for SARS-CoV-2 determined with cobas® SARS-CoV-2 & Influenza A/B

Viral Strain	Kit lot	Panel	95% Probit [TCID ₅₀ /mL]	95% CI of Probit [TCID ₅₀ /mL]	Hit rate ≥95% [TCID ₅₀ /mL]	Mean Ct at ≥95% Hit rate
USA-WA1/2020 heat-inactivated Cat No 0810587CFHI Lot 324045	Lot 1	single-formulated	0.068	0.044– 0.15	0.058	36.9
	Lot 1	co-formulated	0.14	0.086– 0.35	0.12	36.3
	Lot 2	co-formulated	0.13	0.083– 0.26	0.12	36.4
	Lot 3	co-formulated	0.10	0.065– 0.25	0.12	35.9
	Lot 1-3	co-formulated	0.13	0.094– 0.19	0.12	36.2
USA-WA1/2020 infectious culture Cat No NR52281 Lot 70033175*	Lot 1	co-formulated	0.0081	0.0041– 0.049	0.0079	36.2
	Lot 2	co-formulated	0.0071	0.0044– 0.018	0.0079	36.2
	Lot 3	co-formulated	0.0052	0.0032– 0.013	0.0079	35.9
	Lot 1-3	co-formulated	0.0063	0.0046– 0.010	0.0079	36.1

* LoD equivalency was demonstrated via performance studies with cobas® SARS-CoV-2 & Influenza A/B2, refer to table 17.

**Based on the information provided in the Certificate of Analysis from the Vendor, 1 mL is equal to 7,393 genome equivalents per mL.

Table 16 Summary of LoD for pan-Sarbecovirus determined with cobas® SARS-CoV-2 & Influenza A/B*

Viral Strain	Kit lot	Panel	95% Probit [TCID ₅₀ /mL]	95% CI of Probit [TCID ₅₀ /mL]	Hit rate ≥95% [TCID ₅₀ /mL]	Mean Ct at ≥95% Hit rate
USA-WA1/2020 heat-inactivated Cat No 0810587CFHI Lot 324045	Lot 1	single-formulated	0.14	0.08 – 0.37	0.12	35.6
	Lot 1	co-formulated	0.28	0.17 – 0.67	0.55	34.5
	Lot 2	co-formulated	0.23	0.14 – 0.49	0.23	35.1
	Lot 3	co-formulated	0.18	0.11 – 0.37	0.23	34.8
	Lot 1-3	co-formulated	0.23	0.17 – 0.34	0.55	34.2
USA-WA1/2020 infectious culture Cat No NR-52281 Lot 70033175**	Lot 1	co-formulated	0.0090	0.0057 – 0.020	0.016	34.6
	Lot 2	co-formulated	0.0076	0.0049 – 0.016	0.016	34.7
	Lot 3	co-formulated	0.0080	0.0053 – 0.017	0.0079	35.3
	Lot 1-3	co-formulated	0.0082	0.0062 – 0.012	0.016	34.7

* LoD equivalency was demonstrated via performance with cobas® SARS-CoV-2 & Influenza A/Bv2 refer to Table 17

**Based on the information provided in the Certificate of Analysis from the vendor, the LoD is equal to 7,393 genome equivalents by ddPCR.

LoD Equivalency

In order to show LoD equivalency of cobas® SARS-CoV-2 & Influenza A/B v2 compared to the previous version, the cobas® SARS-CoV-2 & Influenza A/B co-formulated panel, consisting of influenza A, influenza B and SARS-CoV-2, was tested with 21 replicates per panel member assay version to compare the LoD in simulated clinical matrix in Universal Transport Media (UTM™). One reagent lot was used of each assay version. This summary of the results is shown in Table 17.

Table 17 LoD Equivalency between cobas® SARS-CoV-2 & Influenza A/B and cobas® SARS-CoV-2 & Influenza A/Bv2

Target	Viral Strain	cobas® SARS-CoV-2 & Influenza A/B 95% Probit (TCID ₅₀ /mL)	cobas® SARS-CoV-2 & Influenza A/B v2 95% Probit Predicted 95% (TCID ₅₀ /mL)	cobas® SARS- CoV-2 & Influenza A/B Lower 95% to Upper 95% Confidence Interval (TCID ₅₀ /mL)	cobas® SARS-CoV-2 & Influenza A/B v2 Lower 95% to Upper 95% Confidence Interval (TCID ₅₀ /mL)
Influenza A	A/Kansas/14/2017 (H3N2) Cat No 0810586CF Lot 323540	0.077	0.040	0.047- 0.190	0.026- 0.088
SARSCoV-2	USA-WA1/2020 heat-inactivated Cat No 0810587CFHI	0.124	0.053	0.077- 0.296	0.036- 0.116
Pan Sarbecovirus	Lot 324045	0.173	0.078	0.112- 0.356	0.005- 0.156
Influenza B	Phuket/3073/2013 (Yamagata lineage) Cat No 0810515CF Lot 320436	0.023	0.011	0.014- 0.057	0.007- 0.026

Inclusivity

The inclusivity for the detection of influenza A was confirmed by testing thirteen influenza A strains with cobas® SARS-CoV-2 & Influenza A/B v2. All strains tested showed 100% rate at the concentrations indicated in Table 18.

Table 18 Summary of inclusivity for influenza A tested with cobas® SARS-CoV-2 & Influenza A/B v2

Viral Target	Strain	Catalog Number	Concentration
Influenza A	A/Canada/6294/09 (H1N1)	0810109CFJ	1.60E02 TCID ₅₀ /mL
	A/California/07/09 (H1N1)	0810165CF	2.31E01 TCID ₅₀ /mL
	A/Mexico/4108/09 (H1N1)	0810166CF	8.71E02 TCID ₅₀ /mL
	A/Singapore/63/04 (H1N1)	0810246CF	1.69E02 TCID ₅₀ /mL
	A/Michigan/45/15 (H1N1)	0810538CF	9.28E02 TCID ₅₀ /mL
	A/California/04/09 (H1N1)	VR-1805	6.41E+00 TCID ₅₀ /mL
	A/England/224020815/2020 (H1N1)*	n/a	4.52E+02 cp/mL***
	A/England/221740513/2020 (H1N1)**	n/a	4.52E+02 cp/mL***
	A/Perth/16/09 (H3N2)	0810251CF	6.30E02 TCID ₅₀ /mL
	A/Wisconsin/67/05 (H3N2)	0810252CF	4.37E01 TCID ₅₀ /mL
	A/Switzerland/9715293/13 (H3N2)	0810511CF	3.63E02 TCID ₅₀ /mL
	A/HongKong/4801/14 (H3N2)	0810526CF	3.57E01 TCID ₅₀ /mL
	A/Texas/50/12 (H3N2)	0810238CF	2.7E+00 TCID ₅₀ /mL

* GISAID IDEPI_ISL_15803829 containing the C124A and G141A mutations in gene 5 (Strain is not commercially available)

** GISAID ID EPI_ISL_14387944 containing the C124A mutation in the M2 gene (Strain is not commercially available)

*** Copy number was determined using Abbott™ Influenza A H5N1 Stock (Catalog # ATFLUAH1-STQ, Lot: 331080) material

The inclusivity for the detection of influenza B and SARS-CoV-2 was confirmed by testing five influenza B and three SARS-CoV-2 strains with cobas® SARS-CoV-2 & Influenza A/B. The lowest target analyte which all four tested replicates were positive are reported (Table 19 and Table 20).

Table 19 Summary of inclusivity influenza B tested with cobas® SARS-CoV-2 & Influenza A/B

Viral Target	Strain	Catalog Number	Lot Number	Lowest Concentration Detected
Influenza B	B/Brisbane/60/2008 (Victoria lineage)	0810254CF	313257 (sublot: 513438)	0.002 TCID ₅₀ /mL
	B/Utah/9/14 (Yamagata lineage)	0810516CF	317295 (sublot: 527062)	0.017 TCID ₅₀ /mL
	B/Alabama/2/17 (Victoria lineage)	0810572CF	322548	0.0064 TCID ₅₀ /mL
	B/Florida/78/2015 (Victoria Lineage)	VR-1931	70020870	0.076 TCID ₅₀ /mL
	B/Wisconsin/1/2010 (Yamagata Lineage)	VR-1883	70012127	0.070 CEID ₅₀ /mL

Table 20 Summary of inclusivity for SARS-CoV-2 tested with cobas® SARS-CoV-2 & Influenza A/B

Viral Target	Strain	Catalog Number	Lot Number	Lowest Concentration Detected
SARS-CoV-2	BetaCoV/France/IDF0372/2020	014V-03890	Not available	0.038 PFU/mL
	BetaCoV/Munich/BavPat1/2020	026V-03883	Not available	0.0036 PFU/mL
	2019-nCoV/Italy-INMII	008V-03893	Not available	0.062 TCID ₅₀ /mL

Further, cobas® SARS-CoV-2 & Influenza A/B was shown to be inclusive for the CDC Human Influenza Virus Panel (2020) (Cat. Number VP2020, Lot Number 2003300). The lowest concentration where at least one out of five replicates was positive is reported as the minimum reactive concentration in Table 21.

Table 21 Summary of CDC Human Influenza Virus Panel (2020) tested with cobas® SARS-CoV-2 & Influenza A/B

Virus	Strain	Minimum Reactive Concentration [EID ₅₀ /mL]
Influenza A	A/Perth/16/2009 (H3N2)	2.62E+00
	A/Hong Kong/2671/2019 (H3N2)	8.29E02
	A/Christ Church/16/2010 (H1N1 pdm)	2.08E+01
	A/Guangdong-maonan/1536/2019 (H1N1 pdm)	6.00E01
Influenza B	B/Michigan/09/2011	1.30E02
	B/Washington/02/2019	2.08E+00
	B/Texas/81/2016	6.54E02
	B/Phuket/3073/2013	2.08E+01

Precision (repeatability)

Within -laboratory precision was examined using a panel composed of spiked influenza A (A/Kansas/14/2017), influenza B (B/Phuket/3073/2013) and SARS-CoV-2 (USA-WA1/2020, heat inactivated) cultures diluted in simulated clinical matrix in UTM-RT®. Sources of variability were examined with a panel consisting of three concentration levels, using three lots of cobas® SARS-CoV-2 & Influenza A/B reagents and two instruments over a time course of 15 days for a total of 30 runs. A description of the precision panel and the observed positivity rates are shown in Table 22. All negative panel members tested negative throughout the study. Analysis of standard deviation and percent coefficient of variation (CV) of the Ct values from tests performed on positive panel members (n=20) yielded overall CV percentages ranging from 12% to 5.1% for influenza A, influenza B, and SARS-CoV-2.

Table 22 Summary of within laboratory precision

Target Concentration	N Tested	N Positive	Positivity Rate	95% Confidence Interval	
				Lower Limit	Upper Limit
Influenza A					
Negative	90	0	0%	0%	4.1%
Weak Positive ~ 0.3 x LoD (0.043TCID ₅₀ /mL)	90	87	96.7%	90.7%	98.9%
Low Positive ~ 1 x LoD (0.14TCID ₅₀ /mL)	90	90	100%	95.9%	100%
Moderate Positive ~ 3 x LoD (0.43TCID ₅₀ /mL)	90	90	100%	95.9%	100%
Influenza B					
Negative	90	0	0%	0%	4.1%
Weak Positive ~ 0.3 x LoD (0.010TCID ₅₀ /mL)	90	81	90.0%	82.1%	94.7%
Low Positive ~ 1 x LoD (0.034TCID ₅₀ /mL)	90	90	100%	95.9%	100%
Moderate Positive ~ 3 x LoD (0.10TCID ₅₀ /mL)	90	90	100%	95.9%	100%
SARS-CoV-2					
Negative	90	0	0%	0%	4.1%
Weak Positive ~ 0.3 x LoD (0.035 TCID ₅₀ /mL)	90	83	92.2%	84.8%	96.2%
Low Positive ~ 1 x LoD (0.12TCID ₅₀ /mL)	90	87	96.7%	90.7%	98.9%
Moderate Positive ~ 3 x LoD (0.35TCID ₅₀ /mL)	90	90	100%	95.9%	100%
pan-Sarbecovirus					
Negative	90	0	0%	0%	4.1%
Weak Positive ~ 0.06 x LoD (0.035TCID ₅₀ /mL)	90	73	81.1%	71.8%	87.9%
Low Positive ~ 0.2x LoD (0.12TCID ₅₀ /mL)	90	87	96.7%	90.7%	98.9%
Moderate Positive ~ 0.6x LoD (0.35TCID ₅₀ /mL)	90	90	100%	95.9%	100%

Table 23 Overall mean, standard deviation, and percent coefficient of variation for Ct values by positive panel member

Target Concentration	Positivity Rate	Mean Ct	Between instrument		Between lot		Between day		Between run		Within run		Total	
			SD	CV%	SD	CV%	SD	CV%	SD	CV%	SD	CV%	SD	CV%
Influenza A														
Weak Positive ~ 0.3 x LoD (0.042 TCID ₅₀ /mL)	96.7%	383	0.00	0.0	0.29	0.8	0.43	1.1	0.00	0.0	1.90	5.0	1.97	5.1
Low Positive ~ 1 x LoD (0.14 TCID ₅₀ /mL)	100%	35.7	0.00	0.0	0.00	0.0	0.19	0.5	0.15	0.4	0.90	2.5	0.98	2.6
Moderate Positive ~ 3 x LoD (0.42 TCID ₅₀ /mL)	100%	34.4	0.11	0.3	0.00	0.0	0.11	0.3	0.00	0.0	0.43	1.2	0.46	1.3
Influenza B														
Weak Positive ~ 0.3 x LoD (0.010 TCID ₅₀ /mL)	90.0%	356	0.11	0.3	0.00	0.0	0.236	0.6	0.09	0.3	0.62	1.7	0.67	1.9
Low Positive ~ 1 x LoD (0.034 TCID ₅₀ /mL)	100%	34.7	0.00	0.0	0.00	0.0	0.19	0.5	0.21	0.6	0.51	1.5	0.58	1.7
Moderate Positive ~ 3 x LoD (0.10 TCID ₅₀ /mL)	100%	33.8	0.07	0.2	0.00	0.0	0.17	0.5	0.00	0.0	0.82	2.4	0.84	2.5
SARS-CoV-2														
Weak Positive ~ 0.3 x LoD (0.035 TCID ₅₀ /mL)	92.2%	366	0.00	0.0	0.00	0.0	0.32	0.9	0.07	0.2	0.60	1.6	0.68	1.9
Low Positive ~ 1 x LoD (0.12 TCID ₅₀ /mL)	96.7%	35.7	0.06	0.2	0.07	0.2	0.00	0.0	0.05	0.1	0.40	1.1	0.42	1.2
Moderate Positive ~ 3 x LoD (0.35 TCID ₅₀ /mL)	100%	34.6	0.17	0.5	0.00	0.0	0.19	0.6	0.00	0.0	0.57	1.7	0.63	1.8
pan-Sarbecovirus														
Weak Positive ~ 0.06 x LoD (0.035 TCID ₅₀ /mL)	81.1%	35.8	0.00	0.0	0.00	0.0	0.16	0.4	0.11	0.3	0.63	1.8	0.66	1.82.0
Low Positive ~ 0.2 x LoD (0.12 TCID ₅₀ /mL)	96.7%	34.9	0.00	0.0	0.04	0.2	0.00	0.0	0.00	0.0	0.52	1.5	0.52	1.5
Moderate Positive ~ 0.6 x LoD (0.35 TCID ₅₀ /mL)	100%	33.9	0.13	0.4	0.00	0.0	0.10	0.3	0.00	0.0	0.54	1.6	0.57	1.7

Analytical specificity (cross-reactivity and microbial interference)

A panel of 40 viruses, bacteria, and fungi (including those commonly found in respiratory tract) plus pooled human nasal wash was tested with cobas® SARS-CoV-2 & Influenza A/B v2 to assess analytical specificity. The organisms listed in Table 24 were spiked at concentrations of 1×10^6 units/mL for viruses and 1×10^8 units/mL for other organisms, unless otherwise noted. Testing was performed with each potential interfering organism in the absence and presence of influenza A, influenza B, and SARS-CoV-2 target (spiked at 420.10 and 0.36 TCID₅₀/mL, respectively). None of the organisms interfered with the test performance by generating false positive results. Testing of SARS-CoV-1 generated an expected parainfluenza virus positive result. Detection of influenza A, influenza B, and SARS-CoV-2 targets was not affected in the presence of the organisms tested. Potential cross-reactivity of influenza C, *Leptospira interrogans*, *Pneumocystis jirovecii*, *Chlamydia psittaci*, *Bacillus anthracis*, and *Coxiella burnetii* was evaluated in silico. Based on the in silico analyses, selected organisms are highly unlikely to interfere with the performance of cobas® SARS-CoV-2 & Influenza A/B v2.

Table 24 Microorganisms tested for analytical specificity/cross reactivity

Microorganism	Concentration
Adenovirus (AdV-1)	1.0E+05 TCID ₅₀ /mL
<i>Bordetella pertussis</i>	1.0E+06 CFU/mL
<i>Candida albicans</i>	1.0E+06 CFU/mL
<i>Chlamydia pneumoniae</i>	7.9E+04 TCID ₅₀ /mL
<i>Corynebacterium diphtheriae</i>	1.0E+06 CFU/mL
Cytomegalovirus	1.0E+05 IU/mL
Enterovirus (EV68)	1.0E+05 TCID ₅₀ /mL
Epstein Barvirus	1.0E+05 cp/mL
<i>Escherichia coli</i>	1.0E+06 CFU/mL
<i>Haemophilus influenzae</i>	1.0E+06 CFU/mL
Human coronavirus 229E	1.0E+05 TCID ₅₀ /mL
Human coronavirus HKU1	6.9E+04 genome cp/mL
Human coronavirus NL63	7.0E+03 TCID ₅₀ /mL
Human coronavirus OC43	1.0E+05 TCID ₅₀ /mL
Human Metapneumovirus	1.0E+05 TCID ₅₀ /mL
<i>Lactobacillus acidophilus</i>	5.0E+05 CFU/mL
<i>Legionella pneumophila</i>	1.0E+06 CFU/mL
<i>Legionella longbeachae</i>	1.0E+06 CFU/mL
Measles virus	1.0E+05 TCID ₅₀ /mL
MERS coronavirus	1.0E+05 cp/mL
<i>Moraxella catarrhalis</i>	1.0E+06 CFU/mL
Mumps Virus	1.0E+05 U/mL
<i>Mycobacterium bovis</i>	1.0E+05 CFU/mL
<i>Mycoplasma pneumoniae</i>	1.0E+06 CCU/mL
<i>Neisseria elongata</i>	1.0E+06 CFU/mL
<i>Neisseria meningitidis</i>	1.0E+06 CFU/mL
Parainfluenza virus 1	1.0E+05 TCID ₅₀ /mL
Parainfluenza virus 2	1.0E+05 TCID ₅₀ /mL
Parainfluenza virus 3	1.0E+05 TCID ₅₀ /mL
Parainfluenza virus 4	1.0E+05 TCID ₅₀ /mL
Parechovirus	1.0E+05 U/mL
<i>Pseudomonas aeruginosa</i>	1.0E+06 CFU/mL
Respiratory syncytial virus	1.0E+05 PFU/mL
Human Rhinovirus	1.0E+05 PFU/mL
SARS coronavirus (SARSCoV-1)	1.0E+07 PFU/mL
<i>Staphylococcus aureus</i>	1.0E+06 CFU/mL
<i>Staphylococcus epidermidis</i>	1.0E+06 CFU/mL
<i>Streptococcus salivarius</i>	1.0E+06 CFU/mL
<i>Streptococcus pneumoniae</i>	1.0E+06 CFU/mL
<i>Streptococcus pyogenes</i>	1.0E+06 CFU/mL

Co-infection (competitive interference)

To assess potential competitive interference between influenza A, influenza B and SARS-CoV-2, samples were tested using cobas® SARS-CoV-2 & Influenza A/B v2 in replicates of 4 where low concentrations (0.42, 0.10 and 0.36 TCID₅₀/mL for influenza A, influenza B, and SARS-CoV-2, respectively) of any two targets were mixed with very high (1.0E+05 Units/mL) concentrations of the third target. None of the targets present at very high concentration interfered with the detection of low levels of the other two targets

Collection media equivalence

Equivalence between different collection media (UTM-RT® cobas® CR Media and saline) was evaluated using one strain each for influenza A (A/Kansas/14/2017 (H3N2)), influenza B (B/Phuket/3073/13 (Yamagata lineage)) and SARS-CoV-2 (USA-WA1/2020, heat inactivated culture). Testing was performed using cobas® SARS-CoV-2 & Influenza A/B. Virus cultures were formulated to a target concentration of approximately 1.2D into simulated clinical matrix formulated either in Universal Transport Media (UTM-RT® cobas® CR Media (CPM) or in 0.9% physiological saline. A total of 21 replicates were tested for each collection media type. All replicates tested were positive in all simulated matrices for influenza A and influenza B. For SARS-CoV-2, positivity rates were 100% for both UTM-RT® and CPM and 95.2% for saline.

Clinical performance evaluation

First, the clinical performance was evaluated at one external site using archived nasopharyngeal swab (NPS) samples from patients with signs and symptoms of a respiratory infection, collected by UTM-RT® or UVT between 2014 and 2020 for the SARS-CoV-2 and influenza B component with cobas® SARS-CoV-2 & Influenza A/B. Clinical samples were collected by qualified personnel according to the package insert of the collection device.

This clinical evaluation study included a total of 349 NPS samples, 57 of which were longitudinal samples from COVID-19 patients. Two FDA cleared molecular tests were utilized as the comparators for assessment of performance of cobas® SARS-CoV-2 & Influenza A/B: one for SARS-CoV-2 and one for influenza B. One of the 349 NPS samples did not have a valid comparator SARS-CoV-2 result and five of the 349 NPS samples did not have valid comparator influenza A/B results, therefore, were excluded from the performance calculations for SARS-CoV-2 and influenza B, respectively.

In a subsequent clinical evaluation, the performance of cobas® SARS-CoV-2 & Influenza A/B v2 was assessed for the influenza A component at one internal site using archived NPS and NS samples from patients with signs and symptoms of a respiratory infection, collected by UTM-RT® or UVT in 2022-2023. Clinical samples were collected by qualified personnel according to the package insert of the collection device. This clinical evaluation study included a total of 75 NPS and 75 nasal swab (NS) evaluable samples. An FDA cleared molecular assay was utilized as the comparator test for assessment of performance of the assay for the influenza A component.

As shown in Table 25, cobas® SARS-CoV-2 & Influenza A/B and cobas® SARS-CoV-2 & Influenza A/B v2 demonstrated high percent agreement with the comparators for the detection of SARS-CoV-2, influenza A and influenza B.

Table 25 Comparison of cobas® SARS-CoV-2 & Influenza A/B and cobas® SARS-CoV-2 & Influenza A/B v2 with the comparator tests.

Virus	Number of Samples	Test Results				Agreement Statistics		
		Concordant Positive (N)	Discordant Positive (N)	Concordant Negative (N)	Discordant Negative (N)	Agreement Parameter	Percent Agreement (%)	95% CI (LCL, UCL)*
SARSCoV-2 [#]	348	53	6	287	2	PPA	96.4%	(87.7%, 99.0%)
						NPA	98.0%	(95.6%, 99.1%)
Influenza A [†]	150 ^a	50	0	100	0	PPA	100.0%	(92.9%, 100.0%)
						NPA	100.0%	(96.3%, 100.0%)
Influenza B	344	37	1	306	0	PPA	100.0%	(90.6%, 100.0%)
						NPA	99.7%	(98.2%, 99.9%)

PPA = Positive Percent Agreement

NPA = Negative Percent Agreement

CI = confidence interval; LCL = Lower confidence Limit; UCL = Upper confidence Limit

*Confidence interval is calculated using Wilson's Score method

[#]A positive result is defined as detection of either of the two SARS-CoV-2 or pan-Sarbecovirus target of the assay

[†]Including six H1N1pdm09 positive samples containing the C124A and G141A mutations in the M gene

^aNPS and NS samples combined

Discordant results between the cobas® SARS-CoV-2 & Influenza A/B assay and the comparator methods were observed for 9 samples. Of these, 8 were longitudinal samples with discordant results for SARS-CoV-2 that showed late Ct values (between 35-43). The candidate test detected an additional influenza B virus positive sample compared to the comparator. Post-PCR analysis of the amplicon from the discordant samples confirmed the presence of SARS-CoV-2 but not influenza B.

cobas® SARS-CoV-2 & Influenza A/B was further evaluated in a prospective clinical study comparison with FDA cleared tests with fresh NPS and NS clinical samples. A total of 604 NPS and a total of 604 self-collected and 304 health care worker collected valid results by the comparator method were evaluated. The PPA for influenza B was 100% (95% Score CI of 99% to 100%) for NPS samples and 99.8% (95% Score CI of 99.1% to 100%) for NS samples. The PPA for NPS samples was not calculable for influenza B as there were no positive samples by the comparator method. The PPA for NS samples for influenza B was 0% (0/1) with a 95% Score CI of 0% to 4.9%, as there was one positive sample by the comparator method and negative on cobas® SARS-CoV-2 & Influenza A/B.

Discordant results between the cobas® SARS-CoV-2 & Influenza A/B assay and the comparator method were observed for 2 samples. One NS sample tested positive for influenza B by the comparator method and was negative on cobas® SARS-CoV-2 & Influenza A/B. Sequencing analysis of the sample did not confirm the presence of influenza B. One NS sample tested positive for influenza B on cobas® SARS-CoV-2 & Influenza A/B and was negative on the comparator method. PCR analysis of the amplicon of this discordant NS sample did not confirm the presence of influenza B.

Additional information

Key test features





















































Sample type	Nasopharyngeal swab samples collected in the Copan UTM-RT® System or the BD™UVT System Nasal swab samples collected in the Copan UTM-RT® System, the BD™UVT System, the cobas ® PCR Media, and 0.9% physiological saline
Minimum amount of sample required	0.6 mL or 1.0 mL*
Sample processing volume	0.4 mL
Test duration	Results are available within less than 3.5 hours after loading the sample on the system.

*Dead volume of 0.2 mL is identified for **cobas**® **omni** Secondary tubes. Dead volume of 0.6 mL is identified for the **cobas**® PCR Media primary tubes. Other tubes compatible with **cobas**® 6800/8800 System (consult User Assistance Guide) may have different dead volume and require more or less minimum volume.

Symbols

The following symbols are used in labeling for Roche PCR diagnostic products.

Table 26 Symbols used in labeling for Roche PCR diagnostic products

 Age/DOB	Age or Date of Birth		Device not for near-patient testing		QS IU per PCR reaction, use the QS International Units (IU) per PCR reaction in calculation of the results.
	Ancillary Software		Device not for self-testing		Serial number
	Assigned Range (copies/mL)		Distributor <i>(Note: The applicable country/region may be designated beneath the symbol)</i>		Site
	Assigned Range (IU/mL)		Do not re-use		Standard Procedure
	Authorized representative in the European Community		Female		Sterilized using ethylene oxide
	Barcode Data Sheet		For IVD performance evaluation only		Store in dark
	Batch code		Global Trade Item Number		Temperature limit
	Biological risks		Importer		Test Definition File
	Catalogue number		In vitro diagnostic medical device		This way up
	CE marking of conformity; this device is in conformity with the applicable requirements for CE marking of an in vitro diagnostic medical device		Lower Limit of Assigned Range		Ultrasensitive Procedure
	Male		Manufacturer		Unique Device Identifier
	Collect date		Negative control		Upper Limit of Assigned Range
	Consult instructions for use		Non-sterile		Urine Fill Line
	Contains sufficient for <n> tests		Patient Name		US Only: Federal law restricts this device to sale by or on the order of a physician.
	Content of kit		Patient number		Use-by date
	Control		Peel here		Positive control
	Date of manufacture		QS copies per PCR reaction, use the QS copies per PCR reaction in calculation of the results.		
	Device for near-patient testing				
	Device for self-testing				

Technical support

For technical support (assistance) please ~~contact~~ reach out to your local affiliate:
https://www.roche.com/about/business/roche_worldwide.htm

Manufacturer and distributor

Table 27 Manufacturer and distributor



RocheMolecular Systems, Inc.
1080 US Highway 202 South
Branchburg, NJ 08876 USA
www.roche.com

Made in USA

Distributed by Roche Diagnostics
9115 Hague Road
Indianapolis, IN 46250-0457 USA
(For Technical Assistance call the
Roche Response Center
toll-free: 1-800-526-1247)

Trademarks and patents

See <https://diagnostics.roche.com/us/en/about/patents>

Copyright

©2023 Roche Molecular Systems, Inc.



References

1. Center for Disease Control and Prevention. Biosafety in Microbiological and Biomedical Laboratories, 5th ed. U.S. Department of Health and Human Services, Public Health Service, Centers for Disease Control and Prevention, National Institutes of Health Publication No. (CDC) 2-1112, revised December 2009.
2. Clinical and Laboratory Standards Institute (CLSI). Protection of laboratory workers from occupationally acquired infections. Approved Guideline-Fourth Edition. CLSI Document M29-A4:Wayne, PA;CLSI, 2014.

Document revision

Document Revision Information	
Doc Rev. 1.0 10/2023	First Publishing

cobas[®] SARS-CoV-2 & Influenza A/B v2



Rx Only

KIT **LOT**



For USA: Emergency Use Authorization only

USA



cobas[®] SCoV2-FluA/B ASAP Version 12.3.0 or higher

**cobas[®] 6800/8800 System Software
Version 1.4 or higher**

This product has not been FDA cleared or approved, but has been authorized for emergency use by FDA under an EUA for use by authorized laboratories.

This product has been authorized only for the detection and differentiation of nucleic acid from SARS-CoV-2, influenza A, and influenza B, not for any other viruses or pathogens.

The emergency use of this product is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of in vitro diagnostics for detection and/or diagnosis of COVID-19 under Section 564(b)(1) of the Federal Food, Drug, and Cosmetic Act, 21 U.S.C. § 360bbb-3(b)(1), unless the declaration is terminated or authorization is revoked sooner.

Non-USA



**cobas[®] 5800 SCoV2-FluA/B ASAP
Version 1.2.0 or higher**

cobas[®] 5800 Software Version 1.0 or higher



cobas[®] SCoV2-FluA/B ASAP Version 12.3.0 or higher

**cobas[®] 6800/8800 System Software
Version 1.4 or higher**

USA



website: <http://e-labdoc.roche.com>

Product No.: 10033401190

10033479001-01 Doc Rev. 1.0

Please contact your local Roche representative at 1-800-526-1247 if you require a printed copy free of charge or need technical support to access the package insert.

Non-USA



website: <http://e-labdoc.roche.com>

Method Sheet Catalog No.: 10033401190 Doc Rev 1.0

Please contact your local Roche representative if you require a printed copy free of charge or need technical support to access the package insert. / Bei Ihrer zuständigen Roche-Vertretung erhalten Sie einen kostenfreien Ausdruck oder technische Unterstützung für den Zugriff auf die Packungsbeilage. / Veuillez contacter votre représentant Roche local pour obtenir un exemplaire papier gratuit ou une assistance technique pour accéder à la notice. / Contattare il rappresentante Roche locale per ottenere gratuitamente una copia stampata o richiedere istruzioni per reperire il foglio illustrativo. / Póngase en contacto con su representante local de Roche si necesita una copia impresa gratuita o ayuda del servicio técnico para acceder al boletín técnico. / Se desejar uma cópia impressa gratuita ou necessitar de assistência técnica para aceder ao folheto informativo, entre em contacto com o representante local da Roche. / Kontakt den lokale Roche-repræsentant, hvis du ønsker en gratis skriftlig kopi eller har brug for teknisk support for at få adgang til indlægssedlen. / Kontakta din Roche-representant om du vill ha en pappersversion kostnadsfritt eller om du behöver teknisk support för att komma åt bipacksedeln.



Roche Molecular Systems, Inc.
1080 US Highway 202 South
Branchburg, NJ 08876 USA

10017244001-01