

LUMISIGHT™ and Lumicell™ Direct Visualization System (DVS) as Adjunct to Standard of Care to Identify Residual Cancer Within the Lumpectomy Cavity

March 5, 2024

Lumicell

Medical Imaging Drugs Advisory Committee (MIDAC)



Introduction

Jorge Ferrer, PhD

Chief Scientific Officer
Lumicell

LUM System Is Real-Time, Intracavity, Fluorescence-Guided Imaging as Adjunct to Standard of Care (SoC)

LUM System

LUMISIGHT



Find

Optical imaging agent that produces fluorescence signal at sites of residual cancer

Lumicell Direct Visualization System (DVS)



+



Detect

Hand-held imaging probe inserted into breast cavity to identify residual cancer

Guide

Real-time cancer detection software guides surgeon to remove residual cancer

LUM System Developed to Fill Important Unmet Clinical Need


Lumpectomy meant to be minimally invasive alternative to mastectomy

Surgeons have no way to see full extent of cancer inside cavity

19% of negative margins have residual cancer remaining¹
9% – 36% second surgeries^{2,3}

Clear unmet need for real-time, intracavity assessment to more effectively determine extent of tumor for more complete resection


LUM System in Action



LUM System Clinical Development Program in Breast Cancer and Cardiovascular Safety

Study ID	Phase / Design	Patients Injected with LUMISIGHT
DUK1-12-137	<p>Phase 1</p> <p>Single site, nonrandomized, open label trial</p>	<p>15</p> <p>(3 breast, 12 sarcoma)</p> <p>(6 at 0.5 mg/kg; 6 at 1 mg/kg; 3 at 1.5 mg/kg)</p>
CLP00201	<p>Phase 1 Cardiovascular Safety Trial, healthy volunteers</p> <p>Single site, randomized, double-blind, placebo controlled, dose-escalation</p>	<p>24</p> <p>32 enrolled</p> <p>(8 in placebo-controlled arm)</p>
LUM-015/2.6-001	<p>Feasibility Phase A</p> <p>Single site, nonrandomized, open label trial</p>	<p>10</p> <p>(5 at 0.5 mg/kg; 5 at 1 mg/kg)</p>
LUM-015/2.6-001	<p>Feasibility Phase B</p> <p>Single site, nonrandomized, open label trial</p>	<p>45</p> <p>(1 mg/kg)</p>
CL0006	<p>Feasibility Phase C</p> <p>Multicenter, nonrandomized, open label trial</p>	<p>234</p> <p>(1 mg/kg)</p>
CLP0008	<p>Feasibility in Patients Receiving Neo-Adjuvant Therapy</p> <p>Multicenter, randomized, blinded trial</p>	<p>12</p> <p>(1 mg/kg)</p>
CL0007	<p>Pivotal Study</p> <p>Multicenter, 2-arm, randomized, blinded trial</p>	<p>406</p> <p>(1 mg/kg)</p>

Data Supporting Efficacy and Safety of LUM System Published in NEJM Evidence

Study ID	Study Title	Study Design
DUK1-12-137	 <p data-bbox="1709 435 2028 468">Published April 27, 2023</p> <p data-bbox="1635 525 2028 558">DOI: 10.1056/EVIDoa2200333</p>	Breast Cancer (Sarcoma) at 1.0 mg/kg; g/kg)
CLP00201	<p data-bbox="471 646 715 672">ORIGINAL ARTICLE</p> <p data-bbox="471 701 1939 843">Intraoperative Fluorescence Guidance for Breast Cancer Lumpectomy Surgery</p>	led controlled arm)
LUM-015/2.6-001	<p data-bbox="471 882 2028 1125">Barbara L. Smith, M.D., Ph.D.,¹ Kelly K. Hunt, M.D.,² David Carr, M.D.,³ Peter W. Blumencranz, M.D.,⁴ E. Shelley Hwang, M.D., M.P.H.,⁵ Michele A. Gadd, M.D.,¹ Kimberly Stone, M.D.,⁶ Donna L. Dyess, M.D.,⁷ Daleela Dodge, M.D.,⁸ Stephanie Valente, D.O.,⁹ Nayana Dekhne, M.D.,¹⁰ Patricia Clark, M.D.,¹¹ Marie Catherine Lee, M.D.,¹² Laila Samiian, M.D.,¹³ Beth-Anne Lesnikoski, M.D.,¹³ Lynne Clark, M.D.,¹⁴ Kate Porta Smith, M.P.H., C.C.R.P.,¹⁵ Manna Chang, Ph.D.,¹⁵ Daniel K. Harris, Ph.D.,¹⁵ Brian Schlossberg, Ph.D.,¹⁵ Jorge Ferrer, Ph.D.,¹⁵ Irene L. Wapnir, M.D.,⁶ for the INSITE Study Team*</p>	at 1.0 mg/kg)
LUM-015/2.6-001		g)
CL0006		g)
CLP0008		g)
CL0007	<p data-bbox="978 1232 1212 1265">Pivotal Study</p> <p data-bbox="705 1275 1485 1308">Multicenter, 2-arm, randomized, blinded trial</p>	<p data-bbox="2135 1232 2198 1265">406</p> <p data-bbox="2084 1275 2249 1308">(1 mg/kg)</p>

Proposed Indication and Dosing

Proposed Indication

- LUMISIGHT is indicated for fluorescence imaging in adults with breast cancer as an adjunct for the intraoperative detection of cancerous tissue within the resection cavity following removal of the primary specimen during lumpectomy surgery (also known as breast-conserving surgery)

Proposed Dosing

- For use only as a single dose of 1 mg/kg, 2-6 hours prior to imaging

Benefits of LUM System Outweigh Risks

Benefits

- Enables removal of residual cancer missed by SoC surgery and pathology
- Converts positive margins to negative, sparing potential second surgeries
- Does not appear to worsen cosmesis
- Provides real-time *in vivo* imaging to surgeons
- As adjunct to SoC, improves surgical outcomes

Risks

- False positives can lead to unnecessary tissue removal
- 0.6% serious hypersensitivity or anaphylaxis risk (4 / 726 patients)

Benefits of removing residual cancer left behind by SoC surgery outweigh any identified risks that can be managed in preoperative setting and with labeling

Sponsor-Proposed Risk Mitigation Strategies for LUMISIGHT Administration

Clear Labeling

Training Program

Enhanced Pharmacovigilance

Postmarket Study

Agenda

Unmet Need	Kelly Hunt, MD, FACS, FSSO Professor and Chair, Department of Breast Surgical Oncology, Division of Surgery MD Anderson Cancer Center President, Society of Surgical Oncologists
Efficacy	Shelley Hwang, MD, MPH Mary and Deryl Hart Distinguished Professor of Surgery Vice Chair of Research Department of Surgery Leader, Breast Cancer Disease Group Duke University and Duke Cancer Institute
Safety	Peter Blumencranz, MD, FACS Medical Director, BayCare Oncology Service Line Medical Director, The Comprehensive Breast Care Center of Tampa Bay
Allergic Reactions / Hypersensitivity	Tanya Laidlaw, MD, FAAAAI Director of Translational Research, Division of Allergy and Clinical Immunology Chief, Section of Clinical and Translational Sciences, Div. of Allergy and Clinical Immunology Brigham and Women's Hospital, Associate Professor, Harvard Medical School
Risk Mitigation Strategies	Jorge Ferrer, PhD Chief Scientific Officer Lumicell
Clinical Perspective	Barbara Smith, MD, PhD Director, Breast Program, Massachusetts General Hospital Massachusetts General Hospital Trustees Chair in Breast Surgery Professor of Surgery, Harvard Medical School

Additional Experts

Gheorghe Doros, PhD, MBA

Professor of Biostatistics
Boston University, School of Public Health
Director of Statistical Consulting
Baim Institute for Clinical Research

Michael Whitworth, MD

American Board of Anesthesiologists
Managing Partner, Prn Anesthesia

Simona Shaitelman, MD

Professor of Breast Radiation Oncology
Director, Division of Radiation Oncology Biomarker
Strategic Initiative Laboratory
UT MD Anderson Cancer Center
Vice-Chair, ASTRO Partial Breast Irradiation Clinical
Practice Guideline

Dorothy Wong, MD

Chair of Pathology, Reg Med Center of San Jose, CA
Medical Director/Staff Pathologist, Dignity Health



Unmet Need


Kelly Hunt, MD, FACS, FSSO

Professor and Chair, Department of Breast Surgical
Oncology, Division of Surgery

MD Anderson Cancer Center

President, Society of Surgical Oncology

Breast Cancer Is Most Common Cancer in Women



women in United States will develop breast cancer in their lifetime¹

Breast Cancer Is Life-Threatening Disease

> 300,000

women estimated to be diagnosed with breast cancer in US in 2023¹

~ 43,000

patients will die from breast cancer each year in US¹


2nd

leading cause of cancer death in women in US¹


> 180,000

patients undergo lumpectomy each year in US²

Caring for Patients with Breast Cancer Is Complex



Standard of Care Can Fail to Achieve Complete Resection



Potential Negative Consequences from Incomplete Resections and Second Surgeries



Anxiety



Morbidity



Cosmesis




Delay in adjuvant therapy



Surgical Complications

Limitations of SoC Pathology Margin Assessment



Pathology Assessment of Excised Tissue Can Lead to False Positives and False Negatives




~ 19%

of negative margins have residual tumor remaining¹



~ 65%

of positive margins have no actual tumor left behind¹



~ 9 – 36%

re-excision rates^{2,3}

Summary of Unmet Need

1

Current tools limited and do not identify extent of tumor accurately enough, making it challenging to achieve complete tumor resection

2

Limitations lead to second surgeries

3

Clear unmet need for visualization tool that looks inside breast cavity for residual cancer during surgery to enable more complete resection

Adjunctive to SoC, LUM System enables *in vivo* cavity assessment in real-time for more effective resection




Pivotal Study CL0007 Efficacy Results


E. Shelley Hwang, MD, MPH

Mary and Deryl Hart Distinguished Professor of Surgery
Vice Chair of Research Department of Surgery
Leader, Breast Cancer Disease Group
Duke University and Duke Cancer Institute

Pivotal Study Included 14 Medical Settings – Academic and Community Hospitals




Study CL0007 Randomized, Blinded Clinical Study



Study Procedures for Patients Randomized to Treatment Arm

Initial
Lumpectomy





Second
Surgery

LUM-Image Results Compared with Pathology Findings

True Positive



Positive LUM-signal



LUM-guided shave **contains tumor**

False Positive

Positive LUM-signal






LUM-guided shave does **not contain tumor**

False Negative

Second surgery finds tumor

No additional tissue excised

Negative LUM-signal






True Negative

Second surgery finds no tumor

No additional tissue excised

Negative LUM-signal



Pivotal Study Endpoints

Primary Endpoints

- **Removal of Residual Cancer:** % pts. with residual cancer in LUM-shave
- **Tissue-level Sensitivity:** true positive rate of LUM-positive signal
- **Tissue-level Specificity:** true negative rate of LUM-negative signal

Clinically Relevant Secondary Endpoints

- Patients converted from positive margins after SoC to final negative margins
- Average volume of LUM-guided shaves; contribution to total excision volume

Exploratory Endpoint

- Impact of LUM-guided shave volume on patient perceived cosmesis

Performance Goals for Co-Primary Endpoints

- **Removal of Residual Cancer: lower bound of CI > 3%**

- Based on published estimates of 5.3% local recurrence after whole breast radiation¹

- **Sensitivity: lower bound of CI > 40%**

- Based on previous study, SoC margin pathology showed 38% sensitivity

- **Specificity: lower bound of CI > 60%**

- Based on previous study showing ~ 1 LUM-guided shave/patient of 68%

Key Enrollment Criteria

Inclusion Criteria

- Female
- Age \geq 18 years
- Histologically or cytologically confirmed primary invasive breast cancer, DCIS, or primary invasive breast cancer with DCIS component
- ECOG 0 or 1

Exclusion Criteria

- Bilateral breast cancer and undergoing bilateral resection procedure
- Received neoadjuvant therapies
- Administration of blue dyes for sentinel lymph node mapping prior to LUM imaging
- History of allergic reaction to polyethylene glycol or any oral or IV contrast agent


Baseline Demographic Characteristics

Characteristic, %	LUMISIGHT / Lumicell DVS N = 357
Age, mean (SD)	62.4 (9.6)
Race	
White	83%
Black	6%
Asian	6%
Other, unknown or not reported	5%
Hispanic or Latino	3%
BMI, mean (SD)	29.8 (6.7)
Menopausal status	
Postmenopausal	84%
Pre/perimenopausal	16%

Baseline Tumor Histology Characteristics

Characteristic, %	LUMISIGHT / Lumicell DVS N = 357
Largest dimension of tumor in main specimen (cm), mean (SD)	1.7 (1.3)
Tumor histology (preoperative)	
DCIS only	20%
IDC	70%
ILC	10%
IDC + ILC	1%
Node positive disease	15%
No lymph node resection	19%

Co-Primary Efficacy Endpoint Met: LUM-Guided Shaves with Residual Cancer Removed in 7.6% of Patients



Aggressive and Extensive Residual Cancer Found in LUM-Guided Shaves


Characteristic, n (%)	LUM System N = 27
Tumor grade	
1	2 (7%)
2	12 (44%)
3	13 (48%)
Residual cancer size 1 – 13mm	20 (74%)
Residual cancer removed after negative margin	19 (70%)

**LUM System enabled removal of aggressive, sizable,
and undetected cancerous tissue**

Tissue-Level Sensitivity Not Met

Missed Lower Bound of PG by 3.6 Percentage Points


		Hierarchy Truth Standard	
		Positive	Negative
LUM Signal	Positive	TP = 34	FP = 337
	Negative	FN = 35	TN = 1,940



Tissue-Level Specificity Co-Primary Endpoint Met

Exceeded Lower Bound of PG

		Hierarchy Truth Standard	
		Positive	Negative
LUM Signal	Positive	TP = 34	FP = 337
	Negative	FN = 35	TN = 1,940




Diagnostic Accuracy of 84%

Exceeds 50% Expected for Random Outcome Test

		Hierarchy Truth Standard	
		Positive	Negative
LUM Signal	Positive	TP = 34	FP = 337
	Negative	FN = 35	TN = 1,940

$$\text{Accuracy} = \frac{TP + TN}{TP + TN + FP + FN} = 84\% \text{ (95\% CI: 82.6, 85.6)}$$


Diagnostic Performance Demonstrates Effectiveness



Secondary Endpoint: LUM-Guided Conversion of Positive Margins to Final Negative Margins

Clinically Relevant Pre-Defined Secondary Endpoint, n (%) (95% CI)	Efficacy Population N = 357
Patients having positive margins after SoC lumpectomy procedure n (%)	62 (17%) (13.6, 21.7)
Percent of patients converted from positive margins after SoC lumpectomy procedure to final negative margins by excising LUM-guided shaves n/n (%) <ul style="list-style-type: none"> ▪ 8 patients avoided second surgery by removing LUM-guided shaves ▪ 1 patient still underwent second surgery despite final negative margins 	9/62 (15%) (6.9, 25.8)
Of remaining 53 patients with SoC pathology-determined positive margins <ul style="list-style-type: none"> ▪ 45 patients proceeded to second surgery <ul style="list-style-type: none"> – 28 patients had no residual cancer found 	28/45 (62%) (48.1, 76.4)

LUMISIGHT Activation in Areas Adjacent to Tumor



**LUMISIGHT's MoA generates elevated fluorescence adjacent to tumor;
reasonable to attribute conversion to negative margins to drug effect**


Secondary Endpoint: Contribution to Excision Volume and (Exploratory) Its Impact to Patient Perceived Cosmesis

Clinically Relevant Secondary Endpoint, % (\pm SD)

N = 166

Contribution of LUM-guided shaves to total excision volume when ≥ 1 removed

20% (\pm 15%)



35 Patients Had Improved Surgical Outcomes by Removing LUM-Guided Shaves



Surgical benefits from LUM-guided shaves, n (%)

**Efficacy Population
N = 357**

Patients with improved surgical outcomes

35* (9.8%)

Residual cancer removed

27

Converted to final negative margins

9

Summary of Efficacy

1

Removal of residual cancer co-primary endpoint met performance goal; LUM System enabled residual cancer removal in 27 (8%) patients

2

Tissue-level sensitivity endpoint missed; tissue-level specificity endpoint met; 84% diagnostic accuracy

3

Converted 15% of positive margins to negative, sparing 8 patients second surgeries

4

Results do not suggest that additional tissue resection driven by LUMISIGHT worsened cosmetic outcomes

5

Provided real-time, *in vivo* examination of lumpectomy cavity



Safety

Peter Blumencranz, MD, FACS

Medical Director, BayCare Oncology Service Line
Health System

Medical Director, The Comprehensive Breast Care
Center of Tampa Bay

LUMISIGHT Safety Profile Well Characterized from 726 Patients

Population	Number of Patients
Overall safety population	726
Breast cancer safety population	703
Other solid tumors	23

Pivotal Study CL0007 includes > 50% of safety population (N = 406)

Administration of LUMISIGHT in Preoperative Area Under Medical Supervision

Administered 2-6 hours prior to imaging at 1 mg/kg dose by IV injection over 3 minutes

Performed in preoperative area under medical supervision

All serious events were managed immediately with standard interventions

Premedication at discretion of physician



Adverse Events Were Infrequent and Mostly Unrelated Other Than Chromaturia as Expected

Preferred Term, n (%)	Overall Safety Population N = 726
AEs	633 (87%)
AEs related to LUMISIGHT	615 (85%)
Chromaturia (discolored urine)	613 (84%)
Hypersensitivity (includes 4 SAEs in next slide)	9 (1%)
Extravasation	4 (0.6%)
Blood creatinine decreased	4 (0.6%)
AEs not related to LUMISIGHT	151 (21%)
Seroma	31 (4%)
Breast Pain	22 (3%)
Nausea	15 (2%)

Few Patients Experienced Serious Adverse Event

Preferred Term, n (%)	Overall Safety Population N = 726
SAEs	7 (1%)
SAEs related to LUMISIGHT	4 (0.6%)
Anaphylactic reaction	3 (0.4%)
Hypersensitivity	1 (0.1%)
SAEs not related to LUMISIGHT	3 (0.4%)
Breast cellulitis	1 (0.1%)
Vascular pseudoaneurysm	1 (0.1%)
Somnolence	1 (0.1%)
Acute kidney injury	1 (0.1%)
Acute respiratory failure	1 (0.1%)

Related AEs Leading to Study Discontinuation Were Infrequent and All Events Resolved

Preferred Term, n (%)	Overall Safety Population N = 726
AEs leading to discontinuation	8 (1%)
AEs related to LUMISIGHT leading to discontinuation	8 (1%)
Hypersensitivity reaction	3 (0.4%)
Anaphylactic reaction	2 (0.3%)
Extravasation event	2 (0.3%)
Nausea	1 (0.1%)
Skin discoloration	1 (0.1%)

All events resolved, most on same day

No Deaths

Category, n	Overall Safety Population N = 726
Deaths	0

Summary of Safety

1

LUMISIGHT safety profile at 1 mg/kg characterized in 726 patients; dose well tolerated

2

All patients with AEs and SAEs recovered and continued to receive SoC lumpectomy procedure

3

Personally enrolled and used LUM System in > 65 patients; comfortable using LUMISIGHT



Allergic Reactions and Hypersensitivity

Tanya Laidlaw, MD, FAAAAI

Director of Translational Research, Division of Allergy and Clinical Immunology

Chief, Section of Clinical and Translational Sciences, Division of Allergy and Clinical Immunology, Brigham and Women's Hospital
Associate Professor, Harvard Medical School

Expert Allergists Involved in Review of Allergic Reaction Events



Tanya Laidlaw, MD, FAACAI

Director of Translational Research, Division of Allergy and Clinical Immunology
Chief, Section of Clinical and Translational Sciences, Division of Allergy and Clinical Immunology,
Brigham and Women's Hospital
Associate Professor, Harvard Medical School



Jamie Waldron, MD

Allergist and Immunologist, Massachusetts General Hospital
Instructor of Medicine, Harvard Medical School



Anna Wolfson, MD, FAACAI

Chair of Quality and Safety, Allergy and Immunology and Assistant Clinical Director, Allergy and
Immunology, Massachusetts General Hospital
Assistant Professor, Harvard Medical School

Trial Reported SAEs of Anaphylaxis or Hypersensitivity

Patient	Reported in Trial
Patient #1	Anaphylaxis Life-threatening
Patient #2	Hypersensitivity Severe
Patient #3	Anaphylaxis Severe
Patient #4	Anaphylaxis Severe

- Goal of post-hoc analysis was to further characterize allergic reactions and suggest appropriate mitigations
- Anaphylaxis guidelines: CTCAE (used in trial), EAACI, NIAID, WAO, USDAR, Ring and Messmer, Brown and NAP6

Patient #1 Summary

- Cefazolin IV given 6 minutes prior to LUMISIGHT
- 1.5-2 minutes into LUMISIGHT administration (received 30mg of 104mg): patient reported chest tightness, dyspnea, upper body pain, noted to have a red face
- Administration of LUMISIGHT stopped
- Anesthesiologist reported patient as nauseous, diaphoretic, dyspneic, appearing cyanotic + apneic, having a weak pulse with generalized rash
- Treatment: 10L oxygen, epinephrine, Pepcid, Solumedrol IV, Benadryl IV; transferred to MICU
- Symptoms resolved in < 12 hours; discharged following day; lumpectomy performed 17 days later

Allergists Post-Hoc Review

- Probably related, life-threatening, anaphylaxis
- Etiology could have been cefazolin or LUMISIGHT; LUMISIGHT more likely given timing
- Patient had history of hives to iodinated contrast media

Patient #2 Summary

- Nuclear medicine injection and image-guided wire insertion 75 minutes prior to LUMISIGHT
- Tylenol 1000mg and gabapentin 300 mg given 32 minutes prior to LUMISIGHT
- 2 minutes into LUMISIGHT administration (received 27 mg of 61 mg): patient reported nausea, vomiting, headache, and lightheadedness; found to have profuse erythema, heart rate in 50s and BP 60/30 mmHg
- Infusion stopped
- Treatment: reclined; 500mL IV normal saline, Zofran 4 mg IV, Benadryl 25 mg IV
- Symptoms resolved within 13 minutes; lumpectomy occurred next day
- Allergy-related lab work: histamine (52 → 22), tryptase (11.5 → 12.6)

Allergists Post-Hoc Review

- Probably related, severe, anaphylaxis
- Tylenol and gabapentin near LUMISIGHT administration → uncommon causes of allergic reactions
- Elevated histamine and less-so tryptase suggest mast cell-mediated mechanism of reaction

Patient #3 Summary

- 1.5 minutes into LUMISIGHT administration (received 22 mg of 91 mg): patient reported dyspnea, tingling in tongue / hands / feet, nausea, swollen lip, eye redness, seeing “black spots”
- Vital signs normal with heart rate 88 and BP 110/89 mmHg
- Treatment: Benadryl 50 mg IV, hydrocortisone 100 mg, Zofran 4 mg, Pepcid 20 mg
- Patient recovered within 20 – 30 minutes; lumpectomy occurred same day
- Allergy-related lab work: histamine (55 → 11), tryptase (3.6 → 4.3)

Allergists Post-Hoc Review

- Probably related, moderate, possible allergic reaction
- Lab results reassuring after quick improvement
- Patient symptoms mostly subjective; documentation shows absence of tachypnea / hypoxia
- Briefly elevated histamine suggests mast cell-mediated mechanism of reaction

Patient #4 Summary

- During 3-minute LUMISIGHT administration (full dose at 61 mg completed), patient reported feeling “funny” with itching in hands, feet and lips; BP 125/76 mmHg immediately after injection, 10 minutes later had BP 98/51 mmHg, 5 minutes after that BP 64/38 mmHg
- Treatment: 1L lactated ringers; reverse Trendelenburg
- Blood pressure normalized within 30 minutes; symptoms resolved within 70 minutes
 - Felt well during subsequent needle localization procedure
 - Brought to PACU in wheelchair felt lightheaded: experienced vasovagal event
 - Treated with 10 mg IV ephedrine, symptoms resolved; Lumpectomy occurred same day
- Allergy-related lab work: histamine (< 8 → < 8), tryptase (4.2 → 4.6)

Allergists Post-Hoc Review

- Possibly related, moderate, vasovagal reaction
- Heart rate remained stable with subsequent vasovagal reaction ~ 3 hours after LUMISIGHT
- Diagnosis of hypersensitivity reaction unlikely due to symptoms resolving with IV fluids alone, and completely normal blood histamine and tryptase levels

Summary of SAEs of Anaphylaxis or Hypersensitivity

Patient	Reported in Trial	Allergist Review
Patient #1	Anaphylaxis Life-threatening	Anaphylaxis Life-threatening
Patient #2	Hypersensitivity Severe	Anaphylaxis Severe
Patient #3	Anaphylaxis Severe	Possible allergic reaction Moderate
Patient #4	Anaphylaxis Severe	Vasovagal reaction Moderate

All 4 patients had reactions that were well identified, well managed, and did not prevent continuing with SoC lumpectomy

Preoperative Settings Well Equipped to Manage Anaphylaxis Risk

- Rate of preoperative mortality due to anaphylaxis expected to be very low
 - Patient is verbal, monitored for reaction, and skin is visible
 - No deaths in any Lumicell clinical trials due to anaphylaxis or any other AEs
- Preoperative and operating rooms already well equipped and well trained to manage anaphylaxis due to commonly used perioperative agents
 - 0.5% of new exposures to cefazolin report an allergic reaction¹
(50% of Lumicell population administered cefazolin prior to surgery)
 - ~ 2% of exposures to isosulfan blue report an allergic reaction²

Summary of Allergic Reactions and Hypersensitivity

1

Events infrequent with rate of 0.6% (4 / 726 patients); study protocol updated after first anaphylactic reaction

2

All events occurred at healthcare setting, treated by trained personnel, fully recovered, and proceeded to SoC lumpectomy

3

Risk of mortality expected to be extremely low in preoperative setting

4

Observed rates of anaphylaxis and hypersensitivity infrequent and acceptable in context of perioperative procedures

5

Mitigations reasonable to manage rate of reactions



Risk Mitigation Strategies


Jorge Ferrer, PhD

Chief Scientific Officer
Lumicell

Risk Mitigation Strategies – Labeling

- Proposed additional warnings and details in Prescribing Information
 - Clearly indicate risk of "life-threatening anaphylaxis" in Highlights and Warnings and Precautions section
 - Advise healthcare providers that before LUMISIGHT administration, obtain history of allergy and prior hypersensitivity reactions
 - Indicate that patients with history of multiple food or drug allergies or other hypersensitivities may be at increased risk
 - Specify to always administer LUMISIGHT in healthcare settings and have emergency resuscitation drugs, equipment, and trained personnel available
 - Instruct to interrupt injection if hypersensitivity reaction is suspected
 - Monitor patients for 15 minutes after injection

Risk Mitigation Strategies – Training



Risk Mitigation Strategies – Enhanced Pharmacovigilance (PV)

- Partnered with PV vendor with experience in combination products
- Provide clear, accurate, and timely medical information
- Collect, evaluate, and report adverse events
- Implement an Adverse Events of Special Interest (AESI) program
- Train users on Lumicell's PV program
- Standardize collection of additional data to help us learn more about etiology of reactions

Risk Mitigation Strategies – Postmarket Study

Design

- Prospective, observational study

Objectives

- Primary objective: evaluate incidence of anaphylactic reactions after administration of LUMISIGHT
- Secondary objective: evaluate incidence of other hypersensitivity symptoms after administration of LUMISIGHT

Data collection

- Baseline and post-injection: vital signs; tryptase and histamine
- Complete medical histories regarding allergies
- Details on patient status and concomitant medications preceding LUMISIGHT injection
- Adverse events or symptoms related to hypersensitivity
- Treatment and outcome

Sponsor-Proposed Risk Mitigation Strategies for LUMISIGHT Administration

Clear Labeling

Training Program

Enhanced Pharmacovigilance

Postmarket Study




Clinical Perspective

Barbara Smith, MD, PhD

Director, Breast Program, Massachusetts General Hospital
Massachusetts General Hospital Trustees Chair in Breast
Surgery

Professor of Surgery, Harvard Medical School

We Believe Change Is Needed Now



Urgent Need for Improved Tool for More Complete Resection

Intraoperative Tools

- Limited to *ex vivo* specimen analysis
- Predict specimen margin status
- Do not directly assess cavity

Pathology Margins

- < 1% of surface¹
- Examines deformed specimen
- 1 to 2 weeks
- Positive margins require 2nd surgery

Second Surgeries

- Healing has deformed cavity during 2nd surgeries
- 65% of time no tumor in positive margin patients²

Benefits of LUM System Outweigh Risks

Benefits

- Enables removal of residual cancer missed by SoC surgery and pathology
- Converts positive margins to negative, sparing potential second surgeries
- Does not appear to worsen cosmesis
- Provides real-time *in vivo* imaging to surgeons
- As adjunct to SoC, improves surgical outcomes

Risks

- False positives can lead to unnecessary tissue removal
- 0.6% serious hypersensitivity or anaphylaxis risk (4 / 726 patients)

Benefits of removing residual cancer left behind by SoC surgery outweigh any identified risks that can be managed in preoperative setting and with labeling



Moderator for Q&A

Jorge Ferrer, PhD

Chief Scientific Officer
Lumicell

LUMISIGHT™ and Lumicell™ Direct Visualization System (DVS) as Adjunct to Standard of Care to Identify Residual Cancer Within the Lumpectomy Cavity

March 5, 2024

















Lumicell

Medical Imaging Drugs Advisory Committee (MIDAC)

Performance by Tumor Histology


Tumor Histology Based on Pre-op Biopsy	Patients	Residual Cancer Removal Rate	Sensitivity	Specificity	Margin Conversion
		% (95% CI)	% (95% CI)	% (95% CI)	% (95% CI)
All	357				
IDC/ILC ± DCIS	287				
DCIS only	70				

Performance by Tumor Grade

Tumor Grade	Patients n	Residual Cancer % (95% CI)	Sensitivity % (95% CI)	Specificity % (95% CI)	Margin Conversion % (95% CI)
All	357				
Grade 1	59				
Grade 2	182				
Grade 3	107				

Performance by SoC Procedure:



Comprehensive, Selective, and No Shaves



LUM-Image Results Compared with Pathology Findings

True Positive



Positive
LUM-signal



n = 34 shaves

False Positive


Positive
LUM-signal



n = 337 shaves

False Negative

Negative
LUM-signal




Second surgery
finds tumor



n = 24 shaves


No additional
tissue excised



n = 11 shaves

True Negative

Negative
LUM-signal




Second surgery
finds **no** tumor



n = 43 shaves

No additional
tissue excised



n = 1,897 shaves