



IceCure ProSense™ Cryoablation System

**General and Plastic Surgery Devices
Panel Meeting**

November 7, 2024

Introduction and Company Overview

Shay Levav

Vice President of Quality Assurance,
Regulatory Affairs and Clinical Application
IceCure Medical

Introducing IceCure Medical




**Bringing advanced
minimally invasive
cryoablation solutions**

**IceCure's flagship
product, ProSense™, is
a cryosurgical tool for
treatment of tumors
in women's health and
interventional oncology
fields**

**ProSense™ utilizes
Liquid Nitrogen (LN2)
for optimal tumor
destruction**


**Freezing tumors
quickly with minimal
pain, rapid recovery,
and great cosmetic
results**

ProSense™ Cryoablation System Approvals



*Under Private License Importation

ProSense™ Cryoablation System Approvals



*Under Private License Importation

ProSense™ Cryoablation System

- FDA cleared since 2007 for tumor ablation (including breast fibroadenomas)
- ICE3 clinical study, initiated in 2014, is the largest clinical trial to evaluate cryoablation without excision for treatment of breast cancer
- FDA granted ProSense™ Breakthrough Device Designation in March 2021



Expands patient choice with minimally invasive alternative to standard-of-care lumpectomy

Proposed Indication

The ProSense™ is indicated for use in the treatment of patients with early stage, low-risk breast cancer* for the treatment of breast cancer with adjuvant endocrine therapy

*Patients ≥ 60 years of age with prognostic stage 1A defined as unifocal tumor size ≤ 1.5 cm, ER+/PR+/-, HER2-, histological grade 1-2 infiltrating ductal carcinoma (excluding lobular carcinoma, extensive intraductal component, or evidence of lymphovascular invasion), and clinically negative lymph node (N0)

ICE3 Study Conclusions

ICE3 Study Shows Benefits of Cryoablation Outweigh Risks

IBTR-free patients

**>95% through 5 years
follow-up**

**>96% through 5 years
when treated with
adjuvant endocrine therapy**

- Near immediate recovery to normal activity – median 1 day
- >99% of patients were satisfied with cosmetic results at 5-years
- Risks are sufficiently mitigated through training and real-time visualization of tumor ablation during treatment and routine annual mammography

ICE3 demonstrates ProSense™ is safe and effective for the treatment of early stage, low-risk breast cancer

Agenda

Patient Selection & Treatment Options

Nathalie Johnson, MD, DHL, FACS

Sr. Medical Director
Legacy Cancer Institute – Portland, Oregon

Principles of Cryoablation & Prior Studies

Robert Carlton Ward, MD

Associate Professor of Diagnostic Imaging
Warren Alpert Medical School of Brown University

ICE3 Study Design & Results

Richard E. Fine, MD, FACS

Breast Program Director & Director of Education & Research
Margaret West Comprehensive Breast Center, WCC & RI

FDA-Requested Post-Hoc Analyses

Margeaux Rogers, MS, RAC

Vice President, Regulatory Affairs
MCRA, LLC, an IQVIA Business

Clinical Perspective

Richard E. Fine, MD, FACS

Patient Selection & Treatment Options

Nathalie Johnson, MD, DHL, FACS

Senior Medical Director, Legacy Cancer Institute

Portland, Oregon

Breast Cancer Background

>300,000
NEW CASES
will be diagnosed in
women in 2024¹

**Heterogeneous complex of diseases,
a spectrum of many subtypes with
distinct biological features**

**Biological features determine response
patterns to various treatment
modalities and clinical outcomes**

Early Stage, Low-risk Breast Cancer Tumor Characteristics

Lowest Risk



Highest Risk

Tumor Size (T)	T1: Tumor size ≤ 2 cm	T2: Tumor size 2-5 cm	T3: Tumor size > 5 cm	T4: Tumor extends to skin or chest wall
Lymph Nodes (N)	N0: No lymph node metastasis	N1: Metastasis to ipsilateral, movable axillary LNs	N2: Metastasis to ipsilateral fixed axillary, or IM LNs	N3: Metastasis to infraclavicular/supraclavicular LN, or to axially and IM LNs
Metastasis (M)	M0: No distant metastasis	M1: Distant metastasis		

Early stage, low-risk breast cancer patients are most appropriate for de-escalation of care

Early Stage, Low-risk Breast Cancer Tumor Characteristics

Lowest Risk ←

→ Highest Risk

Tumor Size (T)	T1: Tumor size ≤2 cm	T2: Tumor size 2-5 cm	T3: Tumor size >5 cm	T4: Tumor extends to skin or chest wall
Lymph Nodes (N)	N0: No lymph node metastasis	N1: Metastasis to ipsilateral, movable axillary LNs	N2: Metastasis to ipsilateral fixed axillary, or IM LNs	N3: Metastasis to infraclavicular/supraclavicular LN, or to axially and IM LNs
Metastasis (M)	M0: No distant metastasis	M1: Distant metastasis		
Subtype	Luminal A			
Nottingham Grade	1 or 2			
Receptor	ER+, PR+			
Her2 Status	Her2neu-			

Early stage, low-risk breast cancer patients are most appropriate for de-escalation of care

Specific Patient Selection Based on Tumor Risk

Risk Factor	ICE3 Selected / Indicated Population
Prognostic Stage	1A
Size	≤2 cm (further reduced to ≤1.5 cm)
Lymph Node Involvement	No lymph node metastases
Metastases	No distant metastases
Age	≥60 years
Receptor Status	ER+/PR+/-
HER2/neu	HER2-
Histological Grade (Nottingham) Subcomponents: nuclear, mitotic, tubule	Composite Grade 1/2

De-escalation Trends in Breast Cancer “Removal” Surgery



Radical Mastectomy
Halstead 1890's

De-escalation Trends in Breast Cancer “Removal” Surgery



Radical Mastectomy
Halstead 1890's



Modified Radical Mastectomy
Dyson & Patey 1948

**Total Mastectomy
+ Radiation**
McWhirter 1948

De-escalation Trends in Breast Cancer “Removal” Surgery



Radical Mastectomy
Halstead 1890's



Modified Radical Mastectomy
Dyson & Patey 1948

**Total Mastectomy
+ Radiation**
McWhirter 1948



**Lumpectomy,
Axillary Dissection,
Radiation**
late 1980's/1990's

De-escalation Trends in Breast Cancer “Removal” Surgery



Radical Mastectomy
Halstead 1890's



Modified Radical Mastectomy
Dyson & Patey 1948

**Total Mastectomy
+ Radiation**
McWhirter 1948



**Lumpectomy,
Axillary Dissection,
Radiation**
late 1980's/1990's



**Lumpectomy/ Sentinel
Lymph Node Biopsy**

De-escalation Trends in Breast Cancer “Removal” Surgery



Radical Mastectomy
Halstead 1890's



Modified Radical Mastectomy
Dyson & Patey 1948

**Total Mastectomy
+ Radiation**
McWhirter 1948



**Lumpectomy,
Axillary Dissection,
Radiation**
late 1980's/1990's



**Lumpectomy/ Sentinel
Lymph Node Biopsy**

CAN WE DO LESS?

Cryoablation of Breast Tumors



Cryoablation
1985 to present

De-escalation of Care – Omission of Radiotherapy

CALGB 9343:

- RT vs. No RT: 3% difference at 5 years and 8% at 10 years, in IBTR
- 1% difference in overall survival

PRIME II:

- RT vs. No RT: 3% difference at 5 years and 9% at 10 years, in IBTR
- No differences in distant metastases, contralateral breast cancer, overall survival, new breast cancers

NCCN CLINICAL PRACTICE GUIDELINES IN ONCOLOGY

Breast Cancer, Version 3.2022

William J. Gradishar, MD^{1*}; Meena S. Moran, MD^{2,†}; Jame Abraham, MD^{3*}; Rebecca Aft, MD, PhD^{4*}; Doreen Agnese, MD⁵; Kimberly H. Allison, MD⁶; Bethany Anderson, MD^{7*}; Harold J. Burstein, MD, PhD⁸; Helen Chew, MD^{9*}; Chau Dang, MD¹⁰; Anthony D. Elias, MD¹¹; Sharon H. Giordano, MD, MPH¹²; Matthew P. Goetz, MD¹³; Lori J. Goldstein, MD¹⁴; Sara A. Hurvitz, MD¹⁵; Steven J. Isakoff, MD, PhD¹⁶; Rachel C. Jankowitz, MD¹⁷; Sara H. Javid, MD¹⁸; Jairam Krishnamurthy, MD¹⁹; Marilyn Leitch, MD²⁰; Janice Lyons, MD^{3*}; Joanne Mortimer, MD²¹; Sameer A. Patel, MD¹⁴; Lori J. Pierce, MD²²; Laura H. Rosenberger, MD, MS²³; Hope S. Rugo, MD²⁴; Amy Sitapati, MD^{25*}; Karen Lisa Smith, MD, MPH^{26*}; Mary Lou Smith, JD, MBA²⁷; Hatem Soliman, MD²⁸; Erica M. Stringer-Reasor, MD²⁹; Melinda L. Telli, MD⁶; John H. Ward, MD³⁰; Kari B. Wisinski, MD⁷; Jessica S. Young, MD^{31*}; Jennifer Burns, BS^{32*}; and Rashmi Kumar, PhD^{32*}

Despite differences in local recurrence rate, no meaningful difference in overall survival

2022/2024 National Comprehensive Cancer Network (NCCN) guidelines updated in 2024:

Consider omitting RT following BCS in women aged 65 years and older with Stage I, ER+ breast cancer who receive adjuvant endocrine therapy

Despite differences in local recurrence rate, the clinical community is recommending de-escalation of care

De-escalation of Care – Less Axillary Surgery, Less Radiation




5 QUESTIONS to Ask Your Doctor Before You Get Any Test, Treatment, or Procedure

- 1** Do I really need this test or procedure?
- 2** What are the risks and side effects?
- 3** Are there simpler, safer options?
- 4** What happens if I don't do anything?
- 5** How much does it cost, and will my insurance pay for it?


CAN WE DO LESS?

Are there even simpler, safer options?

Patient Experience: Lumpectomy



Patient Experience: Cryoablation




Principles of Cryoablation & Prior Studies

Robert C. Ward, MD


Associate Professor of Diagnostic Imaging

Warren Alpert Medical School of Brown University


Cryoablation, a Basic Ultrasound-guided Procedure



Real-time Ultrasound Needle Placement




Pre





Needle Placement

Real-time Ultrasound Procedural Monitoring




Real-time Ultrasound Margin Assessment



Follow-up Imaging: Ultrasound




Pre




6 Months

Follow-up Imaging: Ultrasound





Pre




6 Months

Follow-up Imaging: Mammogram




Pre




6 Months

Follow-up Imaging: MRI




Pre




6 Months

Decades of Research and Prior Studies on Cryoablation in Treatment of Breast Cancer




Decades of Research and Prior Studies on Cryoablation in Treatment of Breast Cancer



Roubidoux 2004

Study of **9 patients**
treated with US-guided
cryo and exision



No residual invasive cancer
in tumors 1.7 cm or smaller or
in cancers without spiculated
margins at US

Littrup 2009

Study of **11 patients** treated with US or US+CT-guided cryoablation



100% procedural success with mean tumor size 1.7 cm

No significant complications, retraction, or scarring were noted

Safely achieved 1 cm visible ice beyond tumor margins with **minimal discomfort, good cosmesis, no short-term tumor recurrences**

Manenti 2013


Study of **80 patients** comparing **RFA vs. cryoablation** in treatment of early breast cancer



Both resulted in **good clinical and cosmetic outcome**

Cryotherapy is the **preferred method** because of analgesic effect of freezing with better patient experience

Decades of Research and Prior Studies on Cryoablation in Treatment of Breast Cancer



2016 ACOSOG Study

Ann Surg Oncol (2016) 23:2438–2445
DOI 10.1245/s10434-016-5275-3

Annals of
SURGICAL ONCOLOGY
OFFICIAL JOURNAL OF THE SOCIETY OF SURGICAL ONCOLOGY

A Phase II Trial Exploring the Success of Cryoablation Therapy in the Treatment of Invasive Breast Carcinoma: Results from ACOSOG (Alliance) Z1072


Rache M. Simmons, MD¹, Karla V. Ballman, PhD², Charles Cox, MD³, Ned Carp, MD⁴, Jennifer Sabol, MD⁴, Rosa F. Hwang, MD⁵, Deanna Attai, MD⁶, Michael Sabel, MD⁷, David Nathanson, MD⁸, Andrew Kenler, MD⁹, Linsey Gold, MD¹⁰, Cary Kaufman, MD¹¹, Linda Han, MD¹², Aaron Bleznak, MD¹³, J. Stanley Smith, MD¹⁴, Dennis Holmes, MD¹⁵, Bruno Fornage, MD¹⁶, Carisa Le-Petross, MD¹⁶, Syed Hoda, MD¹⁷, Linda McCall, MS¹⁸, Kelly K. Hunt, MD⁵, and on behalf of the ACOSOG investigators

KEY CONCLUSIONS

- If 1 cm below, **100% success**
- Excluding multifocal disease **92% success** in tumors ≤ 2 cm
- No cases of ineffective ablation

**Prospective *ablate and resect*, proved effectiveness of cryoablation
Led directly to the ICE3 trial without need for resection**

Decades of Research and Prior Studies on Cryoablation in Treatment of Breast Cancer



Cost of Cryoablation vs. Lumpectomy

Ann Surg Oncol

<https://doi.org/10.1245/s10434-022-12570-5>

Annals of
SURGICAL ONCOLOGY
OFFICIAL JOURNAL OF THE SOCIETY OF SURGICAL ONCOLOGY

The Role of Cryoablation in Breast Cancer Beyond the Oncologic Control: COST and Breast-Q Patient-Reported Outcomes


Sonia Y. Khan, BS^{1,2}, Annie Snitman, MD¹, Zaina Habrawi, MD^{1,2}, Sybil Crawford, PhD³, Michael W. Melkus, PhD^{1,2}, and Rakhshanda Layeequr Rahman, MD^{1,2}

KEY CONCLUSIONS

- Better psychosocial well-being
- Cryoablation significantly associated with better physical, sexual, and cosmetic satisfaction outcomes
- Lower financial toxicity

Better short-term psychosocial well-being, lower financial toxicity with cryoablation vs. lumpectomy procedures

Decades of Research and Prior Studies on Cryoablation in Treatment of Breast Cancer



ICE3 Study Design & Results

Richard E. Fine, MD, FACS

Breast Program Director & Director of Education & Research
Margaret West Comprehensive Breast Center, WCC & RI
Germantown, TN

Overall Study Design

IRB Approved, multi-centered (19 U.S. sites), single-arm, prospective study

Study Population

Early stage, low-risk breast cancer

Women ≥ 60 with unifocal histologic grade 1 or 2, invasive ductal carcinoma, ≤ 1.5 cm, ER+, PR +/-, Her2-, clinically lymph node neg

Intervention


Cryoablation

US-guided precise tumor engulfment with an ice ball

Adjuvant treatment per physician discretion

Follow-up

Mammography (MRI optional*)



Months

Outcome

Ipsilateral Breast Tumor Recurrence (IBTR)

5 years post-cryo

All patients enrolled in the ICE3 study had early stage, low-risk breast cancer

Study Design Considerations

Single-arm design

- **Unable to conduct a blinded study**
 - Cryoablation vs. surgical excision
- **Inherent difficulties in conducting a randomized controlled study**
 - Performed by both surgeons and radiologists

Literature performance goal

- **Primary outcome, IBTR, is objective**
- **Surgical lumpectomy is well-understood and outcomes well-documented in the published literature**
- **Published outcomes on thousands of patients more robust than a prospective control**

Sample size

- **ICE3 builds on a large body of prior work**
 - Data on ablate and resect demonstrate effectiveness in tumor destruction
- **Statistically justified sample size of 200 patients**

Key Inclusion/Exclusion Criteria

Key Inclusion Criteria

- 1. Diagnosis of invasive ductal breast carcinoma by core needle biopsy, meeting the following criteria:**
 - a. Unifocal primary disease**
 - b. Tumor size ≤ 1.5 cm in greatest diameter**
 - c. Histologic (Nottingham) grade 1-2**
 - d. Estrogen receptor positive, and/or progesterone receptor positive, HER2 negative**
- 2. Age ≥ 60 (WCG IRB), age ≥ 50 (Local IRB)**

Key Exclusion Criteria

- 1. Lobular carcinoma**
- 2. Luminal B pathology**
- 3. Microinvasion or invasive breast carcinoma with extensive intraductal component (EIC)**
- 4. Multifocal and/or multicentric in breast cancer, multifocal calcifications**
- 5. Neoadjuvant chemotherapy for breast cancer**
- 6. Prior en bloc open surgical biopsy and/or lumpectomy for the index breast cancer**

Primary Endpoint

- **Protocol Definition of Local Recurrence (IBTR):**
 - Evidence of invasive or in situ breast cancer in the ipsilateral breast or chest wall, including new primary confirmed by biopsy
 - Suspicious findings do not constitute criteria for breast cancer recurrence. Any recurrence of malignant disease should be proven by biopsy or excision
- **Independent DSMB Consideration of Local Recurrence (IBTR):**
 - Per clinical practice, local recurrence does not include new ipsilateral tumor in a different quadrant or at least 5 cm from original tumor. Instead, considered second primary breast cancer
 1. Smith, T., et al. *Int. J. Radiat. Oncol. Biol. Phys.*, (2000).
 2. Yi, M, et al. *Ann. Surg.*, (2011).
 3. Huang, E., et al. *Cancer*, (2002).
 4. Panet-Raymond, V, et al. *Int. J. Radiat. Oncol. Biol. Phys.*, (2011).
 5. Belkacemi Y, et al. *Front. Oncol.*, (2018).

Protocol specified all recurrence events be confirmed by biopsy

Prespecified Literature-based Endpoint

Lumpectomy outcomes are well-established

Systematic literature review using PubMed and EBSCO yielded >1000 hits and resulted in selection of 4 papers reporting IBTR following breast conserving surgery at 2-5 years post-op

Performance goal determined to be <10%

- Literature supports 5-year reference rate of 5%
- Reference margin of 5% was clinically justified
- Performance goal determined by adding the reference margin to the reference rate

If the upper limit of the 95% CI at the 5-year time point is <10%, study will be considered successful

Key Secondary Endpoints: Effectiveness

Regional Occurrence

Tumor in ipsilateral internal mammary, ipsilateral supraclavicular, ipsilateral infraclavicular, and/or ipsilateral axillary nodes or soft tissue of ipsilateral axilla

Distant Metastases

Tumor in any area of the body except those defined as local or regional

Disease Free Survival (Protocol Definition)

Freedom from disease events including local (DCIS or invasive), regional or distant breast cancer recurrence, second primary breast cancer, DCIS or invasive contralateral breast cancer, **second primary non-breast cancer and death due to any cause**

Disease Free Survival (NCI Definition)

Freedom from disease events including local (DCIS or invasive), regional or distant breast cancer recurrence, second primary breast cancer, DCIS or invasive contralateral breast cancer

Overall Survival

Freedom from death due to any cause

Breast Cancer Survival

Freedom from death due to breast cancer or unknown cause

Protocol specified all recurrence events be confirmed by biopsy

Patient Reported Outcomes & Safety



Patient and physician satisfaction
with breast cosmetic outcome




Safety: Adverse events
related to study device or procedure



Time to recovery
and resumption of normal activities

ICE3 Study Results

Analysis Population




- 3 screen fail
- 3 withdrawal before procedure



- 12 excluded by DSMB due to inclusion/exclusion deviations (9) or incomplete treatment (3)

Demographics

Primary Analysis Set (N=194)



Age:

74.9

(Range 55-94)


BMI:

28.8

(Range 15.0-47.6)

Gender:

100% female



Disease Characteristics

Primary Analysis Set

Baseline Characteristics	ICE3 N=194 n (%)
ER (Estrogen Receptor)	
Positive	194 (100)
Negative	0
PR (Progesterone Receptor)	
Positive	180 (93)
Negative	14 (7)
HER2/neu	
Positive	0
Negative	194 (100)
Histologic Grade	
1	96 (49)
2	98 (51)
3	0

Disease Characteristics

Primary Analysis Set

Baseline Characteristics	ICE3 N=194 n (%)
ER (Estrogen Receptor)	
Positive	194 (100)
Negative	0
PR (Progesterone Receptor)	
Positive	180 (93)
Negative	14 (7)
HER2/neu	
Positive	0
Negative	194 (100)
Histologic Grade	
1	96 (49)
2	98 (51)
3	0

Disease Characteristics

Primary Analysis Set

Baseline Characteristics	ICE3 N=194 n (%)
ER (Estrogen Receptor)	
Positive	194 (100)
Negative	0
PR (Progesterone Receptor)	
Positive	180 (93)
Negative	14 (7)
HER2/neu	
Positive	0
Negative	194 (100)
Histologic Grade	
1	96 (49)
2	98 (51)
3	0

Disease Characteristics

Primary Analysis Set

Baseline Characteristics	ICE3 N=194 n (%)
ER (Estrogen Receptor)	
Positive	194 (100)
Negative	0
PR (Progesterone Receptor)	
Positive	180 (93)
Negative	14 (7)
HER2/neu	
Positive	0
Negative	194 (100)
Histologic Grade	
1	96 (49)
2	98 (51)
3	0

Cryoablation Procedure: Intra/periprocedural Data

Mean tumor size:
0.73 cm (0.1–1.49)

Procedure time:
30-40 mins

Anesthesia type:
100% local

**Median time to resume
normal activities:**
1 day (0–8 days)


Procedure-related adverse events:

- Bruising
- Pain
- Localized Edema
- Injection Site Reaction
- Hematoma
- Frost Injury

Procedure-related adverse events were generally mild, occurred acutely, and resolved promptly without intervention

Adjuvant Treatment

Type of adjuvant treatment



*3 patients with unknown adjuvant therapy counted as none

Effectiveness

Primary Endpoint Met

Local IBTR Rate – Primary Analysis Set (N=194)

Time	At Risk*	Cumulative IBTR	IBTR Estimate†	2-sided 95% CI
Operative	194	-	-	-
Month 6	192	0	0.0%	-
Year 1	190	0	0.0%	-
Year 2	184	0	0.0%	-
Year 3	173	1	0.6%	0.1 - 3.9
Year 4	162	3	1.7%	0.6 - 5.3
Year 5	133	7	4.3%	2.1 - 8.7


ICE3 population met pre-specified <10% performance goal

*At risk: number of patients that completed the follow-up time interval with no IBTR event

†Kaplan-Meier (product-limit) estimate with 1-sided 95% CI upper bound (UB) and 2-sided 95% lower and upper bounds (LB and UB); IBTR=Ipsilateral Breast Tumor Recurrence

ICE3 Primary Endpoint


Primary Analysis Population (N=194)



Error bars=95% confidence intervals

ICE3 Primary Endpoint

Primary Analysis Population (N=194)



Error bars=95% confidence intervals


Secondary Kaplan-Meier Analyses

Primary Analysis Population (N=194)

Secondary Effectiveness Endpoint	Events	Year 5 KM Rate	2-sided 95% CI
Regional Recurrence Estimate	0	100.0%	-
Distant Metastases Estimate*	6	96.4%	92.2 - 98.4
Disease Free Survival - NCI Definition*	12	92.8%	87.6 - 95.8
Overall Survival Estimate	20	88.6%	82.9 - 92.5
Breast Cancer Survival Estimate	5	96.7%	92.2 - 98.6

99% of patients free from confirmed breast cancer-related death

Deaths Observed in ICE3 vs. Actuarial Survival



Fewer deaths observed in ICE3 than expected for 74-year-old population

Safety

Safety Overview

Primary Analysis Population (N=194)

	IceCure N=194		
	Events	Patients (n)	% (n/N)
Adverse Events			
All	517	140	72.2%
Procedure Related Adverse Events			
All	180	93	47.9%
Procedure Related Adverse Events by Severity			
Mild	158	83	42.8%
Moderate	18	15	7.7%
Severe	4	4	2.1%
Serious Adverse Events			
All	127	59	30.4%
Non-procedure Related	123	56	28.8%
Procedure Related	4	3	1.5%

Serious Adverse Events

Primary Analysis Population (N=194)

	IceCure N=194		
	Events	Patients (n)	% (n/N)
Serious Adverse Events			
All	127	59	30.4%
Non-procedure Related	123	56	28.8%
Procedure Related	4	3	1.5%

- **3 patients with SAEs possibly related to study procedure or device**
- **All determined to be due to physician error**
 - Probe mispositioning (1)
 - Suboptimal treatment (2)

Procedure-related Events Occurring in >2% of Patients

Primary Analysis Population (N=194)

	IceCure N=194		
	Events	Patients (n)	% (n/N)
Bruising	57	57	29.4
Pain	38	36	18.6
Localized Edema	37	35	18.0
Injection Site Reaction	11	10	5.2
Hematoma	9	9	4.6
Frost Injury	4	4	2.1

**ICE3 study demonstrated no new device-related risks,
all procedure-related events are known risks of cryoablation**

Patient Reported Outcomes & Probable Benefits

Near immediate recovery

- **83%** of the ICE3 study population returned to **full daily activities 48 hours after procedure**
- **Days to resume activity on average:** median 1 day (0 to 8 days)

High percentage of patients and physicians satisfied with cosmetic results at 5 years

- **>99%** of patients 'satisfied' or 'very satisfied' (110/111)
- **97%** of physicians 'satisfied' or 'very satisfied' (99/102)

ICE3 Summary

Robust study with ~80% patient accountability at 5-years

- ICE3 clinical study met the pre-specified primary endpoint
- >95% of patients are IBTR-free through 5-years follow-up
- Near immediate recovery to normal activity (median 1 day recovery time)
- >99% of patients were satisfied with cosmetic results at 5-years

ProSense™ demonstrated to have low risk safety profile

- Non-serious procedure-related risks are common to all cryoablation systems, including ProSense™ System when used per the cleared indications
- All serious, procedure-related events are common to breast cancer treatment and can be mitigated by treatment according to proposed labeling

The totality of evidence demonstrates safety, effectiveness, and positive benefit/risk profile of the ProSense™ for treatment of early stage, low-risk breast cancer

FDA-Requested Post-Hoc Analyses

Margeaux Rogers, MS, RAC

Vice President, Regulatory Affairs

MCRA, LLC, an IQVIA business


FDA-requested PRISMA Systematic Literature Review and Meta-analysis

- **FDA requested PRISMA systematic literature review and meta-analysis**
- **Designed to evaluate 5-year IBTR following lumpectomy without radiation**
- **Literature selection criteria were intentionally aligned with ICE3**
 - Anticipated differences in literature populations
 - Excluded unlike literature populations, excluded populations <100 patients
 - Pre-defined downweighting of data in cases where the literature were similar but not perfectly aligned as described in executive summary
- **Literature review protocol was submitted to FDA for input**
- **Search returned ~800 citations**
- **11 unique articles identified for data extraction**

FDA-requested PRISMA Systematic Review and Meta-analysis

- Overall PRISMA meta-analysis 5-year IBTR rate: 3.52%
- Range: 0% - 15.4%

- PRISMA meta-analysis 5-year IBTR rate with adjuvant endocrine therapy: 2.82%
- Range: 0% - 10.5%



FDA-requested PRISMA Systematic Review and Meta-analysis

- Overall PRISMA meta-analysis 5-year IBTR rate: 3.52%


- Range: 0% - 15.4%

ICE3 5-year recurrence rate
(N=194): 4.3%

- PRISMA meta-analysis 5-year IBTR rate with adjuvant endocrine therapy: 2.82%

- Range: 0% - 10.5%

ICE3 Indicated Population
(N=147): 3.1%



FDA-requested Comparison to LUMINA Study

LUMINA Study: Explored omission of radiotherapy following treatment of breast-conserving surgery and endocrine therapy in patients with luminal A breast cancer

LUMINA Study Population (N=500)	ICE3 Study Population (N=194)
Luminal A subtype	2 patients determined to have luminal B type cancer, 2 unknown
Median age: 67.1 years	Median age: 74.9 years
Selected patients with good surgical outcomes	Prospective study design
All patients received adjuvant endocrine therapy	Only a subset of patients in ICE3 received adjuvant endocrine therapy

ICE3 expected to be worse case due to differences in adjuvant treatment

ICE3 Subpopulation for Comparison to LUMINA




**ICE3 Primary Analysis Population
(N=194)**


**Subpopulation with Adjuvant Endocrine Therapy
without Radiation
(N=124)**

ICE3 Primary Endpoint

Adjuvant Endocrine Therapy Only (N=124) vs. LUMINA (N=500)




ICE3 Subpopulation for Comparison to LUMINA




ICE3 Primary Endpoint

Adjuvant Endocrine Therapy + Biological Characteristics (N=56) vs. LUMINA (N=500)




ICE3 Subpopulation for Comparison to LUMINA



ICE3 Primary Endpoint

+ Endocrine Therapy + Biological Chars. No Nuclear 3/Unknown (N=48) vs. LUMINA (N=500)



Proposed Indication – Updated per FDA Recommendation

The ProSense™ is indicated for use in the treatment of patients with early stage, low-risk breast cancer* for the treatment of breast cancer **with adjuvant endocrine therapy**

*Patients ≥60 years of age with prognostic stage 1A defined as unifocal tumor size ≤1.5 cm, ER+/PR+/-, HER2-, histological grade 1-2 infiltrating ductal carcinoma (excluding lobular carcinoma, extensive intraductal component, or evidence of lymphovascular invasion), and clinically negative lymph node (N0)

FDA-requested Subpopulation Analysis

Aligned with Proposed Indications




**ICE3 Primary Analysis Population
(N=194)**

**Subpopulation of Patients Aligned with
Proposed Indication
(N=147)**

ICE3 Primary Endpoint

Indicated Subpopulation (N=147) vs. PRISMA with Adjuvant Endocrine Therapy



FDA-requested Sub-population Analysis

Aligned with Proposed Indications + Nuclear Grade Restriction




**ICE3 Primary Analysis Population
(N=194)**

**Subpopulation of Patients Aligned with
Proposed Indication
(N=147)**

**Subpopulation without Nottingham Component Scores
OR with Nuclear Grade of 3
(N=120)**

ICE3 Primary Endpoint

Indicated Subpopulation without Nuclear 3/Unknown (N=120) vs. PRISMA




FDA Performed PRISMA Systematic Literature Review and Meta-Analysis

- **25 articles identified with IBTR rates from 0% to 12%**
- **5 articles selected for meta-analysis (IBTR: 0% to 2.3%)**
 - 2 studies included <25 subjects, both with IBTR rates of 0%
- **FDA SLR resulted in overall IBTR rate of 0.61% (CI UB: 3.5%)**
 - IBTR in FDA SLR sub-population with radiation: 0.68% (CI UB: 2.94%)
 - IBTR in FDA SLR population without radiation: 0.47% (CI UB: 35.91%)
- **IceCure concerns:**
 - Unclear if article selection process was applied uniformly
 - IceCure identified articles were excluded, despite relevance
 - Outcome weighting method to derive extremely low IBTR rate

SLR Key Differences – Summary

	IceCure	FDA
Adjuvant Radiotherapy	Excluded	Adjuvant radiotherapy was employed in 21 out of the 25 studies
Sample Size	Excluded samples sizes <100 patients	No exclusion criteria
Adherence to Selection Criteria	Included all included studies in the meta-analysis	Included only 5 of the 25 included studies in the meta-analysis
Weighting Criteria	Weighted based on population alignment – downweighting studies with variability	Used inverse variance weighting – resulted in upweighting of studies with no recurrence reported
Risk of Bias	All studies received “low” risk of bias judgment	Two of the five studies included in the meta-analysis received “serious” risk of bias judgement

FDA Meta-analysis of 5-Year IBTR Rates from Selected Studies



1. Article noted to have moderate or serious risk of bias due to missing data

2. HER2 status was not reported in two studies selected for the quantitative analysis

3. Includes regional recurrences in the reported locoregional recurrence rate

*Indicates corrected values where imputation and zero correction were applied for missing or zero event rates. Corrected values are based on SAS output where missing data for number at risk or events were imputed or corrected as described in methods

FDA Systematic Literature Review (SLR): Weight Compared to Sample Size (ICE3 Study)

Study	Events (Rate)	Number at Risk	Percentage of Total Subjects		Weight
Ciervide, 2018	0 (0.0%)	16	2%	1%	42.1%
Soyder, 2013	0 (0.0%)	11		1%	29.3%
Offersen, 2022 (1)	3 (0.7%)	396	29%		3.8%
Offersen, 2022 (2)	5 (1.2%)	379	28%		6.3%
Kunkler, 2015 (2)	5 (1.3%)	324	24%		6.3%
Whelen, 2023	10 (2.3%)	246	18%		12.3%
	23 (0.61%)	1372	100%		100%

FDA's SLR heavily weights the two smallest populations – introduces bias in the resulting meta-analysis rate

Subpopulation Treated with Adjuvant Radiotherapy



**ICE3 Primary Analysis Population
(N=194)**

**Subpopulation treated with Adjuvant Radiation
(N=29)**

- **ICE3 subpopulation with radiation 0% IBTR rate**
- **FDA SLR rate with radiation is 0.68% (95% CI: 0.15-2.94, Range 0-1.3%)**

ICE3 subpopulation with radiation has better outcomes than FDA SLR

Comparison of IBTR Results by Censoring Approach

Primary Analysis Population (N=194)


- FDA recommended censoring method (CIF) and IceCure censoring method agree that the primary endpoint was met

Number of IBTR events	KM Rate (95% CI) Censoring at IBTR event or study discontinuation (IceCure method)	KM Rate (95% CI) Death censored at time of death, or final follow-up for LTFU	KM Rate (95% CI) Death and LTFU censored at final follow-up (FDA method)	CIF Rate (95% CI) (FDA recommended method)
7 including >Month 60	4.3% (2.1 - 8.7)	4.6% (2.2 - 9.5)	5.0% (2.4 - 10.2)*	4.6% (2.0 - 8.9)
6 through Month 60	3.6% (1.6 - 7.9)	4.0% (1.8 - 8.6)	4.1% (1.9 - 8.9)	4.0% (1.6 - 8.0)

IceCure IBTR rate is conservative given inclusion of >M60 event

*FDA cites 5.2% (2.5-10.7); however, IceCure statisticians are unable to reproduce FDA results
LTFU=lost-to-follow-up or withdrawal from study

Summary of ICE3 Analysis Populations and Respective Literature Comparators



ICE3 IBTR rate similar to lumpectomy in all subpopulation analyses

ICE3 Confirms ProSense™ is Successful in the Treatment of Breast Cancer

Ann Surg Oncol
https://doi.org/10.1245/s10434-024-16181-0

Annals of
SURGICAL ONCOLOGY
OFFICIAL JOURNAL OF THE SOCIETY OF SURGICAL ONCOLOGY



ORIGINAL ARTICLE – BREAST ONCOLOGY

Cryoablation Without Excision for Early-Stage Breast Cancer: ICE3 Trial 5-Year Follow-Up on Ipsilateral Breast Tumor Recurrence

Richard E. Fine, MD¹, Richard C. Gilmore, MD¹, Kenneth R. Tomkovich, MD², Jill R. Dietz, MD^{3,4}, Michael P. Berry, MD¹, Lydia E. Hernandez, MD⁵, Karen S. Columbus, MD⁵, Susan A. Seedman, MD⁶, Carla S. Fisher, MD⁷, Linda K. Han, MD⁷, Eric R. Manahan, MD⁸, Randy D. Hicks, MD⁹, Rashmi P. Vaidya, MD¹⁰, Lisa D. Curcio, MD¹¹, Alexander B. Sevrakov, MD¹², Andrew S. Kenler, MD¹³, Bret Taback, MD¹⁴, Margaret Chen, MD¹⁴, Megan E. Miller, MD³, Linsey Gold, MD¹⁵, Beth V. Anglin, MD¹⁶, Hussein D. Aoun, MD¹⁷, Rache M. Simmons, MD¹⁸, Sheldon M. Feldman, MD¹⁹, and Susan K. Boolbol, MD^{20,21}

¹Margaret West Comprehensive Breast Center, West Cancer Center and Research Institute, Germantown, TN; ²CentraState Medical Center, Freehold, NJ; ³University Hospital Cleveland Medical Center, Cleveland, OH; ⁴Cleveland, OH; ⁵Cincinnati Breast Surgeons Inc, Cincinnati, OH; ⁶Presbyterian Cancer Center, Rio Rancho, NM; ⁷Indiana University IU Health Hospital-University Hospital, Indianapolis, IN; ⁸Dalton Surgical Group, Dalton, GA; ⁹Regional Medical Imaging, Flint, MI; ¹⁰Ironwood Cancer and Research Centers, Glendale, AZ; ¹¹BreastLink Laguna Hills, Laguna Hills, CA; ¹²Thomas Jefferson University Hospital, Philadelphia, PA; ¹³Bridgeport Hospital, Bridgeport, CT; ¹⁴Columbia University Medical Center, New York, NY; ¹⁵Comprehensive Breast Care, Troy, MI; ¹⁶Complete Breast Care, Plano, TX; ¹⁷Karmanos Cancer Institute, Detroit, MI; ¹⁸Weill Cornell Weill Medical College, Cornell University, New York, NY; ¹⁹Montefiore-Einstein Center for Cancer Care, New York, NY; ²⁰Mount Sinai Beth Israel, New York, NY; ²¹Nuvance Health System, Poughkeepsie, NY



Clinical Perspective

Richard E. Fine, MD, FACS

Breast Program Director & Director of Education & Research
Margaret West Comprehensive Breast Center, WCC & RI
Germantown, TN

Patient Experience



So, how was the experience? *It was fine. Great.*

No pain or anything? *No pain.*

Okay, good. Good.

Physician/Patient Discussion

Discuss benefits and risks of all options

- Ask the **Choosing Wisely** questions
- **Provide range of options** including de-escalated alternatives

Consider tumor characteristics and patient history to identify best option for the specific patient

Counsel patient that with any treatment for breast cancer there is a risk of:

- Incomplete treatment
- Recurrent tumors

Procedural risks common to all cryoablation systems

- Edema, bruising, hematoma, hypothermic frost injury, and postoperative pain

De-escalated care involves shared decision making

ProSense™: A Needed Minimally Invasive Option

Patients
deserve
to choose

Less morbid
treatment

More financially
accessible treatment

Better
cosmetic outcome

Greater
quality of life

- Indicated population with early stage, low-risk breast cancer are ideal candidates for de-escalation
- The ICE3 clinical study provides sufficient evidence to conclude that the device benefits outweigh risks for the treatment of early stage, low-risk breast cancer
- IceCure is eager to work with FDA to design appropriate special controls, labeling and training to ensure good patient selection and safe & effective use

Patients should have the choice of all available options



IceCure ProSense™ Cryoablation System

**General and Plastic Surgery Devices
Panel Meeting**

November 7, 2024

Sponsor Backup Slides Shown



What Data is Available as Real-world Data?

- **The ProSense™ System or an equivalent cryoablation device was used in over 1,600 procedures for the treatment of breast cancer in independent clinical studies**
- **Across the clinical literature, 18 cohort studies conducted in United States, Japan, Spain, Germany, Romania, Italy, and the Netherlands with follow-up to 16 years**
- **10 studies reported 22 recurrences in up to 1,366 patients (1.61%)**
- **No serious device or treatment related adverse events were reported**
- **Post market surveillance (PMS) data was collected in 150 BSC patients**
 - In total, 4/140 (2.9%) of patients reported adverse events (1 hematoma, 2 minor skin burns, and 1 not specified)
 - Physician satisfaction was graded as “Excellent” in 124/141 (88%) cases
 - Physician satisfaction was graded as “Good” in 15/141 (11%) cases and “Medium” in the remaining 2/141 (1%) cases

What Happens to the Destroyed Cells after Cryoablation?

30 minutes: Immediate post-thaw phase

Cells swell & irreversible damage

8 hours to 1 week: Hemorrhagic and inflammatory phase

8-12 hours

Apoptosis progressively increases at the peripheral zone

2 days

Coagulation necrosis, characterized by hemorrhage, edema, and inflammation

7 week

Infiltration of inflammatory cells, fibrin and collagen stranding, and capillary ingrowth sharply demarcate the periphery of the lesion

1 to 3 months: Replacement phase

1 month

Necrotic tissue is largely cleared with new blood vessels at the periphery

3 months

A small fibrotic scar with an intact endothelial layer

In weeks to months, the necrotic tissue is slowly removed by phagocytotic activity of the immune system and replaced by a fibrous collagenous scar

Kameda Medical Center Clinical Trial

465 Patients Over 16 Years - Long Term FU

Prof. Eisuke Fukuma, pioneer in cryotherapy for breast cancer

Group	Year	Tumor Size	# of Patients
IceCure	2014-2022	≤15 mm	560 (465 with long term FU)

Local recurrence (IBTR) - 0.65% (3/465)

Patient selection › Cryoablation › Adjuvant Rx › Follow-up



Age: mean 57 (range 31-83)

ER/PR+, HER2-
Sentinel Node Imaging ± Bx Negative
Whole Breast Radiotherapy
Endocrine/Chemo Rx per Biology



All post-cryo annual follow up

PE, Mammo, US, MRI

IBTR Sensitivity Analysis – Ki67

ICE3 Study: Primary Population N=194

Ki67 Stratification	N	Events	5-Year IBTR	95% LB	95% UB
<14%	93	3	3.93%	1.28%	11.76%
≥14%	36	1	3.33%	0.48%	21.39%
Unknown	65	3	5.23%	1.72%	15.36%

No difference in 5-year IBTR rate by Ki67

IBTR Sensitivity Analysis – Age

ICE3 Study (Primary Analysis Population N=194)

Primary Population (N=194)			
Age	Subjects n (%)	Recurrence Rate	95% CI
55 to 60 years	4 (2.1)	0%	–
61 to 70 years	47 (24.2)	5.3%	1.4 – 19.8
71 to 80 years	100 (51.5)	4.5%	1.7 – 11.5
81 to 90 years	41 (21.1)	3.0%	0.4 – 19.6
91 to 94 years	2 (1.0)	0%	–

In ICE3 elevated risk not seen in patients under 70

LUMINA Trial Patient Accounting

- **740 registered patients**
 - 224 patients Ki67 \geq 13.25%
 - 11 patients had insufficient specimens
 - 4 patients were identified by central monitoring as ineligible
 - 1 patient withdraw
- **500 patients enrolled**
- **246 at risk at 5-years (49.2%)**

DSMB Exclusion – Patient 12

- Reason for single freeze: machine malfunctioned. Was only able to do 1st freeze, thaw and then only 2 minutes of the 2nd freeze

Procedure Date	Recurrence Date
18-Jan-2019	No recurrence

DSMB Exclusion – Incomplete Treatment (n=3)

Pt #	Procedure Date	DSMB Meeting Date	Detailed Reason	Outcomes
1	03 Oct 2016	01-03 May 2019	Extremely short treatment that could not cover the lesion. On procedure day, tumor dimensions were 0.7*0.6*0.7 cm, final ice ball dimensions were 1.31 *4.52 cm, treatment cycle times were 1:48F1-8:00T-1:22F2 (mm:sec)	No Recurrence
2	14 Sep 2016	01-03 May 2019	Extremely short treatment that could not cover the lesion. On procedure day, tumor dimensions were 0.86*0.1*0.84 cm, final ice ball dimensions were 0.4 * 2.47 cm, treatment cycle times were 1:22F1-1:51T-1:59F2 (mm:sec)	No Recurrence
3	18 Jan 2019	03-05 Sep 2019	Insufficient treatment due to machine malfunction. On procedure day, tumor dimensions were 0.6*0.7*0.3 cm, final ice ball dimensions were 3.33 *5.3 cm, treatment cycle times were 9:15F1-8:01T- 2:24F2 (mm:sec)	No Recurrence