

Section 13, Product Insert

1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16

Lyra™ Influenza A Subtype H7N9 Assay Instructions for Use

For the qualitative detection and identification of avian influenza A (H7N9) virus (detected in China in 2013) viral RNA extracted from nasal swab and nasopharyngeal swab specimens.

**For Use under an
Emergency Use Authorization Only**

17

18	Contents	
19	For Emergency Use Only	1
20	Intended Use.....	4
21	Summary and Explanation.....	4
22	Principle of the Procedure.....	5
23	Materials Provided.....	6
24	Materials Required But Not Provided.....	6
25	Warnings and Precautions.....	6
26	Storage and Handling of Kit Reagents.....	7
27	Specimen Collection, Storage and Handling.....	7
28	Nucleic Acid Extracts Storage.....	8
29	bioMérieux NucliSENS easyMAG Nucleic Acid Extraction Programming Instructions.....	8
30	Applied Biosystems 7500 Fast Dx Programming Instructions.....	10
31	Assay Procedure.....	14
32	Applied Biosystems® 7500 Fast Dx Thermocycler Test Procedure.....	14
33	Quality Control.....	15
34	Limitations.....	17
35	Clinical Performance.....	17
36	Analytical Performance.....	18
37	Limit of Detection (LoD).....	18
38	Analytical Reactivity – <i>In silico</i> Prediction.....	20
39	Analytical Specificity – Non-A/H7 Influenza A Viruses.....	20
40	Analytical Specificity Study Results – Non-A/H7 Influenza A Viruses.....	20
41	Analytical Specificity/Cross-Reactivity – Common Respiratory Bacteria, Yeast and Viruses.....	22

42	Cross-Reactivity Study Results –Viruses	22
43	Cross-Reactivity Study Results –Bacteria and Yeast	24
44	Microbial Interference – Common Respiratory Bacteria, Yeast and Viruses	25
45	Microbial Interference Study Results –Viruses.....	26
46	Microbial Interference Study –Bacteria and Yeast	28
47	Customer and Technical Assistance.....	30
48	Trademarks	30
49		
50		

51 Intended Use

52 The Lyra™ Influenza A Subtype H7N9 Assay is intended for the *in vitro* qualitative detection and
53 identification of viral RNA from the Avian Influenza A (H7N9) Virus (detected in China in 2013) in nasal and
54 nasopharyngeal swabs from patients with signs and symptoms of respiratory infection. The assay is
55 performed on the Applied Biosystems® 7500 Fast Dx instrument. It is indicated for the presumptive
56 identification of virus in patients who may be infected with the Avian Influenza A (H7N9) Virus (detected
57 in China in 2013) from nasal and nasopharyngeal swab specimens in conjunction with clinical and
58 epidemiological risk factors.

59 The Lyra™ Influenza A Subtype H7N9 Assay targets the hemagglutinin (HA) gene of the Avian Influenza A
60 (H7N9) Virus (detected in China in 2013), which may react to other influenza A/H7 viruses of the Eurasian
61 Lineage.

62 The presence of influenza A viral RNA in a nasal or nasopharyngeal swab specimen must first be
63 established using an FDA-cleared influenza A device and must also be determined as “un-subtypable” by
64 FDA-cleared influenza device(s) with subtyping capabilities for all currently circulating influenza A viruses
65 in the United States (i.e., seasonal A/H3 and A/H1 pandemic) prior to testing with the Lyra™ Influenza A
66 Subtype H7N9 Assay.

67 Testing with the Lyra™ Influenza A Subtype H7N9 Assay should not be performed unless the patient
68 meets clinical and epidemiologic criteria for testing suspect specimens.

69 Negative results do not preclude influenza virus infection and should not be used as the sole basis for
70 treatment or other patient management decisions.

71 If infection with a novel influenza A virus is suspected based on current clinical and epidemiological
72 screening criteria recommended by public health authorities, specimens should be collected with
73 appropriate infection control precautions for novel virulent influenza viruses and sent to state or local
74 health department for testing. Viral culture should not be attempted in these cases unless a BSL 3+ facility
75 is available to receive and culture specimens.

76 The Lyra™ Influenza A Subtype H7N9 Assay is for use under the Food and Drug Administration’s
77 Emergency Use Authorization only.

78

79 Summary and Explanation

80 Human infections with a novel avian influenza A (H7N9) virus as well as poultry infections continue to be
81 reported in China. While some mild illnesses in human H7N9 cases have been observed, most patients
82 have had severe respiratory illness, with about one-third resulting in death.. No cases of this novel virus
83 outside of Asia have been reported, and it has not been detected in people or birds in the United States.

84 Chinese authorities continue to investigate novel avian influenza A (H7N9) cases. Although many patients
85 infected with this virus are reported to have had contact with or exposure to poultry, some cases
86 reportedly have not had such contact or exposure. Close contacts of confirmed H7N9 patients are being
87 followed to determine whether any human-to-human transmission of H7N9 is occurring. No sustained
88 person-to-person transmission of this virus has been found at this time.

89

90 It’s likely that sporadic cases of H7N9 associated with poultry exposure will continue to occur in China.
91 Cases associated with poultry exposure also may be detected in neighboring countries, and it is possible
92 that this virus may be detected in the United States at some point, likely in a traveler returning from an
93 affected area. Most concerning is the pandemic potential of this virus. Influenza viruses constantly

94 change and it is possible that this novel virus could gain the ability to easily and sustainably spread
 95 between people, triggering a pandemic.

96

97 **Principle of the Procedure**

98 The Lyra™ Influenza A Subtype H7N9 Assay detects and identifies the Avian Influenza A (H7N9) Virus
 99 (detected in China in 2013) viral RNA that has been extracted from a patient sample using the NucliSENS®
 100 easyMAG® automated extraction platform. A multiplex real-time RT-PCR reaction is carried out under
 101 optimized conditions in a single tube generating amplicons for the targeted virus (if present) and the
 102 Process Control (PRC) present in the sample. This reaction is performed utilizing the Applied Biosystems
 103 7500 Fast Dx platform. Identification of the Avian Influenza A (H7N9) Virus (detected in China in 2013)
 104 occurs by the use of target specific primers and a fluorescent-labeled probe that hybridizes to a conserved
 105 region of the hemagglutinin gene of the Avian Influenza A (H7N9) Virus (detected in China in 2013).
 106

Lyra™ Probe Labels	
Target	Dye
Influenza A/H7	FAM
Process Control (PRC)	Quasar® 670

107

108 The following is a summary of the procedure:

109

- 110 1. **Sample Collection:** Obtain nasal swabs or nasopharyngeal swabs using standard techniques from
 111 symptomatic patients. These specimens are transported, stored, and processed according to established
 112 laboratory procedures.ⁱ
 113
- 114 2. **Nucleic Acid Extraction:** Extract nucleic acids from the specimens with the NucliSENS® easyMAG® System
 115 following the manufacturer’s instructions and using the appropriate reagents (See **Materials Required but**
 116 **Not Provided**).

117

118 Prior to the extraction procedure add 20 µL of the Process Control (PRC) to each 180 µL aliquot of
 119 specimen. The PRC serves to monitor inhibitors in the extracted specimen, assures that adequate
 120 amplification has taken place and confirms that the nucleic acid extraction was sufficient.

121

- 122 3. **Rehydration of Master Mix:** Rehydrate the lyophilized Master Mix using 135µL of Rehydration Solution.
 123 The Master Mix contains oligonucleotide primers, fluorophore and quencher-labeled probes targeting
 124 conserved region of the hemagglutinin gene of the Avian Influenza A (H7N9) Virus (detected in China in
 125 2013), as well as the process control sequence. The probes are dual labeled with a reporter dye attached
 126 to the 5’ end and a quencher attached to the 3’ end. The rehydrated Master Mix is sufficient for eight
 127 reactions.

128

- 129 4. **Nucleic Acid Amplification and Detection:** Add 15 µL of the rehydrated Master Mix to each plate well. 5
 130 µL of extracted nucleic acids (specimen with PRC) is then added to the plate well. Place the plate into the
 131 Applied Biosystems® 7500 Fast Dx instrument.

132

133 Once the reaction plate is added to the instrument, the assay protocol is initiated. This protocol initiates
 134 reverse transcription of the RNA targets generating complementary DNA, and the subsequent
 135 amplification of the target sequences occurs. The Lyra™ Influenza A Subtype H7N9 Assay is based on
 136 TaqMan® chemistry, and uses an enzyme with reverse transcriptase, DNA polymerase, and 5’-3’
 137 exonuclease activities. During DNA amplification, this enzyme cleaves the probe bound to the
 138 complementary DNA sequence, separating the quencher dye from the reporter dye. This step generates
 139 an increase in fluorescent signal upon excitation by a light source of the appropriate wavelength. With

140 each cycle, additional dye molecules are separated from their quenchers resulting in additional signal. If
 141 sufficient fluorescence is achieved the sample is reported as positive for the detected target sequence.
 142

143 **Materials Provided**

144 SKU # EUA-M100-S

145 Detection Kit (96 Reactions) – Store at 2° to 8°C
 146

#	Component	Quantity
①	Rehydration Solution Part M5003	1 vial/kit 1.9 mL
②	Lyra™ Influenza A Subtype H7N9 Master Mix Part M5104 Lyophilized Contents: DNA polymerase enzyme with reverse transcriptase activity Oligonucleotide primer pairs; Oligonucleotide probes dNTPs (dATP, dCTP, dGTP, dUTP, dTTP) Stabilizers	12 vials/kit, 8 reactions/vial
CONTROL	Process Control Part M5005	1 vial/kit 2.0 mL
Positive Control	Avian Influenza A (H7N9) Synthetic DNA Part M5133	1 vial/kit 1.0 mL

147
 148 • Lyra™ Influenza A Subtype H7N9 Assay Instructions for Use

149
 150 **Materials Required But Not Provided**

- 151 • Micropipettors (range between 1 to 10 µL and 100 to 1000 µL)
- 152 • Non-aerosol pipette tips
- 153 • Applied Biosystems®7500Fast Dx software version 1.4
- 154 • Applied Biosystems®7500Fast Dx 96 well PCR plate
- 155 • Applied Biosystems®optical plate films
- 156 • Plate centrifuge for Applied Biosystems® 96 well plate
- 157 • bioMérieux NucliSENS easyMAG software version 2.0
- 158 • bioMérieux NucliSENS easyMAG Buffers 1, 2, 3
- 159 • bioMérieux NucliSENS easyMAG Lysis Buffer
- 160 • bioMérieux NucliSENS easyMAG Silica Magnetic Beads
- 161 • bioMérieux NucliSENS easyMAG disposables
- 162 • Biohit pipettor

163
 164 **Warnings and Precautions**

- 165 • The presence of influenza A viral RNA in the nasal or nasopharyngeal swab specimen must first be
 166 established using an FDA-cleared influenza A device and must also be determined as “un-subtypable” by
 167 FDA-cleared influenza detection device(s) with subtyping capabilities for all currently circulating influenza
 168 A viruses in the United States (i.e., seasonal A/H3 and A/H1 pandemic) prior to testing with the Lyra™
 169 Influenza A Subtype H7N9 Assay.

- 170 • For *In Vitro* Diagnostic Use
- 171 • The assay has been validated using bioMérieux NucliSENS easyMAG software version 2.0. Please contact
- 172 Quidel Technical Support prior to modifying or upgrading beyond this version of software.
- 173 • The assay has been validated using Applied Biosystems 7500Fast Dx software version 1.4. Please contact
- 174 Quidel Technical Support prior to modifying or upgrading beyond this version of software.
- 175 • Performance characteristics of this test have been established with the specimen types listed in the
- 176 **Intended Use Section** only. The performance of this assay with other specimen types or samples has not
- 177 been evaluated.
- 178 • Use of this product should be limited to personnel with sufficient training in PCR and RT-PCR techniques.
- 179 • Treat all specimen/samples as potentially infectious. Follow universal precautions when handling
- 180 samples, this kit and its contents.
- 181 • Proper sample collection, storage and transport are essential for correct results.
- 182 • Store assay reagents as indicated on their individual labels.
- 183 • Wear suitable protective clothing, gloves, eye and face protection when using this kit.
- 184 • For accurate results, pipette carefully using only calibrated equipment.
- 185 • Thoroughly clean and disinfect all surfaces with a 10% bleach solution followed by molecular grade water.
- 186 • Use micropipettes with an aerosol barrier or positive displacement tips for all procedures.
- 187 • Avoid microbial and cross contamination of the kit reagents. Follow Good Laboratory Procedures.
- 188 • Do not mix reagents from kits with different lot numbers.
- 189 • Do not use reagents from other manufacturers with this kit.
- 190 • Do not use product after its expiration date.
- 191 • Proper workflow planning is essential to minimize contamination risk. Always plan laboratory workflow in
- 192 a uni-directional manner, beginning with pre-amplification and moving through amplification and
- 193 detection.
- 194 • Use dedicated supplies and equipment in pre-amplification and amplification areas.
- 195 • Do not allow cross movement of personnel or equipment between areas.
- 196 • Keep amplification supplies separate from pre-amplification supplies at all times.
- 197 • Do not open sample tubes or unseal plates post amplification.
- 198 • Dispose of amplified material carefully and in accordance with local laws and regulations in order to
- 199 minimize the risk of amplicon contamination.
- 200 • Do not use supplies dedicated for reagent or sample preparation for processing target nucleic acid.
- 201 • MSDS is available upon request or can be accessed on the product website.

202

203 **Storage and Handling of Kit Reagents**

- 204 • Store the unopened kit at 2° to 8°C until the expiration date listed on the outer kit box.
- 205 • The rehydrated Master Mix may be stored at room temperature (20° to 25°C) for up to 24 hours. For
- 206 longer storage the rehydrated Master Mix should be recapped, sealed with parafilm and stored in an
- 207 upright position at ≤-20°C for up to 14 days. Protect the Master Mix from light during storage.

208

209 **Indications of Instability or Deterioration of Reagents:** Cloudiness of the Rehydration Solution, when within
210 expiration, may indicate deterioration of this reagent. Contact Quidel Technical Assistance for a replacement.

211

212 **Specimen Collection, Storage and Handling**

213 Nasal and nasopharyngeal specimens should be collected, transported, stored, and processed according to
214 CLSI M41-A. Specimens should be stored at 2° to 8°C until tested. If specimens cannot be tested within 72
215 hours of collection, they should be frozen at -70°C or colder until tested.

216

217 The following viral transport media (M4, M4-RT, M5, M6, and UTM) (1 mL and 3 mL) are compatible with the
218 Lyra™ Influenza A+B assay.

219

220 **Nucleic Acid Extracts Storage**



221 Eluates from the NucliSENS easyMAG can be stored at room temperature (20° to 25°C) for 2 hours, at 2° to 8°C
 222 for 8 hours and 1 month at –20° to –70°C.
 223



224 **bioMérieux NucliSENS easyMAG Nucleic Acid Extraction Programming Instructions**

225


226 Note: The Avian Influenza A (H7N9) Synthetic DNA Positive Control (i.e. Lyra™ Influenza A Subtype H7N9
 227 Assay Positive Control #M5133), and a negative process control (i.e., viral transport media or previously
 228 characterized influenza A and influenza B negative specimen) should be included in each extraction run.


229 1. Turn on the instrument and wait for instrument light to appear orange. Then switch on the
 230 computer/launch easyMAG software. Do not log into software until the light on the instrument has
 231 turned green.

232 2. Barcode reagents after pressing the 'Instrument'  and 'Reagent Inventory'  buttons.

233 3. To enter samples, press the 'Daily Use'  button, which will default to the 'Define Request' 
 234 screen. Select the following settings:

- 235 a. Sample ID: Enter the **sample name** using the keyboard.
- 236 b. Matrix: Select **Other** from the drop-down menu
- 237 c. Request: Select **Generic** from the drop-down menu
- 238 d. Volume (mL): Select **0.200** from the drop-down menu
- 239 e. Eluate (µL): Select **50** from the drop-down menu
- 240 f. Type: Primary
- 241 g. Priority: Normal


242 4. Upon pressing the 'Save'  button, the sample will appear in the 'Unassigned Sample' window on


243 the left side of the screen. Press the 'Enter New Extraction Request' 
 244 the process for additional samples. Alternatively multiple samples can be entered by pressing the 'Auto

245 Create New Extraction Requests'  button.

246 5. Once all samples are created, go to 'Organize Runs' by clicking on the  icon near the top of the

247 page. Create a run by pressing the 'Create Run'  button. Enter a run name or use the default.

248 6. Add samples to the run by using the 'Auto Fill Run' 
 249 'Unassigned Sample list' on the left hand side of the screen). Alternatively, individual samples can be

250 moved into and out of the run by using the left and right 'Positioning icons'  after

251 selecting the appropriate sample. The sample order within the run can be changed using the ‘Move




252 Extraction Request Up/Down’ buttons

253 7. Obtain 1 to 3 (for 8 to 24 samples, respectively) sample vessel(s), and add 20 µL of Process Control to each
254 sample well used.

255 8. Add 180 µL of each sample to the appropriate well as designated.



256 9. Go to ‘Load Run’ by pressing the  button near the top of the screen. Insert tips and sample
257 vessel(s) into the instrument

258 10. Enter the barcode(s) of the sample vessel(s)

259 11. Enter the barcode(s) of silica beads to be used


260 12. Close the instrument lid.

261 13. Assign silica beads to samples as follows:

262 a. Click the reagents symbol below number 1 in the picture below. The lot number of the silica
263 beads should appear below the Silica tab at number 2 in the picture below.

264 b. Highlight and select the samples in the run for which beads need to be assigned (in the box
265 containing number 3 in the picture below)





266 c. Click the  positioning icon (below number 4 in the picture below) to assign the silica lot
267 number to the selected samples

268 d. If the bead symbol to the right of number 5 in the picture below is selected, the silica bead lot
269 number should be displayed for each sample



270

271 14. Print work list by touching ‘Load Run’ icon followed by pressing the ‘Print Work List’ icon .

272 15. Press the ‘Dispense Lysis’  button. The on-board lysis will take approximately 12 minutes to
273 complete.

274 16. For each sample vessel, prepare magnetic particles using the Biohit pipettor and tips for up to eight
275 reactions as follows:

276 a. Using 1 tip and Program 1, aspirate 550 µL nuclease-free water and dispense into a 1.5 mL DNase
277 / RNase free microfuge tube.

- 278 b. Vortex the magnetic silica. Using 1 tip and Program 1, aspirate 550 µL of magnetic silica, dispense
- 279 into the water and mix by vortexing.
- 280 c. Using 1 tip and Program 2, aspirate 1050 µL of the magnetic silica mixture and dispense 25 µL
- 281 back into the same tube.
- 282 d. Dispense 125 µL magnetic silica mixture each into 8 wells of an ELISA strip plate. Discard tip.
- 283 e. After Lysis is complete (NB: the 'Instrument Status' at the bottom of the screen must be 'IDLE!'),
- 284 using 8 tips and Program 3, aspirate 100 µL of magnetic silica mixture in strip wells, dispense 100
- 285 µL of magnetic silica mixture in strip wells, and aspirate 100 µL of magnetic silica mixture in strip
- 286 wells.
- 287 f. Insert tips into liquid within the sample vessels. Aspirate 800 µL then dispense 900 µL of
- 288 magnetic silica mixture back into vessel. Aspirate 1000 µL of magnetic silica mixture from vessel
- 289 and dispense 1000 µL of magnetic silica back into vessel. Repeat aspiration / dispensing of 1000
- 290 µL two more times.



- 291 17. Close the instrument and press the 'Start' button to begin the run.
- 292 18. Upon completion of run, transfer purified nucleic acid to nuclease-free tubes. Eluates from the easyMAG
- 293 can be stored at room temperature (20° to 25°C) for 2 hours, at 2° to 8°C for 8 hours and 1 month at -20°
- 294 to -70°C.

Applied Biosystems 7500 Fast Dx Programming Instructions

- 295 1. Launch the 7500 Fast Dx software package.
- 296
- 297 2. The **Quick Startup document** dialog window will open. Select the **Create New Document** button to start
- 298 the **New Document Wizard**. Follow each step to initiate the Lyra™ Influenza A Subtype H7N9 Assay
- 299 protocol.
- 300 a. Define Document: Most of the following should be the default setting. If not, change accordingly.
- 301 i. Confirm or enter the following information.

Assay:	Standard Curve (Absolute Quantitation)
Container:	96-Well Clear
Template:	Blank Document
Run Mode:	Fast 7500
Operator:	<i>your operator name</i>
Comments:	SDS v1.4
Plate Name:	'Lyra Influenza A H7N9'

- 302 ii. Select the **Next** button.
- 303
- 304 b. Select Detectors: New detectors for Influenza A, and the process control (PRC) must be added.
- 305 For each target, select the **New Detector** button to open the **New Detector** pop-up window.
- 306 Alternatively, use the **Create Another** button from within the **New Detector** pop-up window for
- 307 the last two detectors.
- 308

309

- i. Enter the following information for each detector.

Name	Reporter Dye	Quencher Dye	Color
Influenza A	FAM	(none)	(Select)
PRC	Cy5	(none)	(Select)

310

- ii. Select a unique color to represent each detector.
- iii. Highlight the new detectors and add to the **Detectors in Document** column using the **Add** button.
- iv. Select **(none)** from the **Passive Reference** drop-down menu.
- v. Select the **Next** button.
- vi. Select the **Finish** button without setting any wells.

311

312

313

314

315

316

317

318

319

320


321

322

323

324

- c. The wizard will close and the software will open, starting with the **Setup** tab. This will show the sample plate that was set up during the quick start. For the initial set up, nothing needs to be changed here.
- d. Defining the Thermocycler Protocol: Select the **Instrument** tab to set up the Lyra™ Influenza A Subtype H7N9 Assay RT-PCR cycling times and temperatures. Under **Thermal Profile** there should be a default 2-stage protocol. Each stage will have 3 user-editable text boxes. The top box value represents the number of reps or cycles for that stage. The middle box value represents the temperature (°C) and the lowest box value represents the time (minutes: seconds).



325

326

327

328

329

330

331

332

333


334

- i. Make the following changes to the default **Thermal Cycler Protocol**:
 1. Stage 1
 - a. Reps: 1
 - b. Temp: 55
 - c. Time: 5:00
 2. Select the bar between Stage 1 and Stage 2. Select the **Add Hold** button to add another stage.
 3. Stage 2

- 335 a. Reps: 1
 336 b. Temp: 60
 337 c. Time: 5:00
 338 4. Select the bar between Stage 2 and Stage 3. Select the **Add Hold** button to add
 339 another stage.
 340 5. Stage 3
 341 a. Reps: 1
 342 b. Temp: 65
 343 c. Time: 5:00
 344 6. Stage 4 (2-Step Dissociation Stage)
 345 a. Reps: 10
 346 b. Step 1
 347 i. Temp: 92
 348 ii. Time: 0:05
 349 c. Step 2
 350 i. Temp: 57
 351 ii. Time: 0:40
 352 7. Select the bar to the right of Stage 4. Select the **Add Cycle** button to add
 353 another stage.
 354 8. Stage 5 (2-Step Dissociation Stage)
 355 a. Reps: 35
 356 b. Step 1
 357 i. Temp: 92
 358 ii. Time: 0:05
 359 c. Step 2
 360 i. Temp: 57
 361 ii. Time: 0:40
 362 9. If a wrong stage is added the stage can be removed by pressing the **Delete**
 363 button after highlighting the stage between the vertical lines
 364 ii. Under **Settings** enter the following:

Sample Volume (µL):	20 (default)
Run Mode:	7500 Fast (default)
Data Collection:	Stage 5, Step 2(57.0 @ 0:40)
NOTE: Do not check the check box next to 'Expert Mode'.	

- 365
 366
 367
 368
 369
 370 iii. Final protocol



371
372
373
374
375
376
377
378
379

- e. Set threshold for each analyte.
 - i. Select the **Results** tab.
 - ii. Select the **Amplification Plot** tab.
 - iii. Select Influenza A from the Detector tab in the top right corner.
 - iv. In the **Analysis Settings** block, set the **Threshold** to **1.5e5**.
 - v. Select the **Auto Baseline** radio button.
 - vi. Repeat iii-v for PRC setting the **Threshold** to **2.7e4**.

Data: Data:

Detector: Detector:

Line Color: Line Color:

Analysis Settings

Auto Ct
 Manual Ct

Threshold:

Auto Baseline
 Manual Baseline:

Start (cycle):

End (cycle):

Analyze

Help

Analysis Settings

Auto Ct
 Manual Ct

Threshold:

Auto Baseline
 Manual Baseline:

Start (cycle):

End (cycle):

Analyze

Help

380
381
382

- f. Save the new protocol as a template for future use.

- 383 i. At the top of the screen select **File** and then **Save As**.
- 384 ii. **Save In:** D:\Applied Biosystems\7500 Fast System\Templates\
385 iii. **File name:** 'Lyra Influenza A H7N9'
- 386 iv. **Save as type:** 'SDS Templates (*.sdt)'
- 387 g. Exit the software.
- 388

389 **Assay Procedure**

390 Run the following procedures at controlled room temperature of 20° to 25°C.

391

392 **Master Mix Rehydration Procedure**

- 393 1. Determine the number of specimens extracted to be tested and obtain the correct number of eight-
394 test lyophilized Master Mix vials for testing.
- 395 2. Return unused reagents to the appropriate storage conditions.
- 396 3. Open Master Mix carefully to avoid disruption of the pellet.
- 397 4. Add 135 µL of Rehydration Solution to the Master Mix.
- 398 5. Place vial at room temperature for 1 to 2 minutes to allow rehydration of pellet.
- 399 6. Gently pipette up and down 2 to 3 times avoiding the formation of bubbles prior to dispensing into
400 the first PCR tube.

401 **Note:** The rehydrated Master Mix is sufficient for 8 reactions.

402 **Note:** The rehydrated Master Mix may be stored at room temperature (20° to 25°C) for up to 24
403 hours. For longer storage the rehydrated Master Mix should be recapped, sealed with parafilm and
404 stored in an upright position at ≤-20°C for up to 14 days. Protect the Master Mix from light during
405 storage.

406

407

408 **RT-PCR Set-up Procedure:**

- 409 1. Add 15 µL of the rehydrated Master Mix to each plate well.
- 410 2. Add 5 µL of extracted nucleic acid (specimen with the process control) into the plate well. Mixing of
411 reagents is not required.
- 412 **Note:** Use a new barrier micropipettor tip with each extracted specimen.
- 413 3. Seal the plate.
- 414 4. Centrifuge the reaction tube or plate for a minimum of 15 seconds. Ensure that all liquid is at the
415 bottom of the tube or plate wells.
- 416 5. Insert tube or plate into the appropriate thermocycler.

417 **Note:** Previously characterized positive influenza A subtype H7N9 specimens serve as an external
418 processing and extraction control and should be treated as a patient specimen and be included in every
419 extraction and PCR run.

420

421 **Applied Biosystems® 7500 Fast Dx Thermocycler Test Procedure**

- 422 1. Switch on Applied Biosystems® 7500 Fast Dx.
- 423 2. Launch the Applied Biosystems® 7500 Fast Dx software v1.4 package.
- 424 3. The **Quick Startup** document dialog window will open.
- 425 4. Click on **Create a new document**.
- 426 5. Most of the following should be the default setting. If not, change accordingly.

427

428

Assay:	Standard Curve (Absolute Quantitation)
Container:	96-Well Clear
Template:	Lyra Influenza A H7N9
Run Mode:	Fast 7500
Operator:	<i>your operator name</i>
Comments:	SDS v1.4
Plate Name:	YYMMDD- Lyra Influenza A H7N9

429

430

6. Set Up Sample Plate

431

a. Under the **Setup** and **Plate** tabs the plate setup will appear.

432

b. Select all wells that will contain sample, right-click and select the **Well Inspector** from the drop-down menu. When the **Well Inspector** pop-up window opens, select the detectors for influenza A and PRC.

433

434

435

c. Use the **Well Inspector** to enter the sample names. Patient IDs can be entered in the Well Inspector window. However it is recommended that this is done prior to re-suspending the lyophilized master mix, post run or using the import function to minimize the time the PCR reactions will sit at room temperature prior to starting the run.

436

437

438

d. Save the run as **YYMMDD- Lyra Influenza A H7N9.sds**.

439

440

e. A window will open asking for the "Reason for change of entry". Enter "**Setup**" and any other comments relevant to the run.

441

442

7. Starting the PCR

443

a. Select the **Instrument** tab.

444

b. Insert the 96 well PCR plate into the machine.

445

c. Under **Instrument Control**, select the **Start** button to initiate the run.

446

8. Post PCR

447

IMPORTANT: When the run is finished press OK.

448

449

a. Analyze the data by pressing the "**Analyze**" button in the top menu and save the file.

450

b. Save the file by pressing **Save Document** in the task bar. A window will open asking for the "Reason for change of entry".

451

452

c. Enter "**Data analysis post run**" and any other comments relevant to the run

453 **Quality Control**

454

The Lyra™ Influenza A Subtype H7N9 Assay incorporates several controls to monitor assay performance.

455

456

1. The **Process Control (PRC)** consists of an inactivated and stabilized MS2 Bacteriophage that contains an RNA genome. It must be used during extraction and amplification in the assay. This control should be added to each sample aliquot prior to extraction. The PRC serves to monitor inhibitors in the

457

458

- 459 extracted specimen, assures that adequate amplification has taken place and confirms that the
 460 nucleic acid extraction was sufficient.
 461
 462 2. The **Positive Control** consists of synthetic Avian Influenza A (H7N9) Virus (detected in China in 2013)
 463 target amplicon specific double-stranded DNA, and must be treated as a patient specimen and be
 464 included in every extraction and RT-PCR run.
 465
 466 3. Viral transport media or previously characterized negative specimen may be used as an external
 467 **Negative Control**. This must be treated as a patient specimen and be included in every extraction and
 468 PCR run.
 469
 470 4. Failure of either the **Positive Control** or the Negative Control invalidates the RT-PCR run and results
 471 should not be reported. The RT-PCR run should be repeated with the extracted controls and
 472 specimens first. Re-extract and retest another aliquot of the controls and the specimens or obtain
 473 new samples and retest is necessary if the controls fail again.
 474

475 **Expected Results from Controls:**

Control Type/ Name	Used to Monitor	Influenza A/H7	Expected Ct Values	PRC	Expected Ct Values
Positive Control	Substantial reagent failure including primer and probe integrity	+	5.0 ≤ Ct ≤ 35.0	+/-	NA
Negative Control	Reagent and/or environmental contamination	-	None detected	+	5.0 ≤ Ct ≤ 35.0

476

477 **Interpretation of Results from Patient Specimens**

478 **Interpretation of the Lyra™ Influenza A Subtype H7N9 Assay Results on the Applied Biosystems®**

479 **7500 Fast Dx Thermocycler**

Assay Result	Detector: Influenza A/H7	Detector: Process Control	Interpretation of Results	Notes and Special Guidance
Negative	Ct < 5.0 or Ct > 35.0	5.0 ≤ Ct ≤ 35.0	No Avian Influenza A (H7N9) Virus (detected in China in 2013) viral RNA detected; PRC Detected.	
Influenza A/H7 Positive	5.0 ≤ Ct ≤ 35.0	NA*	Avian Influenza A (H7N9) Virus (detected in China in 2013)	Contact CDC or qualified Public Health laboratories immediately for

			viral RNA detected.	coordination of additional testing and for further guidance.
Invalid	Ct<5.0 or Ct>35.0	Ct<5.0 or Ct>35.0	No Avian Influenza A (H7N9) Virus (detected in China in 2013) viral RNA and no PRC RNA detected.	Invalid test. Retest the same purified sample. If the test is also invalid, re-extract and retest another aliquot of the same specimen or obtain a new specimen and retest.

480 *No Ct value is required for the Process Control to make a positive call.

481

482 **Limitations**

- 483 • Negative results do not preclude infection with influenza virus and should not be the sole basis of a
- 484 patient treatment decision.
- 485 • Improper collection, storage or transport of specimens may lead to false negative results.
- 486 • Inhibitors present in the sample and/or errors in following the assay procedure may lead to false negative
- 487 results.
- 488 • A trained health care professional should interpret assay results in conjunction with the patient’s medical
- 489 history, clinical signs and symptoms, and the results of other diagnostic tests.
- 490 • Analyte targets (viral sequences) may persist *in vivo*, independent of virus viability. Detection of analyte
- 491 target(s) does not imply that the corresponding virus(es) are infectious, nor are the causative agents for
- 492 clinical symptoms.
- 493 • There is a risk of false positive values resulting from cross-contamination by target organisms, their
- 494 nucleic acids or amplified product, or from non-specific signals in the assay.
- 495 • There is a risk of false negative values due to the presence of sequence variants in the viral targets of the
- 496 assay.
- 497 • The assay performance was not established in immunocompromised patients.

498

499 **Clinical Performance**

500 The Lyra™ Influenza A Subtype H7N9 Assay clinical performance characteristics were estimated

501 using clinical specimens from the 2013 respiratory season and contrived specimens. Due to the

502 lack of available clinical specimens containing the Avian Influenza A (H7N9) Virus (detected in

503 China in 2013), evaluation of the performance of the Lyra™ Influenza A Subtype H7N9 Assay was

504 carried out using an alternative approach. Avian Influenza A (H7N9) Virus (detected in China in

505 2013) positive samples were prepared by spiking the CDC provided BPL inactivated

506 A/Anhui/1/2013 (H7N9) virus at a concentration of approximately 2.5 x LoD in individual

507 nasopharyngeal specimens that were characterized as influenza A negative by the FDA-cleared

508 Lyra™ Influenza A + B Assay (k112172). The Avian Influenza A (H7N9) Virus (detected in China in

509 2013) negative clinical specimen were selected from characterized clinical respiratory specimens

510 from the 2013 influenza season. A total of 87 avian influenza A (H7N9) negative specimens were

511 selected. Twenty-six (26) of the 87 specimens were negative for influenza A as determined by

512 the FDA-cleared Lyra™ Influenza A + B Assay (k112172), 13 of the 87 specimens were negative

513 for respiratory viruses as determined by viral culture, 25 of the 87 specimens were positive for

514 influenza A as determined by viral culture and/or the FDA-cleared Lyra™ Influenza A + B Assay

515 (k112172), 12 of the 87 specimens were positive for influenza B as determined by viral culture

516 and/or the FDA-cleared Lyra™ Influenza A + B Assay (k112172), nine of the 87 specimens were

517 positive for respiratory syncytial virus (RSV) as determined by the FDA-cleared Quidel Molecular
 518 RSV + hMPV Assay (k131813), one of the 87 specimens was positive for human
 519 metapneumovirus (hMPV) as determined by the FDA-cleared Quidel Molecular RSV + hMPV
 520 Assay (k131813), and one of the 87 specimens was positive for parainfluenza virus type 1 as
 521 determined by the FDA-cleared Gen-probe Prodesse ProParaflu + (k132238).

522 Overall, a total of 26 contrived avian influenza A (H7N9) positive specimens and 87 avian
 523 influenza A (H7N9) negative clinical specimens were tested in the study in a randomized and
 524 blinded fashion. Results of the study are summarized in the table below:

525 **Performance Summary**

Lyra™ Influenza A Subtype H7N9 Assay Result	# of Positives ¹	% Positive Agreement (95% CI)	# of Negatives ¹	% Negative Agreement (95% CI)
	26	100.0 (87.1 – 100.0)	87	100.0 (95.8 -100.0)

526 ¹Proportion of true positives or true negatives correctly identified

527

528 **Analytical Performance**

529

530 **Limit of Detection (LoD)**

531

532 A preliminary Limit of Detection (LoD) study was performed initially using a synthetic DNA
 533 template. The synthetic double-stranded DNA template was synthesized based on the published
 534 HA gene sequence of the A/Hangzhou/1/2013 (H7N9) virus strain, and represents a 177-base
 535 pairs sequence that matches the amplicon sequence generated using the Lyra™ Influenza A
 536 Subtype H7N9 Assay primer pair. Serial dilutions of this synthetic DNA template from 5x10⁵ to
 537 5x10⁰ copies per reaction were tested by the Lyra™ Influenza A Subtype H7N9 Assay. The Lyra™
 538 Influenza A Subtype H7N9 Assay detected five copies of the synthetic DNA template per
 539 reaction.

540 The Limit of Detection (LoD) of the Lyra™ Influenza A Subtype H7N9 Assay was further
 541 determined using a β-propiolactone (BPL) inactivated egg cultured A/Anhui/1/2013 (H7N9) virus
 542 strain. The strain was obtained upon request from the Centers for Disease Control and
 543 Prevention (CDC) after authorization by the World Health Organization (WHO) Pandemic
 544 Influenza Preparedness (PIP) Framework. The virus was titered by the CDC using EID₅₀
 545 methodology and subsequently inactivated with BPL. The virus stock was serially diluted in a
 546 negative nasopharyngeal matrix. Each dilution was extracted using the NucliSENS® easyMAG®
 547 system and tested in replicates of 20 per concentration of virus on the Applied Biosystems® 7500
 548 Fast Dx instrument according to the Lyra™ Influenza A Subtype H7N9 Assay protocol. LoD is
 549 defined as the lowest concentration at which 95% of all replicates tested positive.

550 LoD Determination Study Results (Inactivated Virus)

Targets	Lyra™ Influenza A Subtype H7N9 Assay A/H7 (Ct)			Lyra™ Influenza A Subtype H7N9 Assay PRC (Ct)		
	Virus Concentration (EID ₅₀ /mL)	3.95 x 10 ³	1.98 x 10 ³	9.88 x 10 ²	3.95 x 10 ³	1.98 x 10 ³
1	29.4	32.7	Negative	17.9	18.4	18.2
2	29.5	33.9	Negative	18.1	18.5	18.2
3	33.7	Negative	Negative	18.2	18.4	18.2
4	32.6	30.5	31.3	18.3	18.2	18.3
5	29.3	32.5	Negative	18.2	18.6	18.2
6	29.5	Negative	33.6	17.9	18.3	18.2
7	31.4	32.9	34.7	18.1	18.3	18.0
8	30.1	32.7	31.7	18.0	18.4	18.1
9	32.9	32.3	33.1	18.1	18.2	18.2
10	30.9	30.9	32.7	18.2	18.2	18.3
11	28.1	Negative	33.0	18.2	18.3	17.7
12	31.6	34.9	Negative	18.2	18.5	17.9
13	33.8	30.7	Negative	18.7	18.5	17.9
14	34.3	Negative	Negative	18.4	18.4	18.1
15	31.0	32.9	Negative	18.4	18.4	18.1
16	30.2	31.9	32.0	18.3	18.5	18.1
17	28.5	32.5	Negative	18.1	18.3	18.0
18	29.7	31.7	Negative	18.1	18.4	18.0
19	31.4	32.5	Negative	18.4	18.3	18.0
20	30.4	30.8	Negative	18.1	18.4	18.3
Positive Replicates/Total	20/20	16/20	8/20	20/20	16/20	8/20

Replicates						
------------	--	--	--	--	--	--

551

552 The LoD for the Lyra™ Influenza A Subtype H7N9 Assay testing the BPL inactivated
553 A/Anhui/1/2013 (H7N9) virus was determined to be **3.95 x 10³ EID₅₀/mL**.

554 **Analytical Reactivity – *In silico* Prediction**

555

556 Additional isolates of the Avian Influenza A (H7N9) Virus (detected in China in 2013) are not
557 available for testing in the U.S. *In silico* analysis was carried out to predict analytical reactivity to
558 other isolates of the Avian Influenza A (H7N9) Virus (detected in China in 2013). A total of 13
559 full-length HA gene sequences of different isolates of the Avian Influenza A (H7N9) Virus
560 (detected in China in 2013) were obtained from the National Center for Biotechnology
561 Information (NCBI) for this analysis.

562 The Lyra™ Influenza A Subtype H7N9 Assay primer and probe set demonstrated 100% sequence
563 homology to all the 10 Avian Influenza A (H7N9) Virus (detected in China in 2013) isolates
564 recovered from human patients and the three Avian Influenza A (H7N9) Virus (Detected in China
565 in 2013) isolates recovered from the environment and chicken included in this analysis.

566 **Analytical Specificity – Non-A/H7 Influenza A Viruses**

567

568 Analytical specificity of the Lyra™ Influenza A Subtype H7N9 Assay was evaluated with a panel of
569 22 non-A/H7 influenza A viruses. The FDA-cleared Lyra™ Influenza A + B Assay (k112172) was
570 also performed as a quality control measure for the viruses. All of the influenza A viruses tested
571 were positive in the Lyra™ Influenza A + B Assay and negative in the Lyra™ Influenza A Subtype
572 H7N9 Assay at the concentrations tested. The study results are shown in the table below

573 **Analytical Specificity Study Results – Non-A/H7 Influenza A Viruses**

Virus	Virus Titer in TCID ₅₀ /mL	Lyra™ Influenza A Subtype H7N9 Assay A/H7 (Ct)	Lyra™ Influenza A Subtype H7N9 Assay PRC (Ct)	Lyra™ Influenza A + B Assay Result
A/Mexico/4108/2009 (pdm H1N1)	2.80x10 ⁶	Negative	18.9	Influenza A
A/Perth/16/09 (H1N1)	8.80x10 ⁴	Negative	19.2	Influenza A
A/West Virginia/06/2011 (H3N2v)	1.00x10 ⁴	Negative	18.5	Influenza A
A/Pennsylvania/14/2010 (H3N2v)	1.00x10 ⁴	Negative	18.0	Influenza A

Virus	Virus Titer in TCID ₅₀ /mL	Lyra™ Influenza A Subtype H7N9 Assay A/H7 (Ct)	Lyra™ Influenza A Subtype H7N9 Assay PRC (Ct)	Lyra™ Influenza A + B Assay Result
A/Minnesota/11/2010 (H3N2v)	1.00x10 ⁴	Negative	18.1	Influenza A
A/Kansas/13/2009 (H3N2v)	1.00x10 ⁴	Negative	18.1	Influenza A
A/Indiana/08/2011 (H3N2v)	1.00x10 ⁴	Negative	18.5	Influenza A
A/Indiana/10/2011 (H3N2v)	1.00x10 ⁴	Negative	18.3	Influenza A
A/Wisconsin/07/09 (H3N2)	2.00x10 ⁶	Negative	18.5	Influenza A
A/California/07/09 (H3N2)	2.00x10 ⁶	Negative	22.7	Influenza A
A/Port Chalmers/1/73 (H3N2)	4.50x10 ⁶	Negative	21.3	Influenza A
A/Uruguay/7/16/2007 (H3N2)	1.00x10 ⁶	Negative	21.9	Influenza A
A/New Caledonia/20/1999 (H1N1)	1.00x10 ⁶	Negative	20.3	Influenza A
A/Victoria/3/75 (H3N2)	2.00x10 ⁶	Negative	24.0	Influenza A
A/Denver/1/57 (H1N1)	1.00x10 ⁶	Negative	21.3	Influenza A
A1/Mal/302/54 (H1N1)	1.00x10 ⁶	Negative	21.8	Influenza A
A/Perth/16/2009 (H3N2)	1.00x10 ⁸	Negative	19.7	Influenza A
A/Aichi/2/68 (H3N2)	1.00x10 ⁵	Negative	21.7	Influenza A
A/Wisconsin/67/2005 (H3N2)	1.00x10 ⁵	Negative	21.8	Influenza A
A/NWS/33 (H1N1)	1.00x10 ⁵	Negative	22.4	Influenza A
A/Hong Kong/8/68 (H3N2)	1.00x10 ⁷	Negative	20.3	Influenza A
A/New Jersey/8/76 (H1N1)	1.00x10 ⁶	Negative	22.8	Influenza A

574

575

576 **Analytical Specificity/Cross-Reactivity – Common Respiratory Bacteria, Yeast and Viruses**
 577 Potential cross-reactivity of the Lyra™ Influenza A Subtype H7N9 Assay was evaluated by testing
 578 an extensive list of common respiratory bacteria, yeast, and viruses (a total of 30 viruses, 26
 579 bacteria, and one yeast). No cross-reactivity was observed with any of the 57 organisms at the
 580 concentrations tested. The study results are shown in the tables below:

581 **Cross-Reactivity Study Results –Viruses**

Organism	Concentrations Tested (TCID ₅₀ /mL)	Lyra™ Influenza A Subtype H7N9 Assay A/H7 (Ct)	Lyra™ Influenza A Subtype H7N9 Assay PRC (Ct)
Adenovirus type 1	1.51x10 ⁵	Negative	19.9
Coronavirus 229E	2.46x10 ⁷	Negative	18.4
Coronavirus NL63	1.41x10 ⁵	Negative	21.1
Coronavirus OC43	2.42x10 ⁷	Negative	18.4
Coxsackievirus B4:ODH-42385	2.00x10 ⁷	Negative	20.3
Cytomegalovirus	2.14 x10 ⁶	Negative	18.2
Echovirus 6	1.52 x10 ⁹	Negative	18.6
Echovirus 7	4.58 x10 ⁶	Negative	19.6
Echovirus 9	2.17 x10 ⁷	Negative	19.8
Echovirus 11	2.17x10 ⁶	Negative	20.8
Enterovirus 70	2.41x10 ⁶	Negative	18.0
Enterovirus 71	2.03x10 ⁵	Negative	18.1
Epstein Barr Virus	9.27x10 ⁸ genome equivalents/mL*	Negative	19.2
HSV Type 1 MacIntyre Strain	5.89x10 ⁷	Negative	18.3
HSV Type 2 G strain	1.96x10 ⁶	Negative	18.3
Human Metapneumovirus (A1)	3.66x10 ⁵	Negative	18.6

Organism	Concentrations Tested (TCID ₅₀ /mL)	Lyra™ Influenza A Subtype H7N9 Assay A/H7 (Ct)	Lyra™ Influenza A Subtype H7N9 Assay PRC (Ct)
Human Rhinovirus 45	5.87x10 ⁴	Negative	18.2
Human Rhinovirus 52	5.25x10 ⁴	Negative	18.2
Influenza A/Mexico/4108/2009	4.08x10 ⁵	Negative	18.4
Influenza A/Port Chalmers	3.55x10 ⁸	Negative	18.4
Influenza B/Florida/04/2006	1.54x10 ⁶	Negative	18.5
Measles	1.95x10 ⁷	Negative	18.1
Mumps Virus	2.75x10 ⁹	Negative	18.5
Parainfluenza Type 1	1.58x10 ⁵	Negative	18.8
Parainfluenza Type 2	3.15x10 ⁸	Negative	18.9
Parainfluenza Type 3	2.56x10 ⁷	Negative	20.8
Parainfluenza Type 4A	1.04x10 ⁵	Negative	19.0
RSV A (Long)	4.36x10 ⁴	Negative	18.5
RSV B Strain (Wash/18537/62)	3.43x10 ⁵	Negative	18.2
Varicella Zoster Virus	1.11x10 ⁴	Negative	18.3

582 *Quantified by a molecular method

583

584 **Cross-Reactivity Study Results –Bacteria and Yeast**

Organism	Concentrations Tested (CFU/mL)	Lyra™ Influenza A Subtype H7N9 Assay A/H7 (Ct)	Lyra™ Influenza A Subtype H7N9 Assay PRC (Ct)
<i>Bordetella pertussis</i>	9.08 x10 ⁸	Negative	19.7
<i>Bordetella bronchiseptica</i>	5.40 x10 ⁸	Negative	19.5
<i>Chlamydomphila pneumoniae</i>	2.2 ug/mL (DNA)	Negative	18.1
<i>Chlamydomphila trachomatis</i>	2.10x10 ⁶	Negative	19.1
<i>Legionella pneumophila</i>	1.42x10 ⁹	Negative	19.8
<i>Mycobacterium intracellulare</i> (ATCC 95-06)	1.53x10 ⁹	Negative	19.1
<i>Mycobacterium tuberculosis</i>	9.30x10 ⁶	Negative	19.1
<i>Mycobacterium avium</i> (ATCC 25291)	3.18x10 ⁹	Negative	20.2
<i>Mycoplasma pneumoniae</i>	3.16x10 ⁷	Negative	18.6
<i>Haemophilus influenzae</i>	4.00x10 ⁸	Negative	19.4
<i>Pseudomonas aeruginosa</i>	1.32x10 ⁹	Negative	20.9
<i>Proteus vulgaris</i>	6.53x10 ⁸	Negative	19.8
<i>Proteus mirabilis</i>	1.19x10 ⁹	Negative	20.7
<i>Neisseria gonorrhoeae</i>	1.40x10 ⁹	Negative	21.3
<i>Neisseria meningitidis</i>	1.29x10 ⁸	Negative	19.7

Organism	Concentrations Tested (CFU/mL)	Lyra™ Influenza A Subtype H7N9 Assay A/H7 (Ct)	Lyra™ Influenza A Subtype H7N9 Assay PRC (Ct)
<i>Neisseria mucosa</i>	1.61x10 ⁹	Negative	21.5
<i>Klebsiella pneumoniae</i>	9.75x10 ⁸	Negative	20.8
<i>Escherichia coli</i> (ATCC 43895)	1.13x10 ⁹	Negative	21.2
<i>Moraxella catarrhalis</i> (ATCC 8176)	1.26x10 ⁹	Negative	20.8
<i>Corynebacterium diphtheriae</i> (ATCC 19409)	3.44x10 ⁸	Negative	18.5
<i>Lactobacillus plantarum</i>	3.18x10 ⁸	Negative	18.4
<i>Streptococcus pneumoniae</i> (ATCC 6305)	1.43x10 ⁸	Negative	19.3
<i>Streptococcus pyogenes</i> (ATCC 9898)	6.38x10 ⁸	Negative	18.5
<i>Streptococcus salivarius</i>	5.40x10 ⁸	Negative	18.3
<i>Staphylococcus epidermidis</i>	9.23x10 ⁸	Negative	18.3
<i>Staphylococcus aureus</i> (ATCC 12598)	6.08x10 ⁸	Negative	18.2
<i>Candida albicans</i>	9.70x10 ⁷	Negative	18.2

585

586 **Microbial Interference – Common Respiratory Bacteria, Yeast and Viruses**

587

588 The performance of the Lyra™ Influenza A Subtype H7N9 Assay was also evaluated with
 589 potentially interfering common respiratory organisms. The same organisms that were tested in
 590 the Cross-Reactivity Study (described above) were used in this Microbial Interference Study. The
 591 potentially interfering organisms were evaluated with the BPL inactivated A/Anhui/1/2013

592 (H7N9) virus (as described in the LoD section) at a concentration of approximately 1x LoD using
 593 the Lyra™ Influenza A Subtype H7N9 Assay. No interference was observed with any of the 57
 594 organisms at the concentrations tested. The study results are shown in the tables below:

595 **Microbial Interference Study Results –Viruses**

Organism	Concentrations Tested (TCID ₅₀ /mL)	Lyra™ Influenza A Subtype H7N9 Assay A/H7 (Ct)	Lyra™ Influenza A Subtype H7N9 Assay PRC (Ct)
Adenovirus type 1	1.36x10 ⁵	28.7	18.4
Coronavirus 229E	2.21x10 ⁷	28.2	18.1
Coronavirus NL63	1.27x10 ⁵	28.9	18.6
Coronavirus OC43	2.18x10 ⁷	29.8	18.3
Coxsackievirus B4:ODH-42385	1.80x10 ⁷	31.4	19.2
Cytomegalovirus	1.93x10 ⁶	27.6	18.2
Echovirus 6	1.37x10 ⁹	31.2	18.2
Echovirus 7	4.12x10 ⁶	32.0	19.2
Echovirus 9	1.95x10 ⁷	32.6	18.4
Echovirus 11	1.95x10 ⁶	32.8	19.3
Enterovirus 70	2.17x10 ⁶	31.1	18.0
Enterovirus 71	1.83x10 ⁵	29.9	18.0
Epstein Barr Virus	8.34x10 ⁸ genome equivalents/mL*	29.9	18.2
HSV Type 1 MacIntyre Strain	5.30x10 ⁷	28.8	18.2
HSV Type 2 G strain	1.76x10 ⁶	29.7	18.2
Human Metapneumovirus (A1)	3.29x10 ⁵	29.6	18.1
Human Rhinovirus 45	5.28x10 ⁴	32.4	18.2

Organism	Concentrations Tested (TCID ₅₀ /mL)	Lyra™ Influenza A Subtype H7N9 Assay A/H7 (Ct)	Lyra™ Influenza A Subtype H7N9 Assay PRC (Ct)
Human Rhinovirus 52	4.73x10 ⁴	30.5	18.2
Influenza A/Mexico/4108/2009	3.67x10 ⁵	31.1	18.3
Influenza A/Port Chalmers	3.20x10 ⁸	34.8	18.2
Influenza B/Florida/04/2006	1.39x10 ⁶	29.6	18.1
Measles	1.76x10 ⁷	28.8	18.2
Mumps Virus	2.48x10 ⁹	27.0	17.9
Parainfluenza Type 1	1.42x10 ⁵	29.9	18.2
Parainfluenza Type 2	2.84x10 ⁸	31.3	18.3
Parainfluenza Type 3	2.30x10 ⁷	34.0	18.8
Parainfluenza Type 4A	9.36x10 ⁴	29.8	18.1
RSV A (Long)	3.92x10 ⁴	30.0	18.2
RSV B Strain (Wash/18537/62)	3.09x10 ⁵	29.2	18.3
Varicella Zoster Virus	9.99x10 ³	27.0	18.3

596 *Quantified by a molecular method

597

598

599

600

601

602 **Microbial Interference Study Results –Bacteria and Yeast**

Organism	Concentrations Tested (CFU/mL)	Lyra™ Influenza A Subtype H7N9 Assay A/H7 (Ct)	Lyra™ Influenza A Subtype H7N9 Assay PRC (Ct)
<i>Bordetella pertussis</i>	8.17x10 ⁸	28.0	18.9
<i>Bordetella bronchiseptica</i>	4.86x10 ⁸	30.6	18.9
<i>Chlamydomphila pneumoniae</i>	1.98 ug/mL (DNA)	28.8	18.3
<i>Chlamydomphila trachomatis</i>	1.89x10 ⁶	33.8	18.5
<i>Legionella pneumophila</i>	1.28x10 ⁹	30.9	19.0
<i>Mycobacterium intracellulare</i> (ATCC 95-06)	1.38x10 ⁹	29.5	18.3
<i>Mycobacterium tuberculosis</i>	8.37x10 ⁶	29.4	18.5
<i>Mycobacterium avium</i> (ATCC 25291)	2.86x10 ⁹	29.5	18.7
<i>Mycoplasma pneumoniae</i>	2.84x10 ⁷	31.4	18.4
<i>Haemophilus influenzae</i>	3.60x10 ⁸	31.3	19.0
<i>Pseudomonas aeruginosa</i>	1.19x10 ⁹	29.7	19.5
<i>Proteus vulgaris</i>	5.88x10 ⁸	30.6	18.4
<i>Proteus mirabilis</i>	1.07x10 ⁹	29.7	19.3
<i>Neisseria gonorrhoeae</i>	1.26x10 ⁹	30.1	19.5
<i>Neisseria meningitidis</i>	1.16x10 ⁸	28.4	18.7

Organism	Concentrations Tested (CFU/mL)	Lyra™ Influenza A Subtype H7N9 Assay A/H7 (Ct)	Lyra™ Influenza A Subtype H7N9 Assay PRC (Ct)
<i>Neisseria mucosa</i>	1.45x10 ⁹	32.2	19.8
<i>Klebsiella pneumoniae</i>	8.78x10 ⁸	27.6	19.1
<i>Escherichia coli</i> (ATCC 43895)	1.02x10 ⁹	30.3	19.1
<i>Moraxella catarrhalis</i> (ATCC 8176)	1.13x10 ⁹	30.0	18.7
<i>Corynebacterium diphtheriae</i> (ATCC 19409)	3.10x10 ⁸	27.9	18.5
<i>Lactobacillus plantarum</i>	2.86x10 ⁸	30.0	18.5
<i>Streptococcus pneumoniae</i> (ATCC 6305)	1.29x10 ⁸	30.2	18.8
<i>Streptococcus pyogenes</i> (ATCC 9898)	5.74x10 ⁸	29.1	18.6
<i>Streptococcus salivarius</i>	4.86x10 ⁸	30.6	18.4
<i>Staphylococcus epidermidis</i>	8.31x10 ⁸	29.2	18.3
<i>Staphylococcus aureus</i> (ATCC 12598)	5.47x10 ⁸	27.7	18.5
<i>Candida albicans</i>	8.73x10 ⁷	27.5	18.3

603

604

605

606

607 Customer and Technical Assistance

608 To place an order or for technical support, please contact a Quidel Representative at (800) 874-1517 (toll-free
609 in the U.S.A.) or (858) 552-1100, Monday through Friday, between 8:00 a.m. and 5:00 p.m., Eastern Time.
610 Orders may also be placed by fax at (740) 592-9820. For e-mail support contact: customer
611 service@dhiusa.com or technical_services@dhiusa.com. For services outside the U.S., please contact your
612 local distributor. Additional information about Quidel, our products, and our distributors can be found on our
613 website quidel.com.

614

615 Trademarks

616 NucliSENS and easyMAG are registered trademarks of bioMérieux, Inc. TaqMan is a registered trademark of
617 Roche. Applied Biosystems® is a registered trademark of Life Technologies. Dye compounds in this product
618 are sold under license from BioSearch Technologies, Inc. and protected by U.S. and world-wide patents either
619 issued or under application. The license covers R&D use and human in vitro diagnostic (IVD) applications.

ⁱ Clinical and Laboratory Standards Institute. Viral Culture; Approved Guidelines. CLSI document M41-A [ISBN 1562386239] Clinical and Laboratory Standards Institute, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898, USA 2006.