

Contains Nonbinding Recommendations

Technical Performance Assessment of Digital Pathology Whole Slide Imaging Devices

Guidance for Industry and Food and Drug Administration Staff

Document issued on: April 20, 2016

The draft of this document was issued on February 25, 2015

For questions about this document, contact the Division of Molecular Genetics and Pathology at 301-796-6179 and Nicholas Anderson at 301-796-4310 or nicholas.anderson@fda.hhs.gov or Aldo Badano at 301-796-2534 or aldo.badano@fda.hhs.gov.



U.S. Department of Health and Human Services
Food and Drug Administration
Center for Devices and Radiological Health

Office of *In Vitro* Diagnostics and Radiological Health
Division of Molecular Genetics and Pathology
Molecular Pathology and Cytology Branch

Preface

Public Comment

You may submit electronic comments and suggestions at any time for Agency consideration to <http://www.regulations.gov> . Submit written comments to the Division of Dockets Management, Food and Drug Administration, 5630 Fishers Lane, Room 1061, (HFA-305), Rockville, MD 20852. Identify all comments with the docket number [FDA-2015-D-0230]. Comments may not be acted upon by the Agency until the document is next revised or updated.

Additional Copies

Additional copies are available from the Internet. You may also send an e-mail request to CDRH-Guidance@fda.hhs.gov to receive a copy of the guidance. Please use the document number 1400053 to identify the guidance you are requesting.

Table of Contents

I.	Introduction	1
II.	Background	2
III.	Scope	2
IV.	Policy.....	3
IV(A).	Description and Test Methods for Each Component.....	3
IV(A)(1).	Slide Feeder	5
IV(A)(1)(a).	Description	5
IV(A)(2).	Light Source	5
IV(A)(2)(a).	Description	5
IV(A)(2)(b).	Test Method	6
IV(A)(3).	Imaging Optics	6
IV(A)(3)(a).	Description.....	6
IV(A)(3)(b).	Test Methods	7
IV(A)(4).	Mechanical Scanner Movement	7
IV(A)(4)(a).	Description	7
IV(A)(4)(b).	Test Method	8
IV(A)(5).	Digital Imaging Sensor	8
IV(A)(5)(a).	Description.....	8
IV(A)(5)(b).	Test Methods	8
IV(A)(6).	Image Processing Software	9
IV(A)(6)(a).	Description.....	9
IV(A)(6)(b).	Resources.....	9
IV(A)(7).	Image Composition	9
IV(A)(7)(a).	Description.....	9
IV(A)(7)(b).	Test Methods	10
IV(A)(8).	Image Files Formats	10
IV(A)(8)(a).	Description.....	10
IV(A)(9).	Image Review Manipulation Software	11
IV(A)(9)(a).	Description.....	11
IV(A)(9)(b).	Resources.....	11

Contains Nonbinding Recommendations

IV(A)(10). Computer Environment	11
IV(A)(10)(a). Description	11
IV(A)(11). Display	12
IV(A)(11)(a). Description	12
IV(A)(11)(b). Test Methods	12
IV(A)(11)(c). Resources	13
IV(B). System-level Assessment.....	14
IV(B)(1). Color Reproducibility	15
IV(B)(1)(a). Description	15
IV(B)(1)(b). Test Methods.....	15
IV(B)(1)(c). Resources	16
IV(B)(2). Spatial Resolution.....	16
IV(B)(2)(a). Description	16
IV(B)(2)(b). Test Methods.....	16
IV(B)(3). Focusing Test.....	16
IV(B)(4). Whole Slide Tissue Coverage.....	17
IV(B)(4)(a). Description	17
IV(B)(4)(b). Test Method	17
IV(B)(5). Stitching Error	18
IV(B)(5)(a). Description	18
IV(B)(5)(b). Test Methods.....	18
IV(B)(6). Turnaround Time.....	19
IV(B)(6)(a). Description	19
IV(C). User Interface.....	19
IV(C)(1). Description.....	19
IV(C)(2). Test Methods	19
IV(C)(3). Resources.....	22
IV(D). Labeling	22
IV(D)(1). Test Methods	23
IV(D)(2). Resources.....	23
IV(E). Quality Control	23

1 **Technical Performance Assessment** 2 **of Digital Pathology Whole Slide** 3 **Imaging Devices**

4 **Guidance for Industry and Food** 5 **and Drug Administration Staff**

8 *This guidance represents the current thinking of the Food and Drug Administration*
9 *(FDA or Agency) on this topic. It does not establish any rights for any person and is*
10 *not binding on FDA or the public. You can use an alternative approach if it satisfies*
11 *the requirements of the applicable statutes and regulations. To discuss an alternative*
12 *approach, contact the FDA staff or Office responsible for this guidance as listed on the*
13 *title page.*

14 **I. Introduction**

15 FDA is issuing this guidance to provide industry and agency staff with recommendations
16 regarding the technical performance assessment data that should be provided for
17 regulatory evaluation of a digital whole slide imaging (WSI) system. This document
18 does not cover the clinical submission data that may be necessary to support approval or
19 clearance. This document provides our suggestions on how to best characterize the
20 technical aspects that are relevant to WSI performance for their intended use and
21 determine any possible limitations that might affect their safety and effectiveness.
22 23

24 Recent technological advances in digital microscopy, in particular the development of
25 whole slide scanning systems, have accelerated the adoption of digital imaging in
26 pathology, similar to the digital transformation that radiology departments have
27 experienced over the last decade. FDA regulates WSI system manufacturers to help
28 ensure that the images intended for clinical uses are reasonably safe and effective for
29 such purposes. Essential to the regulation of these systems is the understanding of the
30 technical performance of the WSI system and the components in the imaging chain, from
31 image acquisition to image display and their effect on pathologist's diagnostic
32 performance and workflow. Prior to performing non-technical analytical studies (i.e.,
33 those using clinical samples) and clinical studies to evaluate a digital imaging system's
34 performance, the manufacturer should first determine the technical characteristics that are
35 relevant to such performance for its intended use and determine any possible limitations
36

Contains Nonbinding Recommendations

37 that might affect its safety and effectiveness. This guidance provides recommendations
38 for the assessment of technical characteristics of a WSI device.

39
40 FDA's guidance documents, including this guidance, do not establish legally enforceable
41 responsibilities. Instead, guidances describe the Agency's current thinking on a topic and
42 should be viewed only as recommendations, unless specific regulatory or statutory
43 requirements are cited. The use of the word *should* in Agency guidance means that
44 something is suggested or recommended, but not required.

45

46 **II. Background**

47

48 For over a hundred years, the reference method for the diagnosis of cancer and many
49 other critical clinical conditions has been histopathological examination of tissues using
50 conventional light microscopy. This process is known as surgical pathology in the
51 United States.

52

53 In surgical pathology, patient tissue from surgery, biopsy or autopsy goes through a
54 process that includes dissection, fixation, embedding, and cutting of tissue into very thin
55 slices which are then stained, for example by the hematoxylin and eosin (H&E) protocol,
56 and permanently mounted onto glass slides. The slides are examined by a pathologist
57 under a light microscope by dynamically adjusting the focus and using different
58 magnifications. By integrating their interpretations obtained by microscopic examination
59 of the tissue from all slides pertaining to a case, pathologists arrive at a diagnosis of the
60 case.

61

62 WSI refers to the digitization of the stained entire tissue specimen on a glass slide. The
63 glass slide is still prepared and stained just as for conventional light microscopy.

64 Depending on the system used, various magnifications, scanning methodologies,
65 hardware, and software are employed to convert the optical image of the slide into a
66 digital whole slide image. With WSI, the pathologist views the image on a computer
67 monitor rather than through the microscope oculars.

68

69 **III. Scope**

70

71 This document provides guidance regarding only the technical performance assessment
72 of WSI systems for regulatory evaluation. WSI systems are defined here as those
73 consisting of (a) an image acquisition subsystem that converts the content of a glass slide
74 into a digital image file, and (b) a workstation environment for viewing the digital
75 images. If not otherwise specified, the term “image” in the context of whole slide
76 imaging refers to a pyramid structure consisting of multiple images at different
77 resolutions. The baseline image has the highest resolution. This guidance is applicable
78 for surgical pathology tasks performed in the anatomic pathology laboratory. It is
79 intended to provide recommendations to industry and FDA staff regarding only the
80 technical performance assessment data needed for the regulatory evaluation of a WSI
81 device. This document is not meant to provide guidance for special stain techniques or

Contains Nonbinding Recommendations

82 fluorescence imaging or for the non-technical analytical studies (utilizing clinical
83 samples) or pivotal clinical studies necessary to support safety and effectiveness, nor
84 does this guidance alone suffice to demonstrate safety and effectiveness of WSI systems.
85 Interpretation of WSI images on mobile platforms is beyond the scope of this guidance.
86

IV. Policy

87
88
89 The following subsections of this section describe the technical performance assessment
90 data FDA believes will facilitate the regulatory evaluation of a WSI device.
91


IV(A). Description and Test Methods for Each Component

92
93
94 This subsection details the descriptions and the test methods at the component level that
95 should be included in the technical performance assessment of a WSI device. For
96 purposes of this guidance only, a component is a piece of hardware, software, or a
97 combination of hardware and software that processes the image signals flowing through
98 the imaging chain. The concept of a component is based on the transformation of the
99 image signals. For example, the digital imaging sensor is a hardware device that converts
100 optical signals into digital signals. The image composition component is a software
101 program that stitches sub-images together to form a whole slide image. A component
102 and a physical device need not be in close physical proximity. For example, the light
103 source component and the image optics component are usually tightly coupled within the
104 same device, while the display calibration data is often distributed in both the color
105 profile in the computer environment component and the on-screen display settings in the
106 display component.
107

108 The components in a WSI device can be grouped in two subsystems: image acquisition
109 and image display. The image acquisition subsystem digitizes the tissue slide as a digital
110 image file. The image display subsystem converts the digital image file into optical
111 signals for the human reader. In the paradigm of telemedicine, the digital image file can
112 be electronically sent to a remote site for reading, so the image acquisition subsystem and
113 the image display subsystem do not need to be physically coupled. Methods for
114 independently testing the image acquisition and display subsystems are described in
115 Section IV(B).
116

117 Sponsors should provide a block diagram of the components found in the WSI system in
118 the premarket submission. A chart indicating the relationship among the components and
119 the test methods utilized for the specific system characterization should also be provided.
120 Diagram 1 on the following page is offered as an example block diagram of typical
121 components found in current WSI systems. The components of a particular WSI system
122 might not include all of those listed in the diagram or may include additional
123 components. Sponsors are encouraged to provide additional diagrams, illustrations, and
124 photographs of their devices as part of their submissions.
125
126

Diagram 1: Example block diagram of typical components found in current WSI systems



127
128
129
130
131
132
133
134
135
136
137
138
139
140
141
142
143
144
145
146
147
148
149
150
151
152
153
154
155
156
157
158
159
160
161
162
163
164
165
166
167

Contains Nonbinding Recommendations

168 **IV(A)(1). Slide Feeder**

169

170 **IV(A)(1)(a). Description**

171

172 The slide feeder is the mechanism(s) used to introduce the slide(s) to the scanner. For the
173 slide feeder, sponsors should provide the following information, if applicable:

- 174 • Configuration of the slide feed mechanism (a physical description of the
175 equipment)
- 176 ○ Slide configuration (physical description of the slide (i.e., custom or
177 commercial off-the-shelf))
 - 178 ○ Number of slides in queue (carrier)
 - 179 ○ Class of automation (e.g., robotics, pneumatics, etc.)
- 180 • User interaction
- 181 ○ Hardware (e.g., loading of slides into carrier)
 - 182 ○ Software (e.g., does the system recognize the number of slides or is this
183 specified by the user)
 - 184 ○ Feedback (e.g., alarms, notifications, etc.)
 - 185 ○ Failure Mode and Effects Analysis (FMEA) (including severity,
186 likelihood, mitigations, etc.)
- 187

188 **IV(A)(2). Light Source**

189

190 **IV(A)(2)(a). Description**

191

192 The light source, including the light guide, generates and delivers light to the slide being
193 imaged. The two major components are the lamp and condenser. For the light source,
194 sponsors should provide the following information and specifications, if applicable:

- 195 • Lamp
- 196 ○ Bulb type (e.g., halogen, xenon arc, LED)
 - 197 ○ Manufacturer and model
 - 198 ○ Wattage
 - 199 ○ Spectral power distribution
 - 200 ○ Expected lifetime
 - 201 ○ Output adjustment control (electrical/electronic/mechanical)
 - 202 ○ Optical filter(s)
 - 203 ▪ Type (e.g., heat blocking, polarization, neutral density, diffusing)
 - 204 ○ Manufacturer and model
 - 205 ○ Expected intensity variation (coefficient of variation)
 - 206 ▪ Over the duration of scanning a single slide
 - 207 ▪ Over the course of a single workday
 - 208 ▪ Over the lifetime of the device
 - 209 ○ Expected spectral variation
 - 210 ▪ Over the duration of scanning a single slide
 - 211 ▪ Over the course of a single workday
 - 212 ▪ Over the lifetime of the device
 - 213 ○ Capability of tracking intensity and spectral degradation with lifetime

Contains Nonbinding Recommendations

- 214 • Condenser
- 215 ○ Illumination format (e.g., Kohler, critical)
- 216 ○ Manufacturer and model
- 217 ○ Numerical aperture
- 218 ○ Focal length
- 219 ○ Working distance

220

IV(A)(2)(b). Test Method

222

223 The following steps should be used to measure the spectral distribution of light incident
224 on the slide. Position the input of a calibrated spectrometer or monochromator at the
225 plane where the slide would be placed, centered on the illumination spot from the
226 condenser. If desired, the light can be coupled into the spectrometer via light guide (e.g.,
227 fiber optic cable) or an integrating sphere. The measurement aperture should be at least
228 as large as the anticipated field of view on the slide at the lowest magnification of the
229 imaging optics. The wavelength accuracy and relative spectral efficiency of the
230 spectrometer or monochromator in the wavelength range of 360-830 nm should be
231 calibrated prior to measurements and reported. Plots of the measured spectrum with at
232 least 10 nm spectral resolution should be provided, using radiometric units (e.g., spectral
233 irradiance in W/cm²/nm, spectral radiance in W/sr/cm²/nm).

234

IV(A)(3). Imaging Optics

236

IV(A)(3)(a). Description

238

239 The imaging optics comprises the microscope objective and auxiliary lens(es) (e.g., tube
240 lens), which optically transmit an image of the tissue from the slide to the digital image
241 sensor. Sponsors should provide the following information and specifications, if
242 applicable:

- 243 • Optical schematic with all optical elements identified from slide (object plane) to
244 digital image sensor (image plane)
- 245 • Microscope objective
 - 246 ○ Manufacturer
 - 247 ○ Type
 - 248 ○ Magnification
 - 249 ○ Numerical aperture (NA)
 - 250 ○ Focal length
 - 251 ○ Working distance
- 252 • Auxiliary lens(es)
 - 253 ○ Manufacturer
 - 254 ○ Lens type
 - 255 ○ Focal length
- 256 • Magnification of imaging optics: ISO 8039:2014 *Optics and optical instruments*
257 — *Microscopes — Magnification*

258

259

Contains Nonbinding Recommendations

260 **IV(A)(3)(b). Test Methods**

261

262 Sponsors should conduct the following tests in conformance with the International

263 Standards, if applicable:

264 • Relative irradiance of imaging optics at image plane per ISO 13653:1996 *Optics*
265 *and optical instruments – General optical test methods - Measurement of relative*
266 *irradiance in the image field*

267 • Distortion per ISO 9039:2008 *Optics and photonics — Quality evaluation of*
268 *optical systems —Determination of distortion*

269 • Chromatic aberrations per ISO 15795:2002 *Optics and optical instruments —*
270 *Quality evaluation of optical systems — Assessing the image quality degradation*
271 *due to chromatic aberrations*

272

273 **IV(A)(4). Mechanical Scanner Movement**

274

275 **IV(A)(4)(a). Description**

276

277 The mechanical scanner addresses the physical characteristics of the stage upon which
278 the glass slide is affixed. The key components include stage configuration, movement,
279 and control. This information is relevant whether it is only the stage that is moving and
280 the optics are stationary, or if there is movement on all axes. For the mechanical scanner,
281 sponsors should provide the following information and specifications, if applicable:

282 • Configuration of the stage (a physical description of the stage)

283

○ Stage size

284

○ Stage manufacturer and model number

285

○ Stage material (e.g., anodized aluminum)

286

○ Single multi-axis or multiple stacked linear stages (manufacturer and
287 model number)

288

○ Type of guides or ways (e.g., bearings)

289

○ Sample retention mechanism (slide holder)

290

• Method of movement of the stage (e.g., stepper motor, servomotor, piezomotor,
291 etc., coupled with belt, ball-screw, lead-screw, etc.)

292

○ Movement resolution for XY-axes

293

○ Movement in Z-axis

294

○ Speed range

295

○ Travel distance

296

○ Maximum scanning area

297

○ Localization and reading of bar code labels

298

• Control of movement of the stage

299

○ Open or closed loop operation

300

○ Positional accuracy (calibration) and repeatability

301

▪ Lost motion compensation (e.g., backlash)

302

○ Physical control (e.g., joystick) for single-slide, non-batch mode

303

○ Selection of area to be scanned (in accordance to image composition
304 software)

305

▪ whole slide

Contains Nonbinding Recommendations

- 306 ▪ automatically determined area with tissue content
307 • Failure Mode and Effects Analysis (FMEA) (including severity, likelihood,
308 mitigations, etc.)
309

IV(A)(4)(b). Test Method

310
311
312 Sponsors should demonstrate the mechanical performance of the stage with respect to
313 positional repeatability and accuracy on all relevant axes, in accordance with ISO 230-
314 2:2014 Test code for machine tools—Part 2: *Determination of accuracy and*
315 *repeatability of positioning numerically controlled axes.*
316

IV(A)(5). Digital Imaging Sensor

IV(A)(5)(a). Description

317
318
319
320
321 The digital image sensor is an array of photosensitive elements (pixels) that convert the
322 optical signals of the slide to digital signals, which consist of a set of values
323 corresponding to the brightness and color at each point in the optical image. Please
324 provide the following information and specifications:

- 325 • Sensor type (e.g., CMOS, CCD) and manufacturer
326 • Pixel information/specifications
327 ○ Number and dimensions of pixels
328 ○ Design of color filter array
329 ▪ Configuration of color filter array
330 ▪ Spectral transmittance of color filter mask
331 • Responsivity specifications
332 ○ Relative response versus wavelength
333 ○ Linearity
334 ○ Spatial uniformity
335 • Noise specifications
336 ○ Dark current level (electrons per second)
337 ○ Read noise (electrons)
338 • Readout rate (e.g., pixels per second, frames per second)
339 • Digital output format (e.g., bits per pixel, bits per color channel)
340

IV(A)(5)(b). Test Methods

341
342
343 Sponsors should conduct the following tests in conformance with the corresponding
344 International Standards, if applicable:
345

- 346 • Opto-electronic conversion function per ISO 14524:2009 *Photography —*
347 *Electronic still-picture cameras — Methods for measuring optoelectronic*
348 *conversion functions (OECFs)*
349 • Noise measurements per ISO 15739:2013 *Photography — Electronic still-picture*
350 *imaging — Noise measurements*
351

352 **IV(A)(6). Image Processing Software**

353

354 **IV(A)(6)(a). Description**

355

356 Image processing software refers to the embedded software components of the image
357 acquisition device. It typically includes control algorithms for image capture and
358 processing algorithms for raw data conversion into the digital image file. Sponsors
359 should provide the following information and specifications, if applicable:

360

- Exposure control
- White balance
- Color correction
- Sub-sampling
- Pixel-offset correction
- Pixel-gain or flat-field correction
- Pixel-defect correction

361

362

363

364

365

366

367

368 **IV(A)(6)(b). Resources**

369

370 See the guidance entitled “*Guidance for the Content of Premarket Submissions for*
371 *Software Contained in Medical Devices*”

372 ([http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocument](http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/ucm089543.htm)
373 [s/ucm089543.htm](http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/ucm089543.htm)) for the information that should be provided.

374

375 **IV(A)(7). Image Composition**

376

377 **IV(A)(7)(a). Description**

378

379 Image composition is a step present in systems that produce whole slide images as
380 opposed to individual fields of view. Whole slide scanning is typically performed in
381 accordance with the positioning of a stage that moves in submicron steps. At each
382 location of the stage movement, an image of the field of view is acquired. Images can be
383 acquired with a degree of overlapping (redundancy) between them to avoid gaps in data
384 collection. Images can also be acquired at different depths of focus followed by the
385 application of focusing algorithms. At the end of this process, all acquired images are
386 combined (stitched) together to create a composite high resolution image. There are a
387 number of features that can affect this process, and they are listed below. Sponsors
388 should provide a description of these features, if applicable:

389

- Scanning method
 - Single objective or multiple miniature objectives in an array pattern
 - Scanning pattern: square matrix acquisition (tiling), line scanning, etc.
 - Overlap between scanned regions
 - Merging algorithms that stitch the aligned images together into a composite image file. Such algorithms may employ functions to align adjacent fields of view in accordance to the scanning pattern, overlap, etc.

390

391

392

393

394

395

Contains Nonbinding Recommendations

- 396 ○ Automatic background correction functions to eliminate the effect of non-
- 397 uniformities in the microscope's illumination and image merging
- 398 procedure. These non-uniformities if not corrected might create visible
- 399 borders (seams and stitch lines) between the adjacent fields of view.
- 400 ● Scanning speed: time to scan the whole slide. This time is dependent on selected
- 401 magnification, and the amount of tissue on the glass slide.
- 402 ● Number of planes at the Z-axis to be digitized (stack depth)

IV(A)(7)(b). Test Methods

406 Testing for image composition can be performed on a system level using special
407 calibration slides (such as grid patterns) that can test for line uniformity and focus
408 quality. Sponsors should provide the following outputs for these tests, if applicable:

- 409 ● Images of digitized calibration slides
- 410 ● Analysis of focus quality metrics
- 411 ● Analysis of coverage of the image acquisition for the entire tissue slide

IV(A)(8). Image Files Formats

IV(A)(8)(a). Description

417 The final result from image acquisition can be a whole slide image consisting of a stack
418 of all acquired fields of view and magnifications during WSI. The complete digitized
419 image file usually occupies between 1-20 gigabytes of storage space depending on the
420 sample and the magnification of the objective lens used. Images can then be stored in a
421 number of ways and formats. Sponsors should provide the following information:

- 422 ● Compression method (e.g., the wavelet-based JPEG2000 compression standard or
423 TIFF)
- 424 ● Compression ratio: ratio of uncompressed to compressed file size. This metric
425 should be provided along with descriptive information on the data it was
426 measured from, since compression ratio is dependent on the content of the data
427 applied to.
- 428 ● Compression type: lossless or lossy compression
- 429 ● File format: can be formats easily accessible with public domain software such as
430 JPEG or TIFF, or can be proprietary formats only accessible with specific vendor
431 viewers. The file format depends on the file organization and related use.
- 432 ● For systems that interact with DICOM-compliant software and hardware,
433 sponsors should provide a DICOM compatibility report.
- 434 ● File organization:
 - 435 ○ Single file with multi-resolution information (pyramidal organization)
 - 436 ○ Stack of files at different magnifications

Contains Nonbinding Recommendations

440 **IV(A)(9). Image Review Manipulation Software**

441

442

IV(A)(9)(a). Description

443

444 For the image review manipulation software, sponsors should provide the following

445 information, describing software features, if applicable.

446 • Continuous panning (moving in x-y space) and pre-fetching (buffering adjacent
447 images to speed up panning time)

448 • Continuous zooming (magnification)

449 • Discrete Z-axis displacement

450 • Ability to compare multiple slides simultaneously on multiple windows

451 • Ability to perform annotations

452 • Image enhancement such as sharpening functions

453 • Color manipulation, including color profile, white balance, color histogram
454 manipulation, and color filters

455 • Annotation tools

456 • Tracking of visited areas and annotations

457 • Digital bookmarks (revisit selected regions of interest)

458 • Virtual “multihead microscope” (this is when multiple pathologists
459 simultaneously review the same areas remotely)

460

461

IV(A)(9)(b). Resources

462

463 See the guidance entitled “*Guidance for the Content of Premarket Submissions for*
464 *Software Contained in Medical Devices*”

465 ([http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocument](http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/ucm089543.htm)
466 [s/ucm089543.htm](http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/ucm089543.htm)) for additional information on this subject.

467

468 **IV(A)(10). Computer Environment**

469

470

IV(A)(10)(a). Description

471

472 Computer environment refers to the workstation, including both hardware and software
473 components, that retrieves the digital image file and drives the display for the user to

474 review the images. Sponsors should provide the following information and

475 specifications, if applicable:

476 • Computer hardware

477 • Operating system

478 • Graphics card

479 • Graphics card driver

480 • Color management settings

481 • Color profile

482 • Display interface (e.g., DVI or DisplayPort)

483

Contains Nonbinding Recommendations

484 **IV(A)(11). Display**

485

486

IV(A)(11)(a). Description

487

488 The final stage of a WSI system is the display component that presents the scanned image
489 to the pathologists for reading. Technically, display refers to the optoelectronic device
490 that converts the digital image signals in the RGB space into optical image signals. For
491 the display, sponsors should provide the following information and specifications, if
492 applicable:

493

- 494 • Technological characteristics of the display device (e.g., in-plane switching LCD
panel with TFT active-matrix array with fluorescent backlight)
- 495 • Physical size of the viewable area and aspect ratio
- 496 • For transmissive displays, backlight type and properties including temporal,
497 spatial, and spectral characteristics
- 498 • Frame rate and refresh rate
- 499 • Pixel array, pitch, pixel aperture ratio and subpixel matrix scheme (e.g., chevron,
500 RGBW)
- 501 • Subpixel driving to improve grayscale resolution (e.g., spatial and temporal
502 dithering)
- 503 • Supported color spaces
- 504 • Display Interface
- 505 • User controls of brightness, contrast, gamma, color space, power-saving options,
506 etc. via the on-screen display (OSD) menu
- 507 • Ambient light adaptation including the ambient light sensing method,
508 instrumentation, and software tool description
- 509 • Touch screen technology including method, functionality, and any calibration or
510 periodical re-tuning requirements
- 511 • Color calibration tools (sensor hardware and associated software), color profile,
512 and method for color management
- 513 • Frequency and nature of quality-control tests to be performed by the user and/or
514 the physicist with associated action limits.

515

516

IV(A)(11)(b). Test Methods

517

518

- 519 • **User controls:** Modes and settings of the display undergoing testing should be
520 specified, including brightness, contrast, gamma, white point, color space, etc.
See *2.1 Modified-Performance Modes, IDMS 1.03*.
- 521 • **Spatial resolution:** Measurements of the transfer of information from the image
522 data to the luminance fields at different spatial frequencies of interest typically
523 done by reporting the modulation transfer function. Non-isotropic resolution
524 properties should be characterized properly by providing two-dimensional
525 measurements or measurements along at least two representative axes. See *7.7*
526 *Effective Resolution, IDMS 1.03*.

Contains Nonbinding Recommendations

- 527
- 528
- 529
- 530
- 531
- 532
- 533
- 534
- 535
- 536
- 537
- 538
- 539
- 540
- 541
- 542
- 543
- 544
- 545
- 546
- 547
- 548
- 549
- 550
- 551
- 552
- 553
- 554
- 555
- 556
- 557
- 558
- 559
- 560
- **Pixel defects (count and map):** Measurements (counts) and location of pixel defects. This is typically provided as a tolerance limit. Pixel defects can interfere with the visibility of small details in medical images. See *7.6 Defective Pixels, IDMS 1.03*.
 - **Artifacts:** Evaluate for image artifacts such as ghosting and/or image sticking from displaying a fixed test pattern for a period of time. See *4.6 Artifacts and Irregularities, IDMS 1.03*.
 - **Temporal response:** Measurements of the temporal behavior of the display in responding to changes in image values from frame to frame. Since these transitions are typically not symmetric, rise and fall time constants are needed to characterize the system. See *10.2.3 Gray-to-Gray Response Time, IDMS 1.03*.
 - **Maximum and minimum luminance (achievable and recommended):** Measurements of the maximum and minimum luminance that the device outputs as used in the application under recommended conditions and the achievable values if the device is set to expand the range to the limit. See *2.4 Vantage-Point Suite of Measurement, IDMS 1.03*.
 - **Grayscale:** Measurements of the mapping between image values and the luminance. See *6.1 Grayscale, IDMS 1.03*.
 - **Luminance uniformity and Mura test:** Measurements of the uniformity of the luminance across the display screen. See *8.1.2 Sampled Vantage-Point Uniformity and 8.2.3 Mura Analysis, IDMS 1.03*.
 - **Stability of luminance and chromaticity response with temperature and lifetime**
 - **Bidirectional reflection distribution function:** Measurements of the reflection coefficients of the display device. Specular and diffuse reflection coefficients can be used as surrogates for the full bidirectional reflection distribution function. See *11.12 Diagnostic: Characterizing Hemisphere Uniformity, IDMS 1.03*.
 - **Gray Tracking:** Chromaticity at different luminance levels as indicated by the color coordinates in an appropriate units system (e.g., CIE $u'v'$). See *AAPM Task Group 196 Report*.
 - **Color scale:** Color coordinates of primary and secondary colors as a function of the digital driving level and their additivity. See *6. Gray- and Color-Scale Measurement and 5.4 Color-Signal White, IDMS 1.03*.
 - **Color gamut volume:** See *5.31 Volume-Color-Reproduction Capability, IDMS 1.03*.

IV(A)(11)(c). Resources

564 Those interested in learning more about these types of display considerations should
565 consider reading:

- 566
- 567
- 568
- 569
- 570
- 571
- 572
- *IDMS 1.03 - Information Display Measurements Standard Version 1.03, International Committee for Display Metrology, Society for Information Display, www.icdm-sid.org*
 - E. Samei, A. Badano, D. Chakraborty, K. Compton, C. Cornelius, K. Corrigan, M. J. Flynn, B. Hemming, N. Hangiandreou, J. Johnson, M. Moxley, W.

Contains Nonbinding Recommendations

- 573 Pavlicek, H. Roehrig, L. Rutz, J. Shepard, R. Uzenoff, J. Wang, and C. Willis,
574 *Assessment of display performance for medical imaging systems, Report of the*
575 *American Association of Physicists in Medicine (AAPM) Task Group 18,*
576 *Technical Report, AAPM (April 2005).*
577
- 578 • IEC 62563-1:2009, *Medical electrical equipment – Medical image display*
579 *systems – Part 1: Evaluation methods*
580
 - 581 • Amendment 1 to IEC 62563-1: *Medical image display systems – Part 1:*
582 *Evaluation methods*
583
 - 584 • The guidance entitled “*Guidance for Industry and FDA Staff: Display Accessories*
585 *for Full-Field Digital Mammography Systems-Premarket Notification (510(k))*
586 *Submissions*”
587 ([http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceD](http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/ucm107549.htm)
588 [ocuments/ucm107549.htm](http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/ucm107549.htm)).
589


IV(B). System-level Assessment

590
591
592 This subsection details the test methods at the system level that should be included in the
593 technical performance assessment of a WSI device. In this guidance, *system* refers to a
594 series of consecutive components in the imaging chain with clearly defined, measureable
595 input and output. For example, a system-level test can be designed for the image
596 acquisition subsystem, the image display subsystem, or a combination of both. The goal
597 of system-level tests is to assess the composite performance of a series of consecutive
598 components in the imaging chain. System-level tests should be conducted when the
599 component-level tests are either unfeasible or unable to capture the interplay between
600 components.

601
602 The common framework of the system-level tests described in this section is to compare
603 the system under test with an ideal system based on the same input, and then report the
604 difference between their outputs quantitatively. Designing such a system-level test
605 typically involves the following steps: (1) define the scope of the system and its input and
606 output, (2) define the input, which in most cases is a test target or phantom, (3) measure
607 the input to establish the ground truth that would be generated by an ideal system, (4)
608 measure the output of the system under test, and (5) calculate the errors between the truth
609 and the output with a quantitative metric. The framework of a typical system-level test is
610 shown in Diagram 2. Notice that the *ideal system* is a hypothetical device that generates
611 the perfect output with respect to the objective of the test such as color or focus. The
612 purpose of the ideal system is to define the intended behavior of the system under test.
613 The ideal system does not need to be implemented. Instead, the ideal system should be
614 simulated by a test method that establishes the truth of the input phantom.

615
616
617
618

Diagram 2: Framework of a typical system-level test.



IV(B)(1). Color Reproducibility

IV(B)(1)(a). Description

Color reproducibility is one of the key characteristics of a WSI system. The color characteristics are determined by every component in the imaging chain. Therefore, the color characteristics might be best evaluated at the system level. Color reproducibility indicates the accuracy and precision of the color transformation from the tissue sample on the slide to the image on the display. The colors of the tissue specimen should be accurately and precisely reproduced on the display based on the color reproduction intent, which should be clearly defined and justified by the sponsor.

IV(B)(1)(b). Test Methods

The WSI system should be tested with a target slide. The target slide should contain a set of measurable and representative color patches. Ideally the color patches should have similar spectral characteristics to stained tissue. The color patches should include a grayscale ramp for evaluating the grayscale response. The truth of the color patches should be measured with proper apparatuses separately.


For each color patch, the intended color (i.e., the expected output color based on the color reproduction intent defined by the Sponsor) should be calculated based on the truth of the color patches.

The target slide should be scanned and displayed by the WSI system. The output color of each color patch should be measured from the display.

The three datasets – truth, intended color, and output color – should be compared and analyzed. The sponsor should provide a rationale if the intended color is different from the truth.

Contains Nonbinding Recommendations

Diagram 3: Framework of the system-level color reproducibility test.



IV(B)(1)(c). Resources

Useful references on the subject of color reproducibility can be found at the International Color Consortium website <http://www.color.org>.

IV(B)(2). Spatial Resolution

IV(B)(2)(a). Description

Spatial resolution is another key characteristic of a WSI system. The goal of this system-level test is to evaluate the composite optical performance of all components in the image acquisition phase (i.e., from slide to digital image file).

IV(B)(2)(b). Test Methods

The following test is recommended for assessing spatial resolution of the image acquisition phase:

- Resolution and spatial frequency response: ISO 12233:2014(E) — Photography — Electronic still picture imaging — Resolution and spatial frequency responses.

IV(B)(3). Focusing Test


- The quality of focus in WSI can be affected by a number of inter-related factors, including the scanning method and approaches for constructing a focus map. Due to a trade-off between the number of focus points and the overall speed of the scanning process, focusing is typically based on a sample of focus points, determined automatically (auto-focus) or manually by the user. Since tissue can have uneven depth, auto-focus algorithms are needed to detect and adjust for different depths of focus.
- Data demonstrating that the focus quality is acceptable, even in the presence of uneven tissue, should be provided. Such data with proper justification could be derived from a phantom study, from clinical data, or both in a complementary fashion. The technology of phantom construction for testing focus is under development and this guidance will be updated as such technologies become available. Sponsors could attempt to build their own phantoms for testing depth

Contains Nonbinding Recommendations

of focus for their device. Alternatively, sponsors could provide experimental data using clinical tissue slides. Sampling of cases for such an experiment should be enriched for uneven tissue cases within a range representative of typical laboratory output. Alternative approaches for assessing the focus quality of a WSI will be considered along with proper justification. In addition, the following specifications should be provided, if applicable:

- Focus method: auto-focus for high-throughput or user-operated focus points
- Instructions for the selection of manual focus points (if applicable), including number of focus points and location in relation to a tissue sample
- Metrics used to evaluate focusing and description of methods to extract them
- Methods for constructing focus map from sample focus points

Diagram 4: Framework of the system-level focusing test.



IV(B)(4). Whole Slide Tissue Coverage

IV(B)(4)(a). Description

During the scan phase, WSI systems usually skip blank areas where tissue is absent in order to reduce scan time and file size. The purpose of the whole slide tissue coverage test is to demonstrate that all of the tissue specimen on the glass slide is included in the digital image file.


IV(B)(4)(b). Test Method

Sponsors should include a test that demonstrates the completeness of the tissue coverage. Sponsors should describe the test method and include the following items:

- Selection of the input tissue slide
- How to determine the complete coverage of the input tissue slide
- How to measure the actual coverage of the WSI output
- Calculate the ratio of the actual to complete coverage

Contains Nonbinding Recommendations

Diagram 5: Framework of the system-level whole slide tissue coverage test



IV(B)(5). Stitching Error

IV(B)(5)(a). Description


Stitching is the technique that enables a WSI system to combine thousands of sub-images into a single whole-slide image. Although during the scanning process a certain amount of overlapping between adjacent sub-images is maintained for alignment purposes, successful stitching relies on the texture present in the overlapped area. When the stitching algorithm fails to align two sub-images seamlessly, the error may or may not be perceivable by the human reader depending on whether noticeable stitching artifacts are generated. Therefore, a system-level test should be conducted when assessing the stitching quality of the WSI system.

IV(B)(5)(b). Test Methods

Sponsors should include a test that evaluates the stitching errors and include the following items:

- Selection of the input test slide
- Method for sampling of the stitching boundaries where stitching errors might occur
- How to determine the ideal stitching as the ground truth
 - For example, the region of the stitching boundaries can be re-imaged in one shot such that there is no stitching artifact.
- How to evaluate quality of the actual stitching based on the perfect stitching
 - For example, compare the image of stitching boundaries with the perfect one that does not have stitching artifact. The difference between these two images can be used as a figure of merit of the stitching quality.

Diagram 6: Framework of the system-level stitching error test



801

IV(B)(6). Turnaround Time

802

803

804

IV(B)(6)(a). Description

805

806 Turnaround time is the time required by the WSI system to execute a particular user
807 operation such as panning/zooming where the software and I/O (input/output) devices
808 retrieve image data, execute the computation, and refresh the image on the display. The
809 turnaround time starts when the user enters a command via a keyboard stroke or a mouse
810 click/movement and finishes when the image is completely updated on the display.

811 Turnaround time is important for a WSI system when fast and repetitive panning
812 operations are performed during a search task, which is delay-free in an optical
813 microscope. Prolonged, unpredictable turnaround time may impact the user's diagnostic
814 performance. The user interface should properly prompt the user when the operation is
815 incomplete and the requested image is not available. The turnaround time may vary
816 greatly depending on the user-requested operation, image content, data size/location,
817 computer workload, display size, etc. The sponsor should report the typical turnaround
818 time as well as the test method and test conditions.

819

IV(C). User Interface

820

821

IV(C)(1). Description

822

823

824 The user interface covers all components and accessories of the WSI system with which
825 users interact while loading the slides and acquiring, manipulating, and reviewing the
826 images. It also includes preparing the system for use (e.g., unpacking, set up,
827 calibration), and performing maintenance. Elements of the user interface have been
828 noted in many of the preceding sections and include two broad categories:

- 829 • Options through which the user operates the WSI system, such as:
 - 830 ○ Software menu options (e.g., scanning parameters)
 - 831 ○ Physical controls (e.g., clips on the slide feeder)
 - 832 ○ Connectors and connections (e.g., cables connecting system components)
- 833 • Information presented to the user through
 - 834 ○ Visual displays (e.g., scanned image, software menus)
 - 835 ○ Sounds (e.g., tone played when scanning completed)
 - 836 ○ Instructions (e.g., software users' manual)
 - 837 ○ Labels

838

IV(C)(2). Test Methods

839

840

841 It is recommended that the analysis to identify the use-related hazards of the WSI system
842 include the consideration of use errors involving failure to acquire, perceive, read,
843 interpret, and act on information from the WSI system correctly or at all and the harm
844 that could be caused by such errors. A human factors/usability validation test should be
845 performed to demonstrate that representative users of the WSI system can perform
846 essential tasks and those critical to safety under simulated use conditions.

Contains Nonbinding Recommendations

847

848 When selecting participants for validation testing, sponsors should carefully consider user
849 capabilities and expectations that could potentially impact the safe and effective use of
850 the WSI system. Examples of items that should be considered, if applicable, include
851 visual acuity and type of vision correction and the impact of expectations formed from
852 prior experience with other systems (e.g., optical microscope).

853

854 When selecting the critical tasks to be evaluated, sponsors should incorporate all known
855 use related errors and problems from similar devices (devices having similar
856 technological characteristics and indications for use) into the validation testing.

857 Consideration also should be given to whether task performance changes over time, and
858 if test duration needs to account for user fatigue. Examples might include a user altering
859 a task sequence in response to fatigue from repetitive image selection and manipulation
860 with mouse or keyboard.

861

862 When creating the simulated use conditions for validation testing, special consideration
863 should be given to the location of the WSI system primary workstation, its components,
864 their arrangement and how their locations affect user performance. Examples of location
865 considerations might include multiple monitors, a monitor with sub-optimal display
866 settings, or glare on a monitor from indoor lighting.

867

868 A human factors/usability validation test report should generally include the information
869 found in Table 1.

870

871 **Table 1: Items a Human Factors/Usability Validation Test Report Should Include**

872

Section	Contents
1	Intended device users, uses, use environments, and training <ul style="list-style-type: none">• Intended user population(s) and critical differences in capabilities between multiple user populations• Intended uses and operational contexts of use• Use environments and key considerations• Training intended for users and provided to test participants
2	Device user interface <ul style="list-style-type: none">• Graphical depiction (drawing or photograph) of device user interface• Verbal description of device user interface
3	Summary of known use problems <ul style="list-style-type: none">• Known problems with previous models• Known problems with similar devices

Contains Nonbinding Recommendations

	<ul style="list-style-type: none">• Design modifications implemented in response to user difficulties
4	User task selection, characterization and prioritization <ul style="list-style-type: none">• Risk analysis methods• Use-related hazardous situation and risk summary• Critical tasks identified and included in HFE/UE validation tests
5	Summary of formative evaluations <ul style="list-style-type: none">• Evaluation methods• Key results and design modifications implemented• Key findings that informed the HFE/UE validation testing protocol
6	Validation testing <ul style="list-style-type: none">• Rationale for test type selected (i.e., simulated use or clinical evaluation)• Number and type of test participants and rationale for how they represent the intended user populations• Test goals, critical tasks and use scenarios studied• Technique for capturing unanticipated use errors• Definition of performance failures• Test results: Number of device uses, success and failure occurrences• Subjective assessment by test participants of any critical task failures and difficulties• Description and analysis of all task failures, implications for additional risk mitigation
7	Conclusion <p>A statement to the effect that “The <device name/model> has been found to be reasonably safe and effective for the intended users, uses and use environments” should be included under the following conditions:</p> <ul style="list-style-type: none">• The methods and results described in the preceding sections support this conclusion.• Any residual risk that remains after the validation testing would not be further reduced by modifications of design of the user interface (including any accessories and the Instructions for Use (IFU)), is not needed, and is outweighed by the benefits that

Contains Nonbinding Recommendations

may be derived from the device's use.

873

874

Recommended methods for performing a human factors/usability validation test are described in the resources listed in section IV(C)(3) entitled “Resources” directly below.

875

876

The goal of testing is to assure that users can operate the WSI system successfully for the intended uses without negative clinical consequences to the patient and that potential use errors or failures have been eliminated or reduced.

877

878

879

880

IV(C)(3). Resources

881

882

FDA recognizes standards published by national and international organizations that apply human factors engineering/usability engineering (HFE/UE) principles to device design and testing. The recognized standards listed below provide suggestions on conducting an analysis of use-related hazards and a human factors/usability validation test to assess the safety and effectiveness of the final device design.

883

884

885

886

887

888

- ISO 14971:2007, *Medical Devices – Application of Risk Management to Medical Devices*: Provides systematic process to manage the risks associated with the use of medical devices.

889

890

891

- AAMI/ANSI HE75:2009, *Human Factors Engineering – Design of Medical Devices*: Comprehensive reference of recommended practices related to human factors design principles for medical devices.

892

893

894

- IEC 62366-1:2015, *Medical devices – Application of usability engineering to medical devices*: Describes the process to conduct medical device usability testing and incorporate results into a risk management plan.

895

896

897

In addition, FDA has published guidance with human factors related recommendations to assist manufacturers and facilitate premarket review. The guidance entitled “Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices” (<http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/ucm089543.htm>). This guidance document provides recommendations to industry regarding premarket submissions for software devices, including stand-alone software applications and hardware-based devices that incorporate software. It includes test methods to assure that the software conforms to the needs of the user and to check for proper operation of the software in its actual or simulated use environment.

898

899

900

901

902

903

904

905

906

907

IV(D). Labeling

908

909

The premarket application must include labeling in sufficient detail to satisfy the requirements of 21 CFR Part 801 and 21 CFR 809.10. The labeling includes supplementary information necessary to use and care for the WSI system such as instruction books or direction sheets and software user manuals.

910

911

912

913

Contains Nonbinding Recommendations

914 Although instructions, labeling, and training can influence users to use devices safely and
915 effectively, they should not be the primary strategy used to control risk. Modification of
916 the user interface design is a more effective approach to mitigate use-related hazards.

917

IV(D)(1). Test Methods

918

919 It is recommended that studies on labeling and training be conducted separately from
920 other human factors/usability validation testing. Human factors/usability validation
921 testing should be conducted with the final version of the labeling and related materials.
922 Timing and content of training should be consistent with that expected of actual users.

923

IV(D)(2). Resources

924

925 FDA has published several guidance documents on labeling to facilitate premarket
926 review and assist manufacturers.

927

- 928 • The guidance entitled “Labeling - Regulatory Requirements for Medical Devices”
929 ([http://www.fda.gov/downloads/MedicalDevices/DeviceRegulationandGuidance/
930 GuidanceDocuments/UCM095308.pdf](http://www.fda.gov/downloads/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/UCM095308.pdf)).
931
 - 932 ○ This publication covers labeling issues that device manufacturers,
933 reconditioners, repackers, and relabelers should consider when a product
934 requires labeling. Labeling includes adequate instructions for use,
935 servicing instructions, adequate warnings against uses that may be
936 dangerous to health, or information that may be necessary for the
937 protection of users.
- 938 • The guidance entitled “Device Labeling Guidance #G91-1 (blue book memo)”
939 ([http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceD
940 ocuments/ucm081368.htm](http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/ucm081368.htm)).
941
 - 942 ○ This guidance is intended to ensure the adequacy of, and consistency in
943 device labeling information. It is intended for use by industry in preparing
944 device labeling.
- 945 • The guidance entitled “Human Factors Principles for Medical Device Labeling”
946 ([http://www.fda.gov/downloads/MedicalDevices/DeviceRegulationandGuidance/
947 GuidanceDocuments/UCM095300.pdf](http://www.fda.gov/downloads/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/UCM095300.pdf)).
948
 - 949 ○ This report presents the principles of instruction, human factors, and
950 cognitive psychology that are involved in designing effective labeling for
951 medical devices.

952

IV(E). Quality Control

953

954 Sponsors should provide information on the quality control procedures, including
955 frequency and testing methods to be performed by the laboratory technologists and/or
956 field engineers with associated quantitative action limits. Discussions of tests for
957 constancy should include discussions of the slide feeder and scanning mechanisms,
958 coverage of the entire tissue slide, the bar code reader, the light source, the imaging
959 sensor device, and the calibrations at the component and system level. A detailed quality
960 control manual should be provided.