



St. Jude Medical Nanostim Leadless Pacemaker

Presentation to the
Circulatory System Devices Panel
February 18, 2016

Introduction

Mark Carlson, MD

Chief Medical Officer and VP Global Clinical
Affairs, St. Jude Medical

Adjunct Professor

Case Western Reserve University

Leadless Pacemaker Rationale for Development: Eliminate Issues with Pacemaker Pockets


- Discomfort (1.9%)¹
- Cosmetic concerns
- Hematomas (3.0%)¹
- Infections (2.7%)¹



¹ Udo et al, *Heart Rhythm* 9:728–735 (2012)


Leadless Pacemaker Rationale for Development: Eliminate Issues with Pacemaker Leads

- Mechanical failures (1.5%)¹
- Infections (0.2%)¹
- Mobility restrictions



¹ Udo et al, *Heart Rhythm* 9:728–735 (2012)

Substantial Incidence of Acute and Chronic Complications with Standard Pacemakers



Udo et al, *Heart Rhythm* 9:728-735 (2012)

Note: Includes both single and dual chamber pacemakers

Description of Device and Procedure

Today's Leadless Pacemaker System

The Nanostim Device

- 42 mm (~1 2/3") long
- 6 mm (~1/4") wide
- Percutaneous femoral vein delivery
 - 18F introducer
 - Steerable catheter
- Self-contained in ventricle
 - No lead or surgical pocket
- Provides traditional single chamber pacing therapy in patients clinically indicated for VVI(R) pacemaker therapy




FDA Question
Q4

Today's Leadless Pacemaker System

The Nanostim Device

- Single-turn helix and short stabilizing nylon tines secure fixation
- Steroid eluting electrode
- Temperature-based rate response
- Long battery life (8-18 years)
- Catheter-based retrieval
- Magnet Mode
- MRI Compatible*



Leadless Pacemaker System Implantation Procedure



**THIS IS A 40 SECOND VIDEO OF THE
IMPLANT PROCEDURE**

Agenda

**Safety and Effectiveness
for Leadless II Study**

Vivek Reddy, MD

Professor of Medicine and Cardiology
Mount Sinai Hospital, New York

**Nanostim Leadless
EU Post Market Study
US Post Approval Study
Training Program**

Mark Carlson, MD

Chief Medical Officer
St. Jude Medical
Adjunct Professor
Case Western Reserve University

Additional Experts

Dr. Paul Friedman, MD

Director – Implantable Device Lab
Professor, Medicine
Mayo Clinic, Minnesota

Dr. Joshua Cooper, MD

Director – Cardiac Electrophysiology
Professor, Medicine
Temple University, Pennsylvania

Barathi Sethuraman, Ph.D

Vice President – Clinical Science
St. Jude Medical

Chris Hubbard

Vice President – Nanostim Technology
St. Jude Medical

Safety and Effectiveness of a Leadless Pacemaker: Leadless II Clinical Trial Results


Vivek Reddy, MD

Professor of Medicine and Cardiology
Mount Sinai Hospital, New York

Leadless II Clinical Trial Overview

- Prospective, non-randomized
- Single chamber right ventricular pacing in clinically-indicated patients for traditional systems
- 56 Centers in US, Canada and Australia
 - 100 Operators
- N=667
- N=300 for pre-specified primary analysis

Leadless II Clinical Trial Patient Disposition



Demographics Reflect Elderly Population With Significant Comorbidities

| Demographic Variable | Primary Analysis Cohort (N=300) | Total Cohort (N=526) |
|--------------------------------|--|-----------------------------|
| Mean Age (years)± SD | 75.7 ± 11.6 | 75.8 ± 12.1 |
| Sex - Female | 35.7 | 38.2 |
| Coronary Artery Disease | 40.3 | 38.2 |
| Hypertension | 84.0 | 79.8 |
| Diabetes Mellitus | 27.3 | 27.3 |
| Anticoagulants | 60.0 | 58.9 |
| Antiplatelets | 47.7 | 47.0 |

Key Procedural Characteristics

| Procedural Characteristics | Primary Analysis Cohort (N=300) | Total Cohort (N=526) |
|---|--|-----------------------------|
| Successful implantation - n (%) | 289 (96.3%) | 504 (95.8%) |
| Device Repositioning | | |
| None | 68.9% | 70.2% |
| 1 | 18.3% | 17.7% |
| 2 | 8.3% | 7.7% |
| >2 | 4.5% | 4.4% |
| Final Device Position in Right Ventricle | | |
| Apical | 48.4% | 38.1% |
| Septum | 51.6% | 60.7% |
| Other | 0 | 1.2% |

Primary Effectiveness and Safety Endpoints Achieved

| | Population | P-Value |
|--|------------|---------|
| Effectiveness: | | |
| Acceptable pacing capture threshold | ITT | 0.007 |
| AND | | |
| Therapeutically acceptable sensing amplitude at 6 months | Implanted | <0.001 |
| Safety: | | |
| Freedom from Serious Adverse Device Effects through 6 months | ITT | <0.001 |

Serious Adverse Device Effects

| Serious Adverse Device Effect | Primary Analysis Cohort (N=300) | | Total Cohort (N=526) | |
|-------------------------------|---------------------------------|------------|----------------------|------------|
| | n | % | n | % |
| Total Patients | 20 | 6.7 | 34 | 6.5 |
| Cardiac perforation | 4 | 1.3 | 8 | 1.5 |
| Vascular complications | 4 | 1.3 | 6 | 1.1 |
| Device dislodgement | 5 | 1.7 | 6 | 1.1 |
| Pacing threshold elevation | 4 | 1.3 | 4 | 0.8 |
| Other | 4 | 1.3 | 13 | 2.5 |

Other Events included: Arrhythmia during device implantation, Intra-procedural device migration, Orthostatic hypotension with weakness, Pericarditis, presumed Pulmonary embolism, Hemothorax, Angina pectoris, Acute confusion and expressive aphasia, Dysarthria and lethargy after implantation, Contrast-induced nephropathy, Left-leg weakness during implantation, Ischemic stroke


Events of Interest: Cardiac Perforation and Vascular Complications

- Cardiac perforation (N=8)
 - 3 with surgical intervention
 - 2 with percutaneous intervention
 - 3 no intervention
 - 1 received traditional pacemaker
- Vascular complications (N=6)
 - 2 Access Site Hematoma
 - 2 Pseudoaneurysms
 - 1 AV Fistula
 - 1 Vascular Closure Malfunction

Events of Interest: Dislodgement and Retrieval

- Device dislodgement (N=6)
 - All reported in early post-op period (1 - 14 days)
 - All devices retrieved without issue
- Pacing threshold elevation with retrieval and new implantation (N=4)
 - All devices retrieved without issue


Freedom from SADEs: SADEs Occurred Within First Few Weeks of Procedure. No late SADEs.



| | | | | | | | | |
|----------|-------|-------|-------|-------|-------|-------|-------|-------|
| At Risk | 526 | 479 | 430 | 381 | 334 | 308 | 279 | 70 |
| Event | 23 | 34 | 34 | 34 | 34 | 34 | 34 | 34 |
| Success | 95.6% | 93.5% | 93.5% | 93.5% | 93.5% | 93.5% | 93.5% | 93.5% |
| 97.5 LCB | 93.5% | 91.1 | 91.1 | 91.1 | 91.1 | 91.1 | 91.1% | 91.1% |

Total Cohort (N=526)

Comparison of SADE Rates Nanostim vs. Traditional Pacemakers



Total Cohort (N=526)


¹ Udo et al, *Heart Rhythm* 9:728–735 (2012)

² Perforation and dislodgement for VVI pacemakers only; other data include single and dual chamber pacemakers

Procedure Related Mortality Adjudicated by Independent CEC

- No intraprocedural deaths
- 3 deaths (0.6%) adjudicated as procedure related
 1. 71 y.o. with cancer, respiratory arrest during implant (abandoned), required tracheotomy, made DNR, expired ~2 weeks later
 2. 89 y.o., successful implant, right groin hematoma, discharged home, expired ~ 2 weeks later
 3. 74 y.o., right atrial perforation, implant abandoned, large MCA stroke 2 days later and expired

No Significant Predictors of SADEs



Retrieval Animation



**THIS IS A 30 SECOND VIDEO OF THE
DEVICE RETRIEVAL**

Retrieval of Implanted Devices: 7 Retrievals, 100% Success Without SADEs

- Retrieval an important capability
- Time from implant to retrieval
 - Average 160 ± 180 days (Median = 100)
 - Range 1 - 413 days
- Reasons for retrieval
 - Elevated pacing thresholds (n=4)
 - CRT implantation (n=2)
 - Elective explant (n=1)

Leadless II Clinical Trial Summary

- Successfully implanted in ~96% of patients
- Trial met pre-specified Safety and Effectiveness endpoints
- Complication rate similar to conventional pacemakers
- Device is retrievable

Nanostim Leadless EU Post Market Study US Post Approval Study Training Program

Mark Carlson, MD

Chief Medical Officer and VP Global Clinical
Affairs, St. Jude Medical

Adjunct Professor

Case Western Reserve University

Learnings and Enhancements from the EU Post Market Study

- Enhanced patient selection criteria
- Required high resolution fluoroscopy
- Recommended septal rather than apical implants
- Enhanced training program

EU SADE Rate Decreased After Changes From Key Learnings Were Implemented

| | EU Post Market Study | | | |
|---|--------------------------|-----|--------------------------|-----|
| | Pre-Learnings (N=147) | | Post-Learnings (N=93) | |
| | % | % | n | % |
| Cardiac perforation or pericardial effusion | 6 | 4.1 | 2 | 2.2 |
| Device dislodgement | 2 | 1.4 | 0 | 0.0 |

U.S. Post Approval Study

U.S. Post Approval Study Overview (1 of 2)


- Prospective
- Non-randomized
- Multi-center
- Acute and long term safety including
 - Complications and success rate of removal/extraction
- Primary endpoint:
 - Freedom from Complication

FDA Question
Q2.C.iii


FDA Question
Q2.B.ii

U.S. Post Approval Study Overview (2 of 2)

- Data collected at:
 - Implant
 - Pre-discharge
 - Two weeks and
 - Semi-annually for 7 years
- Patient management at time of device replacement or deactivation
 - 30 day post replacement with traditional pacemaker
 - Continued follow up if replaced with Nanostim



FDA Questions
Q2.A.i
Q2.B.iv



FDA Questions
Q2.B.i
Q2.C.i
Q2.C.ii

Post Approval Study Overview

Sample Size

FDA Questions

Q2.A.i

Q2.B.i

Q2.B.ii

Q2.B.iii

Q2.C.i

34

- 1,700 patients
- Design allows for early and late AEs to be estimated to within a 90% CI width of 1%
- Study to include Leadless II and newly enrolled patients clinically indicated for single chamber pacing therapy

Mandatory Nanostim Physician Training Program

Prerequisite Requirements

- Qualified for pacemaker implantation
- An established practice affiliation with institution that has:
 - Resources to support implantation
 - High resolution fluoroscopy equipment
 - Proper emergency facilities for cardioversion, defibrillation, pericardiocentesis and cardio-pulmonary resuscitation

7-Module Training Program

Comprehensive Content

- Didactic Training / Patient Selection (Module 1)
- Hands-on Training
 - Implant Demonstration (Module 2)
 - Animal Lab Training (Module 3) or Virtual Reality Training (Module 5)
- Video Compendium Review (Module 4)
- Site-Training and onboarding, Case Observation, Technical and Implant Support and In-case Training provided by SJM certified personnel (Modules 6 and 7)

FDA Question
Q3

Virtual Reality

Reinforcing Correct Technique

- Benefits over animal lab
- Virtual reality demonstrates:
 - Catheter handle operations
 - Procedural steps
 - Best and worst practices
 - How to avoid complications
- Provides real-time critical warning messages and feedback

7-Module Training Program

Comprehensive Content

- Didactic Training / Patient Selection (Module 1)
- Hands-on Training
 - Implant Demonstration (Module 2)
 - Animal Lab Training (Module 3) or Virtual Reality Training (Module 5)
- Video Compendium Review (Module 4)
- Site-Training and onboarding, Case Observation, Technical and Implant Support and In-case Training provided by SJM certified personnel (Modules 6 and 7)

Physician Certification Contingent on Completion of Training Program

- Physician Certification received after successful completion of
 - All modules
 - 10 procedures with technical and implant support and in-case training provided by SJM certified personnel

Summary

- Complication rates similar to alternative therapies
- Absence of longer-term SADEs
- Absence of certain complications associated with standard pacemakers
- Robust training program will support safe use
- Event rates will continue to be monitored in post-approval study




St. Jude Medical Nanostim Leadless Pacemaker

Presentation to the
Circulatory System Devices Panel
February 18, 2016

Changes to Transfer Nanostim Learnings


| Lessons Learned | Action Taken |
|---|--|
| Nanostim used as device of last resort | Align inclusion/exclusion criteria with IDE study. Stress care during patient selection by SJM field personnel |
| All perforations were associated with RV apical implants | Placement in lower septum |
| Quick rotation may cause catheter to torque and over-rotate | Slowly rotate catheter with pauses, 1-1 ¼ turns rather than 1¼ |
| Pressure on the endocardium increases if protective sleeve is not fully retracted and if the catheter buckles | Fully pull back protective sleeve before engaging endocardium and apply forward pressure gently so the device is moving with the cardiac cycle |
| COI associated with active fixation leads, associated with higher initial thresholds | Wait up to 20 mins for COI to resolve |
| Suboptimal imaging equipment contributed to at least one cardiac perforation | Sites were required to use high resolution fluoroscopy equipment for implantation |
| Presence of an existing perforation before device implant was observed in at least one case | IFU warning added to not implant device in presence of an existing perforation |

Primary Effectiveness Endpoint Surpassed the Performance Goal



SADE Rate

First 10 Cases vs >10 Cases



Threshold Elevation (N=4)

- Patient #1
 - Elevated Pacing Threshold at Implant
 - 2-week visit (Device reprogrammed)
 - 100 days post implant-LP retrieved and replaced with another LP
- Patient #2
 - Elevated Pacing Threshold at implant
 - The next day- LP retrieved and replaced with another LP
- Patient #3
 - Elevated Pacing Threshold at implant
 - The next day- LP retrieved and replaced with another LP
- Patient #4
 - Elevated Pacing Threshold- 72 hrs. post implant
 - Device Reprogrammed/Temporary pacer placed the following day
 - 23 days post implant-LP retrieved and replaced with transvenous ppm

Table 3-10: Deaths Classified by CEC Adjudication in Total Cohort

| Cause of Death | Number of Patients | Relation to Device or Procedure | Number of Days Post-Implant |
|--|--------------------|----------------------------------|-----------------------------|
| Cardiac | | | |
| Arrhythmic | 2 | Not Related (1); Procedure (1) | 18, 100 |
| Heart failure | 1 | Not Related (1) | 99 |
| Unknown | 1 | Procedure/Introducer (1) | 14 |
| Non-cardiac | | | |
| Accidental gunshot wound | 1 | Not Related (1) | 47 |
| Renal or liver failure | 5 | Not Related (5) | 73, 82, 89, 135, 320 |
| Respiratory failure | 3 | Procedure (1) Not Related (2) | 10, 103, 182 |
| Multiple organ failure | 2 | Not Related (2) | 34, 38 |
| Ischemic bowel/small bowel obstruction | 2 | Not Related (2) | 185, 270 |
| Mixed respiratory and metabolic acidosis | 1 | Not Related (1) | 176 |
| Unknown* | | | |
| Death- Sudden with antecedent worsening heart failure | 1 | Not Related (1) | 267 |
| Death- Sudden without antecedent worsening heart failure | 1 | Not Related (1) | 274 |
| Death- Non-sudden with antecedent worsening heart failure | 2 | Not Related (2) | 18, 42 |
| Death- Non-sudden with antecedent worsening heart failure status unknown | 1 | Not Related (1) | 281 |
| Death-Unknown (presumed sudden) with no antecedent worsening heart failure | 3 | Not Related (2) Unknown (1) | 5, 69, 126 |
| Death-Unknown (presumed sudden) with antecedent worsening heart failure status unknown | 1 | Not Related (1) | 219 |
| Death- Unknown temporal cause and antecedent worsening heart failure status unknown | 1 | Not Related (1) | 409 |
| Total | 28 | | |


* Sudden death: death \leq 1 hour after onset of symptoms

Non-sudden death: death > 1 hour after onset of symptoms

Death Unknown (presumed sudden): documentation of patient's condition by a witness within 24 hours


Death Unknown: death where onset of symptoms cannot be determined

Enrollment By Operator



Total Cohort N=526

PAS - Projection of Patient Enrollment and Device End of Service



Assumes 15% attrition per year and 2.7% end of service per year