

MCH Water Path & Airflow Design

Cleanability, Compatibility, Compliance, Communication





Presenters

✦ Douglas E. Platt, Partner

- * 25+ years experience servicing & refurbishing perfusion equipment
- * Two patents in cooler-heater design

✦ Jon L. Gardner, Chief R&D Engineer


- * 8 years experience with perfusion equipment
- * One patent in cooler-heater design
- * Fifth-generation medical professional

CardioQuip MCH-1000



- ✦ Same hardware, software, & design
- ✦ Different tank sizes for different cooling needs

Open vs. Closed Water Path



Open Water Path Design



- No inaccessible airspace to harbor microbes
- Water-borne disinfectants contact all wetted components
- Tank is easy to inspect & clean

Cleaning Protocol

- ✦ Regular inspection
- ✦ Add disinfectant & circulate
- ✦ Wipe/scrub the tank
- ✦ Flush & fill twice
- ✦ Wipe down external surfaces



Cleaning Protocol Validation

- ✦ Fill device with contaminated water
- ✦ Sample water for HPC testing (pre-cleaning)
- ✦ Perform cleaning procedure
- ✦ Sample water for HPC testing (post-cleaning)
- ✦ Evaluate “cleanability”
 - * Objective: time required, inspection results
 - * Subjective: ease of inspection & cleaning


Validation Results

- ✦ The MCH is easy to clean
- ✦ The cleaning protocol is effective

Disinfectant	CFU/mL (Pre)	CFU/mL (Post)
PineSol™	15,100	4.0
Chlorine Bleach (5%)	2,570	60.0
Hydrogen Peroxide (3%)	1,770,000	9.3
MinnCare™ (22% H ₂ O ₂ , 5% Peracetic Acid)	195,000	<200

Why Use Water?

- ✦ Regulatory history
- ✦ Materials compatibility
- ✦ Similarity to blood
- ✦ Superior performance



Water Risk Mitigations

- ✦ Valved connectors
- ✦ Dripless hose kits
- ✦ Anti-microbials
- ✦ Education




The Other Factor to Mitigate:

Airflow



Surgical Ventilation System



Factors for Airflow Validation



- Analyze with respect to Surgical Ventilation System standards (ASHRAE 17-2008, ISO 14644-3 et al)
- Consider airflow & history of similar devices

Airflow Validation

- Standard unit: ~50 CFM fan, bottom intake/exhaust
- Optional Refrigeration Module: front/rear intake, side exhaust, ~300 CFM (~150 CFM per side), smoke tested



Airflow Validation

- Compact unit: ~50 CFM fan, bottom intake, side/rear exhaust
- Optional Thermo-Electric Cooling Module, front/rear intake, side exhaust, ~100 CFM (~50 CFM per side)



Summary

Water Quality
Airflow

“Cleanability”
Effectiveness

Risk Assessment
Mitigation

Communication
Education



Contact Info

- ✦ CardioQuip, LP

- * ☎ 979-691-0202

- ✦ Douglas E. Platt, Partner

- * dplatt@cardioquip.com

- ✦ Jon L. Gardner, Chief R&D Engineer

- * jon@cardioquip.com

www.cardioquip.com