

From: Tull, Lori
Sent: Tuesday, February 02, 2016 2:26 PM
To: Wouter Van't Hof (wvanthof@clevelandcordblood.org)
Subject: FW: IR for STN 125594

Importance: High

Hi Wouter,

The CMC reviewer has the following request for information. Would you please respond by email and follow with an official amendment to the BLA by 2/16/16?

1. We note that although you have cross-referenced the [REDACTED] you have used [REDACTED] gauge needles to inoculate the [REDACTED] bottles for your validation studies (the [REDACTED] media package insert recommends that you use [REDACTED] gauge needles for the inoculation of [REDACTED] bottles). Please clarify your rationale for using [REDACTED] gauge needles (as using a larger gauge needle may compromise the seal of the [REDACTED]T bottles and the sterility test for the strict anaerobic microorganisms). Also, please clarify if you will be switching to [REDACTED] gauge needles for the [REDACTED] bottle inoculation for the actual sterility test (if applicable please update and submit the respective SOP) of HPC, Cord Blood.

2. Reference - Qualification of new lots of test media: Please change (based on the data from your Sterility Test method validation) the "expected outcomes" to "All tests Positive for Growth in [REDACTED] days" for [REDACTED] and "All tests Positive for Growth in [REDACTED] days" for [REDACTED].

Regards,
Lori Tull
Team Lead
Regulatory Management Staff
Office of Cellular, Tissue, and Gene Therapies
Center for Biologicals Evaluation and Research
(240) 402-8361

THIS MESSAGE IS INTENDED ONLY FOR THE USE OF THE PARTY TO WHOM IT IS ADDRESSED AND MAY CONTAIN INFORMATION THAT IS PRIVILEGED, CONFIDENTIAL, AND PROTECTED FROM DISCLOSURE UNDER LAW. If you are not the addressee, or a person authorized to deliver the document to the addressee, you are hereby notified that any review, disclosure, dissemination, copying, or other action based on the content of this communication is not authorized. If you have received this document in error, please immediately notify the sender immediately by e-mail or phone.