

InteliSwab®
COVID-19 RAPID TEST • Rx

Instrucciones de use en español
ubicadas dentro de la caja.

The InteliSwab® COVID-19 Rapid Test detects active infection in individuals 2 years and older. See Instructions for Use for full intended use.

QUICK AND EASY RESULTS

- Use at home or anywhere there is a flat surface
- No professional supervision or video consultation required
- No phone or other equipment needed to interpret results

CONTENTS

- One (1) pouch, with:
 - Single-use test device
 - Tube with developer fluid
- Instructions for Use in English and Spanish
- Test stand

You will need a way to time the test for 30 minutes (a watch, clock or cellphone). To report results to public health, visit www.InteliSwab.com/app.

Determining a negative result requires multiple tests. You may need to purchase additional tests to perform serial (repeat) testing. This test is more likely to give you a false negative result when you have COVID-19 than a lab-based molecular test.

Store the InteliSwab® COVID-19 Rapid Test in a dry location between 35°-86°F (2°-30°C).

For more information on expiration dating for this product, please refer to www.InteliSwab.com

QUESTIONS? Go to www.InteliSwab.com or call 1-833-601-0127

In the USA, this product has not been FDA cleared or approved, but has been authorized by FDA under an EUA.

The emergency use of this product is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of in vitro diagnostics for detection and/or diagnosis of COVID-19 under Section 564(b)(1) of the Federal Food, Drug, and Cosmetic Act, 21 U.S.C. § 360bbb-3(b)(1), unless the declaration is terminated or authorization is revoked sooner.

This test does NOT determine if you had COVID-19 in the past or if you have immunity.

This product has been authorized only for the detection of proteins from SARS-CoV-2, not for any other viruses or pathogens.

■ OraSure Technologies, Inc.
220 East First Street | Bethlehem, PA 18015 USA
(+1) 610-882-1820 | www.OraSure.com

© 2022 OraSure Technologies, Inc.
3001-3456 rev. 01/23



86°F



Swab
both
nostrils



Swirl
in the
tube



See
your
results

InteliSwab®
COVID-19 RAPID TEST • Rx

1 TEST

1 TEST

InteliSwab®
COVID-19 RAPID TEST • Rx

HOME TEST
RESULTS IN
MINUTES



actual size

Rx
ONLY

InteliSwab®
COVID-19 RAPID TEST • Rx

**SIMPLE &
CONVENIENT
TESTING**

No assembly
required
— ready to use

No waiting for
mail or lab
processing

**SIMPLY
SWAB, SWIRL
AND SEE**

**SMART SCIENCE
MADE SIMPLE™**

**ANYTIME,
ANYWHERE**

**LEARN
MORE**



**GENTLE SHALLOW
NASAL SWAB**

For Emergency Use Authorization (EUA) only
For *in vitro* diagnostic use

REF
LOT

YYYY-MM-DD
XXXXXX
XXXXXX



6 08337 00163 1

See side panel for more information on expiration dates.