

In-Depth Impurity Assessment of Synthetic Oligonucleotides Enabled by HRMS

*SBIA 2022: Advancing Generic Drug Development:
Translating Science to Approval*

Day 1, Session 1B: Oligonucleotide Active Pharmaceutical Ingredient (API) Sameness and Impurity Assessment Considerations

Kui Yang, Ph.D.

Sr. Research Scientist

Division of Complex Drug Analysis

Office of Testing and Research, Office of Pharmaceutical Quality

CDER | US FDA

9/20/2022

Learning Objectives

- Identify unique scientific and regulatory challenges in synthetic oligonucleotides
- Understand in-depth assessment of product-related oligonucleotide impurities
- Resolve complex impurity mixtures by high resolution mass spectrometry (HRMS)


Disclaimer: This presentation reflects the views of the authors and should not be construed to represent FDA's views or policies.

Challenges in Synthetic Oligonucleotides




- *Target a broad range of mRNAs* that encode critical cellular proteins (“undruggable”; great variety)
- *Pose unique scientific and regulatory challenges* (“Big small molecules”; currently no ICH regulatory guidelines or FDA general CMC guidances)
 - ✓ 1st PSG – Nusinersen, draft guidance in Feb 2022

Challenges in Impurity Assessment



Challenges in Impurity Assessment



Outline


- In-depth impurity assessment case study
 - HRMS-enabled identity and quantitation
 - Isobaric impurity case: deletion sequence variants

- High vs low resolution MS


HRMS-enabled quantitation



HILIC-HRMS: Extracted ion chromatogram (EIC) of isotopic peaks



RSD (%) of EIC peak area < 15%



RSD (%) of EIC peak area < 10%

A challenging scenario ...

Case study: resolve isobaric sequence deletion impurities

Isobaric: equal nominal mass

Full-length product (FLP):

2'-O-MOE, PS modified RNA 18-mer (nusinersen)

U-C-A-C-U-U-U-C-A-U-A-A-U-G-C-U-G-G (C and U are methylated)

Exact mass: 7122.2763; m/z (-4 charge state): 1779.5618

Isobaric sequence deletion impurities

Impurities	Exact mass (Da)	m/z (-4 charge state)
<i>n-U</i>	6728.2163	1681.0468
<i>n-C</i>	6729.2003	1681.2928




Individual n-U and n-C






Accurate mass

Isotopic envelope / distribution

n-U



n-C




Individual n-U and n-C






Accurate mass

Isotopic envelope / distribution


n-U



n-C



<i>n-U</i>		<i>n-C</i>		Δ
M	1681.0468			
M+1	1681.2975	M	1681.2928	0.0047
M+2	1681.5478	M+1	1681.5435	0.0043
M+3	1681.7981	M+2	1681.7938	0.0043
M+4	1682.0483	M+3	1682.0441	0.0042
M+5	1682.2984	M+4	1682.2943	0.0041
M+6	1682.5486	M+5	1682.5445	0.0041
M+7	1682.7987	M+6	1682.7946	0.0041
M+8	1683.0488	M+7	1683.0447	0.0041
M+9	1683.2989	M+8	1683.2948	0.0041



Coexisting *n-U* / *n-C* in a mixture:

- Not chromatographically separated
- Not mass differentiated
- Not individually reported in QC testing or specification

n-p(MOE^{Me}C) / n-p(MOE^{Me}U)


@ m/z 1682.2

n-U / n-C mixtures



<i>n-U</i>		<i>n-C</i>	
M	1681.0468		
M+1	1681.2975	M	1681.2928
M+2	1681.5478	M+1	1681.5435
M+3	1681.7981	M+2	1681.7938

- Monoisotopic peak at 1681.0 indicates **the presence of n-U**




n-U / n-C mixtures



<i>n-U</i>		<i>n-C</i>	
M	1681.0468		
M+1	1681.2975	M	1681.2928
M+2	1681.5478	M+1	1681.5435
M+3	1681.7981	M+2	1681.7938


- Monoisotopic peak at 1681.0 indicates **the presence of n-U**
- Mismatching isotopic pattern indicates **the coexistence** of n-U with n-C in the mixture




Calculated (predicted) vs Measured isotopic distributions



Individual
n-U and n-C






+




Linear combination

n-U / n-C
mixtures



Quantification of isobaric composition of n-U/n-C mixtures




Workflow

Best match:


Defined by: the difference between the predicted and the targeted is minimum.

Quantification of isobaric composition of n-U/n-C mixtures (covering 5%/95% n-U/n-C to 95%/5% n-U/n-C)




Unpublished data: A.M. Abdullah, C. Sommers, J. Rodriguez, D. Zhang, D. Kozak, J. Hawes, M. Sapru, and K. Yang. Manuscript under revision.

Quantification of isobaric composition of the n-U/n-C standard *spiked* in FLP



Unpublished data: A.M. Abdullah, C. Sommers, J. Rodriguez, D. Zhang, D. Kozak, J. Hawes, M. Sapru, and K. Yang. Manuscript under revision.

Quantification of isobaric composition of the n-U/n-C impurity **present** in FLP




FLP:

Custom synthesized by vendor:

- Synthesis – Solid phase phosphoramidite synthesis
- Purification method – RP HPLC

On-column n-U/n-C
< 0.02 pmol > 0.02 pmol



Total n-U/n-C in FLP quantified at
 $0.27 \pm 0.02\%$ of FLP

High vs low resolution MS

	HRMS	Unit- or low-resolution MS
Isotopic peak-resolved envelope	Yes	No
Identification of coexisting n-U and n-C (isobaric)	Yes	No (n-U/n-C)
Quantitation of coexisting n-U/n-C ratio	Yes	No (sum of n-U/n-C)

Summary



- Synthetic oligonucleotides, as an evolving class of therapeutics, pose unique scientific and regulatory challenges.
- In-depth impurity assessment offers comprehensive impurity evaluation and is crucial for comparative impurity analysis in generic drug development.
- High resolution MS demonstrates advantages in resolving complex oligonucleotide impurities, such as isobaric sequence variants.

Acknowledgement



- **Office of Testing and Research (OTR)**

- ✓ Oligo research team
A.M. Abdullah
- ✓ Division of Complex Drug Analysis (DCDA)
- ✓ OTR Management
Cindy Sommers
Jason Rodriguez
David Keire

- **Office of Generic Drugs (OGD)**

Deyi Zhang
Darby Kozak

- **FDA Critical-Path Research Grants (FY20, FY22)**

- ✓ PI: Kui Yang
- ✓ Collaboration offices:
OPQ, OGD, OND, NCTR

- **OTR-led Oligo SME team / working group**

- ✓ Participating offices:
OPQ, OGD, OND, NCTR, OTS



Challenge Question #1

Which of the following statements is NOT correct for synthetic oligonucleotide therapeutics (ONTs)?

- A. ONTs target a broad range of mRNAs.
- B. Like mRNAs, ONTs are regulated as biologics by FDA.
- C. Currently no ICH regulatory guidelines specifically address quality aspects for ONTs.
- D. FDA published the first PSG draft guidance for ONTs in 2022.



Challenge Question #2

True or False?

In-depth impurity assessment by HRMS is recommended in PSG draft guidance for Nusinersen published in Feb 2022.

- A. True
- B. False

Questions?

Kui Yang

Sr. Research Scientist, Division of Complex Drug Analysis
Office of Testing and Research, Office of Pharmaceutical Quality
CDER | US FDA

kui.yang@fda.hhs.gov