

MOBILE RADIOGRAPHIC SYSTEMS

FIELD TEST RECORD EDIT CHECKS

(Test Procedure MRA - Form FDA 2783)

Verify that:

1. The kVp is in the above 70 kVp range in item 4.
2. The certification status of each component is indicated in item 2.
3. If only four values are entered for reproducibility, no two exposures differ by more than ten percent of the highest value.
4. The value of tube current in item 5 does not differ from the value in item 34 by more than a factor of two.
5. If item 3 is marked "C," data is entered at items 46 and 47.
6. If a value of tube current is not entered in item 34, a value of mAs is entered at item 7 or the system is identified as fixed tube current in Remarks.
7. If the control was manufactured on or after May 1994 then data entered at item 34 is mAs product (mA x s) and not just mA.